**OPERATING INSTRUCTIONS** 



# TYPE 1611-B CAPACITANCE TEST BRIDGE

GENERAL RADIO COMPANY

**OPERATING INSTRUCTIONS** 

# TYPE 1611-B CAPACITANCE TEST BRIDGE

Form 1611-0100-H January, 1965

GENERAL RADIO COMPANY WEST CONCORD, MASSACHUSETTS, USA

### **SPECIFICATIONS**

CAPACITANCE RANGE:	0 to 11,000 $\mu$ f in eight ranges.
TEST-VOLTAGE FREQUENCY:	60 cps with internal supply and filter; 120 cps with ex- ternal supply on four higher multipliers; other audio fre- quencies to 1000 cps with reduced sensitivity, with ex- ternal supply and filter.
CAPACITANCE ACCURACY:	$\pm(1\%$ + $1\mu\mu$ f) over the entire range.
DISSIPATION-FACTOR RANGE:	0 to 60% at 60 cps (x $\frac{f}{60}$ at other frequencies).
DISSIPATION-FACTOR ACCURACY:	$\pm(2\% \text{ of dial reading } \pm 0.05\% \times \frac{f}{60}$ dissipation factor).
	above 100 pf. Power factor equals $\frac{D}{\sqrt{1+D^2}}$ where D = dissipation factor.
VOLTAGE APPLIED TO SAMPLE:	Varies from 125 v at 100 $\mu\mu$ f to less than 3 v at 10,000 $\mu$ f with the internal 60-cycle supply.
POLARIZING VOLTAGE:	On the four highest multipliers (1 to 1000) a d-c polar- izing voltage can be applied.
POWER SUPPLY VOLTAGE:	105 to 125 (or 210 to 250) volts, 60 cps.
ACCESSORIES SUPPLIED:	Power cord and spare fuses.
DIMENSIONS:	Width $14\frac{1}{2}$ in., height 10 in., depth 16 in., over-all, in- cluding cover and handles. (370 x 255 x 410 mm).
WEIGHT:	30½ lb. (14 kg).

# TABLE OF CONTENTS

Section 1 INTRODUCTION	:
1.1 Purpose	
1.2 Description	I
1.3 Accessories	\$
Section 2 OPERATING PROCEDURE	}
2.1 Power Supply	
2.2 External Generator and Filter	
2.3 Zero Adjustment	
2.4 Lead Capacitance	•
2.5 Measurement of Unknown Capacitance and Dissipation Factor .	
2.6 Measurement of Resistance	
2.7 Measurement of Electrolytic Capacitor with Polarizing Voltage .	
2.8 Measurements of Three-Terminal Capacitance	
2.9 Inductance Measurement	
2.10 Resistance Measurement	
2.11 Measurement of Dielectric Samples	
2.12 Measurement of Liquid Insulation	
2.13 Measurement of Insulators	
2.14 Transformer Measurements	
2.15 Measurement of Parallel Capacitance	4
2.16 Use as a Limit Bridge	4
Section 3 PRINCIPLES OF OPERATION	5
3.1 Bridge Circuit	5
3.2 Detector	5
3.3 Bridge Voltage Supply	5
3.4 Compensating Circuits	6
Section 4 SERVICE AND MAINTENANCE	7
4.1 General	7
4.2 Routine Maintenance	7
	7
4.4 Trouble-Shooting Procedure	7
	2



Figure 1. Type 1611-B Capacitance Test Bridge.

# Section 1 INTRODUCTION

1.1 PURPOSE. The Type 1611-B Capacitance Test Bridge (Figure 1) is a direct-reading capacitance bridge designed for wide-range capacitance and dissipation-factor measurements at 60 and 120 cps. (At 120 cps, an external generator is required; the Type 1214-D is recommended.) Measurements can be made at other audio frequencies up to about 1000 cps with reduced sensitivity. For such measurements an external generator and filter are required. The bridge is suitable for laboratory and shop testing of paper, mica, and electrolytic capacitors. It also meets the needs of the electric-power industry for shop testing of insulators, particularly for dissipationfactor measurement of bushings, insulators, transformer insulation, rotating machinery, and cables. For the wire and cable manufacturer, the bridge offers a convenient and rapid means of locating breaks in cable and of performing laboratory and production tests of dissipation factor and capacitance on all kinds of cable. The bridge's many applications in the communications industry include checking capacitance to ground of transformer windings, shields, and circuit elements, as well as routine capacitance and dissipation-factor tests on electrolytic capacitors and other components.

#### 1.2 DESCRIPTION.

1.2.1 GENERAL. The circuit used is the series resistance-capacitance bridge (Figure 2), in which one ratio arm is continuously variable and calibrated to read directly in capacitance. The other ratio arm is variable in decade steps and serves as a multiplier for the directreading dial. The variable resistors in series with the standard capacitors are calibrated directly in dissipation factor. A visual null detector is used, consisting of a tuned amplifier and an electron-ray tube.



Figure 2. (a, left) Bridge Circuit for Four Lower Multiplier Positions. (b, right) Bridge Circuit for Four Higher Multiplier Positions.

# 1.2.2 CONTROLS. The following controls are on the front panel of the Type 1611-B Capacitance Test Bridge:

Name	Туре	Function
POWER	2-pos toggle switch	Main power switch
FILTER	3-pos rotary switch	Selects either internal 60-cps, 120-cps, or external filter.
DIR-STANDBY-REV	3-pos rotary switch	Provides opposite polarities of voltage applied to unknown or removes test voltage on four lower multipliers.
SENSITIVITY	Continuous rotary control	Controls detector sensitivity.
GEN	2-pos toggle switch	Adapts instrument for internal or external audio supply.
MULTIPLY CAPACITANCE BY	9-pos rotary switch	Selects capacitance range, or permits detector to be used independently of bridge.
DISSIPATION FACTOR	Continuous rotary control	Dissipation-factor balancing controls
CAPACITANCE	Continuous rotary control	Capacitance balancing controls
D	2-pos toggle switch	Adds 30% dissipation factor to continuous range.

1.2.3 CONNECTORS. The following connectors are on the Type 1611-B Capacitance Test Bridge.

Name	Type	Function
SHIELD	Jack-top binding post	Used in three-terminal measurements (paragraph 2.8)
EXT GEN	Binding-post pair	Connection for external audio generator
EXT FILTER	Phone jack	Connection for external filter
DC	Binding-post pair	Connection for dc polarizing voltage

#### 1.2.4 MECHANICAL DETAILS.

1.2.4.1 General. The bridge is housed in a shielded airplane-luggage-type cabinet, with a cover that protects the panel and controls when the instrument is carried about and that keeps dust out when the bridge is not in use. Simplified operating instructions are mounted in the cover, as is the power cord.

1.2.4.2 Standard Capacitors. Two standard capacitors are used. One is a  $0.01-\mu f$  mica capacitor, Type 505, mounted on a subpanel (see Figure 15). The other is a special 1.0- $\mu f$  capacitor, made in two sections, with polystyrene tape used as the dielectric. After special heat treatment, aging, adjustment, and impregnation, it is mounted and sealed in the cylindrical metal containers located under the main shelf, as shown in Figure 15. 1.2.4.3 Dissipation-Factor Rheostats. The two rheostats (one for each capacitor) used to balance for dissipation factor are ganged to a common shaft and mounted on a subpanel (see Figure 14).

Each rheostat winding consists of two tapered sections. The resistances of the two sections are chosen so that the resulting scale permits precise readings at low values of dissipation factor, while at the same time retaining the convenience of a range up to 30% on a single scale. An etched dial scale is used, and four adjustable shunt resistors are provided, one across each section of each rheostat, to bring the actual resistance into agreement with the value required by the scale.

Two resistors controlled by a double-pole switch are connected one to each rheostat to increase the dissipation-factor range from 30% to 60%.

1.2.4.4 "A" Arm. The variable resistor  $R_A$ , by means of which the capacitance balance is obtained, is a tapered rheostat with a total resistance of about 11,000 ohms (see Figure 14). The taper is such that the dial scale is approximately logarithmic. An adjusting plate and cam built into the unit permit adjustment of the arm with respect to the dial at eight points. As adjusted at the factory, the resistance in kilohms corresponds to the dial reading within 0.5% over the main decade from 1.0 to 10.

1.2.4.5 Bridge Transformer. A shielded transformer couples the voltage source to the bridge (Figure 2a) or the bridge to the detector (Figure 2b). An electrostatic shield surrounds each winding, insuring a very low electrostatic

coupling between windings. The primary shield is grounded and the secondary shield is connected to the bridge. The secondary shield-to-ground capacitance is part of the standard capacitance arm. For this reason the transformer is wound with polystyrene tape as insulation to insure a low and constant dielectric loss.

The transformer is mounted in a high-permeability case to minimize magnetic coupling to the rest of the circuit and to any external fields (see Figure 13). The entire assembly is treated and impregnated to minimize the possible effects of moisture. 1.2.4.6 The "Magic-Eye" Tube. The electron-ray tube (Type 6E5) is so mounted in a slotted tube that its position with respect to the panel can readily be changed. It is held in position by a thumbscrew, which is accessible when the instrument is removed from the case. When the bridge is used in brightly lit locations, or when the brilliance of the eye has been reduced with age, the tube can be slid back in its mounting to provide additional light shielding.

1.3 ACCESSORIES. Supplied with the Type 1611-B Capacitance Test Bridge are a power cord and spare fuses.

## Section 2 OPERATING PROCEDURE

2.1 POWER SUPPLY. The instrument is shipped wired for either 105-125 or 210-250 volts, 60 cycles, as indicated on the nameplate near the power receptacle. The power transformer can be reconnected for either supply voltage in accordance with the information given on the schematic diagram, Figure 16.

For best results, connect the bridge to a good external ground. Use the ground terminal of either the UN-KNOWN or DC binding-post pair.

The seven-foot line connector cord is mounted in the instrument cover. Plug it into the power source and snap the POWER switch on.

2.2 EXTERNAL GENERATOR FILTER. For measurements at 60 cps, both generator and filter are contained in the Type 1611-B.

2.2.1 FREQUENCIES OTHER THAN 60 CPS — FOUR HIGHER MULTIPLIERS. For measurements at 120 cps on the four higher multipliers, connect a suitable external generator (refer to paragraph 3.3.2) to the EXT GEN terminals, set the GEN switch to EXT and the FILTER switch to 120 c.

For frequencies other than 60 or 120 cps, connect a suitable external generator (refer to paragraph 3.3.2) to the EXT GEN terminals, connect the proper filter at the EXT FILTER terminals, set the GEN switch to EXT and the FILTER switch to EXT.

A transformer is normally required between the voltage source and the EXT GEN terminals, both for impedance matching and for isolation. since the dc polarizing voltage is connected on the ground side. The Type 1214-D Unit Oscillator is equipped with such a transformer. (Refer to paragraph 3.3.2.) A connection diagram for external audio oscillator, matching transformer, and limiting resistors is shown in Figure 3.

2.2.2 120 CPS — FOUR LOWER MULTIPLIERS. Place GEN switch in EXT position. Connect the external 120-cycle supply to the power socket of the bridge and remove the rectifier tube (V3) to reduce power consumption. Connect the external detector to the high EXT GEN binding post through a coupling capacitor, and to ground. An external polarizing voltage may be connected to the d-c binding posts, with a high-impedance choke between the EXT GEN binding posts.

2.3 ZERO ADJUSTMENT. The following adjustments have been made at the factory and in general need be checked only occasionally. For the best possible accuracy, especially for the measurement of small capaci-



Figure 3. Connection Diagram for External Audio Oscillators.

tance, or to balance out lead capacitance, it is well to check the zero adjustment fairly often.

a. Set the MULTIPLY CAPACITANCE BY switch to 0.0001, the CAPACITANCE dial to zero, and the SEN-SITIVITY control fully clockwise. (Although capacitance less than 1000  $\mu\mu$ f should normally be measured on the 0.0001 multiplier, it is sometimes more convenient to use the 0.001 range.)

b. Set the DIR-STANDBY-REV switch to DIR, and adjust the two screw-driver controls under the snap buttons at the upper left of the panel for a null (maximum angle of the dark sector in the null detector). The left-hand control is for dissipation factor, the right-hand control for capacitance.

c. Set the DIR-STANDBY-REV switch to REV, and adjust similarly the two screw-driver controls under the snap buttons at the upper right of the panel. The bridge is now ready for measurements.

2.4 LEAD CAPACITANCE. Lead capacitance may be included in the zero adjustment if desired. The procedure is as follows:

a. Before making the zero adjustments described in paragraph 2.3, connect the unknown capacitance to the UNKNOWN terminals with the high lead disconnected at the capacitor terminal. A capacitance of about 15  $\mu\mu$ f can be compensated in this manner.

b. If the lead capacitance is greater than 15  $\mu\mu$ f, as for long leads or a cable, its effects on the measured values of capacitance and dissipation factor can be eliminated by measurement of cable or lead alone. Let the capacitance and dissipation factor of the lead or cable alone be C<sub>E</sub> and D<sub>E</sub>. Then the unknown capacitance and dissipation factor can be expressed by the following formulas:

$$C_{XS} = C - C_E$$
(1)

$$D_{X} = \frac{DC - D_{E}C_{E}}{C - C_{E}}$$
(2)

where C and D are the values with the unknown connected. 2.5 MEASUREMENT OF UNKNOWN CAPACITANCE AND DISSIPATION FACTOR.

a. Set the FILTER switch to the position corresponding to the test frequency (EXT for frequencies other than 60 or 120 cps).

b. Set the GEN switch to INT or EXT depending on whether the internal or an external test voltage supply is used.

c. Set the DIR-STANDBY-REV switch to STANDBY.

d. Connect the unknown capacitance to the UNKNOWN terminals, with the low or ground terminal of the unknown connected to the bridge G terminal.

e. Set the DIR-STANDBY-REV switch to DIR.

f. Set the SENSITIVITY control so that the eye just closes without overlap.

g. If the nominal capacitance value is known:

(1) Set the MULTIPLY CAPACITANCE BY to a range that will include the nominal value.

(2) Adjust the CAPACITANCE and DISSIPATION FACTOR controls for a null, changing the SENSITIVITY setting as necessary.

h. If the nominal capacitance value is not known:

(1) Set the MULTIPLY CAPACITANCE BY switch to 1000.

(2) Adjust the CAPACITANCE control for a null, changing the SENSITIVITY setting as necessary.

(3) If the null point is less than 1.0 on the dial, turn the MULTIPLY CAPACITANCE BY switch to 100, 10, etc, till the null is within the main decade (1.0 to 10) of the dial.

i. To eliminate the effect of extraneous voltages that may be induced (usually electrostatically) from adjacent electromagnetic fields at the power-line frequency:

(1) Obtain readings with the DIR-STANDBY-REV switch in both DIR and REV positions. The true result is the average of the two readings if the error is less than 10 percent.

(2) For strong interference causing large differences in readings, apply the following formulas:

$$D_{X} = \frac{D_{1}C_{1} + D_{2}C_{2}}{C_{1} + C_{2}}$$
(3)

$$C_{\rm XS} = \frac{C_1 + C_2}{2} \tag{4}$$

where  $C_1$  and  $D_1$  are measured at one position of the STANDBY switch;  $C_2$  and  $D_2$  are measured in the other switch position.

j. Capacitance is indicated directly by the CAPACI-TANCE dial and multiplier switch to an accuracy of  $\pm(1\%$ +1  $\mu\mu$ f) over the entire range from 1  $\mu\mu$ f to 11,000  $\mu$ f.

k. Dissipation factor is indicated in percent on the DISSIPATION FACTOR dial, and this reading should be multiplied by  $\frac{f}{60}$  where f is the frequency of the test voltage. In all unmodified formulas the dissipation factor must be expressed as a ratio (divide reading by 100). The accuracy of dissipation-factor measurement is  $\pm(2\%)$  of dial reading  $\pm\frac{f}{60} \times .05\%$  dissipation factor) over the range from 0.05% to 60% for capacitances greater than

100  $\mu\mu f$ . Below 100  $\mu\mu f$ , the decreasing sensitivity of balance increases the error inversely with capacitance.

1. Values of dissipation factor beyond the directreading range (60%) will occasionally be encountered, especially for electrolytic capacitors of very high value. The most convenient way of extending the range is by the use of an external parallel variable resistor between the SHIELD terminal and ground, with the DISSIPATION FACTOR control set at zero and the bridge balanced by the added resistance and CAPACITANCE dial. The dissipation factor of the unknown capacitor is given by:

$$D_{X} = \frac{1}{R_{NP}\omega C_{N}}$$
(5)

where  $D_X$  is in absolute value

R<sub>ND</sub> is the added parallel resistance, in ohms

 $\omega$  is  $2\pi$  times test frequency, in cycles per second C<sub>N</sub> is standard capacitance from Table I, in farads

For a test frequency of 60 cycles, formula (5) becomes:

$$D_{\rm X} = \frac{265.3}{R_{\rm NP}C_{\rm N}} \tag{6}$$

with  $D_x$  in %,  $R_{NP}$  in kilohms, and  $C_N$  in microfarads.

For a test frequency of 120 cycles, formula (5) becomes:

$$D_{\rm X} = \frac{132.6}{R_{\rm NP}C_{\rm N}} \tag{7}$$

with D<sub>x</sub> in %, R<sub>NP</sub> in kilohms, and C<sub>N</sub> in microfarads.

The measurement of a dissipation factor of 60% requires 442 kilohms on the four lower multipliers and 4.42 kilohms on the four upper multipliers. The capacitance thus measured is the parallel capacitance. The relation between series and parallel capacitance is given in paragraph 2.8.3 and in Figure 4.

2.6 MEASUREMENT OF RESISTANCE. Since the bridge measures series capacitance, the series resistance of the unknown capacitors is given by the formula:

$$R_{XS} = \frac{D_X}{\omega C_{XS}}$$
(8)

where R<sub>vs</sub> is series resistance, in ohms

- $D_X$  is dissipation factor, expressed as a ratio (or absolute value)
- $\omega$  is  $2\pi$  times test frequency, in cycles per second
- Cys is series capacitance, in farads.

For a test frequency of 60 cycles, formula (8) becomes:

$$R_{XS} = \frac{26.53 D_X}{C_{XS}}$$
(9)

with D<sub>x</sub> in % and C<sub>xs</sub> in microfarads.

For a test frequency of 120 cycles, formula (8) becomes:

$$R_{XS} = \frac{13.26 D_X}{C_{XS}}$$
(10)

with D<sub>y</sub>in % and C<sub>ys</sub> in microfarads.

Figure 6 can be used to determine series resistance for any pair of values of series capacitance and dissipation factor. Use the lower chart to determine the order of magnitude, and the upper chart to obtain significant figures.

Series resistance can be converted to parallel resistance by means of Figure 5.

# 2.7 MEASUREMENT OF ELECTROLYTIC CAPACITOR WITH POLARIZING VOLTAGE.

2.7.1 GENERAL. An external polarizing voltage up to a maximum of 500 volts can be applied to the capacitor under test on the four upper multipliers (1 $\mu$ f to 11,000 $\mu$ f). Connect the polarizing voltage to the D-C binding posts, observing the proper polarity. The impedance of the power source is placed in series with the a-c voltage supply of the bridge, and reduces both applied voltage and sensitivity. An adequate bypass capacitor across the D-C binding posts will restore sensitivity.

Measurement procedure is the same as that described in paragraph 2.4.

#### CAUTION

To protect against the short-circuit current drawn by a defective capacitor, use a protective resistor of about 50 ohms per volt in series with the d-c supply.

2.7.2 LEAKAGE CURRENT. The measurement of leakage current of electrolytic capacitors is often required in conjunction with the capacitance measurement. The addition of a milliammeter in series with the polarizing voltage source permits this leakage current to be measured at the same time the bridge measurement is made.

2.7.3 CONTROL OF TEST VOLTAGE. In certain measurements, the maximum allowable a-c test voltage may be specified. The applied voltage can be controlled by a rheostat inserted across the D-C terminals or in series with the external d-c power supply, if one is used.

2.7.4 CONVERSION OF DISSIPATION FACTOR TO P. In Electronic Industries Association (formerly RETMA) Standard RS-154A (December 1957) dealing with polarized











a Function of Series Capacitance and Dissipation Factor. Use lower chart to determine order of magnitude, and upper chart for significant figures.



1000 10000

10

μf

I

100

100

10

ł DX

.

01

1

10

րու

100

.001

.01

.

dry electrolytic capacitors for general use, maximum values of P are prescribed for various rated d-c voltages, where P is defined as rated capacitance in microfarads times equivalent series resistance in ohms, at 120 cycles. To determine the value P, simply multiply the dissipation factor (in percent) of the capacitor at 120 cps by 13.26.

#### 2.8 MEASUREMENTS OF THREE-TERMINAL CAPACI-TANCE.

2.8.1 THREE-TERMINAL NETWORKS. Figure 7 shows two methods by which the bridge can be connected to measure one or all the capacitances of a three-terminal network. In Figure 7a, the third terminal is connected to the bridge SHIELD terminal to measure the capacitance or dissipation factor of  $C_{13}$ . In this connection, the capacitance  $C_{12}$  is placed diagonally across the bridge where it has no effect on the measurement,  $C_{23}$ , on the other hand, goes directly across the standard arm, and will cause both capacitance and dissipation factor to read low by an amount determined by the ratios of  $C_{23}$  and  $D_{23}$  to the capacitance and dissipation factor of the standard arm. The errors are usually negligible if  $C_{23}$  is less than 100  $\mu\mu$ f.

An alternate method, shown in Figure 7b, takes advantage of the fact that the fourth corner of the bridge is also available on the panel, at the high EXT GEN terminal. With the third terminal of the unknown network connected to the EXT GEN terminal as shown in Figure 7b, the capacitance  $C_{23}$  is placed diagonally across the bridge and does not affect the measurement. Now, however, the capacitance  $C_{12}$  is across the B arm of the bridge. In this position it has little effect on the capacitance balance, but it causes the dissipation factor to read high by an amount  $R_B\omega C_{12}$ .

2.8.2 THE THREE-MEASUREMENT METHOD. To determine accurately all three capacitances of the networks of Figure 7, as well as associated dissipation factors, it is necessary to make three measurements, as follows:

a. Measure  $C_{12}$  and  $C_{23}$  in parallel, by grounding point 1 and measuring from point 2 to ground.



b. Measure  $C_{12}$  and  $C_{13}$  in parallel, by grounding point 2 and measuring from point 1 to ground.

c. Measure  $C_{13}$  and  $C_{23}$  in parallel, by connecting points 1 and 2 together.

d. Using the data from step a and the following formulas, compute  $C_A$  and  $D_A$ .

$$C_{A} = C_{12} + C_{23}$$
(11)

$$D_{A} = \frac{D_{12}C_{12} + D_{23}C_{23}}{C_{12} + C_{23}}$$
(12)

e. Using the data from step b and the following formulas, compute  $C_B$  and  $D_B$ .

$$C_{\rm B} = C_{12} + C_{13} \tag{13}$$

$$D_{B} = \frac{D_{12}C_{12} + D_{13}C_{13}}{C_{12} + C_{13}}$$
(14)

f. Using the data from step c and the following formulas, compute  $C_C$  and  $D_C$ .

$$C_{C} = C_{13} + C_{23}$$
 (15)

$$D_{C} = \frac{D_{13}C_{13} + D_{23}C_{23}}{C_{13} + C_{23}}$$
(16)

In solving the above simultaneous equations, the following expressions are useful:

$$C = \frac{1}{2}(C_{A} + C_{B} + C_{C})$$
(17)

$$DC = \frac{1}{2}(D_{A}C_{A} + D_{B}C_{B} + D_{C}C_{C})$$
(18)

In terms of these expressions, the three unknown capacitances and dissipation factors can be written as follows, in a form convenient for calculation:



Figure 7. Equivalent Representation of a Three-Terminal Capacitance, and Two Methods of Connecting to Bridge.

$$C_{13} = C - C_{A}$$

$$C_{23} = C - C_{B}$$

$$C_{12} = C - C_{C}$$

$$D_{13} = \frac{DC - D_{A}C_{A}}{C_{13}}$$

$$D_{23} = \frac{DC - D_{B}C_{B}}{C_{23}}$$

$$D_{12} = \frac{DC - D_{C}C_{C}}{C_{12}}$$
(19)

A convenient tabular form of computation can be used, as illustrated by the following numerical example:

	D	С	DC
CA	3.0	600	1800
C <sub>B</sub>	4.0	500	2000
c	3.5	700	2450
Sum		1800	6250
<sup>1</sup> / <sub>2</sub> Sum		900	3125
C <sub>13</sub>	4.41	300	1325
C <sub>23</sub>	2.81	400	1125
C <sub>12</sub>	3.38	200	675

2.8.3 EFFECTS OF LARGE DISSIPATION FACTORS. The formulas of paragraph 2.8.2 are based on the assumption that the measured values are parallel values of capacitance. For values of dissipation factor below 0.10(10%) the difference between series and parallel components is insignificant. For larger values of dissipation factor, the measured series values of capacitance must be converted to parallel values before the formulas can be used. The conversion formula is:

$$C_{XP} = \frac{C_{XS}}{1 + D^2}$$
(20)

where  $C_{XP}$  is parallel capacitance

C<sub>XS</sub> is series capacitance

D is dissipation factor (absolute value)

A chart by means of which parallel capacitance can be determined for any pair of values of series capacitance and dissipation factor is given in Figure 4.

#### 2.9 INDUCTANCE MEASUREMENT.

2.9.1 SERIES-SUBSTITUTION METHOD. Inductance can be determined from the measurement of the net effective

capacitance of the unknown inductor in series with a known capacitor of suitable value. The series capacitor must be small enough so that the net reactance of the combination is capacitive. At the same time, it must be large enough so that a significant change in effective capacitance results. The procedure is as follows:

a. Connect the unknown inductor and series capacitor as shown in Figure 8.

b. Short-circuit the inductor and measure the capacitor alone. Use  $C_1$  and  $D_1$  to denote the bridge readings.

c. Remove the short circuit and rebalance the bridge. Use  $C_2$  and  $D_2$  to denote the new readings.

d. Compute series inductance as follows:

$$L_{\rm XS} = \frac{\Delta C}{\omega^2 C_1 C_2} \tag{21}$$

where  $L_{XS}$  is the unknown series inductance in henrys

 $\omega$  is  $2\pi$  times frequency in cycles per second

 $C_1$ ,  $C_2$  are as given in b and c above, in farads  $\Delta C = C_2 - C_1$ 

For 60 cycles:

$$L_{\rm XS} = 7.04 \frac{\Delta C}{C_1 C_2}$$
 (22)

where  $\triangle C$ ,  $C_1$ , and  $C_2$  are in  $\mu f$ .

For 120 cycles:

$$L_{\rm XS} = 1.76 \frac{\Delta C}{C_1 C_2}$$
 (23)

where  $\Delta C$ ,  $C_1$ , and  $C_2$  are in  $\mu f$ .

e. Compute series resistance as follows:

$$R_{XS} = \frac{D_2 C_1 - D_1 C_2}{\omega C_1 C_2}$$
(24)



Figure 8. Connection Diagram for Series-Substitution Measurement of Inductance.

where R<sub>vs</sub> is series resistance, in ohms

 $\omega$  is  $2\pi$  times frequency in cycles per second

 $C_1$ ,  $C_2$ ,  $D_1$ ,  $D_2$  are from b and c above, in farads and absolute value (or ratio).

For 60 cycles:

$$R_{xs} = 26.53 \frac{D_2 C_1 - D_1 C_2}{C_1 C_2}$$
(25)

where D<sub>1</sub> and D<sub>2</sub> are in percent

 $C_1$  and  $C_2$  are in  $\mu f$ .

For 120 cycles:

$$R_{XS} = 13.27 \frac{D_2 C_1 - D_1 C_2}{C_1 C_2}$$
(26)

where D<sub>1</sub> and D<sub>2</sub> are in percent

C<sub>1</sub> and C<sub>2</sub> are in  $\mu$ f.

f. Compute storage factor, Q, as follows:

$$Q_{X} = \frac{100 \Delta C}{D_{2}C_{1} - D_{1}C_{2}}$$
(27)

where Q<sub>X</sub> is storage factor

 $\Delta C = C_2 - C_1$ C<sub>1</sub>, C<sub>2</sub>, D<sub>1</sub>, D<sub>2</sub> are from b and c above, in farads and percent.

2.9.2 PARALLEL-SUBSTITUTION METHOD. Inductance can also be measured by a parallel-substitution method, with the conditions regarding the size of the parallel capacitor the reverse of those given for the series capacitor in paragraph 2.9.1. However, since the bridge measures series components, the parallel-substitution formulas are rather complicated. The easiest way of handling the computation is probably to convert bridge readings to parallel capacitance by means of equation (20). The procedure is as follows:

a. Connect the unknown inductor and parallel capacitor as shown in Figure 9.



Figure 9. Connection Diagram for Parallel-Substitution Measurement of Inductance.

b. Disconnect the high lead of the inductor and balance the bridge. Use  $C_1$  and  $D_1$  to denote the capacitance and dissipation factor,

c. Reconnect the inductor and rebalance the bridge. Use  $C_2$  and  $D_2$  to denote the new readings.

d. Convert  $C_1$  and  $C_2$  to effective parallel values by using equation (20), and designate the new values  $C'_1$  and  $C'_2$ .

e. Compute parallel inductance as follows:

$$L_{XP} = \frac{1}{\omega^2 \Delta C'}$$
(28)

where L<sub>XP</sub> is parallel inductance, in henrys

 $\omega$  is  $2\pi$  times frequency in cycles per second

$$\triangle C'$$
 is  $C'_2 - C'_1$ , in farads.

For 60 cycles:

$$L_{XP} = \frac{7.04}{\Delta C'}$$
(29)

where  $\Delta C'$  is in  $\mu f$ .

For 120 cycles:

$$L_{XP} = \frac{1.76}{\Delta C'}$$
(30)

where  $\Delta C'$  is in  $\mu f$ .

f. Compute parallel resistance as follows:

$$R_{XP} = \frac{1}{\omega(D_2C_2' - D_1C_1')}$$
(31)

where  $R_{XP}$  is parallel resistance, in ohms

 $\omega$  is  $2\pi$  times frequency in cycles per second  $C'_1$  and  $C'_2$  are as given in d above, in farads  $D_1$  and  $D_2$  are as given in b and c above, in absolute value.

For 60 cycles:

$$R_{XP} = \frac{265,300}{D_2C_2' - D_1C_1'}$$
(32)

where D<sub>1</sub> and D<sub>2</sub> are in percent

 $C_1'$  and  $C_2'$  are in  $\mu f$ .

For 120 cycles:

$$R_{XP} = \frac{132,700}{D_2 C_2' - D_1 C_1'}$$
(33)

where  $D_1$  and  $D_2$  are in percent  $C_1'$  and  $C_2'$  are in  $\mu$ f.

g. Compute storage factor, Q, as follows:

$$Q_{X} = \frac{1}{D_{X}} = \frac{100\Delta C'}{D_{X}(D_{2}C_{2}' - D_{1}C_{1}')}$$
(34)

where  $Q_x$  is storage factor

 $\Delta C'$  is  $C'_{2} - C'_{1}$ , in farads

C<sub>1</sub> and C<sub>2</sub> are as given in d above, in farads

D<sub>1</sub> and D<sub>2</sub> are as given in b and c above, in absolute value.

h. Compute series inductance as follows:

$$L_{XS} = \frac{L_{XP}}{1 + D_{y}^{2}}$$
(35)

where  $L_{YP}$  is series inductance

 $L_{YP}$  is parallel inductance

$$D_X \text{ is } \frac{1}{Q_X}$$
 from (34).

i. Compute series resistance as follows:

$$R_{XS} = \frac{R_{XP}}{1 + Q_X^2}$$
(36)

where R<sub>XS</sub> is series resistance

 $R_{\chi p}$  is parallel resistance

 $Q_x$  is storage factor from (34).

#### 2.10 RESISTANCE MEASUREMENT.

2.10.1 SERIES-SUBSTITUTION METHOD. To measure the effective a-cresistance of a resistor, connect a capacitor of suitable size in series with the resistor and measure the dissipation factor. The difference between this value and the dissipation factor of the capacitor alone can be used to compute resistance. The procedure is as follows:

a. Connect the capacitor and resistor as shown in Figure 8, replacing the inductor with a resistor.

b. Short-circuit the resistor and measure the capacitor alone. Use  $C_1$  and  $D_1$  to denote the capacitance and dissipation factor.

c. Remove the short circuit and rebalance the bridge. Use  $C_2$  and  $D_2$  to denote the new values.  $C_1$  and  $C_2$  will be equal unless the resistor has a reactive component.

d. Compute the unknown resistance as follows:

$$R_{XS} = \frac{\Delta D}{\omega C}$$
(37)

where R<sub>XS</sub> is series resistance, in ohms

D is D<sub>2</sub> - D<sub>1</sub>, absolute value

 $\omega$  is  $2\pi$  times frequency in cycles per second

 $C = C_1 = C_2$  in farads.

For 60 cycles:

$$R_{XS} = 26.53 \frac{\Delta D}{C}$$
(38)

where  $\Delta D$  is in percent

C is in  $\mu f$ .

For 120 cycles:

$$R_{xs} = 13.27 \frac{\Delta D}{C}$$
(39)

where  $\Delta D$  is in percent

C is in  $\mu f$ .

Figure 5 shows the relation between series and parallel resistance as a function of dissipation factor.

2.10.2 PARALLEL-SUBSTITUTION METHOD. The procedure described in paragraph 2.9.2 may be used with the resistor in parallel with the known capacitor. The procedure is as follows:

a. Connect the resistor and capacitor as shown in Figure 9, replacing the inductor with a resistor.

b. Disconnect the high lead of the resistor and balance the bridge. Use  $C_1$  and  $D_1$  to denote the capacitance and dissipation factor.

c. Reconnect the resistor and rebalance the bridge, using  $C_2$  and  $D_2$  to denote the new values.

d. Convert  $C_1$  and  $C_2$  to effective parallel values by using equation (20) or Figure 4.  $C_1$  and  $C_2$  will be equal unless the resistor has a reactive component. Designate the new values  $C'_1$  and  $C'_2$ .

e. Compute the unknown resistance as follows:

$$R_{XP} = \frac{1}{\omega C' \Delta D}$$
(40)

where  $R_{XP}$  is parallel resistance, in ohms

 $\omega$  is  $2\pi$  times frequency in cycles per second

C' is as given in d above, in farads

 $\Delta D$  is  $D_2 - D_1$ , absolute value.

For 60 cycles:

$$R_{XP} = \frac{265,300}{C' \Delta D}$$
(41)

where C' is in  $\mu f$ 

 $\Delta D$  is in percent.

For 120 cycles:

$$R_{XP} = \frac{132,700}{C'\Delta D}$$
(42)

where C' is in  $\mu f$ 

 $\Delta D$  is in percent.

Figure 5 shows the relation between series and parallel resistance as a function of dissipation factor.

#### 2.11 MEASUREMENT OF DIELECTRIC SAMPLES.

2.11.1 GENERAL. The dielectric constant and dissipation factor of an insulating material can be determined from the measurement of the capacitance and dissipation factor of an elementary capacitor, with the material used as the insulating medium between metallic electrodes of suitable dimensions.

2.11.2 TWO-ELECTRODE METHOD. A simple two-electrode method is sufficient for most purposes. The procedure is as follows:

a. If possible, choose a sample of such shape and dimension as to yield a capacitance of 100  $\mu\mu$ f or more. The calculation of dielectric constant is simplified if the thickness and area are easily measured and calculated, such as a disk or rectangle. If measurements are to be made at various frequencies, it is best to use sizes and shapes as specified in ASTMD-150 (available from American Society for Testing Materials, 260 Race Street, Philadelphia, Pennsylvania).

b. Measure and record the dimensions of the sample, and clean it thoroughly. (A mixture of half grain alcohol and half ether is recommended unless it is a solvent for the material.)

c. When the sample is dry, apply a very thin film of refined petrolatum to one surface. Place a thin metal foil electrode, preferably less than 1 mil thick, and larger than the sample, on this surface.

d. Press the electrode in place with a pad of cloth or squeegee roller and rub out any bubbles, so that the foil is in intimate contact with the surface. Then trim the foil to the same size as the sample.

e. Apply the other electrode to the sample as described in steps c and d.

#### NOTE

An alternate method of forming electrodes is to brush a good silver paint (such as Dupont No. 4132 Silver Paste) on the sample and to dry it overnight at 60°C. Such an electrode is porous to moisture, so that the dielectric can be conditioned at any desired relative humidity without removing the electrode.

f. Measure capacitance as described in paragraph 2.4.

g. Compute dielectric constant (to a first approximation) as follows:

$$K = \frac{4.45 t C}{A}$$
(43)

where K is dielectric constant

t is thickness of the sample, in inches

C is measured capacitance, in  $\mu\mu f$ 

A is area of the electrodes, in square inches.

For a complete discussion of the effects of stray electric field at the edges of the electrodes, and the effect of the capacitance of the high electrode to ground, refer to ASTM D-150.

2.11.3 THREE-ELECTRODE (GUARD-ELECTRODE) METHOD. The guard arrangement, shown in Figure 10, is electrically equivalent to the three-terminal capacitance discussed in paragraph 2.7.1 and shown in Figure 7. To measure the characteristics of different parts of a nonhomogenous sample, place the electrodes over the part to be measured and block off the remainder by means of the guard electrode. If the guard electrode is connected to the SHIELD terminal, as described in paragraph 2.7.1, the guard-to-ground capacitance is placed across the capacitance arm of the bridge, where it will generally have a negligible effect on the measurement.

2.12 MEASUREMENT OF LIQUID INSULATION. Liquid insulation, such as transformer oils and askarels, requires some type of cell for measurement of capacitance and dissipation factor. The cell in its simplest form can be a multiple-plate air capacitor immersed in the liquid, or a grounded cylindrical can with a slightly smaller insulated cylindrical electrode. Such cells do not allow the accurate calculation of dielectric constant, nor do they maintain a constant voltage gradient on the liquid. These difficulties are overcome by the use of a threeelectrode cell, such as described in ASTM D-150. Such a cell is electrically equivalent to Figure 10, with a ter-



Figure 10. Guard-Electrode Arrangement for Measuring Dielectric Samples.

minal capacitance small enough to permit use of the SHIELD terminal as described in paragraph 2.7.1.

2.13 MEASUREMENT OF INSULATORS. Electrical insulators, such as bushings, may be measured as described in paragraph 2.5, since such insulators have one terminal grounded. To eliminate the effect of leakage over the surface of the insulation, place an electrode around the insulator surface and connect this electrode to the SHIELD terminal, as described in paragraph 2.8.1.

2.14 TRANSFORMER MEASUREMENTS. The insulation in a transformer, together with the primary and secondary windings and the transformer case, form a three-terminal network (refer to paragraph 2.8). All the direct capacitances are generally so large that the three methods described in paragraph 2.8.2 must be used to separate the three insulations: primary to case, secondary to case, and primary to secondary.

A sample of the oil in a transformer can be drawn off and tested by the method described in paragraph 2.12. 2.15 MEASUREMENT OF PARALLEL CAPACITANCE. Effective parallel capacitance can be measured by the use of an external decade-resistance box or other variable resistor to balance for dissipation factor. Connect the external resistor from SHIELD to ground and set the DISSIPATION FACTOR control at zero. Balance the bridge by means of the CAPACITANCE dial and the external resistance. Used in this way, the CAPACITANCE dial reads directly the parallel capacitance of the unknown. Compute parallel resistance as follows:

$$R_{XP} = R_{NP} \frac{R_B}{R_A}$$
(44)

$$R_{XP} = R_{NP} \frac{C_N}{C_{XP}}$$
(45)

where  $R_{yp}$  is the effective parallel resistance, in ohms

R<sub>NP</sub> is the external resistance, in ohms

 $R_{A}$  in k $\Omega$  = the CAPACITANCE DIAL reading in  $\mu f$ 

 $R_B$  is the resistance of  $R_B$  from Table 1 in k $\Omega$ 

 $\boldsymbol{C}_{N}$  is the value of the standard capacitor, from Table 1

 $C_{yp}$  equals bridge reading in  $\mu f$ .

2.16 USE AS A LIMIT BRIDGE. Because of the limiting characteristics of the null-detector response, the bridge is not as readily suitable for limit work as one having a linear detector response. Such linear response can be obtained by the removal of the Thyrite resistor across the grid of the 6E5 tube.

T	A	В	L	Ε	1
---	---	---	---	---	---

MULTIPLIER	R <sub>B</sub>	I <sub>M</sub>	C <sub>N</sub>	VOLTS APPLIED TO UNKNOWN WITH C DIAL AT	
	RESISTANCE	MAX CURRENT		MIN	MAX
0.0001	100,000 ohms	1.5 ma	0.01 <i>µ</i> f	128	121
0.001	10,000 ohms	5 ma	$0.01 \mu f$	57	43
0.01	1,000 ohms	16 ma	$0.01 \mu f$	13	5
0.1	100 ohms	55 ma	$0.01 \mu f$	1.5	0.5
1	1,000 ohms	16 ma	1.0 $\mu$ f	22	5
10	100 ohms	55 ma	$1.0 \ \mu f$	12	1.5
100	10 ohms	170 ma	$1.0 \ \mu f$	2.4	0.5
1000	1 ohm	550 ma	1.0 $\mu$ f	1.4	0.15

 $R_A$  resistance in k $\Omega$  equals dial reading in  $\mu f$ 

# Section 3 PRINCIPLES OF OPERATION

3.1 BRIDGE CIRCUIT. In order to cover an extremely wide capacitance range effectively, two series-resistance bridge circuits are used, differing in the value of standard capacitance used and in the method of connecting the generator and detector. All necessary changes in circuit connections are made by the panel controls.

Figure 2a shows the basic bridge circuit used for the four lower multiplier positions (0.0001, 0.001, 0.01, 0.1). Figure 2b shows the arrangement used for the four high multipliers (1, 10, 100, 1000), where generator and detector are interchanged from the circuit of Figure 2a and the capacitance and resistance of the standard arm have been changed by a factor of 100.

The balance conditions are independent of generator and detector connections. They are given by:

$$C_{XS} = C_N \frac{R_A}{R_B}$$
(46)

$$D_{x} = R_{x}\omega C_{x} = R_{N}\omega C_{N} = D_{N}$$
(47)

where  $C_{XS}$  and  $D_X$  are the series capacitance and dissipation factor of the unknown capacitor and  $C_N$  and  $D_N$  are the series capacitance and dissipation factor of the standard arm.

With "high multiplier" connection of Figure 2b, there is provision for connecting a d-c voltage in series with the a-c bridge voltage for polarizing capacitors. One terminal of the d-c source is at ground potential, so that almost any available d-c source can be used. This arrangement also avoids the necessity for a high-impedance choke.

3.2 DETECTOR. The detector system consists of a single stage of amplification and an electron-ray tube ("magic eye"), which is used as a visual indicator. As shown in



Figure 11. Diagram of Detector System, Showing Resonant Circuit Used for Frequency Selectivity.

Figure 11, the amplifier is made selective to the operating frequency by means of a parallel-resonant circuit in the plate circuit. The detector system is designed to be very sensitive when the bridge is at or near balance, but relatively insensitive when the bridge is out of balance. This is accomplished by means of a Type 6E5 remote-cutoff electron-ray tube, in conjunction with a voltage-sensitive resistance element (Thyrite) shunting its input. The resulting characteristics greatly simplify the process of locating the balance, and the necessity of using the manual sensitivity control is greatly reduced.

To disconnect the detector from the bridge, it is necessary only to turn the multiplier switch to DET. The bridge detector can then be used as a null indicator for any other bridge. The input to the detector is between the high EXT GEN binding post and a G terminal.

The FILTER switch provides a choice of: 60-cycle internal filter, 120-cycle internal, or untuned. With the switch in the EXT (untuned) position, an external filter (parallel-resonant tuned circuit) can be connected to the EXT FILTER jack for measurements at frequencies other than 60 or 120 cps. An external detector may be connected to these terminals if desired.

#### 3.3 BRIDGE VOLTAGE SUPPLY.

3.3.1 INTERNAL 60-CYCLE SUPPLY. Because of the wide impedance range presented by the bridge, no single source of voltage can efficiently supply the test voltage for all ranges. Four separate sources are used, one for each pair of the eight multiplier positions (see Figure 12).



Figure 12. Diagram of Internal Bridge-Voltage-Supply Connections.

The voltages and series resistances are chosen to deliver approximately the maximum safe power to the bridge under all conditions. The proper source is selected by a switch mechanically connected to the switch that controls the bridge ratio arm. The approximate voltages applied to the unknown are given in Table 1.

3.3.2 EXTERNAL SUPPLIES. At 120 cycles, the impedance presented by the bridge at the EXT GEN terminals is approximately 1000 ohms for 1  $\mu$ f and approximately 1 ohm for 10,000  $\mu$ f. To provide reasonable impedance match, an arrangement similar to that used in the internal 60-cycle supply is recommended. (Refer to paragraph 3.3.1.)

With the voltages and resistances as shown in Figure 3, the bridge cannot be overloaded at any setting. If a low-power audio oscillator is used as a source, the limiting resistors may not be needed. The maximum power that should be delivered to the bridge is  $\frac{1}{2}$  watt, the rating of the resistors  $R_B$ .

Most R-C oscillators commercially available are not able to deliver the desired power to the bridge on lowimpedance, high-capacitance multipliers. A Type 1206-B Unit Amplifier (with a Type 1203-B Unit Power Supply) can be used to raise the voltage to a satisfactory level.

A Type 1304-B Beat-Frequency Oscillator can be used directly without an external transformer (one is included in its output circuit). Connect the oscillator to the EXT GEN terminals and set the output to 30 volts on X1 and X10 and to 3 volts on X100 and X1000 multipliers with the unknown disconnected. Use the proper limiting resistors as shown in Figure 3.

The Type 1214-D Unit Oscillator is designed specifically as a 120-cycle source for the Type 1611-B Bridge, and requires no additional accessories. The transformer of this Unit Oscillator may be used with any oscillator to provide impedance match.

#### 3.4 COMPENSATING CIRCUITS.

3.4.1 RATIO - ARM COMPENSATION. Any loss in the standard capacitor causes the bridge to read low in dissipation factor by an amount equal to the dissipation factor of the standard capacitor. Since this value may be as much as 0.0012, it must be compensated in order to realize the rated accuracy of the bridge. Compensation is accomplished by fixed capacitors connected across the resistors in the opposite arm. These capacitors are of such values that  $R_B\omega C_B$  approximates  $D_O$ , the dissipation factor of the standard capacitor. No compensation is used on the 1- and 10-ohm resistors corresponding to the 100 and 1000 multipliers, since the dissipation factor of the standard capacitor associated with these multipliers is only 0.0003.

# Section 4 SERVICE AND MAINTENANCE

4.1 GENERAL. The two-year warranty given with every General Radio instrument attests the quality of materials and workmanship in our products. When difficulties do occur, our service engineers will assist in any way possible.

In case of difficulties that cannot be eliminated by the use of these service instructions, please write or phone our Service Department, giving full information of the trouble and of steps taken to remedy it. Be sure to mention the serial and type numbers of the instrument.

Before returning an instrument to General Radio for service, please write to our Service Department or nearest district office (see back cover), requesting a Returned Material Tag. Use of this tag will insure proper handling and identification. For instruments not covered by the warranty, a purchase order should be forwarded to avoid unnecessary delay.

4.2 ROUTINE MAINTENANCE. Switch blades, contacts, and contact surfaces of the wire-wound rheostats should periodically be cleaned and lubricated. A recommended cleaning method is to use a mixture of half ether and half alcohol, followed by a few light strokes with crocus cloth. A very thin smear of petrolatum should then be applied on the wire-wound rheostats.

4.3 DIAL CALIBRATION PROCEDURES.

4.3.1 GENERAL. The following calibration procedures require the use of a Wheatstone bridge accurate to at least  $\pm 0.25\%$ .

4.3.2 CAPACITANCE DIAL CALIBRATION.

a. Make certain that the rheostat R19 is clean, that the setscrews are tight, and that the cam follower mechanism is operating freely.

b. Remove one connection to R19, and connect the Wheatstone bridge to the rheostat.

c. Measure the resistance in ohms, using the Wheatstone bridge. The resistance should be 1000 times the setting of the CAPACITANCE dial. If it is not, proceed to step d.

d. Turn the CAPACITANCE dial to a point near its center. This will bring the cam follower to a point opposite one of the cam plate screws. Adjust this screw until the resistance in ohms is 1000 times the dial reading.

e. Progressing up and then down the scale, adjust each cam plate screw in the same manner. As a final check, see that each main point on the dial as well as points 0.5 and 0.1 check within 0.5% of the correct resistance value. If it is found that the whole cam plate is adjusted too high or too low, loosen the set screws in the dial, shift its position on the shaft, and start the cam plate adjustment over again.

4.3.3 DISSIPATION FACTOR DIAL CALIBRATION.

a. Make certain that the blade on each rheostat passes from the last wire to the terminal within one-half division of the dial ZERO. If necessary, shift the dial and/or blades to meet this condition. It is especially important that the ZERO of the high-resistance rheostat (R31A, B) match the dial.

b. Connect the Wheatstone bridge between ground and the junction of R26A, R30, C13, and C14. Set the DISSI-PATION FACTOR dial to 5% and adjust R30 to obtain exactly 132.6 ohms.

c. Set the dial to 30% and adjust R27 to obtain exactly 795.8 ohms.

d. Connect the Wheatstone bridge between ground and the junction of R28, R31A, and C12.

e. Set the DISSIPATION FACTOR dial to 5% and adjust R28 to obtain 13,500 ohms.

f. Set the dial to 30% and adjust R29 to obtain 81,100 ohms.

#### 4.4 TROUBLE - SHOOTING PROCEDURE.

4.4.1 GENERAL. If the bridge fails to function, check the positions of the bridge controls and the condition of the sample being measured before looking for trouble in the bridge. Failure to obtain a balance may well be caused by reasons outside the bridge, e.g. a short-circuited sample. Make certain, if possible, that the impedance of the sample under test does not exceed in magnitude the limits of the bridge.

Before removing the bridge from its case, check that the voltage and frequency of the power line are correct, test the power cord for open circuits or for poor contacts in the power outlets, and check the fuses mounted on top of the panel both for continuity and for tightness.

4.4.2 DETAILED TROUBLE - SHOOTING PROCEDURE.

4.4.2.1 Indicating Eye Does Not Light.

- a. Check tube V1 and voltages per Table 2.
- b. Check operation of switch S1.

c. Check transformer T1 for open windings and proper connections.

d. Check resistor R1 for open or short circuit and for proper value.

e. Check capacitor C1 for open or short circuit and for leakage.

- 4.4.2.2 Eye Does Not Deflect Properly.
  - a. Refer to paragraph 4.4.2.1.
  - b. Check V2 and voltages per Table 2.

c. Check R3 for open or short circuit and for proper operation.

- d. Check C2 for open or short circuit and for leakage.
- e. Check inductor L1 for open or short circuit.

f. Check contacts and operation of S2, second and third decks from panel.

- g. Check detector circuit gain:
  - Turn MULTIPLY CAPACITANCE BY switch to DET.
  - (2) Apply 0.5 v, 60 cps between DET terminal and chassis.
  - (3) Check that voltage from pin 3, V1 to chassis is approximately 5.0 v.

h. Check Thyrite resistor R2:

- (1) Apply between 1.1 and 3.3 v across R2.
- (2) Check that current is appoximately 1  $\mu$ a.
- (3) Apply between 3.8 and 5.2 v across R2.
- (4) Check that current is approximately 10  $\mu$ a.

4.4.2.3 No Balance On Any Position of MULTIPLY CA-PACITANCE BY Switch.

a. Check that the sample under test is not defective.

b. Check T2 for open windings and proper connections. Check that inner conductors of shielded leads, as well as sheath to SHIELD terminal and S2, are not grounded.

c. Check R19 through R31, and R34 for open or short circuit and for proper values.

d. Check capacitors C8 and C10 through C15 for open or short circuit and for leakage.

e. Check contacts and operation of S2, fourth and fifth decks from panel.

f. Check that UNKNOWN + terminal is not shorted to ground.

g. Refer to paragraphs 4.4.2.1 and 4.4.2.2.

4.4.2.4	No	Balance	on	One	or	More	Positions	of	MUL	TI-
PLY CA	PA	CITANC	ΕB	Y Sv	vit	ch.				_

a. Check T1 for open windings and proper connections.

b. Check contacts and operation of S2, deck nearest panel.

c. Check R7 through R10 for open or short circuit and proper values for capacitance range involved.

#### 4.4.2.5 Zero Adjustments Impossible.

a. Check R11 through R18, and R32 for open or short circuit and for proper values.

- b. Check C7 for open or short circuit and for leakage.
- c. Check contacts and operation of S3.
- d. Refer to paragraphs 4.4.2.1 through 4.4.2.4.

4.4.3 TEST VOLTAGES. Table 2 lists tube socket voltages from pin to ground, as measured with a 20,000-ohmper-volt meter. Line voltage should be 115 volts, 60 cps, MULTIPLY CAPACITANCE BY switch set to 1, and DIS-SIPATION FACTOR, CAPACITANCE, and SENSITIVITY controls all set to maximum. Deviations up to 20% in d-c voltages should not be considered abnormal.

TABLE 2

TUBE (TYPE)	PIN	VOLTS	TUBE (TYPE)	PIN	VOLTS
V1 (6E5)	1-6 2 3 4	6.3 ac 155* -6.6** 155*	V2 (Cont.)	5 6 8	0 22* 150
	5	0	V3 (6X5GT/G)	1	0
V2 (6SJ7)	1 2-7 3 4	0 6.3 ac 0 -3.0	(04)01/0)	2-7 3 5 8	6.3 ac 260 ac 260 ac 300

\*Use VTVM

\*\*May vary up to ±50%







Figure 14. Bottom Interior View.



Figure 15. Left Side Interior View.

### Section 5

## PARTS LIST

RESISTORS

Ref No

Part No Ref No

CAPACITORS

Part No

R2Thy POTENTIOMETER, $1 M\Omega \pm 20\%$ 6000-1300 R4C740-0300 6000-1300C2Mica, $0.035 \ \mu\text{F} \pm 10\%$ (C3A33 (C4 (C5)A33 (C4 (C5)C3C0 (C4 (C5)C3C0 (C4 (C5)C3C0 (C4 (C5)C3C0 (C4 (C5)C4 (C5)Electrolytic Block, 20 \ \mu\text{F}} (C6)A433 (C6)A433 (C6)C3C0 (C4 (C6)C3C3 (C4 (C6)C3 (C4 (C6)C3 (C4 (C6)C3 (C4 (C6)C3 (C4 (C6)C3 (C4 (C6)C3 (C4 (C6)C3 (C4 (C6)C3 (C4 (C6)C3 (C4 (C6)C3 (C4 (C6)C3 (C4 (C6)C3 (C4 (C6)C3 (C4 (C6)C3 (C4 (C6)C3 (C4 (C6)C4 (C6) (C6)C3 (C4 (C6)C3 (C4 (C6)C3 (C4 (C6)C3 (C4 (C6)C3 (C4 (C6)C4 (C6) (C4) (C6)C3 (C4 (C6)C3 (C4 (C6)C3 (C4) (C6)C4 (C6) (C6)C3 (C4) (C6)C4 (C6) (C6)C3 (C4) (C6)C4 (C6) (C6)C3 (C4) (C6)C4 (C6) (C6)C3 (C4) (C6)C4 (C6) (C6)C3 (C4) (C6)C4 (C4) (C6)C4 (C4) (C6)C4 (C4) (C6)C4 (C4) (C6)C4 (C4) (C6)C4 (C4) (C10)C4 (C10)C4 (C10)C4 (C10)C4 (C10)C4 (C10) (C11)C4 (C12) (C12)C4 (C12) (C12)C4 (C12) (C12)C4 (C12) (C12) (C13)C4 (C14) (C14) (C14) <br< th=""><th></th><th></th><th></th><th></th><th></th><th></th></br<>						
R3POTENTIOMETER, $1 M\Omega \pm 20\%$ 6000-1300C2Mica, $0.053 \mu^{2} \pm 10\%$ and $483$ R4Composition, $27 k\Omega \pm 10\%$ $6100-3275$ C3 $20 \mu^{F}$ $450 V$ $446$ R5Composition, $27 k\Omega \pm 10\%$ $6100-4565$ C4Electrolytic Elock, $20 \mu^{F}$ $450 V$ $446$ R7Composition, $2.2 k\Omega \pm 10\%$ $6120-3229$ C6 $20 \mu^{F}$ $450 V$ $446$ R6Omposition, $3.3 k\Omega \pm 10\%$ $6120-3229$ C6 $20 \mu^{F}$ $450 V$ $466$ R10Wire-wound, $3.3 \Omega \pm 10\%$ $6120-3229$ C7Mica, $0.001 \mu^{F} \pm 10\%$ $466$ R11Composition, $10 M\Omega \pm 10\%$ $6100-3275$ C1Mica, $0.005 \mu^{F} \pm 10\%$ $466$ R12Composition, $27 k\Omega \pm 10\%$ $6100-3275$ C1Mica, $0.005 \mu^{F} \pm 10\%$ $466$ R12Composition, $10 R\Omega \pm 10\%$ $6100-3275$ C13Plastic, $0.5 \mu^{F} \pm 1.5\%$ $480$ R14Composition, $10 k\Omega \pm 10\%$ $6100-3275$ C13Plastic, $0.5 \mu^{F} \pm 1.5\%$ $486$ R15POTENTIOMETER, $20 k\Omega \pm 5\%$ $973-4080$ C14Plastic, $0.5 \mu^{F} \pm 1.5\%$ $486$ R16POTENTIOMETER, $11.2 - 11.7 k\Omega$ $443-4020$ MisceLLANEOUS $710 - 75\%$ $8073-4080$ R18POTENTIOMETER, $10.2 + 01.5\%$ part of 1611-0280F1FUSE, $115 V$ , $0.2 A$ 533R21Precision, $10 k\Omega \pm 0.1\%$ part of 1611-0280F2FUSE, $115 V$ , $0.2 A$ 533R24Wire-wound, $9.0 \pm 0.1\%$ part of 1611-0280F1INDUCTOR <t< th=""><th></th><th></th><th>6100-5335</th><th>C1</th><th>Mica, 0.02 μF ±10%</th><th>4830-0500</th></t<>			6100-5335	C1	Mica, 0.02 μF ±10%	4830-0500
R4Composition, 27 k0 ±10% $6000-1300$ C320 µF433R5Composition, 27 k0 ±10% $6100-3275$ C4Electrolytic Elock, 20 µF $450 V$ $446$ R6Composition, 22 k0 ±10% $6120-2339$ C6 $20 µF$ $450 V$ $446$ R7Composition, 3.3 k0 ±10% $6120-2339$ C7Mica, 0.001 µF ±10% $466$ R8Composition, 3.3 k0 ±10% $6760-1339$ C7Mica, 0.001 µF ±10% $466$ R10Wire-wound, 330 0 ±10% $6760-1339$ C1Mica, 0.005 µF ±10% $466$ R11Composition, 10 M0 ±10% $6100-6105$ C1Mica, 0.005 µF ±10% $466$ R12Composition, 27 k0 ±10% $6100-3275$ C12Mica, 0.009 µF ±0.25% $650$ $650$ R13Composition, 27 k0 ±10% $6100-3105$ C14Plastic, 0.5 µF ±1.5% $486$ R14Composition, 10 k0 ±10% $6100-3105$ C14Plastic, 0.5 µF ±1.5% $486$ R15POTENTIOMETER, 20 k0 ±5% $0973-4080$ C15Mica, 0.0088 µF ±10% $476$ R16POTENTIOMETER, 20 k0 ±5% $0973-4080$ C16Mica, 0.0088 µF ±10% $480$ R19POTENTIOMETER, 11.2 - 11.7 k0 $0433-4020$ $115 V, 0.2 A$ $533$ R21Precision, 100 k0 ±0.1%part of 1611-0280 $F1$ FUSE, 115 V, 0.2 A $533$ R24Wire-wound, 90 $\Omega$ ±0.1%part of 1611-0280 $F2$ FUSE, 115 V, 0.2 A $533$ R24Wire-wound, 9.98 $\Omega$ ±0.2%part of 1611-0280 $F1$ PLUG, Power </td <td></td> <td></td> <td>6740-0300</td> <td>C2</td> <td>Mice 0.035 uE +10% and</td> <td>4830-0500</td>			6740-0300	C2	Mice 0.035 uE +10% and	4830-0500
R5Composition, 27 kΩ ±10%6100-3275C4C5Electrolytic Block, 20 μF450 V446R6Composition, 22 kΩ ±10%6120-3229C620 μFC620 μF466R8Composition, 3.3 kΩ ±10%6120-3229C620 μF466R9Wire-wound, 330 Ω ±10%6760-1339C8Mica, 50 pF ±10%466R10Wire-wound, 3.3 Ω ±10%6760-9339C10Mica, 0.005 μF ±10%480R11Composition, 10 MΩ ±10%6100-6105C11Mica, 0.0005 μF ±10%480R12Composition, 27 kΩ ±10%6100-3275C13Plastic, 0.5 μF ±1.5%486R14Composition, 10 kΩ ±10%6100-3105C14Plastic, 0.5 μF ±1.5%486R15POTENTIOMETER, 20 kΩ ±5%0973-4080C15Mica, 0.0088 μF ±10%476R16POTENTIOMETER, 20 kΩ ±5%0973-4080C16Mica, 0.0088 μF ±10%480R17POTENTIOMETER, 20 kΩ ±5%0973-4080C16Mica, 0.0088 μF ±10%480R18POTENTIOMETER, 10 kΩ ±0.25%0602-3052F1FUSE, 230 V, 0.1 A533R21Precision, 10 kΩ ±0.025%0602-3052F1FUSE, 230 V, 0.1 A533R21Precision, 10 kΩ ±0.7%part of 1611-0280J1JACK, EXT FILTER426R25Wire-wound, 90 Ω ±0.1%part of 1611-0280J1IJACK, EXT FILTER789R24Wire-wound, 90 Ω ±0.1%6010-2600SWITCH, Rotary Wafer,SWITCH, Rotary Wafer,533R24 <td></td> <td></td> <td>6000-1300</td> <td></td> <td><math>10^{-10}</math> and <math>10^{-10}</math></td> <td>4830-0300</td>			6000-1300		$10^{-10}$ and $10^{-10}$	4830-0300
R6Composition, 0.56 MΩ ±10%6100-4565C5Electrolytic kick, 20 µF450 V440R7Composition, 22 kΩ ±10%6120-3329C620 µF20 µF20 µFR8Composition, 3.3 kΩ ±10%6120-3339C7Mica, 0.001 µF ±10%466R9Wire-wound, 3.3 Ω ±10%6760-9339C7Mica, 0.005 µF ±10%466R10Wire-wound, 3.3 Ω ±10%6100-6105C1Mica, 0.005 µF ±10%480R11Composition, 27 kΩ ±10%6100-3275C13Plastic, 0.5 µF ±1.5%486R14Composition, 27 kΩ ±10%6100-3105C13Plastic, 0.5 µF ±1.5%486R15POTENTIOMETER, 20 kΩ ±5%0973-4080C16Mica, 0.008 µF ±10%476R16POTENTIOMETER, 20 kΩ ±5%0973-4080C16Mica, 0.008 µF ±10%476R18POTENTIOMETER, 20 kΩ ±5%0973-4080C16Mica, 0.008 µF ±10%480R19POTENTIOMETER, 11.2 - 11.7 kΩ0433-4020MiSCELLANEOUSF1FUSE, 115 V, 0.2 A533R21Precision, 10 kΩ ±0.25%0602-3052G10-1600MiSCELLANEOUSF2FUSE, 115 V, 0.2 A533R22Wire-wound, 90 Ω ±0.1%part of 1611-0280F1FUSE, 115 V, 0.2 A533R23Wire-wound, 90 Ω ±0.1%part of 1611-0280F1FUSE, 230 V, 0.1 A533R24Wire-wound, 90 Ω ±0.1%6010-2600SWITCH, Rotary Wafer,MULTIPLY CAP BY789R25POTENTIOMETER, 30 kΩ ±10%6010-2600SWITCH, Rotary		Composition, 27 k $\Omega$ ±10%	6100-3275			
R7Composition, 22 kΩ ±10%6120-3239C6 $20 \mu F$ R8Composition, 33 Ω ±10%6120-2339C7Mica, 0.001 μF ±10%466R9Wire-wound, 33 Ω ±10%6760-9339C8Mica, 50 pF ±10%466R10Wire-wound, 33 Ω ±10%6760-9339C10Mica, 0.005 μF ±10%480R11Composition, 27 kΩ ±10%6100-3275C12Mica, 0.009 μF ±0.25%050R13Composition, 27 kΩ ±10%6100-3275C13Plastic, 0.5 μF ±1.5%486R14Composition, 10 kΩ ±10%6100-3105C14Plastic, 0.5 μF ±1.5%486R15POTENTIOMETER, 20 kΩ ±5%0973-4080C15Mica, 0.0088 μF ±10%480R17POTENTIOMETER, 20 kΩ ±5%0973-4080C16Mica, 0.0088 μF ±10%480R17POTENTIOMETER, 20 kΩ ±5%0973-4080C16Mica, 0.0088 μF ±10%480R17POTENTIOMETER, 20 kΩ ±5%0973-4080C16Mica, 0.0088 μF ±10%480R18POTENTIOMETER, 11.2 - 11.7 kΩ0433-4020MisCeLLANEOUSF1FUSE, 230 V, 0.1 A533R21Precision, 10 kΩ ±0.1%part of 1611-0280J1JACK, EXT FILTER426R22Wire-wound, 90 Ω ±0.1%part of 1611-0280J1JACK, EXT FILTER426R24Wire-wound, 90 Ω ±0.1%part of 1611-0280J1JACK, EXT FILTER426R25Wire-wound, 90 R ±0.1%6010-1200S3SWITCH, Rotary Wafer,S3R24POTENTIOMETER, 30 kΩ ±10%601	R5	Composition, 27 k $\Omega$ ±10%	6100-3275		Electrolytic Block 20 µF 450 V	4460-0900
R8Composition, $3.3 k\Omega \pm 10\%$ $6120-2339$ C7Mica, $0.001 \ \mu\text{F} \pm 10\%$ $466$ R9Wire-wound, $3.3 \Omega \pm 10\%$ $6760-1339$ C8Mica, $50 \ p\text{F} \pm 10\%$ $466$ R10Wire-wound, $3.3 \Omega \pm 10\%$ $6760-9339$ C10Mica, $0.005 \ \mu\text{F} \pm 10\%$ $460$ R11Composition, $10 \ M\Omega \pm 10\%$ $6100-3275$ C12Mica, $0.005 \ \mu\text{F} \pm 10\%$ $480$ R12Composition, $27 \ k\Omega \pm 10\%$ $6100-3275$ C12Mica, $0.009 \ \mu\text{F} \pm 1.2\%$ $456$ R14Composition, $10 \ k\Omega \pm 10\%$ $6100-3275$ C12Mica, $0.009 \ \mu\text{F} \pm 1.2\%$ $456$ R14Composition, $10 \ k\Omega \pm 10\%$ $6100-3275$ C12Mica, $0.009 \ \mu\text{F} \pm 1.2\%$ $456$ R15POTENTIOMETER, $20 \ k\Omega \pm 5\%$ $0973-4080$ C14Plastic, $0.5 \ \mu\text{F} \pm 1.5\%$ $486$ R16POTENTIOMETER, $20 \ k\Omega \pm 5\%$ $0973-4080$ C15Mica, $0.0088 \ \mu\text{F} \pm 10\%$ $476$ R17POTENTIOMETER, $20 \ k\Omega \pm 5\%$ $0973-4080$ C16Mica, $0.0088 \ \mu\text{F} \pm 10\%$ $480$ R19POTENTIOMETER, $20 \ k\Omega \pm 5\%$ $0973-4080$ R15FUFUSE, $230 \ V$ , $0.1 \ A$ $533$ R21Precision, $10 \ k\Omega \pm 0.1\%$ part of 1611-0280F1FUSE, $230 \ V$ , $0.1 \ A$ $533$ R22Wire-wound, $90 \ \Omega \pm 0.1\%$ part of 1611-0280F1FUSE, $230 \ V$ , $0.1 \ A$ $533$ R23Wire-wound, $90 \ \Omega \pm 0.1\%$ $6010-1200$ S2SWITCH, Rotary Wafer, $791$ R27POTENTIOMETER, $30 \ M \pm 20\%$ $6010-2600$ S3	R6	Composition, $0.56 M\Omega \pm 10\%$	6100-4565		Electronytic mock, $20 \mu\text{F}$	4400-0900
R9Wire-wound, 330 Ω ±10%6760-1339 6760-9339C8Mica, 50 pF ±10%466R10Wire-wound, 3.3 Ω ±10%6760-9339 6760-9339C10Mica, 0.005 µF ±10%480R11Composition, 10 MΩ ±10%6100-3275 6100-3275C11Mica, 0.009 µF ±0.25%050R13Composition, 27 kΩ ±10%6100-3275 6100-3105C12Mica, 0.009 µF ±0.25%050R14Composition, 10 kΩ ±10%6100-3275 6100-3105C13Plastic, 0.5 µF ±1.5%486R15POTENTIOMETER, 20 kΩ ±5%0973-4080 6710-1600C15Mica, 0.01 µF ±10%476R17POTENTIOMETER, 20 kΩ ±5%0973-4080 6710-1600C16Mica, 0.008 µF ±10%480R17POTENTIOMETER, 11.2 - 11.7 kΩ0433-4020 433-4020MiscelLaneous480R18POTENTIOMETER, 11.2 - 11.7 kΩ0433-4020 433-4020F1FUSE, 230 V, 0.1 A533R21Precision, 100 kΩ ±0.1% 700 Ω ±0.1%part of 1611-0280 710 - 750 ΩF1FUSE, 230 V, 0.1 A533R24Wire-wound, 90 Ω ±0.1% 710 - 750 Ωpart of 1611-0280 710 - 750 ΩJ1JACK, EXT FILTER426R25Wire-wound, 9.98 Ω ±0.25% 72part of 1611-0280 710 - 750 ΩSWITCH, Rotary Wafer, MULTIPLY CAP BY789R27POTENTIOMETER, 30 Ω ±10% 73,000 - 77,000 ΩG010-2600 73,000 - 77,000 ΩS3SWITCH, Rotary Wafer, 73,000 - 77,00 ΩR31A R31BPOTENTIOMETER, 10 Ω ±10% 73,000 - 77,000 Ω6100-3105 73,000 - 77,000 Ω6100-3105 73,000 - 77	R7	Composition, 22 k $\Omega$ ±10%	6120-3229		20 µF	
R10Wire-wound, $3.3 \Omega \pm 10\%$ 6760-9339C10Mica, $0.005 \mu F \pm 10\%$ 480R11Composition, $10 M\Omega \pm 10\%$ 6100-6105C11Mica, $0.005 \mu F \pm 10\%$ 480R12Composition, $27 k\Omega \pm 10\%$ 6100-3275C12Mica, $0.009 \mu F \pm 0.25\%$ 050R13Composition, $10 k\Omega \pm 10\%$ 6100-3105C14Plastic, $0.5 \mu F \pm 1.5\%$ 486R14Composition, $10 k\Omega \pm 10\%$ 6100-3105C14Plastic, $0.5 \mu F \pm 1.5\%$ 486R15POTENTIOMETER, $20 k\Omega \pm 5\%$ 0973-4080C15Mica, $0.0088 \mu F \pm 10\%$ 476R16POTENTIOMETER, $20 k\Omega \pm 5\%$ 0973-4080C16Mica, $0.0088 \mu F \pm 10\%$ 480R17POTENTIOMETER, $11.2 - 11.7 k\Omega$ 0433-4020MisceLLANEOUS480R18POTENTIOMETER, $11.2 - 11.7 k\Omega$ 0433-4020MisceLLANEOUS533R21Precision, $10 k\Omega \pm 0.025\%$ 0602-3052F1FUSE, $230 V$ , $0.1 A$ 533R22Wire-wound, $90 \Omega \pm 0.1\%$ part of 1611-0280F2FUSE, $230 V$ , $0.1 A$ 533R24Wire-wound, $0.988 \Omega \pm 0.25\%$ part of 1611-0280J1JACK, EXT FILTER426R25Wire-wound, $0.988 \Omega \pm 10\%$ 6010-1200S1SWITCH, Toggle, POWER791R27POTENTIOMETER, $3 M\Omega \pm 10\%$ 6010-2600MULTIPLY CAP BY789R31APOTENTIOMETER, $13,980 - 14,760 \Omega$ 0371-4150S3SWITCH, Rotary Wafer,791R31BPOTENTIOMETER, $10\%$ 6100-3105S6SWITCH, Rotary Wafer, FILTER789 <td>R8</td> <td>Composition, <math>3.3 \text{ k}\Omega \pm 10\%</math></td> <td>6120-2339</td> <td>C7</td> <td>Mica, 0.001 μF ±10%</td> <td>4660-6400</td>	R8	Composition, $3.3 \text{ k}\Omega \pm 10\%$	6120-2339	C7	Mica, 0.001 μF ±10%	4660-6400
R11Composition, 10 MΩ ±10%6100-6105 6100-3275C11Mica, 0.005 $\mu$ F ±10%480R12Composition, 27 kΩ ±10%6100-3275 (C12C12Mica, 0.0099 $\mu$ F ±0.25%050R13Composition, 10 kΩ ±10%6100-3105C14Plastic, 0.5 $\mu$ F ±1.5%486R15POTENTIOMETER, 20 kΩ ±5%0973-4080C15Mica, 0.01 $\mu$ F ±1.5%486R16POTENTIOMETER, 20 kΩ ±5%0973-4080C16Mica, 0.0088 $\mu$ F ±10%476R18POTENTIOMETER, 11.2 - 11.7 kΩ0433-4020MiscelLaneous480R19POTENTIOMETER, 11.2 - 11.7 kΩ0433-4020F1FUSE, 230 V, 0.1 A533R21Precision, 10 kΩ ±0.025%0602-3052F1FUSE, 230 V, 0.1 A533R22Wire-wound, 90 Ω ±0.1%part of 1611-0280F2FUSE, 230 V, 0.1 A533R24Wire-wound, 9.0 $\alpha$ ±0.1%part of 1611-0280F2FUSE, 230 V, 0.1 A533R24Wire-wound, 9.0 $\alpha$ ±0.1%0371-4140F1PLUG, Power426R25Wire-wound, 9.0 $\alpha$ ±0.1%6010-1200S2SWITCH, Rotary Wafer,719R27POTENTIOMETER, 3 MΩ ±10%6010-2600S3SWITCH, Rotary Wafer,789R31APOTENTIOMETER, 5 kΩ ±10%6010-2600S3SWITCH, Rotary Wafer,789R31BPOTENTIOMETER, 13,980 - 14,760 Ω73,000 - 77,000 ΩS3SWITCH, Toggle, ORMAL - ADD 30%791R31BComposition, 10 k1 ±10%6100-5335T1TRANSFORMER734 </td <td>R9</td> <td>Wire-wound, 330 <math>\Omega</math> ±10%</td> <td>6760-1339</td> <td></td> <td>Mica, 50 pF ±10%</td> <td>4660-2000</td>	R9	Wire-wound, 330 $\Omega$ ±10%	6760-1339		Mica, 50 pF ±10%	4660-2000
R12Composition, 27 kΩ ±10%6100-3275C12Mica, 0.0099 μF ±0.25%050R13Composition, 27 kΩ ±10%6100-3275C13Plastic, 0.5 μF ±1.5%486R14Composition, 10 kΩ ±10%6100-3105C14Plastic, 0.5 μF ±1.5%486R15POTENTIOMETER, 20 kΩ ±5%0973-4080C15Mica, 0.01 μF ±10%476R16POTENTIOMETER, 20 kΩ ±5%0973-4080C16Mica, 0.0088 μF ±10%480R17POTENTIOMETER, 20 kΩ ±5%0973-4080C16Mica, 0.0088 μF ±10%480R18POTENTIOMETER, 11.2 - 11.7 kΩ0433-4020MiSCELLANEOUS533R20Precision, 10 kΩ ±0.1%part of 1611-0280F1FUSE, 230 V, 0.1 A533R21Precision, 100 kΩ ±0.1%part of 1611-0280F2FUSE, 230 V, 0.1 A533R22Wire-wound, 90 Ω ±0.1%part of 1611-0280F2FUSE, 230 V, 0.1 A533R23Wire-wound, 9.0 ±0.1%part of 1611-0280F2FUSE, 230 V, 0.1 A533R24Wire-wound, 9.0 ±0.1%part of 1611-0280F2FUSE, 230 V, 0.1 A533R25Wire-wound, 9.0 ±0.1%part of 1611-0280F1INDUCTOR034R26APOTENTIOMETER, 136.6 - 144.6 Ω0371-4140S1SWITCH, Rotary Wafer,791R27POTENTIOMETER, 0.5 MΩ ±10%6010-2600S3SWITCH, Rotary Wafer,791R31APOTENTIOMETER, 13,980 - 14,760 Ω0371-4150S4SWITCH, Rotary Wafer,791R31BPOTENT	R10	Wire-wound, $3.3 \Omega \pm 10\%$	6760-9339			4800-1600
R13Composition, 27 kΩ ±10%6100-3275C13Plastic, 0.5 μF ±1.5%486R14Composition, 10 kΩ ±10%6100-3105C14Plastic, 0.5 μF ±1.5%486R15POTENTIOMETER, 20 kΩ ±5%0973-4080C15Mica, 0.01 μF ±10%476R16POTENTIOMETER, 20 kΩ ±5%0973-4080C16Mica, 0.0088 μF ±10%486R17POTENTIOMETER, 20 kΩ ±5%0973-4080C16Mica, 0.0088 μF ±10%486R18POTENTIOMETER, 11.2 - 11.7 kΩ0433-4020C16MiscelLaneous480R20Precision, 10 kΩ ±0.025%0602-3052F1FUSE, 115 V, 0.2 A533R21Precision, 100 kΩ ±0.1%part of 1611-0280F2FUSE, 115 V, 0.2 A533R23Wire-wound, 90 Ω ±0.1%part of 1611-0280F2FUSE, 230 V, 0.1 A533R24Wire-wound, 0.988 Ω ±0.25%part of 1611-0280J1JACK, EXT FILTER426R25Wire-wound, 0.988 Ω ±0.25%0371-4140J1INDUCTOR034R26APOTENTIOMETER, 30 Ω ±10%6010-2600K11INDUCTOR034R27POTENTIOMETER, 30 Ω ±10%6010-2600S3SWITCH, Rotary Wafer,53R31APOTENTIOMETER, 5 kΩ ±10%6010-2600S3SWITCH, Rotary Wafer,54R31BPOTENTIOMETER, 13,980 - 14,760 Ω6100-3105S4SWITCH, Toggle, NORMAL - ADD 30%791R31BPOTENTIOMETER, 13,980 - 14,760 Ω6100-3105S6SWITCH, Rotary Wafer, FILTER789R33 <t< td=""><td>R11</td><td>Composition, <math>10 M\Omega \pm 10\%</math></td><td>6100-6105</td><td></td><td></td><td>4800-0800</td></t<>	R11	Composition, $10 M\Omega \pm 10\%$	6100-6105			4800-0800
R14Composition, 10 kΩ ±10%6100-3105C14Plastic, 0.5 μF ±1.5%486R15POTENTIOMETER, 20 kΩ ±5%0973-4080C15Mica, 0.01 μF ±10%476R16POTENTIOMETER, 20 kΩ ±5%0973-4080C16Mica, 0.0088 μF ±10%480R17POTENTIOMETER, 20 kΩ ±5%0973-4080C16Mica, 0.0088 μF ±10%480R18POTENTIOMETER, 11.2 - 11.7 kΩ0433-4020MiscellaneousMiscellaneous480R20Precision, 100 kΩ ±0.1%602-3052G10-1600Miscellaneous533R21Precision, 100 kΩ ±0.1%part of 1611-0280F1FUSE, 230 V, 0.1 A533R22Wire-wound, 90 Ω ±0.1%part of 1611-0280F2FUSE, 230 V, 0.1 A533R24Wire-wound, 0.988 Ω ±0.25%part of 1611-0280J1JACK, EXT FILTER426R26APOTENTIOMETER, 136.6 - 144.6 Ω0371-4140J1JACK, EXT FILTER426R26APOTENTIOMETER, 10%6010-1200S2SWITCH, Rotary Wafer,791R27POTENTIOMETER, 0.5 MΩ ±10%6010-2600S3SWITCH, Rotary Wafer,789R28POTENTIOMETER, 13,980 - 14,760 ΩG0371-4150S4SWITCH, Rotary Wafer,789R31BPOTENTIOMETER, 13,980 - 14,760 ΩG0371-4150S4SWITCH, Toggle, NORMAL - ADD 30%791R31BPOTENTIOMETER, 13,980 - 14,760 ΩG0371-4150S4SWITCH, Rotary Wafer, FILTER789R32Composition, 10 kΩ ±10%6100-5335T1TRANSFORMER	R12	Composition, 27 k $\Omega$ ±10%	6100-3275	C12	Mica, 0.0099 µF ±0.25%	0505-4630
R15POTENTIOMETER, $20 \ k\Omega \pm 5\%$ 0973-4080C15Mica, $0.01 \ \mu\text{F} \pm 10\%$ 476R16POTENTIOMETER, $20 \ k\Omega \pm 5\%$ 0973-4080C16Mica, $0.0088 \ \mu\text{F} \pm 10\%$ 480R17POTENTIOMETER, $20 \ k\Omega \pm 5\%$ 0973-4080C16Mica, $0.0088 \ \mu\text{F} \pm 10\%$ 480R18POTENTIOMETER, $20 \ k\Omega \pm 5\%$ 0973-4080Miscellaneous480R19POTENTIOMETER, $11.2 - 11.7 \ k\Omega$ 0433-4020Miscellaneous533R20Precision, $10 \ k\Omega \pm 0.025\%$ 0602-3052F1FUSE, $230 \ V$ , $0.1 \ A$ 533R21Precision, $100 \ k\Omega \pm 0.1\%$ part of 1611-0280F2FUSE, $230 \ V$ , $0.1 \ A$ 533R23Wire-wound, $90 \ \Omega \pm 0.1\%$ part of 1611-0280F2FUSE, $230 \ V$ , $0.1 \ A$ 533R24Wire-wound, $0.988 \ \Omega \pm 0.25\%$ part of 1611-0280J1JACK, EXT FILTER426R26APOTENTIOMETER, $136.6 - 144.6 \ \Omega$ 0371-4140J1INDUCTOR034R26BPOTENTIOMETER, $30 \ k\Omega \pm 10\%$ 6010-2600S3SWITCH, Rotary Wafer,7910R27POTENTIOMETER, $3 \ M\Omega \pm 20\%$ 6010-2600S3SWITCH, Rotary Wafer,7890R31APOTENTIOMETER, $13,980 - 14,760 \ \Omega 30371-4150$ S4SWITCH, Rotary Wafer,7910R31BPOTENTIOMETER, $13,980 - 14,760 \ \Omega 30371-4150$ S4SWITCH, Rotary Wafer,7910R31APOTENTIOMETER, $73,000 - 77,000 \ \Omega 3071-4150$ S4SWITCH, Rotary Wafer,7910R33Composition, $10 \ k\Omega \pm 10\%$ 6100-31	R13	Composition, 27 k $\Omega$ ±10%	6100-3275			4860-4602
R16POTENTIOMETER, 20 kΩ ±5%0973-4080C16Mica, 0.0088 μF ±10%480R17POTENTIOMETER, 20 kΩ ±5%0973-4080MiscellaneousMiscellaneous480R18POTENTIOMETER, 20 kΩ ±5%0973-4080MiscellaneousMiscellaneous480R19POTENTIOMETER, 11.2 - 11.7 kΩ0433-4020Miscellaneous533R20Precision, 10 kΩ ±0.025%0602-3052F1FUSE, 115 V, 0.2 A533R21Precision, 100 kΩ ±0.1%part of 1611-0280F2FUSE, 230 V, 0.1 A533R23Wire-wound, 90 Ω ±0.1%part of 1611-0280F2FUSE, 230 V, 0.1 A533R24Wire-wound, 0.988 Ω ±0.25% part of 1611-0280J1JACK, EXT FILTER426R25Wire-wound, 0.988 Ω ±0.25% part of 1611-0280J1INDUCTOR034R26APOTENTIOMETER, 30 kΩ ±10%6010-1200S2SWITCH, Toggle, POWER791R27POTENTIOMETER, 30 kΩ ±10%6010-2600S3SWITCH, Rotary Wafer,13,980 - 14,760 ΩR30POTENTIOMETER, 3 MΩ ±20%6010-2600S3SWITCH, Rotary Wafer,DIR-STANDBY-REV7890R31APOTENTIOMETER, 73,000 - 77,000 Ω0371-4150S4SWITCH, Toggle, NORMAL - ADD 30%7910R312Composition, 10 kΩ ±10%6100-3105S6SWITCH, Rotary Wafer, FILTER7890R33Composition, 3.3 MΩ ±10%6100-5335T1TRANSFORMER034	R14	Composition, $10 \text{ k}\Omega \pm 10\%$	6100-3105			4860-4602
R17POTENTIOMETER, 20 kΩ ±5%0973-4080R18POTENTIOMETER, 20 kΩ ±5%0973-4080R19POTENTIOMETER, 11.2 - 11.7 kΩ0433-4020R20Precision, 10 kΩ ±0.025%0602-3052R21Precision, 100 kΩ ±0.1%6710-1600R22Wire-wound, 900 Ω ±0.1%part of 1611-0280R23Wire-wound, 90 Ω ±0.1%part of 1611-0280R24Wire-wound, 90 Ω ±0.1%part of 1611-0280R25Wire-wound, 0.988 Ω ±0.25%part of 1611-0280R26APOTENTIOMETER, 136.6 - 144.6 Ω0371-4140R26BPOTENTIOMETER, 30 kΩ ±10%6010-1200R27POTENTIOMETER, 30 kΩ ±10%6010-2600R28POTENTIOMETER, 30 kΩ ±10%6010-2600R31APOTENTIOMETER, 13,980 - 14,760 Ω0371-4150R32Composition, 10 kΩ ±10%6100-3105R33Composition, 3.3 MΩ ±10%6100-3105R33Composition, 3.3 MΩ ±10%6100-5335R33Composition, 3.3 MΩ ±10%6100-5335	R15	POTENTIOMETER, 20 k $\Omega$ ±5%	0973-4080	C15	Mica, 0.01 μF ±10%	4760-0100
R18POTENTIOMETER, 20 kΩ ±5%0973-4080R19POTENTIOMETER, 11.2 - 11.7 kΩ0433-4020R20Precision, 10 kΩ ±0.025%0602-3052R21Precision, 100 kΩ ±0.1%6710-1600R22Wire-wound, 900 Ω ±0.1%part of 1611-0280R23Wire-wound, 90 Ω ±0.1%part of 1611-0280R24Wire-wound, 9 Ω ±0.1%part of 1611-0280R25Wire-wound, 0.988 Ω ±0.25%part of 1611-0280R26APOTENTIOMETER, 136.6 - 144.6 Ω0371-4140R26BPOTENTIOMETER, 30 kΩ ±10%6010-1200R27POTENTIOMETER, 0.5 MΩ ±10%6010-2600R31APOTENTIOMETER, 5 kΩ ±10%6010-2600R31APOTENTIOMETER, 73,000 - 77,000 ΩR31AR31BPOTENTIOMETER, 13,980 - 14,760 Ω0371-4150R32Composition, 10 kΩ ±10%6100-3105R33Composition, 3.3 MΩ ±10%6100-5335R33Composition, 3.3 MΩ ±10%6100-5335	R16	POTENTIOMETER, $20 \text{ k}\Omega \pm 5\%$	0973-4080	C16	Mica, 0.0088 μF ±10%	4800-1900
R19POTENTIOMETER, $11.2 - 11.7 k\Omega$ $0433-4020$ R20Precision, $10 k\Omega \pm 0.025\%$ $0602-3052$ R21Precision, $100 k\Omega \pm 0.1\%$ $6710-1600$ R22Wire-wound, $90\Omega \pm 0.1\%$ part of 1611-0280R23Wire-wound, $90\Omega \pm 0.1\%$ part of 1611-0280R24Wire-wound, $9\Omega \pm 0.1\%$ part of 1611-0280R25Wire-wound, $0.988 \Omega \pm 0.25\%$ part of 1611-0280R26APOTENTIOMETER, $136.6 - 144.6 \Omega$ $0371-4140$ R26BPOTENTIOMETER, $136.6 - 144.6 \Omega$ $0371-4140$ R27POTENTIOMETER, $0.5 M\Omega \pm 10\%$ $6010-1200$ R28POTENTIOMETER, $0.5 M\Omega \pm 10\%$ $6010-2600$ R30POTENTIOMETER, $13,980 - 14,760 \Omega$ $0371-4150$ R31APOTENTIOMETER, $13,980 - 14,760 \Omega$ $0371-4150$ R31APOTENTIOMETER, $13,980 - 14,760 \Omega$ $0371-4150$ R32Composition, $10 k\Omega \pm 10\%$ $6100-3105$ R33Composition, $3.3 M\Omega \pm 10\%$ $6100-5335$ R33Composition, $3.3 M\Omega \pm 10\%$ $6100-5335$	R17	POTENTIOMETER, 20 k $\Omega$ ±5%	0973-4080			
R19POTENTIOMETER, $11.2 - 11.7 k\Omega$ $0433-4020$ R20Precision, $10 k\Omega \pm 0.025\%$ $0602-3052$ R21Precision, $100 k\Omega \pm 0.1\%$ eart of 1611-0280R22Wire-wound, $900 \Omega \pm 0.1\%$ part of 1611-0280R23Wire-wound, $90 \Omega \pm 0.1\%$ part of 1611-0280R24Wire-wound, $9 \Omega \pm 0.1\%$ part of 1611-0280R25Wire-wound, $0.988 \Omega \pm 0.25\%$ part of 1611-0280R26APOTENTIOMETER, $136.6 - 144.6 \Omega$ $0371-4140$ R26BPOTENTIOMETER, $136.6 - 144.6 \Omega$ $0371-4140$ R27POTENTIOMETER, $10\%$ $6010-1200$ R28POTENTIOMETER, $3 0 k\Omega \pm 10\%$ $6010-2600$ R30POTENTIOMETER, $13,980 - 14,760 \Omega$ $0371-4150$ R31APOTENTIOMETER, $13,980 - 14,760 \Omega$ $0371-4150$ R32Composition, $10 k\Omega \pm 10\%$ $6100-3105$ R33Composition, $10 k\Omega \pm 10\%$ $6100-5335$ R33Composition, $3.3 M\Omega \pm 10\%$ $6100-5335$	R18	POTENTIOMETER, $20 \text{ k}\Omega \pm 5\%$	0973-4080		MISCELLANEOUS	
R21Precision, $100 \ k\Omega \pm 0.1\%$ $6710-1600$ F1FUSE, $230 \ V$ , $0.1 \ A$ $533$ R22Wire-wound, $900 \ \Omega \pm 0.1\%$ part of $1611-0280$ F2 $115 \ V$ , $0.2 \ A$ $533$ R23Wire-wound, $90 \ \Omega \pm 0.1\%$ part of $1611-0280$ F2 $FUSE$ , $230 \ V$ , $0.1 \ A$ $533$ R24Wire-wound, $9 \ \Omega \pm 0.1\%$ part of $1611-0280$ F2 $FUSE$ , $230 \ V$ , $0.1 \ A$ $533$ R24Wire-wound, $9 \ \Omega \pm 0.1\%$ part of $1611-0280$ F2 $FUSE$ , $230 \ V$ , $0.1 \ A$ $533$ R25Wire-wound, $0.988 \ \Omega \pm 0.25\%$ part of $1611-0280$ J1JACK, EXT FILTER $426$ R26APOTENTIOMETER, $136.6 - 144.6 \ \Omega$ $710 - 750 \ \Omega$ $0371-4140$ S1SWITCH, Toggle, POWER $791$ R27POTENTIOMETER, $30 \ k\Omega \pm 10\%$ $6010-1200$ S2SWITCH, Rotary Wafer, DIR-STANDBY-REV $789$ R28POTENTIOMETER, $3 \ M\Omega \pm 20\%$ $6010-2600$ S3SWITCH, Rotary Wafer, DIR-STANDBY-REV $789$ R31APOTENTIOMETER, $5 \ k\Omega \pm 10\%$ $6010-300$ $S4$ SWITCH, Toggle, GEN $791$ R31BPOTENTIOMETER, $13,980 - 14,760 \ \Omega \\ 73,000 - 77,000 \ \Omega$ $6100-3105$ S6SWITCH, Rotary Wafer, FILTER $789$ R32Composition, $10 \ k\Omega \pm 10\%$ $6100-3105$ S6SWITCH, Rotary Wafer, FILTER $789$ R33Composition, $3.3 \ M\Omega \pm 10\%$ $6100-5335$ T1TRANSFORMER $034$	R19	POTENTIOMETER, $11.2 - 11.7 \text{ k}\Omega$	0433-4020			
R21Precision, 100 km $\pm 0.1\%$ $0.010^{-1000}$ $230^{-1}$ , $0.14^{-1}$ $533^{-100}$ R22Wire-wound, 900 $\pm \pm 0.1\%$ part of $1611-0280$ F2FUSE, $115^{-1}$ V, 0.2 A $533^{-1}$ R23Wire-wound, 90 $\pm \pm 0.1\%$ part of $1611-0280$ F2FUSE, $230^{-1}$ V, 0.1 A $533^{-1}$ R24Wire-wound, $90 \pm \pm 0.1\%$ part of $1611-0280$ J1JACK, EXT FILTER $426^{-1}$ R25Wire-wound, $0.988^{-1} \pm 0.25\%$ part of $1611-0280$ J1JACK, EXT FILTER $426^{-1}$ R26BPOTENTIOMETER, $136.6 - 144.6^{-1} - 0.280^{-1}$ $0371-4140^{-1}$ S1SWITCH, Toggle, POWER $791^{-1}$ R26BPOTENTIOMETER, $30^{-1} \times 10\%$ $6010-1200^{-1}$ S2SWITCH, Rotary Wafer, $791^{-1}$ R27POTENTIOMETER, $30^{-1} \pm 10\%$ $6010-2600^{-2}$ S3SWITCH, Rotary Wafer, $791^{-1}$ R29POTENTIOMETER, $30^{-1} \pm 20\%$ $6010-2600^{-2}$ S3SWITCH, Rotary Wafer, $10^{-1} \times 10\%$ R31APOTENTIOMETER, $13,980 - 14,760^{-1} \Omega$ $0371-4150^{-1} \times 10\%$ S4SWITCH, Toggle, NORMAL - ADD $30\%$ $791^{-1} \times 10\%$ R31APOTENTIOMETER, $13,980 - 14,760^{-1} \Omega$ $6100-3105^{-1} \times 10\%$ S4SWITCH, Rotary Wafer, FILTER $791^{-1} \times 10\%$ R32Composition, $10^{-1} \times 10\%$ $6100-3105^{-1} \times 10\%$ S6SWITCH, Rotary Wafer, FILTER $791^{-1} \times 10\%$ R33Composition, $3.3^{-1} \times 10\%$ $6100-5335^{-1} \times 10\%$ $510^{-1} \times 10\%$ $510^{-1} \times 10\%$ R33Composition, $3.3^{-1} \times 10\%$ $6100-5335^{-1}$	R20	Precision, $10 \text{ k}\Omega \pm 0.025\%$	0602-3052	<b>F1</b>		5330-0600
R23Wire-wound, $90 \Omega \pm 0.1\%$ Part of $1611-0280$ R24Part of $1611-0280$ Part of $1611-0280$ R25Part of $1611-0280$ Part of $1611-0280$ R26Part of $1611-0280$ PARPart of $1610-0280$ PARPart of $1610-0280$ PARPart of $1610-02800$ PAR <td>R21</td> <td></td> <td>6710-1600</td> <td>1.1</td> <td><sup>103</sup><sup>103</sup><sup>103</sup><sup>10</sup>, 230 V, 0.1 A</td> <td>5330-0400</td>	R21		6710-1600	1.1	<sup>103</sup> <sup>103</sup> <sup>103</sup> <sup>10</sup> , 230 V, 0.1 A	5330-0400
R23Wire-wound, 90 $\Omega \pm 0.1\%$ part of 1611-0280230 V, 0.1 A533R24Wire-wound, 9 $\Omega \pm 0.1\%$ part of 1611-0280J1JACK, EXT FILTER426R25Wire-wound, 0.988 $\Omega \pm 0.25\%$ part of 1611-0280J1JACK, EXT FILTER426R26APOTENTIOMETER, 136.6 - 144.6 $\Omega$ 0371-4140N1PL1PLUG, Power424R26BPOTENTIOMETER, 30 k $\Omega \pm 10\%$ 6010-1200S2SWITCH, Toggle, POWER791R27POTENTIOMETER, 0.5 M $\Omega \pm 10\%$ 6010-2600S3SWITCH, Rotary Wafer,791R28POTENTIOMETER, 3 M $\Omega \pm 20\%$ 6010-2600S3SWITCH, Rotary Wafer,789R30POTENTIOMETER, 5 k $\Omega \pm 10\%$ 6010-2600S3SWITCH, Rotary Wafer,789R31APOTENTIOMETER, 13,980 - 14,760 $\Omega$ 0371-4150S4SWITCH, Toggle, ORMAL - ADD 30\%7910R31BPOTENTIOMETER, 13,980 - 14,760 $\Omega$ 6100-3105S4SWITCH, Rotary Wafer,7910R32Composition, 10 k $\Omega \pm 10\%$ 6100-3105S6SWITCH, Rotary Wafer, FILTER789R33Composition, 3.3 M $\Omega \pm 10\%$ 6100-5335T1TRANSFORMER034	R22	Wire-wound, 900 $\Omega$ ±0.1% part of	1611-0280	- F2		5330-0600
R25Wire-wound, $0.988 \Omega \pm 0.25\%$ part of $1611-0280$ L1INDUCTOR $034$ R26APOTENTIOMETER, $136.6 - 144.6 \Omega$ $0371-4140$ PL1PLUG, Power $424$ R26BPOTENTIOMETER, $30 k\Omega \pm 10\%$ $6010-1200$ S2SWITCH, Toggle, POWER $7910$ R27POTENTIOMETER, $30 k\Omega \pm 10\%$ $6010-1200$ S2SWITCH, Rotary Wafer, $7910$ R28POTENTIOMETER, $0.5 M\Omega \pm 10\%$ $6010-2600$ S2SWITCH, Rotary Wafer, $7910$ R29POTENTIOMETER, $3 M\Omega \pm 20\%$ $6010-2600$ S3SWITCH, Rotary Wafer, $7910$ R30POTENTIOMETER, $5 k\Omega \pm 10\%$ $6010-2600$ S3SWITCH, Rotary Wafer, $7910$ R31APOTENTIOMETER, $13,980 - 14,760 \Omega$ $0371-4150$ S4SWITCH, Toggle, NORMAL - ADD 30\% $7910$ R318POTENTIOMETER, $73,000 - 77,000 \Omega$ $6100-3105$ S6SWITCH, Rotary Wafer, FILTER $7890$ R32Composition, $10 k\Omega \pm 10\%$ $6100-5335$ T1TRANSFORMER $0344$	R23	Wire-wound, 90 $\Omega$ ±0.1% part of	1611-0280	1.7	230 V, 0.1 A	5330-0400
R26A R26B POTENTIOMETER, $136.6 - 144.6 \Omega$ $710 - 750 \Omega$ $0371 - 4140$ PL1PLUG, Power $424$ R26 	R24	Wire-wound, 9 $\Omega$ ±0.1% part of	1611-0280	J1	JACK, EXT FILTER	4260-1500
R26BPOTENTIOMETER, POTENTIOMETER,710 - 750 $\Omega$ 03/1-4140S1SWITCH, Toggle, POWER791R27POTENTIOMETER, POTENTIOMETER, $30 k\Omega \pm 10\%$ $6010-1200$ S2SWITCH, Rotary Wafer, MULTIPLY CAP BY789R28POTENTIOMETER, POTENTIOMETER, $3M\Omega \pm 20\%$ $6010-2600$ S3SWITCH, Rotary Wafer, DIR-STANDBY-REV789R30POTENTIOMETER, POTENTIOMETER, $5k\Omega \pm 10\%$ $6010-2600$ S3SWITCH, Rotary Wafer, DIR-STANDBY-REV789R31A R31BPOTENTIOMETER, 73,000 - 77,000 $\Omega$ $30371-4150$ S4SWITCH, Toggle, NORMAL - ADD 30%7910R32Composition, $10 k\Omega \pm 10\%$ $6100-3105$ S6SWITCH, Rotary Wafer, FILTER7890R33Composition, $3.3 M\Omega \pm 10\%$ $6100-5335$ T1TRANSFORMER $0344$	R25	Wire-wound, 0.988 $\Omega$ ±0.25% part of	1611-0280	L1	INDUCTOR	0345-4450
R27POTENTIOMETER, $30 k\Omega \pm 10\%$ $6010-1200$ S1SWITCH, $10ggle, 10WHC$ $7910$ R28POTENTIOMETER, $0.5 M\Omega \pm 10\%$ $6010-1200$ S2SWITCH, Rotary Wafer,R29POTENTIOMETER, $3 M\Omega \pm 20\%$ $6010-2600$ S3SWITCH, Rotary Wafer,R30POTENTIOMETER, $5 k\Omega \pm 10\%$ $6010-2600$ S3SWITCH, Rotary Wafer,R31APOTENTIOMETER, $13,980 - 14,760 \Omega$ $0371-4150$ S4SWITCH, Toggle, NORMAL - ADD 30\%R31BPOTENTIOMETER, $13,980 - 14,760 \Omega$ $0371-4150$ S4SWITCH, Toggle, GEN7910R32Composition, $10 k\Omega \pm 10\%$ $6100-3105$ S6SWITCH, Rotary Wafer, FILTER7890R33Composition, $3.3 M\Omega \pm 10\%$ $6100-5335$ T1TRANSFORMER $0344$	R26A		0371-4140		PLUG, Power	4240-0700
R28POTENTIOMETER, 0.5 M $\Omega$ ±10%6010-2600MULTIPLY CAP BY789R29POTENTIOMETER, 3 M $\Omega$ ±20%6010-2600S3SWITCH, Rotary Wafer,R30POTENTIOMETER, 5 k $\Omega$ ±10%6010-0800DIR-STANDBY-REV789R31APOTENTIOMETER, 73,000 - 77,000 $\Omega$ S4SWITCH, Toggle, NORMAL - ADD 30%791R32Composition, 10 k $\Omega$ ±10%6100-3105S6SWITCH, Rotary Wafer, FILTER789R33Composition, 3.3 M $\Omega$ ±10%6100-5335T1TRANSFORMER034	R26B	710 - 750 Ω	0371-4140	S1	SWITCH, Toggle, POWER	7910-1300
R29       POTENTIOMETER, $3 M\Omega \pm 20\%$ 6010-2600       S3       SWITCH, Rotary Wafer,         R30       POTENTIOMETER, $5 k\Omega \pm 10\%$ 6010-0800       DIR-STANDBY-REV       7890         R31A       POTENTIOMETER, $73,000 - 77,000 \Omega$ $0371-4150$ S4       SWITCH, Toggle, NORMAL - ADD 30%       7910         R32       Composition, $10 k\Omega \pm 10\%$ 6100-3105       S6       SWITCH, Rotary Wafer, FILTER       7890         R33       Composition, $3.3 M\Omega \pm 10\%$ 6100-5335       T1       TRANSFORMER       0341	R27	POTENTIOMETER, $30 \text{ k}\Omega \pm 10\%$	6010-1200	S2	SWITCH, Rotary Wafer,	
R30       POTENTIOMETER, $5 k\Omega \pm 10\%$ 6010-0800       DIR-STANDBY-REV       7890         R31A       POTENTIOMETER, $13,980 - 14,760 \Omega$ $\Omega$ (371-4150       S4       SWITCH, Toggle, NORMAL - ADD 30%       7910         R31B       POTENTIOMETER, $73,000 - 77,000 \Omega$ $G$ (3071-4150       S5       SWITCH, Toggle, GEN       7910         R32       Composition, $10 k\Omega \pm 10\%$ $6100-3105$ S6       SWITCH, Rotary Wafer, FILTER       7890         R33       Composition, $3.3 M\Omega \pm 10\%$ $6100-5335$ T1       TRANSFORMER       0345	R28	POTENTIOMETER, $0.5 M\Omega \pm 10\%$	6010-2600		MULTIPLY CAP BY	7890-0240
R31A R31B R31BPOTENTIOMETER, 73,000 - 77,000 $\Omega$ 13,980 - 14,760 $\Omega$ 0371-4150S4SWITCH, Toggle, NORMAL - ADD 30% 7910 S5R32Composition, 10 k $\Omega$ ±10%6100-3105 6100-5335S4SWITCH, Toggle, GEN7910 7910R33Composition, 3.3 M $\Omega$ ±10%6100-5335S6SWITCH, Rotary Wafer, FILTER7890 7910	R29	POTENTIOMETER, $3 M\Omega \pm 20\%$	6010-2600	S3	SWITCH, Rotary Wafer,	
R31BPOTENTIOMETER, 73,000 - 77,000 $\Omega$ 0.071-4130S5SWITCH, Toggle, GEN7910R32Composition, 10 k $\Omega$ ±10%6100-3105S6SWITCH, Rotary Wafer, FILTER7890R33Composition, 3.3 M $\Omega$ ±10%6100-5335T1TRANSFORMER0349	R30	POTENTIOMETER, $5 k\Omega \pm 10\%$	6010-0800		DIR-STANDBY-REV	7890-0252
R31B73,000 - 77,000 $\Omega$ S5SWITCH, Toggle, GEN7910R32Composition, 10 kΩ ±10%6100-3105S6SWITCH, Rotary Wafer, FILTER7890R33Composition, 3.3 MΩ ±10%6100-5335T1TRANSFORMER0349	R31A	DOTENTIONETED 13,980 - 14,760 S			SWITCH, Toggle, NORMAL - ADD 30%	7910-1300
R32         Composition, 10 kΩ ±10%         6100-3105         S6         SWITCH, Rotary Wafer, FILTER         789           R33         Composition, 3.3 MΩ ±10%         6100-5335         T1         TRANSFORMER         034	R31B	POTENTIONETER, 73,000 - 77,000 Ω	20371-4130	S5	SWITCH, Toggle, GEN	7910-1500
R33 Composition, $3.3 \text{ M}\Omega \pm 10\%$ 6100-5335 T1 TRANSFORMER 0344				S6		7890-1680
			6100-5335	T1		0345-4430
K34 Composition, 1 MM ±10% 0100-5105   12 TRANSFORMER 0343	R34	Composition, $1 M\Omega \pm 10\%$	6100-5105	T2	TRANSFORMER	0345-4430
R35 Precision, 796 Ω ±0.5% 6700-0400 V1 TUBE, Type 6E5 8366	R35	Precision, 796 $\Omega$ ±0.5%	6700-0400	V1	TUBE, Type 6E5	8360-4600
		Precision, 81,100 $\Omega$ ±0.5%	6700-3500		TUBE, Type 6SJ7	8360-6600
			6110-4279	<b>V</b> 3		8360-8300



\*<u>1 S-2 five wafer switch arm mechanically connected as viewed fra</u> #2 Twin Twisted Conductor Shielded lead from S-2 to S-3



WEST CONCORD, MASSACHUSETTS,\* 01781

617 369-4400

617 646-7400

#### SALES ENGINEERING OFFICES

#### METROPOLITAN NEW YORK\*

Broad Avenue at Linden Ridgefield, New Jersey, 07657 Telephone N.Y. 212 964-2722 N.J. 201 943-3140

#### SYRACUSE

Pickard Building East Molloy Road Syracuse, New York, 13211 Telephone 315 454-9323

#### PHILADELPHIA

General Radio Company Fort Washington Industrial Park Fort Washington, Pennsylvania 19034 Telephone 215 646-8030

#### WASHINGTON\* and BALTIMORE

Rockville Pike at Wall Lane Rockville, Maryland 20852 Telephone 301 946-1600

#### ORLANDO

113 East Colonial Drive Orlando, Florida, 32801 Telephone 305 425-4671

\* Repair services are available at these offices.

#### CHICAGO\*

6605 West North Avenue Oak Park, Illinois, 60302 Telephone 312 848-9400

#### CLEVELAND

5579 Pearl Road Cleveland, Ohio, 44129 Telephone 216 886-0150

#### LOS ANGELES\*

1000 North Seward Street Los Angeles, California, 90038 Telephone 213 469-6201

#### SAN FRANCISCO

1186 Los Altos Avenue Los Altos, California, 94022 Telephone 415 948-8233

#### DALLAS

2501-A West Mockingbird Lane Dallas, Texas, 75235 Telephone 214 Fleetwood 7-4031

#### **TORONTO\***

99 Floral Parkway Toronto 15, Ontario, Canada Telephone 416 247-2171

#### MONTREAL

Office 395, 1255 Laird Boulevard Town of Mount Royal, Quebec, Canada Telephone 514 737-3673

General Radio Company (Overseas), 8008 Zurich, Switzerland General Radio Company (U.K.) Limited, Bourne End, Buckinghamshire, England Representatives in Principal Overseas Countries