

"The emulsion used on this material is the only one with which I have ever worked which reproduces everything there is in the negative."

*—Opinion of a leading
commercial photographer
in Chicago.*

Referring to

Raygraphs

A high speed emulsion on an opaque fibre base.

Fast enough for direct exposure.

Ten times faster than printing paper.

Extreme latitude for over-exposure.

Use Raygraphs on negatives usually considered failures and note results—contact print or enlargement.

Every laboratory doing any amount of photography will find RAYGRAPHS invaluable for various phases of their work.

Send for samples and prove them out.

VICTOR X-RAY CORPORATION




Supply Sales Department

2012 Jackson Blvd.,

Chicago



LIGHT SOURCES



For
Optical Research
and
Demonstration

L568. Colored Flame Burner. A Bunsen burner, a vessel for solution and an atomizer are mounted together. The atomizer when attached to an air supply blows a fine spray of the solution in the gas, giving a bright, broad flame with the color characteristic of the salt used. The flame if viewed edgewise gives considerable intensity. All metal parts exposed to the solution are heavy porcelain enameled. \$8.50

L575. Arc Lamp with carbons in line. Distance adjustable with rack and pinion motion which dis-places both carbons symmetrically \$14.50

L522. Support for Spark Terminals. The clips holding the metal electrodes are readily adjustable and are well insulated \$15.00

THE GAERTNER SCIENTIFIC CORPORATION

1201 Wrightwood Ave.

Successor to Wm. Gaertner & Co.

Chicago, U. S. A.

NEW MICROSCOPE SLIDE CATALOG

Several hundred new preparations are described in the new and enlarged 1926 edition of the Turtor Microscope Slide Catalog

Write for your copy today



The Sign of the Turtor Pledges Absolute Satisfaction

General Biological Supply House
761-8 E. 69th Place Chicago, Illinois

Alkali Metal Hydride PHOTOELECTRIC CELLS

of extreme sensitivity

for

Laboratory and stellar photometry
Measurement of reflecting power
Sorting of materials according to color
Measurement of illumination
Operation of relays by changes of light

Write for circular

THE G-M SCIENTIFIC INSTRUMENT CO.
202 Broadway North Urbana, Illinois



VACUUM PUMPS

28" vacuum guaranteed. 1½ cu. ft. free air per minute. Piston type for continuous duty. Complete with motor gauge, etc.

Write for circular.


ALSO ENGINEERING CO.
47 W. 68rd St., New York

PHONELESCOPE

TEACHES THROUGH THE
The Movies for all
SOUND

WAVE MOTION

ELECTRICITY
HERBERT GROVE DORSEY



"MEYLAN" STOP WATCHES

For all Scientific Purposes

Approved and used by prominent colleges and universities as well as leading industrial laboratories. Accurate, durable, dependable and guaranteed. Ask for free illustrated catalogue of various models of "MEYLAN" stop watches.

A. R. & J. E. MEYLAN
552 Seventh Ave.,
New York, N. Y.



WE REPAIR ALL MAKES OF STOP WATCHES