

**ELECTRO-MEDICAL
LABORATORY, INC.**
HOLLISTON, MASSACHUSETTS, U.S.A.

PRECISION INSTRUMENTS


FOR

ELECTRICAL STIMULATION

**VARIABLE FREQUENCY STIMULATOR
CHRONAXIMETRE
SIXTY-CYCLE STIMULATOR**

Write for Catalogue B

**MAKERS OF THE
ELECTROENCEPHALOGRAPH**



No. 10 Analytical Balance


HENRY TROEMNER

The Standard of Excellence for 99 years.

Manufacturer of the highest grade Analytical, Assay and Pulp Balances and Weights of Precision.

Catalog 1929S

911 Arch Street
Philadelphia, Pa.



Instructor's Table No. F-1111

*Your
Wishes
Come
True!*

Science Instructors often marvel at the thoroughness of understanding of Kewaunee Engineers and Designers in providing the many wished for conveniences. They say, "How could Kewaunee understand so perfectly the instructor's problems?" "How could they foresee so perfectly the needs of my school, my students and myself?" It's very simple when you know the secret. For 35 years Kewaunee has lived with your problems, invited your counsel, and acted on your suggestions. No wonder "Your Wishes Come True" when you install Kewaunee Furniture.

If you do not have a copy of the Kewaunee Catalog of Laboratory and Vocational Furniture, write today.



Ever-Hold
Automatic
Adjustable
Stool
No. E-1824

Kewaunee Mfg. Co.
LABORATORY FURNITURE EXPERTS

C. G. CAMPBELL, Pres. and Gen. Mgr.
115 Lincoln St., Kewaunee, Wis.
Eastern Branch: 220 E. 42nd St.,
New York, N. Y.
Mid-West Office: 1208 Madison St.,
Evanston, Ill.
Representatives in Principal Cities

JENA



Fritted Glass Filters for Bacteria Filtration



The structure of the fritted glass in these filters is denser than diatomaceous earth filters. Users report complete sterilization of broth cultures, and the successful use of these filters for the preparation of sterile solutions of carbohydrates to be added to culture media.

Jena Fritted Glass Filters have been found particularly advantageous when it is necessary to work with small amounts, and where it is desirable to use a filter free from antigenic substances. They are easily cleaned after use and withstand repeated sterilization.

Made with a disc of #5 porosity having an average pore diameter of between 1 and 1.5 microns, over a disc of #3 porosity. Now also available in micro filter type.

In following capacities: micro filter—2 cc.; pressure filter—35 cc.; crucible—30 cc.; funnels in 30, 60, 140 and 450 cc.



Available at leading laboratory supply dealers

Write for Catalogue J.L. 270

FISH-SCHURMAN CORPORATION

260 East 43rd Street, New York City
U. S. Agents, Jena Glass Works, Schott & Gen.

Fish-Schurman

