


BAUSCH & LOMB PHOTOMICROGRAPHIC EQUIPMENT



CERTAINTY PLUS EASE IN OPERATION

In many laboratories, large and small, where a photomicrographic record is a required part of every microscopic observation, Bausch & Lomb Photomicrographic Equipment is used because of its speed, accuracy and convenience in operation. Here are some of the major reasons why B&L Equipment GBVP, one of the most universally used, simplifies the making of high quality photomicrographs:

1. This one equipment provides all magnifications from the lowest Micro Tessars to the highest oil immersion objectives and high power eyepieces.
2. Adjustable microscope support accommodates practically any standard type of microscope.
3. Camera has 40" bellows extension and all-metal back which takes up to 8" x 10" plates.
4. Extension device permits focusing microscope from operating position back of camera.
5. Illuminating unit and microscope support are mounted on the same base.
6. Light source, microscope and camera are rigidly aligned.
7. Concealed, adjustable spring suspension shock absorber counteracts vibration.
8. Heavy rigid supporting stand has convenient plate holder and accessory cabinet.

In the complete B&L line of Photomicrographic apparatus, you will find an equipment best suited to your individual requirements. For complete details, write to Bausch & Lomb Optical Co., 642 St. Paul Street, Rochester, N. Y.

For The Metallographer

The B&L Large Metallographic Equipment is used in leading industrial, commercial and educational laboratories. It is one of a complete line of instruments designed especially for metallographic work. Catalog E-225, "Metallographic Equipment and Accessories" will be sent on request.

BAUSCH & LOMB OPTICAL COMPANY

EST.  1853

FOR YOUR EYES, INSIST ON BAUSCH & LOMB EYEWEAR, MADE FROM BAUSCH & LOMB GLASS TO BAUSCH & LOMB HIGH STANDARDS OF PRECISION