


Savage Hall, new headquarters of Cornell University's School of Nutrition, being dedicated today at a convocation entitled "Nutrition—A Key to Human Welfare." The new \$500,000 structure, named in honor of the late Elmer Seth Savage, a pioneer in nutrition at Cornell and a member of its faculty for 35 years, was financed by farmers of the Northeast; the equipment was provided by the State of New York.





Speedomax for Thermostat Test Temperatures—Vernay Laboratories, Inc.

Micromax for Concrete Temps. Dewey & Almy Chem. Co.

**Micromax for Hydrogenation** Temperatures—Univ. N. Carolina

## **Recorders and Controllers for** "EVERYBODY"

Did you ever try this way of selecting the correct Recorders or Controllers for laboratory use?

First step is an analysis (with our help, if you wish) of such factors as:

Type of instrument—recording or non-recording; singlepoint or multiple; round-chart or strip; speed of response. In automatic control, whether to use a recorder or a nonrecorder.

Type of Control action if any-on-and-off; fully proportioning; manual regulation based on instrument readings.

Characteristics of the problem such as size and construction of vessels, inherent lags, etc. which influence instrumentation.

## Meeting Science's Need for Stable Control

For any combination of requirements which may appear, Micromax and/or Speedomax instruments meet the need -exactly and completely. And they do more-they contribute a highly stable continuous operation.


Stability means "no interference" by weather, temperature, vibration, dust, electric fields, etc. Stability enables Micromax and/or Speedomax to meet the needs of almost "everybody". For details, write Leeds & Northrup Co., 4926 Stenton Ave., Phila. 44, Pa.



HEAT-TREATING FURNACES



**Speedomax in Plane Logs** Surveys-Aero Service Corp.



Speedomax Shows Spectrum Densities—Amer. Hammered Piston Ring Div.