

VOL. 107 • NO. 2783 • PAGES 433-454


April 30, 1948

Science



Presentation of the Viking Fund Medal

(See page 441)



Twice every minute the **TECHNICON** Pipette Washer completes its syphon-cycle . . . seven surging gallons of water to loosen and dislodge clinging residues. Wastes are then literally swept out of the pipettes by the violent suction set up during the brisk ten-second emptying phase.

Only the Technicon washer offers such rapidity and thoroughness, delivering in six cycles a full basket-load of immaculate pipettes (capacity 200 of 1 ml in 100th's, up to 18" long). The stainless steel

basket, proof against corrosive reagents, carries its load of pipettes undisturbed through every stage: cleansing, washing and drying. Requires no individual handling, no shifting from container to container, practically eliminates breakage and chippage. A *must* for the busy laboratory: let us send you Bulletin #6150 giving particulars.

THE TECHNICON COMPANY
215 E. 149 St., New York 51, N. Y.

No fuss or muss . . .
pipette washing can
be as easy as this

Technicon
automatic **PIPETTE WASHER**

