


July 9, 1948

Science



Advisory Committee of the
First International Poliomyelitis Conference

(See page 34)



**you know
your tissues
will spread
e-v-e-n-l-y**

The water temperature will be just right for smooth, even spreading: never too hot or too cold. No melting paraffin to distort tissues; no folds or wrinkles. No bubbles due to frequent water change . . . temperature stays constant day in, day out. Bulletin 6320 gives particulars: send for it.

TECHNICON COMPANY, 215 E. 149th St., New York 51, N. Y.

.... it's a



Technicon
constant temperature
water bath

