


KLETT ELECTROPHORESIS



CUSTOM MADE

TOOL FOR THE ANALYSIS
OF COMPLEX COLLOID SYSTEMS, AND FOR
THE CONTROL OF PRODUCTION OF
PURIFIED PROTEINS, ENZYMES, HORMONES

KLETT MANUFACTURING CO.
179 EAST 87TH STREET
NEW YORK, N. Y.



***Immediately Available
for Research Investigators***

GELATIN GLUTATHIONE

- For intravenous use as a blood extender in shock burns, radiation injury, severe depletions of fluid and toxic infections.
- As a new field for investigation in rheumatic fever, arthritis and in ACTH induced diabetes.

Gelatin in hydrolysis is well known to yield a non-antigenic product. Glutathione is the tripeptide of the three amino acids glycine, cysteine and glutamine used in biologic detoxifications. Glutathione increases the osmotic pressure of gelatin 20% and inhibits syneresis (ageing of colloids) thus advancing gelatin usefulness.

Supplied in 500 cc vials, sterile, ready for use intravenously containing 6% gelatin (knox) and 0.5% glutathione S H (reduced form) or S S (oxidized form) as desired.

Bibliography on Request

PHYSIOLOGICAL CHEMICALS CO., INC.
3 LE FEVRE LANE, NEW ROCHELLE, N. Y.



Proven

DAYLIGHT
FLUORESCENCE
TESTER

Hundreds of users in laboratory, field and factory now make fluorescent checks quickly and easily—in broad daylight—with the Menlo Fluoretor (blacklight).

Exclusive dark chamber holds comparison test tubes, solid samples, powders or fluids. Neoprene viewing cone excludes surrounding light for examination of large areas.

Choice of filtered 2537 or 3650-63 Angstrom frequencies; interchangeable head for alternate frequency. Transmission curve with every instrument. Standardized mercury-vapor discharge tubes. Instant-starting, troublefree. Portable—weighs less than 2 pounds, ready to use. Models for operation from either two standard flashlight batteries or 110-120 v AC. Write for brochure, data on U-V use in your field.

MENLO RESEARCH LABORATORY
Drawer A-9, Box 522, Menlo Park, Calif.
DEALER INQUIRIES INVITED