


## THE MICRO-FLOW REGULATOR

Here is a new clamp for regulation of flow through plastic or rubber tubing, 0"— $\frac{3}{8}$ ".



Adjusts by smooth-acting ratchet which permits minute changes of flow for infusions etc., with ready return to a former rate. One "click" equals .00078" movement of the pressure plate.

Constructed entirely of stainless steel, with the traditionally fine workmanship of the Harvard Apparatus Company, Inc. The cut shows actual size.

From stock, \$5.00

*Catalogue and price-list sent on request.*

THE HARVARD APPARATUS COMPANY, INC.  
Dover Massachusetts



## DETERMINATIVE BACTERIOLOGY LABORATORY MANUAL

by Thomas H. Lord, Kansas State College

This new manual was written for a course in bacterial identification. It is compiled to meet the needs for a step-by-step procedure for identifying the common members of the Eubacteriineae encountered as unknown organisms.

The manual leads the student through the various progressive steps of isolation, purification, morphologic study, cultural study, and biochemical study of an organism. Designed for one semester of three two-hour classes per week.

PRICE \$2.25

Send for an On-Approval copy



**BURGESS PUBLISHING COMPANY**  
428 SOUTH SIXTH STREET • MINNEAPOLIS 15, MINNESOTA



## Light-Scattering INDICATING and RECORDING MICROPHOTOMETER

FOR THE MEASUREMENT OF:

Droplet size of aerosols . . . Purity of proteins and other colloids . . . Rate of change of molecule size in solution . . . Particle size and concentration in dispersed systems . . . Concentration and size of undissolved particles . . . Opacity of almost-opaque materials . . . Molecular aggregation in antigen-antibody reactions . . . Microfluorescence in dilute solutions . . . Depolarization, turbidity, and dissymmetry of high-molecular-weight solutions . . . Light reflectance of dark surfaces . . . Minute precipitation, flocculation, and crystallization . . . Molecular weight of high polymers in solution . . . Scintillation efficiency of nuclear-radiation-sensitive crystals and phosphors . . . Microchemical turbidity analyses.

Excellent sensitivity.  
Completely a-c. operated.  
Has outlet for recorder operation.  
Internal electronic stabilization.

### FEATURES

Compact . . . accessible.  
Fully accessible interior.  
Comparisons of light intensities indicated as decimal percentages.

Built-in stable d-c. amplifier.  
Maladjustment of the controls cannot damage the unit.  
Removable slit and optical systems.



Measures 20  
Micro-microlumens  
to 20 Lumens  
of  
Scattered Light



WRITE FOR BULLETIN 2182-K

**AMERICAN INSTRUMENT COMPANY, INC.**

SILVER SPRING, MARYLAND — IN METROPOLITAN WASHINGTON, D. C.