

## PROTHROMBIN TIME

## NOW DETERMINED SIMPLY, RAPIDLY, ACCURATELY



An important, new scientifically designed apparatus that determines PROTHROMBIN TIME accurately and simply. It is a completely integrated unit that eliminates variables in temperature and results. The apparatus is thermostatically controlled at 37°-38°C. and includes an electronically operated stop watch actuated by a foot switch.

For full details write to:

## J. BEEBER CO., Inc.

838 Broadway, N.Y.C. AL. 4-3510 1109 Walnut Street, Philadelphia

