

keysort-punched to give you quick answers

No available reference material or system can answer your questions faster on the properties of the nuclides than the new Raychronix NUCLIDENTIFIER Set.

You can identify unknown materials, or select isotopes for a particular experiment by sorting in nine major classifications. This valuable research tool comes in a handy metal tray, complete with two probes.

Do you know about these Radio-Steroids?



The Raychronix Carbon 14 compound program has resulted in a stock of labelled Cholesterol, Cholestenone, Progesterone and Testosterone in specific activities as high as 10 uc/mg. Why not send today for full technical data on these?

RAYOHRO	Please send chemical NIX and Nuclidentifier data.
	e Products, Inc.
443 W. Co	ngress. Detroit 26. Mich.
	ngress, Detroit 26, Mich.
Name	
Name Title or Dept.	
Name	

<section-header><text><text><text><text>

black-and-white and color photomicrography Write for Bulletin #810 to Price \$72.---

PHOTOVOLT CORP. 95 Madison Ave. New York 16, N. Y.

BACTO CONCENTRATION DISKS

Bacto-Concentration Disks are for determination of penicillin level in milk, serum, urine and other body fluids.

The concentration disks contain measured amounts of penicillin.

The following standardized materials are recommended for this procedure:

Bacto-Penicillin Disks Bacto-Penase Disks Bacto-Subtilis Spore Suspension Bacto-Whey Agar

DIFCO LABORATORIES Detroit 1, Michigan

Complete details of this procedure available upon request.