


## PHOTOVOLT Line-Operated Multiplier FLUORESCENCE METER Mod. 540



- High-sensitivity for measurement of low concentrations (full-scale setting for 0.001 microgram quinine sulphate)
- Micro-fluorimetry with liquid volumes as low as 1 ml
- Low blank readings, strict linearity of instrument response
- Universally applicable due to great variety of available filters, sample holders, adapters and other accessories
- Interference filters for high specificity of results and for determining spectral distribution of the fluorescent light
- High-sensitivity nephelometry for low degrees of turbidities
- Fluorescence evaluation of powders, pastes, slurries, and solids, also for spot-tests on filter paper without elution

Write for Bulletin #392 to

### PHOTOVOLT CORP.

95 Madison Ave.

New York 16, N. Y.

## DIFCO

## Dehydrated Culture Media DIFCO

The utmost in efficiency and economy in the bacteriological laboratory is realized through use of Dehydrated Culture Media, Difco.

- Convenience—any medium can be instantly prepared.
- Stability—media can be kept without deterioration.
- Availability—each medium is instantly available.
- Completeness—no additional ingredients are required.
- Uniformity—successive lots are identical.
- Comparability—dependable, comparative studies are possible in widely separated laboratories over long periods of time when standardized Difco Products are employed.

Specify "DIFCO"

DIFCO LABORATORIES  
DETROIT 1, MICHIGAN



When MYOCARDIAL INFARCTION OR LIVER DAMAGE IS SUSPECTED  
**DETERMINE SERUM TRANSAMINASE**  
using the inexpensive Reagents from Sigma —

- MALIC DEHYDROGENASE
- REDUCED FORM DPN
- $\alpha$  KETO GLUTARIC ACID
- ASPARTIC ACID

Requires only a few minutes per the exciting procedure outlined by Karmen, Wroblewski, & LaDue, J. Clin. Invest. 34, 126, 1955.

Early reports indicate a remarkable correlation between an infarction and greatly elevated *Transaminase* which reaches

a maximum in approximately 24 hours. The simple test can be performed at anytime up to 24 or more hours following the infarction. In liver damage *Transaminase* levels may be several times higher than the maximum following an infarction and may remain constant for an extended period of time.

Certainly the **PROCEDURE OF CHOICE**  
For **SERUM PHOSPHATASE**  
SIGMA 104 PHOSPHATASE SUBSTRATE

For **UREA NITROGEN**,  
we urge you to use  
SIGMA TYPE II UREASE, Powder or Tablets  
Highest activity, surprising solubility, and low blank.

BE SURE TO GET YOUR SPRING 1956 CATALOG  
LISTING PREVIOUSLY UNAVAILABLE COMPOUNDS

PHONE COLLECT: DAN BROIDA, ST. LOUIS, MISSOURI  
Office: PRospect 1-5750 Home: WYdown 3-6418



## SIGMA CHEMICAL COMPANY

3500 DeKalb STREET, ST. LOUIS 18, MO., U.S.A.

MANUFACTURERS OF THE FINEST BIOCHEMICALS AVAILABLE