


SCIENCE


21 February 1958
Volume 127, Number 3295

Editorial	Education and Strategy	381
Articles	Measurement and Man: <i>S. S. Stevens</i>	383
News of Science	News Articles and Briefs; Scientists in the News; Recent Deaths	389
Association Affairs	Paul E. Klopsteg, President-Elect: <i>G. R. Harrison</i>	395
	AAAS Council Meeting, 1957: <i>D. Wolfe</i>	397
	AAAS Officers, Committees, and Representatives for 1958	400
	Indianapolis Meeting in Retrospect: <i>R. L. Taylor</i>	403
	Public Information Service: <i>S. S. Negus</i>	409
	Reports of Sections and Societies	410
Meetings and Societies	Forthcoming Events	431
	Equipment News	436



New personality for hard-working labs

...J-M Colorlith®
handsome, tough, colorful



Colorlith is available in Surf Green, Charcoal Gray and Cameo Brown.

Here's the new look in labs that's just as practical as it is attractive. It's Johns-Manville's chemical-resistant work-surface material—Colorlith—and it's ready to help you bring new color and design into the school and industrial laboratory.


Made of asbestos and cement—Colorlith offers all these important advantages: **COLOR**—makes any lab a showplace. **WORKABILITY**—offers new design freedom. (Large 4' x 8' sheets can be cut in any shape). **HIGH UNIFORM STRENGTH**—resists years of strenuous service—can be used in thicknesses

as low as 1/4" for resurfacing and fume hoods. **HEAT RESISTANCE**—asbestos and cement offer unmatched protection against thermal shock and fire. And Colorlith offers unusually good chemical resistance, *plus* a smooth-writing surface that lab workers consider so important.

For free Colorlith specification sheet plus a coast-to-coast list of J-M Colorlith laboratory furniture manufacturers, see your local J-M representative. Or write Johns-Manville, Box 14, New York 16, N. Y. In Canada: Port Credit, Ontario.



JOHNS-MANVILLE



The place of the Particle Accelerator in Basic Research...

Particle Accelerators in Biological Research—III

High-energy charged particles are of great value in the study of biological matter. The effects of ionizing radiation can provide information that is invaluable in the development of models to describe the structure and organization of living matter.

Energy-transfer Mechanism




The basic mechanisms by which incident energy is transferred to the target material are of importance. These mechanisms include the "direct effect", which can be described by the target theory of Lea and others, and the "indirect effect" where the energy is absorbed by the solvent and subsequently transferred by radiation-produced intermediates to the solute. Both of these effects take place within the living cell, as shown by irradiation in the wet and dry states. A knowledge of the kinetics of the indirect effect can be important to consideration of the biological effects of ionizing radiation.

Effects of Light and Heavy Particles

The physical structure of microorganisms such as bacteria and viruses has been studied extensively with ionizing radiation. Light particles, such as electrons or x-rays, are easily scattered in matter and may strike the organism from any direction. Analysis of the inactivation from this type of radiation can give a radio-sensitive volume.

Heavy particles, on the other hand, are not easily scattered and travel in straight lines until they are stopped. Analysis of the in-

Characteristics of particles obtainable from Van de Graaff accelerators

	electron	x-ray	proton	deuteron	neutron
Mev	μa	r/min—100cm	μa	μa	n/sec.
1	250	20	50	50	2×10^{10}
2	250	75	50	50	10^{11}
3	1000	1200	200	200	1.2×10^{12}
Typical particle trajectory in target material					

activation data will give a radio-sensitive cross-section for the organism. From data of this type it has been possible to formulate models that agree with structures developed by other experimental methods.

Effects on Man

In addition to fundamental studies on biological materials, there is a great need for better understanding of the effects of ionizing radiation on man. Extensive studies under way on the effects of radiation on animals are providing data that, extrapolated to man, will help greatly in our current problem of providing protection for persons working in radiation environments such as radiation clinics, reactor facilities, and research laboratories. Knowledge of the tolerance dose and the relative biological effects of all the types of ionizing radiation is essential for our safety and for the fulfillment of the ultimate promise of the atomic age.

Necessary Equipment

The radiation equipment necessary to gain this knowledge depends on the major field of interest. The equipment should be variable in energy, should produce all of the particles which are of interest, and should have sufficiently high outputs to produce data in a reasonable period of time. The Van de Graaff® is such a tool. The energies and types of particles it can provide are shown in the table above, with diagrams of their typical modes of energy transfer in matter. The particle energy can be varied from the maximum shown to approximately one tenth of this value. The current is continuously variable from a remote control station. Installation is relatively simple, and the ease of operation and flexibility make it an ideal source of ionizing radiation for biological research.

Write for informative reprints on the effects of charged particles on bacteria.

HIGH VOLTAGE ENGINEERING CORPORATION

BURLINGTON



MASSACHUSETTS

Just Published: Volume 3 (approx. 502 pages)

ANNUAL REVIEW OF ENTOMOLOGY

The Nervous System
Chemoreception in Arthropods
Internal Symbiosis in Insects
Nutritional Requirements of Phytophagous Insects
Recent Advances in Silkworm Nutrition
Uses of Sounds by Insects
Dynamics of Insect Populations
Ovarian Structure and Vitellogenesis in Insects
Genetics and Breeding of the Honey Bee
The Phylogeny of the Panorpid Orders
Zoogeography of Insects
Hybridization and Speciation in Mosquitoes
The Feeding Habits of Biting Flies and Their Significance in Classification
Resistance of Plants to Insects
Biological Control of Insect Pests
Biology of Scarabaeidae
Insect Eradication Programs
The Chemistry and Action of Acaricides
Organic Phosphorus Insecticides for Control of Field Crop Insects
Insecticides for Control of Adult Diptera
Forage Insects and Their Control
Control of Forest Insects
Transmission of Plant Viruses by Arthropods

\$7.00 postpaid (U.S.A.); \$7.50 postpaid (elsewhere)

ANNUAL REVIEWS, INC.
231 Grant Avenue, Palo Alto, California

Choose this...

RELIABLE

source of supply for animal serum,
bloods, complement, cells,
hemolysin, plasma

COLORADO SERUM CO.



- ★ Over 34 years experience producing veterinary biologicals
- ★ Constant research helps fulfill your exacting needs
- ★ Devotion to quality ensures product purity
- ★ Finest animals from high altitude pastures

Write for **FREE** catalogs today
No salesman will call.




**COLORADO
SERUM CO.**



PEAK OF QUALITY

Laboratory and General Office
4950 YORK ST. • DENVER 16, COLORADO
Phone MAin 3-5373



**FARRAND
OPTICAL CO. INC.**

Bronx Blvd. and E. 238th St., New York 70, N. Y.

Engineering
Research • Development
Design • Manufacture



Precision Optics
Electronic and
Scientific Instruments

THE FARRAND® Photoelectric FLUOROMETER

FOR FLUOROMETRIC METHODS of ANALYSIS

Designed for stability and linear response over a wide range of sensitivities, the Farrand Fluorometer provides precise, reliable measurements, especially for extremely low concentrations in micro or macro volumes.


Descriptive Bulletin No. 803R, including list of users, sent upon request.

Bulletin No. 820, describing the FARRAND Spectrofluorometer, also available.

the HITACHI HU-10 ELECTRON MICROSCOPE

— a superb instrument

for metallurgical investigations



The Hitachi HU-10 Universal Electron Microscope is extremely useful for metallurgical investigations through X-Ray micrography, electron diffraction, heating and cooling electron diffraction, high resolution diffraction and selected area diffraction—all done directly on the specimen material. With this instrument, X-Rays at 50, 75 and 100 KV intensities are routinely obtained. The intensity of these rays allows even greater penetrating power than is obtainable from microscopes designed for X-Ray work exclusively. Hitachi Electron Microscopes are unique in this adaptability.

The Hitachi HU-10 Electron Microscope provides continuous magnification from 400x to 100,000x, and its high resolving power of 10 Angstrom units allows photographic enlargements up to 1,000,000x!

For further information, contact

ERB & GRAY SCIENTIFIC

Sole Hitachi distributor for the United States

854 S. FIGUEROA ST., LOS ANGELES 17, CALIF.

HARSHAW *Synthetic* OPTICAL CRYSTALS

● The rapid development and extensive application of instrumental analysis has been made possible to a large extent through the infra-red and ultra-violet transmitting optics supplied by The Harshaw Chemical Co. The development of the process for growing large size crystals and the production of these synthetic crystals commercially are among the Harshaw Laboratory's important contributions to science.


Rigid control of all steps in the manufacturing process assures a uniform product of lasting quality.

Sodium Chloride	Silver Chloride	Cesium Bromide
Potassium Bromide	Calcium Fluoride	Cesium Iodide
Potassium Chloride	Lithium Fluoride	Thallium Bromide
Potassium Iodide	Barium Fluoride	Iodide


Your Requests for Technical Assistance and
Price Quotations are Welcomed

THE HARSHAW CHEMICAL COMPANY

1945 East 97th Street, Cleveland 6, Ohio




• Write for our free 32-page booklet "Harshaw Synthetic Optical Crystals"



NEW LABLINE SEPARATORS

CUT TESTING TIME
FROM HOURS TO MINUTES!



LABLINE Chromatofuge

(Accelerated Chromatography)


HELPS YOU TO:

- Make your tests in 5 to 30 minutes.
- Run 25 different separations at same time.
- Identify samples immediately.
- Make two dimensional separations.

In the new Labline Chromatofuge, the centrifugal developing force is controllable within a wide range, is simple to select, and will not change with voltage or load fluctuations. Rapidity of separation minimizes decomposition and temperature effect, sharpens definition, reduces drying. R_f values are readily determined.

Built-in, one hour timer turns off at end of preset time . . . additional equipment available for continuous, automatically controlled, large volume separations. Strong, stain-free construction, built for long, dependable service.

- Catalog No. 22976—
\$1,095.00 complete.



LABLINE Electro- Chromatofuge

(Continuous Planar
Electrophoresis)

HELPS YOU TO:


- Run one or two samples at same time.
- Make chromatographic and electrophoretic separations—independently, in sequence, or simultaneously.
- Preset centrifugal force (including zero) and electrical force desired for test.
- Perform discontinuous or continuous operations, as desired.

The Labline Electro-Chromatofuge offers greater centrifugal and electrical forces at higher collection speeds. Rectangular head permits separate collection of electrolysis products. Has chromatofuge features plus dependable electrical field and contact components.

- Catalog No. 22977—
\$1,495.00 complete.
- Catalog No. 22978 (Power
Supply, 100-1000 Volts DC
for use with Cat. No. 22977)
\$295.00

• Write for detailed information

KEY TO ALL
LABORATORY SUPPLIES




NEW YORK
Laboratory Supply
COMPANY, INC.

76 Varick Street, New York 13, N. Y.
Telephone: CAnal 6-6504

POLARAD
PROVEN RELIABILITY

ELECTRON RESONANCE SPECTROMETER



Custom-designed for diverse and specific
spectroscopic analysis


Utilizing the selective absorption of electromagnetic radiation by a substance when under the influence of a powerful, steady magnetic field, the unit supplies the following kinds of data:

- Free radical determinations
- Substance identification through its 'g' value
- Molecular structure
- Fine and Hyperfine structure
- Dipole broadening
- Exchange narrowing
- Qualitative prediction concerning the nature and extent of radiation damage to biological materials
- Metals, semiconductors and transitional elements.

Scientific Instruments Division

POLARAD ELECTRONICS CORPORATION

43-20 34th Street, Long Island City 1, N. Y. • EXeter 2-4500



Cambridge pH Meters

The three portable models illustrated above incorporate new and important advantages. A.C. line operation eliminates battery nuisance. Laboratory and Research models use electron-ray tube for precise null-point indication. Industrial model is direct-reading, continuous-indicating; ruggedly built for plant use.

Accuracy: Research .02 pH, Laboratory .05 pH, Industrial .10 pH. Cambridge pH equipment also includes single- and multi-point indicators and recorders. Send for Bulletin 910-M.

CAMBRIDGE INSTRUMENT COMPANY, INC.
3556 GRAND CENTRAL TERMINAL, NEW YORK 17, N.Y.

Pioneer Manufacturers of
PRECISION INSTRUMENTS

FOR THOROUGH CLEANING
of Laboratory Glassware and Equipment

NON-TOXIC FOR TISSUE CULTURE




for Hand Washing



for Mechanical and Automatic Washers

DISPOSO-TRAY

for
VIRUS TITRATIONS




Use as Test Tubes and Culture Tubes. So Low in Price You Can Use Them Once and Throw Them Away!
Sold by Dealers the World Over

LINBRO CHEMICAL CO., INC.

681 DIXWELL AVE., NEW HAVEN 11, CONN.

AUTOMATIC BLOOD CELL & PARTICLE COUNTER

- 1 Lamp
- 2 Filter
- 3 Reservoir
- 4 Counting Chamber
- 5 Photo Cell Stabiliser
- 6 Photo Cell
- 7 Vacuum Gauge
- 8 Level Gauge
- 9 Sample Platform
- 10 Oscilloscope
- 11 Flow Gauge
- 12 Operating Lever
- 13 Pump Selector Switch
- 14 Displacement Pumps
- 15 Counter
- 16 Discrimination Control



ACCURATE: Accuracy better than 3%.

FAST: Complete count in less than 1 minute.

SIMPLE: No preparations of blood samples necessary.

No special training necessary.

COMPLETELY SELF-CONTAINED

Only a touch of a lever is necessary to get complete count.

Each cell is counted individually, therefore extremely accurate.

Developed for the British Atomic Energy Commission.

TECHNICAL INSTRUMENT CO.


98 Golden Gate Ave.

San Francisco, Calif.

GLASS ABSORPTION CELLS

made by

KLETT



Makers of Complete Electrophoresis Apparatus

SCIENTIFIC APPARATUS

Klett-Summerson Photoelectric Colorimeters—Colorimeters—Nephelometers—Fluorimeters—Bio-Colorimeters—Comparators—Glass Standards—Klett Reagents.

Klett Manufacturing Co.

179 East 87 Street, New York, New York


Self Illuminating 7 POWER MAGNIFIER

Thousands of satisfied users among:

- HOBBYISTS
- STAMP COLLECTORS
- MODEL BUILDERS
- MECHANICS
- JEWELERS
- STUDENTS
- LAB. TECHNICIANS
- ARTISTS
- TEACHERS
- INSTRUCTORS
- MED. PROFESSION



Life Size Enlargement



A regular \$9.95 value now available at less than 1/2 former price.

This quality built POWER MAGNIFIER uses not one but TWO PRECISION-GROUND OPTICALLY-PERFECT COATED LENSES. Magnifies with startling sharpness and clarity. Built-in self-illuminator uses two standard flashlight batteries. Opening in housing permits using with other instruments. Built-in gauge measures in inches and centimeters.

Does not include Batteries

Ruggedly constructed of chromed steel and bakelite for years of satisfactory service.

Larger 10 Power Deluxe Model \$3.98
Sorry no C.O.D.

NOVEL MFG. CO.

33 2nd Ave., Dept. M-2309
New York 3, N.Y.


MULTI-BLOCK HEATER

For digestions, extractions,
and other reactions requiring
close temperature control.


Interchangeable heat-transfer
blocks mountable on a single
heating base: hot plate,
tube heater, liquid bath...
sold individually or in
combination.

Constant temperature control
 $\pm 1^\circ\text{C}$.


Variable from ambient to
300°C.



Hot Plate



Tube Heater



Liquid Bath



SEND FOR FREE CATALOG 2125-A

RESEARCH SPECIALTIES CO.
2005 Hopkins St. Berkeley 7, Calif.

Dajac LABORATORIES

OVER 150 NEW
REAGENTS FOR
MEDICAL RESEARCH
AND BIOCHEMICAL
TESTING

Compound

Acetyl Choline Chloride
p-Bromomandelic Acid
p-Bromoaniline
Benzidine (purified)
2-Hydroxy-3-Naphthaldehyde
AET Hydrobromide (2-amino-ethyl isothiuronium)
Beta-Naphthyl Acetate, and
6-Bromo-2-Naphthyl acetate
Sodium-6-Benzoyl-2-naphthyl
phosphate
2-(4-Hydroxybenzeneazo)
benzoic acid
Carbonaphthoxy choline iodide
6-Bromo-2-naphthyl-beta-D-
glucuronide
Indoxyl acetate
N-Benzoyl-dl-phenylalanine beta
naphthol ester
L-Leucyl-beta-naphthylamide
hydrochloride
Reinecke Salt
Thioacetamide
Neutral Buffered formaldehyde
solution
5-Nitrobenzotriazole
Myristoyl choline chloride
2,4,7-Trinitrofluorenone
2,3,5-Triphenyltetrazolium
chloride

Use

Acetyl choline esterase
Zirconium
Pentose detection
Peroxidase
Protein bound amino
groups
Antiradiation
Acetyl choline esterase
Phosphatase
Albumin in blood plasma
Serum choline esterase
Beta-glucuronidase
Esterase and acetyl
choline esterase
Chymotrypsin
Leucine amino
peptidase
Choline
Precipitant for metal ions
Fixation
Nitroreductase
Choline esterase
Polynuclear complexes
Reducing enzymes,
carbohydrates


Custom Syntheses Invited

The Borden Company
5000 Langdon Street
P.O. Box 9522
Philadelphia 24, Pa.

THE *Borden* COMPANY
CHEMICAL DIVISION

BIRD

Kymograph Recording Camera



This instrument makes possible the recording of many
physiological reactions by the photographic process. It is
modestly priced and is suitable for student demonstration
purposes and student use. The kymograph movement is
from our standard 70-140 continuous feed kymograph.

The speeds in cms per minute are approximately:

Used in conjunction with
a light source (including
focusing lens and vertical
exit slit) and a mirror gal-
vanometer, or other type
of optical lever, the phys-
iological action under ob-
servation is recorded on the paper as the vertical line of
light travels back and forth along the adjustable hori-
zontal slit in the front panel of the camera.

	With change gears in nor- mal position	In reverse
1st	2.2	4.3
2nd	11.0	21.5
3rd	55.0	107.5
4th	275.	537.5

The supply roll quill accommodates a roll of Linagraph No.
483 paper 100 feet long by 6 inches wide.

The camera may be detached from the driving unit, by
turning a single knurled screw, and removed to a dark room
for loading or unloading. The driving unit and base need not
be disturbed during this procedure.

**Kymograph Recording Camera with one roll of
Linagraph No. 483 paper**

**Linagraph Paper No. 483 in rolls 6 inches wide, 100
feet long**

PHIPPS & BIRD, INC.

Manufacturers & Distributors of Scientific Equipment



6th & Byrd Streets - Richmond, Va.

**AT LAST! a stabilized
human whole blood
hemoglobin standard and control**

Hemotrol*
Trade Mark

Hemotrol is a natural whole blood internal standard which goes through each step of the procedure in the identical manner used by the laboratory in its normal routine. In this way, Hemotrol is employed to evaluate the entire method and equipment, i.e., technical precision at each step, pipette volume, dilution volume, reagents, hemoglobinometer calibration, and reading. When used as a control, Hemotrol provides a simple and direct check of each day's tests.

Hemotrol has been precisely standardized by independent laboratories employing Total Iron, Oxyhemoglobin, and Cyanmethemoglobin methods. The data are statistically evaluated to insure an accurate value within $\pm 1.0\%$ of the true value, based on 99% probability. Hemotrol carries an 18 month expiration date when stored under refrigeration; no other precautions are necessary. Available at you dealer in vials containing 2.4 ml.; sufficient for 40 control tests and 4 complete standardizations.

Vial (2.4 ml.) \$5.00

Write for brochure SC describing our products.

CLINTON LABORATORIES


6010 Wilshire Boulevard, Los Angeles 36, California

* PATENTS APPLIED FOR

COPYRIGHT 1957, CLINTON LABORATORIES

VARICORD
VARIABLE-RESPONSE RECORDER
for densitometry in paper electrophoresis

Recording function electrically
adjustable from linear to
logarithmic and trans-
logarithmic



Sensitivity
adjustable from
10 to 200 millivolts

Designed specifically to be used with scanning densitometers for correctly-compensated quantitative evaluation of electrophoretic patterns on filter paper

Write for Bulletin #1100

Also : Densitometers pH Meters Colorimeters Fluorescence
Meters Electronic Photometers Multiplier Photometers

PHOTOVOLT CORP.

95 Madison Ave.

New York 16, N. Y.

for chromatography . . .

NEW RINCO

2-Way

FRACTION COLLECTOR

**Automatically Collects by Either
"Timed-Flow" or "Volumetric" Method**


This precision apparatus is a valuable aid to accurate fraction measurement and collection in chromatographic analysis. Operated by a heavy-duty capacitor motor, the RINCO Fraction Collector provides completely automatic operation for any desired period. NO MERCURY, CHEMICALS OR ELECTRIC CURRENT COME IN CONTACT WITH YOUR SAMPLE. When utmost accuracy is required, fractions are measured with manometer assembly, as in diagram.

**3 Receiver Tables Available —
Interchangeable with One Indexer**

24" Diameter—holds 237 tubes (18 mm. x 150 mm.)
18" Diameter—holds 60 tubes (25 and 29 mm.)
14" Diameter—holds 120 tubes (1 to 5 ml, 10 ml, 20 ml, and 50 ml.)

Complete Units, with a variety of attachments and supplementary equipment,

from \$390.00 to \$533.50



For Routine Work or
Critical Determinations

SCHAAR AND COMPANY

Complete Laboratory Equipment

7300 W. Montrose Avenue • Chicago 34, Illinois