



GLASS ABSORPTION CELLS made by **KLETT**



SCIENTIFIC APPARATUS
Klett-Summerson Photoelectric Colorimeters—
Colorimeters—Nephelometers—Fluorimeters—
Bio-Colorimeters—Comparators—Glass Stand-
ards—Klett Reagents.

Klett Manufacturing Co., Inc.
179 East 87 Street, New York, New York

CONFERENCE ON NEWER RESPIRATORY DISEASE VIRUSES



Held at the
National Institutes
of Health

Over 100 outstanding virologists examine the diagnosis, treatment, and present status of the new respiratory diseases in this collection of 42 papers presented at a three-day conference in Bethesda, Maryland, in October, 1962. This exposition of scientific advances in the field of acute respiratory diseases, as reported by participants from seven countries, is all in readable English. Furnishing up-to-date material nowhere else available, it was published by the National Tuberculosis Association as a Supplement to the American Review of Respiratory Diseases. September, 1963, 432 pages.

- Paperbound, \$3.00 per copy.

Mail your order to:

American Review of Respiratory Diseases
Dept. S, 1790 Broadway, New York, N. Y. 10019

for **CIRCULATING**
or **STIRRING**
while **HEATING**

to **150°C**

**HAAKE
MODEL ED**



The Haake constant temperature circulator Model ED may be used either to maintain accurate temperature control of any open bath (stirring) or for control by pumping through externally mounted appliances (circulating) such as spectrophotometers, refractometers, polarimeters, titration vessels, etc.

Unique features include magnetic setting thermostats, pumping to 5 gal/min., accuracy to $\pm 0.01^\circ\text{C}$, built-in cooling coil and automatic "water level" safety shut-off device. For control near or below ambient temperature, the ED may be used in conjunction with heat exchangers down to -60°C .

For complete descriptive literature, please contact:

BRINKMANN

(NEW ADDRESS)
BRINKMANN INSTRUMENTS INC.,
CANTAGUE ROAD, WESTBURY, L. I., NEW YORK

PHILADELPHIA • CLEVELAND • HOUSTON • MIAMI • MENLO PARK, CAL. • ST. LOUIS

LABELED NUCLEOSIDES, NUCLEOTIDES AND RELATED COMPOUNDS

CARBON-14 LABELED

		SP. ACT.*	PRICES				
		mc/mM	10 μ C	50 μ C	100 μ C	0.5 mc	1.0 mc
NEC-5	Adenine-8-C ¹⁴	11.3	\$ —	\$ 25	\$ 35	\$175	\$ 350
NEC-347	Adenosine-5'-Monophosphate-C ¹⁴ (u.l.)	231	40	100	200	†	†
NEC	Bromodeoxyuridine-2-C ¹⁴	15	In process				
NEC-348	Cytidine-5'-Monophosphate-C ¹⁴ (u.l.)	225	55	200	400	†	†
NEC-355	Deoxyuridine-2-C ¹⁴	32.4	—	80	140	700	1400
NEC-51	Guanine-8-C ¹⁴ Sulfate	12.7	—	30	40	200	400
NEC-354	Guanosine-5'-Monophosphate-C ¹⁴ (u.l.)	196	55	200	400	†	†
NEC-149	Hypoxanthine-8-C ¹⁴	4.1	—	35	53	265	530
NEC-185	6-Mercaptopurine-8-C ¹⁴ Hydrate	4.4	—	40	66	330	660
NEC-73	Orotic-6-C ¹⁴ Acid Hydrate	6.5	—	25	34	170	340
NEC-176	Orotic-carboxyl-C ¹⁴ Acid Hydrate	8.4	—	20	30	150	300
NEC-106	Thiothymine-2-C ¹⁴	1.1	—	30	49	245	490
NEC-123	Thiouracil-2-C ¹⁴	2.0	—	—	—	170	340
NEC-156	Thymidine-2-C ¹⁴	30.0	—	70	130	650	1300
NEC-107	Thymine-2-C ¹⁴	8.3	—	35	50	250	500
NEC-107H	Thymine-2-C ¹⁴	30	—	35	50	250	500
NEC-177	Thymine-methyl-C ¹⁴	1.8	—	30	42	210	420
NEC-109	Uracil-2-C ¹⁴	6.6	—	25	31	155	310
NEC-109H	Uracil-2-C ¹⁴	46.7	—	25	31	155	310
NEC-167	Uridine-2-C ¹⁴	28	—	80	140	700	1400
NEC-	Uridine-5'-Monophosphate-C ¹⁴ (u.l.)	200	In process				

*Present lot

†Quote on request

TRITIUM LABELED

		SP. ACT.*	PRICES			
		curies/mM	250 μ C	1 mc	5 mc	25 mc
NET-63	Adenine-H ³	1.3	\$ 20	\$ 40	\$120	\$360
NET-64	Adenosine-H ³	0.3	30	70	210	630
NET-29	Cytidine-H ³	1.2	30	70	210	630
NET-107	Cytosine-H ³	1.1	20	40	120	360
NET-123	Deoxyadenosine-H ³	0.5	40	100	300	900
NET-41	Deoxyuridine-H ³	0.1	30	80	240	720
NET-27	Thymidine-methyl-H ³	6.7	15	40	150	600
NET-27A	Thymidine-methyl-H ³	2.0	15	40	150	600
NET-27X	Thymidine-methyl-H ³	10.0	—	40	150	600
NET-42	Thymine-methyl-H ³	6.6	20	40	120	360
NET-49	Uracil-H ³	3.2	20	40	120	360
NET-28	Uridine-H ³	1.2	30	70	210	630

*Present lot

Write for Schedule K to get complete listing of other NEN labeled chemicals.



NEW ENGLAND NUCLEAR CORP.

575 ALBANY STREET, BOSTON 18, MASSACHUSETTS

TEL. 426-7311 | AREA CODE 617

New Products

Spectrodiffractometer is a dual-function instrument that can be switched from x-ray fluorescence analysis to x-ray diffraction without requiring changing or realignment of parts. The instrument has two permanently mounted x-ray tubes, one for each technique. Either can be switched on without affecting the other. For fluorescence analysis, constant-potential radiation of high intensity is supplied by a vertically mounted end-window tube that is available with tungsten, molybdenum, chromium, or platinum target. For diffraction studies, a horizontally mounted constant potential diffraction tube provides a variety of focal spots ranging from 0.4 to 1.5 mm. The specimen chamber used for emission analysis is water cooled below room temperature. Individual holders allow rapid exchange of analyzing crystal. Located at the base of the instrument are a slow motor and encoder for the manufacturer's single-angle programmer. Once set for a desired sequence, the programmer will automatically analyze up to ten elements by fluorescence analysis or up to five pairs of Bragg angles by x-ray diffraction.—J.s. (Picker X-Ray Corp., Dept. S13, Cleveland, Ohio)

Digital voltmeter (model 2660) is an automatic instrument for the measurement of d-c voltage from 1 μ v to 1000 volts. Resolution is extended by a front-panel meter that verifies the null balance condition. Automatic

ranging is provided down to 3 volts full scale with manual switching for ranges to 30 mv full scale. Automatic polarity switching and decade sensitivity control are standard features. Bridge switching is effected by dry reed relays operated under zero-current conditions. Input resistance is greater than 100 meg to 3 volts and greater than 1000 megohms on higher ranges. Accuracy in the 30-mv range is said to be $\pm 3 \mu$ v; in the 300 mv range, 0.01 percent or 2 digits; in higher ranges, 0.01 percent or 1 digit. Input is isolated from ground and may be operated as much as 500 volts above ground. Common mode rejection is specified greater than 100 db at 60 cy/sec.—J.s. (Houston Instrument Corp., Dept. S61, 4950 Terminal Ave., Belaire 101, Tex.)

Spectrum analyzer (model 2112) covers the frequency range from 22 to 45,000 cy/sec and provides 1/3-octave band analysis as well as full-octave band analysis. The 1/3-octave band analysis is said to be accurate within ± 0.5 db throughout the entire frequency spectrum of the instrument. The instrument uses contiguous filtering with all data provided in terms of 33 standardized frequency bands. The analyzer provides measurement of root mean square, absolute average, or instantaneous peak value of the input signal. The r.m.s. reading is said to be accurate for sinusoidal, periodic, or random signals.—J.s. (B & K Instruments, Inc., Dept. S62, 3044 W. 106th St., Cleveland 11, Ohio)

Vector voltmeter (type 248) requires no adjustment for frequency change between 40 cy and 40 kcy/sec when it is used for direct reading of phase angle between an unknown and a reference signal, or for when measuring the in-phase component of an unknown signal with respect to a reference signal, or for reading the total amplitude of an unknown signal. When

used with the manufacturer's type 209 phase shifter, measurement of the quadrature component with respect to a reference signal can also be performed over the same frequency range. The instrument consists of an input attenuator, a bandpass filter, and a feedback amplifier in the unknown channel, and a phase inverter, phase shifting network, and a push-pull amplifier in the reference channel. The outputs of the two channels are applied to a phase-sensitive detector that rejects all even harmonics and attenuates the odd harmonics in inverse proportion to their orders. The panel meter has its center at 90° or 270° and provides full-scale deflections of $\mp 90^\circ$, $\mp 30^\circ$, and $\pm 10^\circ$. Accuracy is said to be $\frac{1}{4}^\circ$ or 1 percent relative and 1° or 2 percent absolute.—J.s. (Ad-Yu Electronics, Inc., Dept. S63, 249 Terhune Ave., Passaic, N.J.)

Combination spectrograph and scanning monochromator is a one-meter instrument of the Czerny-Turner type. The instrument is furnished with a 1200 line/mm grating, 4- by 4-inch square, mounted kinematically. With this grating, the dispersion is 8.3 Å/mm, resolution is better than 0.15 Å with a 10- μ slit, and wavelength range is 2000 to 10,000 Å. Effective aperture ratio is approximately f/9. A light-tight housing behind the focal curve is designed to hold up to 18 photomultiplier tubes to adapt the instrument for direct reading. Provision is made for purging with helium or nitrogen to improve efficiency at wavelengths below 2200 Å. Twelve scanning speeds are provided for manual and electrical operation.—J.s. (McPherson Instrument Corp., Dept. S64, 530 Main St., Acton, Mass.)

Fluid amplifiers and flip-flops fabricated of photosensitive glass are available in sample quantities. More complex devices with as many as six layers contained in dimensions 5- by 4- by 2 inches are made to purchaser designs. The design is transferred photographically to a flat piece of glass. The portions of the glass exposed to light are more soluble in acid than the unexposed portions, so that passages and ports can be produced by etching away the image partially or completely through the glass. A number of plates made in this way are fused under heat to make a solid, monolithic, three-dimensional system. The system, when connected to a source of fluid, can

The material in this section is prepared by the following contributing writers:

Robert L. Bowman (R.L.B.), with the assistance of Denis J. Prager (D.J.P.), Laboratory of Technical Development, National Heart Institute, Bethesda 14, Md. (medical electronics and biomedical laboratory equipment).


Joshua Stern (J.S.), Basic Instrumentation Section, National Bureau of Standards, Washington 25, D.C. (physics, computing, electronics, and nuclear equipment).

The information reported here is obtained from manufacturers and from other sources considered to be reliable. Neither *Science* nor the writers assume responsibility for the accuracy of the information.

Address inquiries to the manufacturer, mentioning *Science* and the department number.

Super Conducting MAGNET DEWARS

FAST DELIVERY



This Sulfrian 30 liter stainless steel dewar is a breakthrough in cryogenics. Liquid helium loss rate is less than 70 cc per hour. Send for engineering drawings.


Sulfrian also makes liquid nitrogen/oxygen containers, liquid hydrogen/helium containers, maser dewars, free radical dewars and a large line of low temperature equipment. Send for 64 page catalog and price list.

SULFRIAN
Cryogenics, Inc.
1290 Central Ave., Hillside, N. J.
Phone: ELizabeth 5-1975

perform control or logic functions. Fabrication accuracy is said to be better than ± 1 mil and operating parameters of input, control, and output pressures are said to remain within 5 percent for all devices within any one family. The glass fluid amplifiers can be converted by heat treatment into ceramic with three times greater strength, temperature capability up to 500°C, and increased radiation resistance, with some forfeit of dimensional accuracy.—J.S. (Corning Glass Works, Dept. S54, Corning, N.Y.)

Frequency multiplier chain (model 101c) is a crystal-controlled source for generating microwave harmonics. Frequency control is provided by a 5-Mcy/sec crystal. Four multiplier stages are used to increase the frequency from 5 to 450 Mcy/sec. Harmonics at 50 and 150 Mcy/sec are also present in the output at lower amplitudes and can be used to provide additional frequencies. Either fixed or variable crystal-controlled oscillators are available. The fixed-frequency oscillator is said to have a long-term stability of 0.3×10^{-6} and the variable-frequency oscillator, adjustable over a 0.5-percent range of frequencies, a long-term stability of 1×10^{-6} . Crystal current is adjustable up to 50 ma for 1N23 diodes. Power available at 450 Mcy/sec is 75 mw. With addition of the manufacturer's Model 301 frequency multiplier, harmonics up to 40 gcy/sec are said to be obtainable.—J.S. (Micro-Now Instrument Co., Inc., Dept. S53, 6124 N. Pulaski Rd., Chicago 46, Ill.)

Frequency synthesizer furnishes approximately 700,000 discrete frequencies in the 2 to 34 Mcy/sec range with stability said to be 1 part in 10^8 per day with spurious frequencies 120 db down. Design of the instrument is based on the indirect method of frequency synthesis and employs a double superheterodyne circuit that electrically locks four free-running oscillators to a high-stability 1 Mcy/sec reference. The total range of frequencies is covered in four bands: 2 to 4 Mcy/sec in steps of 12.5 cy/sec; 4 to 8 Mcy/sec in steps of 25 cy/sec; 8 to 16 Mcy/sec in steps of 50 cy/sec; and 16 to 34 Mcy/sec in steps of 100 cy/sec. Output level is adjustable from 1 to 2.5 volts between 2 and 34 Mcy/sec. Output impedance is 50 ohms nominal.—J.S. (Manson Laboratories, Inc., Dept. S66, 375 Fairfield Ave., Stamford, Conn.)



ESCAPING RADIOACTIVITY

Problem: A convoy has been attacked by a thermonuclear weapon. What maneuvers will permit the ships to best avoid the radioactivity that will follow?

This is an example of the challenging tasks assigned to the Center for Naval Analyses of The Franklin Institute.

Several possible tactics have been evaluated by CNA analysts. Each has its advantages and disadvantages. If the convoy stays together and maintains course, possibility of collision is minimized, but the swiftest escape from contamination may not be attained. While other maneuvers may reduce the possibility of contamination, they may lead to confusion, minimizing the convoy's over-all progress, and increasing the danger of repeat attack. The conclusion is that one of the intermediate tactics is best.

CAREER OPPORTUNITIES in this and other problem areas are now available with CNA for Operations Analysts, Mathematicians, Physicists, and Engineers. For additional information, write:

Director
CENTER FOR NAVAL ANALYSES
Dept. 5
1710 H St., N. W., Washington, D. C.

CNA
CENTER FOR NAVAL ANALYSES
OF THE FRANKLIN INSTITUTE

OEG • OPERATIONS EVALUATION GROUP
INS • INSTITUTE FOR NAVAL STUDIES
NAVWAG • NAVAL WARFARE ANALYSIS GROUP

An equal opportunity employer

SUPERVISORY POSITIONS AT MELPAR, INC.

BIOPHYSICS GROUP LEADER

Melpar's Research Division has an immediate need for a Biophysical Chemist to head a group in physico-chemical optics. The applicant must have strong academic training and research experience in properties of solutions of biological macromolecules.

BRANCH SUPERVISOR

To be responsible for the technical and administrative guidance of a new group of scientists to work in

- IMMUNOCHEMISTRY
- ENZYMOLOGY
- ANALYTICAL BIOCHEMISTRY
- BIOPHYSICAL CHEMISTRY

Requirements include an advanced degree in one of the above disciplines.

For further details, write in strictest confidence to:

JOHN A. HAVERFIELD

Manager—Professional Placement

MELPAR INC.

A Subsidiary of Westinghouse Air Brake Company

3422 Arlington Blvd., Falls Church, Virginia

(a suburb of Washington, D.C.)

an equal opportunity employer

PHARMACOLOGISTS

*for new research
programs with*

I I T RESEARCH INSTITUTE

The increasing number of investigations in the Life Sciences Division has created select opportunities for pharmacologists at the Ph.D. level. Desired qualifications for these appointments are a background and/or experience in the effects of new drugs on physiological systems and on behavior response.

The advantages of a research career at IITRI are many, including the use of advanced facilities and equipment, opportunity to pursue investigations of personal interest, salaries in competition with industry and opportunities for personal and professional growth.

Your inquiries are invited. Please send complete résumé in confidence to Mr. George Zima. U. S. citizenship required.

IITRI

**IIT
Research
Institute**

*Formerly Armour Research Foundation
of Illinois Institute of Technology*

10 West 35th St., Chicago 16, Illinois

An Equal Opportunity Employer

PERSONNEL PLACEMENT

Replies to blind ads should be addressed as follows:

Box (give number)
SCIENCE
1515 Massachusetts Ave., NW
Washington, D.C. 20005

POSITIONS WANTED

(a) **Ph.D. Zoologist** (embryologist); experienced in, and prefers, teaching and research. (b) **M.S. Pharmaceutical Chemist**, strong industrial and administrative experience; prefers Northeast. (c) **Ph.D. Zoologist-Physiologist**, biochemical, medical and animal parasitology, high-altitude physiology background; prefers research. (Please write for information regarding these and other scientists, senior and junior, in all fields.) Science Division, The Medical Bureau Inc., Burneice Larson, Chairman, 900 North Michigan Avenue, Chicago 11, Illinois. X

POSITIONS OPEN

CLINICAL INVESTIGATION DEPARTMENT

Large midwest pharmaceutical company seeks a man with clinical research experience to become a part of our select group assigned to clinical investigation. Age to 45. The right man will have an opportunity to grow with one of the nation's leading growth companies. For full information, in absolute confidence, write. Box 254, SCIENCE.

ELECTRON MICROSCOPIST, Ph.D., M.S. or equivalent, to supervise Electron Microscope Laboratory and Participate in investigating Project, in Hospital connected Research Institute located in New York Area, Opportunity Available to Start Own Research Project. Salary Open. Box 243, SCIENCE.

POSITIONS OPEN

Research Associate in Physiological Chemistry Ph.D. at post doctoral level for advanced research and study in rapidly growing health center. Investigations vary from enzyme mechanism to nutritional biochemistry. Salary open. Position available Feb. 1, 1964. **Dr. Harold L. Rosenthal, Department Physiological Chemistry, Washington University, 4559 Scott Ave., St. Louis, Mo.**

'Bacteriologist' The National Research Council of Canada, Saskatoon, Sask. requires a Bacteriologist, Ph.D. to study the metabolism of micro-organisms. Applicant should be proficient in the techniques of isolation and identification of bacteria and have a knowledge of genetics, or should be well acquainted with biochemical methods and enzymes.

The salary will be commensurate with qualifications. Address enquiries to Dr. H. R. Salans, Director, Prairie Regional Laboratory, National Research Council, Saskatoon, Saskatchewan.

McGill University, Montreal, Canada Institute of Parasitology

Owing to the approaching retirement of the present Director, the University invites applications for this post. The appointment will be at the rank of Professor and in the first instance for a period of five years. The salary, not less than 12,000 dollars Canadian, will be in accordance with the experience of the person appointed, with retirement, hospitalization and life assurance benefits. Letters of enquiry and of application should be addressed to Dean S. B. Frost, The Faculty of Graduate Studies and Research, McGill University, Montreal 2.

POSITIONS OPEN

ANALYTICAL CHEMIST

Excellent opportunity for chemist at the B.S. or M.S. level with experience in analytical instrumentation & clinical chemistry. Duties will relate to the development of new, and the modification of existing diagnostic procedures. Send résumé to: Employment Manager ORTHO PHARMACEUTICAL CORP., Route 202, Raritan, N.J. An equal opportunity employer.

PROTEIN CHEMIST (Ph.D.) To collaborate in physicochemical research on proteins in relation to chronic diseases.

Starting Salary—\$10,000

Send complete résumé to:—A. Saifer, Ph.D. Chief, Dept. of Biochemistry, Isaac Albert Research Institute, 86 E. 49 St., Bklyn 3, N.Y.

(a) **Physical Chemist, Ph.D./M.S.**, strong instrumentation background; design, test medical, industrial instruments; \$12,000; Central. (b) **Ph.D. Pharmacologist**, cardiovascular research; eastern pharmaceutical house; to \$10,800. (c) **Ph.D. Microbiologist**; strong antibacterial antifungal experience; section chief eastern drug firm. (d) **M.S. Pharmaceutical Chemist**, some industrial experience; supervisory duties; midwestern pharmaceutical house. (e) **Ph.D./M.S. Clinical Biochemist**; 200-bed New England hospital; to \$15,000. (f) **Research Biologist, Ph.D.** or equivalent; biological waste treatment processes; outside U.S. (g) **Virologist**; biological R&D; western laboratory; \$10,000. *Faculty Appointments:* (h) **Ph.D. Zoologist**, ecology training, research interests; Central university. (i) **Microbiology, Biology**; junior college; East; to \$10,230. (j) **M.S. Physical Science**, small eastern college. Many other opportunities available for both Junior and Senior scientists. Please write to Science Division, The Medical Bureau, Inc., Burneice Larson, Chairman, 900 North Michigan Avenue, Chicago 11, Illinois. X

basic research

Biological Sciences with WARNER-LAMBERT RESEARCH INSTITUTE

Fundamental studies here cover a wide field—infectious, cardiovascular, mental, allergic and inflammatory diseases, and endocrinology. The biological science division has the following openings for experienced scientists in independent and collaborative research with chemists and members of the medical profession.

HISTOLOGY/HISTOCHEMISTRY

Experienced Ph.D. or equivalent required for basic research on experimental cardiovascular and reproductive physiological and histopathological problems. Independent as well as collaborative research with endocrinologists, biochemists and pharmacologists.

CARDIOVASCULAR PHARMACOLOGY/ PHYSIOLOGY

Experienced Ph.D. or M.D. required to direct and supervise cardiovascular research group involved in drug evaluation and basic studies. Independent research experience in cardiovascular physiology or pharmacology essential.

BIOCHEMISTRY

Recent Ph.D. with experience in metabolic studies and an interest in host defense mechanisms required for basic studies in these areas.

PSYCHOPHARMACOLOGY—Assistant

B.S. with an academic background in biology and psychology plus 1-3 years of experience in pharmacological evaluation of compounds affecting the central nervous system or experience in experimental animal psychology required to perform preliminary and advanced CNS evaluations of new compounds, including use of operant behaviour techniques. In addition to participation in the screening program will be assigned projects for independent investigation.

ENDOCRINOLOGY—Assistant

B.S. in Biology with 1-3 years of experience in endocrinological or endocrine biochemical work required to assist in related research. Opportunity for an experienced man to assist materially in initial planning for group.

Suburban location between New York and Philadelphia. Professional associations and publication of work encouraged. Library, technical information and abstracting service. Complete benefit program including an attractive education assistance plan. Retraction expenses paid. Send resume to Employment Department.

WARNER-LAMBERT RESEARCH INSTITUTE

Tabor Road, Morris Plains, N.J.

An Equal Opportunity Employer

POSITIONS OPEN

(a) **Assistant Director, Biological Sciences**; req's Ph.D. to supv plann'g, execution rsrch proj animal sci, incl animal disease, parasitol, virol, nutrition; share admin duties w/director R&D; MW affil noted pharm co; to \$17,500, incl fee, mov'g expenses. (b) **Biochemist**; M.S. hsp exper to supv dept, lge gen hsp; 3 pathol; \$9000 up; Los Angeles area. (c) **Biochemist**; M.S. for modern co hlth lab serv'g 2 hsp; to \$8000; coll tw; NE. (d) **Microbiologist**; Ph.D. w/appr 10 yrs lab, supv exper to dir human microbiol sect, prominent pharm co; emphasis develop antibacterial agents. Science Division, Woodward Medical Personnel Bureau, Ann Woodward, Director, 185 North Wabash Avenue, Chicago. X

Clinical Biochemist

Experienced, with supervisory ability, Ph.D. preferred; for 200 bed hospital. Excellent salary and benefits. Send résumé to W. A. Towle, Jr., Administrator, Bristol Hospital, Bristol, Conn.

VIROLOGIST

Melpar's Research Division has an immediate need for a virologist to engage in:

The investigation of the biological, physical and chemical changes occurring during a virus cell interaction. Requirements include a Ph.D. and a background in physics or biochemistry. Experience in instrumentation is also desired.

For further details, write in strictest confidence to:

JOHN A. HAVERFIELD
Manager—Professional Placement

MELPAR INC

A Subsidiary of
Westinghouse Air Brake Company

3419 Arlington Blvd.
Falls Church, Virginia
(a suburb of Washington, D. C.)
an equal opportunity employer

FELLOWSHIPS

UNIVERSITY OF BIRMINGHAM DEPARTMENT OF ANATOMY Research Fellow in Endocrinology

Applications are invited for a Research Fellowship in the Department of Anatomy (Professor Sir Solly Zuckerman, C.B., F.R.S.) to study the influence of environmental and other factors on seasonal changes in the ferret.

The salary will depend on qualifications and experience and will include F.S.S.U. and family allowance benefits.

Further particulars may be obtained from Dr. W. P. Dallas Ross, The Department of Anatomy, The Medical School, Birmingham, 15, to whom applications, naming two referees, should be sent as soon as possible. 254, SCIENCE.

FELLOWSHIPS

UNIVERSITY OF BIRMINGHAM DEPARTMENT OF ANATOMY Research Fellowship in the Physiology of Reproduction

Applications are invited for a Research Fellowship in the Department of Anatomy (Professor Sir Solly Zuckerman, C. B., F.R.S.) to study the physiology of reproduction and in particular to investigate the control of fertility in primates and lower mammals.

The salary will be determined according to qualifications and experience. F.S.S.U. and family allowance will be paid. Initially the appointment will be for a period of two or three years, which may be extended.

Further particulars may be obtained from the Assistant Registrar, the Medical School, Birmingham 15, to whom applications, naming two referees, should be sent as soon as possible.

Graduate Fellowships in Anatomical Sciences are available at the University of Tennessee Medical Units for predoctoral training in the research areas of cytology, histology, cytogenetics, histochemistry, cytochemistry, electron microscopy, tissue culture, neuroanatomy, embryology, teratology, fetal physiology, immunomorphology and endocrinology. Stipends range from \$2800. to \$4200. with exemption from tuition and course fees. For further details write to **Chairman, Department of Anatomy, University of Tennessee Medical Units, Memphis 3, Tennessee.**

Postdoctoral Research Fellowship in Virology

A Postdoctoral Research Fellowship is available to a virologist who has recently completed a Ph.D. Current research in department includes biophysical characteristics of human and animal enteroviruses and myxoviruses; antigenic analysis of reoviruses; rubella virus; certain epidemiologic problems in virology; and electron microscopy.

The present Fellowship is offered to those interested in studying the interference phenomenon caused by rubella and other viruses. Candidates will be encouraged to develop their own research along these lines. In addition, some graduate teaching will be expected. Excellent prospects to right person in well equipped and growing graduate department.

Inquiries, giving academic background, publications, references, and salary expected may be addressed to: Dr. A. J. Rhodes, Director, School of Hygiene, and Head, Department of Microbiology, University of Toronto, Toronto 5, Ontario.

The Market Place

BOOKS AND MAGAZINES

Your sets and files of scientific journals

are needed by our library and institutional customers. Please send us lists and description of periodical files you are willing to sell at high market prices. Write Dept., A35, CANNERS, Inc. Boston 20, Massachusetts

SUPPLIES AND EQUIPMENT



tissue grinders


Homogenizers

- Macro • 8 sizes • 4 models
McShane • Potter-Elvehjem
Parallel • Ten Broeck
- Micro • 3 sizes • 3 models

Available from local dealer
Write for Bulletin AGR-H2

MACALASTER
SCIENTIFIC CORPORATION
243 BROADWAY • CAMBRIDGE, MASS.

SPRAGUE-DAWLEY, INC.
Pioneers in the development of the
STANDARD LABORATORY RAT.
P.O. Box 4220
Madison, Wisconsin
CE 3-5318



The Case of the Meandering ML's

In the vicinity of many research laboratories, small organic solutions are being mysteriously lost, daily. They are placed in small flasks, heated, and then . . . they vanish. Chemists patiently wait for their meandering ml's to show up. Sometimes a few straggle into fraction cutters, but most of them escape. The M.O. is always the same. Our job . . . find them.

Acting on a tip from Dr. Raymond Firestone at Merck & Co., Rahway, we went to work on our new K-28480 Bantam-ware Short Path Fractionator. Make the column about an inch long and let the thermometer bulb and stem serve as a packing, he said. Construct the condenser and vacuum connection in two

inches, total length. Use only one joint on the fraction cutter but design it so that pipettes can be inserted to the bottom of any of the four receivers. Then, tilt the whole assembly so that drops fall from the drip tip $\frac{1}{4}$ 14/20 joint into one of the four receivers without unnecessary wetting. But nothing is vertically aligned, we pleaded. I know, he said. In desperation, we tried it.

We're happy to report that most of the meandering ml's now show up as product. The case is closed. Sequel: K-28480 has been added to the line because it fits in with the whole Bantam-ware concept. Very compact, low hold-up . . . functional.

Incidentally, that little rectangular box in the photo is our new Un-A-Watt voltage controller. It occupies $\frac{1}{15}$ th the space and weighs only $\frac{1}{10}$ th as much as a conventional coil transformer. It is custom made for Bantam-ware heating mantles with 0-60 Volt range, sparkless operation, removable support rod and high sensitivity. It fits in with the Bantam-ware concept also.

Bantam-ware Kits and hundreds of individual pieces are described in Catalog BW-2. Write for it today.

KONTES
GLASS COMPANY
Vineland, New Jersey



Thomas announces the new...

Whirlpool Thermoelectric IMMERSION COOLER FOR LABORATORY BATHS



COOLER, IMMERSION, THERMOELECTRIC, Whirlpool. A Peltier-principle thermal pump for lowering the temperature of laboratory baths. Capable of producing sub-zero temperatures (to -25°C from 25°C ambient) in bath volumes up to 4 liters in insulated vessels. Can be used in larger baths to provide cooling required to maintain control near ambient temperature. Operates on a 115-volt, 50 or 60 cycle a.c. line, on which the unit draws 165 watts. A slow-flowing source of cooling water at normal temperature is also required.

Cooling Unit. Made in the form of a thick nickel-plated copper blade, the cooling surface is 6 inches long \times $2\frac{3}{4}$ wide \times $1\frac{1}{8}$ inches thick, providing an effective area of approximately 45 square inches when immersed to recommended depth. Three-pin receptacle for power input and tapped hole for horizontal attachment of $5\frac{1}{2} \times \frac{3}{8}$ -inch diameter mounting rod are provided in polypropylene shank. Two $\frac{1}{4}$ -inch o.d. tubulations emerge from top for connection of tubing for cooling water input and discharge. Thermoelectric modules within the outer shell are bonded to a cast bronze heat sink which is cooled by the passage of tap water. A thermal cut-off switch attached to the heat sink closes at 71°C if cooling water is inadvertently shut off (see *Power Supply*).

Power Supply. Suitable for unattended, continuous operation to provide 3-volt, 25-ampere d.c. operating power. Contains transformer, two silicon diode rectifiers, choke for reducing output ripple, fuse, pilot lamp and thermal cut-off power switch. Heater in thermal cut-off switch is activated by closure of a bimetal switch attached to the heat sink inside the cooling unit. Although bimetal switch in cooling unit reopens when heat sink temperature drops to 55°C , power will not be restored until power supply switch is reset manually. Plug on d.c. power cable connects the power supply to cooling unit. A.c. input is supplied through a 3-wire cord, with 3-prong plug and adapter.

Cooling Rate. Non-linear; depends upon bath volume and temperature of heat sink coolant. In test with 3.5 liters of ethylene glycol, with coolant and initial bath temperature of 24°C , hourly temperatures were recorded as follows: 14, 4, -3, -8, -12, -14, -16, -17, -18, -19, -20, -21, -22, -23, -24 degrees Centigrade.

4049. THERMOELECTRIC COOLER, as above described, complete with d.c. power supply unit. **250.00**

Bulletin 143, giving detailed information, sent upon request



ARTHUR H. THOMAS COMPANY

Scientific Apparatus and Reagents

VINE STREET AT 3RD • P. O. BOX 779 • PHILADELPHIA 5, PA., U. S. A.