


9 February 1979 • Vol. 203 • No. 4380

\$1.50

SCIENCE

AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE





Announcing the Beckman J2-21... a centrifuge you'll love.

The J2-21 is an exciting new high-speed centrifuge. It keeps the components that made the 21,000-rpm J-21C so reliable—and adds proven features from other Beckman centrifuges. The result: by far the finest high-speed centrifuge ever built.

The J2-21 features: a fast-accelerating, high-torque drive; variable power-braking system; superior sound suppression; easy temperature setting; new, automatic vacuum system; rotor imbalance detector; easy-opening counterbalanced door; compact cabinet; built-in speed calibration; 2°



to 40°C operation; electronic door interlock; and it is UL listed and CSA approved.

There are eleven J2-21 rotors to choose from, including the new JA-17 which spins fourteen 50-ml tubes, and the new JA-20.1 which spins thirty-two 15-ml tubes.

The J2-21: it gives you reliability, quiet efficiency unmatched by any other centrifuge, and it's backed by Beckman sales and service—worldwide.



Send for brochure SB-366 to Beckman Instruments, Inc., Spinco Division, 1117 California Ave., Palo Alto, CA 94304.

BECKMAN®

Circle No. 133 on Readers' Service Card

OLYMPUS VANOX

Modular, universal research microscope system.



This is the most versatile research microscope system ever developed. Naturally, by Olympus. Modern design, rugged squareline stand with the broadest choice of application-oriented options that allow you to customize the microscope to your needs.

Choose transmitted and/or reflected light, fluorescence, phase contrast, differential interference contrast, Hoffman modulation contrast. Select a tungsten, halogen or mercury arc light source, to name a few.

Use any of a wide range of color-corrected, hard coated achromats, plan achromats, fluorites, apochromats or plan apochromats. Combine with our automatic camera system, the PM-10, to obtain flawlessly sharp photomicrographs, 35mm or large format.

It's a virtually unlimited system designed with you and your applications in mind. And all this at unprecedented economy. For details, write Olympus Corporation of America, 4 Nevada Dr., New Hyde Park, NY 11040.

OLYMPUS
SEEING BEYOND MAN'S VISION

In Canada: W. Carsen Co., Ltd., Ontario.

LETTERS	EPA Air Pollution Standard: <i>D. B. Menzel</i> ; Federal Regulation: What Role for Universities?: <i>B. Wolfe</i> ; Problems: <i>G. Green</i> ; <i>G. B. Kolata</i>	500
EDITORIAL	Scientific Relations with China	503
ARTICLES	China in Transition Education, Science, and Technology in China: <i>P. H. Abelson</i> The Chinese Scene: <i>W. D. Carey</i> China: Objectives, Contradictions, and Social Currents: <i>E. E. David, Jr.</i> The Major Histocompatibility Complex of the Mouse: <i>J. Klein</i>	505 509 512 516
NEWS AND COMMENT	Albert-Szent-Gyorgyi, Electrons, and Cancer Yearly Report on Carcinogens Could Be a Potent Weapon in the War on Cancer . . . Congress Seeks New Approach to Arms Control <i>Briefing</i> : Laetrile's Day in Court; Meteorites and Nuclear Power; The Public Interest versus Lollipops; Smog's Not So Bad, EPA Decides	522 525 526 528
RESEARCH NEWS	Mantle Geochemistry: Probing the Source of the Earth's Crust The Scrapie Agent: Is It a Viroid?	530 532
AAAS NEWS	AAAS Board Visit to China—A Brief Report: <i>W. D. Carey</i> ; Warren Weaver, 1894–1978; Fourth Annual R & D Policy Colloquium; AAAS Introduces New Life Insurance Program; AAAS Membership List; Climate Project Organizing Workshop; Edison Centennial Symposium Planned; Chautauqua Short Course Program; For the Library	533

BOARD OF DIRECTORS

EDWARD E. DAVID, JR.
Retiring President, Chairman

KENNETH E. BOULDING
President

FREDERICK MOSTELLER
President-Elect

ELOISE E. CLARK
MARTIN M. CUMMINGS

RENÉE C. FOX
ANNA J. HARRISON

CHAIRMEN AND SECRETARIES OF AAAS SECTIONS

MATHEMATICS (A)
Garrett Birkhoff
Ronald Graham

PHYSICS (B)
Arthur L. Schawlow
Rolf M. Sinclair

CHEMISTRY (C)
Fred Basolo
William L. Jolly

ASTRONOMY (D)
Peter S. Conti
Donat G. Wentzel

PSYCHOLOGY (J)
Frances K. Graham
Meredith P. Crawford

SOCIAL AND ECONOMIC SCIENCES (K)
David L. Sills
Gillian Lindt

HISTORY AND PHILOSOPHY OF SCIENCE (L)
Melvin Kranzberg
Diana L. Hall

ENGINEERING (M)
Daniel C. Drucker
Donald E. Marlowe

EDUCATION (Q)
Fletcher G. Watson
James T. Robinson

DENTISTRY (R)
Carl J. Witkop, Jr.
Harold M. Fullmer

PHARMACEUTICAL SCIENCES (S)
Samuel Elkin
Robert A. Wiley

INFORMATION, COMPUTING, AND COMMUNICATION (T)
Mary E. Corning
Madeline M. Henderson

DIVISIONS

ALASKA DIVISION

Daniel B. Hawkins
President

Keith B. Mather
Executive Secretary

PACIFIC DIVISION

Glenn C. Lewis
President

Alan E. Leviton
Secretary-Treasurer

SOUTHWESTERN AND ROCKY MOUNTAIN DIVISION

James W. O'Leary
President

Lora M. Shields
Executive Officer

SCIENCE is published weekly on Friday, except the last week in December, by the American Association for the Advancement of Science, 1515 Massachusetts Avenue, NW, Washington, D.C. 20005. Second-class postage (publication No. 484460) paid at Washington, D.C., and at an additional entry. Now combined with *The Scientific Monthly*. Copyright © 1979 by the American Association for the Advancement of Science. Domestic individual membership and subscription (51 issues): \$31. Domestic institutional subscription (51 issues): \$65. Foreign postage extra: Canada \$12, other (surface mail) \$15, air-surface via Amsterdam \$35. First class, airmail, school-year, and student rates on request. Single copies \$1.50 (\$2 by mail); back issues \$2.50 (\$3 by mail); classroom rates on request. **Change of address:** allow 6 weeks, giving old and new addresses and seven-digit account number. **Postmaster:** Send Form 3579 to *Science*, 1515 Massachusetts Avenue, NW, Washington, D.C. 20005. *Science* is indexed in the *Reader's Guide to Periodical Literature* and in several specialized indexes.

AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE

BOOK REVIEWS	R. A. Fisher, reviewed by A. P. Dempster; Social Change and Scientific Organization, J. L. Heilbron; Thermophilic Microorganisms and Life at High Temperature, R. W. Castenholz; The Biology and Chemistry of the Compositae, P. H. Raven; Perspectives in Urban Entomology, M. J. Tauber; Books Received and Book Order Service	537
REPORTS	Cloning of Herpes Simplex Type 1 DNA Fragments in a Bacteriophage Lambda Vector: L. W. Enquist et al.	541
	Radioimmunoassay of the Insulin Receptor: A New Probe of Receptor Structure and Function: L. C. Harrison et al.	544
	Transfer of Experimental Allergic Encephalomyelitis with Guinea Pig Peritoneal Exudate Cells: B. F. Driscoll, M. W. Kies, E. C. Alvord, Jr.	547
	Rotaviral Immunity in Gnotobiotic Calves: Heterologous Resistance to Human Virus Induced by Bovine Virus: R. G. Wyatt et al.	548
	Eye Malformations in Rats: Induction by Prenatal Exposure to Nickel Carbonyl: F. W. Sunderman, Jr., et al.	550
	Fertilization: A Uterine Glycosaminoglycan Stimulates the Conversion of Sperm Proacrosin to Acrosin: T. J. Wincek, R. F. Parrish, K. L. Polakoski . . .	553
	A Regulatory Locus for Mouse β -Glucuronidase Induction, <i>Gur</i> , Controls Messenger RNA Activity: K. Paigen, C. Labarca, G. Watson	554
	Thalamic Projections to Layer I of Striate Cortex Shown by Retrograde Transport of Horseradish Peroxidase: R. G. Carey, D. Fitzpatrick, I. T. Diamond	556
	Sex Pheromone in the Dog: M. Goodwin, K. M. Gooding, F. Regnier	559
	Postmating Reproductive Isolation in <i>Pseudophryne</i> and the Evolutionary Significance of Hybrid Zones: D. S. Woodruff	561
	<i>Technical Comments</i> : The <i>Salmonella</i> Mutagenicity Assay: Recommendations: F. J. de Serres and M. D. Shelby; Manganese Nodules on the Sea Floor: Are Economic Mining Operations Feasible?: F. J. Lampietti and L. F. Marcus .	563
PRODUCTS AND MATERIALS	Spectrophotometer; Chromatogram Solvent Equilibration; Ultrasonic Processor; Calculator; Pressure Meter for HPLC; Digital Image Analyzer; Chromatography Column; Literature	566

MIKE MC CORMACK RUSSELL W. PETERSON	JOHN C. SAWHILL CHEN NING YANG	WILLIAM T. GOLDEN Treasurer	WILLIAM D. CAREY Executive Officer
GEOLOGY AND GEOGRAPHY (E) Linn Hoover Ramon E. Bisque	BIOLOGICAL SCIENCES (G) Donald S. Farnier Walter Chavin	ANTHROPOLOGY (H) James B. Watson Priscilla Reining	
MEDICAL SCIENCES (N) Theodore Cooper Leah M. Lowenstein	AGRICULTURE (O) Election in progress Coyt T. Wilson	INDUSTRIAL SCIENCE (P) Herbert I. Fustfeld Robert L. Stern	
STATISTICS (U) Richard L. Anderson Ezra Glaser	ATMOSPHERIC AND HYDROSPHERIC SCIENCES (W) Eugene W. Bierly Glenn R. Hilst	GENERAL (X) Ruth B. Pitt S. Fred Singer	

COVER

Coal-burning locomotive on main line between Peking and Inner Mongolia. China has abundant reserves of coal and is discovering more. This fuel, supplemented by hydropower, will serve as the nation's principal energy source during the next decades, thus sparing oil for high value petrochemicals and for export. See page 505. [William Wight, AAAS]

The American Association for the Advancement of Science was founded in 1848 and incorporated in 1874. Its objects are to further the work of scientists, to facilitate cooperation among them, to foster scientific freedom and responsibility, to improve the effectiveness of science in the promotion of human welfare, and to increase public understanding and appreciation of the importance and promise of the methods of science in human progress.

CORNING

FRESH

PEEL AN ORANGE

Our new zip bag makes it easier to get to our dishes. On the outside you'll find a stacking ring that prevents scratching. All get a special CORNING surface treatment to improve plating efficiency, speed up confluency, and optimize cell recovery.

ORANGES BY THE CRATE

The free orange rack makes our tissue culture tubes easier to store and use. Plug seal cap is airtight, obviates pH shift. All tubes receive special CORNING surface treatment.



THE NEW ORANGE

The latest in a series of innovations, our new 850 cm² disposable roller bottle, increases surface area over 70%...the greatest surface area of any bottle designed for two-roller use. Like all our tissue culture products, the new orange is packed sterile, ready to use. Our little orange - 490 cm² is also available.

ORANGES


Orange means CORNING™ Cellware and that means innovations in product and packaging... a tradition since we originated the borosilicate PYREX® Erlenmeyer flask in 1915.

A SACK OF ORANGES

Our disposable tissue culture flasks come in three sizes – 25 cm², 75 cm², and 150 cm². All feature canted necks. The 75 cm² size also comes with a straight neck. Both necks in the 75 cm² size have a 14% bigger mouth for easier accessibility. The 75 cm² straight neck and the 150 cm² come with an orange plug seal or hard phenolic cap. All flasks feature accurate graduations, marking spots, and heavy walls.



JUICY OFFER



SE-I-2

TOTAL VALUE

We try diligently to see that you get the total of Corning's many capabilities in every product we make. We know that's how we earn your trust and retain your business. We know that you would never settle for less than total value from Corning.

THE JUICY OFFER

We'll send you a new catalog listing our CORNING Cellware products and information on how to obtain free samples of the products mentioned above. Just clip this coupon and attach it to your letterhead.


Corning Glass Works
Science Products Division
Mail Station 18, ATTN: LOC
Corning, NY 14830

Circle No. 161 on Readers' Service Card

CORNING™
cellware



1. Mobile stainless steel rack with suspended polycarbonate cages.



2. Disposable formed plenum filters shown above cages. Flat filter sheets also available.

See-Through™ suspended cage systems **The most versatile small animal housing systems available.**

(So many options...you "build" it to your own specs.)

Basically, this system involves a mobile stainless steel rack with suspended polycarbonate cages.

Two rack sizes are available: the large rack accommodates either 60 small cages, or 30 large cages, or *combinations* of both size cages. The *new* smaller rack size takes either 50 small cages, or 25 large cages, or *combinations* of both.

This See-Through™ system is ideal for housing rats, mice, hamsters, and guinea pigs.

The many other options

Cage sizes: 10½" x 9½" x 8" deep or 19" x 10½" x 8" deep.


Cage bottoms: solid or stainless steel wire mesh.

Walkfloor mesh: 2 x 2 (for rats) or 3 x 3 (for mice).

Food hoppers: Four available (for rats, mice, and guinea pigs).

Watering systems: Bottle or automatic (many variations available).

Filtering systems: New disposable formed plenum polyester spun-bonded filter that provides at least 100% more filter area than flat sheet filters. Effective in reducing cross-contamination and in minimizing the stress of a high ammonia atmosphere. Also available: autoclavable filters.



- 3.** Rack accommodates small or large cages or combinations of both.



- 4.** One available automatic watering system shown; other automatic systems or bottle watering also available.

Other important features/benefits

Of fundamental importance regardless of how you "build" your See-Through system initially to meet your current needs, its flexibility assures easy modification to your *future* requirements. And in an era where change is inevitable and budgets tight, this kind of foresight can substantially reduce "future shock."

No separate cage lids! Each cage is suspended from the shelf above, and that shelf is *perforated*. Accordingly, the *perforated* shelf itself serves as the "lid" for the cage immediately below. No initial lid cost. No subsequent handling. No replacement.

Visibility. As the See-Through name implies and the photos demonstrate, animal observation is completely unimpaired.

Excreta trays. Easily removable.

New Product!

See-Through II: meets ILAR's new housing specifications (see "Guide for the Care and Use of Laboratory Animals"). 8" deep cage with 200 sq. in. floor area at dead bottom houses 5 rats. 6" deep cage with 75 sq. in. floor at dead bottom houses 5 mice.

Circle No. 80 on Readers' Service Card


For more information about the most versatile small animal housing available, write or call Lab Products Inc., 365 W. Passaic St., Rochelle Park, N.J. 07662 (phone: 201/843-4600)

Why not also request our complete catalog?

lab products
inc a **bioMedic** company

Lab Products...not just plastic cages, metal cages, custom fabrication, laminar flow systems, bedding, automatic watering systems, accessories...

Copyright © BioMedic Corporation 1976



**"I have used Revco freezers
for the past 10 years. Nothing can compare
to their dependability."**

Rockville, Maryland

**"We've used Revco's in our department for years.
We have had excellent service."**

Iowa City, Iowa

"Reliability, high performance. . ."

Republic of Singapore

**"Completely satisfied with performance of
the first unit. It has been in
service over 10 years!"**

Windsor, Ontario

At Revco, we know we're doing our job when we receive comments like these. For over 40 years Revco has built its reputation around the world by designing and manufacturing freezers that are not only dependable but economical to operate.

Our concern for you, the customer, does not end with selling you what we feel to be the finest freezer obtainable. Revco goes far beyond normal industry standards with features such as listing by Underwriters Laboratories, Inc.; a nation-wide network of factory-trained people to assist and serve you; and a warranty that's second to none.

**When you need
dependable ULTra-Low® temperature,
you need Revco.**

(The comments appearing above were unsolicited and are on file at Revco for inspection.)

REVCO

The world's leader in ULTra-Low® temperature equipment.

Rheem Refrigeration Products Division

1100 Memorial Drive, West Columbia, S.C. 29169 Telephone: 803/796-1700 TWX: 810-666-2103 Cable: Revco


For additional information and a copy of our new color catalog, drop us a line or call. Or, contact one of our Regional Offices:

Northeast: (215) 547-0150

Midwest: (414) 781-3639

Western: (714) 661-6588

Circle No. 181 on Readers' Service Card



Choose Your Topic...

CA SELECTS Will Keep You Informed

You can't afford to be unaware of recently published literature in your field. But you can afford your own subscription to one or more publications of the CA SELECTS series — a current-awareness service from the publishers of CHEMICAL ABSTRACTS.

CA SELECTS contains informative abstracts of chemistry-related articles, patent literature, reports, and books from around the world. Delivered every two weeks, CA SELECTS provides you with a thorough, up-to-date survey of research in one well-defined area of chemistry — the area (or areas) you specify from the 76 choices now available.

HERE'S AN OPPORTUNITY TO EVALUATE AN ISSUE OF CA SELECTS IN YOUR AREA OF INTEREST — AT NO COST TO YOU! THE ISSUE YOU GET IS YOURS TO KEEP.

To receive a FREE issue of CA SELECTS, indicate the topic(s) pertaining to your interest from the list below. We'll also send complete ordering instructions, though you will not be obligated to order a subscription.

Clip This Form (or Mail a Photocopy) Today



NOW — 76 TOPICS

Please send me my free issue(s) of CA SELECTS, covering the topic(s) I have indicated below. Also send complete ordering instructions.

I understand this request does not obligate me to order a subscription.

- ☐ Adhesives
- ☐ Analytical Electrochemistry
- ☐ Animal Longevity & Aging
- ☐ Anti-Inflammatory Agents & Arthritis
- ☐ Antitumor Agents
- ☐ Atherosclerosis & Heart Disease
- ☐ Batteries & Fuel Cells
- ☐ Biogenic Amines & the Nervous System
- ☐ Biological Information Transfer
- ☐ Blood Coagulation
- ☐ Carcinogens, Mutagens, & Teratogens
- ☐ Catalysis (Applied & Physical Aspects)
- ☐ Catalysis (Organic Reactions)
- ☐ Chemical Hazards
- ☐ Coal Science & Process Chemistry
- ☐ Colloids (Applied Aspects)
- ☐ Colloids (Physicochemical Aspects)
- ☐ Computers in Chemistry
- ☐ Corrosion
- ☐ Crystal Growth
- ☐ Detergents, Soaps & Surfactants
- ☐ Drug & Cosmetic Toxicity
- ☐ Electrochemical Reactions
- ☐ Electrodeposition
- ☐ Electron & Auger Spectroscopy
- ☐ Electron Spin Resonance (Chemical Aspects)
- ☐ Energy Reviews & Books
- ☐ Engine Exhaust
- ☐ Environmental Pollution
- ☐ Flammability
- ☐ Flavors & Fragrances
- ☐ Food Toxicity
- ☐ Forensic Chemistry
- ☐ Fungicides
- ☐ Gas Chromatography
- ☐ Gaseous Waste Treatment
- ☐ Gel Permeation Chromatography

- ☐ Heat-Resistant & Ablative Polymers
- ☐ Herbicides
- ☐ High Speed Liquid Chromatography
- ☐ Infrared Spectroscopy
- ☐ Insecticides
- ☐ Ion Exchange
- ☐ β -Lactam Antibiotics
- ☐ Liquid Crystals
- ☐ Liquid Waste Treatment
- ☐ Mass Spectrometry
- ☐ Metallo Enzymes & Metallo Coenzymes
- ☐ Moessbauer Spectroscopy
- ☐ New Books in Chemistry
- ☐ Nuclear Magnetic Resonance (Chemical Aspects)
- ☐ Organic Reaction Mechanism
- ☐ Organofluorine Chemistry
- ☐ Organophosphorus Chemistry
- ☐ Organosilicon Chemistry
- ☐ Organo-Transition Metal Complexes
- ☐ Paper & Thin-Layer Chromatography
- ☐ Photobiochemistry
- ☐ Photochemistry
- ☐ Pollution Monitoring
- ☐ Prostaglandins
- ☐ Psychobiochemistry
- ☐ Radiation Chemistry
- ☐ Raman Spectroscopy
- ☐ Recovery & Recycling of Wastes
- ☐ Silver Chemistry
- ☐ Solar Energy
- ☐ Solid & Radioactive Waste Treatment
- ☐ Solvent Extraction
- ☐ Steroids (Biochemical Aspects)
- ☐ Steroids (Chemical Aspects)
- ☐ Substituent Effects & Linear Free Energy Relationships
- ☐ Surface Chemistry (Physicochemical Aspects)
- ☐ Thermochemistry
- ☐ Trace Element Analysis
- ☐ X-Ray Analysis & Spectroscopy

Mail Coupon To:

Chemical Abstracts Service
Marketing Department — SMA
P. O. Box 3012
Columbus, Ohio 43210 U.S.A.



Name _____

Title _____

Organization _____

Address _____

How much solar energy falls on the roof of a typical residence?

THE SUN

Our Future Energy Source

David K. McDaniels, *University of Oregon*

Now you can give your students a clear introductory understanding of the significant impact solar heating and cooling are having on our society. With this outstanding new volume, you can give them detailed, up-to-date treatment of solar radiation, flat-plate collectors, active and passive solar heating systems, solar thermal electric generating schemes, and photovoltaic devices, with a simplified examination of the overall energy problem, such as the energy crisis, fossil fuel reserves, and nuclear fission. You and your students will also appreciate the numerous illustrations and the annotated bibliographies at the end of every chapter. What's more, this new book gives your wide spectrum of nonscience majors the informative, comprehensive introduction they need—all in a nonmathematical approach.

(0 471 04500-4) 1979

approx. 304 pp. \$11.95 (tent.)

Why can you hear the sea in a conch shell?

THE FLYING CIRCUS OF PHYSICS WITH ANSWERS

Jearl Walker, *Cleveland State University*

At last—a popular physics text complete with answers to over 600 stimulating questions. Walker covers “naked-eye” physics by focusing on the everyday world of Frisbees, sounds of thunder, rainbows, sand dunes, soap bubbles...using such familiar objects as rubber bands, eggs, teapots, and Coke bottles. With the inclusion of answers, students get immediate feedback and see whether or not they understand the underlying principles. The many references to other, in-depth sources encourage your students to read more about the topics that interest them.

(0 471 02984-X) 1977 295 pp. \$8.95

WILEY TEXTS HAVE ANSWERS THAT WILL HAVE YOUR STUDENTS COMING BACK FOR MORE

Is music really just good vibrations?

PHYSICS AND THE SOUND OF MUSIC

John S. Rigden, *University of Missouri, St. Louis*

A clearly written and organized new book about the production, propagation, and perception of that sound we call music. Rigden develops the conceptual basis of the pure tone before moving on to discussions in which the musical instruments themselves become the sources of interest...the concert hall is treated as the environment affecting the transmitted sound...and the perceptual problems of consonance and dissonance are fully explored. The final part covers the subject of sound reproduction—the components of sound systems and those systems in their room environment.

(0 471 02433-3) 1977 286 pp. \$13.50

Is there one source for a solid introduction to physical science?

ESSENTIALS OF PHYSICAL SCIENCE

Kenneth R. Atkins, *University of Pennsylvania*,
John R. Holum, *Augsburg College*, &
Arthur N. Strahler

Here it is...the first triumvirate approach to teaching all of physical science from *one* text! Stressing the *phenomena* of science rather than its structural unity, three of Wiley's best known authors have combined their individual talents to give your students a thorough knowledge of the different fields of science. Each of the authors deals with a different aspect of nature, following a straight course from the physicist's world of the small and very small, through the small and medium world of the chemist, on through the immense worlds of the geologist and the astronomer.

Supplements available: Student Study Guide
Teacher's Manual

(0 471 03617-X) 1978 546 pp. \$17.95

To be considered for complimentary copies, write to Merrill Floyd, Dept. 3525. Please include course name, enrollment and title of present text.

Prices subject to change without notice.



JOHN WILEY & SONS, Inc.

605 Third Avenue
New York, N.Y. 10016

In Canada: 22 Worcester Road, Rexdale, Ontario

A 3525-11

Get horizontal...

- For electrophoresis of DNA and DNA fragments
- For DNA sequencing

Get horizontal...

with the Cancalco HSG Chamber for Horizontal Slab Gel Electrophoresis in agarose or polyacrylamide gel. Designed expressly for DNA research, the horizontal-format apparatus features a removable gel tray which provides support for low-mechanical-strength gels (0.18% agarose) during electrophoresis, staining, and photography. Gels are 19.5 cm wide by 40 cm long—long enough to resolve very complex mixtures in a single run. The chamber is capable of simultaneous analysis of 96 specimens; and it may be used preparatively, with gels up to 2.5 cm thick.



Circle No. 122 on Readers' Service Card

© 1978 MILES LABORATORIES, INC.

MRP-5002

☐ Please send more information about the new Cancalco® HSG Chamber for Horizontal Slab Gel Electrophoresis.

Name _____ (please print)
Institution _____
Address _____
City _____ State _____ Zip _____

Miles Laboratories, Inc.
Research Products Division
P.O. Box 2000
Elkhart, IN 46515
Tel: 219-264-8804

Miles Laboratories, Ltd.
P.O. Box 37
Stoke Court, Stoke Poges
Slough SL2 4LY
England
Tel: Farnham Common 2151


MRP-5002

The Cancalco HSG Chamber has numerous other features designed to help you perform electrophoresis effectively, including: Provision for buffer recirculation • Reusable gel wicks • Temperature control (up to +60°C for heat denaturation) • Safety interlock • Optional UV-

transparent gel tray (for use with ethidium bromide or other fluorescent stains)


Cancalco® HSG Chamber (#43-001-1) \$695
UVT Gel Tray (#43-002-1) \$135

New Cancalco® HSG Chamber



Get horizontal.

The Beckman gamma choice. Big enough to let you be choosy.



Beckman now offers the research and clinical laboratory a complete choice of gamma-counting instruments. Check the features below to find the one instrument that can most effectively satisfy your needs.

The low cost Gamma 4000 is a high performance, bench-top unit featuring ease of operation, 200 or 400 sample capacity, preset and variable windows, and a choice of print-out devices. The optional DP-5000 computer/printer provides on-line data reduction capability—for Radio-Immunoassay final answer calculations.

The Gamma 7000 is micro-processor controlled and offers high

performance features with push-button simplicity. Permanently stored counting programs eliminate all knobs, manual controls and switches to make operation fast, easy and more reliable. Simply select the desired program, load samples and depress the AUTO pushbutton—a completely simple and thorough approach to counting.

The Gamma 8000 takes micro-processor-control one step further and allows the user to design special counting programs, while maintaining the simplicity of operation unique to Beckman Instruments.

Finally, there's the Gamma 9000 with the ultimate in data-handling

capability. Built right into the 9000 is a complete computer that both controls the operation of the instrument and takes your data through to final answer.

Of course, if you're really going to be choosy about your next gamma instrument, you'll need complete technical information. Why not get it now?

Contact Scientific Instruments Division, Beckman Instruments, Inc., P.O. Box C-19600, Irvine, CA 92713.

Innovation you can count on

Circle No. 105 on Readers' Service Card

BECKMAN®



WITH OUR ARRAY PROCESSING MINIS, ALBERT WOULD HAVE DISCOVERED RELATIVITY AND STILL HAD TIME TO EXPLAIN IT.

Extra time is what our two integrated array processing minicomputers, AP/130 and ECLIPSE S/250, give you. They're a lot faster than plain minicomputers. Either Data General system can perform 1024 point complex floating point FFTs in less than 10 milliseconds. The ECLIPSE S/250 computer goes even further. It's capable of performing double precision adds in 600 nanoseconds and double precision SINE in 25.3 microseconds. Which makes it the world's fastest FORTRAN minicomputer.

Data General integrated array processing minis save money as well as time. The AP/130 actually costs less than some separate AP units. Less, in fact, than many comparable general purpose minicomputers with-

out array processing. And either system saves you the expense and trouble of designing and programming a special interface.

System integration has other advantages. For example, our CPU can address AP memory directly. Software development is also completely integrated, with the array processing capability fully supported under RDOS and FORTRAN 5. And of course you get Data General design, engineering, testing and support for the whole system. When you need the speed of an array processor built together with the flexibility of a mini for a wide range of scientific and research applications, such as digital signal processing, structural design and simulation, and image processing,

we can make you look like a genius. Return the coupon below for AP/130, S/250 details, or call (617) 366-8911, Ext. 4735.

Send information on			<input type="checkbox"/> Array Processors	<input type="checkbox"/> S/250
<input type="checkbox"/> Scientific Computers			<input type="checkbox"/> Have someone contact me.	
Name _____				
Title _____				
Company _____				
Tel. No. _____				
Address _____				
City _____		State _____	Zip _____	
<small>Data General Corporation, Westboro, MA 01581, (617) 366-8911 Data General (Canada) Ltd., Ontario, Data General Europe, 61 rue de Courcelles, Paris, France, 766.51.78. Data General Australia (03) 89-0633. Data General Ltda. Sao Paulo, Brazil 543-0138. Data General Middle East, Athens, Greece, 952-0557. © Data General Corporation, 1978. ECLIPSE is a registered trademark of Data General Corporation. S2979</small>				

Data General

We make computers that make sense.

Circle No. 180 on Readers' Service Card

Seek peaks at 206 nm
and get up to 200x
the sensitivity of
monitoring at 280 nm...



...with the new Uvicord S UV monitor

Sensitivity is increased up to 200x for proteins when you monitor at 206 nm with LKB's new Uvicord® S UV-monitor. This unique instrument will detect non-aromatic peptides, polysaccharides, nucleotides, lipids and steroids as well as proteins. And, naturally, you can also monitor at 254 or 280 nm.

Enhanced versatility has required no compromise in stability. Quite the contrary. Sophisticated optics and solid state circuitry provide outstanding linearity. And you can monitor simultaneously at high and low sensitivities.

Unlike others, the new Uvicord S UV-monitor is contained in a single small case which mounts easily on a fraction collector or ring stand. And its low price matches its small size.

Contact LKB today for full details.

LKB

LKB Instruments Inc.
12221 Parklawn Drive, Rockville, Maryland 20852
301: 881-2510

Circle No. 165 on Readers' Service Card

3A-300

LETTERS

EPA Air Pollution Standard

Eliot Marshall's article (News and Comment, 1 Dec. 1978, p. 949) on the proposed Environmental Protection Agency (EPA) standard for photochemical oxidants raises several serious questions, some in a misleading light.

Marshall suggests EPA is proceeding in the face of opposition from the scientific community and without consultation. Having been a scientific reviewer of the "Photochemical oxidant criteria document," an author of the "Oxides of nitrogen document," and one of those presently preparing the "Sulfur oxides and particulate document," I find this contention incredulous. The former two documents were circulated in at least two drafts, and criticisms were incorporated from both EPA's science advisory board (SAB) and outside reviewers. True, there were and are differences of opinion but, to my knowledge, the photochemical oxidant document was never "rejected" by the SAB. Marshall does not distinguish between a request for revision and rejection. The criteria documents are immense efforts by a broad spectrum of scientists. EPA has taken the proper position that, because they are public documents, comments by the public should be incorporated as far as possible. To my knowledge, this has happened. To suggest that the contributions of nongovernment scientists to these documents have been discarded is a slur on the efforts of many.

The issue of using ozone in place of total oxidants as a pollution index has been, and continues to be, controversial. As Marshall points out, some scientists, such as James Pitts, Jr., regard this move as impractical. I both agree with Pitts' contention that ozone alone is not a sufficient standard for polluted air and contend that ozone is a reasonable interim approach to regulation as proposed by EPA. Ozone has been better studied than other photochemical pollutants and its toxicity is better understood. The toxicity of the complex mixture of polluted air is poorly understood, and there is simply not enough data available to make a reasonable judgment or defense of a more general standard.

While no one can fault obtaining a better inventory of the toxicants in air, the Clean Air Act amendments require action before added data can be obtained. I am not so pessimistic about support for added research on the composition of

SCIENCE, VOL. 203

polluted air or its health effects. The storm clouds of the new budget may prove that I am incorrect.

A critical issue is the safety of a higher air quality standard than the present 0.08 part per million (ppm) of pollutant. Present data suggest that chronic toxicity will occur at a level of at least 0.5 ppm. To ignore the animal toxicity data supporting this contention is to ignore the basic concept of prediction of safety from experimental studies as a whole. Human pulmonary function tests are much less sensitive than biochemical, physiological, and morphological studies in animals and thus tend to underestimate, short of lifetime exposures, the hazards from ozone. A standard of 0.2 ppm will provide little, or no, safety margin. More important, the effects of cyclical, short-term exposure regimens, mimicking human exposures in urban areas, are lacking.

The gaps in critical data in air pollutant toxicology are not for want of trying. Rather, they reflect the overall sophistication of toxicology and the improved sensitivity of pathophysiological measurements. What is needed is stronger support for intensive, long-term scientific studies of the highest quality. If the initiative to clean up the environment is to survive, EPA and other agencies regulating the outpouring of pollutants need the support and contribution of the scientific community as a whole.

DANIEL B. MENZEL

Departments of Pharmacology and Medicine, Duke University Medical Center, Durham, North Carolina 27710

Federal Regulation:

What Role for Universities?

Eliot Marshall (News and Comment, 1 Dec. 1978, p. 955) reports on Massachusetts Institute of Technology president Jerome B. Wiesner's concern over federal intrusion into internal university operation and management. The suggested solution appears to be the establishment, or reestablishment, of a special "federal-academic relationship" that would exempt universities from the enforcement of government regulations applied to others.

One concerned about the welfare of our universities can only sympathize with the problem Wiesner raises and applaud him and others who speak out. But one concerned about the overall welfare of our society may also ask whether it is responsible for the universities to plead

for special treatment. After all, government regulations are intended to help achieve ideals to which all of us ascribe: equality of opportunity regardless of race, religion, sex, or physical handicap; prevention of cruelty to animals; control or elimination of activities with risk to health or safety; fair labor practices; sound financial management—especially of public funds; decent treatment of the aged, and so forth, and so forth.

Do universities have a special role that should exempt them from trying to reach these ideals? Might one not argue, instead, that universities should show exemplary conformance with government regulations intended to help achieve these goals?

Let me suggest that it is not the goals, but unbalanced and intrusive federal actions connected with the achievement of those goals, that are the problem. If this is so, shouldn't the universities, on the basis of their own experience, be seeking to change this federal approach wherever it is applied?

Is it responsible for those in universities—so long as they are not directly affected—to ignore, or even support, the same regulatory approach they find so "repressive"?

BERTRAM WOLFE

*15453 Via Vaquero,
Monte Sereno, California 95030*

Problems

I am frankly baffled by Soviet mathematician Grigori Freiman's reported bafflement. Gina Bari Kolata, in her article on alleged anti-Semitism in Soviet mathematics (News and Comment, 15 Dec. 1978, p. 1167), describes Freiman as being amazed that a Soviet Jewish student could solve the "difficult" problem, "Which is greater, $\sin 8/7$ or $8\pi/7$?" Either the problem is garbled, or I am missing something, because the relationship sought is elementary (easier, I'd judge, than proving the problem alleged to be less difficult—that " $\log_2 3$ " is irrational). So my curiosity is aroused. What was the problem?

GEORGE GREENE

*524 Oakdale Avenue,
Chicago, Illinois 60657*

Both problems were given incorrectly in my article. The difficult problem was to decide, without the use of tables, which is greater, $\sin 8/7$ or $8\pi/27$. The other problem was to prove that $\log_2 3$ is irrational.—GINA BARI KOLATA

New!

A continuously variable length chromatography column



For gel filtration, ion exchange, affinity and adsorption chromatography—the new LKB 2137 column adjusts to any length from 5 to 95 cm. No longer do you need different length columns for each technique. Now, simplify inventory, cut costs dramatically!

In minutes you can have a completely jacketed high-resolution column of optimum length. For either descending or ascending flow. For sharper, more symmetrical peaks than ever before.

This versatile column comes complete with flow adaptors. The price? Less than you pay for other suppliers' flow adaptors alone.


*Contact LKB today
for full details.*

LKB

LKB Instruments Inc.
12221 Parklawn Drive
Rockville, Maryland 20852
301: 881-2510

13A-301

Circle No. 166 on Readers' Service Card



3367R

The first inverted microscope designed as an inverted microscope.

The new Leitz Diavert® was engineered from the base up to meet the special requirements of inverted microscopy. It is exceptionally stable with a large rectangular base area. The stand will support large or small vessels and photomicrographic equipment, including the Leitz Orthomat® fully automatic camera and the Leicina® Super 8 movie camera for time lapse studies.

The large object stage is positioned low so you can see what you are doing while manipulating the specimen.

The image moves in the same direction as the specimen. Reason: image-erecting optics. The transmitted light source moves with the stage so it doesn't have to be refocused when

you change the specimen plane.

Special long-working-distance, brightfield as well as phase objectives, and also long-working-distance condensers (up to 64mm) are important features.

The Diavert is a building block system so you can readily employ all the techniques of light microscopy, including new interference contrast and incident fluorescence accessories. If you own a Leitz Ortholux II research microscope, you can interchange all Ortholux accessories.

For information call 201-767-1100 or write E. Leitz, Inc., Dept. SC 29 Rockleigh, N.J. 07647.

Leitz

Where most new developments start.

Circle No. 129 on Readers' Service Card

AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE

Science serves its readers as a forum for the presentation and discussion of important issues related to the advancement of science, including the presentation of minority or conflicting points of view, rather than by publishing only material on which a consensus has been reached. Accordingly, all articles published in *Science*—including editorials, news and comment, and book reviews—are signed and reflect the individual views of the authors and not official points of view adopted by the AAAS or the institutions with which the authors are affiliated.

Editorial Board

1979: E. PETER GEIDUSCHEK, WARD GOODENOUGH, N. BRUCE HANNAY, MARTIN J. KLEIN, FRANKLIN A. LONG, NEAL E. MILLER, JEFFREY J. WINE
1980: RICHARD E. BALZHISER, WALLACE S. BROECKER, CLEMENT L. MARKERT, FRANK W. PUTNAM, BRYANT W. ROSSITER, VERA C. RUBIN, MAXINE F. SINGER, PAUL E. WAGGONER, F. KARL WILLENBROCK

Publisher

WILLIAM D. CAREY

Editor

PHILIP H. ABELSON

Editorial Staff

Managing Editor ROBERT V. ORMES
Business Manager HANS NUSSBAUM

Assistant Managing Editor JOHN E. RINGLE
Production Editor ELLEN E. MURPHY

News Editor: BARBARA J. CULLITON

News and Comment: WILLIAM J. BROAD, LUTHER J. CARTER, CONSTANCE HOLDEN, ELIOT MARSHALL, DEBORAH SHAPLEY, R. JEFFREY SMITH, NICHOLAS WADE, JOHN WALSH. **Editorial Assistant,** SCHERRAINE MACK

Research News: BEVERLY KARPLUS HARTLINE, FREDERICK F. HARTLINE, RICHARD A. KERR, GINA BARI KOLATA, JEAN L. MARX, THOMAS H. MAUGH II, ARTHUR L. ROBINSON. **Editorial Assistant,** FANNIE GROOM

Consulting Editor: ALLEN L. HAMMOND

Associate Editors: ELEANORE BUTZ, MARY DORFMAN, SYLVIA EBERHART, JUDITH GOTTLIEB, RUTH KULSTAD

Assistant Editors: CAITILIN GORDON, LOIS SCHMITT, DIANE TURKIN

Book Reviews: KATHERINE LIVINGSTON, **Editor:** LINDA HEISERMAN, JANET KEGG

Letters: CHRISTINE KARLIK

Copy Editor: ISABELLA BOULDIN

Production: NANCY HARTNAGEL, JOHN BAKER; YA LI SWIGART, HOLLY BISHOP, ELEANOR WARNER; JEAN ROCKWOOD, LEAH RYAN, SHARON RYAN

Covers, Reprints, and Permissions: GRAYCE FINGER, **Editor:** CORRINE HARRIS, MARGARET LLOYD

Guide to Scientific Instruments: RICHARD SOMMER

Assistant to the Editors: RICHARD SEMIKLOSE

Membership Recruitment: GWENDOLYN HUDDLE

Member and Subscription Records: ANN RAGLAND
EDITORIAL CORRESPONDENCE: 1515 Massachusetts Ave., NW, Washington, D.C. 20005. Area code 202. General Editorial Office, 467-4350; Book Reviews, 467-4367; Guide to Scientific Instruments, 467-4480; News and Comment, 467-4430; Reprints and Permissions, 467-4483; Research News, 467-4321. Cable: *Advancesci*, Washington. For "Instructions for Contributors," write the editorial office or see page xi, *Science*, 29 September 1978.

BUSINESS CORRESPONDENCE: Area Code 202. Business Office, 467-4411; Circulation, 467-4417.

Advertising Representatives

Director: EARL J. SCHERAGO

Production Manager: MARGARET STERLING

Advertising Sales Manager: RICHARD L. CHARLES

Marketing Manager: HERBERT L. BURKLUND

Sales: NEW YORK, N.Y. 10036: Steve Hamburger, 1515 Broadway (212-730-1050); SCOTCH PLAINS, N.J. 07076: C. Richard Callis, 12 Unami Lane (201-889-4873); CHICAGO, ILL. 60611: Jack Ryan, Room 2107, 919 N. Michigan Ave. (312-DE-7-4973); BEVERLY HILLS, CALIF. 90211: Winn Nance, 111 N. La Cienega Blvd. (213-657-2772); DORSET, VT. 05251: Fred W. Dieffenbach, Kent Hill Rd. (802-867-5581)

ADVERTISING CORRESPONDENCE: Tenth floor, 1515 Broadway, New York, N.Y. 10036. Phone: 212-730-1050.

Scientific Relations with China

President Carter's decision to normalize relations with China has met a mixed public response. However, in scientific and educational circles the reaction is more favorable. There one finds awareness and appreciation of the intellectual potentials of the Chinese, of the great inventive contributions during past millennia, and a feeling that the United States should be on friendly terms with a nation that includes one-fourth of the world's population. U.S. scientists are disposed to give Chinese visiting scholars a warm welcome and many are eager to spend some time in China. The extent of interchange and cooperation will be determined by the Chinese rather than by reluctance here.

Current official attitudes of the Chinese toward interaction are in sharp contrast to those that prevailed only a few years ago. From 1948 to 1971 there were few contacts. Then a gradual change occurred in which science and technology had a major role. A pivotal moment was the Nixon visit to China in 1972. But even before that, scholarly groups here were trying to foster relationships. In 1966 a Committee on Scholarly Communication with Mainland China was formed to explore and pursue opportunities to facilitate scientific and other scholarly communication. The Committee included four major sectors of American scholarly activities—the American Council of Learned Societies, the National Academy of Engineering, the National Academy of Sciences, and the Social Science Research Council.

Early efforts to achieve a working relationship with China were not fruitful, but in 1972 the Committee, renamed the Committee on Scholarly Communication with the People's Republic of China, became a major factor in exchanges. From 1972 through 1978 the Committee and its staff served as hosts to 37 Chinese scientific delegations to the United States. Since the spring of 1973 the Committee has sent 30 American scholarly delegations to China, selecting the exchange topics and participants.

This organized effort was supplemented by private groups and individuals. Notably effective were Chinese-Americans, some of whom prepared trip reports that were published. The various visits made clear to the Chinese that a good and friendly climate for interchange existed.


In July 1978 President Carter sent his science adviser, Frank Press, and a delegation of heads of science and technology-related government agencies to China. Their mission was to investigate the possibility of cooperation going beyond the survey trips. Press had earlier served as chairman of the Committee on Scholarly Communication with the People's Republic of China and he had been a member of an American visiting delegation. This background was helpful in discussions leading to opening the way for serious scholarly cooperation. The Press mission probably had a significance that went beyond scholarship. When two major powers make an agreement, many factors are involved. But in the normalization of relations, considerations of scientific and technological interchange surely were an important factor.

With the resumption of diplomatic relations between the United States and China, interchange will occur through other mechanisms such as commercial relationships and exchanges between universities. To obtain information to better formulate a cooperative program, most of the members of the Board of Directors of AAAS visited China from 14 November to 3 December 1978. The delegation was able to examine a broad cross section of Chinese life, education, and science and technology. In the course of the visit, many hundreds of photographs were taken both on film and in the mind. In this issue there are a few pictures and impressions of changing scenes in China in late 1978. What is presented is somewhat different, less euphoric than the views of earlier visitors. In turn, it is to be expected that what will be reported a few years hence will differ from the present portrayal, for the vast energy of China is being channeled in new directions and change for better or worse is inevitable.—PHILIP H. ABELSON

15,000 RPM (12,800 xG) REGARDLESS OF LOAD

**With some micro centrifuges,
constant speed is a constant problem,
but not with an Eppendorf®.**

Within ten seconds, the Eppendorf 5412 attains 15,000 rpm, generating a force of 12,800xG, regardless of load. Performance like this means rapid sample separation, in most cases within 60 seconds or less. An angled rotor accommodates twelve disposable Eppendorf 1.5ml micro test tubes, or twelve 500 μ l, 400 μ l or 250 μ l tubes using adapters. (For higher capacity requirements, the Eppendorf Model 5413 accepts forty 1.5ml, 400 μ l or 250 μ l disposable test tubes in four carriers, but operates at lower speeds.) Eppendorf Micro Centrifuges are equipped with automatic 15 min. timer, safety switch (prevents operation with lid open) and safety lid lock (lid stays locked while rotor is spinning).




Eppendorf Micro Test Tubes have attached caps and are ideal for centrifuging, mixing, or storing reagents. Economically priced, they are available in the following sizes: polypropylene — 1.5ml, 500 μ l, 400 μ l; polyethylene — 400 μ l, 250 μ l.

For complete literature, write: Eppendorf Division, Brinkmann Instruments, Inc., Cantiague Road, Westbury, N.Y. 11590.
In Canada: Brinkmann Instruments (Canada), Ltd.



Eppendorf Micro Centrifuges





Spectrophotometer

The Spectromatic system includes a spectrophotometer, a thermal flow cell, a pump, and a controller. A thermal printer identifies each sample, prints data, and alerts the operator to errors. A quality-control feature, called Slope Check, assures stability of electronic function and of the reagents. The range of coverage is 325 to 650 nanometers (± 2 nanometers). The operator selects readout in transmittance, absorbance, or concentration. Fisher Scientific. Circle 778.

Chromatogram Solvent Equilibration

A solvent vapor equilibration tank holds up to six plates, 0.125-inch thick, or up to 12 thinner plates. The apparatus includes a rectangular glass developing jar, a special glass cover with two Teflon collars and glass hooks, and anodized aluminum multiplate rack. Distortion of the chromatogram may be avoided if the layers are suspended above the solvent inside the closed tank during solvent equilibration. The retractable glass hooks are used to raise and lower the plate rack inside the closed tank. After the tank's atmosphere is saturated with vapor, the racks and plates are lowered into the solvent. Scientific Manufacturing Industries. Circle 779.

Ultrasonic Processor

Model 190 is a 900-watt system for processing dispersions, emulsions, or solutions. Electrical energy is transformed from 60 to 20,000 hertz and converted to ultrasonic vibration. This mechanical en-

ergy is coupled to a horn that transmits the energy to the sample. Controls include an adjustable power output, a timer, and an output meter. The converter operates without a fan to eliminate the possibility of cross contamination. Applications include dispersion, emulsification, extraction, separation, scrubbing, leaching, mixing, cell disruption, disintegration, and others. Branson Sonic Power. Circle 780.

Calculator

Model EL-502 is designed for scientific applications. It handles inverse trigonometric and logarithmic functions and has an independently accessible memory. It features an eight-digit display or a five-digit mantissa and two-digit exponent. Its light weight and small size make it highly suitable for bench or field calculations. Sharp Electronics. Circle 782.

Pressure Meter for HPLC

An electronic pressure meter functions from a remote transducer of small volume and flowthrough design. It has facilities for setting upper and lower limits. When these limits are reached and the transducer is tripped, the power is isolated from two outlets on the back of the meter. An analog outlet is also provided. Pressure is read on a digital display in the range of 0 to 5 pounds per square inch. Phase Separations. Circle 784.

Digital Image Analyzer

The MOP system allows the operator to interact with the high-speed calculation and memory of its microcomputer to reduce significantly the systematic errors of fully automatic image analysis. The system is modular and the three basic models offered are suited to quantitative analysis of length, area, diame-

ter, orientation, count, form, projected length, and center of mass coordinates. Images from photographs, x-rays, drawings, and other sources may be analyzed or the system may be used directly as an adjunct to microscopy. The system consists of a compact hybrid microcomputer and a multipurpose measuring tablet. It accurately measures up to eight parameters simultaneously and can summarize and store data in up to 20 channels. Carl Zeiss. Circle 781.

Chromatography Column

Model 2137 offers continuous variability of length. Bench workers may select lengths from 5 to 95 cm. Three interchangeable borosilicate glass tubes, 35, 65, and 95 centimeters long, are available in diameters of 16 or 26 millimeters. A compartmented storage box also houses three polymethylmethacrylate jacket sections with locking rings and seals, two flow adaptors with snap-on filters, two special end fittings, a three-way valve, sample syringe, and polyethylene tubing. Assembly is rapid and, by changing tubes and adjusting flow adaptors, one converts lengths easily. Ascending or descending flows are permitted and the filters minimize dead volume. LKB Instruments. Circle 783.

Literature

Liquid Enzyme Reagents describes the Liquid-STAT line of stable ultraviolet reagents that are ready for use without reconstitution for a variety of clinical assays. Beckman Instruments. Circle 785.

Computers for Science is devoted to the Eclipse S/250, a system of software, hardware, and peripherals for research applications in the laboratory. Data General. Circle 788.


Cryogenic Controls catalogs a product line for low-temperature applications including control of liquid level, temperature alarms, temperature controls, freezing chambers, and more. Cryo-Med. Circle 789.

Bio-Chemical Research Instruments includes pH meters, titrators, chromatography apparatus, Kjeldahl analyzers, homogenizers, peristaltic pumps, and colorimeters. Brinkmann Instruments. Circle 790.

Lectins and Related Biochemicals describes a new line of biochemicals that depend on the high-affinity interaction of biotin and avidin. Vector Laboratories. Circle 791.

Newly offered instrumentation, apparatus, and laboratory materials of interest to researchers in all disciplines in academic, industrial, and government organizations are featured in this space. Emphasis is given to purpose, chief characteristics, and availability of products and materials. Endorsement by *Science* or AAAS is not implied. Additional information may be obtained from the manufacturers or suppliers named by circling the appropriate number on the Readers' Service Card (on pages 502A and 566A) and placing it in the mailbox. Postage is free.

—RICHARD G. SOMMER



You are not alone.

And we are happy to share the experience we've accumulated over the years. Request our helpful technical bulletin on the uses of our high-specific-activity Bolton-Hunter Reagent.

Bolton-Hunter Reagent [125 I]

~2000Ci/mmol NEX-120
Monoiodinated: Prepared weekly
~4000Ci/mmol NEX-120H
Diiodinated: Prepared biweekly

Not for use in humans or clinical diagnosis.



New England Nuclear

549 Albany Street, Boston, Mass. 02118

Call toll-free: 800-225-1572

(In Massachusetts and International: 617-482-9595)

NEN Chemicals GmbH, Dreieich, W. Germany; NEN Canada Ltd., Lachine, Quebec

Circle No. 164 on Readers' Service Card



Hellma—the largest assortment of highest precision glass and quartz cells.

Standard · Flow-through · Constant-temperature Anaerobic · Special Designs

Also available—ULTRAVIOLET LIGHT SOURCES
Deuterium Lamps · Mercury Vapor Lamps
Hollow Cathode Lamps · Power Supplies



Write for literature
Box 544

Borough Hall Station
Jamaica, New York 11424
Phone (212) 544-9534

Circle No. 149 on Readers' Service Card

ELECTRODE IMPEDANCE TESTER

For EEG, EOG, EMG, EKG, evoked potentials and wherever surface electrodes are used to record bioelectric events.



Compact, portable instrument quickly and conveniently checks surface electrode impedances. Electrodes are simply plugged in and a button pushed. Color-coded meter shows at a glance whether impedance is in the "safe", "doubtful" or "unacceptable" range. Alternating test current prevents electrode polarization. Multiple safety features guard against test current exceeding medical safety standards.



Life-Tech Instruments, Inc.

BOX 36221 • HOUSTON, TEXAS 77036 • (713) 783-6490

Circle No. 220 on Readers' Service Card



Why Sol[®] is the right small computer for scientists and engineers.

Sol[®] small computer systems offer the flexibility you need to analyze a wide range of data and solve complex problems.

- Language Flexibility. Sol small computers use Extended BASIC, FORTRAN, FOCAL, PILOT and Assembler on both disk and cassette.

- Precision Flexibility. Extended Disk BASIC is available with 6, 8, 10, 12, 14 or 16 digits of precision.


- Memory Storage Flexibility. Up to 1.5 million bytes of disk and 64K bytes of onboard memory fully supported by an outstanding disk operating system.

For a full report on why Sol stands up to the minis, address Box RD, Processor Technology Corporation, 7100 Johnson Drive, Pleasanton, CA 94566, or see your nearest dealer.

Toll-free Dealer Locator Hotline, (800) 227-1241. In California call toll-free (800) 972-5951.

Processor Technology

Circle No. 217 on Readers' Service Card



Little things mean a lot in HPLC

Inlet filters. Rheodyne inlet filters can be connected between the sample injection valve and the column to protect the column from plugging. The 2 micron filter element prevents plugging caused by particles in the samples or by injection valve wear particles. Pressure rating of the filter assembly is 7000 psi (500 bar). Only type 316 stainless steel and PTFE contact the stream.

U.S. prices are \$45 for the Model 7302 Column Inlet Filter and \$20 for a package of 5 2 μ m filter elements and gaskets.


Pressure Relief Valve. Rheodyne's Model 7037 Pressure Relief Valve protects your equipment against damage from over-pressure. You can set it anywhere in the range of 2000 to 7000 psi (140 to 500 bar). U.S. price is only \$270. So, when you consider the fact that you can blow a pressure gauge in less than 1 second, it's a real bargain.

Teflon Rotary Valves. Type 50 Rheodyne Valves are real workhorse accessories for your LC equipment.

Use them for sample injection, column switching, recycling, reagent switching, fraction collection, stream sampling and quantitative reagent injection. Available in four different versions, these valves are chemically inert with zero dead volume, operate at 300 psi. They are offered in 0.8 or 1.5 mm bore and in either manual or automatic versions.

U.S. price of 0.8 mm bore 3 and 4-way valves is \$80. The 6-position and sample injection valves are priced at \$95. Cost of 1.5 mm bore valves is \$2 more.

Write for more data. For full information, please address Rheodyne, Inc., 2809 Tenth Street, Berkeley, CA 94710. Phone (415) 548-5374.


RHEODYNE
THE LC CONNECTION COMPANY
Circle No. 168 on Readers' Service Card

BOOKS RECEIVED AND

BOOK ORDER SERVICE

(Continued from page 540)

Operating Handbooks, vol. 2. *To order this book circle No. 580 on Readers' Service Card*

Coniferen im westlichen malayischen Archipel. Studien zu ihrer Verbreitung und Ökologie innerhalb der vollhumiden südostasiatischen Tropen. N. Stein. Junk, The Hague, 1978 (U.S. distributor, Kluwer Boston, Hingham, Mass.). xii, 168 pp., illus. + plates. \$29.50. Biogeographica, vol. 11.

Consciousness and Self-Regulation. Advances in Research and Theory. Vol. 2. Gary E. Schwartz and David Shapiro, Eds. Plenum, New York, 1978. xxiv, 446 pp. \$22.50.

Contrary to Nature. Being an Illustrated Commentary on Some Persons and Events of Historical Importance in the Development of Knowledge Concerning Cancer. Michael B. Shimkin. National Institutes of Health, Bethesda, Md., 1977 (available from the Superintendent of Documents, Washington, D.C.). xx, 498 pp., illus. \$12.75. DHEW Publication No. (NIH) 76-720.

Desert Tortoise Council Proceedings of 1977 Symposium. Las Vegas, Mar. 1977. Desert Tortoise Council, San Diego, Calif., 1978. viii, 188 pp., illus. Paper, \$5.

Developmental Neurobiology. Marcus Jacobson. Plenum, New York, ed. 2, 1978. xii, 562 pp., illus. \$37.50.

Dynamic Morphology of Leukemia Cells. A Comparative Study by Scanning Electron Microscopy and Microcinematography. H. Felis, G. Haemmerli, and P. Sträuli. Springer-Verlag, New York, 1978. xiv, 192 pp., illus. \$49. *To order this book circle No. 581 on Readers' Service Card*

The Ecology of Social Boundaries. Agta Foragers of the Philippines. Jean Treloggen Peterson. University of Illinois Press, Urbana, 1978. xiv, 142 pp., illus. + plates. \$8.95. Illinois Studies in Anthropology, No. 11.

Economic Microbiology. Vol. 2, Primary Products of Metabolism. A. H. Rose, Ed. Academic Press, New York, 1978. xviii, 470 pp., illus. \$44.

Epidemiology of Cancer in Texas. Incidence Analyzed by Type, Ethnic Group, and Geographic Location. Eleanor J. Macdonald and Evelyn B. Heinze. Raven, New York, 1978. xx, 650 pp. \$45.

Experimental Psychology. A Methodological Approach. F. J. McGuigan. Prentice-Hall, Englewood Cliffs, N.J., ed. 3, 1978. xvi, 496 pp., illus. \$13.95.

Faanakkar: Cultural Values in a Micro-nesian Society. John L. Caughey. University of Pennsylvania Department of Anthropology, Philadelphia, 1977. x, 180 pp., illus. Paper, \$17. University of Pennsylvania Publications in Anthropology, No. 2.

Fatty Acids and Glycerides. Anis Kuksis, Ed. Plenum, New York, 1978. xviii, 470 pp., illus. \$35. Handbook of Lipid Research, vol. 1.

Handbook of Sex Therapy. Joseph LoPiccolo and Leslie LoPiccolo, Eds. Plenum, New York, 1978. xx, 532 pp., illus. \$29.50. Perspectives in Sexuality.

Highly Hazardous Materials Spills and Emergency Planning. J. E. Zajic and W. A. Himmelman. Dekker, New York, 1978. xii, 226 pp. \$24.50. Hazardous and Toxic Substances, 1.

Historical Archaeology. A Guide to Substantive and Theoretical Contributions. Rob-

ert L. Schuyler, Ed. Baywood Publishing Company, Farmingdale, N.Y., 1978. xiv, 286 pp., illus. \$15.

Historical Studies in the Physical Sciences. Ninth Annual Volume. Russell McCormmach, Lewis Pyenson, and Roy Steven Turner, Eds. Johns Hopkins University Press, Baltimore, 1978. xii, 362 pp. \$18.95.

Humanscape. Environments for People. Stephen Kaplan and Rachel Kaplan, Eds. Duxbury Press (Wadsworth), North Scituate, Mass., 1978. xvi, 480 pp., illus. Paper, \$9.95.

Hyperspaces of Sets. A Text with Research Questions. Sam B. Nadler, Jr. Dekker, New York, 1978. xvi, 708 pp. \$45. Monographs and Textbooks in Pure and Applied Mathematics, 49.

Image and Pilgrimage in Christian Culture. Anthropological Perspectives. Victor Turner and Edith Turner. Columbia University Press, New York, 1978. xx, 282 pp. + plates. \$15. Lectures on the History of Religions.

In the Center of Immensities. Bernard Lovell. Harper and Row, New York, 1978. xx, 172 pp. \$10.95. World Perspectives, vol. 53.

Indicators of Viruses in Water and Food. Gerald Berg, Ed. Ann Arbor Science, Ann Arbor, Mich., 1978. viii, 424 pp. \$39.95.

Inelastic Electron Tunneling Spectroscopy. Proceedings of a conference, Columbia, Mo., May 1977. T. Wolfram, Ed. Springer-Verlag, New York, 1978. viii, 242 pp., illus. \$29. Springer Series in Solid-State Sciences, vol. 4. *To order this book circle No. 582 on Readers' Service Card*

The Integrated Mind. Michael S. Gazzaniga and Joseph E. LeDoux. Plenum, New York, 1978. xii, 168 pp., illus. \$15.95.

International Health. Paul F. Basch. Oxford University Press, New York, 1978. viii, 380 pp., illus. Cloth, \$16.95; paper, \$9.95.

An Introduction to Medical Genetics. J. A. Fraser Roberts and Marcus E. Pembrey. Oxford University Press, New York, ed. 7, 1978. xviii, 324 pp., illus. Paper, \$13.95.

An Introduction to the Aquatic Insects of North America. Richard W. Merritt and Kenneth W. Cummins, Eds. Kendall/Hunt, Dubuque, Iowa, 1978. xiv, 442 pp., illus. \$18.95.

Introduction to X-Ray Spectrometric Analysis. Eugene P. Bertin. Plenum, New York, 1978. xiv, 486 pp., illus. \$23.50.

Laboratory Manual for Invertebrate Zoology. Emphasizing Marine Forms. Ju-Shey Ho. Illustrated by David M. Marschall. Hwong Publishing Co., Los Alamitos, Calif., 1978. viii, 152 pp. Paper, \$7.95.

Linguistic Theory and Psychological Reality. Morris Halle, Joan Bresnan, and George A. Miller, Eds. MIT Press, Cambridge, Mass., 1978. xx, 330 pp., illus. \$17.50. MIT Bicentennial Studies, 3.

The Mammals of the Palearctic Region. A Taxonomic Review. C. B. Corbet. British Museum (Natural History), London, and Cornell University Press, Ithaca, New York, 1978. viii, 314 pp., illus. \$38.50.

Man-Made Radio Noise. Edward N. Skomal. Van Nostrand Reinhold, New York, 1978. x, 342 pp., illus. \$19.95. *To order this book circle No. 583 on Readers' Service Card*

The Man Who Rode Sharks. William R. Royal with Robert F. Burgess. Dodd, Mead, New York, 1978. xii, 254 pp. + plates. \$8.95.

Methods in Pharmacology. Vol. 4B, Renal Pharmacology. Manuel Martinez-Maldonado, Ed. Plenum, New York, 1978. x, 404 pp., illus. \$39.50.

The Mind Unfolded. Essays on Psychology's Historic Texts. Daniel N. Robinson. Uni-

versity Publications of America, Washington, D.C., 1978. xii, 544 pp. Cloth, \$24; paper, \$8.

Mineral Economics and Basic Industries in Asia. K. P. Wang and E. Chin. Westview Press, Boulder, Colo., 1978. xxvi, 358 pp., illus. \$23.70. Westview Special Studies on China and East Asia/South and Southeast Asia.

Molecular and Cellular Aspects of Vascular Smooth Muscle in Health and Disease. Proceedings of a symposium, Honolulu, Jan. 1977. D. F. Bohr and F. Takenaka, Eds. Karger, Basel, 1978. 216 pp., illus. Paper, \$47.50. *Blood Vessels*, vol. 15, Nos. 1-3.

The Natural History of the Whale. L. Harrison Matthews. Columbia University Press, New York, 1978. xiv, 220 pp., illus. + plates. \$19.95.

Neurochemical and Immunologic Components in Schizophrenia. Proceedings of a conference, Galveston, Tex., Oct. 1976. Daniel Bergsma, Allan L. Goldstein, Bernard Haber, Geraldine H. Cohen, Natalie W. Paul, and Jeanne Mudge, Eds. Liss, New York, 1978. xvi, 432 pp., illus. \$46. Birth Defects: Original Article Series, vol. 14, No. 5. *To order this book circle No. 584 on Readers' Service Card*

The Ocean Realm. National Geographic Society, Washington, D.C., 1978. 200 pp., illus. \$5.75.

Origin and Relationships of the California Flora. Peter H. Raven and Daniel I. Axelrod. University of California Press, Berkeley, 1978. viii, 134 pp., illus. Paper, \$6. University of California Publications in Botany, vol. 72.

Oscar Zariski: Collected Papers. Vol. 3, Topology of Curves and Surfaces, and Special Topics in the Theory of Algebraic Varieties. M. Artin and B. Mazur, Eds. MIT Press, Cambridge, Mass., 1978. xxvi, 480 pp. \$35. *Mathematicians of Our Time*.

Paper Basics. Forestry, Manufacture, Selection, Purchasing, Mathematics and Metrics, Recycling. David Saltman. Van Nostrand Reinhold, New York, 1978. xiv, 226 pp., illus. \$10.95.

Photovoltaic Power Generation. David L. Pulfrey. Van Nostrand Reinhold, New York, 1978. xvi, 218 pp., illus. \$16.95. *To order this book circle No. 577 on Readers' Service Card*

Physical Chemistry. Principles and Applications in Biological Sciences. Ignacio Tinoco, Jr., Kenneth Sauer, and James C. Wang. Prentice-Hall, Englewood Cliffs, N.J., 1978. xvi, 624 pp., illus. \$19.95.

Physiological Aspects of Crop Nutrition and Resistance. U. S. Gupta, Ed. Atma Ram, Delhi, 1977. viii, 384 pp., illus. Rs. 75.

Placodermi. Robert Denison. Gustav Fischer, New York, 1978. vi, 128 pp., illus. \$65.70. *Handbook of Paleoichthyology*, vol. 2.

Plant Tissue Culture and Its Biotechnological Application. Proceedings of a congress, Munich, Sept. 1976. W. Barz, E. Reinhard, and M. H. Zenk, Eds. Springer-Verlag, New York, 1977. xvi, 422 pp., illus. \$38.80. *Proceedings in Life Sciences. To order this book circle No. 578 on Readers' Service Card*

PL 1: A Structured Approach. Don Cassel. Reston (Prentice-Hall), Reston, Va., ed. 2, 1978. x, 262 pp. Paper, \$9.95.

Polar Research. To the Present, and the Future. Papers from a AAAS meeting, Denver, Feb. 1977. Mary A. McWhinnie, Ed. Published for the American Association for the Advancement of Science by Westview Press, Boulder, Colo., 1978. xxx, 310 pp., illus. \$20. AAAS Selected Symposia Series, 7.

Pollen Flora of Argentina. Modern Spore



Alfa's new 1979-80 Catalog NOW WITH ORGANICS.

Alfa, the recognized leader in inorganic and organometallic research chemicals and materials, has now added organic chemicals to its new, 1979-80 catalog, to better supply all of your needs from one convenient source.

In addition to its new organic line, the Alfa Catalog includes:

- Inorganics & Organometallics • Ultrapure Chemicals • Pure Elements
- Puratronic® Products • Analytical Standards and • Research & Safety Equipment...at the lowest industry prices. Plus, attractive quantity discounts available to all customers. Compare for yourself.

We hope that you will call on us, not only with your requirements for specific products, but also for applications and new product assistance or to exchange your thoughts and ideas with our staff chemists.

This is the one valuable chemicals and materials handbook which you can't afford to be without. Send today for your free Alfa Catalog.

Alfa DIVISION

Ventron Corporation

Dept.15, 152 Andover St., Danvers, MA 01923, (617) 777-1970, (415) 836-1777

Name: _____

Title: _____

Company: _____

Street: _____

City: _____

State: _____ Zip: _____

Circle No. 72 on Readers' Service Card

L/I's new REPIPET[®] II dispenser Virtually unbreakable!

Better than any \$85 plastic dispenser! But only \$52.50!*

Here it is, the best value ever offered in reagent dispensers. The tough REPIPET II reflects more than 14 years of L/I's experience in liquid dispensing—your best assurance of quality and performance. All REPIPET II dispensers feature the use of fluorocarbon and extra-thick borosilicate glass for strength and durability. They're tough lab workhorses, virtually impossible to break.

Note the 100-division scales, printed where they should be—on the glass barrel, not on an elastic plastic! L/I's exclusive sturdy magnifying indicator is a further aid to precise volume settings.

L/I guarantees accuracy of 1.5% full scale and 0.2% reproducibility for all REPIPET II dispensers. They are suitable for use with all laboratory reagents (except HF).


(Although TEFLON[®] plastic and TEFLON[®] plastic are chemically inert to all laboratory reagents, slight distortions may occur at autoclaving temperatures and with some strong solvents.)

*REPIPET II dispensers are offered in three popular sizes, 5, 10, and 20 ml. Price of 5 and 10 ml units is \$52.50; 20 ml, \$64.50.

To order, contact your distributor. For literature, write, call or circle the reader service number below.


LABINDUSTRIES

620 Hearst Avenue, Berkeley, CA 94710
Phone (415) 843-0220



Circle No. 147 on Readers' Service Card

get the most out of your gradients



ISCO density gradient fractionators produce a continuous absorbance profile as the gradient is fractionated into a built-in collector. Zones undetectable by other methods are easily resolved and their exact locations plotted. Any size centrifuge tube can be quantitatively scanned at a choice of wavelengths, and fractionated with or without piercing.

The Model 570 gradient former sequentially produces up to 22 absolutely identical linear or curved gradients of any size from 3 to 80 ml. Gradient size and shape are precisely and reproducibly determined by differential pumping.

For more information, send for your free catalog today. Or dial direct, toll free (800) 228-4250 (continental U.S.A. except Nebraska). Instrumentation Specialties Company, P.O. Box 5347, Lincoln, Nebraska 68505.



Instruments with a difference

Circle No. 119 on Readers' Service Card