

14 September 1979 • Vol. 205 • No. 4411

\$1.50

SCIENCE

AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE



A photograph of a Beckman Model TJ-6 Centrifuge. The centrifuge is a teal-colored, boxy machine with a large circular lid that is open, revealing the rotor inside. The rotor has four sample tubes: an orange one at the top, a blue one on the left, a green one on the right, and a red one at the bottom. The control panel on the front of the machine includes a large speed dial on the right, a smaller timer dial, and a power switch. The text "BECKMAN Model TJ-6 Centrifuge" is visible on the front panel. The background is a collage of various text snippets, some of which are repeated, such as "The TJ-6 is an excellent product" and "We're delighted we have one".

No Wonder the Beckman TJ-6 is the Most Popular Centrifuge Ever Introduced!

Users appreciate how quiet it is. They like the removable bowl that makes cleaning a snap. They like the colorful adapters which bring a new practicality and ease of handling to table-top centrifuges.

are hundreds of Beckman service representatives nationwide—at least one in every major city in the U.S. and Canada.

Beckman Instruments,
Inc., Spinco Division
1117 California Avenue
Palo Alto, CA 94304.

Circle No. 365 on Readers' Service Card

The Touch Meter

originated by CORNING

You have to touch it to believe it. You have to touch it to see its speed. Touch it and see how smart it is.

No analog dials. No thumb-wheels or switches. Just a smooth-surface data entry pad you control with the touch of a finger.

TOUCH here to choose your mode. A variety of Ion Selective Electrode operations. pH. mV. °C.

TOUCH here to use its memory. Instruct, store values, compute.

TOUCH here to set your reference. The Touch Meter responds at once, and keeps watch to hold these values true.

TOUCH here for instant readouts on the bright, clear LED display. Concentration reads directly in the units you have chosen.

All the features of the TALL meter. Space saving, fluid movement sidearm, up front operation. Three year warranty! Plus the speed, wit, and simplicity of micro-processor control.



It even does scientific notation. Touch. And see.

GET IN TOUCH

Get to handle a Touch Meter, up close in your own lab. See for yourself what great sense it makes.

Attach this coupon to your letterhead today for our literature. We'll contact you for a demonstration.

Science Products Division
Mail Station 18
Attention: TSK
Corning Glass Works
Corning, New York 14830

CORNING

SE-T-9

Suggested list price \$1299



CORNING
pHware

Circle No. 310 on Readers' Service Card

SCIENCE

LETTERS	Lunar Studies: <i>P. J. Coleman, Jr.</i> , et al.; Soviet Anti-Semitism: Reply by Pontryagin: <i>L. S. Pontryagin</i> ; Academic Ties with the Middle East: <i>V. Gregorian</i> ; Nutrition Program: <i>N. Kretchmer</i> and <i>D. M. Hegsted</i> ; Interferon: <i>J. Vilček</i>	1082
EDITORIAL	Education for the 21st Century	1087
ARTICLES	The Colorado Plateaus: Cultural Dynamics and Paleoenvironment: <i>R. C. Euler</i> et al.	1089
	Economic Benefits from Research: An Example from Agriculture: <i>R. E. Evenson, P. E. Waggoner, V. W. Ruttan</i>	1101
NEWS AND COMMENT	Agricultural Network Fights Unwelcome Gift	1108
	"Truth-in-Testing" Attracts Diverse Support	1110
	<i>Briefing: Scientific American</i> To Be Published in China; Surgeon General Says Get Healthy, Eat Less Meat; Pop Nutrition Books Face Legal Hurdles; Dollar-a-Gallon Gas Fuels Metric Conversion	1112
	Test Coaching Dispute Lingers	1114
RESEARCH NEWS	How Is New Ocean Crust Formed?	1115
BOOK REVIEWS	The Eighth Day of Creation, reviewed by <i>E. J. Yoxen</i> ; The SALT Experience and Arms Control and SALT II, <i>G. H. Quester</i> ; The Outer Shores, <i>M. T. Ghiselin</i> ; The Nutrition of Aborigines in Relation to the Ecosystem of Central Australia, <i>N. Howell</i> ; Fundamental Problems in Statistical Mechanics IV, <i>R. F. Fox</i> ; Books Received.	1119
REPORTS	Nitrous Oxide Emissions from an Irrigated Cornfield: <i>G. L. Hutchinson</i> and <i>A. R. Mosier</i>	1125

BOARD OF DIRECTORS					
EDWARD E. DAVID, JR. Retiring President, Chairman	KENNETH E. BOULDING President	FREDERICK MOSTELLER President-Elect	ELOISE E. CLARK MARTIN M. CUMMINGS	RENÉE C. FOX ANNA J. HARRISON	
CHAIRMEN AND SECRETARIES OF AAAS SECTIONS					
MATHEMATICS (A) Garrett Birkhoff Ronald Graham	PHYSICS (B) Arthur L. Schawlow Rolf M. Sinclair	CHEMISTRY (C) Fred Basolo William L. Jolly	ASTRONOMY (D) Peter S. Conti Donat G. Wentzel		
PSYCHOLOGY (J) Frances K. Graham Meredith P. Crawford	SOCIAL AND ECONOMIC SCIENCES (K) David L. Sills Gillian Lindt	HISTORY AND PHILOSOPHY OF SCIENCE (L) Melvin Kranzberg Diana L. Hall		ENGINEERING (M) Daniel C. Drucker Donald E. Marlowe	
EDUCATION (Q) Fletcher G. Watson James T. Robinson	DENTISTRY (R) Carl J. Witkop, Jr. Harold M. Fullmer	PHARMACEUTICAL SCIENCES (S) Samuel Elkin Robert A. Wiley	INFORMATION, COMPUTING, AND COMMUNICATION (T) Mary E. Corning Madeline M. Henderson		
DIVISIONS					
ALASKA DIVISION		PACIFIC DIVISION		SOUTHWESTERN AND ROCKY MOUNTAIN DIVISION	
Daniel B. Hawkins President	T. Neil Davis Executive Secretary	William L. Sims President	Alan E. Leviton Secretary-Treasurer	Harry J. Shine President	M. Michelle Balcomb Executive Officer
SCIENCE is published weekly on Friday, except the last week in December, by the American Association for the Advancement of Science, 1515 Massachusetts Avenue, NW, Washington, D.C. 20005. Second-class postage (publication No. 484460) paid at Washington, D.C., and at an additional entry. Now combined with The Scientific Monthly®. Copyright © 1979 by the American Association for the Advancement of Science. Domestic individual membership and subscription (51 issues): \$34. Domestic institutional subscription (51 issues): \$70. Foreign postage extra: Canada \$12, other (surface mail) \$15, air-surface via Amsterdam \$40. First class, airmail, school-year, and student rates on request. Single copies \$1.50 (\$2 by mail); back issues \$2.50 (\$3 by mail); classroom rates on request. Change of address: allow 6 weeks, giving old and new addresses and seven-digit account number. Postmaster: Send Form 3579 to Science, 1515 Massachusetts Avenue, NW, Washington, D.C. 20005. Science is indexed in the Reader's Guide to Periodical Literature and in several specialized indexes.					

AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE


Carbon-13 in Tree Rings Indicates No Record of Sea-Surface Temperature: <i>F. MacIntyre</i>	1127
Giant Larvacean Houses: Observations from Deep Submersibles: <i>E. G. Barham</i> . .	1129
The Viscous Sublayer at the Sea Floor: <i>D. R. Caldwell</i> and <i>T. M. Chriss</i>	1131
Gravitational Lens of the Sun: Its Potential for Observations and Communications over Interstellar Distances: <i>V. R. Eshleman</i>	1133
Enhanced 5-Fluorouracil Nucleotide Formation After Methotrexate Administration: Explanation for Drug Synergism: <i>E. Cadman, R. Heimer, L. Davis</i>	1135
Iodine-125-Labeled Estradiol: A Gamma-Emitting Analog of Estradiol That Binds to the Estrogen Receptor: <i>R. B. Hochberg</i>	1138
Molecular Cloning of Polyoma Virus DNA in <i>Escherichia coli</i> : Oncogenicity Testing in Hamsters: <i>M. A. Israel</i> et al.	1140
Sporozoite-Induced Malaria: Therapeutic Effects of Glycolipids in Liposomes: <i>C. R. Alving</i> et al.	1142
Determination of the Pore Size of Cell Walls of Living Plant Cells: <i>N. Carpita</i> et al. .	1144
[³ H] GABA Binding in Brains from Huntington's Chorea Patients: Altered Regulation by Phospholipids?: <i>K. G. Lloyd</i> and <i>L. Davidson</i>	1147
Vibrations: Their Signal Function for a Spider Kleptoparasite: <i>F. Vollrath</i>	1149
Pharmacologic Effects in Man of a Potent, Long-Acting Dopamine Receptor Agonist: <i>L. Lemberger</i> and <i>R. E. Crabtree</i>	1151
Genetic Variance in Nonverbal Intelligence: Data from the Kinships of Identical Twins: <i>R. J. Rose</i> et al.	1153
The Cheetah: Native American: <i>D. B. Adams</i>	1155
Regrowth of Severed Axons in the Neonatal Central Nervous System: Establishment of Normal Connections: <i>K. Kalil</i> and <i>T. Reh</i>	1158
High-Speed Cinematographic Evidence for Ultrafast Feeding in Antennariid Anglerfishes: <i>D. B. Grobecker</i> and <i>T. W. Pietsch</i>	1161
Plasma β -Endorphin Immunoreactivity in Schizophrenia: <i>M. Ross, P. A. Berger,</i> <i>A. Goldstein</i>	1163
PRODUCTS AND MATERIALS	
Growth Factor; Water Purification; Sampler for Respiratory and Metabolic Studies; Spectrophotometer; Total Exhaust Hood; Satellite- Controlled Clock; Fraction Collector; Literature	1166

MIKE MC CORMACK RUSSELL W. PETERSON	JOHN C. SAWHILL CHEN NING YANG	WILLIAM T. GOLDEN Treasurer	WILLIAM D. CAREY Executive Officer
GEOLOGY AND GEOGRAPHY (E) Linn Hoover Ramon E. Bisque	BIOLOGICAL SCIENCES (G) Donald S. Farner Walter Chavin	ANTHROPOLOGY (H) James B. Watson Priscilla Reining	
MEDICAL SCIENCES (N) Theodore Cooper Leah M. Lowenstein	AGRICULTURE (O) J. Lawrence Apple Coyt T. Wilson	INDUSTRIAL SCIENCE (P) Herbert I. Fustfeld Robert L. Stern	
STATISTICS (U) Richard L. Anderson Ezra Glaser	ATMOSPHERIC AND HYDROSPHERIC SCIENCES (W) Eugene W. Bierly Glenn R. Hilst	GENERAL (X) Ruth B. Pitt S. Fred Singer	

COVER

North American Pleistocene cheetah *Acinonyx trumani*. Reconstruction based on the living cheetah, *Acinonyx jubatus*, and skeletons of *A. trumani* from Natural Trap Cave, Wyoming. See page 1155. [Drawing by Dawn A. Adams]

The American Association for the Advancement of Science was founded in 1848 and incorporated in 1874. Its objects are to further the work of scientists, to facilitate cooperation among them, to foster scientific freedom and responsibility, to improve the effectiveness of science in the promotion of human welfare, and to increase public understanding and appreciation of the importance and promise of the methods of science in human progress.



We pay these scientists to make our production lines obsolete.

As soon as there's a better way to do something, the old way becomes obsolete.

And these scientists do nothing *but* design new and better ways for Western Electric to manufacture products for the Bell System. In the last two years, they and their colleagues at Western Electric's Engineering Research Center have helped cut costs by over \$156 million.

They are a vital part of Western Electric's commitment to helping the Bell System keep the cost of communications service down, while keeping the quality up.

The Factory's Laboratory

The people at our Engineering Research Center work full time on manufacturing research. Over 50% of them are PhD's.

At home in the factory as well as the lab, they turn theoretical

experiments into real savings. For example, while exploring ways to improve production of telephone circuit boards, they invented a totally new material with which to make the boards and a whole new way to mass-produce them.

Another breakthrough is already saving more than \$2 million a year. It's a new soldering technique that makes 20,000 connections in less than 60 seconds.

The list of money-saving developments goes on. A new computer program that makes technical drawings in one-fourth the time of manual methods. An automatic system for sorting telephone parts by color. A process that measures the production of optical fiber 1000 times per second, keeping its diameter accurate within 30-millionths of an inch. And there's a new system that takes only 30 seconds to help align tiny integrated circuits that are

200 times thinner than the thinnest human hair.

These are just a few examples of what the Engineering Research Center does throughout *all* of Western Electric's manufacturing operations — developing new ideas in technology to improve telephone products and produce them at lower cost.

Anticipating the Future


Much of the work of the Engineering Research Center is designed to keep Western Electric ahead of technology, by anticipating what new processes will be needed for tomorrow's products.

It's a total and ongoing commitment. A commitment to provide you with the highest quality telephone service at the lowest possible cost.

Keeping your communications system the best in the world.



Western Electric



**If you can't believe the price tag,
pretend they're \$500 more.**

The no-nonsense balances from Fisher.

In Spring 1976 we introduced our first **it-can't-be-a-balance** balance.


Today there's a whole family of Fisher/Ainsworth electronic toploaders, including a compact printing balance and the lab world's first battery-pack portable. Plus our two newsmaking microprocessor-

controlled **M** models. Each a remarkably rugged lightweight instrument. With single touchbar operation (instant tare, instant sample readout on LED display). And so **modestly priced** they've made a whole generation of lab people re-think their concept of what a balance can be!

Fisher/Ainsworth toploaders. Solid examples of American electronic technology, innovative manufacturing, and sound value.

Say **No!** to inflation. Say **Yes!** to savings. Up to 50% per balance. Send for product data today.


Model	A-50	A-200	A-200BP battery-operated	A-200BCD printing balance	M-300	M-300DR dual-range	A-2000
Capacity	50g	200g	200g	200g	300g	30/300g	2000g
Readability	0.001g	0.01g	0.01g	0.01g	0.01g	0.001/.01g	0.1g
Catalog no.	1-913-265	1-913-260	1-913-261	1-913-262	1-913-252	1-913-253	1-913-267
U.S. price	\$1395	\$895	\$1195	\$1695	\$995	\$1495	\$1195



We Know Laboratories.



Fisher Scientific Company
711 Forbes Avenue
Pittsburgh, Pennsylvania 15219
(412) 562-8543



Why are the world's most advanced LS instruments tested with these nuclear supplies?

The Beckman total nuclear supply capability evolved along with our LS counting instruments to give us a constant source of

fully quality-controlled cocktails, solubilizers and counting vials.

Today, the same supplies we used in our own state-of-the-art instrument R&D programs are available to LS researchers nationwide.

And whether or not you own a Beckman Instrument you'll benefit from optimized counting efficiencies, superior reliability, and the convenience of one expert source for all your nuclear needs.

Right now, orders of premixed cocktails, cocktail ingredients, and counting vials are available with 24-hour response times from a national network of


shipping depots. And major cost reduction programs have been designed to meet the individual requirements of research facilities both small and large.

So contact your local Beckman representative now and ask about our special nuclear supplies introductory programs.

Or contact Scientific Instruments Division, Beckman Instruments, Inc., P.O. Box C-19600, Irvine, CA 92713.


Innovation you can count on.

BECKMAN®



Guinness Book of World Records, please note:
**This is probably the largest selection of
plastic laboratory animal cages in the world.**

29 cages that accommodate mice.
22 cages that accommodate rats.
18 cages that accommodate hamsters.
12 cages that accommodate guinea pigs.




Plus: Many types of covers, filter caps, bedding, racks, accessories, watering systems, metabolism cages, etc., etc., etc.


Plus: Our housing systems and accessories adhere to all NAS, NIH, and Public Law 89-544 requirements.

Plus: We manufacture all our cages and covers and virtually everything else listed above. We control the quality of everything we sell.

Please write or call for details.
Lab Products Inc., 365 West Passaic
Street, Rochelle Park, New Jersey 07662
Phone: 201 / 843-4600

lab products inc
a **Medic** company





From the beginnings of biomedical research, instruments have played an important role in the drama of discovery. Crude though they were, early microscopes served to extend the vision of the determined researchers who, in turn, extended our understanding of life itself.

Today, researchers are still characterized by the rare insight and unusual skill of their predecessors. But their vision of the microscopic world is considerably sharper; thanks in part to Nikon's tireless dedication to the creation of superior optics.

Now, with a single stunning innovation, Nikon has produced an unprecedented advance in microscope technology. Called CF optics, it's a development so revolutionary that it represents far more than a refinement. It stands as nothing less than an entirely new standard of light microscope performance. A standard against which all others will be measured.

To achieve this dramatic breakthrough, Nikon engineers combined a new, extra-low dispersion glass with an ingenious concept in optical design to virtually eliminate chromatic aberration. The Nikon CF system delivers unparalleled resolution, color fidelity, flatness of field and image brightness.

In fact, CF optics bring to Nikon resolving power closely approaching theoretical limits.

Biomedical research is a matter of life itself. It demands Nikon. The difference between seeing and not seeing. Knowing and not knowing.

For further information on the entire line of Nikon instruments for biomedical research, contact: Nikon Incorporated, Instrument Division, Ehrenreich Photo-Optical Industries Inc. 623 Stewart Avenue, Garden City, NY 11530 (516) 222-0200.

Look to Nikon


ニコンとご用命下さい

Blicken Sie auf Nikon


Circle No. 290 on Readers' Service Card

IN BIOMEDICAL RESEARCH, LOOK TO NIKON.

I BIOPHOT RESEARCH PHOTOMICROSCOPE WITH NOMARSKI DIC — Everything an advanced research microscope should be, plus CF optics for unsurpassed resolution and contrast. Complete capabilities for Nomarski differential interference contrast, phase contrast, and other techniques. Complete photomicrographic capability.




II OPTIPHOT RESEARCH MICROSCOPE WITH HFM — Versatile, rugged microscope features CF optics for superb field flatness, contrast and resolution. Shown here with trinocular bodytube and the advanced, computer controlled HFM camera system for consistently correct exposures — automatically.




III SMZ-10 ZOOM STEREO MICROSCOPE — High-resolution zoom stereo microscope with the ability to photograph actual stereo pairs by the simple touch of a lever. Features a continuous zoom over a wide magnification range without the need to change focus or working distance.




IV LABOPHOT LABORATORY MICROSCOPE — Developed with special attention to the sturdiness and simplicity of operation required in laboratory and clinical work. Employs CF optics for unsurpassed optical performance. The microscope accepts the full line of Nikon accessories, including photomicrography.




V FLUOPHOT RESEARCH FLUORESCENCE MICROSCOPE — Provides transmitted-light or reflected-light fluorescence with convenient one-step selector knob. Built-in turret with four reflected-light fluorescence excitation modes. Brighter-than-ever images with 200-watt DC mercury illumination and high-transmission trinocular head.




Research for impatient



a. CK Isoenzymes, separated by electrophoresis on agarose. Type 665 Positive/Negative film.



b. Micrograph of longitudinal section of rat incisors, 450x. Type 667 Coaterless film.



c. Green peach aphids (*myzus persicae*), 8x. Type 57 High Speed 4 x 5 film.



d. Copy of chart. Type 51 High Contrast 4 x 5 film.

e. Hook and dentition, fresh water clam. 520x, scanning electron microscope. Type 52 Fine Grain 4 x 5 film.

photography people.

If you need photographs of your work, why sit around waiting for them? With Polaroid instant photographic equipment you get immediate results.

Our wide range of equipment lets you record everything from microbes to metal stress tests. With Polaroid self-developing films you can have professional quality results in color in 1 minute, or in black and white (with or without a usable negative) in seconds. And our equipment is simple to use. So you can take the photographs you need without needing to know a lot about photography.

The MP-4 Multipurpose camera (1) is a versatile, self-contained photo studio anyone can operate. It uses 14 different Polaroid instant films to keep you out of the darkroom. And it copies, delivers close-ups, reductions, macrophotographs and photomicrographs, to bring your answers to light.

Our CU-5 Close-up camera (2) is a lightweight, hand-held system you can take almost anywhere and get instant photos from ¼ to 3 times life size. Exposure is easy to set,

lighting is built in and framing is automatic. So all you have to think about is the picture. You can even use the CU-5 to capture a transient image on a cathode ray tube.



Many cameras and instruments can be adapted quickly and easily for instant photography with Polaroid film holders (3). They come in 3 models to handle 3 different sizes of Polaroid Land film (3¼ x 4¼, 4 x 5, and 8 x 10 in.), so you can get instant results in almost any format.

Many other manufacturers of cameras and instruments also supply Polaroid Land camera backs that adapt their equipment to instant photography. These backs use 9 different convenient Polaroid pack films, so you can see your project in a new light.

To find out how you can get instant results, mail the coupon below. Or if you're a very impatient person, call us toll-free from the continental U.S. at 800-225-1618 (from Massachusetts, call collect 617-547-5177).

And stop tapping your foot while you wait for your photos to come back from the lab.

Circle No. 347 on Readers' Service Card



Polaroid

Instant Laboratory Pictures

Polaroid Corporation
Industrial Marketing, Dept. A434
575 Technology Square
Cambridge, Mass. 02139

I'm impatient. Please send me more information on how I can use Polaroid instant photography in my work.

Name _____
Title _____
Company _____
Address _____
City _____ State _____ Zip _____
Current camera or instrument _____
Application or need _____

pH Electrode Quiz

QUESTION

The most common cause of pH electrode failure is:

- ☐ Other people break them.
- ☐ Other people let the reference dry out.
- ☐ Cosmic forces de-energize them.

ANSWER

Eventually all pH electrodes age; response becomes slow and span is short. Unfortunately, most electrodes are broken or let dry out before old age takes its toll.

PRIZE

Sensorex Combination pH Electrodes are prize-winners. They have:

- Epoxy Bodies with Recessed Bulb/Safeguard Tips to minimize breakage.
- Sealed, gel-filled references that never need refilling.
- Fast response over the full pH range.

**STANDARD
SIZE
MODEL
S200C**

\$36

Semi-micro, flask, flat surface, other models available. Send for Bulletin 200A




9713 Bolsa Ave.
Westminster, CA
USA 92683
Phone: 714/554-7090

Circle No. 35 on Readers' Service Card

HIGHLY PURIFIED

PHORBOL ESTERS

AND RELATED DITERPENES



Partial Listing:

		1g	300mg	100mg	25mg	5mg
PMA	Phorbol Myristate Acetate	\$2450	\$ 895	\$325	\$ 95	\$28
MPMA	4-O-Methylphorbol Myristate Acetate		2400	895	275	65
PDBu	Phorbol Dibutyrate	3400	1300	530	150	38
PDBuAL	Phorbol Dibutyrate Aldehyde	4500	1800	690	205	55
PDA	Phorbol Diacetate	2950	1080	390	110	28
4α-PDD	4α-Phorbol Didecanoate		2400	895	275	65
PHR	Phorbol	1450	580	235	70	19
MEZ	Mezerein		1800	690	205	55

- Complete Selection of Related Compounds
- Custom Synthesis (Labeled or Unlabeled)
- Bulk Quantities

For Complete Catalog
Call or Write:



Lifesystems Company

1254 Chestnut Street
Newton, Massachusetts 02115
617-244-8800

Circle No. 367 on Readers' Service Card

GUTH HANDI-MAT® Belongs In Your Lab

- Clear durable, heavy weight polyethylene
- Easily cut to shape with shears or knife
- Protects benchtops and floors
- Cushions equipment
- Excellent splash protection
- Lines drawers and shelves



Thickness
through ribs:
.028" approx.

- Available From Your Leading Laboratory Supplier!

GUTH

GUTH PRODUCTS COMPANY
332 S. Center St., Hillside, IL 60162 USA
(312)547-9050 Cable Address: "GUTHCO"

Circle No. 330 on Readers' Service Card

COMBATING THE #1 KILLER

The **SCIENCE**
Report on
Heart
Research



**JEAN L. MARX and
GINA BARI KOLATA**

\$17.00 casebound

\$ 7.50 paperbound

10% discount to AAAS members

Send name, address and remittance to

**AMERICAN ASSOCIATION
FOR THE ADVANCEMENT
OF SCIENCE**

Department B-6
1516 Massachusetts Avenue, NW
Washington, D.C. 20006

Ohaus Brainweigh 300

Never before, for so little, could you get so much.

The Ohaus Brainweigh 300 costs only \$1,195 and gives you a capacity of 300g with a readability of 0.01 g. (If you want to find out how much that price is going to save you, just check what our competitors are charging for balances with similar capabilities.)

And we didn't stop with just economy. We went on to design in features that would be costly "extras" on any other balance on the market—if you could even get them.

Like extended visibility



Stand up; sit down; go to the back of the room. The display on the Ohaus 300 is still easy to read. That's

because the no-glare window with the big, easy-to-read digits is deliberately angled for maximum visibility wherever you are. And that can save a lot of headaches whether you're in a lab, or on the line or in a classroom.


It tells you when it's ready—and when you've goofed



No more guessing whether the balance has settled down; no more annoying digit flutter. When the "g" appears beside the digits—in 3 seconds or less—it's your guarantee that the reading is stable.

And if the pan is overloaded or dislodged, you'll get an error sign in the window.

The 300 is a new generation balance



This new, low-profile balance is only a little larger than an ordinary desk phone. But packed inside is enough advanced technology to make the Ohaus 300 a generation ahead of any other balance on the market. Its state-of-the-art electronic components are almost entirely plug-in for easy service. And there are fewer of them—for better reliability and more dependable performance.


Want a BCD interface now—or later?



You don't have to worry about it. Unlike other electronic bal-



ances, every Ohaus Brainweigh comes equipped to handle a BCD option. You can add a BCD interface right away—or whenever you come to need it.



Parts and service are as near as your dealer

The Ohaus 300 is made right here in the U.S. And although its built-to-take-it die-cast case and solid state components put maintenance at a minimum, when the balance finally needs a little attention, it's nice to know your local dealer can do the job.



Get all the facts

If you'd like to see the rest of the reasons why the Ohaus Brainweigh 300 outclasses the competition, and the many ways it can make your operation faster, easier and more accurate, just return the coupon.

**THE OHAUS
BRAINWEIGHS™**
The weight of the future.

©1979 Ohaus Scale Corporation

OHAUS

© Ohaus Scale Corporation, Dept. 14-099
29 Hanover Road, Florham Park, N.J. 07932
(201) 377-9000

- ☐ Yes, show me the weight of the future. Send me the full-color catalog on the Ohaus Brainweighs.
- ☐ I'm ready to see a demonstration. Please have a dealer sales rep call me.

Name _____

Title _____

Organization _____

Address _____


City _____

State _____

Zip _____

Phone Number _____

Circle No. 311 on Readers' Service Card



Want the only dispenser with a Teflon[®]-coated plunger?

Ask for the 'Teflon Dispenser'

Why does a Brinkmann Dispensette seem to work more smoothly than other bottle-top dispensers?

The secret is its Teflon-coated plunger.

Dispensers with glass or polypropylene plastic plungers often 'freeze' or 'stick' when used with alkaline reagents. On a Dispensette, the Teflon coating on the plunger insures smooth operation with any reactive chemical (except HF), even with concentrated acids and bases.

Of course, it takes more than a smoothly-working plunger to make a dependable dispenser. Instead of external glass tubing that could break or chip, Dispensettes are equipped with flexible Teflon filling and discharge tubes that can easily be cut to any desired length. Simply pull them off and a Dispensette is ready for autoclaving at 120°C without further disassembly.

Dispensette offers a wide selection of models for fast, accurate dispensing of volumes from 0.1 to 50ml, with better than $\pm 1.0\%$ accuracy and $\pm 0.1\%$ reproducibility. There is hardly a standard size screw-neck reagent bottle, can or container (even STJ 24/40 and 29/42 glassware) which a Dispensette will not fit, either directly or with optional screw-in adapters.


Some Dispensette design features have been copied by other dispensers, but none has its smooth, foolproof operation. That's because only Dispensette has a Teflon-coated plunger. Write for literature: Brinkmann Instruments, Inc., Cantiague Rd., Westbury, N.Y. 11590. In Canada: Brinkmann Instruments (Canada), Ltd.

Brinkmann Dispensette[®]


Available from: Ace Scientific
Beckman Instruments, Science Essentials Operations
Bio-Rad Laboratories/Cole-Parmer Instrument Co.
Curtin Matheson Scientific/Fisher Scientific
General Scientific/Markson Scientific/Preiser Scientific
Sargent Welch/Scientific Products/SGA Scientific
Arthur H. Thomas/VWR Scientific/Wilkins-Anderson Co.

Dispensette[®] is a registered trademark of R. Brand Co.
Teflon[®] is a Du Pont trademark.

Circle No. 250 on Readers' Service Card



The Illumitran eliminates Original Sins.



The Bogen Illumitran 3C, with its exclusive contrast control, can create copies that exactly reproduce or improve your original slide. And it can do it using the widest choice of films.

How? With two unusual features: Contrast variation, and induced shadow variation (that lets



you vary the color of shadow areas without affecting highlight reproduction).

Already have an Illumitran, but no contrast control? No problem. Yours can be easily updated.

For all the facts on the Bogen Illumitran please complete the coupon at right.

Bogen Photo Corp., 100 So. Van Brunt St., P.O. Box 448, Englewood, NJ 07631.

Please send me complete information on Bogen's Illumitran.

Name

Street

City State Zip 9/79 S

Lester Bogen doesn't sell anything he wouldn't buy himself.

bogen

MOP+YOU= instant image analysis.

For a multitude of geometric dimensions

The new MOP quantitative system for morphological and stereological image analysis is unique: Only MOP lets the operator's discriminating judgment interact with the high-speed calculation ability and memory of a microcomputer. This mode of operation is highly reliable in applications where fully automatic systems fail.

Single key selects a multitude of functions


MOP is portable, requires no installation and practically no operator training. Yet it can measure 8 parameters simultaneously, store data in 20 channels, and perform statistical analyses such as distribution histograms. Easily interfaced to computers for data-banking and high-level statistics. **Nationwide service.**

The great name in optics



Circle No. 213 on Readers' Service Card

Carl Zeiss, Inc., 444 5th Avenue, New York, N.Y. 10018 (212) 730-4400. Branches: Atlanta, Boston, Chicago, Houston, Los Angeles, San Francisco, Washington, D.C. In Canada: 45 Valleybrook Drive, Don Mills, Ontario, M3B 2S6. Or call (416) 449-4660.



The MOP measuring tablet can be used with images directly generated through a microscope, by projection, or with photomicrographs, electron micrographs, X-Rays, drawings, etc.

MORE THAN EVER, CA SELECTS Will Keep You Informed

We've added more new topics to the CA SELECTS series. Now there are 80 reasons to look into CA SELECTS, the current-awareness service from the publishers of CHEMICAL ABSTRACTS.

CA SELECTS brings you...

- Informative abstracts of current chemistry-related literature
- Your choice of subject coverage, from the topics listed below
- Delivery every two weeks

All this at a price you can afford.

See For Yourself,
Get a FREE Issue
of CA SELECTS

To receive a FREE issue of CA SELECTS, indicate the topic(s) pertaining to your interest from the list below. We'll also send complete ordering instructions but you won't be obligated to order a subscription.

CA SELECTS

Clip This Form (or Mail a Photocopy) Today

CA SELECTS

Now — 80 Topics

Please send me my free issue(s) of CA SELECTS, covering the topic(s) I have indicated below. Also send complete ordering instructions.

I understand this request does not obligate me to order a subscription.

- ☐ Adhesives
- ☐ Analytical Electrochemistry
- ☐ Animal Longevity & Aging
- ☐ Anti-Inflammatory Agents & Arthritis
- ☐ Antitumor Agents
- ☐ Atherosclerosis & Heart Disease
- ☐ Batteries & Fuel Cells
- ☐ Biogenic Amines & the Nervous System
- ☐ Biological Information Transfer
- ☐ Blood Coagulation
- ☐ Carbon & Heteroatom NMR ★
- ☐ Carcinogens, Mutagens, & Teratogens
- ☐ Catalysis (Applied & Physical Aspects)
- ☐ Catalysis (Organic Reactions)
- ☐ Chemical Hazards, Health, & Safety
- ☐ Coal Science & Process Chemistry
- ☐ Colloids (Applied Aspects)
- ☐ Colloids (Macromolecular Aspects) ★
- ☐ Colloids (Physicochemical Aspects)
- ☐ Computers in Chemistry
- ☐ Corrosion
- ☐ Crystal Growth
- ☐ Detergents, Soaps & Surfactants
- ☐ Drug & Cosmetic Toxicity
- ☐ Electrochemical Reactions
- ☐ Electrodeposition
- ☐ Electron & Auger Spectroscopy
- ☐ Electron Spin Resonance (Chemical Aspects)
- ☐ Energy Reviews & Books
- ☐ Engine Exhaust
- ☐ Environmental Pollution
- ☐ Flammability
- ☐ Flavors & Fragrances
- ☐ Food Toxicity
- ☐ Forensic Chemistry
- ☐ Fungicides
- ☐ Gas Chromatography
- ☐ Gaseous Waste Treatment
- ☐ Gel Permeation Chromatography
- ☐ Heat-Resistant & Ablative Polymers

- ☐ Herbicides
- ☐ High Speed Liquid Chromatography
- ☐ Infrared Spectroscopy (Organic Aspects) ★
- ☐ Infrared Spectroscopy (Physicochemical Aspects) ★
- ☐ Insecticides
- ☐ Ion Exchange
- ☐ β -Lactam Antibiotics
- ☐ Liquid Crystals
- ☐ Liquid Waste Treatment
- ☐ Mass Spectrometry
- ☐ Metallo Enzymes & Metallo Coenzymes
- ☐ Moessbauer Spectroscopy
- ☐ New Books in Chemistry
- ☐ Organic Reaction Mechanisms
- ☐ Organofluorine Chemistry
- ☐ Organophosphorus Chemistry
- ☐ Organosilicon Chemistry
- ☐ Organo-Transition Metal Complexes
- ☐ Paper & Thin-Layer Chromatography
- ☐ Photobiology
- ☐ Photochemistry
- ☐ Pollution Monitoring
- ☐ Prostaglandins
- ☐ Proton Magnetic Resonance ★

- ☐ Psychobiology
- ☐ Radiation Chemistry
- ☐ Raman Spectroscopy
- ☐ Recovery & Recycling of Wastes
- ☐ Silver Chemistry
- ☐ Solar Energy
- ☐ Solid & Radioactive Waste Treatment
- ☐ Solvent Extraction
- ☐ Steroids (Biochemical Aspects)
- ☐ Steroids (Chemical Aspects)
- ☐ Substituent Effects & Linear Free Energy Relationships
- ☐ Surface Chemistry (Physicochemical Aspects)
- ☐ Thermochemistry
- ☐ Trace Element Analysis
- ☐ X-Ray Analysis & Spectroscopy
- ☐ Zeolites ★

★ New as of July 1979

Mail Coupon To:

Chemical Abstracts Service
Marketing Department — SMA
P. O. Box 3012
Columbus, Ohio 43210 U.S.A.

Name _____

Title _____

Organization _____

Address _____



AAAS BOOKS AAAS BOOKS

From heart disease to solar energy, from food to electronics—the American Association for the Advancement of Science publishes books which cover these topics and many more. Books that will keep you informed and abreast of important scientific advancements. Books that are designed for the researcher, scientist, engineer, student, and general public. Books that you need to read.

AAAS SCIENCE REPORT SERIES

Combating the #1 Killer: The SCIENCE Report on Heart Research, by Jean L. Marx and Gina Bari Kolata. Based on the American Heart Association's Blakeslee Award-winning articles in SCIENCE, this book presents an overview of the status of current research on all aspects of cardiovascular diseases. xi + 205 pp. \$17.00 (case); \$7.50 (paper).

Solar Energy in America, by William D. Metz and Allen L. Hammond. Chapters deal with research strategies of both government and industry: "Power Towers," photovoltaic cells and other collector strategies, fuel photosynthesis, wind and ocean thermal power, energy storage, and more. Fall 1978. ca. \$17.00 (case); ca. \$7.50 (paper).

Two other AAAS series, *Compendia* and *R&D Reports*, provide authoritative information on some of the principal issues that concern us all today. These publications can serve not only as valuable reference sources for scientists and researchers, but also as supplementary reading material in courses devoted to these issues.

AAAS SCIENCE COMPENDIA SERIES

ENERGY: Use, Conservation and Supply; Volume II. vi + 202 pp. \$14.00 (case); \$6.00 (paper).

ELECTRONICS: The Continuing Revolution. iv + 217 pp. \$12.00 (case); \$4.50 (paper).

MATERIALS: Renewable and Nonrenewable Resources. x + 198 pp. \$12.00 (case); \$4.50 (paper).

POPULATION: Dynamics, Ethics and Policy. viii + 184 pp. \$12.00 (case); \$4.50 (paper).

FOOD: Politics, Economics, Nutrition and Research. vi + 202 pp. \$12.00 (case); \$4.50 (paper).

ENERGY: Use, Conservation and Supply; Volume I. vi + 154 pp. \$10.00 (case); \$3.50 (paper).

HEALTH CARE: Regulation, Economics, Ethics, Practice. iv + 256 pp. Fall 1978. \$14.00 (case); \$6.00 (paper).

AAAS RESEARCH & DEVELOPMENT REPORT SERIES

Research & Development: AAAS Report III. xiv + 122 pp. \$6.00 (paper).

Research & Development in the Federal Budget: FY 1978. xii + 148 pp. \$5.00 (paper).

Research & Development in the Federal Budget: FY 1977. ix + 100 pp. \$5.00 (paper).

Now . . . it's your turn to let us know how we can help you. . . .


To order your individual copies of the above titles, please send your name, address, and list of book titles (indicate case or paper) to AAAS, Department J-1, 1515 Massachusetts Avenue, NW, Washington, DC 20005. Remittance must accompany all orders under \$10.00. Please allow 6-8 weeks for delivery.

^3H Cum Laude



Each year we make it tougher for our graduates. For three years our most promising tritiated, multi-labeled candidates have had to pass TNMR—T for tritium, NMR for Nuclear Magnetic Resonance. It's a very precise analytical technique that ends up assuring you of more confidence in tracing metabolic and enzymatic processes with our ^3H labeled compounds.

Now, to characterize our ^3H chemicals even more rigorously, we have purchased a super-conducting NMR spectrometer. It was quite an investment, but we consider it a necessity to serve the ever more sophisticated requirements of biochemical science.

 **New England Nuclear**
Label of Excellence

Circle No. 315 on Readers' Service Card

Now you can weigh with analytical precision (0.1mg) on a fully automatic electronic top-loader, or on a fully automatic electronic suspended-pan balance, as the weighing situation or your personal preference dictates.

The new top-loader that weighs with analytical precision is our Model 1201MP. Compact and very affordable, this addition to our popular 1200 Series has a capacity of 30g and weighs to 0.1mg at the touch of a button. Its many advanced

features include a built-in micro-processor, instant electronic taring, memory and weight recall, BCD output for printer connection, and a large, bright digital display with reading stability indicator. Pan access is provided from the top and both sides of the glass enclosed weighing chamber.

For 0.1mg precision in an electronic analytical balance, we offer our new Model 2003MP. Fully automatic push-button operation and other features are similar to our

electronic top-loaders. In addition, a door-activated switch blocks display of the last digit except when both chamber doors are closed, thereby preventing air current-induced instability of the readout. With a capacity of 166g, larger pan and larger weighing chamber, the 2003MP is more versatile and somewhat more costly.

For literature write: Sartorius Balances, Division of Brinkmann Instruments, Inc., Cantiague Road, Westbury, N.Y. 11590.

Sartorius introduces an electronic top-loader with the same precision (0.1mg) as our new electronic analytical balance.



FULLY AUTOMATIC
ELECTRONIC
TOP-LOADING
BALANCE WITH
0.1mg PRECISION,
MODEL 1201MP



FULLY AUTOMATIC
ELECTRONIC
ANALYTICAL
BALANCE WITH
0.1mg PRECISION,
MODEL 2003MP



sartorius

Circle No. 251 on Readers' Service Card

WHY DO SO MANY PEOPLE BUY THE KIND OF PROTECTION THAT'S WRONG FOR THEM?



Because they don't clearly understand the difference between Term and Whole Life insurance ... and that a Term policy can typically provide five to ten times as much immediate family protection for the same money!

As shown in the table below, a first year net outlay of only \$150 buys a \$142,000 TIAA 20-Year Decreasing Term policy or a \$97,000 5-Year Renewable Term policy, for a man aged 30 or a woman aged 35. The same outlay used to purchase a Whole

Life policy would give their family only \$13,000 of family protection.

If you're like most of your colleagues in education and research, your first priority is to give your family the most financial security for the least amount of money—and that's what TIAA Term is all about! Whatever *your* age, the difference in the level of protection available for just \$150 is dramatic, as the table illustrates.

Immediate Protection Available For \$150 First Year Net Cost *

For a Man Aged	For a Woman Aged	20-Year Decreasing Term Plan	5-Year Renewable Term Plan	Whole Life Plan
25	30	\$206,000	\$112,000	\$16,000
30	35	142,000	97,000	13,000
35	40	88,000	74,000	10,000
40	45	52,000	48,000	8,000
45	50	31,000	31,000	7,000

*Annual premium less cash dividend payable at the end of first policy year, based on 1979 dividend scales. While not guaranteed, dividends have been paid every year since TIAA's founding in 1918.

Why such big differences? With Term you pay only for protection, while a Whole Life policy combines protection with a savings (cash value) element. Remember, the protection part of life insurance is relatively inexpensive; it's when savings are included as well that a high outlay per \$1,000 of coverage is required.

If you'd like to know more, ask us to send two highly informative articles reprinted from *The New York Times* and *Business Week* that discuss the choice between Term and Whole Life. We'll gladly mail them, along with detailed personal information about TIAA policies for you.

Mail coupon or call collect

Just complete and mail the coupon at right, or call the TIAA Life Insurance Advisory Center (collect) 212-490-9000. If you wish, an Insurance Counselor will review your insurance needs with you and help you select the plan and amount that's right for you. No obligation, of course.

Eligibility to apply for a life insurance policy from Teachers Insurance and Annuity Association is limited to employees of colleges, universities, private schools and certain other nonprofit educational or research institutions. The employee's spouse is also eligible, provided more than half of the combined earned income of husband and wife is from a qualifying institution.



Established as a nonprofit service organization by the Carnegie Foundation for the Advancement of Teaching

Life Insurance Advisory Center Teachers Insurance and Annuity Association 730 Third Avenue, New York, NY 10017	SP7909K1
Please send me the reprinted articles contrasting Term and Whole Life, and personal illustrations of these TIAA policies:	
<input type="checkbox"/> Decreasing Term <input type="checkbox"/> 5-Year Renewable Term <input type="checkbox"/> Whole Life	
Please print	
Name _____	
Job Title _____ Birthdate _____	
Home Address _____	
City _____ State _____ Zip _____	
Phone #s (optional) Home () Office ()	
Your nonprofit employer (college, univ., private school, research institution, etc.) _____	
If your spouse is also eligible under the rules given at left, please provide:	
Spouse's name _____ Birthdate _____	

Bon apatite!



The youngest member of the growing LKB family of Ultrogel® chromatography and gel filtration media is a unique beaded form of hydroxyapatite microcrystals trapped in a highly cross-linked agarose network. Tradenamed HA-Ultrogel®, the medium is recommended for adsorption chromatography when it is important to have high flow rates, resistance to pH as high as 13 and to SDS, urea and other denaturing agents.

Another new member is Ultrogel AcA-202, an agarose-polyacrylamide type very high in acrylamide, very low in exclusion range. It is ideal for desalting and for separation of peptides, oligosaccharides and other low MW compounds. Uniform, rigid beads provide much sharper peaks than previously available commercial products.

Still three other new additions to the LKB line are Ultrogels A-2, A-4 and A-6. Straight agarose types, they give excellent results with affinity chromatography as well as gel filtration. The very high exclusion limit of Ultrogel A-2 makes it a fine tool for gel filtration of large nuclear fragments and viruses.

Each of the new Ultrogel grades — and each of the four existing members of the LKB family — is characterized by a narrow, well-defined particle size distribution curve. You have no fines to contend with. And each grade comes pre-swollen, ready to use direct from the bottle. You can run your column the same day you prepare it.

If you have not examined LKB Ultrogel media lately, you are in for a pleasant surprise. Contact LKB today for full information on this growing family.

Ultrogel	% Acrylamide	% Agarose	MW fractionation range (linear)	MW exclusion limit (globular proteins)
AcA-202	20	2	1,000 to 15,000	22,000
AcA-54	5	4	5,000 to 70,000	90,000
AcA-44	4	4	10,000 to 130,000	200,000
AcA-34	3	4	20,000 to 350,000	750,000
AcA-22	2	2	100,000 to 1,200,000	3,000,000
A-6	0	6	25,000 to 2,400,000	4,000,000
A-4	0	4	55,000 to 9,000,000	20,000,000
A-2	0	2	120,000 to 25,000,000	50,000,000
HA	*	*	<5,000,000	5,000,000

*hydroxyapatite

LKB

LKB Instruments Inc.
12221 Parklawn Drive Rockville, MD 20852
301: 881-2510 Telex: 8-9682

Circle No. 356 on Readers' Service Card

LETTERS

Lunar Studies

Richard A. Kerr's carefully constructed article "Ten years later: Whence the moon?" (Research News, 20 July, p. 292) treats the scientific study of the moon—discovery, solution, and enigma—as a matter of geology and geochemistry. To a substantial degree, our present understanding of the moon arises as well from the extensive program of physical measurements that has always been an integral and substantial part of lunar science. In the space program, the work in lunar physics began even before Apollo, and included Explorer 35 as well as earlier attempts, although unsuccessful, extending back to the late 1950's (Pioneers 1 and 2 and Atlas-Able). In the Apollo program the physical exploration of the moon included automated observatories at fixed locations on the lunar surface (ALSEP), instruments carried by the astronauts, experiments set up by the astronauts and operated subsequently by remote control, the Apollo subsatellites, and, as part of the study of the moon's gravitational field, the Apollo command and service modules themselves. Since the emplacement of hardware and the acquisition of data depended upon the capabilities and expertise of the astronauts, lunar physics shared the manned experience with lunar geology. As with the geochemistry of the moon, the verdict is not in as to the meaning of all the information on the physics of the moon gathered in the course of this remarkable program. Thus, lunar physics shares some of the enigmas of lunar chemistry and geology, as it has shared in the discoveries. Clearly, any viable hypothesis of lunar origin and evolution must account for the relic magnetism of the moon. Even among those physical measurements which are not yet well understood conceptually, important connections have already been made with the other terrestrial planets. These include the absence of a so-far-detectable dipole moment, the seismicity of the crustal scattering layer, the structural support for the mascons, the heat flow from the interior, the rate of retreat of the moon from Earth, the interaction of solar wind with the moon and the associated downstream cavity and reimplantation of argon into lunar rocks by the associated electric field, and the gravitational field of the moon. The last of these has now provided meaningful constraints on the value of the mean moment of inertia, eagerly awaited by aficionados of the moon for some 200 years.

By way of emphasizing the interdisci-

plinary nature of lunar science, we note that the dating of samples places the remanent magnetic field in a chronological sequence, showing that the field was acquired very early (3.2 to 3.95 billion years ago) in lunar history. If the source of the background field is eventually found to have been the result of a regenerative dynamo operating within the moon, then, to a considerable degree, the very early thermal regime leading to the iron fractionation required for core formation, and the subsequent source of power for this dynamo, become fundamental problems in selenophysics and selenochemistry, jointly associated with the very young moon. Their answers would place serious constraints upon models of early lunar evolution and, in turn, of the evolution of the early solar system.

P. J. COLEMAN, JR.

*Institute of Geophysics and
Planetary Physics, University of
California, Los Angeles 90024*

M. FULLER

*Department of Geological Sciences,
University of California,
Santa Barbara 93106*

R. M. HOUSLEY

*Rockwell International Science Center,
Thousand Oaks, California 91360*

G. V. LATHAM

*Geophysics Laboratory, Marine
Science Institute, University of
Texas, Galveston 77550*

S. K. RUNCORN

*Lunar and Planetary Institute, School of
Physics, University of Newcastle upon
Tyne, Newcastle, United Kingdom*

W. J. SJOGREN

*Jet Propulsion Laboratory,
Pasadena, California 91103*

SEAN C. SOLOMON

*Department of Earth and Planetary
Sciences, Massachusetts Institute of
Technology, Cambridge 02139*

C. P. SONETT

*Department of Planetary Sciences,
Lunar and Planetary Laboratory,
University of Arizona, Tucson 85721*

Soviet Anti-Semitism: Reply by

Pontryagin

My name was mentioned several times in the article "Anti-Semitism alleged in Soviet mathematics" published in *Science* (News and Comment, 15 Dec. 1978, p. 1167). Except for the unfounded statement that I am an anti-Semite based on the evidence of persons unknown to me, the concrete accusations consist of the following:

1) "Pontryagin represents the Soviet

Union in the International Mathematical Union."

2) "He leads the editorial board that makes the final decision on every book in mathematics proposed for publication."

3) "He is editor of the prestigious journal *Matematicheski Sbornik*."

4) "He controls a voting bloc on the National Committee of Soviet Mathematicians and thereby determines the international contacts of Russian mathematicians."

5) "Two independent sources say that Soviet authorities were embarrassed by Pontryagin's denial of a visa to Margoulis, whose absence at the International Congress highlighted Pontryagin's policies."

6) "Pontryagin . . . was called to account for Margoulis' absence. Pontryagin angrily said that Margoulis could not be nominated for a Fields Medal because he was not a Russian nominee."

7) "The sources say he then got himself in deeper trouble by writing to every Soviet satellite, telling them they must support him in his threat."

8) ". . . Pontryagin invited a writer, Ivan Shetsov, . . . to speak at the prestigious Steklov Institute."

9) ". . . Soviet émigrés . . . report that since 1975 the number of published papers by Jewish authors in *Matematicheski Sbornik* has declined to zero."

10) "Particularly offensive to the Jewish mathematicians, sources say, is the habit of Pontryagin and his supporters of referring privately to this and several other journals as 'judenfrei'—the German expression for 'free of Jews.' "

I answer point by point:


1) I actually represented the Soviet Union on the executive committee of the International Mathematical Union from 1971 till 1978. But I have never done anything which could be considered as anti-Semitism. I believe my colleagues from the executive committees of 1971 to 1974 and 1975 to 1978 will confirm this.

2) I lead the mathematical group of one of the sections of the publishing house Nauka. This group inspects only a small part of the mathematical books printed in the Soviet Union and only prepares all necessary documentation for the mentioned section. I believe that this group has actually improved the publication activity, but it does not undertake any actions against Jews.

3, 9) I have actually been the editor-in-chief of the journal *Matematicheski Sbornik* since 1975. It may be seen from the table in *Science* that from 1970 till 1974 the number of Jewish papers in the journal constituted 36 percent of the total amount and from 1975 till 1978, 9 percent. I cannot consider that evidence of edito-

New!

A continuously variable length chromatography column



For gel filtration, ion exchange, affinity and adsorption chromatography—the new LKB 2137 column adjusts to any length from 5 to 95 cm. No longer do you need different length columns for each technique. Now, simplify inventory, cut costs dramatically!

In minutes you can have a completely jacketed high-resolution column of optimum length. For either descending or ascending flow. For sharper, more symmetrical peaks than ever before.

This versatile column comes complete with flow adaptors. The price? Less than you pay for other suppliers' flow adaptors alone.

*Contact LKB today
for full details.*


LKB

LKB Instruments Inc.
12221 Parklawn Drive
Rockville, Maryland 20852
301: 881-2510

13A-301

Circle No. 357 on Readers' Service Card

High, wise and handsome



The MultiRac fraction collector

LKB's new MultiRac™ fraction collector is a real space saver. On the opposite page you can see that reservoir, column, pump, monitor and recorder all mount on it easily. And you can keep adding decks to take as much more equipment as you want.

The new LKB fraction collector is bright too. It will collect from microliters to liters, adjust fraction size automatically according to OD, stop all flow as the head traverses, and on command will channel all void volume to waste.

Its good looks go well beneath the surface: solid state electronics, rugged materials of construction and a clear, bright, unambiguous LED display all make for an instrument that's safe, dependable and easy to use.

Contact LKB today for full details.

LKB

LKB Instruments Inc.
12221 Parklawn Drive Rockville, MD 20852
301: 881-2510 Telex: 8-9682

58A-311

Circle No. 358 on Readers' Service Card

rial anti-Semitism. The statement contained in 9) that the number of the Jewish papers has diminished to zero is false, as may be seen from those data. I emphasize that I cannot check the truthfulness of the information in the table, since we do not consider our authors from the point of view of their nationality.

4) I am actually a member of the National Committee of Soviet Mathematicians, which includes 25 persons; as one of them I really have an influence 1/25 the size of the committee.

5) The negative decision about Margoulis' participation in the International Congress was unanimously made at a sitting of the National Committee of Soviet Mathematicians.

6) I have not been asked any questions about Margoulis' absence from the congress, and I have not given any answers about that subject.

7) I have not written any letters to the socialist countries ("satellites") concerning the organization of the congress.

8) I did not invite the writer Shetsov to speak at the Steklov Institute.

10) I never use the expression "judenfrei," which is foreign to my native language; nor do I use any Russian expression equivalent to that.

I consider the article in *Science* to be founded on false information and to be a personal, malicious slander against me.

L. S. PONTRYAGIN

*Steklov Institute, Vavilov Street 33,
Moscow, Soviet Union*

Academic Ties with the Middle East

R. Jeffrey Smith's article on the controversies over Middle East investments in American universities (News and Comment, 2 Feb., p. 421) has recently come to my attention. As the article notes, a proposed agreement between the University of Pennsylvania and the Arab Development Institute in Libya was under consideration when the above institute withdrew the proposal. It was a controversial issue on our campus, and opinions were divided. However, I wish to state publicly my confidence in the fairness of the review process, which was taking place under procedures developed by our faculty, and my confidence in the judgment of the members of our faculty who were involved in the process, both the sponsors of the proposal and those who were reviewing it. There were no improper political pressures or discriminatory actions.

The University of Pennsylvania will maintain a strong interest in the Middle

East, as in many other areas of the world, and will welcome research and training programs that are consistent with its mission, standards, and policies.

VARTAN GREGORIAN

*Office of the Provost,
University of Pennsylvania,
Philadelphia 19104*

Nutrition Program

In a recent article, William J. Broad (News and Comment, 8 June, p. 1060) cites the infant nutrition program of the U.S. Department of Agriculture (USDA) at the Baylor College of Medicine and that of the National Institute of Child Health and Human Development (NICHD) as examples of needless duplication of research and lack of coordination. Unfortunately, the article does not acknowledge the efforts we have made to avoid these problems.

The language of the Conference Report of the House of Representatives No. 95-1579, 18 September 1978, directs the Children's Nutrition Laboratory (CNL) of USDA at Baylor "to insure that there is close cooperation and coordination" with NICHD. The Conference Report further stipulates "that the USDA and the NIH develop appropriate interagency mechanisms to insure coordination of this program" and that our two agencies make the "research and program plan of [the CNL] consistent with the purposes of the [NICHD]." In order to carry out these directives, we and our staffs have met six times within the past year and have established a firm working relationship. In November 1978 we held a series of intensive discussions with the staff of the CNL in Houston. During these meetings we worked out a cooperative agreement among USDA, NICHD, and the Baylor College of Medicine. This agreement establishes an interagency advisory committee for the CNL. Both of us belong to this advisory committee, which is designed to foster cooperation and avoid unnecessary duplication of research on infant nutrition.

In addition to this arrangement, we agreed that both agencies must concur in the selection of a director of the CNL. We also agreed which areas of infant nutrition research each agency will emphasize. NICHD will continue to foster research on the clinical aspects of infant nutrition, especially those relating to prevention, diagnosis, and treatment. The CNL will devote its research efforts to ascertaining nutritional requirements of normal infants and to defining stan-

dards for nutrient intake during infancy. The CNL will also work to develop precise methodology for such determinations.

Our agencies share a keen interest in research on lactation, human milk, and colostrum. In fact, NICHD has contracted with the Baylor College of Medicine to develop procedures for the collection, storage, processing, and distribution of milk for clinical research. Some of the staff of the CNL work on this project, and the staff of the Human Nutrition Center of USDA follow the progress of this contract with interest. Cooperation, not competition, best describes this venture.

We are now in the midst of planning a conference on optimal nutrition during pregnancy that our agencies will sponsor jointly. We anticipate more such cooperative enterprises in the near future.

NORMAN KRETCHMER

*National Institute of Child Health
and Human Development,
Bethesda, Maryland 20205*

D. MARK HEGSTED

*Human Nutrition Center,
Science and Education Administration,
U.S. Department of Agriculture,
Washington, D.C. 20250*

Interferon

In her two articles on interferon research (Research News, 15 June, p. 1183; 22 June, p. 1293), Jean L. Marx does an admirable job summarizing the many aspects of work in a most active field. I realize it is impossible to mention all investigators who made significant contributions. However, I would like to point out the omission of the name of Kurt Paucker of the Medical College of Pennsylvania, who has made many important original observations in interferon research. One of his milestone contributions (with K. Cantell and W. Henle) was the demonstration that interferon can inhibit cell division.

In Marx's first article, Edward Havell and I are given credit for finding "that the principal interferon made by leukocytes is not the same as the major form of interferon produced by fibroblasts." Although Havell and I did make contributions in this area, I believe the major portion of work elucidating antigenic differences between leukocyte and fibroblast interferon was done in Paucker's laboratory.

JAN VILČEK

*New York University School of
Medicine, New York 10016*

14 SEPTEMBER 1979

The first microprocessor-controlled fraction collector... LKB's MultiRac^{T.M.}



LKB's MultiRac fraction collector is unlike any other in versatility of function, reliability of operation and simplicity of use.

Just push a few buttons and you can program the unit to collect in any size vessel from test tubes to carboys — by time, drop or LKB's precise volume. The 250 msec switch-over between tubes minimizes dribbling and a "stop" feature can eliminate it entirely.

Microprocessor control conserves test tubes by diverting column void volume to a built-in waste container. Still more tubes are conserved by varying fraction sizes according to OD. During peak elution, small fractions are collected for optimal resolution... in the valleys, tubes are filled to capacity. You can start your run and leave for an early dinner, knowing that when you come back the next day, not a single tube will have been wasted.

Non-tipping spring-loaded racks accommodate everything from 8 mm test tubes to 28 mm scintillation vials to funnels feeding carboys. Racks fit in a tray that is removed in one easy movement.

Keeping wet operations below the circuitry and using tough, inert materials make the MultiRac safe and easy to use with all liquids, even aggressive solvents and radioactive solutions. Power interruptions are also tolerated — a memory retains data up to 16 hours. And in the unlikely event that something *should* go wrong, a self-diagnosis function helps you locate the problem.

By the way — fraction collector, pump, column, monitor and recorder all occupy less than two square feet of bench space.

LKB

LKB Instruments Inc.

12221 Parklawn Drive Rockville, MD 20852
301: 881-2510 Telex: 8-9682

Circle No. 359 on Readers' Service Card

57A-310

This is the microprocessor meter that simplifies electrode analyses.

Here's how:



The operating mode now selected gives you ion or dissolved gas concentrations directly by the standard calibration method. Other modes let you read millivolts, pH, or concentrations by the various incremental methods without the complication of having to enter analate and standard solution volumes.

The concentration of your standardizing solution (or buffer) is entered in the unit system of your choice. The digital switches let you see what value has been selected at all times.

The electrode slope value is entered (and adjusted as required) on these digital switches. Whenever you use the 901, you can see at a glance if your analytical parameters are correctly entered – without the complicated recall procedures of other microprocessor meters.



You standardize the 901 at the push of a button – and you're reading concentrations directly on the LED display. If you'd like to see for yourself how simple it is to analyze your samples with ORION electrodes and the 901, call our toll-free number and we'll arrange for one of our technical specialists to set up a full-fledged demonstration in your laboratory.

ORION RESEARCH

(800) 225-1480 in U.S.A. • (617) 864-5400 in MA • (800) 363-9270 in Canada

AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE

Science serves its readers as a forum for the presentation and discussion of important issues related to the advancement of science, including the presentation of minority or conflicting points of view, rather than by publishing only material on which a consensus has been reached. Accordingly, all articles published in *Science*—including editorials, news and comment, and book reviews—are signed and reflect the individual views of the authors and not official points of view adopted by the AAAS or the institutions with which the authors are affiliated.

Editorial Board

1979: E. PETER GEIDUSCHEK, WARD GOODENOUGH, N. BRUCE HANNAY, MARTIN J. KLEIN, FRANKLIN A. LONG, NEAL E. MILLER, JEFFREY J. WINE
1980: RICHARD E. BALZHISER, WALLACE S. BROECKER, CLEMENT L. MARKERT, FRANK W. PUTNAM, BRYANT W. ROSSITER, VERA C. RUBIN, MAXINE F. SINGER, PAUL E. WAGGONER, F. KARL WILLENBROCK

Publisher

WILLIAM D. CAREY

Editor

PHILIP H. ABELSON

Editorial Staff

<i>Managing Editor</i> ROBERT V. ORMES	<i>Business Manager</i> HANS NUSSBAUM
<i>Assistant Managing Editor</i> JOHN E. RINGLE	<i>Production Editor</i> ELLEN E. MURPHY

News Editor: BARBARA J. CULLITON
News and Comment: WILLIAM J. BROAD, LUTHER J. CARTER, CONSTANCE HOLDEN, ELIOT MARSHALL, DEBORAH SHAPLEY, R. JEFFREY SMITH, NICHOLAS WADE, JOHN WALSH. *Editorial Assistant:* SCHERRAINE MACK

Research News: BEVERLY KARPLUS HARTLINE, FREDERICK F. HARTLINE, RICHARD A. KERR, GINA BARI KOLATA, JEAN L. MARX, THOMAS H. MAUGH II, ARTHUR L. ROBINSON. *Editorial Assistant:* FANNIE GROOM

Consulting Editor: ALLEN L. HAMMOND

Associate Editors: ELEANORE BUTZ, MARY DORFMAN, SYLVIA EBERHART, JUDITH GOTTLIEB, RUTH KULSTAD

Assistant Editors: CAITILIN GORDON, LOIS SCHMITT

Book Reviews: KATHERINE LIVINGSTON, *Editor:* LINDA HEISERMAN, JANET KEGG

Letters: CHRISTINE KARLIK

Copy Editor: ISABELLA BOULDIN

Production: NANCY HARTNAGEL, JOHN BAKER; YA LI SWIGART, HOLLY BISHOP, ELEANOR WARNER; MARY MCDANIEL, JEAN ROCKWOOD, LEAH RYAN, SHARON RYAN

Covers, Reprints, and Permissions: GRAYCE FINGER, *Editor:* CORRIE HARRIS, MARGARET LLOYD

Guide to Scientific Instruments: RICHARD SOMMER

Assistant to the Editors: RICHARD SEMIKLOSE

Membership Recruitment: GWENDOLYN HUDDLE

Member and Subscription Records: ANN RAGLAND
EDITORIAL CORRESPONDENCE: 1515 Massachusetts Ave., NW, Washington, D.C. 20005. Area code 202. General Editorial Office, 467-4350; Book Reviews, 467-4367; Guide to Scientific Instruments, 467-4480; News and Comment, 467-4430; Reprints and Permissions, 467-4483; Research News, 467-4321. Cable: Advancisci, Washington. For "Instructions for Contributors," write the editorial office or see page xi, *Science*, 29 June 1979.

BUSINESS CORRESPONDENCE: Area Code 202. Membership and Subscriptions: 467-4417.

Advertising Representatives

Director: EARL J. SCHERAGO

Production Manager: MARGARET STERLING

Advertising Sales Manager: RICHARD L. CHARLES

Marketing Manager: HERBERT L. BURKLUND

Sales: NEW YORK, N.Y. 10036: Steve Hamburger, 1515 Broadway (212-730-1050); SCOTCH PLAINS, N.J. 07076: C. Richard Callis, 12 Unami Lane (201-889-4873); CHICAGO, ILL. 60611: Jack Ryan, Room 2107, 919 N. Michigan Ave. (312-DE-7-4973); BEVERLY HILLS, CALIF. 90211: Winn Nance, 111 N. La Cienega Blvd. (213-657-2772); DORSET, VT. 05251: Fred W. Dieffenbach, Kent Hill Rd. (802-867-5581)

ADVERTISING CORRESPONDENCE: Tenth floor, 1515 Broadway, New York, N.Y. 10036. Phone: 212-730-1050.

Education for the 21st Century

Education stands out as the best basis for hope that this country and others will somehow manage to avoid enormous trauma during the transitions that lie ahead. But are young people being well counseled in preparation for the 21st century? What will be the shape of the future and the corresponding demands for trained people?

We have lived in an era of conspicuous consumption, rising expectations, and exponential growth. We have entered a period of uncertainty, conflicts in value systems, and possible decline in living standards. We face the necessity of creating an economy based on less oil, less energy, and changed raw materials. Society will not and cannot return to the Stone Age or even to a universal simple agrarian existence. It will not abandon knowledge nor the ability to harness it. There will be no shortage of people ready to tell the politicians what should be done. There will be no shortage of instant solutions. But there may be a shortage of trained people capable of meeting society's physical needs.

In many ways the educational system is excellent, but in at least one respect it can be faulted. Its overall performance in counseling the young has been mediocre. Universities have practiced a policy of conducting an intellectual smorgasbord.

Since the students have not received adequate guidance from the educational institutions, their decisions have been based on other sources of information, notably the media. In the past this has resulted in some tragic wastes of talent. The hoopla of the "space age" led to career decisions that brought later disappointments. A lesser example was a vogue in oceanography. Currently it is fashionable to attempt to go to medical school. The pressures are easing somewhat but an enormous wastage continues. Stories appearing in the media can now attenuate the flow. A notable example a decade ago was the impact of a few stories about engineers driving taxicabs. Only now have enrollments in engineering recovered from the slump.

A basic problem in providing career guidance is that no one knows precisely what will be needed decades hence. In addition, at the universities campus politics gets in the way. At a time of budget squeezes and of departmental support based on enrollments, the fight for survival makes it difficult for faculty to be objective in advice. Another problem is that neither individual students nor their potential counselors can accurately gauge aptitudes, talents, drive, and judgment at the time of entry into the university. Few students are aware of their own potentials or are in command of themselves. An additional complication for potential scientists and engineers is that many curricular decisions are, in practice, irreversible. Basic courses in mathematics, chemistry, and physics must be taken in a timely fashion or options are foreclosed.

Whatever the changing shape of society, scientists and engineers will have essential roles. The uncertainties, though, make it advisable to caution against excessive specialization. In contrast, it seems desirable to adopt policies of maximum flexibility, of preservation of options, of being prepared to pursue lifelong learning.


For those who are capable of handling mathematics and abstract reasoning, this means building the necessary foundation in mathematics, chemistry, physics, and engineering design as well as obtaining adequate competence in verbal and written communications. Those with the necessary aptitudes should also be exposed to biology, earth sciences, behavioral sciences, and the humanities.

One cannot insist that students with no aptitude for the hard sciences take the fundamental courses, but it is difficult to visualize how they could be comfortable in the world of tomorrow without some inkling of the forces shaping their lives.

For their part, the universities should have a searching look at their counseling policies. There must be better ways than entrusting young lives to a hit-or-miss system.—PHILIP H. ABELSON

Philips EM-400 Twin System

Imaging and analysis with no compromise



The EM-400 Electron Microscope with Twin Lens objective provides atomic level resolution combined with nano-PROBE capability for analysis of specimen areas as small as a billionth of a meter. For the first time, these capabilities are available to you in TEM, S(T)EM and SEM modes with no compromise in performance.

No compromise means TEM lattice resolution to 0.14 nanometres with $\pm 60^\circ$ of specimen tilt. Single switch operation to S(T)EM, SEM or micro-micro diffraction modes with probe sizes to 2.0 nanometres. X-ray detector collection angle of 0.13 sr — more than double currently available sensitivity.

No compromise in signal handling means providing six different signals. And the ability of mixing these in any sum or difference combination and display them as mono/dual, dual magnification and stereo images.


There's more. The EM-400 was specifically designed to provide an ultra-clean, high vacuum specimen environment to insure consistently successful microanalysis.

Don't compromise. Get complete details. Write to Philips Electronic Instruments, Electron Optics Group, 85 McKee Dr., Mahwah, NJ 07430.

Electron Diffraction pattern of CdS from a single grain in a solar cell. Courtesy of Dr. J.W. Edington, Dept. of Aeronautical Engineering, Univ. of Del.




**Electronic
Instruments**



PHILIPS

Circle No. 270 on Readers' Service Card



LUMINESCENCE... What can a lifetime do for you?

PRA's pulsed nanosecond fluorometer systems now make the measurement of luminescence lifetimes routine procedure with highest result confidence.

Uses:

- Elucidate complex photophysical mechanisms
- Identify the presence of more than one fluorescent compound
- Open new ways of probing luminescent systems
- Can be highly sensitive to sample impurities
- Lamp can be used as pulsed light source
- Determine lifetimes as short as 200 picoseconds
- Less than 100 picoseconds time uncertainty

Features:

- High signal to noise ratio
- Low RF interference
- Low scattered light levels
- Complete computer interface and software options available



Competent consulting staff with broad experience in instrumentation and applications to advise you

PRA

**Photochemical Research
Associates Inc.**
45 Meg Drive, London, Ontario
Canada N6E 2V2
(519)686-2950 Telex 064-7597
U. S. OFFICE: ROGER LEACH (615) 483-3433

Circle No. 305 on Readers' Service Card

A new approach in tissue culture




PETRIPERM

a novel culture dish
with gas permeable
membrane as cellular
support


Additional advantages of the PETRIPERM membrane are:

- UV permeable 200 nm
 - Suitable for high power light and fluorescence microscopy
 - The membrane can be cut for cell cloning, electron microscopy, staining and documentation
 - Chemically resistant to acids, bases and organic solvents
- PETRIPERM is available either with a hydrophilic surface for optimal attachment, or with a hydrophobic surface for easy detachment of adherent cells.



In the intact animal, the pH, and the partial pressure of both oxygen and CO₂ of the cell environment are under total biochemical and physiological control. This state is not even approximated by available techniques of monolayer cell culture in conventional tissue culture ware. Continuously changing gradients of these parameters in the micro-environment of cultured cells will obscure the biological significance of culture dependent phenomena. PETRIPERM is the prerequisite for the precise control of pO₂, pCO₂ and pH in the cellular environment as gases reach the cultured cells directly through the highly gas permeable 25 μm thick membrane.

**Heraeus-
CO₂-Incubator
B 5060 EK/CO₂**
We build the
incubator with
appropriate
pCO₂, pH and pO₂
control units



W.C. HERAEUS GMBH
Produktbereich
Elektrowärme
D-6450 HANAU

Heraeus

Circle No. 38 on Readers' Service Card

PRODUCTS and MATERIALS

Growth Factor

CR-ECGS supplements the growth of human vascular endothelial cells in vitro. This makes it possible for investigators to maintain these cells at low densities in monolayers. When it is added in recommended concentration to serum-containing medium, CR-ECGS promotes a six- to tenfold proliferation of human vascular endothelial cells over 10 days. It also supports repeated passage of established cultures of these cells. Collaborative Research. Circle 795.

Water Purification

The Critical Applications Adsorption column removes organic material from deionized or distilled water for high-performance liquid chromatography or gas chromatography-mass spectrometry. It consists of a column packed with an adsorbent medium, a diaphragm control valve, a membrane filter, a holder, and mounting accessories. Filtering is done with activated carbon and a 0.2-micrometer polycarbonate membrane filter. Hydro Service & Supplies. Circle 798.

Sampler for Respiratory and Metabolic Studies

Douglas bags are used for collecting expired air for volume analysis and composition studies. This line is constructed of 0.020-inch polyvinyl chloride with electrically sealed seams. Advantages include light weight, flexibility, durability, and retentive power. Inside diameter of the tubing that connects the mouthpiece to the bag is 30 millimeters. The mouth-


piece features a one-way valve and is removable for sterilization. Bags are offered in 50-, 100-, and 150-liter capacities, and in a 2-liter sampling size. Ealing. Circle 799.

Spectrophotometer

Model 100-80 is a computerized, ultraviolet-visible, double-beam, ratio-recording spectrophotometer. It features flexibility in wavelength scanning and expanded ability to collect, process, and report quantitative data. An x-y recorder and a digital printer are integral. Standard functions include wavelength scanning, repetitive scanning, calculation of derivatives, wavelength programming, time scan, automatic correction for base line, correction of curves, and automatic printing. A complete line of accessories is available. NSI Hitachi Scientific Instruments. Circle 794.

Total Exhaust Hood

The Biofume model is a laminar-flow biological safety cabinet and fume hood of class II type B. No air recirculates in the cabinet. The supply HEPA filter is located above the work area. Additional ambient air enters an opening at the bottom of the viewing panel at sufficient ve-



Newly offered instrumentation, apparatus, and laboratory materials of interest to researchers in all disciplines in academic, industrial, and government organizations are featured in this space. Emphasis is given to purpose, chief characteristics, and availability of products and materials. Endorsement by *Science* or AAAS is not implied. Additional information may be obtained from the manufacturers or suppliers named by circling the appropriate number on the Reader Service Card (on pages 1074A and 1170E) and placing it in the mailbox. Postage is free.

—RICHARD G. SOMMER

locity to contain airborne hazards within. There are 4- and 6-foot models with stainless steel work areas and trim. Germfree Laboratories. Circle 796.

Satellite-Controlled Clock

The SCC is a microprocessor-controlled time code receiver and decoder. It derives its precision (30 to 50 microseconds relative to Coordinated Universal Time) from a signal relayed by the GOES east or west synchronous satellite to virtually any location in the Western Hemisphere. The instrument requires no calibration, correction, or initial setting once latitude and longitude are entered. In addition to a selection of pulse-coded outputs, time is displayed in hours, minutes, and seconds on a digital readout. LaBarge. Circle 797.

Fraction Collector

Model 1850 features a solvent-resistant, stainless steel tray which may be lifted out for sample processing and storage or for washing. Test tubes from 12 to 18 millimeters in diameter are held in 35 small self-standing racks; total capacity is up to 210 tubes. Digital displays continuously indicate fraction intervals in time, drop, or volumetric units. Number of units deposited per tube and number of tubes filled are also indicated. ISCO. Circle 800.

Literature

Ultra-Low Temperature Freezers includes chest, upright, and walk-in models with guaranteed performance to -85°C in a 32.2°C ambient. Forma Scientific. Circle 801.

Gel Filtration Theory and Practice is a handbook of experimental design, procedures, and applications. Pharmacia Fine Chemicals. Circle 802.

Incident Light Microscopy covers theory and principles from sample preparation through optics and photography. Carl Zeiss. Circle 804.

Gas Chromatograph is devoted to the model 560 glass capillary GC's features, applications, and specifications. Tracor Instruments. Circle 805.

Glassware and Equipment is a catalog of scientific apparatus for general laboratories. Reliance Glass Works. Circle 806.


Spectrodensitometer describes the SD 3000 instrument for quantitative analysis of thin-layer chromatograms. Schoeffel Instrument Division, Kratos. Circle 807.

"OPTO-VARIMEX"

Ultimate Animal Activity Meter

MEASURES TRAVELED
DISTANCE, PLOTS PATTERN OF
MOVEMENTS.

Measures
Animal
Respiration
and Tremors



- * Measures horizontal activity by interruption of invisible infrared beams criss-crossing animal cage 1" above the bottom.
- * Measures rearing and jumping activity by means of two adjustable light curtains.
- * Plots pattern of animal movements on X-Y chart recorder by decoding position of animal from interrupted IR beams.
- * *Completely* insensitive to urine bedding or feces as IR beams cross the cage 1" above bottom.
- * Large 15" x 15" transparent plastic cage—not attached to the instrument.
- * Measures tremors and respiration by proximity sensor.



for more information call or write:
COLUMBUS INSTRUMENTS

950 NORTH HAGUE AVENUE
COLUMBUS, OHIO 43204 U.S.A.
PHONE (614) 488-6176

Circle No. 324 on Readers' Service Card

MONOCLONAL ANTI-THY 1.2

Extraordinarily Specific

To identify,
localize, and
selectively isolate
mouse T-cell lymphocytes.
Highly cytotoxic, high titer
antibody.

Superior to conventionally
produced antisera.

Monoclonal Anti-Thy 1.2

~1mg in 1ml buffer system

NEI-001

~0.1mg in 1ml buffer system

NEI-001A

Not for use in humans or clinical diagnosis.



New England Nuclear

549 Albany Street, Boston, Mass. 02118
Call toll-free: 800-225-1572

(In Massachusetts and International: 617-482-9595)

NEN Chemicals GmbH, Dreieich, W. Germany; NEN Canada Ltd., Lachine, Quebec
Circle No. 303 on Readers' Service Card

A photographic base now gives the Swift M100 Tissue Culture Microscope greater-than-ever versatility, with capabilities for standard single lens reflex and/or Polaroid® cameras. It's an unbeatable team for advanced research in living cells and biological specimen. You will be pleased with the sheer practicality of this inverted instrument with its erect, natural image; large, sturdy extension stage; brilliant, variable intensity illumination; large nosepiece and

smooth, precise focusing; angled 45° binocular or monocular heads with eyepoint 335mm above table surface for comfortable viewing; exceptionally fine optics. Easily adapted to phase contrast and/or polarizing microscopy by specific configuration or accessories.

Write or call now
for literature
and name of
nearest dealer for
demonstration.


Now, right on target for tissue culture photography



SWIFT INSTRUMENTS, INC.

Scientific Instrument Division
P.O. BOX 562, SAN JOSE, CA 95106 • 408/293-2380
(MAIN OFFICE: Boston, MA)

SWIFT AGENCIES are located throughout the U.S. and in most foreign countries.



Circle No. 37 on Readers' Service Card

HEALTH CARE



Project Hope

Regulation, economics, ethics, and practice are just a few of the topic areas discussed in this 7th compendium from SCIENCE. This highly acclaimed publication will help the professional and the lay person understand not only recent advances in medical practice but other less obvious but equally vital components of our health-care system.

HEALTH CARE: Regulation, Economics, Ethics, Practice

Edited by Philip H. Abelson

256 pp., 8½" x 11";
0-87168-302-4 (cloth), \$14.00*;
0-87168-240-0 (paper), \$6.00*

November 1978

To order your copy, send your name, address and remittance to



AAAS Department HC-3
1515 Massachusetts Avenue, NW
Washington, D. C. 20005

Please allow 6-8 weeks for delivery.

*AAAS members deduct 10%.

BOOKS RECEIVED

(Continued from page 1124)

Biochemical Thermodynamics. M. N. Jones, Ed. Elsevier, New York, 1979. xii, 390 pp., illus. \$73.25. Studies in Modern Thermodynamics 1.

Biochemistry of Cell Walls and Membranes II. J. C. Metcalfe, Ed. University Park Press, Baltimore, 1978. xii, 288 pp., illus. \$29.50. International Review of Biochemistry, vol. 19.

Carcinogens. Identification and Mechanisms of Action. Papers from a symposium, Houston, 1978. A. Clark Griffin and Charles R. Shaw, Eds. Published for the University of Texas System Cancer Center M. D. Anderson Hospital and Tumor Institute by Raven, New York, 1979. xviii, 488 pp., illus. \$38.

Cardiovascular Physiology III. Arthur C. Guyton and David B. Young Eds. University Park Press, Baltimore, 1979. xii, 368 pp., illus. \$24.95. International Review of Physiology, vol. 18.

Central Nervous System Pharmacology. A Self-Instruction Text. Donald E. McMillan. Little, Brown, Boston, ed. 2, 1979. xiv, 168 pp. Paper, \$12.50.

Characterisation of Protein Conformation and Function. Papers from a school, Wye, England, Sept. 1977. Felix Franks, Ed. Symposium Press, London, 1978. vi, 158 pp., illus. Cloth, \$40; paper, \$19.

Chemical and Biochemical Applications of Lasers. Vol. 4. C. Bradley Moore, Ed. Academic Press, New York, 1979. xii, 414 pp., illus. \$25.

Chemical Bonding. Audrey L. Companion. McGraw-Hill, New York, ed. 2, 1979. xii, 180 pp., illus. Cloth, \$7.95; paper, \$5.95.

Chemical Principles. Richard E. Dickerson, Harry B. Gray, and Gilbert P. Haight, Jr. Benjamin/Cummings, Menlo Park, Calif., ed. 3, 1979. xx, 944 pp., illus. + appendixes. \$18.95.

Chemistry. Karen Timberlake. Harper and Row, New York, ed. 2, 1979. xvi, 548 pp., illus. + appendixes. \$14.95.

Chemistry. Principles and Applications. Michell J. Sienko and Robert A. Plane. McGraw-Hill, New York, 1979. xii, 692 pp., illus. \$19.

Circulatory System Dynamics. Abraham Noordergraaf. Academic Press, New York, 1978. xvi, 352 pp., illus. \$36. Biophysics and Bioengineering Series, 1.

Climber's Guide to Devil's Lake. William Widule and Sven Olof Swartling. University of Wisconsin Press, Madison, 1979. 200 pp., illus. Paper, \$8.95.

Computer Technology. Logic, Memory, and Microprocessors. A Bibliography. A. H. Agajanian, Ed. IFI/Plenum, New York, 1978. xiv, 346 pp. \$95. IFI Data Base Library.

Conditioning. An Image Approach. Donald L. King. Gardner Press, New York, 1979 (distributor, Halsted [Wiley], New York). xii, 502 pp., illus. \$27.95.

Conservation Directory 1979. A List of Organizations, Agencies, and Officials Concerned with Natural Resources Use and Management. Jeannette Bryant, Ed. National Wildlife Federation, Washington, D.C., 1979. x, 272 pp. Paper, \$4.

Contested Terrain. The Transformation of the Workplace in the Twentieth Century. Richard Edwards. Basic, New York, 1979. x, 262 pp. \$12.95.

Controversies about Stuttering Therapy. Hugo H. Gregory, Ed. University Park Press, Baltimore, 1979. x, 326 pp. \$18.95.

The Courage to Fail. A Social View of Organ Transplants and Dialysis. Renée C. Fox and Judith P. Swazey. University of Chicago Press, Chicago, ed. 2, 1978. xx, 438 pp. Paper, \$6.95.

Current Issues in Clinical Geriatrics. Michael B. Miller. Tiresias Press, New York, 1979. xii, 244 pp. \$13.50.

Cytochrome Oxidase. Proceedings of a seminar, Osaka, Japan, Aug. 1978. Tsao E. King, Yutaka Orii, Britton Chance, and Kazuo Okunuki, Eds. Elsevier/North-Holland, New York, 1979. xiv, 426 pp., illus. \$55.50. Developments in Biochemistry, vol. 5.

Dangerous Behavior. A Problem in Law and Mental Health. Calvin J. Frederick, Ed. National Institute of Mental Health, Rockville, Md., 1978 (available from the Superintendent of Documents, Washington, D.C.). viii, 192 pp. Paper, \$3. Crime and Delinquency Issues. DHEW Publication No. (ADM) 78-563.

Death and Dying. Challenge and Change. Robert Fulton, Eric Markusen, Greg Owen, and Jane L. Scheiber, Eds. Addison-Wesley, Reading, Mass., 1978. xviii, 476 pp., illus. Paper, \$9.95. A Course by Newspaper Reader.

Design for Accessibility. Robert James Sorensen. McGraw-Hill, New York, 1979. viii, 264 pp., illus. \$19.50.

The Determination of Dynamic Properties of Polymers and Composites. B. E. Read and G. D. Dean. Halsted (Wiley), New York, 1979. x, 208 pp., illus. \$49.95.

Development of Sensory Systems. Marcus Jacobson, Ed. Springer-Verlag, New York, 1978. xii, 472 pp., illus. \$115. Handbook of Sensory Physiology, vol. 9.

Developmental Biology of Acetabularia. Proceedings of a meeting, Moscow, Aug. 1978. Silvano Bonotto, Valentin Kefeli, and Simone Puiseux-Dao, Eds. Elsevier/North-Holland, New York, 1979. xxiv, 312 pp., illus. \$49. Developments in Cell Biology, vol. 3.

Developments in Cardiovascular Medicine. Proceedings of a symposium, London, July 1978. C. J. Dickinson and J. Marks, Eds. University Park Press, Baltimore, 1978. x, 372 pp., illus. \$24.50.

Diagnosis and Therapy of Porphyrrias and Lead Intoxication. Proceedings of a symposium, Marburg, Germany, June 1977. Manfred Doss, Ed. Springer-Verlag, New York, 1978. xxii, 310 pp., illus. Paper, \$34.

Diagnostic Methods in Viral Hepatitis. Proceedings of a symposium, London, Jan. 1978. P. S. Gardner and Colin R. Howard, Eds. Liss, New York, 1978. xii, 92 pp., illus. \$12. Laboratory and Research Methods in Biology and Medicine, vol. 2. Reprinted from *Journal of Medical Virology*, vol. 3, No. 1.

Dictionary of Biomedical Acronyms. Jacques Dupuyrat. Editions SCM, Paris, 1978. 110 pp. Paper, 40 F.

Dictionary of Scientific Biography. Vol. 15, Supplement I, Roger Adams—Ludwik Zejszner; Topical Essays. Charles Coulston Gillispie, Ed. Scribner, New York, 1978. xiv, 818 pp., illus. \$55.

Environmental Systems. Philosophy, Analysis and Control. R. J. Bennett and R. J. Chorley. Princeton University Press, Princeton, N.J., 1979. xii, 624 pp., illus. \$65.

The Essential Tension. Selected Studies in Scientific Tradition and Change. Thomas S. Kuhn. University of Chicago Press, Chicago, 1979. xxiv, 366 pp. Paper, \$5.95. Reprint of the 1977 edition.

Experimental Methods in Heavy Ion Physics. K. Bethge, Ed. Springer-Verlag, New York, 1978. vi, 254 pp., illus. Paper, \$14.30. Lecture Notes in Physics, vol. 83.

LIGHTWEIGHT AND PORTABLE METALLURGICAL MICROSCOPE WITH FULL BHM SYSTEM CAPABILITY

Newest member of the BHM family of precision microscopes is this metallurgical focusing unit on a lightweight, portable universal post mounting, or for mounting on any machine or instrument you desire.

Especially suited for use in both field and factory, it offers an extreme focusing range, with removable stage plate permitting the BHMJ to be placed directly on a workpiece. Or by raising it to the top of its mounting post, it accommodates specimens up to approximately 55mm.

Available with Neo achromats for brightfield and darkfield, or plan achromats or achromats for brightfield. A wide choice of light sources, eyepieces, observation tubes, phototube, stage micrometers, polarizers and filters and other precision accessories is available for any metallurgical requirement.

A very special instrument for production and inspection equipment, with particularly versatile application for the microelectronics industry. For complete details, write Olympus Corporation of America, 4 Nevada Drive, New Hyde Park, New York 11042.



OLYMPUS BHMJ

Canada: W. Carsen Co., Ltd., Ontario.




Circle No. 33 on Readers' Service Card

Now Available!

TABLE OF MOLECULAR WEIGHTS

A Companion Volume to THE MERCK INDEX



The new *Table of Molecular Weights* is a 300-page listing of more than 9,000 chemical compounds — pharmaceuticals, pesticides, industrial chemicals, solvents, and common reagents — in increasing order of molecular weight. Each chemical is referenced, for full encyclopedic, information, to *The Merck Index*, 9th Edition. In addition to the molecular weights, reported to six decimal places, this volume gives the empirical formula for each compound named. It also contains a table of masses and abundances of naturally occurring isotopes. Page size is 6-7/8" x 9-7/8", in durable, gold-stamped Lexotone® cover. Mass spectroscopists, chemists, biochemists, pharmacists, botanists, physicists, chemical engineers, and other life scientists will find the *Table of Molecular Weights* an important addition to their reference library.

The *Merck Index*, Ninth Edition is an encyclopedia of chemicals, drugs, and biologicals, complete with a cross-index of 50,000 synonyms, formula index, Chemical Abstracts registry numbers, and a variety of valuable information. Of some 10,000 monographs, more than half are illustrated with stereochemical structural formulas, and 7,500 offer information on uses and toxicity.

Merck & Co., Inc. • Publications Department
P.O. Box 2000 • Rahway, New Jersey 07065

SC109

Please send on approval

- ☐ TABLE OF MOLECULAR WEIGHTS . . . \$12.00
☐ THE MERCK INDEX, NINTH EDITION . . . \$18.00
☐ Payment enclosed (Merck pays shipping and any sales tax).
☐ Bill me for book(s) plus shipping and any tax.

Name _____

Address _____

City _____

State _____ Zip _____


Circle No. 273 on Readers' Service Card

Vibration protection for sensitive instruments! MICRO-g[®] Tables and Bench Top Units.

Low cost MICRO-g air-isolation systems can reduce building vibration at least 95 to 98%. Even low frequency (7-25 Hz) vibrations that conventional static systems can't handle. The self-leveling top surface is supported on ultra-sensitive air-pistons that adjust automatically to changing loads. Performance is guaranteed.

Standard MICRO-g Tables and Bench Top units cost approximately \$1000 and will support equipment loads of over 900-lbs. All that is needed for operation is a compressed air supply or bottled air or gas.

If you're using microscopes, microtomes, interferometers, mask alignment tools or other ultra-sensitive equipment you need to know about MICRO-g systems. There are six standard configurations in a wide range of sizes. Special configurations are readily available. Write for a catalog and price list today.



BACKER-LORING CORPORATION

65 Walnut Street, Peabody, Mass. 01960 • 617-531-3600

BL-1

Experimental Organic Chemistry. David Todd. Prentice-Hall, Englewood Cliffs, N.J., 1979. xxxvi, 348 pp., illus. Paper, \$13.95.

Fibrinolytics and Antifibrinolytics. Fritz Markwardt, Ed. Springer-Verlag, New York, 1978. xxii, 656 pp., illus. \$160. Handbook of Experimental Pharmacology, vol. 46.

Fifth Midwest Prairie Conference Proceedings. Ames, Iowa, Aug. 1976. David C. Glenn-Lewin and Roger Q. Landers, Jr., Eds. Extension Courses and Conferences, Iowa State University, Ames, 1978. 230 pp., illus. Paper, \$3.50.

Fine Structure of Parasitic Protozoa. An Atlas of Micrographs, Drawings and Diagrams. E. Scholtyseck. Springer-Verlag, New York, 1979. viii, 208 pp. Paper, \$39.80.

Fish Physiology. Vol. 7, Locomotion. W. S. Hoar and D. J. Randall, Eds. Academic Press, New York, 1978. xx, 576 pp., illus. \$47.50.

Fluid Mechanics. An Introductory Course. Pau-Chang Lu. Iowa State University Press, Ames, 1979. xxvi, 468 pp., illus. \$20.

Fluid Mechanics, Thermodynamics of Turbomachinery. S. L. Dixon. Pergamon, New York, ed. 3, 1978. xvi, 264 pp., illus. Cloth, \$25; paper, \$10. Thermodynamics and Fluid Mechanics Series.

Genetics of Industrial Microorganisms. Proceedings of a symposium, Madison, Wis., June 1978. O. K. Sebek and A. I. Laskin, Eds. American Society for Microbiology, Washington, D.C., 1979. xii, 284 pp., illus. \$12.

Geochemical Processes. Water and Sediment Environments. A. Lerman. Wiley-Interscience, New York, 1979. x, 482 pp., illus. \$29.95.

Geochemistry and the Environment. Vol. 3. Distribution of Trace Elements Related to the Occurrence of Certain Cancers, Cardiovascular Diseases, and Urolithiasis. Papers from a workshop, Captiva Island, Fla., Oct. 1974. National Academy of Sciences, Washington, D.C., 1978. xii, 200 pp., illus. Paper, \$14.

Géologie des Marges Continentales. G. Boillot. Masson, New York, 1979. viii, 140 pp., illus. Paper, \$18.20.

Geothermal Energy. Its Past, Present and Future Contributions to the Energy Needs of Man. H. Christopher, H. Armstead. Spon, London, and Halsted (Wiley), New York, 1979. xxviii, 358 pp., illus. \$27.50.

Geothermal World Directory. 1978/79 Edition. Geothermal World Corporation, Reseda, Calif., 1978. 524 pp., illus. \$45.

GLC and HPLC Determination of Therapeutic Agents. Part 2. Kiyoshi Tsuji, Ed. Dekker, New York, 1978. xiv + pp. 415-936, illus. \$45. Chromatographic Science, vol. 9.

Impact of the Federal Safe Drinking Water Act on Texas. Lyndon B. Johnson School of Public Affairs, University of Texas, Austin, 1978. xvi, 190 pp., illus. Paper, \$5. Policy Research Project Report, No. 21.

L'Induzione e l'Ordine dell'Universo. Alberto Meotti, Ed. Edizioni di Comunità, Milan, 1978. 274 pp. Paper, L 7500. Luoghi Critici, 3.

Infection, Immunity, and Genetics. Papers from a symposium, Philadelphia, June 1977. Herman Friedman, T. Juhani Linna, and James E. Prier, Eds. University Park Press, Baltimore, 1978. xiv, 186 pp., illus. \$16.50.

Infinitesimal Calculus. James M. Henle and Eugene M. Kleinberg. MIT Press, Cambridge, Mass., 1979. x, 136 pp., illus. \$12.50.

International Clay Conference 1978. Oxford, England, July 1978. M. M. Mortland and V. C. Farmer, Eds. Elsevier, New York,

Somatostatin 1-tyr [¹²⁵I]

Tested for binding and displacement
in a specific RIA
Carrier-free, monoiodinated
Fresh monthly

Somatostatin, 1-tyrosine [¹²⁵I]-, monoiodinated
Shipped frozen in 0.05M sodium phosphate
buffer (pH 7.5) containing a stabilizer and
proteolytic enzyme inhibitor.
NEX-129 10μCi 2x10μCi 50μCi 2x50μCi

Not for use in humans or clinical diagnosis.




New England Nuclear

549 Albany Street, Boston, Mass. 02118
Call toll-free: 800-225-1572
(In Massachusetts and International: 617-482-9595)

NEN Chemicals GmbH, Dreieich, W. Germany; NEN Canada Ltd., Lachine, Quebec

Circle No. 345 on Readers' Service Card

SPEED VAC **EVAPORATOR/CONCENTRATOR** *is a natural for Microcentrifuge Tubes.*




The **SPEED VAC** (Pat Pending)
EVAPORATOR/CONCENTRATOR is
now widely used with these
Microcentrifuge Tubes to minimize
sample loss and maximize sample
concentration. A combination of rotary
action and vacuum application account
for the high efficiency volume reduction
of the sample.

A wide variety of rotors is available to handle
a wide range of sample sizes from a few
microliters to ten milliliters and more.

Savant can supply the entire system
including the **SPEED VAC CONCENTRATOR**,
the **Refrigerated Trap** and the **Vacuum**
Pump as well as many accessories which
makes the laboratory completely operational
in a short time. The entire system does not
take up valuable laboratory space. Send for
literature today.

SAVANT INSTRUMENTS, INC.
221 Park Ave., Hicksville, N.Y. 11801
(516) 935-8774

Savant




The **SPEED VAC CONCENTRATOR** will yield virtually
100% sample recovery with
Eluates, Hydrolysates, Extracts,
Dialyzed solutions.

NO BUMPING **NO DILUTION**
NO SMEARING **NO TRANSFER**
NO WASHDOWN **NO LOSSES**

Excellent when preparing
samples by Derivatization,
GLC, HPLC, Amino Acid
Analysis, Electrophoresis and
Chromatography.

Circle No. 34 on Readers' Service Card

W. B. Saunders Announces a LANDMARK new edition!



Bondy & Rosenberg **METABOLIC CONTROL** **AND DISEASE** **8th Edition**

Why landmark? Because it is the only text which covers both human
metabolism and endocrinology in one integrated volume! The
thorough and totally up-to-date revision of this superb work is de-
signed to present the clearest and most current picture of the many
factors involved in producing disease. It is meticulously organized
and carefully written for fluidity and greater understanding.

Invaluable to the clinician or academician, **Bondy & Rosenberg**
contains clinically applicable information, a strong theoretical back-
ground, and abundant current and classic references—so it makes
as excellent an ongoing reference as it does a text. Shouldn't the
new 8th Edition of **Metabolic Control and Disease** be on your shelf?
Order now!

Edited by **Philip K. Bondy**, MD, Prof. of Medicine; Assoc. Chief of Staff for
Research and Development, West Haven Veterans Administration Hospital,
Yale University; and **Leon E. Rosenberg**, MD, Prof. and Chairman, Dept. of
Human Genetics, Prof. of Pediatrics and Medicine, Yale Univ. With contributors.
About 2000 pp., 825 ill. (About 4 pp. in color). About \$75.00. Ready soon.
Order #1844-8.

Please send on 30-day approval: SCIENCE 9/79
Bondy & Rosenberg **\$75.00**

Metabolic Control and Disease, 8th Edition
#1844-8

FULL NAME _____

POSITION & AFFILIATION (IF APPLICABLE) _____

ADDRESS _____

CITY _____

STATE _____

ZIP _____

☐ check enclosed
(send postpaid)

☐ bill me
(plus postage
and handling)

Prices differ outside US
and are subject to change

W.B. Saunders Co.

P.O. Box 207, Philadelphia, PA 19105

1 Goldthorne Avenue, Toronto, Ontario M8Z 5T9, Canada

1 St. Anne's Road, Eastbourne, East Sussex BN21 3UN, England

9 Waltham Street, Artarmon NSW 2064, Australia

Can you find your specialty in this list?

- Animal Surgery
- Pharmacology
- Electrophysiology
- Physiology Teaching
- Neurophysiology
- Physiological Psychology
- Respiratory Physiology

If you work in any of
these fields, you will want
the new Ealing BioSciences
Catalog.



The 216 page catalog
contains more than 500
products made by Ealing. To
receive your free copy call or
write your nearest Ealing
office today.



Pleasant St., So. Natick, MA 01760
(617) 655-7000 • Newport Beach, CA 92660
(714) 833-9826

International Offices in: Watford, England • Douai,
France • Höchst, W. Germany • Novara, Italy •
Montreal, Canada

1979. xiv, 662 pp., illus. \$55.50. Developments in Sedimentology, 27.

The Interpretation of Visual Motion. Shimon Ullman. MIT Press, Cambridge, Mass., 1979. xiv, 230 pp., illus. \$17.50. MIT Press Series in Artificial Intelligence.

Introduction to Animal Biology. Claude A. Villee, Warren F. Walker, Jr., and Robert D. Barnes. Saunders, Philadelphia, 1979. xvi, 616 pp., illus. + plates + appendixes. \$16.95.

Introduction to Educational Gerontology. Ronald H. Sherron and D. Barry Lumsden. Hemisphere, Washington, D.C., 1978. xii, 308 pp. Paper, \$14.95.

Introduction to Solid-State Theory. O. Madelung. Translated from the German by B. C. Taylor. Springer-Verlag, New York, 1978. xii, 488 pp., illus. \$28. Springer Series in Solid-State Sciences, vol. 2.

Mammalian Metabolism of Plant Xenobiotics. Ronald R. Scheline. Academic Press, New York, 1978. xii, 502 pp., illus. \$58.

Management of Low-Level Radioactive Waste. Papers from a symposium, May 1977. Melvin W. Carter, A. Alan Moghissi, and Bernd Kahn, Eds. Pergamon, New York, 1979. Two volumes. xx, 1214 pp., illus. \$125.

Manipulation of the Immune Response in Cancer. Proceedings of a workshop, Oxford, England, Sept. 1977. N. Avrion Mitchison and Maurice Landy, Eds. Academic Press, New York, 1978. xvi, 358 pp., illus. \$23.75. Perspectives in Immunology.

Manual of Mink Anatomy. An Illustrated Guide to Dissection. Allen A. Smith and Candice G. Krasulak. Burgess, Minneapolis, 1979. viii, 118 pp. Paper, \$7.95.

Mare Crisium. The View from Luna 24. Proceedings of a conference, Houston, Dec. 1977. Pergamon, New York, 1978. xiv, 710 pp., illus. + indexes. \$40. *Geochimica et Cosmochimica Acta*, Supplement 9.

Marine Natural Products. Chemical and Biological Perspectives. Vol. 2. Paul J. Scheuer, Ed. Academic Press, New York, 1978. xiv, 392 pp., illus. \$38.

Mass Spectrometry. Part A. Charles Merritt, Jr., and Charles N. McEwen, Eds. Dekker, New York, 1979. xiv, 284 pp., illus. \$35. Practical Spectroscopy, vol. 3.

Maximizing Treatment Gains. Transfer Enhancement in Psychotherapy. Arnold P. Goldstein and Frederick H. Kanfer, Eds. Academic Press, New York, 1979. xii, 488 pp., illus. \$25.

The Maximum Entropy Formalism. Papers from a conference, Cambridge, Mass., May 1978. Raphael D. Levine and Myron Tribus, Eds. MIT Press, Cambridge, Mass., 1979. xiv, 498 pp., illus. \$22.50.

Noise in Physical Systems. Proceedings of a conference, Bad Nauheim, Germany, Mar. 1978. Dietrich Wolf, Ed. Springer-Verlag, New York, 1978. x, 338 pp., illus. \$32.50. Springer Series in Electrophysics, vol. 2.

Nonmetallic Materials and Composites at Low Temperature. Proceedings of a symposium, Munich, July 1978. A. F. Clark, R. P. Reed, and G. Hartwig, Eds. Plenum, New York, 1979. xvi, 440 pp., illus. \$39.50. Cryogenic Materials Series.

Nutritional and Safety Aspects of Food Processing. Steven R. Tannenbaum, Ed. Dekker, New York, 1979. xii, 448 pp., illus. \$45. Food Service, vol. 6.

Ocean Wave Climate. Proceedings of a symposium, Herndon, Va., July 1977. Marshall D. Earle and Alexander Malahoff, Eds. Plenum, New York, 1979. xii, 368 pp., illus. \$36. Marine Science, vol. 8.

On the Theory and Practice of Voice Identifi-

cation. National Academy of Sciences, Washington, D.C., 1979. x, 162 pp. Paper, \$7.

Optical Fibers for Transmission. John E. Midwinter. Wiley, New York, 1979. xviii, 410 pp., illus. \$24.95. Wiley Series in Pure and Applied Optics.

Organization in Plants. W. M. M. Baron. Halsted (Wiley), New York, ed. 3, 1979. viii, 264 pp., illus. Paper, \$16.95.

Organometallic Mechanisms and Catalysis. The Role of Reactive Intermediates in Organic Processes. Jay K. Kochi. Academic Press, New York, 1978. xviii, 624 pp., illus. \$59.50.

Picosecond Phenomena. Proceedings of a conference, Hilton Head, S.C., May 1978. C. V. Shank, E. P. Ippen, and S. L. Shapiro, Eds. Springer-Verlag, New York, 1978. xii, 360 pp., illus. \$35.20. Springer Series in Chemical Physics, vol. 4.

Pioneer Conservationists of Western America. Peter Wild. Mountain Press, Missoula, Mont., 1979. xxvi, 246 pp., illus. \$12.95.

Plant Physiology. R. G. S. Bidwell. Macmillan, New York, and Collier Macmillan, London, ed. 2, 1979. xxii, 726 pp., illus. \$18.95.

The Plant Protection Discipline. Problems and Possible Developmental Strategies. Webster H. Sill, Jr. Allanheld, Osmun, Montclair, N.J., and Halsted (Wiley), New York, 1979. x, 190 pp. \$25.

Polymeric Delivery Systems. Papers from a symposium, Midland, Mich., Aug. 1976. Robert J. Kostelnik, Ed. Gordon and Breach, New York, 1978. viii, 312 pp., illus. \$39.80.

A Portrait of Isaac Newton. Frank E. Manuel. New Republic Books, Washington, D.C., 1979. xviii, 478 pp. + plates. Cloth, \$19.95; paper, \$7.95. Reprint of the 1968 edition.

Potential Industrial Carcinogens and Mutagens. Lawrence Fishbein. Elsevier, New York, 1979. x, 534 pp., illus. \$66.75. Studies in Environmental Science, vol. 4.

Practical Aspects of Memory. Proceedings of a conference, Cardiff, Wales, Sept. 1978. M. M. Gruneberg, P. E. Morris, and R. N. Sykes, Eds. Academic Press, New York, 1978. xiv, 786 pp., illus. \$38.50.

Preserving Our Natural Heritage. Vol. 2, State Activities. Prepared for the National Park Service by the Nature Conservancy, Arlington, Va., 1978 (available from the Superintendent of Documents, Washington, D.C.). vi, 672 pp. \$12.50.

Principles of Organic Synthesis. R. O. C. Norman. Chapman and Hall, London, and Halsted (Wiley), New York, ed. 2, 1979. xiv, 800 pp., illus. Paper, \$19.95.

Psychology and Women. In Transition. Jeanne E. Gullahorn. Winston, Washington, D.C., and Halsted (Wiley), New York, 1979. xiv, 224 pp., illus. \$15.95. Scripta Series in Personality and Social Psychology.

Psychotheatrics. The New Art of Self-Transformation. Robert D. Allen and Nina Krebs. Garland STPM Press, New York, 1979. xviii, 216 pp. \$17.50.

Public Health and Private Gain. The Economics of Licensing Clinical Laboratory Personnel. William D. White. Maaroufa Press, Chicago, 1979. x, 138 pp. \$14.95.

Rational Quadratic Forms. J. W. S. Cassels. Academic Press, New York, 1978. xvi, 414 pp. \$36.25. L.M.S. Monographs, 13.

Reflectors and Localization. Application to Sheaf Theory. F. Van Oystaeyen and A. Verschoren. Dekker, New York, 1979. viii, 166 pp. Paper, \$23.50. Lecture Notes in Pure and Applied Mathematics, vol. 41.

The Regulation of Psychotherapists. Daniel B. Hogan. Ballinger (Harper and Row), Cambridge, Mass., 1979. Four volumes. Vol. 1, A

solar energy
alternative energy technologies
food, energy, society
chemicals from biomass
toxins in the environment
engineering and ethics

world's arid regions and cultures
ethnography through film
community power studies
radical perspectives on economics
computers and society

risk-benefit analysis
patterns of problem solving

exploratory/robust statistics
numerical analysis
mathematical modeling
analysis of form

microcomputers
computers and learning science

cognition and teaching
refreshing calculus
women and mathematics
handicapped student in science
S.T.I.M.
energy and society

aging, family, bureaucracy
death, dying, ethics
holistic health
food and nutrition hassles

genetics and society
behavior programming
culture, body, behavior
behavior and evolution
rhythms
chronobiology

frontiers of neurosciences
psychobiology of pain
drug action
immunobiology
human diseases


flowering plants
animal communication
life in the oceans

how life began
early evolution of life
cosmology
theory of relativity
quarks to quasars
foundations of modern science

industrial organic chemistry
science of sound
holography

ANNOUNCEMENT ACADEMIC YEAR 1979-80

NSF Chautauqua-Type Short Courses for College Teachers



Conducted by
the American Association for the Advancement of Science
and the University of Missouri—Kansas City
and 14 Regional Field Centers
with support from the National Science Foundation

Office of Science Education
American Association for the Advancement of Science
1776 Massachusetts Avenue, N.W., D.C. 20036

For college teachers: Fifty-four short courses will be held throughout the United States. Each course is an intensive four days—two days in the fall (October-November) plus two days in the spring (February-March), with projects during the interim.

For nonacademic scientists and engineers: In the 1979-80 series, twelve short courses—at the University of Hartford and the Oregon Graduate Center—will be open to a mixed audience of college teachers and nonacademic professionals in science and engineering. Fall sessions are scheduled in October-November; spring follow-up sessions in February-March.

For further information (course descriptions, applications, etc.) please write or call:

American Association for the
Advancement of Science
Office of Science Education
1776 Massachusetts Avenue, N.W.
Washington, D.C. 20036
202/467-4465

RTE THE QUALITY REFRIGERATED BATH CIRCULATOR



Reliable temperature control of:

- Electrophoresis Chambers
- Calibration Procedures
- Spectrophotometer Cells
- Enzyme Assays
- External circulation with temperature control


Refrigerated bath circulators engineered for continuous operation within a temperature range of -30°C to +100°C. All stainless steel construction prevents corrosion. Rugged pump design. Low cost, 20-30% greater cooling capacity.



Call the leader — toll free
1-800-258-0830
In N.H. call collect 603-436-9444

NESLAB the name in circulation
NESLAB INSTRUMENTS, INC. PORTSMOUTH, N.H. 03801 U.S.A.

CRYOCOOL OPERATES CONTINUOUSLY. REPLACES DRY ICE & LN₂



- Range from -25°C to -100°C
- Flexible immersion probes — 3 styles
- CC II Series will control to -80°C with ±0.5°C stability

Refrigeration system which eliminates the hazards of handling and storing dry ice and liquid nitrogen. Save time and refrigerant costs.



Call the leader — toll free
1-800-258-0830
In N.H. call collect 603-436-9444

NESLAB the name in circulation
NESLAB INSTRUMENTS, INC. PORTSMOUTH, N.H. 03801 U.S.A.

EXACAL BATH CIRCULATORS UNMATCHED STABILITY



- 3-Point Sensing Probe
- Continuous industrial grade pump
- Range -35°C to +150°C
- Stability of ±0.01°C or better

Engineering achievements not available at twice the price. EXACAL Constant Temperature Baths perform trouble-free. Optional Z-Pump.



Call the leader — toll free
1-800-258-0830
In N.H. call collect 603-436-9444

NESLAB the name in circulation
NESLAB INSTRUMENTS, INC. PORTSMOUTH, N.H. 03801 U.S.A.

Unitron's zoomer

Model ZSB zoom stereo microscope

A big 6:5 to 1 ratio zoom optical system. Excellent resolution and contrast even at high powers. From 7X to 45X with 10X eyepieces. As low as 3.5X or as high as 180X with accessories. And there's a big choice of accessories . . . eyepieces, lens attachments, reticles, transmitted light base, illuminators, polarizing attachment and more.

Ask for free literature:
Unitron Instruments, Inc.,
101 Crossways Park West,
Woodbury, NY 11797, U.S.A.
Phone (516) 364-8046




UNITRON The value line

Subsidiary of Ehrenreich Photo-Optical Industries Inc.

Circle No. 338 on Readers' Service Card

DRINKING MONITOR WITH PRINTER SIMULTANEOUSLY MONITORS 8 ANIMALS



Multichannel Drinking Monitor works on principle of electrical conductivity of extremely small current (below 1 μ A) which is unnoticeable to animal touching water dispenser.

Unit can be installed in ordinary plastic or metal animal cages and with ordinary water dispenser. Monitor automatically prints out in preset time intervals (1 to 99 min.) number of animal licks on the dispenser. The same instrument can be used as printing multichannel event counter.



for more information call or write:
COLUMBUS INSTRUMENTS

950 NORTH HAGUE AVENUE
COLUMBUS, OHIO 43204 U.S.A.
PHONE (614) 488-6176

Circle No. 323 on Readers' Service Card

Looking for a grant?

GRANTS MAGAZINE

The Journal of Sponsored Research and Other Programs

will help you unravel the mysteries of obtaining one.

editor: Virginia White

assistant editor: Patrick Tarantino

Grants Magazine provides a forum for the lively exchange of information between those who seek grants and those charged with the responsibility for allocating funds and formulating programs through which grants are made available. Articles are selected for their relevance to public and private philanthropy from the philosophical, economic, social, and political viewpoint. *Grants Magazine* also publishes information about grant programs; guides to effective grantsmanship; listings of seminars, conferences, and workshops in the grants and development fields.

Grants Magazine is the recipient of the 1978 TSM Award for Outstanding Journal from the Technical, Scientific, and Medical Division of the Association of American Publishers.

Subscription: Volume 2, 1979 (4 issues) \$22.50

Institutional rate \$45.00

Please call for your sample copy: (212) 255-0713 ext 220

Prices slightly higher outside the U.S.



227 West 17th Street, New York, N.Y. 10011