

13 April 1984 • Vol. 224 • No. 4645

\$2.50

SCIENCE


AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE



John P. O'Neill - 1983

THE AIRFUGE® ULTRACENTRIFUGE:

*For high forces
and small volumes,
it's superb!*



The Beckman Airfuge® brings tabletop ease to ultracentrifugation. Compact and simple to use, it's the answer to separating microsample volumes in a hurry.

Rotors for Wide Range of Applications

The Airfuge can sediment any material that can be pellet-

ed in a full size ultracentrifuge. It's the convenient choice for molecular weight determinations, binding studies, cellular separations, lipoprotein fractionation and steroid hormone receptor assays. Fixed angle rotors hold 450- μ L, 240- μ L and 175- μ L tubes and a batch rotor holds volumes to 7 mL. The Electron Microscopy Rotor pellets particles directly onto EM grids, saving time and increasing sensitivity of particle detection 50 to 1000-fold.

Faster Separations


In the Airfuge, separations take only a fraction of the time

that full size ultracentrifuges require. It reaches top speeds and forces in just 30 seconds — to 110,000 rpm, 199,000 *g*.

New Refrigerator Kit for Increased Versatility

The Airfuge can be operated at temperatures as low as 12°C with optional new refrigerator accessories.

The Beckman Airfuge. For microsample volumes it's the ultracentrifuge to use. For details and a bibliography of 500 references, write: Beckman Instruments, Inc., Spinco Division, 1050 Page Mill Road, Palo Alto, CA 94304.



BECKMAN


If you hate waiting you'll love ENLIGHTNING™

**The new one-step
fluorography enhancer
that's ten times faster
than PPO-DMSO.**

ENLIGHTNING™ Rapid
Autoradiography Enhancer
takes only 15 to 30 minutes
to process your gels,
compared to five hours
for PPO-DMSO.

And in only one step:
simply soak the gels in
ENLIGHTNING™. No multi-
step procedure. Quick
processing minimizes
diffusion artifacts and loss
of smaller molecular weight
components, too.

The resulting fluorograms
are fully comparable to
those obtained using other
fluorographic techniques,
with either polyacrylamide
or 2% agarose gels (3 mm
min. thickness).



**And easier to live with
in the lab.**

ENLIGHTNING™ is odorless
and safer to work with than
enhancers containing DMSO.
It requires nothing more
than standard laboratory
gloves. What's more,
refrigeration is unnecessary;
ENLIGHTNING™ has a long
shelf life at room temperature.

Free 16-page Guide

Ask us to send you our
Guide to Autoradiography
Enhancement, with its
helpful tips on the
enhancement, exposure,
and development of
film—complete with a
trouble shooting checklist.
Free for the asking.

Not for use in humans or clinical diagnosis

New England Nuclear
549 Albany Street, Boston, MA 02118
Call toll free: 800-225-1572, Telex: 94-0996
Mass. and Internat'l: 617-482-9595
Europe: NEN Chemicals GmbH D-6072,
Dreieich, W. Germany, Postfach 401240,
Tel. (06103) 803-0, Telex 4-17993 NEN D
NEN Canada: 2453 46th Avenue,
Lachine, Que H8T 3C9
Tel. 514-636-4971, Telex 05-821808
© 1983 NEN

NEN New England Nuclear®



SCIENCE

LETTERS	U.S. Withdrawal from Unesco: <i>V. Rabinowitch</i> ; Juarez Radiation Accident: <i>L. Pagliaro</i> ; <i>E. Marshall</i> ; Nuclear Winter Scenario: <i>G. F. White</i> and <i>J. London</i> ; Illmensee Inquiry: <i>J. Gasser</i>	110
EDITORIAL	Amending the National Science Foundation Act: <i>F. Press</i>	115
ARTICLES	New Laser Source Technology: <i>C. P. Christensen</i>	117
	Fiber-Optic Sensors for Biomedical Applications: <i>J. I. Peterson</i> and <i>G. G. Vurek</i>	123
	Scientists and Congress: <i>H. B. Gonzalez</i>	127
NEWS AND COMMENT	Iraq's Chemical Warfare: Case Proved	130
	Pentagon Names Star Wars Czar	131
	Lawyers Flush Out Toxic Shock Data	132
	Universities Gag on Research Controls	134
	Will Deficits Put a Damper on R & D?	135
	Boom Time for British Biotechnology?	136
	Pressure for Patent Reform	137
	<i>Briefing</i> : Meselson Meets a Shower of Yellow Rain from Bees; Pediatrician May Head ADAMHA; In the Doghouse at OSHA	138
RESEARCH NEWS	Computing in the Language of Science	140
	Ice Cap of 30 Million Years Ago Detected	141

BOARD OF DIRECTORS

E. MARGARET BURBIDGE Retiring President, Chairman	ANNA J. HARRISON President	DAVID A. HAMBURG President-Elect	ROBERT W. BERLINER LAWRENCE BOGORAD	NANCIE L. GONZALEZ WALTER E. MASSEY
MATHEMATICS (A) Lipman Bers Lynn Arthur Steen	PHYSICS (B) James A. Krumhansl Rolf M. Sinclair	CHEMISTRY (C) Murray Goodman William L. Jolly	ASTRONOMY (D) Paul W. Hodge Donat G. Wentzel	
PSYCHOLOGY (J) Janet T. Spence Bert F. Green	SOCIAL, ECONOMIC, AND POLITICAL SCIENCES (K) Kenneth J. Arrow David L. Sills	HISTORY AND PHILOSOPHY OF SCIENCE (L) Daniel J. Kevles David L. Hull	ENGINEERING (M) Eric A. Walker W. Edward Lear	
EDUCATION (Q) Hans O. Andersen Roger G. Olstad	DENTISTRY (R) Erling Johansen Harold M. Fullmer	PHARMACEUTICAL SCIENCES (S) Stanley A. Kaplan David A. Knapp	INFORMATION, COMPUTING, AND COMMUNICATION Robert Lee Chartrand Madeline M. Henderson	

DIVISIONS

ARCTIC DIVISION	PACIFIC DIVISION	SOUTHWESTERN AND ROCKY MOUNTAIN I
John Davies President	Gunter E. Weller Executive Secretary	Barbara Wright President
		Alan E. Leviton Executive Director
		Charles E. Holley, Jr. President
		M. Michelle B. Executive Offi

SCIENCE is published weekly on Friday, except the last week in December, by the American Association for the Advancement of Science, 1515 Massachusetts Avenue, NW, Washington, D.C. 20005. Second-class postage (publication No. 484460) paid at Washington, D.C., and at an additional entry. Now combined with *The Scientific Monthly* © Copyright © 1984 by the American Association for the Advancement of Science. Domestic individual membership and subscription (51 issues): \$56. Domestic institutional subscription (51 issues): \$93. Foreign postage extra: Canada \$24, other (surface mail) \$27, air-surface via Amsterdam \$65. First class, airmail, school-year, and student rates on request. Single copies \$2.50 (\$3 by mail); back issues \$3 (\$3.50 by mail); Biotechnology issue, \$5 (\$5.50 by mail); classroom rates on request. **Change of address:** allow 6 weeks, giving old and new addresses and seven-digit account number. Authorization to photocopy material for internal or personal use under circumstances not falling within the fair use provisions of the Copyright Act is granted by AAAS to libraries and other users registered with the Copyright Clearance Center (CCC) Transactional Reporting Service, provided that the base fee of \$1 per copy plus \$0.10 per page is paid directly to CCC, 21 Congress Street, Salem, Massachusetts 01970. The identification code for *Science* is 0036-8075/83 \$1 + .10. **Postmaster:** Send Form 3579 to *Science*, 1515 Massachusetts Avenue, NW, Washington, D.C. 20005. *Science* is indexed in the *Reader's Guide to Periodical Literature* and in several specialized indexes.

AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE

	PET Scan Controversy Aired	143
	A New Route to Chiral Compounds	145
BOOK REVIEWS	Supernova Remnants and Their X-Ray Emission, <i>reviewed by S. P. Reynolds</i> ; Fundamental Processes in Energetic Atomic Collisions, <i>C.-D. Lin</i> ; Coated Grains, <i>J. D. Milliman</i> ; The Order of Man, <i>B. T. Shea</i> ; Books Received	146
REPORTS	Helium-3 Emission Related to Volcanic Activity: <i>Y. Sano et al.</i>	150
	Early Eocene Vertebrates from Baja California: Evidence for Intracontinental Age Correlations: <i>J. J. Flynn and M. J. Novacek</i>	151
	Topography from Single Radar Images: <i>R. L. Wildey</i>	153
	Gaseous Nitrate Radical: Possible Nighttime Atmospheric Sink for Biogenic Organic Compounds: <i>A. M. Winer, R. Atkinson, J. N. Pitts, Jr.</i>	156
	Cytomegalovirus Replicates in Differentiated but Not in Undifferentiated Human Embryonal Carcinoma Cells: <i>E. Gönczöl, P. W. Andrews, S. A. Plotkin</i>	159
	Single-Copy Inverted Repeats Associated with Regional Genetic Duplications in γ Fibrinogen and Immunoglobulin Genes: <i>A. J. Fornace, Jr., et al.</i>	161
	Phospholipase-Induced Crystallization of Channels in Mitochondrial Outer Membranes: <i>C. A. Mannella</i>	165
	Immunohistochemistry of Adenosine Deaminase: Implications for Adenosine Neurotransmission: <i>J. I. Nagy et al.</i>	166
	Sequence of the Human Somatostatin I Gene: <i>L.-P. Shen and W. J. Rutter</i>	168
	High Incidence of "Leapfrog" Pattern of Geographic Variation in Andean Birds: Implications for the Speciation Process: <i>J. V. Remsen, Jr.</i>	171
	<i>Technical Comments:</i> Arctic Biostratigraphic Heterochroneity: <i>D. V. Kent et al.</i> ; <i>G. Norris and A. D. Miall</i> ; <i>L. J. Hickey, R. M. West, M. R. Dawson</i>	173

Y NELKIN
SAWYER

SHEILA E. WIDNALL
HARRIET ZUCKERMAN

WILLIAM T. GOLDEN
Treasurer

WILLIAM D. CAREY
Executive Officer

Y AND GEOGRAPHY (E)
Merriam
s Dutro, Jr.

SCIENCES (N)
retchmer
E. Rhoads

CS (U)
Moses
Wegman

BIOLOGICAL SCIENCES (G)
Charlotte P. Mangum
Walter Chavin

AGRICULTURE (O)
Leo M. Walsh
Coyt T. Wilson

ATMOSPHERIC AND HYDROSPHERIC (W)
Hans A. Panofsky
Bernice Ackerman

ANTHROPOLOGY (H)
Richard A. Gould
Priscilla Reining


INDUSTRIAL SCIENCE (P)
Nat C. Robertson
Robert L. Stern

GENERAL (X)
Lora M. Shields
Rodney W. Nichols

COVER

Individuals from three populations of *Hemispingus superciliaris*, a small tanager of the Andes that illustrates the "leapfrog" pattern of geographic variation. Yellow populations of the northern and southern Andes are geographically separated by very differently colored gray populations in central Peru. This type of pattern is found in about 20 percent of all Andean birds with three or more subspecies. See page 171. [Painting by John P. O'Neill, Museum of Natural Science, Louisiana State University, Baton Rouge 70803]

**We want you to
know why we
put a handle on
our microbalances.**



In the minds of many, the word *microbalance* conjures up images of white gloves, dust-free chambers, and operators with the nerves and patience of a brain surgeon. The Cahn Series 20 Automatic Electrobalances put those images to rest once and for all.

They're rugged.

Our Series 20 balances are so rugged and compact that we even put a handle on them. You can take them where you need them *for microanalysis applications* such as *weighing expensive reagents, filters, hazardous materials, radioisotopes...* whenever you need to weigh small amounts of material.

They're sensitive and easy to use.

With Cahn microbalances, there are no internal weights or special set-up procedures. Push-button operation ensures *precise weighing to 0.1µg*—even for inexperienced users.

We've been making them for over 25 years.

And many of the first ones we built are still in use today. Which just goes to show, a microbalance doesn't have to be a delicate balance.

Contact Cahn for complete information and a demonstration. Our *toll-free number* is (800) 423-6641.

CAHN

Cahn Instruments, Inc.,
16207 South Carmenita Road, Cerritos, CA 90701.
Phone (213) 926-3378. TWX 910-583-4806.

Circle No. 86 on Readers' Service Card

INVITATION FOR PAPERS

**2nd U.S. - Dutch
International Symposium:
Aerosols**

**April 28 - May 3, 1985
Raleigh, NC**

As a part of cooperative activities in environmental protection under a Memorandum of Understanding between the United States and the Netherlands, an international symposium dealing with particulate matter air pollution is being organized. The first symposium in this series was conducted in Maastricht, The Netherlands in May 1982 addressing problems associated with nitrogen oxides pollution.

This second symposium will focus on issues related to air pollution by aerosols. The subtopic areas to be discussed are:

- Physical and Chemical Characterization
- Sampling, Measurements and Monitoring
- Sources, Emissions and Environmental Levels
- Atmospheric Transport and Transformation
- Biological Disposition of Inhaled Aerosols
- Toxicological Studies
- Controlled Human Exposure Studies
- Epidemiologic Studies
- Identification of Population Groups at Risk
- Ecological Effects of Acidic and Other Deposition
- Welfare Effects Assessment (Vegetation, Soiling, Materials Damage, Visibility, and Climate)
- Control Technology
- Benefit Analysis
- Economic Impact
- Control Policies
- Policy Implementation
- Other Relevant Issues


Abstracts: Persons wishing to make oral or poster presentations should submit a two-hundred word abstract by October 1, 1984, to:

**Dr. Si Duk Lee, U.S. Coordinator
2nd U.S.-Dutch
International Symposium: Aerosols
U.S. Environmental Protection Agency
(MD-52)
Research Triangle Park, NC 27711**

Papers based on oral presentations and selected poster presentations will appear in published conference proceedings.

OLYMPUS
The Image of Quality

Cambridge **Q10** Image Analyzer



The new Quantimet 10 Image Analyzer makes available for the first time, a high speed, automatic, sophisticated desk top analyzer at remarkably low cost. *In fact*, it offers speedy computerized image display and printout of feature measurements to those who regarded automated image analysis as unaffordable.

Teamed with the optically incomparable Olympus BHS microscope and TV camera, the Q10 can analyze images generated by macro equipment or any televised subject including photos and real objects. Output of the Q10 is displayed on its high resolution monitor in numerical terms or as distribution in bar charts, all printable with the Q10's high speed printer.

For a demonstration of this remarkable instrument, write to Olympus Corporation, P.A.I. Division, 4 Nevada Drive, Lake Success, NY 11042-1079 describing your applications. Or phone (516) 488-3880, Ext. 410. In the West: (415) 342-3384; Midwest: (913) 648-8323.

Cambridge **Q10**
Image Analyzer
from
OLYMPUS
Corporation

Circle No. 28 on Readers' Service Card

YOU'VE GOT A BRIGHT IDEA

for new biological or
chemical compounds.

WE HAVE THE CAPITAL to develop new drugs.

We are not a pharmaceutical firm but a Swiss-based, independant financial group which is prepared to let you participate in a capital venture for the exploitation of your project.

Our objective: the development of new biological or chemical products having previous pharmlological screening up to the stage of clinical tolerance and efficacy required to permit licensing negotiations.

To obtain all the information required for initial selection, please write in English, French or German to:
**DEBIOPHARM S.A. Petit-Chêne 38 -
1001 LAUSANNE (Switzerland).**

Circle No. 14 on Readers' Service Card

Announcing the new current-
awareness journal from Cambridge

BIO/ESSAYS

Editor-in-chief: **Dr W.J. Whelan** *University of Miami*

Cambridge University Press is pleased to announce this important new monthly which will be launched in July 1984 and published for the ICSU Press, the publishing house of the International Council of Scientific Unions.

BIO/ESSAYS meets the demand for a current-awareness journal devoted to reviews, comment, debate and progress in the frontier areas of molecular, cellular, and developmental biology and their impact on science and society.

While BIO/ESSAYS will concentrate primarily on advances in knowledge from basic research, it will also pay significant attention to actual or potential applications in medicine, industry and agriculture and is of particular interest to those working, or planning a career, in this field. It will debate speculative or controversial aspects of biological research, and disseminate authoritative views and exchanges of views.

CAMBRIDGE UNIVERSITY PRESS

32 East 57th Street, New York, N.Y. 10022

The Edinburgh Building Shaftesbury Road Cambridge CB2 2RU England

Circle No. 120 on Readers' Service Card

How To Buy An Incubator Without Laying An Egg

Learn how constant temperature happens. Precision™ incubators show you how. We've made more constant °C appliances for the lab than anyone else in the world.

**HOW
TO BUY AN
INCUBATOR
WITHOUT
LAYING AN
EGG**

FREE GUIDE on how to buy an incubator plus literature on general and special purpose, low temp, and CO₂ incubators. Write GCA/Precision Scientific Group, 3737 West Cortland Street, Chicago, Illinois 60647.

Precision™
Incubators

Digital temperature setpoint, 0.1°C precise. No rheostats, no guesswork.

LED continuously shows chamber °C. Never an average.

Solid state controller cycles heat near setpoint for closer uniformity, quicker recovery.

Platinum probe measures heat up top, in the best air flow.

Stainless steel chamber. Thick fiber-glass blanket. Epoxy coated finish.

Model 4L incubator, 37°C ± 0.4°C uniformity.



GCA CORPORATION
Precision Scientific Group

FOR TODAY'S FACULTY AND COLLEGE STAFF MEMBERS* FROM 18 TO 80.

Whether you're thinking retirement or not, review the plan that provides for cash withdrawal and/or a lifetime income.

TIAA-CREF Supplemental Retirement Annuities (SRAs) offer you substantial flexibility including cash withdrawal and/or lifetime retirement income. You can even reduce your income taxes now!

You can begin contributions to an SRA at any age and begin benefits at any age up to age 71 unless you are still employed (then you can delay beginning benefits until age 80). For example, you could start contributions at age 25, and choose to begin benefits or withdraw cash at age 34, 40 or 50, regardless of your employment status.

Get your money at any time.

You can receive benefits as a lifetime income or over a fixed period of from 2 to 10 years. What's more, if you need it (even while employed by your current employer), you can withdraw all the money you have accumulated by surrendering your contracts. Or, you can withdraw \$1,000 or more every six months. There is never a cash surrender charge.

Contributions are tax-deferred, so you pay less income taxes now.

The federal income tax on your contributions is deferred until they are paid to you as benefits. So, you pay less tax now.

Changing employers? Take SRAs with you.


Since you own your Supplemental Retirement Annuities, you take them with you if you leave your current employer. You can make contributions through any institution that makes Supplemental Retirement Annuities available to staff members. Contributions can be as little as \$25 a month.

Full information.

Complete and mail the coupon for an SRA Information Kit today. You'll get full details about all the advantages SRAs have to offer, why this plan suits so many financial situations and age groups and how much you may contribute to the plan.

*TIAA-CREF provides annuities and other services for employees of colleges, universities, private schools and certain other nonprofit tax-exempt educational and research institutions.

HELP YOURSELF TO A BRIGHTER FINANCIAL FUTURE ... SEND FOR A FREE INFORMATION KIT.



Teachers Insurance and
Annuity Association of
America—College Retirement
Equities Fund
730 Third Avenue
New York, N.Y. 10017

Please send me full details about TIAA-CREF Supplemental Retirement Annuities, the flexible tax-deferred annuity plan that offers the opportunity to accumulate funds for additional retirement income and the option for cash withdrawal.

Name _____

Address _____

City _____

State _____ Zip _____ Date of Birth _____

Name of Institution _____

Please let us know if you are participating in a TIAA-CREF retirement plan at your institution. ☐ YES ☐ NO

S41384

INTRODUCING EPHORTEC™
THE FIRST COMPLETE, INTEGRATED ELECTROPHORESIS SYSTEM...

ELECTROPHORESIS WITHOUT COMPROMISE

Now, there's no need to compromise in choosing an electrophoresis system for your laboratory. Just mention one name—EPHORTEC. And you've automatically included everything you could ever need for electrophoresis. It's the first and only complete electrophoresis line.

Ephortec offers all-new levels of consistency and quality—in the form of greater sensitivity, improved resolution and reduced running times. EPHORTEC's combination of instruments, accessories, reagents and technical expertise provide the user with maximum performance. They form one integrated system capable of meeting any application need.

Ephortec also means unprecedented convenience—all supplies and support are now available from a single source, easily obtained through local dealers.

Then, confidence. Ephortec is an already-proven technology, developed and backed by Haake Buchler experience and its tradition of highest quality instrumentation.


The Ephortec line of over 20 application-matched and system compatible instruments features:

DNA 100 cm Sequencer—Cuts typical running time in half. Highest resolution—over 200 nucleotides easily read in a single run. Specially designed thermoregulation allows overnight running, eliminates smearing; no circulator needed. Simple protocols, and uncomplicated gel preparation. Also available is a 45 cm model with the advantage of using standard size X-Ray films.

Ephortec Power Supplies

The Ephortec line features 6 power supplies from 250 to 5000 V—one to match every lab application. For example our 3000 V Power Supply—state-of-the-art for ultrathin, high voltage needs. Automatic crossover regulation in all three modes—constant voltage, current and power—assures optimal speed without overheating. Electronic timing from 15 min. to 18 hours; load detection and auto shutdown.

IEF Cell—Exceptional cooling, with flexibility for preparative (lengthwise) and analytical (horizontal) sample runs. Adjustable platinum electrodes for optimal spacing. Buffer tanks included.



Vertical Slab Cell—Exceptional cooling capacity, yielding running time under an hour for most procedures. Ease guaranteed through a unique clamping system for leak proof seals. Dual gel 40 sample capacity with optional combs and spacers.

Ephortec also makes available from one source the highest quality reagents and chemicals for all electrophoretic procedures.

Discover electrophoresis without compromise. To find out how easily and economically Ephortec can work for you, just call or write for our new free catalog with full details.

Ephortec™

Electrophoresis Instrumentation

Haake Buchler Instruments, Inc.,
244 Saddle River Road, P.O. Box 549,
Saddle Brook, New Jersey 07662
(201) 843-2320 (800) 631-1369

A HAAKEBUCHLER
Product Group




Circle No. 108 on Readers' Service Card

Ask us to solve your most difficult protein separation problems.


Hundreds of researchers have seen proof of how Waters new HPLC systems isolate and analyze proteins faster with higher resolution and recovery.

High speed purification of monoclonal antibodies.




High performance anion exchange chromatography successfully separates and purifies antibodies from cell culture and ascites fluid in less than 30 minutes.

High resolution mapping of cellular proteins.



Subtle differences in cell constituents can be detected by using high resolution ion exchange chromatography for the separation of crude cell lysates. (Courtesy of G.H.R. Rao et al.)

Total recovery of enzymatic activity from bacterial extracts.



The use of non-denaturing buffers with high speed isolation methods allows full activity recovery for many labile enzymes.

Let us show you how Waters HPLC systems can improve your protein separations.

Call ext. 2835 to arrange a technical presentation or a demonstration on your sample.



Waters Associates

MILLIPORE 34 Maple St., Milford, MA 01757 (617) 478-2000

Circle No. 46 on Readers' Service Card

Polaroid introduces instant 35mm slides.

Now you can use Polaroid's Instant Slide System with any 35mm camera or instrument with a standard 35mm back.



Introducing Polaroid's three new instant slide films: Polaroid Polachrome CS 35mm color transparency film (ISO 40/17). Polaroid Polapan CT 35mm continuous tone black and white transparency film (ISO 125/22). And Polaroid Polagraph HC 35mm high contrast black and white transparency film (ISO 400/27).

The new Polaroid 35mm Instant Slide System lets you develop, mount, and project these films in under 5 minutes. And at a cost that's about the same as conventional transparency films with processing.*

These Polaroid films are totally unlike any other 35mm films you've

ever used. And the technology is totally different. Polachrome CS film is based on an additive process; its base density is greater, and the ISO film speed is a little slower. Yet it produces unbelievably accurate, saturated colors with medium grain, sharp detail, and high resolution. The black and white films offer you the world's simplest, most direct way to get high quality, fine grain, continuous tone or high contrast slides. Instantly.

Processing is easy. You load the film into the AutoProcessor. Turn the handle. Wait a minute. And rewind the film into the cartridge. Mounting is simple, and you're ready to pro-

ject. All in under 5 minutes.

The scientific applications are virtually limitless.

Polaroid instant slide films can be used in hand-held 35mm cameras or in the most sophisticated laboratory instruments equipped with standard 35mm backs.

And because the Instant Slide System can give you slides in under five minutes, you can have scientific images ready on-the-spot. Securely, in-house. No darkroom is needed.

And you can be sure you have the results you need for documentation, presentation, or consultation.

You've waited long enough. It's available at dealers now.

See your local dealer now. Or find the name of the Polaroid Dealer nearest you by calling toll-free 800-526-7843 ext. 400 in the continental U.S.; 800-522-4503 ext. 400 in N. J.

*Comparison based on suggested list price of film and processing for 36 exposure roll.

The new 35mm Instant Slide System from POLAROID



"Polaroid", "Polapan", "Polachrome" and "Polagraph"

Circle No. 62 on Readers' Service Card

UV Fluorescence of CHO cells, Hoechst stain.
Nikon Diaphot microscope.
CF 20X fluor objective, 4 second exposure.
450X print magnification.

LIVING PROOF

For inverted microscopes, the proof is in the optical performance. And for years Nikon has outperformed the field. Now, the Diaphot sets even higher standards.

Utilizing Nikon's exclusive CF optics, it produces astonishingly bright, crisp, high-contrast images. Special Nikon objectives provide the unique range of 4X

to 100X phase contrast. A large stationary stage and moveable nosepiece system make micromanipulation sure and precise. Its wide range of accessories will allow the Diaphot to meet all the requirements of your laboratory. Best of all, the Diaphot's striking superiority is available at a very modest price.

Ask to look into a Diaphot. For full information and a reproduction of this photomicrograph, contact: Nikon Inc., Instrument Division, 623 Stewart Avenue, Garden City, NY 11530. (516) 222-0200.

Nikon
Extending Man's Vision

For information circle reader service number 102
For a demonstration circle reader service number 103

Nikon Diaphot



AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE

Science serves its readers as a forum for the presentation and discussion of important issues related to the advancement of science, including the presentation of minority or conflicting points of view, rather than by publishing only material on which a consensus has been reached. Accordingly, all articles published in *Science*—including editorials, news and comment, and book reviews—are signed and reflect the individual views of the authors and not official points of view adopted by the AAAS or the institutions with which the authors are affiliated.

Editorial Board

FREDERICK R. BLATTNER, BERNARD F. BURKE, ARNOLD DEMAIN, CHARLES L. DRAKE, ARTHUR F. FINDEIS, E. PETER GEIDUSCHEK, GLYNN ISAAC, NEAL E. MILLER, FREDERICK MOSTELLER, ALLEN NEWELL, RUTH PATRICK, BRYANT W. ROSSITER, VERA C. RUBIN, WILLIAM P. SLICHTER, SOLOMON H. SNYDER, PAUL E. WAGGONER, JOHN WOOD

Publisher: WILLIAM D. CAREY
Associate Publisher: ROBERT V. ORMES

Editor: PHILIP H. ABELSON

Editorial Staff

Assistant Managing Editor: JOHN E. RINGLE
Production Editor: ELLEN E. MURPHY
Business Manager: HANS NUSSBAUM
News Editor: BARBARA J. CULLITON
News and Comment: COLIN NORMAN (deputy editor), JEFFREY L. FOX, CONSTANCE HOLDEN, ELIOT MARSHALL, R. JEFFREY SMITH, MARJORIE SUN, JOHN WALSH

European Correspondent: DAVID DICKSON
Contributing Writer: LUTHER J. CARTER
Research News: ROGER LEWIN (deputy editor), RICHARD A. KERR, GINA KOLATA, JEAN L. MARK, THOMAS H. MAUGH II, ARTHUR L. ROBINSON, M. MITCHELL WALDROP

Administrative Assistant, News: SCHERRAINE MACK;
Editorial Assistant, News: FANNIE GROOM
Senior Editors: ELEANORE BUTZ, MARY DORFMAN, RUTH KULSTAD

Associate Editors: MARTHA COLLINS, SYLVIA EBERHART, CAITILIN GORDON, LOIS SCHMITT
Assistant Editors: STEPHEN KEPPEL, LISA MCCULLOUGH, EDITH MEYERS

Book Reviews: KATHERINE LIVINGSTON, *Editor*; LINDA HEISERMAN, JANET KEGG

Letters: CHRISTINE GILBERT
Copy Editor: ISABELLA BOULDIN
Production: JOHN BAKER; HOLLY BISHOP, ELEANOR WARNER; JEAN ROCKWOOD, SHARON RYAN, BEVERLY SHIELDS

Covers, Reprints, and Permissions: GRAYCE FINGER, *Editor*; GERALDINE CRUMP, CORRINE HARRIS

Guide to Scientific Instruments: RICHARD G. SOMMER
Editorial Administrator: SUSAN ELLIOTT
Assistant to the Associate Publisher: ROSE LOWERY
Assistant to the Managing Editor: NANCY HARTNAGEL
Membership Recruitment: GWENDOLYN HUDDLE

Member and Subscription Records: ANN RAGLAND
EDITORIAL CORRESPONDENCE: 1515 Massachusetts Avenue, NW, Washington, D.C. 20005. Area code 202. General Editorial Office, 467-4350; Book Reviews, 467-4367; Guide to Scientific Instruments, 467-4480; News and Comment, 467-4430; Reprints and Permissions, 467-4483; Research News, 467-4321. Cable: *Advances*, Washington. For "Information for Contributors," write to the editorial office or see page xi, *Science*, 30 March 1984.

BUSINESS CORRESPONDENCE: Area Code 202. Membership and Subscriptions: 467-4417.

Advertising Representatives

Director: EARL J. SCHERAGO
Production Manager: GINA REILLY
Advertising Sales Manager: RICHARD L. CHARLES
Marketing Manager: HERBERT L. BURKLUND
Sales: NEW YORK, N.Y. 10036: Steve Hamburger, 1515 Broadway (212-730-1050); SCOTCH PLAINS, N.J. 07076: C. Richard Callis, 12 Unami Lane (201-889-4873); CHICAGO, ILL. 60611: Jack Ryan, Room 2107, 919 N. Michigan Ave. (312-337-4973); BEVERLY HILLS, CALIF. 90211: Winn Nance, 111 N. La Cienega Blvd. (213-657-2772); SAN JOSE, CALIF. 95112: Bob Brindley, 310 S. 16 St. (408-998-4690); DORSET, VT. 05251: Fred W. Diefenbach, Kent Hill Rd. (802-867-5581).
ADVERTISING CORRESPONDENCE: Tenth floor, 1515 Broadway, New York 10036 (212-730-1050).

Amending the National Science Foundation Act

The intent of a proposal before Congress to amend the National Science Foundation Statutory Act of 1950 is to provide for a stronger emphasis on engineering in NSF programs, principally by inserting the words "engineering" or "engineers" in almost every place that "science" or "scientist" appears in the text. I fully share the concerns that prompted the suggested changes, and I have spoken to the importance of efforts to ensure our continuing technological strength, both nationally and internationally. Among the many problems we face are faculty shortages; inadequate, often outdated instruments and equipment; insufficient funding for engineering research; inefficient transfer of new scientific and engineering knowledge into public and commercial sectors; and inefficient transfer from industry to the universities of industrial knowledge and needs.

We do have to act. That agreed, there are, in terms of the proposed amendment, several questions before us. What can effectively be done within the existing statutory framework? Are new arrangements needed, and what are their likely effects? What processes should be followed for forming new arrangements and for implementing them?

First, what can be effectively done within existing arrangements? The Administration's proposed fiscal year 1985 budget answers the question. For example, there is a substantial increase for engineering research in the NSF budget—the largest proportionate increase for any of the directorates. In addition, support for engineering is embedded in the budgets of the mission agencies, with individual agencies providing much more support than the \$147 million proposed in fiscal 1985 for engineering research within NSF. In 1984, an estimated \$786 million was provided by the federal government for basic research in engineering, most of that coming from agencies other than NSF.

Let us assume that the new arrangements of the sort suggested in the proposed revisions are adopted. What are the likely outcomes? A major change may be to dilute the fundamental mission for which NSF was created—to support basic research in all the sciences. One needs to remember that, in contrast to engineering, many of the sciences are not related to any agency's mission except NSF. Is there any assurance that, as engineering research grows within NSF, there will be commensurate growth in its overall budget? If not, what impact will competition for funds have on the sciences that depend heavily on NSF support, such as mathematics and chemistry?

Although partitioning NSF between science and engineering is not the intent of the House Committee on Science and Technology, the wording used in the amendment can be given this unfortunate interpretation. Given the very few changes made to date in the 1950 act, can that act not, in its present form, provide for a larger role for engineering within NSF?

The science and engineering communities—academic, industrial, and governmental—should work together in a constructive way to address the problems of engineering. But one must ask whether legislative remedies will be truly effective in promoting this dialogue or may instead lead to corrosive frictions? Amendment of the NSF organic act is a serious undertaking whose consequences are uncertain. Will there be other requests to amend the act? Will we have a "discipline of the month" series of amendments, leading to the ultimate decay of one of the most successful institutions established by the federal government?


It is not clear that the changes are necessary. They may even be counterproductive. What is critical is to continue the dialogue on engineering that has been under way for several years. That dialogue has already produced substantial results. And, if conducted in a collegial manner with all parties involved, then we may fashion more effective remedies that will benefit not only engineering but also science and the federal role in research and development.—FRANK PRESS, *President, National Academy of Sciences, Washington, D.C. 20037*

Serum free cell growth with new ULTROSER G batch, after batch, after batch

ULTROSER G is a serum free synthetic growth support for multi-purpose cell culture—and cuts out the pretesting needed for fetal calf and other alternatives. As a true serum substitute, its well defined composition guarantees you batch uniformity to consistently give you regular cell development and reliable control over tissue culture variability. You'll soon find out how

ULTROSER G saves you time, is easy to use, simple to store and gets you results—batch after batch after batch. With 5 times less protein than FCS to simplify purification and bioassays, and the minimal amounts needed (1–3%), you'll know why so many labs like it even more. And when you discover that it comes lyophilized, ready for use and can be kept outside a freezer, you'll want to check it yourself. Get in touch with LKB, and we'll send you the latest on how to grow cells without serum—and without problems.

Trials performed with Vero cells in HAM'S F12 culture medium supplemented with 2% ULTROSER™ G. Tests carried out by Dr. E. Boschetti at Réactifs IBF, France.



Head office: LKB-Produkter AB, Box 305, S-161 26 Bromma, Sweden
Tel 08-98 00 40, telex 10492

Main US sales office: LKB Instruments, Inc. 9319 Gaither Road, Gaithersburg, Maryland 20877. Tel. 301-963-3200, telex 909870 (dom.), 64634 (intern.)

UK sales office: LKB Instruments Ltd., 232 Addington Road, S. Croydon, Surrey CR2 8YD, England. Tel 01-657 8822, telex 264414

Other sales offices in: Antwerp, Athens (for Middle East), Bangkok, Copenhagen, Lucerne, Madras, Moscow, Munich, Paris, Rome, Turku, Vienna, Zoetermeer

Developed and
manufactured by
IBF