


# SCIENCE

AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE



Trends in Analytical Instrumentation  
Multi-Step Chemistry Without Test Tubes



*The L8M:*

# FOUR REASONS WHY THIS IS THE BEST ULTRACENTRIFUGE YOU CAN BUY.

## Touch-Programmable Operation

With full microprocessor control, continuously monitored operation, digital readouts, built-in memory for delayed starting, and Memory-Pak™ modules, the L8M is the easiest-to-operate ultracentrifuge, ever.

## Quiet, Energy-Efficient Power

The L8M's advanced Ultra-Smooth direct induction drive is vacuum-encased with no high-speed vacuum seal to wear, no brushes to replace. Now even more powerful, the L8M acceler-

ates rotors up to 30% faster than ever before.

## Ultimate Versatility

You get  $\omega^2t$  computation, selectable accel/decel rates, the widest selection of rotors and tubes including exclusive Quick-Seal® Tubes, plus accessories for analytical capabilities, high tem-

perature work, and more. The L8M can even be controlled by your personal or lab computer.

## Reliable Performance

With user and service diagnostics, comprehensive safety systems, and a proven design based on 35 years of experience, the L8M is the most reliable, best-performing ultracentrifuge ever made.

Add our unmatched service and applications support, and it's no wonder that the L8M is the standard of performance, quality and value, worldwide. For more information contact your local Beckman Representative or write: Beckman Instruments, Inc., Spinco Division, 1050 Page Mill Road, Palo Alto, California 94304.



# BECKMAN

Circle No. 54 on Readers' Service Card



# COME TO THE SOURCE

We are the original and only source of true ATCC cultures. While many companies refer to an ATCC number, their cultures are derivatives.

The American Type Culture Collection, a non-profit organization, began acquiring, preserving and distributing microorganisms in 1925.

Today ATCC is *the* recommended source of test organisms cited by U.S. Pharmacopeia, Code of Federal Regulations, Title 21, Food and Drugs; AOAC Official Methods of Analysis; standards established by ASTM, ANSI, NCCLS, and more.

ATCC is the only "one-stop" source of a wide range of biologicals including: • Algae • Animal Viruses and Antisera • Bacteria • Cell Lines • Chlamydiae • rDNA Vectors • Fungi • Human Tumor Cells • Hybridomas • Molecularly Cloned Viruses • Mycoviruses • Oncogenes • Phages • Plant Viruses and Antisera • Plasmids • Protozoa • Rickettsiae • Tumor Immunology Cells • Yeasts

## ATCC Offers

**Reliability** — A culture purchased from the ATCC has been extensively characterized and quality controlled. ATCC's staff of experts make every effort to ensure that the culture you purchase is as unchanged as possible from the original deposit.

**Stability** — Strains are maintained in cryogenic storage rather than by continuous subculturing.

## Service

- Place orders toll-free. Dial 800-638-6597.
- Rapid turnaround — most orders *received* within one week.
- Unsatisfactory cultures promptly checked and replaced.
- Expert staff for technical support.

Reduce your risk of genetic drift, use thoroughly tested, authenticated cultures. Buy direct. ATCC — **the source** for quality cultures since 1925. Contact us for catalogues of available strains.



**American Type Culture Collection**

12301 Parklawn Drive, Rockville, MD 20852 USA  
(301) 881-2600

<b>LETTERS</b>	Space Commerce: <i>B. W. Ferguson</i> ; Stereochemistry: <i>K. Mislow and J. Siegel</i> ; DeBranges' Proof: <i>D. Drasin</i> ; Lamarckian Evolution: <i>J. G. Kaplan</i> ; On Making Comparisons: Reminded Again: <i>I. A. Greaves</i> ; Recombinant RNA Research: <i>R. J. Cedergren</i> . . . . .	240
<b>EDITORIAL</b>	Analytical Instruments. . . . .	249
<b>ARTICLES</b>	Trends in Analytical Instrumentation: <i>F. W. McLafferty</i> . . . . .	251
	Trends in Analytical Scale Separations: <i>J. W. Jorgenson</i> . . . . .	254
	Fourier Transform Mass Spectrometry: <i>M. L. Gross and D. L. Rempel</i> . . . . .	261
	Pyrolysis Mass Spectrometry of Complex Organic Materials: <i>H. L. C. Meuzelaar et al.</i> . . . . .	268
	High-Performance Liquid Chromatography–Mass Spectrometry: <i>J. L. Vestal</i> . . . . .	275
	High-Resolution Nuclear Magnetic Resonance of Solids: <i>G. E. Maciel</i> . . . . .	282
	Nuclear Magnetic Resonance Technology for Medical Studies: <i>T. F. Budinger</i> and <i>P. C. Lauterbur</i> . . . . .	288
	Laser Chemical Analysis: <i>R. N. Zare</i> . . . . .	298
	Contemporary Methodology for Protein Structure Determination: <i>M. W. Hunkapiller, J. E. Strickler, K. J. Wilson</i> . . . . .	304
	Chemical Sensing in Process Analysis: <i>T. Hirschfeld, J. B. Callis,</i> <i>B. R. Kowalski</i> . . . . .	312
<b>NEWS AND COMMENT</b>	A Cloudy Forecast for Solar Cells . . . . .	319
	Rifkin Takes Aim at USDA Animal Research . . . . .	321
	The Mystery of Florida's Citrus Canker . . . . .	322
	<i>Briefing:</i> NIMH Study Finds One in Five Have Disorders; WIPP in Good Shape Says Academy; Union Sets Up Legal Fund for Occupational Health	

## BOARD OF DIRECTORS

## CHAIRMEN AND SECRETARIES OF AAAS SECTIONS

ANNA J. HARRISON Retiring President, Chairman	DAVID A. HAMBURG President	GERARD PIEL President-Elect	ROBERT W. BERLINER LAWRENCE BOGORAD	WALTER E. MASSEY DOROTHY NELKIN
MATHEMATICS (A) Gail S. Young Lynn Arthur Steen	PHYSICS (B) Chen Ning Yang Rolf M. Sinclair	CHEMISTRY (C) Fred W. McLafferty Jean'ne M. Shreeve	ASTRONOMY (D) Patrick Palmer Donat G. Wentzel	
PSYCHOLOGY (J) Gregory A. Kimble William N. Dember	SOCIAL, ECONOMIC, AND POLITICAL SCIENCES (K) Robin M. Williams, Jr. David L. Sills	HISTORY AND PHILOSOPHY OF SCIENCE (L) Wesley C. Salmon David L. Hull	ENGINEERING (M) Raymond L. Bisplinghoff W. Edward Lear	
EDUCATION (Q) Marvin Druger Joseph D. Novak	DENTISTRY (R) Robert J. Fitzgerald Harold M. Fullmer	PHARMACEUTICAL SCIENCES (S) Stuart Feldman David A. Knapp	INFORMATION, COMPUTING, AND COMMUNICATION Joseph Becker Madeline M. Henderson	

## DIVISIONS

ARCTIC DIVISION	PACIFIC DIVISION	SOUTHWESTERN AND ROCKY MOUNTAIN DIVISION
John Davies President	Gunter E. Weller Executive Secretary	Barbara Wright President
		Alan E. Leviton Executive Director
		Charles E. Holley, Jr. President
		M. Michelle Bal Executive Director

SCIENCE is published weekly on Friday, except the last week in December, by the American Association for the Advancement of Science, 1515 Massachusetts Avenue, NW, Washington, D.C. 20005. Second-class postage (publication No. 484460) paid at Washington, D.C., and at an additional entry. Now combined with *The Scientific Monthly*. Copyright © 1984 by the American Association for the Advancement of Science. Domestic individual membership and subscription (51 issues): \$56. Domestic institutional subscription (51 issues): \$93. Foreign postage extra: Canada \$24, other (surface mail) \$27, air-surface via Amsterdam \$65. First class, airmail, school-year, and student rates on request. Single copies \$2.50 (\$3 by mail); back issues \$3 (\$3.50 by mail); Biotechnology issue, \$5 (\$5.50 by mail); classroom rates on request. Change of address: allow 6 weeks, giving old and new addresses and seven-digit account number. Authorization to photocopy material for internal or personal use under circumstances not falling within the fair use provisions of the Copyright Act is granted by AAAS to libraries and other users registered with the Copyright Clearance Center (CCC) Transactional Reporting Service, provided that the base fee of \$1 per copy plus \$0.10 per page is paid directly to CCC, 21 Congress Street, Salem, Massachusetts 01970. The identification code for Science is 0036-8075/83 \$1 + .10. Postmaster: Send Form 3579 to Science, 1515 Massachusetts Avenue, NW, Washington, D.C. 20005. Science is indexed in the *Reader's Guide to Periodical Literature* and in several specialized indexes.



# AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE

	Test-Maker Sues Stanford; Comings and Goings; DOD Buys Secure Telephones.....	324
<b>RESEARCH NEWS</b>	Climate Since the Ice Age Began to Melt .....	326
	Surprise Finding with Insect Globin Genes.....	328
	Femtosecond Laser Annealing of Silicon .....	329
<b>BOOK REVIEWS</b>	From the American System to Mass Production, 1800-1932, <i>reviewed by J. Atack</i> ; The Chemical Evolution of the Atmosphere and Oceans, <i>J. F. Kasting</i> ; The Science of Musical Sound, <i>R. N. Shepard</i> ; Books Received .....	331
<b>REPORTS</b>	Spread of <i>Diadema</i> Mass Mortality Through the Caribbean: <i>H. A. Lessios</i> , <i>D. R. Robertson, J. D. Cubitt</i> .....	335
	The Jovian Nebula: A Post-Voyager Perspective: <i>J. T. Trauger</i> .....	337
	Photocontraction of Liquid Spiropyran-Merocyanine Films: <i>I. Cabrera et al.</i> ....	341
	Temporal Variability of the Antarctic Circumpolar Current Observed from Satellite Altimetry: <i>L.-L. Fu and D. B. Chelton</i> .....	343
	Extensive Volcanism Associated with the Separation of Australia and Antarctica: <i>J. B. Jones and M. J. Fitzgerald</i> .....	346
	Chromosomal Location of Human T-Cell Receptor Gene $T_i\beta$ : <i>P. E. Barker et al.</i> .....	348
	Frequency-Dependent Noradrenergic Modulation of Long-Term Potentiation in the Hippocampus: <i>W. F. Hopkins and D. Johnston</i> .....	350
	<i>Technical Comments</i> : Uncertainty of Histologic Classification of Experimental Tumors: <i>S. J. Galli; A. Konstantinidis, J. B. Smulow, C. Sonnenschein</i> ; Impact Event at the Cretaceous-Tertiary Boundary: A Possible Site: <i>B. M. French</i> .....	352


## COVER

Sequence of mass spectra, obtained using Fourier transform mass spectrometry (FTMS), shows the in situ synthesis of the bare dimer ion  $\text{CoFe}^+$  and its subsequent reactions with 1-pentene in the gas phase. Recent developments in FTMS, involving both software and hardware, now enable selected ions of interest to be stored on the order of seconds during which time they can be mixed with various reagents and their chemistries monitored through a series of more than eight reaction steps (all under computer control). See page 261. [Research by D. B. Jacobson and B. S. Freiser, Department of Chemistry, Purdue University, West Lafayette, Indiana 47907; illustration by J. Christenson and R. Cody, Nicolet Instrument Corp., Madison Wisconsin 53711]


B. SLAUGHTER E. SAWYER	SHEILA E. WIDNALL LINDA S. WILSON	WILLIAM T. GOLDEN Treasurer	WILLIAM D. CAREY Executive Officer
OGY AND GEOGRAPHY (E) W. Hay mas Dutro, Jr.	BIOLOGICAL SCIENCES (G) Dorothy M. Skinner Walter Chavin	ANTHROPOLOGY (H) Priscilla Reining	
AL SCIENCES (N) A. Good an E. Rhoads	AGRICULTURE (O) John Pesek Ralph J. McCracken	INDUSTRIAL SCIENCE (P) J. Kenneth Craver Robert L. Stern	
STICS (U) a A. Bailar J. J. Wegman	ATMOSPHERIC AND HYDROSPHERIC (W) William W. Kellogg Bernice Ackerman	GENERAL (X) George C. Sponsler Rodney W. Nichols	

American Association for the Advancement of Science was founded in 1848 and incorporated in 1874. Its objects further the work of scientists, to facilitate cooperation among them, to foster scientific freedom and responsibility, to promote the effectiveness of science in the promotion of human welfare, and to increase public understanding and appreciation of the importance and promise of the methods of science in human progress.





ASYST multiple windows permit side-by-side comparisons. The two-dimensional FFT routine, shown, is one of the many built-in functions.



With ASYST, data acquisition and analysis routines can be combined to create powerful, menu-driven functions.

# Now acquire, analyze, and graph data all at the same time.


## Introducing ASYST™ Scientific Software.

**ASYST. The most powerful scientific software package available for a personal computer.**

ASYST is the only software that turns an IBM PC, or compatible, into a complete scientific workstation—offering full data acquisition, analysis, and graphics capabilities.

- ASYST puts you in direct control of your data. Its unique interactiveness allows you to start reducing data even as it is being collected—graphing as you proceed.
- Multiple graphics windows permit quick visual comparison and allow you to select segments for further manipulation.

ASYST gives you the speed and precision of a minicomputer—at a fraction of the cost. And because it works on a PC, you




ASYST and your IBM PC give you the power, speed and precision of a minicomputer in your own lab.


can pull information through the analytical process without timesharing delays or software bottlenecks.

- ASYST fully utilizes an IBM PC with the 8087 coprocessor. This unique combination takes all intermediate calculations to the 80-bit precision level, offering precision exceeding that of other micros—as well as many minis and mainframes.
- A 1024-point Fast Fourier Transform, with ASYST and an IBM PC with 8087 chip, takes less than 3 seconds. An optimum performance custom routine tested on a DEC 11/23+ minicomputer with FPF 11™ took 2 seconds—at five times the price.





ASYST automatic graphics are flexible enough to handle the most demanding applications.



ASYST integrates analysis functions with graphics. On-screen cursors allow interactive selection of curve segments.

### Three integrated modules.

ASYST is composed of three separate, but fully-integrated, modules—Module 1: System/Graphics/Statistics; Module 2: Analysis; and Module 3: Data Acquisition.

- Module 1 works alone—or with either of the other two modules—allowing you to tailor the system to your specific needs.
- All functions are always available for your use. You never have to leave the system to access any of its capabilities.

**Module 1: System/Graphics/Statistics** establishes the environment. Among other features, it provides data representation and storage capability; supplies arithmetic, trigonometric, hyperbolic and other mathematical and statistical functions; and provides direct graphics output and display. Array manipulation, control of vectors and matrices, automatic plotting, file manipulation, programming control structures, and a built-in text editor are all included in Module 1.

**Module 2: Analysis** reduces and analyzes data and includes a powerful selection of analytical functions. They include Eigenvalues, Eigenvectors, and polynomials. Least squares approximations, curve fitting, convolutions, integration, differentiation, smoothing, and Fast

Fourier Transform are only a few of the automatic functions provided.

**Module 3: Data Acquisition** allows ASYST to interface with laboratory instruments to capture data directly with a minimum of keyboarding. Standard commands such as "A/D.IN" are all that are needed to carry out communication between ASYST and standard interface boards. A/D and D/A conversions, digital I/O, timing, and triggering are all supported. Commands can be combined in programs to provide customized automatic acquisition and control.

### Built-in routines make ASYST easy to use. Programmability lets you customize.

You don't have to be a computer expert to utilize ASYST. Instead of rigid, unreadable computer syntax, you deal with an assortment of built-in commands that do what their names indicate—XY.DATA.PLOT, ARRAY.EDIT, FFT, etc.

- Interactively use pre-programmed commands to get immediate control of your work.
- Combine and modify commands to extend the capability of the system for your custom applications. ASYST is fully programmable.

### Supported by extensive, top-down documentation.

ASYST layered documentation—with examples, a quick reference card, and a cross-referenced index—lets you go only as deep as your needs require. A comprehensive, on-line "help" system is always ready to supply additional assistance.

**30-Day  
No-Risk Offer.**  
**1-800-348-0033**

In New York State call (212) 702-3241.

Call our ASYST hotline for more information.




**ASYST**

MACMILLAN SOFTWARE CO.  
An Affiliate of  
Macmillan Publishing Company  
866 Third Avenue  
New York, NY 10022

Circle No. 75 on Readers' Service Card





# WE CONFIGURE IT OUT

## Whatever your application requirements, our FPS-5000 Series of Array Processors offer a power-memory-I/O configuration to match.

No two application environments are exactly the same. No single array processing system can meet every need.

That's why, unlike other systems on the market, the FPS-5000 Series from Floating Point Systems gives you a choice:

The FPS-5000 family, with 16 primary configurations, offers a performance range from 8 to 62 million floating-point operations per second; data memory from 0.5M to 1M words; program memory to 32K words; and a variety of I/O options.

By taking this modular approach to system architecture, we're able to work with you in determining the optimum

configuration to match your particular performance requirements. Should your requirements change, cost-efficient field-installable upgrades assure continued compatibility.

## Distributed Processing Architecture

The unequalled configurability and cost efficiency (less than \$2,000 per MFLOP)\* of the FPS-5000 Series is achieved by combining a distributed processing architecture with the latest VLSI technology.

Throughput is maximized by allocating the computational load to a set of independent, self-contained 18-MFLOP-capacity processing units called Arithmetic Coprocessors. I/O processors can be added to interface with bulk storage and real-time devices. Data flow is managed by both

the central control processor and by independent I/O processors.

With distributed processing architecture, you select the power-memory-I/O configuration that is tailored specifically to your needs.

## We work with you

Floating Point Systems pioneered the array processor in 1975, and has led the industry ever since, with product developments and enhancements designed to optimize each customer's system and protect their investment.

The FPS-5000 Series is no exception. It provides a software migration path for our previous 38-bit processors, interfaces with a wide range of host computers, and offers the same impeccable quality, reliability and worldwide support services that have always set Floating Point Systems products above the rest.

For more information, call toll free, 1-800-547-1445.

### FPS-5000 SERIES

### Rated MFLOP

Product Group	8	12	26	30	44	48	62
5100	5105	—	5110	—	—	—	—
5200	—	5205	—	5210	—	—	—
5300	—	—	—	5310	—	5320	—
5400	—	—	5410	—	5420	—	5430
5500	—	5505	—	5510	—	—	—
5600	—	—	—	5610	—	5620	—

\*Based on U.S. Domestic Prices, Model FPS-5430.

Circle No. 112 on Readers' Service Card



**FLOATING POINT  
SYSTEMS, INC.**

P.O. Box 23489  
Portland, OR 97223 (503) 641-3151  
TLX: 360470 FLOATPOIN BEAV



## The Case of the Captive Investment

W

illiam and Judith Trembath needed some of the extra money they were saving for retirement. "No!" they were told, "You haven't met the conditions for its release." "Conditions?" they asked, "What conditions? We didn't realize that there were any conditions." If only, they thought, there were a way to save extra money for retirement without tying it up.

Finding a way to gain access to your extra savings in a tax-deferred retirement account—without a federal penalty tax—needn't be a mystery. There is a simple solution. Supplemental Retirement Annuities from TIAA-CREF. Or, very simply, SRAs. The flexible tax-deferred annuity plan that can increase your personal retirement savings and reduce your taxes while still giving you access to your money when needed.

Contributions to SRAs are made through your institution's payroll system on a *before-tax* basis. You don't pay any federal taxes on that money, or on the interest and dividends it earns, until you decide to withdraw your funds. And you can withdraw cash from your SRAs even *before* retirement—without paying a federal penalty tax.

Should you ever decide to change jobs, you can take your SRAs with you. You can continue contributions through any of

the 2,700 institutions in the country that make SRAs available.

And you're not locked into receiving your income benefits at retirement. You can begin benefits *before* you retire without paying a federal penalty tax. And you can choose from a number of *lifetime* income options for you and your spouse.

You can enjoy all this freedom and still have the benefit of a low-risk, high-security approach to investing. Funds in TIAA grow at the current effective annual interest rate of 11.50%\* and those in CREF share in the investment experience of the Fund's broadly diversified common stock portfolio. And TIAA-CREF have over \$30 billion in assets and 65 years of experience in retirement plan investing for people in education.

So if you want to have access to your extra savings in a tax-deferred retirement account, investigate the advantages of Supplemental Retirement Annuities.



**SRAs—Simple. Reliable. Accessible. Tax-deferred annuities. From TIAA-CREF.**

\*Effective March 1, 1984 through February 28, 1985 for funds credited January 1, 1982 or later. This interest rate is not guaranteed after February 28, 1985. A 1.5% expense charge is deducted from SRA premiums. If you would like to know more about SRAs, simply mail the coupon below or write TIAA-CREF, 730 Third Avenue, New York, NY 10017.



SCI 101984 2

TIAA-CREF  
730 Third Avenue  
New York, NY 10017

**Yes!** I'd like to increase my extra retirement savings and still have access to my money. Please send me your free brochure on solving this and other mysteries of my financial future.

Name \_\_\_\_\_

Name of Institution \_\_\_\_\_

Address \_\_\_\_\_

City \_\_\_\_\_

State \_\_\_\_\_

Zip \_\_\_\_\_



# THE FOURTH ANNUAL CONGRESS FOR HYBRIDOMA RESEARCH

FEBRUARY 3-6

SAN FRANCISCO HILTON HOTEL, SAN FRANCISCO, CALIFORNIA

Zenon Steplewski  
The Wistar Institute  
Philadelphia, PA

*Co-Chairmen:*  
Hilary Koprowski  
The Wistar Institute  
Philadelphia, PA

Joseph Davie  
Washington University  
St. Louis, MO

The Annual Congress for Hybridoma Research, now in its fourth year, is recognized as the premier symposium on the subject, and presents the important developments that impact on current and future work.

## **Opening Remarks: CONGRESS CHAIRMAN**

## **Workshop Topics & Chairmen:**

MONOCLONAL ANTIBODIES AND LYSIS OF TUMOR CELLS

Dolf Adams, Duke University

IMMUNOREGULATION

Fritz Melchers, Basel Institute for Immunology

IDIOTYPES

Hilary Koprowski, The Wistar Institute

GENETIC ENGINEERING OF MONOCLONAL MOLECULES

Donald Capra, University of Texas Southwestern Medical School

DIAGNOSTIC APPROACHES

Noel Warner, Becton-Dickinson Company

WORKSHOP VI - TO BE ANNOUNCED

## **Poster sessions: TECHNOLOGICAL ADVANCES IN HYBRIDOMA RESEARCH**

**Participants are invited to submit abstracts for the poster sessions. These abstracts will be reviewed up until the time of the meeting; however, only those accepted by Dec. 15 will be published in the journal, Hybridoma. Contact Dr. Zenon Steplewski (215) 898-3924.**

**Organized by Scherago Associates, Inc., in conjunction with  
Mary Ann Liebert, Inc.**

## REGISTRATION FEES:

\$400 On-site registration — includes a one year subscription or renewal to the journal, Hybridoma.

\$350 ADVANCE REGISTRATION — (Received by Dec. 1) — Includes a one year subscription or renewal to the journal, Hybridoma.

\$175 STUDENT REGISTRATION — Student status must be confirmed in writing by department chairman.

4-7 registrations received together from same organization \$300 each. Includes 4 journal subscriptions only.

8-10 registrations received together from same organization \$200 each. Includes 4 journal subscriptions only.

Larger group rates available upon request.

Attendance will be limited. Make checks payable to: **Scherago Assoc., Inc., DNA/HYBRIDOMA**

☐ Please reserve \_\_\_\_\_ space(s): Registration Fee of \$ \_\_\_\_\_ enclosed.

☐ Please send abstract form.

Name \_\_\_\_\_

Dept. \_\_\_\_\_

Organization \_\_\_\_\_

Street \_\_\_\_\_

City \_\_\_\_\_ State \_\_\_\_\_ Zip \_\_\_\_\_

Telephone: ( ) \_\_\_\_\_

Return to: **Hybridoma; c/o Scherago Associates, Inc.**  
1515 Broadway, New York, NY 10036 • (212) 730-1050



# MicroGenie: a powerful, easy-to-use sequence analysis program.

Now there's a program for the IBM PC or XT microcomputer that's fast enough to find every site for 100 restriction enzymes on pBR322 in two minutes, and potent enough to analyze sequences up to 60,000 base pairs in length.



It's called MicroGenie, and was developed by Drs. Korn and Queen. Users call it the most versatile and convenient program ever developed for analyzing DNA, RNA and proteins. And you don't need to know anything about computers to use it: MicroGenie contains built-in menus and help messages to guide you at every step.

## MicroGenie gives you the entire National Nucleic Acid Sequence Data Bank at your fingertips.

MicroGenie has all the features needed to collect, analyze and display sequence data. It comes with the more than 3000 sequences in "GenBank™" and can search the whole data bank for sequences homologous to yours. A unique sequence editor helps you enter and change sequences quickly. A sonic digitizer can be used to enter sequences directly from autoradiograms.

MicroGenie determines residue and codon frequencies, translates and reverse translates, locates open reading frames, merges sequences, even predicts RNA and protein secondary structures. And more.

Convince yourself! Send for a free demonstration diskette to sample what MicroGenie can do for your research. Write on your letterhead to Beckman Instruments, Inc., Spinco Division, 1050 Page Mill Road, Palo Alto, CA 94304; or phone 800-DNA-ACGT.



MicroGenie is a Trademark of SciSoft, Inc., and the MicroGenie sequence analysis program is licensed for exclusive distribution by Beckman Instruments, Inc.

Genbank is a Trademark NIH.

# BECKMAN

Circle No. 56 on Readers' Service Card




## The Scientist's Selection

The best microcomputer scientific software for the laboratory in one catalog. Selected by scientists for scientists . . . from Scientific Software Products, Inc.

After an extensive survey of scientists and engineers, we've compiled a unique single-source catalog for Apple II and IBM PC software specifically for the laboratory.

Our own scientists have evaluated and selected the very best software for **computer graphics, scientific word processing, research data management, statistics**



**and more.** We've also added difficult-to-find hardware items, accessories and books on scientific computing.

This second edition of our catalog emphasizes scientific graphics and includes software for: analytical and presentation graphics, creative drawing, computer-aided design (CAD), contour and 3-D plotting, and a variety of graphics utilities.


### SCIENTIFIC SOFTWARE PRODUCTS, INC.

For a FREE copy of our catalog, write or call: 5726 Professional Circle, Suite 105  
Indianapolis, IN 46241  
(317) 244-6163

Graphic design courtesy of Golden Graphics, Inc. Apple® II is a registered trademark of Apple Computer, Inc. IBM PC is a registered trademark of International Business Machines Corporation.

Circle No. 28 on Readers' Service Card

### **Microtex Scientific Imaging Lets Your Computer See The Light**



An integrated scientific imaging workstation from Microtex can solve your image acquisition problems, and provide results in minutes instead of days.

Whether your application is high speed, low light, UV or IR, a Microtex system can capture, analyze and display almost any image an experiment can create.

For expert applications engineering and further information, contact:

Microtex Corporation  
80 Trowbridge Street  
Cambridge, MA 02138  
617-491-2874 Telex 948-014

**MICROTEX**  
DIGITAL IMAGING

Circle No. 171 on Readers' Service Card

## AAAS Report IX

## Research & Development

## FY 1985

## Intersociety Working Group

This timely document analyzes major budget and policy issues relating to R&D in the FY 1985 budget, presents data on federal agency and industry support for R&D, and discusses trends in R&D funding in light of current policy issues.

**284 pp.**  
**Paperback \$10.00**

*To order, please write:*

**AAAS Sales Dept.,  
1515 Massachusetts  
Avenue, NW,  
Washington, DC 20005.**

Please allow 6-8 weeks for delivery. All orders under \$10 must be prepaid.

Visa and MasterCard customers include account number, expiration date, and signature.

**American Association for the  
Advancement of Science**




# We cut 64-bit computation down to size and made it practical...

Mainframe computational speed in a compact package about the size of a microwave oven. That's the new MAP-6420 array processor. Now computation for scientific and engineering tasks doesn't demand three refrigerators' worth of real estate. Now it doesn't require special air-conditioned rooms. And now, for the economical price of \$100,000, you can boost the speed of your host computer into the mainframe class.

**Practical Number-Crunching:** The basic MAP-6420 configuration with 4 mega-



bytes of data memory is packed into a rack-mountable chassis that stands only 21" high. Data memory can be expanded to 64 megabytes to handle even the most demanding applications. MAP was designed with a low



component count that results in high reliability. And low component count offers the benefits of low power, low cooling, and minimal space requirements. In dollar outlay—from purchase price, to installation, to operating costs—MAP-6420 is extremely economical.

**FORTRAN Programmability:** If you know how to program in FORTRAN, you already know how to get high performance from the MAP-6420. We planned MAP's hardware around its FORTRAN compiler so that programs you write will

deliver the performance you paid for. Parallel units process algorithms efficiently, even if they're not vector-oriented. Uniform, versatile performance delivered from FORTRAN—that's MAP-6420.

Call CSPI today to arrange a MAP-6420 demonstration or to get technical specs and price information. You know about the alternatives, now get to know MAP-6420. Discover a great investment in 64-bit computation. **Call toll free 1 800 325-3110.**

SEE US AT SCICOMP—BOOTH #306.



## CSPI

THE ARRAY PROCESSORS

40 Linnell Circle, Billerica, Massachusetts 01821 • 617/272-6020 • TWX: 710-347-0176

Circle No. 123 on Readers' Service Card





# All microscopes should be created upright... even inverted microscopes!


The new Leitz Labovert revolutionizes inverted microscopy. Its unique upright design is an exceptional innovation offering unprecedented convenience and ease of operation in an inverted microscope.

Users no longer have to reach around the microscope's body or up in the air as with conventional inverted microscopes. The controls for focussing and stage travel are exactly where human anatomy dictates they should be. The stage is in full view and easily accessible for loading. Just like an upright microscope.

The Labovert scores many other firsts. The lamp housing, an annoying source of stray light and heat, has been removed from the working space. High performance optics (field of view up to 22.5mm), condensers with long working distance (up to 95mm) and a formidable array of peripheral equipment adapt the Labovert to almost any application.

Innovators expect their ideas to be copied — but you can't improve on a true Leitz original. If you need an inverted microscope, you will want to compare the Labovert with all others.

For more information or a demonstration of the Labovert call or write E. Leitz Inc., Rockleigh, NJ 07647, (201) 767-1100.



*The Reinheimer Principle  
Illuminating and Image Forming  
Beams Cross at a Point in Space*

## Leitz Innovates Others Imitate

**Leitz Means  
Precision.  
World Wide.**



R-3013

Circle No. 81 on Readers' Service Card





# For the small lab, it's often sink or swim. Here's help from Perkin-Elmer.

## **First, we assessed your needs.**

The smaller lab is often a hectic environment. You have oceans of paperwork. And workloads are heavier, with no budget to increase your workforce. With our LIMS (Laboratory Information Management System) for the small lab, you can: provide rapid response to submitter inquiries; produce backlog reports for scheduling and workload analysis; obtain information to construct reports on turnaround time, productivity, and quality; make workload projections; reduce your paperwork and provide automatic report generation. Just for starters.

## **Then, we made an affordable solution.**

With recent advances in computer

architecture and design, we created several LIMS using the Perkin-Elmer 3205 minicomputer. This full 32-bit mini has hardware floating point capabilities. You can add additional memory, disk capacity and hardware peripherals. And best of all, the Perkin-Elmer 3205 is very economical.

## **With great flexibility.**

The 3205 can support as many as 16 channels of simultaneous data collection and four concurrent terminal users in the C/LAS (Chromatography Laboratory Automation System) environment. An integrated LIMS/2000 system which incorporates both information management and chromatographic data handling is also available on the 3205.

## **And all the support you need.**

Our LIMS for the smaller lab is genuinely easy to use. And to expand, if the need arises. It's also very easy to begin. We will work closely with you to choose exactly what you need. For complete information, talk to your Perkin-Elmer representative today. For literature, please call 1-800-762-4000.

Perkin-Elmer Corp., Main Ave. (MS-12),  
Norwalk, CT 06856 U.S.A. Tel: (203) 762-1000.  
Telex 965-954.

Bodenseewerk Perkin-Elmer & Co., GmbH,  
Postfach 1120, 7770 Ueberlingen, Federal  
Republic of Germany. Tel: (07551) 811.

Perkin-Elmer Ltd., Post Office Lane,  
Beaconsfield, Bucks HP9 1QA, England.  
Tel: Beaconsfield (049 46) 6161.

Circle No. 125 on Readers' Service Card

# PERKIN-ELMER

The science and computer company.  
Where solutions come first.



## SCIENCE / SCOPE

A redesigned equipment rack for the F/A-18 strike fighter's radar will allow the U.S. Navy and Marine Corps to save millions of dollars on new radars over the life of the program. The original rack consisted of many honeycombed sections bonded together. The Value-Engineered rack is a single piece fabricated by an automatic, numerically controlled machine center. The changes stem from technology that was not available at the time the original contract was signed. Hughes Aircraft Company will share some of the savings with the Navy through the Department of Defense Value Engineering program. This program is designed to encourage employees to look at the functions of a product and develop alternatives that cost less, perform better, and are more reliable.

The first U.S. facility for making gallium arsenide solar cells on a standard production line is now under construction at Spectrolab, Inc., a Hughes subsidiary. Gallium arsenide cells, which are now being made on a prototype line at Hughes Research Laboratories, will help satellites and spacecraft become more efficient in converting sunlight into electricity. Compared to conventional silicon cells, gallium arsenide cells generate up to 30% more power and operate at much higher temperatures. The first cells are expected to come off the production line in mid-1985. Full-scale mass production at rates to 15,000 cells per year is scheduled for January 1986.

Two high-power communications satellites have been ordered by Advanced Business Communications Inc. (ABCI). The satellites are widebody HS 393 models, now being developed at Hughes for launch from the space shuttle. The drum-shaped spacecraft are designed to fill the width of the space shuttle's cargo bay in order to take advantage of launch pricing policies. Each satellite will operate over the Ku band and carry 16 channels designed for business communications, video distribution, and teleconferencing. The spacecraft's higher power will let users receive signals through small, low-cost earth terminals. Hughes Communications Inc., a Hughes subsidiary, will provide launch, tracking, telemetry, and control services for the Ku band system. Hughes Communications Galaxy Inc. will market 16 transponders.

A microwave/millimeter-wave radar cross-section measurement system, designed to take automated measurements on full-size or scale-model targets, joins the solid-state millimeter-wave product line at Hughes. The new Model 42260H system, originally built for Boeing Military Airplane Co., is the latest in coherent short-pulse instrumentation measurement systems. It is suitable for indoor and outdoor ranges as close as 50 feet and as far away as beyond 4,000 feet. The system uses a modular design concept that accommodates up to six separate radio-frequency transceivers.

An ultramodern facility spanning 1.75 million square feet is the showcase where outstanding Hughes engineering combines advanced manufacturing techniques and production processes. Our complex is complete, so we're looking for experienced and graduating engineers to work on such programs as: infrared thermal imaging systems, laser rangefinders and designators, artificial intelligence, signal processing, VLSI, component qualification, hybrid microcircuitry, and focal plane. Send your resume to Hughes Electro-Optical and Data Systems Group, Professional Employment, P.O. Box 913, E9/W101, Dept. S2, El Segundo, CA 90245. Equal opportunity employer. U.S. citizenship required.

For more information write to: P.O. Box 11205, Marina del Rey, CA 90295






# YOU CAN JUDGE A BOOK BY ITS COVER!

*The latest MSD Isotopes catalogue of stable isotope labeled compounds. Almost two thousand compounds for research in chemistry, biology, physics and related fields. If you would like one—just call us...*

You could say we wrote the book on synthesized stable isotope labeled compounds. And here it is. The widest range of such compounds available anywhere. All of them made to the highest degrees of isotopic and chemical purity.

At MSD Isotopes we pioneered the commercial production of stable, isotope labeled compounds almost thirty years ago and today our catalogue lists almost 2,000 compounds labeled with deuterium, carbon-13, nitrogen-15 alone or in combination. Equally important are another 1,000



products which we do not list. In other words, if you don't find what you are looking for in our catalogue just call us. If your compound is not available our custom synthesis laboratory will be happy to consider making it for you.

If anyone can, we can, if it can be made.

Our latest catalogue contains over 250 new products and over three hundred price reductions as volume of demand for certain products has

reduced our costs.

The MSD Isotopes book of stable isotope labeled compounds. Just ask us for one. Its worth looking into.



#### West Coast

P.O. Box 2951  
Terminal Annex  
Los Angeles, CA 90051  
Los Angeles Area: (213) 723-9521  
Outside State of CA: (800) 423-4977  
State of CA: (800) 372-6454

## MSD ISOTOPES

DIVISION OF MERCK FROSST CANADA INC., Montreal, Canada

#### East Coast & Central

4545 Oleatha Avenue  
St. Louis, MO 63116  
Outside state of MO: (800) 325-9034  
State of Missouri: (314) 353-7000

#### Canada

MSD Isotopes  
P.O. Box 899  
Pointe Claire/Dorval, Quebec  
Canada H9R 4P7  
Telephone: (514) 697-2823

Circle No. 113 on Readers' Service Card



# Why Pay More for High Performance?



## **Chromatopac C-R3A Demands Comparison**

### **Chromatopac C-R3A Data Processor**

Shimadzu's cost- and performance-efficient data processor for chromatography includes:

- 176 Kbyte internal memory for up to 6 hrs general chromatogram storage.
- BASIC programming capability is standard.
- full alphanumeric capability.
- function keys provide true simplicity of operation.
- two-way communication with Shimadzu microprocessor-equipped GC and LC systems.
- two-way communication between another C-R3A or a computer.
- only low cost chromatography data processor with floppy disk expandability.

### **FDD-1A Floppy Disk Drive Accessory**

- single or dual drive
- 5" 1.2 megabyte floppy disk for storage of programs or data.
- up to 100 hrs general chromatogram storage, 20 hrs capillary chromatography storage.
- does not detract from large internal memory of C-R3A.
- allows use of CRT as well as print-copy of CRT display.

### **Chromatopac C-R3A with Accessories**

(Floppy Disk Drive, CRT)



• Write today for more information on these and other Shimadzu instruments.



*Shimadzu*

#### **SHIMADZU SCIENTIFIC INSTRUMENTS, INC.**

7102 Riverwood Drive, Columbia, Maryland 21046, U.S.A. Phone: (301) 997-1227

#### **SHIMADZU (EUROPA) GMBH**

Acker Strasse 111, 4000 Dusseldorf, F.R. Germany. Phone: (0211) 666371 Telex: 08586839

#### **SHIMADZU CORPORATION, INTERNATIONAL MARKETING DIV.**

Shinjuku-Mitsui Building, 1-1, Nishishinjuku 2-chome, Shinjuku-ku, Tokyo 160, Japan.  
Phone: Tokyo 03-346-5641 Telex: 0232-3291 SHMDT J.

Circle No. 149 on Readers' Service Card



# The Hoefer SE 600 Vertical Slab Unit has become the Standard for the Research Laboratory.

## **Multiple gel runs**

The SE 600 runs two 14 x 16cm slabs at a time; with divider plate accessories, four at a time.

## **Reproducible results**

Temperature regulation in the SE 600 really works. You can hold buffer temperature at 1°C. No smiling gels!

## **Easy casting**

Gel sandwiches can be assembled and sealed in the gel casting stand in two minutes. No leaks!

## **Volume gel production**

The new SE 615 Multiple Casting Stand (pictured below) casts ten gels at once, uniform pore or gradient.

## **Isoelectric focusing**

Both sides of the gel are cooled uniformly and protected from surface contact with CO<sub>2</sub>.

## **2-D gels**

Set up is easy using specially molded troughs in the upper buffer chamber.

## **Agarose gels**

Run them in the vertical mode between frosted glass plates.

## **Flexibility**

The SE 600 has more varieties of sample well formers and spacer thicknesses than any competitive unit on the market.

## **Safety**

Every Hoefer unit comes equipped with electrical safety interlocks.

**TRULY THE SLAB UNIT FOR ALL REASONS!**

For more information write or phone;

In California: 415-282-2307,


Outside California toll free: 800-227-4750

Designed, patented and manufactured by

**HOEFER SCIENTIFIC INSTRUMENTS**

654 Minnesota Street, P.O. Box 77387

San Francisco, CA 94107



Circle No. 44 on Readers' Service Card



# The Ventilated Animal Rack:

Some day soon most research animals, personnel, and research programs will be protected this way.\*

## What is a "Ventilated Animal Rack"?

It is a portable, totally enclosed animal rack with four separate, independent, isolated-from-each-other chambers. Low velocity air enters each chamber, makes a single pass over the cages, and is exhausted by negative pressure directly to the main exhaust system. This special rack (VR-1) most effectively isolates the animals from the animal room... and the research personnel from the animals.

## What are the benefits to the animals?

There are many. Cross-contamination is substantially decreased because air from an infected animal goes to the exhaust system with an absolute minimal exposure of the other animals. Animal stress is also significantly reduced: the enclosed environment is quiet; drafts and thermal and humidity fluctuations are greatly minimized; and animals can be easily observed without inducing stress. The success of this environment is attested to by the fact that the total number of animals born to a species that breeds poorly (DBA/2J mice) is increased and the percent survival is also appreciably higher. Additional evidence: judging by acceleration of weight gain, newly arrived animals housed in this system become acclimated more rapidly. Further evidence? Even multiple species can be successfully housed in the same rack.

## What are the benefits to the research workers?

Since the air in the rack is exhausted into the main exhaust system and does *not* re-enter the animal room itself, research workers are effectively isolated from animal dander or other allergens, odor, pheromones, microorganisms, and food and bedding dust. Even with the doors of the unit open, the direction of air flow tends to be *from* the room and *into* the unit which helps to contain contaminated air *within* the unit. Result: virtual elimination of allergic reactions and generally, a cleaner, safer, odor-free work environment for the research people.

## What are the benefits to research programs?

Because this system greatly reduces the chance of cross-contamination, and because it provides a much less stressful environment generally (e.g., it tends to reduce the amount of animal handling required), the chances of jeopardizing expensive research programs are substantially minimized.

## Are there other benefits?

The air velocity is variable and is separately adjustable for *each* shelf. The system offers a choice of bottle watering or a specially designed upfeed serpentine automatic watering configuration that eliminates stagnant water, permits flushing during the day, and significantly minimizes contamination. This rack also permits excellent space utilization since multiple species can be safely housed in the same room. Cleaning is easy; VR-1 can be handled by most standard rack washers. The unit is quiet. And, in summary, it is a most effective isolation system *that can actually divide a room into multiple separate, isolated environments.*

## From Lab Products, Inc.—The leader in environmental control products

Lab Products now offers the widest selection of systems for environmental protection: five Stay-Clean™ laminar flow systems; Isosystem™ housing system consisting of a disposable filter cap, cage cover, and plastic cage; Enviro-Gard™ filter system with permanent filter bonnets; and See-Through™ suspended cage systems with a special filtering system. We are now likely to have at hand solutions to virtually all of your environmental problems.

## For more information

Why wait for "some day soon"... write or call Lab Products, Inc., 255 West Spring Valley Avenue, Maywood, N.J. 07607 or complete the coupon. (phone 201 / 843-4600)

Lab Products, Inc.  
255 West Spring Valley Avenue  
Maywood, N.J. 07607

C

I would appreciate:

- ☐ Additional information on your VR-1 Ventilated Animal Rack.
- ☐ Information on your *other* environmental control products.
- ☐ Your 68-page catalog of animal housing and care systems and accessories.
- ☐ Seeing your local representative—Please call and set up an appointment:

Area Code \_\_\_\_\_ Number \_\_\_\_\_ Ext. \_\_\_\_\_

Name \_\_\_\_\_

Title \_\_\_\_\_

Organization \_\_\_\_\_

Address \_\_\_\_\_

Zip \_\_\_\_\_

\*Many of these systems are already installed in major research institutions... and conversion to these ventilated animal racks is accelerating.





The VR-1 is the subject of one or more pending patents.

**lab products inc**  
a **bioMedic** company

Circle No. 59 on Readers' Service Card

LP81-1667



# SECOND INTERNATIONAL CONGRESS ON COMPUTERS IN SCIENCE

28 October - 1 November 1984

Washington Hilton, Washington DC

Chairman: Stephen R. Heller  
Environmental Protection Agency, Washington DC 20460

*THEME: Workstations for the Scientist*

## OPENING ADDRESS

Computer Education/Literacy for the Scientists – C.D. Martin, George Washington University, Washington, D.C.

## SESSION TOPICS

Scientific Data Bases – L. Gevantman, National Bureau of Standards, Washington, D.C.  
Computer Aided Molecular Design – R. Potenzzone Jr., American Cyanamid Co., Stamford, Conn.  
Artificial Intelligence – P. Friedland, Stanford University, Stanford, CA.  
Biological Workstations – P. Armstrong, Intelligenetics, Inc., Palo Alto, CA.  
Management of the Electronic Laboratory – R. Dessy, Virginia Polytechnic Institute, Blacksburg, VA.  
Laboratory Automation and Robotics – J. Lemanowicz, Exxon Research and Engineering, Annadale, N.J.  
Workshops – Joseph L. Modelevsky, Lilly Research Labs.  
Molecular Biology Software: Conventions, distribution and exchange  
Molecular Biology Computing: Perils, prospects, and payoffs  
\* Other Workshops to be announced

## POSTER SESSIONS: WORKSTATION APPLICATIONS

Participants are invited to submit abstracts for the poster sessions. These abstracts will be reviewed up until the time of the meeting; however, only those accepted by 15 September 1984 will be listed in the conference program. For abstract forms contact address listed below.

## SCICOMP EXHIBITION

A large exhibit of state of the art hardware, software and systems will also be part of the overall conference.

**Organized by Scherago Associates, Inc.  
in conjunction with SCIENCE Magazine.**

**Scherago Associates, Inc., 1515 Broadway, New York, NY 10036, (212) 730-1050**

---

## REGISTRATION FEES: Make checks payable to: Scherago Assoc., Inc., COMPUTERS IN SCIENCE

\$200 On-site registration

\$150 ADVANCE REGISTRATION

\$125 Group Registrations – 4 or more received *together* from same organization.

\$ 50 Student Registrations – *Must* have status confirmed in writing by department chairman.

☐ Please reserve \_\_\_\_\_ space(s): Registration Fee of \$ \_\_\_\_\_ enclosed.

☐ Please send EXHIBIT SPACE information.

*Please Print*

Name \_\_\_\_\_ Street \_\_\_\_\_

Dept. \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip \_\_\_\_\_

Organization \_\_\_\_\_ Telephone: ( \_\_\_\_\_ ) \_\_\_\_\_

**Return to: Computers in Science, Scherago Associates, Inc., 1515 Broadway, New York, NY 10036, Tel: (212) 730-1050**

**All cancellations must be received in writing by Scherago Associates, Inc. no later than September 21, 1984 – no refunds will be given after this date.**

Circle No. 145 on Readers' Service Card



# YOU CAN COUNT ON TM ANALYTIC



For Quality

## $\beta$ and $\gamma$ Systems

When you are in need of accurate, reproducible results and maximum reliability of your liquid scintillation and gamma counting instruments, it makes sense to count on the company that has continually set the technological pace in these fields.

TM Analytic is the same dedicated organization that has designed, produced and serviced innovative, high-performance, nuclear counting systems for over 35 years.



# Tm Analytic

Allow us to tell you about the latest in beta and gamma counting technology soon. Let's discuss your applications, your problems, your present or future equipment requirements. Call or write us for a complete literature portfolio.

- **BetaTrac<sup>®</sup>, Mark III<sup>™</sup>, and Delta<sup>™</sup> Liquid Scintillation Systems.**
- **GammaTrac<sup>®</sup> Gamma Counting Systems.**
- **Counting Accessories.**

## TM Analytic, Inc.

1842 Brummel Drive  
Elk Grove Village, IL 60007  
800-323-5405  
In Illinois: (312) 364-9100

YOU CAN COUNT ON US





# Close your eyes.

Imagine you're creating a weekly current awareness service for scientists and scholars — one that would alert you to just those articles you were interested in. You think to yourself, this is a good idea...but impossible.

Now open your eyes. Such a service is available right now. It's called ASCA®. ASCA is custom-designed...for each individual...anywhere in the world...in any discipline. And your individualized ASCA reports can cost as little as \$3.00 per week. Intrigued? Want more information? Fill out the coupon and mail it. You'll be glad you did.

## Let ASCA Open Your Eyes.

☐ Send me more information about ASCA®.

Name / Title

Address

City

State

Country

Organization / Department

Zip

Phone

**ISI® Institute for  
Scientific Information®**

**Customer Services Department**

3501 Market Street, Philadelphia, PA 19104 U.S.A.

Tel: (215)386-0100, Cable: SCINFO, Telex: 84-5305

**European Office:** 132 High Street, Uxbridge,

Middlesex UB8 1DP, United Kingdom


Tel: 44-895-70016, Telex: 933693 UKISI

29-3549 © 1984 ISI



## UV/VIS PERFORMANCE:

The plot quickens!



Packed into the compact cabinetry of the Beckman DU-7 are a range of features that have established the standards for performance in spectrophotometry.

Our unique Stable Beam optics, together with a response time of just 0.05 second, allows the DU-7 to scan at a phenomenal 1200nm/min. — with no tracking error.

But you don't have to sacrifice sensitivity for speed. The DU-7's signal-to-noise ratio is far superior to any other comparable UV/Vis spectrophotometer. Baseline flatness is  $\pm 0.001A$  for precise readings. And the system has the ability to compensate for highly absorbing blanks up to approximately 2.0A.

Equally advantageous is the DU-7's ability to collect and manipulate large amounts of data easily. For *kinetic analysis* you simply view the data, subtract the blank and choose your time

interval. In five seconds your calculations are complete with statistical analysis. And, because it can gather and store 20 readings per second, it allows you to monitor even the fastest kinetic reactions.

To perform *protein analysis*, your options include both standard calibration curves and non-linear corrections for Lowry and Bradford curves, plus the ability to monitor and flag out-of-tolerance calibrations.

*Gel scanning* is equally comprehensive. In addition to tube gels, the DU-7 can handle dried slab gels and films, performing area and peak position calculations and graphic selection of sloping baseline.

Start plotting to acquire a DU-7 today by contacting your Beckman representative or Beckman Instruments, Inc., P.O. Box C-19600, Irvine, CA 92713.

© 1984 Beckman Instruments, Inc.

# BECKMAN

Circle No. 99 on Readers' Service Card





Early cuneiform tablet.

# 3500 B.C.

**To cope with a need to record information, man learns to write.**

# Today,

**DIALOG® helps scientists find key knowledge in a world overwhelmed with information.**

Thousands of years ago, an information revolution began when writing was invented in Mesopotamia.

Today, the world faces a revolution of a different kind. The problem isn't how to record information. It's how to find it.

The solution is DIALOG®, the world's largest scientific and technical databank. For over a decade we've been providing online access to scientific information.

To millions of journal articles, conference papers, technical reports, books and patents. The information you need to do your job.

With more databases than any other online system and the most powerful software available, DIALOG pinpoints new ideas—emerging technologies, product developments, research policies, social and economic issues.

Instantly, from your desktop terminal or personal computer. With no monthly minimums or initiation fees.

To learn more about this revolutionary tool, contact DIALOG Marketing, 3460 Hillview Ave., Dept. 83, Palo Alto, CA 94304. Call 800/227-1927 (800/982-5838 in Calif.).

**df DIALOG**

**INFORMATION SERVICES, INC.**

A subsidiary of the Lockheed Corporation.

© DIALOG Information Services, Inc. 1984

# Particles of Difference

IDC rigid uniform latex spheres, flown on 7th mission of Space Shuttle.

## Ultraclean latex from IDC.

IDC's uniformly sized, polystyrene latex particles are produced *without* the use of surfactants. As a result, tedious clean-up procedures are eliminated!

These particles can be suspended in buffers or sera for a wide variety of protein, enzyme, or antibody coupling applications.

In one short year, IDC's unique microspheres have become the accepted standard for biomedical applications. We supply highly reproducible batches, prepared under the most precisely controlled manufacturing conditions. It's a small wonder that these particles have become sizing equipment calibration standards.

IDC offers latex with sulfate, carboxyl, and carboxylate-modified surface groups, from 0.03 to 1.05 micron in small or large quantities. If the quality of your latex material makes a particle of difference to you, find out more about IDC—the standard.

Interfacial Dynamics  
P.O. Box 40306  
Portland, OR 97240  
(503) 230-6138



Specialists in colloidal and surface chemistry.  
Circle No. 101 on Readers' Service Card



# More from LKB Wallac with an extended RackBeta range

Features of the new RackBeta.Spectral make it an outstanding instrument for those who demand much from a beta counter. A new much faster CPU speeds up all aspects of your counting and result calculation. Eight count windows give much greater flexibility: in setting different windows, in counting several isotopes, in optimizing count conditions. The fast built-in "Help" function makes operation easier than ever before because it is simple to get immediate information. But above all in the Spectral, the 1024 channel multichannel analyser brings real advantages to LSC users: a completely new method of quench correction using SQP (spectral quench parameter), and the radical new spectrum display that really makes it possible for you to visualise your sample, to see the real spectrum almost immediately, to check for interfering effects such as chemiluminescence, and to select optimum windows.

The Spectral is the latest member of the RackBeta range that with the Universal and Compact models, and the various options, including a wide variety of data reduction programs, gives the most complete selection on the market. You can find out more by contacting:



LKB Instruments Inc., Nuclear Division, 9319 Gaither Road, Gaithersburg, Maryland 20877 USA

Circle No. 41 on Readers' Service Card



# Why settle for one When you can get **312 DETECTORS** Introducing the automated **AMBIS BETA SCANNING SYSTEM**

A unique imaging detection system that utilizes state-of-the-art instrumentation to solve today's tough gel scanning problems.

The Ambis system features a specially designed 2D proportional detector with 312 discrete elements that simultaneously scan gels two-dimensionally. This produces an image with up to 45,000 image points for resolution, accuracy, and sensitivity that replaces most autoradiography.



The Ambis system is controlled by a microcomputer that not only routinely performs the entire scanning procedure but also analyzes and quantitates all data quickly and precisely.

It is the Ambis software which provides the power to display "pictures" of scans, to extract channels of data and to compare patterns until a matching pattern is found. With this sophistication comes the simplicity of menu selection which guides the operator through each step of the software.

Contact us for more information on the Ambis line that is creating a revolution in resolution.




**automated microbiology systems inc.**

3939 Ruffin Road  
San Diego, CA 92123  
619-571-0113

Circle No. 201 on Readers' Service Card

## How To Buy A CO<sub>2</sub> Incubator Without Choking Up

Learn how a constant environment happens. Precision™ CO<sub>2</sub> incubators show you how. We've made more constant °C appliances for the lab than anyone else in the world.




**FREE GUIDE** on how to buy a CO<sub>2</sub> incubator plus literature. Write

GCA/Instruments and Equipment Group, 3737 West Cortland Street, Chicago, Illinois 60647. Or call 1-800-524-9482.

**Precision™**  
CO<sub>2</sub> Incubators

Digital pushbutton control to set °C and CO<sub>2</sub>%.  
Bright LED display for each.


Shelves and supports pop out to easily clean and decontaminate chamber.



GCA CORPORATION  
Instruments and Equipment Group

For a Sales Representative to call circle reader service number 116  
For Literature circle reader service number 120





## Easy separation of peptides — with the new ZORBAX® Bio Series Peptide Column and comprehensive User's Guide from DuPont.


The new ZORBAX® reversed-phase peptide column, the first in the versatile new Bio Series, is specifically designed for the successful separation of small proteins and peptides from protein digests. And to help you achieve the best possible separations, DuPont includes a detailed *User's Guide* with every column shipped.

The unique column configuration and packing characteristics of the ZORBAX® Peptide Column permit rapid separation of protein fragments with maximum resolution at normal flow rates and operating pressures. The exclusive ZORBAX® process provides exceptional uniformity of column packing and performance so you can expect your method to yield consistently repeatable results time after time.

To help you work out an effective method, the *User's Guide* clearly outlines specific procedures for developing and controlling reversed-phase separations of peptides with the new ZORBAX® column. Recommendations for selection

of gradient and other parameters are easy to follow and easy to apply to your particular separation challenge.

A *Separation Development Guide*, abstracted from the *User's Guide* and detailing the formulation of optimal separation conditions, is available to you at no charge. To obtain your copy, and for more information on the new ZORBAX® Bio Series, call DuPont at the following 24-hour toll free-number: (800) 345-8600, Ext. 502, or write: DuPont Company, Room X40540, Wilmington, DE 19898.



**DuPont makes HPLC make sense**  
**Liquid Chromatography Columns**





# NEW BioMag<sup>TM</sup> Magnetic Separations

## Fast and easy for:

- Particulate Systems
- Cell Sorting Applications
  - IgG Complexes
- Double Antibody Assays
- Monoclonal Screening
  - Serum Stripping
- Affinity Chromatography

## Separations are:


- Rapid
- Convenient
- Safe, No Containment Required
- No Centrifugation, No Aerosols
  - No Filters, No Fouling

## BioMag available as amine terminated particles or coupled to:

- Protein A
- Goat Anti-Rabbit IgG
- Goat Anti-Mouse IgG
  - Charcoal

Magnetic Separator Racks for  
Microtiter Plates or Test Tubes  
Available

For further information call or write




**Advanced Magnetics Inc.**

45 Spinelli Place • Cambridge, Massachusetts 02138  
Telephone (617) 497-2070  
Telex: 95-1417 TX Network BSN Ref: BIOC

Circle No. 121 on Readers' Service Card

# HELLMA<sup>®</sup> HPLC



A complete line  
of Flow Cells  
capable of  
turning any  
spectrophotometer  
into an  
HPLC detector.

Deuterium Lamps  
available for every  
instrument.

For 20 years...

**HELLMA**  
tomorrow's designs today.  
Hellma Cells, Inc.  
Box 544, Borough Hall Station  
Jamaica, New York 11424  
Telephone (212) 544-9534

Circle No. 43 on Readers' Service Card

## PRINT BIBLIOGRAPHIES IN THE STYLE OF ANY JOURNAL WITH YOUR MICROCOMPUTER AND Martz-BIBLIOFILE<sup>TM</sup>

A program which doesn't care which text editor you use  
(Wordstar ["Micropro"] is fine).

Martz-BIBLIOFILE<sup>TM</sup> remembers which journals you have used, and lists them by name  
on a menu. Martz-BIBLIOFILE<sup>TM</sup> can sequentially

### Renumber Citations in Your Manuscript,

then get the cited references from your master database, and put the bibliography for  
your manuscript in a standard ASCII file for printing. Also supports author/year  
citation.

### Alphabetize Any Bibliography or Your Entire Collection!

### Search with And-Or-Not Logic

on designated fields (author, year, title, journal/book,  
volume-pages, or keywords).

### Include Personalized Keywords

up to 500 characters (the limit per field; a reference can total up to 2000 characters). A  
SSDD 5 1/4" disk will hold 700 references; a DSDD 8" disk, 4000. Large collections can be  
spread over several floppy disks.

Martz-BIBLIOFILE<sup>TM</sup> is available for most 64K CP/M-80 Vers. 2.2 ("Digital Research")  
microcomputers for \$199 (8" SSDD) or \$219 (5 1/4" inch disks) including Apple II (with  
CP/M card only), DEC Rainbow, Epson QX-10, Heath/Zenith, Kaypro, Northstar,  
Osborne, Superbrain, Televideo, Xerox 820, Zorba, and many others. IBM-PC/MS-  
DOS version also available. Manual only, \$29.

POWER TOOLS, 48 Hunter's Hill Circle, Amherst MA 01002.

For detailed brochure, call 413-256-0751. M/C, Visa, PO.


Overseas air shipment, \$15. Check on foreign bank, \$15. Mass. residents add 5% tax.

Reviewed in *Annals of Internal Medicine* 100:621, *Educational Computer* April 1984, *Im-  
munology Today* 4:271.

Circle No. 177 on Readers' Service Card



# The collector that can't lose a drop, doesn't waste your time and won't cost a fortune.



*You'll find more LKB fraction collectors at work in laboratories throughout the world than any other.*

## The new HeliRac™ Fraction Collector from LKB

This new compact guarantees that no tube will ever change until the drop has been collected. So you won't waste your sample or your time. And the world's fastest time for tube change gives you perfect collection in any mode you choose: time, volume or drop. You'll also get simplified waste collection, with zero dead-volume and hence no contamination. So if you want reliability, HeliRac's got it. And if you want to save time, it's got every feature to automate your entire collection.

With more intelligence than many larger instruments, HeliRac is so simple to program and so easy to use. Let it run your system or be its slave. You'll get 6 time windows for optimal collection, 4 memories for program storage and an escape key to speed up parameter changes. That's the kind of flexibility you've always wanted in a collector so small—and with a price that's so sensible. When you've seen it at work in your own laboratory, you'll know why the performance and value of HeliRac just can't be matched. Try it today and you'll see for yourself.

### NOW FILL THIS IN TODAY

- ☐ Send me the facts.
- ☐ Show me how HeliRac can help my lab.
- ☐ Send me a quotation.

Name \_\_\_\_\_

Dep \_\_\_\_\_

Comp/Instit \_\_\_\_\_

Address \_\_\_\_\_

Tel \_\_\_\_\_

Circle No. 20 on Readers' Service Card



Head office: LKB-Produkter AB, Box 305, S-161 26 Bromma, Sweden  
Tel 08-98 00 40, telex 10492

Main US sales office: LKB Instruments, Inc. 9319 Gaither Road, Gaithersburg, Maryland 20877. Tel 301-963-3200, telex 909870 (dom.), 64634 (intern.)

UK sales office: LKB Instruments Ltd., 232 Addington Road, S. Croydon, Surrey CR2 8YD, England. Tel 01-657 8822, telex 2644 14

Other sales offices in: Antwerp, Athens (for Middle East), Bangkok, Copenhagen, Lucerne, Madras, Moscow, Munich, Paris, Rome, Turku, Vienna, Zoetermeer



NEW UV/VIS:

# Spectrophotometry hasn't been this simple since we invented it.

1941

Dr. Beckman introduces the first practical UV spectrophotometer, the DU. Thousands of these dependable instruments are still in use.



The Beckman Model B offers a low-cost, easy-to-use solution for routine, repetitive testing. This instrument was sold for 25 years.

1949





The DB Series introduces the first simple double-beam optical system utilizing programmed slits for scanning.

1960

The 20 Series with digital readout, scanning and modular accessories continues the tradition of affordable, easy-to-use instruments.

1973



1984

Remember the first Beckman Spectrophotometer you used? The simplicity and dependability?

The tradition continues.

Behind the simple operation of the new DU-50 Series Spectrophotometers are technological advances that give these instruments incredible capability.

They feature full programmability so you have control of all operations. Once you've

established your test parameters, you simply store them on a simple plug-in Memory-Pac™

module for instant reuse.

And that's that. You save time. You eliminate errors.

There are four models available offering you the choice of scanning or programmability or both.

All of them deliver the day in, day out dependability we built an industry on.

A utilitarian, functional instrument so simple anyone can run it; so dependable you'll never give it a second thought.

Now is the perfect time to find out about it. For more information, contact your Beckman representative or Beckman Instruments, Inc., Spectroscopy Instruments Operations, Campus Drive at Jamboree Blvd., P.O. Box C-19600, Irvine, CA 92713; (714) 833-0751.






## Now Available!

### Origin and Evolution of Deserts

Edited by Stephen G. Wells  
and Donald R. Haragan

This is the first study to integrate current geosciences research into an overview of four of the world's major deserts: the Sahara, the Great Basin, the Mohave, and the Sonoran. The nine essays discuss geology, vegetation, climate, and land management issues.

1983 Paper \$12.00



#### Also available

#### The Reclamation of Disturbed Arid Lands

Edited by Robert A. Wright  
1978 Paper \$8.00

#### Energy Resource Recovery in Arid Lands

Edited by Klaus D. Timmerhaus  
1981 Paper \$8.00

All orders must be prepaid.  
Send to **M. Michelle Balcomb,**  
**Exec. Officer, SWARM/**  
**AAAS, Colorado Mountain**  
**College, Spring Valley Cam-**  
**pus, 3000 County Rd. 114,**  
**Glenwood Springs, CO**  
**81601.** Ask for a **free** listing of  
related titles.

quired characteristics, is more precisely the stable inheritance of *adaptive* changes induced in an individual organism by an altered environment. In other words, the acquired characteristics that are inherited must be adaptive with respect to the environmental change that induced them if the inheritance is to be considered Lamarckian. So far as I know, no authenticated evidence for this type of heredity exists, although many false claims have been made.

In the cases described by Marx, the phenotypic changes seem to have had no *adaptive* relation to the altered environments that provoked them and therefore have nothing whatever to do with Lamarckian evolution. They are, to me, not "reminiscent" of Lamarckism, but of a number of misguided attempts to use experimental results showing nonadaptive hereditary change induced by an environmental alteration as evidence for the Lamarckian model of evolution; these experiments, like those cited by Marx, are interesting and important, but do not constitute evidence for, or even bear on, "the Lamarckian concept of evolution."

J. GORDIN KAPLAN

Department of Biochemistry,  
University of Alberta,  
Edmonton, Alberta, Canada T6G 2J9

#### On Making Comparisons: Reminded Again

As scientists we frequently make comparisons: between arbitrarily defined groups, between species, between countries, between ideas, theories, or whatever else takes our fancy. Yet, in making these comparisons, we are frequently guilty of a gross illiteracy. I refer to the common, almost universal, practice of comparing one thing to another. The verb *to compare* may be used with either the preposition *with* or *to*, but the two uses carry almost opposite meanings. If one thing is *compared with* another, the two are set side by side and the degree to which they differ is determined. This is the meaning usually implied in scientific writings. If one thing is *compared to* another, however, the two things are being likened and their similarities emphasized, as, for example, in "Shall I compare thee to a summer's day?"

It could be argued perhaps that the common usage of *to compare to* is sanctioned, by the frequency of its use, as an acceptable alternative to the less common (but correct) usage of *to compare*

*with*. At the risk of being pedantic, I do not think this argument can be accepted even for common American parlance; the latest edition of *Webster's Unabridged Dictionary* clearly reiterates the above distinctions.

Scientists pride themselves on their clarity and precision of thought. We should be no less concerned about the clarity and precision of what we write. This applies not only to authors of scientific works but more important to their editors, who have a particular responsibility to ensure the accuracy of the written word; the latter could surely afford the extra two ems to make comparisons *with* instead of *to*.

IAN A. GREAVES

Occupational Health Program,  
Harvard School of Public Health,  
Boston, Massachusetts 02115

#### Recombinant RNA Research

I would like to offer supplementary information about the article "The birth of recombinant RNA technology" by Roger Lewin (23 Dec. 1983, p. 1313). First, *in vitro* recombinant RNA has been practiced in many laboratories for several years, thanks to the pioneering work of Uhlenbeck, Gumpert, and their co-workers at the University of Illinois. Lewin's article attributes this technology to E. Miele, D. Mills, and F. R. Kramer (*1*). It is clear in Lewin's article that the major impact of the work by Miele *et al.* is the *in vitro* amplification of RNA with biological function by the introduction of foreign RNA into a vector derived from Q $\beta$  RNA. However, in a similar article published in May 1982, Shen Tongjian and Jiang Meiyuan reported the introduction of poly(A) into Q $\beta$  RNA, which was subsequently used to infect and thereby amplify the RNA insert. Although characterization of both the *in vitro* recombinant RNA and the product of replication *in vivo* was not fully described, it appears that the Chinese scientists have achieved the *in vivo* amplification of recombinant RNA.

R. J. CEDERGREN

Department of Biochemistry,  
Université de Montréal,  
Montréal, Québec, Canada H3C 3J7


#### References and Notes

1. E. Miele, D. Mills, F. R. Kraker, *J. Mol. Biol.* 171, 281 (1983).

**Erratum:** In the last paragraph of Arthur L. Robinson's Research News article "High spatial resolution ion microprobe" (14 Sept., p. 1137), J. Ronald Hass' name was misspelled.



# BIOSYSTEMS UPDATE



**Over 350 laboratories worldwide now use Applied Biosystems products. In just two years, these remarkable instrument-reagent systems have become the standards for performance, reliability and ease of operation.**

## **Protein/Peptide Sequencing, Oligonucleotide Synthesis, and Peptide Synthesis**

Our contributions to the chemistries and automation of these techniques will enhance your research, reduce your operating costs and free you for more productive work.

## **A Commitment to Scientific Contributions**

We pioneered automated systems for *gas phase* protein sequencing, *phosphoramidites* in oligonucleotide synthesis and *preactivation of amino acids* in peptide synthesis. Each one of these innovations has resulted in improved performance over conventional approaches.

## **A Commitment to Reliability**

Scientists and researchers in over 350 laboratories will attest to the reliability of our instrument-reagent systems. We design, manufacture and test instruments and reagents

side-by-side to ensure the optimum system performance you expect.

## **A Commitment to Product Support**

The quality of our worldwide product support is unmatched. We install each new system and teach you how to get the best from it. We do not leave your laboratory until your system's performance meets published specifications.

Follow-up service, if needed, is equally thorough and just a telephone call away from our network of product support centers. In addition, our research team is prepared to help you solve problems unique to your own applications—whether in research or in the production of peptides or oligonucleotides.

For more information, ask your local Applied Biosystems representative, contact us at one of the offices listed below, or circle reader service number 93.



APPLIED BIOSYSTEMS, INC., 850 Lincoln Centre Drive, Foster City, CA 94404 • (800) 874-9868 • In California (800) 831-3582 Telex: 470052  
IN EUROPE: APPLIED BIOSYSTEMS GMBH, Bergstrasse 104, D6102 Pfungstadt, West Germany • 06157-6036 Telex: 4191746






## New standards of vibration isolation **MICRO-g<sup>®</sup>** systems

- Standard Post-Mount Isolators combined with Series 76 Honeycomb Tops achieve exceptional levels of attenuation and damping.
- Gimbal Piston isolators attenuate both vertical and horizontal vibration in a single compact module.
- New Series 76 honeycomb core tops offer superior damping, rigidity and structural integrity.
- Standard sizes to 4' x 12'; top thickness 2", 4", 8" and 12"; wide range of load ratings.

### Technical Manufacturing Corporation

15 Centennial Drive, Peabody, Massachusetts 01960  
Telephone: 617-532-6330 • Telex: 951408



TMC-14R

Circle No. 67 on Readers' Service Card

## SCIENCE

### Posters

The following posters of *Science* covers are available:

30 March 1979, Tropical flowering tree;

23 February 1983, Landsat photo of Detroit, Michigan;

29 July 1983, Cheetah;


2 December 1983, Snowshoe hare;

23 December 1983, Cathedral window/DNA molecule.

Combination of space covers in scroll format: 1 June 1979, 23 November 1979, 10 April 1981.


Price is \$5 each (prepaid).

Write to AAAS, Department POST, 1515 Massachusetts Avenue, NW, Washington, D.C. 20005.



## Environmental Control From Percival

Have control of the environment with a biological, insect rearing, dew, or plant growth chamber from Percival. Whether you need a table model or a giant walk-in, Percival has the knowledge, experience, and quality product line to meet your specific needs. Or, Percival will prepare a recommendation to fit your requirements. Write today for more information and the complete Percival catalog or call collect 515/432-6501.



*Percival*

**Percival Manufacturing Company**  
P.O. Box 249 Dept. S-10 Boone, IA 50036

... The name to remember for versatile Biological Incubators, Insect Rearing Chambers, Dew Chambers, and Plant Growth Chambers.

Circle No. 10 on Readers' Service Card



# MPF-66 fluorescence combines superior optical performance and powerful data handling.



The new Perkin-Elmer MPF-66 Fluorescence Spectrophotometer combines the superior optical performance that you have come to expect from the MPF Series with the state-of-the-art Perkin-Elmer Series 7000 Professional Computer. A powerful, yet easy-to-use, software program provides simple instrument operation and data manipulation.

## **The MPF-66 provides excellent optical performance.**

The MPF-66 optimizes both sensitivity and resolution, the two most important considerations in fluorescence. A new optical design which includes continuously variable slits allows you to select the largest slit width that provides the desired resolution, maximizing the sensitivity of your analysis.

## **A versatile microcomputer.**

All instrument commands are entered via the Series 7000 Professional Computer, a unit specifically designed for the analytical laboratory. Soft keys on the computer make instrument operation simple. The Series 7000 Computers also support BASIC and FORTRAN, word processing and electronic spreadsheets.

## **Data manipulation-data you need in the format you want.**

Perkin-Elmer's Computerized Luminescence Spectroscopy Software (PECLS) program, regarded as the most complete and flexible data manipulation program available for fluorescence, can be used with the MPF-66. PECLS includes an interactive viewing feature to let you see the part of

the spectrum that you really need. The "OBEY" programming feature is included so you can automate your analysis and format the output.

## **Find out more!**

For more information about Perkin-Elmer's new Model MPF-66 Fluorescence Spectrophotometer, call (800) 323-7155 (in Illinois, call (312) 887-0770). If you prefer, contact one of the offices below.

Perkin-Elmer Corp., Analytical Instruments,  
Main Ave. (MS-12), Norwalk, CT 06856  
U.S.A. Tel: (203) 762-1000. Telex 965-954

Bodenseewerk Perkin-Elmer & Co., GmbH,  
Postfach 1120, 7770 Ueberlingen, Federal  
Republic of Germany. Tel: (07751) 811

Perkin-Elmer Ltd., Post Office Lane,  
Beaconsfield, Bucks HP9 1QA, England.  
Tel: Beaconsfield (049 46) 6161

# PERKIN-ELMER



# Pure Magic.



You're looking at an electron micrograph of MonoBeads®—the most significant advancement in separation technology since the discovery of Sephadex®. And it's the MonoBeads which give a magical high performance to our new Fast Protein Liquid Chromatography (FPLC) System.

With a uniform particle size of 9.8 µm, this new hydrophilic, monodisperse medium is the secret of high performance.

High resolution of complex biological samples such as proteins, peptides, enzymes, monoclonal antibodies and polynucleotides can now be conjured up with typical separation times of 10-20 minutes.

**A small selection of the substances which have already been studied using the FPLC System.**

Bacitracin.	Tryptic peptides.	L-amino acid
Monoclonal antibodies.	Haemoglobin A <sub>1c</sub> .	oxidase
Serum proteins.	Haemoglobin chains.	isoenzymes.
Urine proteins.	Insulins.	Restriction enzymes.
Synthetic oligonucleotides.	NADP.	Labile citrate synthase.
Viral proteins.	Venoms.	Placental homogenate.
Membrane proteins.	Pollen extracts.	Steroids.
	Ovalbumin.	
	Gc-globulin.	

When we had perfected the separation media, we turned our attention to the instrumentation. And we designed a system which represents the state of the art in protein separation technology.

We've included a gradient programmer, non-stainless steel pumps, programmable

fraction collector, UV monitor and pre-packed separation columns to fulfill all your needs for biological separations. Now, in just one or two steps, electrophoretically pure enzymes can be purified from crude tissue extract in just 10 to 40 minutes, with an overall yield in enzyme activity of greater than 80%.



In fact, we're so confident that we invite you to try out any biological sample, in aqueous or organic solution, so that you can see for yourself just how magical our new FPLC System would be for you.

**Pharmacia Fine Chemicals**

**The Pharmacia FPLC System. Pure Magic. Pure Logic.**

**Pharmacia Fine Chemicals**

800 Centennial Avenue, Piscataway, New Jersey 08854

Telephone: (201) 457-8068 or call your local representative.


Pharmacia Products will be on display at ASCB - Booths 530, 531, 532.

Circle No. 40 on Readers' Service Card



# Introducing RepliPlate™

## Colony Transfer Pads. The precise, convenient method of whole plate *replication*. *Put It to the Test!*



RepliPlate colony transfer pads are designed specifically for replica plating colonies of bacteria. By using a specially coated pad, RepliPlate can easily deliver at least six precise replicates.

RepliPlate pads can produce whole plate colony replicas with more precision, speed, and ease than velvet pads or nitrocellulose membranes\*.

RepliPlate pads are pre-sterilized and ready-to-use so there is no fear of contamination, no special preparation or storage, and convenient biohazard containment is provided.

Put RepliPlate colony transfer pads to the test in your laboratory! Discover for yourself how RepliPlate pads offer more economy than

conventional methods through time savings and accuracy.

For further information and to order, contact FMC Corporation, Marine Colloids Division, BioProducts Department, 5 Maple Street, Rockland, ME 04841.

**1-800-341-1574.** In Maine call (207) 594-4436.

\* Results may depend upon the microorganisms used, their stage of growth, plating density and medium used.


RepliPlate and **FMC** are trademarks of FMC Corporation. © FMC Corporation, 1984

**FMC** Marine Colloids  
BioProducts




# Waters™ HPLC Systems for Amino Acid Analysis.

## Now analyze primary and secondary amino acids with one picomole sensitivity.




**Introducing Waters high-speed PICO-TAG™ method.** This innovative pre-column derivatization method is pre-packaged for guaranteed results. In 15 minutes, you can run a complete compositional analysis from as little as 100 nanograms of protein, with a detection limit of one picomole. Waters provides the complete methodology, plus reagents, buffers, column, Waters HPLC System and PICO-TAG Work Station.

**With Waters Systems, choose the method that best meets your needs.**



**Waters AUTO-TAG™ OPA fully automated derivatization method.** This exclusive method offers fully automated pre-column derivatization, separation and data reduction for high sensitivity analysis of primary amino acids. Perfect for process monitoring of amino acids in fermentation broths and cell cultures.



**Classical Ion Exchange method with better quantitation.** Waters improves this high resolution method with raw chromatographic data storage and re-integration capability for better quantitation, and full documentation of method, chromatography and integration.

**More than offering just a choice of methodologies, Waters helps you get results.**

For more information, contact Waters Life Science Dept., 34 Maple Street, Milford, MA 01757, 617-478-2000 (X 2665).

**MILLIPORE**  
Waters Chromatography Division



## AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE

*Science* serves its readers as a forum for the presentation and discussion of important issues related to the advancement of science, including the presentation of minority or conflicting points of view, rather than by publishing only material on which a consensus has been reached. Accordingly, all articles published in *Science*—including editorials, news and comment, and book reviews—are signed and reflect the individual views of the authors and not official points of view adopted by the AAAS or the institutions with which the authors are affiliated.

### Editorial Board

FREDERICK R. BLATTNER, BERNARD F. BURKE, ARNOLD DEMAIN, CHARLES L. DRAKE, ARTHUR F. FINDEIS, E. PETER GEIDUSCHEK, GLYNN ISAAC, NEAL E. MILLER, FREDERICK MOSTELLER, ALLEN NEWELL, RUTH PATRICK, BRYANT W. ROSSITER, VERA C. RUBIN, WILLIAM P. SLICHTER, SOLOMON H. SNYDER, PAUL E. WAGGONER, JOHN WOOD

**Publisher:** WILLIAM D. CAREY

**Editor:** PHILIP H. ABELSON

### Editorial Staff

*Assistant Managing Editor:* JOHN E. RINGLE  
*Production Editor:* ELLEN E. MURPHY  
*Business Manager:* HANS NUSSBAUM  
*News Editor:* BARBARA J. CULLITON  
*News and Comment:* COLIN NORMAN (deputy editor), JEFFREY L. FOX, CONSTANCE HOLDEN, ELIOT MARSHALL, R. JEFFREY SMITH, MARJORIE SUN, JOHN WALSH  
*European Correspondent:* DAVID DICKSON  
*Contributing Writer:* LUTHER J. CARTER  
*Research News:* ROGER LEWIN (deputy editor), RICHARD A. KERR, GINA KOLATA, JEAN L. MARX, THOMAS H. MAUGH II, ARTHUR L. ROBINSON, M. MITCHELL WALDROP  
*Administrative Assistant, News:* SCHERRAINE MACK;  
*Editorial Assistant, News:* FANNIE GROOM  
*Senior Editors:* ELEANORE BUTZ, RUTH KULSTAD, MARY PRESCOTT  
*Associate Editors:* MARTHA COLLINS, SYLVIA EBERHART, CAITILIN GORDON, WILLIAM GREAVES, LOIS SCHMITT  
*Assistant Editors:* STEPHEN KEPPLER, LISA MCCULLOUGH, EDITH MEYERS  
*Book Reviews:* KATHERINE LIVINGSTON, *Editor:* LINDA HEISERMAN, JANET KEGG  
*Letters:* CHRISTINE GILBERT  
*Copy Editor:* ISABELLA BOULDIN  
*Production:* JOHN BAKER; HOLLY BISHOP, ELEANOR WARNER; JEAN ROCKWOOD, SHARON RYAN, BEVERLY SHIELDS  
*Covers, Reprints, and Permissions:* GRAYCE FINGER, *Editor:* GERALDINE CRUMP, CORRINE HARRIS  
*Guide to Scientific Instruments:* RICHARD G. SOMMER  
*Editorial Administrator:* SUSAN ELLIOTT  
*Assistant to the Associate Publisher:* ROSE LOWERY  
*Assistant to the Managing Editor:* NANCY HARTNAGEL  
*Membership Recruitment:* GWENDOLYN HUDDLE  
*Member and Subscription Records:* ANN RAGLAND  
**EDITORIAL CORRESPONDENCE:** 1515 Massachusetts Avenue, NW, Washington, D.C. 20005. Area code 202. General Editorial Office, 467-4350; Book Reviews, 467-4367; Guide to Scientific Instruments, 467-4480; News and Comment, 467-4430; Reprints and Permissions, 467-4483; Research News, 467-4321. Cable: *Advancesci*, Washington. For "Information for Contributors," write to the editorial office or see page xi, *Science*, 28 September 1984.  
**BUSINESS CORRESPONDENCE:** Area Code 202. Membership and Subscriptions: 467-4417.

### Advertising Representatives

*Director:* EARL J. SCHERAGO  
*Production Manager:* DONNA RIVERA  
*Advertising Sales Manager:* RICHARD L. CHARLES  
*Marketing Manager:* HERBERT L. BURKLUND  
**Sales:** NEW YORK, N.Y. 10036: Steve Hamburger, 1515 Broadway (212-730-1050); SCOTCH PLAINS, N.J. 07076: C. Richard Callis, 12 Unami Lane (201-889-4873); CHICAGO, ILL. 60611: Jack Ryan, Room 2107, 919 N. Michigan Ave. (312-337-4973); BEVERLY HILLS, CALIF. 90211: Winn Nance, 111 N. La Cienega Blvd. (213-657-2772); SAN JOSE, CALIF. 95112: Bob Brindley, 310 S. 16 St. (408-998-4690); DORSET, VT. 05251: Fred W. Diefenbach, Kent Hill Rd. (802-867-5581).  
**ADVERTISING CORRESPONDENCE:** Tenth floor, 1515 Broadway, New York 10036 (212-730-1050).

## Analytical Instruments

Continued improvements in analytical instrumentation are perhaps the most crucial factor in the vitality of scientific research. This issue presents a sampling of further advances in analytical capabilities. In recent years the kinds of substances that can be studied have increased. Sensitivity of detection has been enhanced. Abilities to isolate minor constituents of complex mixtures have been improved. The time required for analytical procedures has diminished.

Gas chromatography coupled with mass spectrometry has been for some time the analytical method of choice for volatile substances. Even so, it continues to be improved, for example, through the introduction of durable, coated capillary columns. With such a device it has been possible to separate a mixture of 100 plant hydrocarbons in only 6 minutes.

Liquid chromatography is especially useful in dealing with hydrophilic substances, including many large molecules of biological importance. Resolution and speed of liquid chromatography columns have been improved. However, the usefulness of liquid chromatography has been less than ideal. The bottleneck has been detector systems. In special cases spectrophotometry, electrophoresis, or fluorescence may be employed, but these methods lack the power and versatility of mass spectrometry. Earlier when mass spectrometry was employed, the usual practice was to collect a peak, evaporate off the eluting liquid, and coat part of the residue on an electrode that was subsequently inserted into the mass spectrometer. The solid was then subjected to a high electric field, fast atom ion bombardment, or other excitation that created ions. New techniques now permit better interfacing of liquid chromatography with mass spectrometry. A thermospray method described in this issue facilitates removal of the eluting liquid and ionization of the eluted substances. Thus the fluid from the liquid chromatography can be injected directly into the mass spectrometer. The molecular weights of undegraded molecules can be determined as well as the weights of fragments.

The mass spectrometer can be used to provide information about solids that ordinarily cannot be analyzed by liquid chromatography. A small portion of the material is subjected to practically instantaneous pyrolysis within the mass spectrometer. The fragments produced give a good representation of the nature of the parent substances.

Mass spectrometers have been undergoing development with, for example, new ion sources and the tandem configuration. An interesting version is the Fourier transform mass spectrometer. The crucial volume in the device is a cube 2.54 centimeters on a side, bounded by electrodes, and maintained in a magnetic field. The device is operated as an ion trap with decay times of as much as half a day. Within the trap a series of operations is conducted, including ion formation, storage, manipulation, and detection. Resolution at a mass to charge ratio of 18 is 100,000,000.

The value of applications of nuclear magnetic resonance continues to increase. Organic chemists find it almost indispensable in the determination of structure of compounds present in a liquid phase. It is being applied increasingly to study of solids, and the earlier emphasis on observations of protons has been supplemented by measurements of other nuclides. One dynamic application is use in noninvasive imaging in medicine. For many tissues nuclear magnetic resonance is superior to the CAT scan.

Lasers are being increasingly used as key components in analytical processes. One use involves the application of tunable dye lasers and takes advantage of resonant multiple photon ionization to attain excellent selectivity. Illustrative of the power of the method is a demonstration that  $10^3$  atoms of krypton-81 can be sorted out and counted in a sample along with  $10^7$  atoms of krypton-80 and of krypton-82 plus  $10^{12}$  atoms of other atoms or molecules. With the same techniques it was possible to detect  $10^{-8}$  to  $10^{-10}$  of a monolayer of anthracene and naphthalene adsorbed on graphite.

—PHILIP H. ABELSON



Now as low as  
\$995

# Sartorius accuracy in 0.1 second flat!



Series 1400 and 1500 electronic toploaders offer unbeatable performance at affordable prices.

Electronic weighing has never been faster—or more accurate—than Sartorius makes it with Series 1400 and 1500 MP8 toploading balances. Display update speed—0.1 second when weighing in—practically eliminates lag time. Advanced electronic technology, including very large scale integration (VLSI), puts routine functions under microprocessor control. And optional built-in RS232C data output lets you interface your Sartorius MP8 balance to a printer, computer, or data keyboard for maximum system flexibility.

Both Series 1400 and 1500 feature compact, space-efficient design and are available in single-, dual-, or quad-range models with readabilities from 0.001 to 0.1 g and full-range electronic taring. Series 1400 toploaders weigh to 6000 g; Series 1500, to 12,000 g. Both include the unique AutoCheck feature, which

automatically reviews and verifies the status of all electronic circuits *before* the balance enters the weighing mode. And both feature new built-in keyboard options to simplify program capabilities such as counting, conversions, animal weighing, and formulation.

Put Sartorius speed and accuracy to work in *your* laboratory. For more information or a demonstration, call or write: Brinkmann Instruments Co., Division of Sybron Corporation, Cantiague Road, Westbury, NY 11590, Tel: 800-645-3050; in New York, 516-334-7500.



Sartorius electronic  
weighing systems  
Brinkmann

SYBRON