

AMERICAN
ASSOCIATION FOR THE
ADVANCEMENT OF
SCIENCE

SCIENCE


16 OCTOBER 1987
VOL. 238 ■ PAGES 249-436

\$2.50

*Bioanalytical
Instrumentation*



ASYST™ 2.0...Superiority Reflected By Its Performance



We Just Made The World's Best Engineering and Scientific Software ... Better!

With ASYST 2.0, we have taken the industry standard in PC based scientific software and expanded its formidable capabilities.

ASYST continues to be the only technical software available for the PC and compatibles which integrates data acquisition, data analysis, statistics and astounding graphics in a fully programmable environment. This provides great flexibility and offers you virtually limitless application opportunities.

ASYST provides its users with minicomputer speed at PC prices. ASYST supports analog and digital I/O in foreground and background operation on many popular boards and its A/D is fast and flexible. The IEEE-488 and RS-232 capabilities interface some 10,000 instruments. The graphics facilitate axonometric, cartesian and contour plots. You can perform waveform, matrix and polynomial analysis and a 1024-pt. complex FFT can be done in less than 2.5 seconds on

an IBM-AT or compatible. In addition, the software accommodates linear and non-linear curve fitting as well as smoothing ... and that just scratches the capabilities surface.

Now ASYST 2.0 has added these expanded performance features:

- A/D rate of 130 kHz (run on an 8 MHz, 286)
- Lotus® 1-2-3® interface
- Supports LIM expanded memory boards
- Overlay architecture for fully customized systems
- Set-up menus for rapid and easy configuration
- Support for new high-speed A/D boards
- Automatic command memory
- Dynamic arrays for programming simplicity
- Intelligent on-screen help

Thousands of users feel that ASYST software has no equal, in speed, power or flexibility. ASYST 2.0 has widened the competitive gap even further.

ASYST 2.0 makes your PC as powerful as a mini while it unlocks your own personal potential and expands productivity.

ASYST 2.0 includes free technical support for 60 days and a comprehensive Extended Support Plan is available.

We would like to include you among our ASYST User Group Members. Please call us today for more detailed information.

1-800-348-0033

ASYST is a trademark of Asyst Software Technologies, Inc.
IBM-AT is a registered trademark of International Business Machines Corp.
Lotus and 1-2-3 are registered trademarks of Lotus Development Corporation

ASYST
SOFTWARE TECHNOLOGIES, INC.

100 Corporate Woods • Rochester, N.Y. 14623

Circle No. 130 on Readers' Service Card

Brandel Introduces The Ultimate In Automatic Sample Preparation

Filter up to 48 samples, deposit filter discs into scintillation vials and dispense scintillation fluid to all vials simultaneously. Elapsed time? Less than two minutes. And that's just the beginning.

Flexibility? Brandel Harvesters can accommodate 96 well plates, test tubes or most any other well plate configuration. Interchangeable harvesting probes offer the ability to change from one test tube or well configuration to another—in seconds. Dispensers and Deposit/Dispensers can accommodate most any vial size

and rack configuration. The modular design permits each component to be used together as a system, or separately.

If you use a Beckman LS 5801 Series Counter, vial racks from the Beckman Counter can be loaded directly into the Brandel Deposit/Dispensing System. There's no need to dry filter paper, punch filter discs

by hand, measure cocktail or reload vial racks. It's automatic.

Accuracy? Consistent filtering, washing, depositing, and dispensing increase accuracy, and make obtaining repeatable results easier. Brandel systems are keyed alike to provide positive sample identification every step of the way.

Cost? Surprisingly affordable.



BRANDEL

8561 Atlas Drive, Gaithersburg, MD 20877 USA
(301) 948-6506 TLX:898083

Circle No. 250 on Readers' Service Card

255 This Week in *Science*

Editorial

257 Instrumentation and Equipment

Letters

259 Chemical Carcinogenesis: M. MEISLER; E. E. MCCONNELL ■ The Candidates' Budgets: P. IVERSEN ■ Human Control: D. H. SHAPIRO, JR., G. W. EVANS, J. SHAPIRO; J. RODIN ■ Comparison of High Schools: A. M. LEE

News & Comment

262 President's AIDS Panel in Disarray
263 Bumps and Falls on the Road to Stockholm
265 Spotlight Brighter on Minorities in Science
266 Reinventing the Space Truck

Research News

269 What Is Worse Than "The Big One"?
271 New Sequencers to Take on the Genome
273 My Close Cousin the Chimpanzee
276 "On the Shelf" AIDS Drug in Clinical Trial

Articles

Bioanalytical Instrumentation

305 Time-Resolved X-ray Diffraction of Biological Materials: S. M. GRUNER
312 New Opportunities in Synchrotron X-ray Crystallography: C. T. PREWITT, P. COPPENS, J. C. PHILLIPS, L. W. FINGER
319 The Role of Protein Structure in Chromatographic Behavior: F. E. REGNIER
323 High-Speed Chromosome Sorting: J. W. GRAY, P. N. DEAN, J. C. FUSCOE, D. C. PETERS, B. J. TRASK, G. J. VAN DEN ENGH, M. A. VAN DILLA
329 Magnetoencephalography and Epilepsy Research: D. F. ROSE, P. D. SMITH, S. SATO

Research Articles

336 A System for Rapid DNA Sequencing with Fluorescent Chain-Terminating Dideoxynucleotides: J. M. PROBER, G. L. TRAINOR, R. J. DAM, F. W. HOBBS, C. W. ROBERTSON, R. J. ZAGURSKY, A. J. COCUZZA, M. A. JENSEN *et al.*

- **SCIENCE** is published weekly on Friday, except the last week in December, and with an extra issue in February by the American Association for the Advancement of Science, 1333 H Street, NW, Washington, DC 20005. Second-class postage (publication No. 484460) paid at Washington, DC, and at an additional entry. Now combined with **The Scientific Monthly**® Copyright © 1987 by the American Association for the Advancement of Science. The title **SCIENCE** is a registered trademark of the AAAS. Domestic individual membership and subscription (51 issues): \$65. Domestic institutional subscription (51 issues): \$98. Foreign postage extra: Canada \$32, other (surface mail) \$27, air-surface via Amsterdam \$65. First class, airmail, school-year, and student rates on request. Single copies \$2.50 (\$3 by mail); back issues \$4 (\$4.50 by mail); Biotechnology issue, \$5.50 (\$6 by mail); classroom rates on request; Guide to Biotechnology Products and Instruments \$16 (\$17 by mail). **Change of address:** allow 6 weeks, giving old and new addresses and seven-digit account number. Authorization to photocopy material for internal or personal use under circumstances not falling within the fair use provisions of the Copyright Act is granted by AAAS to libraries and other users registered with the Copyright Clearance Center (CCC) Transactional Reporting Service, provided that the base fee of \$1 per copy plus \$0.10 per page is paid directly to CCC, 21 Congress Street, Salem, Massachusetts 01970. The identification code for *Science* is 0036-8075/83 \$1 + .10. **Postmaster:** Send Form 3579 to *Science*, 1333 H Street, NW, Washington, DC 20005. *Science* is indexed in the *Reader's Guide to Periodical Literature* and in several specialized indexes.
- The American Association for the Advancement of Science was founded in 1848 and incorporated in 1874. Its objects are to further the work of scientists, to facilitate cooperation among them, to foster scientific freedom and responsibility, to improve the effectiveness of science in the promotion of human welfare, and to increase public understanding and appreciation of the importance and promise of the methods of science in human progress.



COVER Molecular structure can be inferred from spectral analysis of biological samples. Our understanding of structure and function has become increasingly fundamental, in part, because of advances in bioanalytical instrumentation. See page 305. [Cover design by Sharon Wolfgang, Washington, D.C.]

Reports

- 342 The Microstructure of High-Critical Current Superconducting Films: P. CHAUDHARI, F. K. LEGOUES, A. SEGMÜLLER
- 344 Downregulation of L3T4⁺ Cytotoxic T Lymphocytes by Interleukin-2: C. C-Y. SHIH AND R. L. TRUITT
- 347 Conservation of the Duchenne Muscular Dystrophy Gene in Mice and Humans: E. P. HOFFMAN, A. P. MONACO, C. C. FEENER, L. M. KUNKEL
- 350 Depolarization Without Calcium Can Release γ -Aminobutyric Acid from a Retinal Neuron: E. A. SCHWARTZ
- 355 Blockade of "NMDA" Receptors Disrupts Experience-Dependent Plasticity of Kitten Striate Cortex: A. KLEINSCHMIDT, M. F. BEAR, W. SINGER
- 358 Conduction Velocity Variations Minimize Conduction Time Differences Among Retinal Ganglion Cell Axons: L. R. STANFORD
- 360 Heritability at the Species Level: Analysis of Geographic Ranges of Cretaceous Mollusks: D. JABLONSKI
- 363 Apolipoprotein B-48 Is the Product of a Messenger RNA with an Organ-Specific In-Frame Stop Codon: S.-H. CHEN, G. HABIB, C.-Y. YANG, Z.-W. GU, B. R. LEE, S. WENG, S. R. SILBERMAN, S.-J. CAI, J. P. DESLYPERE *et al.*
- 366 Behavioral Recovery Induced by Applied Electric Fields After Spinal Cord Hemisection in Guinea Pig: R. B. BORGES, A. R. BLIGHT, M. E. MCGINNIS
- 369 Phylogenetic Relations of Humans and African Apes from DNA Sequences in the $\psi\eta$ -Globin Region: M. M. MIYAMOTO, J. L. SLIGHTOM, M. GOODMAN
- 373 Linkage of Functional and Structural Heterogeneity in Proteins: Dynamic Hole Burning in Carboxymyoglobin: B. F. CAMPBELL, M. R. CHANCE, J. M. FRIEDMAN

AAAS Meetings

- 402 *Forum '87: Students and Science Learning* ■ Advance Registration and Housing Form

Book Reviews

- 409 A Very Special Relationship, *reviewed by* T. FINDLAY ■ The Shaky Game, D. HOWARD ■ Ecological Theory and Integrated Pest Management Practice, R. J. PROKOPY ■ Books Received

Products & Materials

- 412 Experimental Optimization Software ■ Reagents for Biotin Labeling ■ HPLC Software ■ Separatory Funnel Platforms ■ Culture Dish for High Magnifications ■ Antibacterial Antiviral Wipes ■ Calculator-Emulation Software ■ Beads for ELISA, RIA ■ HPLC-Grade Water Purifier ■ Literature

Board of Directors

Lawrence Bogorad
*Retiring President,
Chairman*

Sheila E. Widnall
President

Walter E. Massey
President-elect

Robert McC. Adams
Floyd E. Bloom
Mary E. Clutter
Mildred S. Dresselhaus
Beatrix A. Hamburg
Donald N. Langenberg
Frank von Hippel
Linda S. Wilson

William T. Golden
Treasurer
Alvin W. Trivelpiece
Executive Officer

Editorial Board

Elizabeth E. Bailey
David Baltimore
William F. Brinkman
Philip E. Converse
Joseph L. Goldstein
James D. Idol, Jr.
Leon Knopoff
Seymour Lipset
Oliver E. Nelson
David V. Ragone
David M. Raup
Vera C. Rubin
Larry L. Smarr
Solomon H. Snyder
Robert M. Solow
James D. Watson

Board of Reviewing Editors

John Abelson
Qais Al-Awqati
James P. Allison
Don L. Anderson
Elizabeth H. Blackburn
Floyd E. Bloom
Charles R. Cantor
Ralph J. Cicerone
James H. Clark
Bruce F. Eldridge
Stanley Falkow
Theodore H. Geballe
Roger I. M. Glass
Stephen P. Goff
Robert B. Goldberg

Corey S. Goodman
Stephen J. Gould
Richard M. Held
Gloria Heppner
Eric F. Johnson
Konrad B. Krauskopf
I. Robert Lehman
Karl L. Magleby
Joseph B. Martin
John C. McGiff
Alton Meister
Mortimer Mishkin
Peter Olson
Gordon H. Orians
Carl O. Pabo
John S. Pearce

Yeshayau Pocker
Jean Paul Revel
James E. Rothman
Thomas C. Schelling
Ronald H. Schwartz
Stephen M. Schwartz
Otto T. Solbrig
Robert T. N. Tjian
Virginia Trimble
Geerat J. Vermeij
Martin G. Weigert
Harold Weintraub
Irving L. Weissman
George M. Whitesides
Owen N. Witte
William B. Wood



Wouldn't you prefer to have the proven leader in peptide chemistry behind your research?


**Our total system for
synthesis & purification
consistently provides
the purest peptides.**

Synthesis


- Flexible, versatile, fully programmable.
- Multiple chemistry approaches. *t*-BOC / Fmoc / NMP.DMSO
- Efficient small scale peptide production with under 30 - minute cycle times.

Purification

- Highest purity of up to gram quantities of peptide.
- Flexibility of preparative and analytical operation with a simple column change.
- Confidence with a complete system including optimized protocols, columns and reagents.



Analytical Screen of Crude Peptide



Preparative Isolation



Final Analysis of Purified Product

Applied Biosystems is committed to total customer support. You get advanced applications, quality instrumentation, reagents and complete after sales support from the proven leader. For information on our complete peptide chemistry system contact us at the nearest office below.

**Applied
Biosystems**
Foster City Biosystems Division

Sales and Service Worldwide: U.S.A. 777 Lincoln Centre Drive, Foster City, California 94404 U.S.A., Tel. (415) 570-6667, (800) 874-9868, in Calif. (800) 831-3582, Telex 470052 U.K. Birchwood Science Park North, Warrington, Cheshire, U.K., WA3 7PB, Tel. 0925-825650, Telex 629611 • Japan 5-13-9, Nakakasai, Edogawa-Ku, Tokyo, Japan 134, Tel. (03) 877-0071 W. Germany Robert Koch Strasse 16, 6108 Weiterstadt, B.R.D., Tel. 06151-87940, Telex 4197318 • Australia Suite 2, 1401 Burke Road, Kew East, Victoria, Australia 3102, Tel. (03) 859-9571 France B.P. 50086, 13 Rue de la Perdrix, Z.A.C. Paris Nord II, 95948 Roissy, Charles de Gaulle, Cedex, R.F., Tel. (1) 48 63 24 44, Telex (1) 230 458 ABIF.

This Week in SCIENCE

Mapping the brain's magnetic fields

ABNORMAL electric discharges from neurons in the cortex of the brain accompany seizures in epileptic patients; these discharges generate magnetic fields that can be detected outside the brain (the skull and scalp do not interfere with the fields) by magnetoencephalography (MEG) (page 329). For some (an estimated 50,000 in the United States) patients with epilepsy, surgical removal of epileptogenic subregions could provide freedom from seizures; this requires exact identification of where the unusual neuronal discharges are originating in the brain. Rose *et al.* discuss the components of MEG technology, how environmental magnetic fields are canceled from measurements (brain fields are about one-billionth the size of the earth's magnetic field), how data are collected and projected into three-dimensional maps, how these data both compare and can be coordinated with data from other noninvasive (like MEG) and invasive techniques, and expected improvements in the technology that should make it increasingly useful in the future. Other papers in this special issue on bioanalytical instrumentation appear on pages 305 to 335 and are discussed in Abelson's editorial (page 257) and Roberts' Research News feature (page 271).

Superconductor microstructure

THE appearance at close range of one of the new superconductors—a film of oxides of three metallic substances (yttrium, copper, and barium) deposited on a strontium titanate crystal (the substrate)—is described by Chaudhari *et al.* (page 342). The microstructure of the superconducting film was studied by x-ray diffraction methods and by transmission electron microscopy. The film was aligned with the substrate in a true epitaxy; that is, the crystal structures of film and substrate are the same and the

interface between them is nearly flat atomically. Within the system there were many precipitates of yttrium oxides distributed in random orientations and fewer precipitates of copper oxides distributed nonrandomly. The superconducting layer contained numerous twin lattice structures, mirror images of each other. Differences in the microstructure of the superconducting film with that of a polycrystalline film provide clues to the features of the superconductor that contribute to its ability to carry large currents.

Cycling cytotoxicity of T cells

ENVIRONMENTAL signals can affect the ability of cytotoxic T lymphocytes (CTLs) to lyse target cells (page 344). CTLs kill target cells after exposure to target-cell antigens. Shih and Truitt show that exposure of lytic cells to interleukin-2 (IL-2) halts the killing; in its absence, CTLs reexpress their lytic activity. The cells can cycle back and forth many times between their lytic and nonlytic phenotypes. Although IL-2 inhibits lytic activity, it continues to promote clonal proliferation. This in vitro model may illustrate what happens in the body when this subset of cells is alternately and intermittently exposed to antigen, IL-2, and perhaps other up- and down-regulating substances. It may also explain why lytic activity sometimes cannot be demonstrated in vitro for mixed cell populations or for cells grown in the presence of IL-2 producers.

Calcium-independent transmitter release


RELEASE of neurotransmitters from nerve cells often is dependent upon the availability of calcium ions (page 350). However, it is also possible for a neurotransmitter to be released without calcium acting as a mediator. Schwartz shows that release of the neurotransmitter gamma-amino-

butyric acid (GABA) from catfish retinal horizontal cells occurs directly in response to depolarization of the cell. The release of GABA in cell cultures was detected by a second cell type, the bipolar cell from the goldfish retina; exposure of this cell to GABA opens chloride ion channels, negative ions flow into the cell, and a large positive outward current is produced. Measurement of the outward currents of the bipolar cell and measurement of the membrane voltage of the horizontal cell showed that the two "communicated" through GABA under calcium-independent conditions. This type of neurotransmitter release may operate in other parts of the brain as well since most synapses in the brain use GABA or glutamate as transmitters.

Humans and African apes

CHIMPANZEES and humans may be more closely related to each other than either is to the gorillas (page 369). Thus, shared attributes of chimpanzees and gorillas (such as their ability to walk on their knuckles) are possibly primitive features lost by humans as they evolved. The phylogenetic relations of these primates were drawn from new sequence data of the psi-eta region of the beta globin gene. To known sequence data of the $\psi\eta$ locus, Miyamoto *et al.* added the sequences of regions flanking this locus, thereby expanding the sequenced stretch to 7.1 kilobase pairs. Overall, there are few differences (1.6 to 2.1%) among all sequences compared for these three primates that are believed to have branched from each other 5 to 10 million years ago. The differences found indicate that nucleotide evolution has not occurred at the same rate in all (the human evolutionary rate has been the slowest), and this further complicates determinations of how long species have been independently evolving. Lewin explains how these and other molecular data lead to a different conclusion than do anatomic data about who is whose closest cousin (page 273).

INTRODUCING AN UNPRECEDENTED CONCEPT IN THE LITERATURE OF SCIENCE:



THE ISI ATLAS OF SCIENCE®

INNOVATIVE. TIMELY. UNIQUE.

That's the *ISI Atlas of Science*, an exciting new series of publications that explores the present and future direction of research in science, medicine, and technology. The *Atlas* combines the power of computerized citation analysis with the expert interpretation of leading researchers to guide you through the landscape of science, helping you understand what is happening in science research today . . . and where it might lead tomorrow.

The *Atlas* series will consist of 12 sections, each focusing on a distinct area of scientific research. The first section, the *ISI Atlas of Science: Pharmacology*,

is available now. Over the next four years, 11 additional *Atlas* sections will begin regular publication.

A yearly subscription to each *Atlas* section consists of four quarterly issues, each containing 20-30 surveys analyzing the latest research in the topic covered. Institutional subscriptions also include an annual cumulation containing all surveys for that year, plus extensive indexes and reference materials.

Flexible . . .

A variety of subscription options are available for the *ISI Atlas of Science*, designed to meet the research and budgetary needs of individuals, departments, and small or large libraries. And, for a limited time, subscribers can take advantage of our Charter Subscription

Program and receive all 12 sections of the *Atlas* as they become available . . . at significantly reduced rates!

Explore your options . . .

To explore the variety of *ISI Atlas of Science* options available, call toll-free, 1-800-523-1850, extension 1405. Outside the United States and Canada, call 215-386-0100, extension 1405. We'll send you a complimentary information package, including a detailed brochure and sample articles from the *ISI Atlas of Science: Pharmacology* . . . at no obligation.

**FROM THE PUBLISHERS OF
CURRENT CONTENTS®
AND THE
SCIENCE CITATION INDEX®.**

isi®

Institute for Scientific Information®

©1987 29-4981

3501 Market Street, Philadelphia, PA 19104 U.S.A., Telephone: (215) 386-0100, ext. 1405, Cable: SCINFO, Telex: 84-5305
European Branch: 132 High Street, Uxbridge, Middlesex UB8 1DP, United Kingdom, Telephone: 44-895-70016, Telex: 933693 UKISI

Circle No. 229 on Readers' Service Card

American Association for the Advancement of Science

Science serves its readers as a forum for the presentation and discussion of important issues related to the advancement of science, including the presentation of minority or conflicting points of view, rather than by publishing only material on which a consensus has been reached. Accordingly, all articles published in *Science*—including editorials, news and comment, and book reviews—are signed and reflect the individual views of the authors and not official points of view adopted by the AAAS or the institutions with which the authors are affiliated.

Publisher: Alvin W. Trivelpiece

Editor: Daniel E. Koshland, Jr.

Deputy Editors: Philip H. Abelson (*Engineering and Applied Sciences*); John I. Brauman (*Physical Sciences*)

EDITORIAL STAFF

Managing Editor: Patricia A. Morgan

Assistant Managing Editor: Nancy J. Hartnagel

Senior Editors: Eleanor Butz, Ruth Kulstad

Associate Editors: Martha Collins, Barbara Jasny, Katrina L. Kelner, Edith Meyers, Phillip D. Szurmi, Kim D. Vandegriff, David F. Voss

Letters Editor: Christine Gilbert

Book Reviews: Katherine Livingston, *editor*; Deborah F. Washburn

This Week in Science: Ruth Levy Guyer

Contributing Editor: Lawrence I. Grossman

Chief Production Editor: Ellen E. Murphy

Editing Department: Lois Schmitt, *head*; Mary McDaniel, Barbara E. Patterson

Copy Desk: Lyle L. Green, Sharon Ryan, Beverly Shields, Anna Victorien

Production Manager: Karen Schools

Assistant Production Manager: James Landry

Graphics and Production: Holly Bishop, James J. Olivari, Eleanor Warner

Covers Editor: Grayce Finger

Manuscript Systems Analyst: William Carter

NEWS STAFF

News Editor: Barbara J. Culliton

News and Comment: Colin Norman, *deputy editor*; William Booth, Mark H. Crawford, Constance Holden, Eliot Marshall, Marjorie Sun, John Walsh

Research News: Roger Lewin, *deputy editor*; Deborah M. Barnes, Richard A. Kerr, Jean L. Marx, Leslie Roberts, M. Mitchell Waldrop

European Correspondent: David Dickson

BUSINESS STAFF

Associate Publisher: William M. Miller, III

Business Staff Manager: Deborah Rivera-Wienhold

Classified Advertising Supervisor: Karen Morgenstern

Membership Recruitment: Gwendolyn Huddle

Member and Subscription Records: Ann Ragland

Guide to Biotechnology Products and Instruments: Shauna S. Roberts

ADVERTISING REPRESENTATIVES

Director: Earl J. Scherago

Traffic Manager: Donna Rivera

Traffic Manager (Recruitment): Gwen Canter

Advertising Sales Manager: Richard L. Charles

Marketing Manager: Herbert L. Burkland

Sales: New York, NY 10036: J. Kevin Henebry, 1515 Broadway (212-730-1050); Scotch Plains, NJ 07076: C. Richard Callis, 12 Unami Lane (201-889-4873); Chicago, IL 60611: Jack Ryan, Room 2107, 919 N. Michigan Ave. (312-337-4973); San Jose, CA 95112: Bob Brindley, 310 S. 16 St. (408-998-4690); Dorset, VT 05251: Fred W. Dieffenbach, Kent Hill Rd. (802-867-5581); Damascus, MD 20872: Rick Sommer, 24808 Shrubbery Hill Ct. (301-972-9270); U.K., Europe: Nick Jones, +44(0647)52918; Telex 42513; FAX (0392) 31645.

Information for contributors appears on page XI of the 25 September 1987 issue. Editorial correspondence, including requests for permission to reprint and reprint orders, should be sent to 1333 H Street, NW, Washington, DC 20005. Telephone: 202-326-6500.

Advertising correspondence should be sent to Tenth Floor, 1515 Broadway, NY 10036. Telephone 212-730-1050 or WU Telex 968082 SCHERAGO.

Instrumentation and Equipment

The scientific enterprise continues to develop new, rewarding frontiers to explore. A key factor is exploitation of opportunities that are created by instrumentation and equipment. This issue of *Science* presents a sample of such opportunities.

Synchrotrons capable of producing ultraviolet radiation and x-rays have been in operation for some years. However, their numbers have increased both here and abroad. In addition their outputs of radiation have been greatly improved by installation of periodic magnetic insertion devices called wigglers and undulators. The flux of energy per unit area of synchrotron radiation exceeds that of older conventional sources by factors of 10^6 to 10^8 . An arbitrary portion of the spectrum can be sharply defined and used for particular experiments. An enormous number of studies that previously were impractical can now be performed. Gruner discusses the feasibility of performing time-resolved x-ray diffraction of biological materials. Many biological processes occur on a millisecond time scale, and these are now close to being accessible to investigation.

Prewitt and colleagues list a large number of experiments that have been performed on both inorganic and organic materials. As crystallographers they are pleased with an ability to determine structure of crystals having dimensions on the order of 10 micrometers. They also mention experiments that determine structure of surfaces both of inorganics and organics. In many—or even most—solids the surface structure differs from that of the bulk. Ability to obtain such information is of increasing potential economic importance. Molecular beam epitaxy is now producing a large number of new types of layered materials. Another application is in macromolecular crystallography. The small size of the x-ray beam and its collimation allow possible resolution of reflection maxima from viruses with unit cells as large as 1000 angstroms (mass of about 10,000,000 daltons).

One of the most exciting opportunities of our times is the sequencing of the human genome. It will be an enormous task. But it will be done. Were the techniques of yesteryear employed, the cost would be many billions of dollars, and many years—perhaps decades—would be required. However, human ingenuity is being applied to make the task easier and less costly. In this issue Gray and colleagues describe their work on chromosome purification and, in a research article, Prober and colleagues present a new system for rapid DNA sequencing. Chromosomes can be isolated from cells, stained with a DNA-specific fluorescent dye, classified by flow cytometry, and purified by flow sorting. With newly developed high-speed sorting, microgram quantities of some individual chromosomes having a purity of 90 percent can be isolated in a day or less. The authors suggest that the speed of the separation process might eventually be improved by an order of magnitude.

In the new DNA sequencing method, fluorescent molecules are added to DNA fragments by enzymatic chain extension reactions. The labeled fragments are separated by polyacrylamide gel electrophoresis and identified as they migrate past a fluorescent detection system. The wavelength of the fluorescent radiation differs for each terminal base. Sensitivity of detection is very good. After an initial period of electrophoresis the sequencer is capable of determining 50 bases per hour per lane. Twelve lanes can be used.

Proteins for therapeutic applications are major crude products of recombinant DNA technology. To be safe for human use, the proteins must be highly purified. A major means of achieving this is the use of liquid chromatography. The crucial factors in such separations are interactions between the proteins and the column packing. Regnier has analyzed the interactions that occur when various types of packing are used. He points out that only a fraction of those amino acid residues at or near the exterior surface of a three-dimensional protein can interact with a particular type of chromatographic matrix.

Electroencephalography is now being supplemented by an increasingly useful magnetoencephalography (MEG). Rose, Smith, and Sato describe use of MEG in determining the location of epileptic discharges. Sensitive superconducting quantum interference device (SQUID) sensors can detect the tiny magnetic fields that are produced when the electric discharges occur. With information derived from a number of detectors the location of the source of the discharge can be determined. The procedure holds promise of becoming a valuable noninvasive diagnostic aid.—PHILIP H. ABELSON


Sample preparation

Try the new Amprep™ minicolumns from Amersham.

● Easy to use


A full range of accessories makes it easy for you to prepare:

24 samples



Use the Super Separator-24

10 samples



Use the Manifold-10


1 sample



Use a syringe and free adaptor in each pack of columns

● The best sorbent for your analyte

You may already be using a C18 column. This is a good general purpose sorbent, but we have shown that other sorbents can produce better results for many analytes. That is why there are 9 sorbents in the Amprep range.



● Fully developed protocols

We have developed protocols for many analytes. They are now available to help you get the best sample preparation using Amprep minicolumns.

For further information contact your nearest Amersham Office

Amersham International plc
Amersham UK

Amersham

UK Sales Aylesbury (0296) 395222 Amersham Australia PTY Ltd Sydney 888-2288 Amersham Belgium SA/NV Brussels (2) 538 9194
Amersham Buchler GmbH & Co KG Braunschweig West Germany (05307) 8080 Amersham Canada Ltd Ontario (416) 842 2720
Amersham Corporation Arlington Heights USA (312) 364 7100 Amersham Denmark ApS Birkerød (02) 82 0222 Amersham France SA Paris (1) 69.28.83.00
Amersham Japan Tokyo (03) 816 6161 Amersham Nederland BV Houten 03403 76660 Amersham Norway Sandvika (02) 54 63 18
Amersham Sweden AB Solna (8) 734 0800

Circle No. 235 on Readers' Service Card

5. T. Holmes and R. Rahe, *J. Psychosom. Res.* 11, 213 (1967); J. G. Rabkin and E. L. Struening, *Science* 194, 1013 (1976).
6. R. S. Lazarus, in *Models for Clinical Psychopathology*, C. Eisendorf, D. Cohen, A. Kleinman, P. Maxim, Eds. (SP Medical and Scientific Books, New York, 1981), pp. 177-214; G. Cohen, G. W. Evans, D. Stokols, D. S. Krantz, Eds., *Behavior, Health, and Environmental Stress* (Plenum, New York, 1986).
7. J. Wolpe, *The Practice of Behavior Therapy* (Pergamon, New York, 1969), M. R. Goldfried, in *Behavior Change Through Self-Control*, M. R. Goldfried and M. Merbaum, Eds. (Holt, Reinhart, Winston, New York, 1973); D. H. Shapiro and S. M. Zifferblatt, *Am. Psychol.* 31, 7 (1976).
8. R. N. Walsh and D. H. Shapiro, Eds., *Beyond Health and Normality: Explorations of Exceptional Psychological Wellbeing* (Van Nostrand/Reinhold, New York, 1983); M. N. Palsane, S. N. Bhavsar, R. P. Goswami, G. W. Evans, *Indian J. Psychol.* 5, 1 (1986).

Response: Shapiro, Evans, and Shapiro point out complexity in the construct of control, correctly noting that there are different types of control, and they argue for conceptual and empirical rigor. While there is undoubtedly need for clarity, I am not confident that the several distinctions among types of control that Shapiro *et al.* propose have been shown to have significant heuristic value in generating research questions. One must be careful not to overemphasize the importance of differentiation of terms and concepts when indeed there may be fewer, rather than more, underlying constructs in this area. This remains an empirical question, however, and one that most urgently needs to be addressed.

Shapiro *et al.* comment on the distinction between control-enhancing interventions offered by the environment and self-control strategies. They discuss impairment of self-control as an essential feature of many clinical problems, for example, obesity, bulimia, and alcoholism, and cite weak results from the use of self-management strategies in these areas. These data are then used to imply that self-control interventions may not work. One must be careful, however, to separate studies of clinical populations from studies of normative samples, for example, the aged or children, who may suffer an impairment in control because of developmental stage and environmental change. In addition, many clinical disorders have been intractable after self-control interventions, not because of problems with self-management as an intervention strategy, but because they often have a large genetic contribution and involve a heavy burden of biological change once the disorder is initiated (1). Control-relevant intervention may not work where biological and genetic factors influence the disorder. Indeed, I have argued that teaching self-management strategies in these domains can convey an implicit message of personal blame for the cause of the disorder, leading to feelings of shame and reduced ability to exercise control (2).

The study of control in human populations is an exciting and timely one, especially with increasing demonstrations of potential health significance. Like Shapiro *et al.*, I believe strongly that the time has come to understand the underlying similarities and dissimilarities among the various constructs that have been used in the control literature. These are not context-free evaluations, however; setting accounts for a substantial portion of the variance when studying the construct of control. As Bandura (3) has suggested, individuals with a high degree of self-efficacy still can recognize when there are no response-outcome contingencies, that is when events in the environment are uncontrollable by anyone, despite the individual's own sense of personal mastery.

JUDITH RODIN
Department of Psychology,
Yale University,
New Haven, CT 06520

REFERENCES

1. J. Rodin, paper presented at the Master Lecture Series, American Psychological Association, Washington, DC, 1977; *Pers. Soc. Psychol. Bull.* 7, 147 (1981); *Health Psychol.* 4, 1 (1985).
2. J. Rodin, in *Human Services in the 1980's: Social Science Perspectives*, A. Johnson, O. Grunsky, B. Raven, Eds. (Auburn House, Boston, MA, 1982).
3. A. Bandura, *Psychol. Rev.* 84, 2 (1977).

Comparison of High Schools

Luther B. Otto's review (11 Sept., p. 1357) of *Public and Private High Schools: The Impact of Communities* by James S. Coleman and Thomas Hoffer does not stress a problem presented by overall comparisons of public with private and Roman Catholic high schools. It is easy to make a case for private and Roman Catholic schools if one does not take into consideration the superb records of selected public schools.


The public high schools welcome students of all sorts to the extent that the students can meet current standards of admission. Public schools are required to be as inclusive as possible. The students in public high schools in such choice areas as Short Hills or Madison, New Jersey, however, have records of accomplishment that go well beyond those of students in the local Roman Catholic schools.

To compare these school systems without reference to class and area is like treating statistics on elephants, dogs, and cats as comparable.

ALFRED MCCLUNG LEE
Graduate School and University Center,
City University of New York,
New York, NY 10036-8099 and
Drew University, Madison, NJ 07940

Tiny Flow Metering Pumps With AC Motor Drive

Those handy little RH pump heads for piddling flow rates now come packaged with their own synchronous AC motor drive. Each has an RH miniature pump head and totally enclosed three-speed motor drive assembly. This provides a choice of 150, 300 or 600 piston strokes per minute by means of a safe and simple belt change arrangement. Motor can be stalled indefinitely without overheating!



Pumps plug into standard 110-volt, 60-Hz three-prong sockets. Are steplessly adjustable even while running from zero to maximum flow rate by means of calibrated control ring. Are reversible while operating. And have the exceptional chemical resistance of alumina ceramic and fluorocarbon in the wetted stream. Compression fittings at the ports accept 1/4 inch O.D. standard laboratory tubing. Pressures to 50 psi at low flows, to 20 psi at any flow.

Repeat accuracy to better than 1%!

Model RHSY0CKC \$400
Max flow rate: 7.5 ml/min at 150 SPM,
15 ml/min at 300 SPM,
30 ml/min at 600 SPM

Model RHSY1CKC \$400
Max flow rate: 15 ml/min at 150 SPM,
30 ml/min at 300 SPM,
60 ml/min at 600 SPM

Available from stock for immediate shipment. Write or call for Catalog RP401—the FMI LAB PUMP STORY.



CALL
TOLL-FREE
800-223-3388
IN NY STATE
516-922-6050

FMI FLUID METERING, INC.
P.O. Box 179, 29 Orchard Street
Oyster Bay, New York 11771
Telex 5101001281 / FAX 516-624-8261

An open & shut case for the compact new HBI Micro Centrifuge. You'll love it.

Open, you'll love how easy it is to load, empty or clean its top-mounted, autoclavable rotor; far easier than rotors in wells, as on other microcentrifuges.

Closed, you'll love how smoothly and quietly it spins, how quickly it accelerates (13,000 max. rpm/13,400 x g in just 10 sec.), and how incredibly compact and light it is (you can easily carry it into a cold room, another lab, or anywhere you need it).

You'll also love its convenient features: dual speed selector, safety lid lock, 10-min. timer, pulse switch and open lid light. Interchangeable adaptor plates accept 12 x 1.5, 14 x 0.75 and 18 x 0.4 mL micro test tubes and B-D Microtainer* Tubes.



Last but not least, you'll love its price. Need more information? Call us toll free: 1-800-631-1369 (in NJ, 201-843-2320) or write: HBI, 244 Saddle River Road, Saddle Brook, NJ 07662-6001.

*Microtainer is a trademark of Becton Dickinson and Company.

For unique laboratory equipment
HAAKE BUCHLER INSTRUMENTS, INC.



Circle No. 163 on Readers' Service Card



Meet System Gold™ The Personal Chromatograph, new from Beckman. It's very special. It's easy to get to know. And it's the one you've been waiting for.

System Gold has all the highly capable modules essential for HPLC. But now they're linked together on a digital communications network. So you control them all from one place. Either the portable NEC lap-top computer. Or the powerful IBM PC-AT or -XT. In fact, with the latter's convenient mouse, you just point and click your way to more productivity.

That's Personal™ Chromatography.

Solvent Delivery. Performance, Precision and Reliability.

From 60 mL/minute to just 1 μ L/minute, System Gold offers you single-piston technology you can count on.

Single and binary modules provide a uniform and pulse-free flow. With reproducibility that's better than $\pm 0.1\%$.

What's more, they offer stand-alone programmability, a unique piston wash, high pressure dynamic mixing, and automated multi-solvent selection. In other words, System Gold solvent delivery modules think of everything.

Detection That's Sensitive and Versatile.

All-digital System Gold detectors give you nothing less than flawless and consistently reliable data. Featuring sharp, undistorted peaks. Low baseline noise. And variable wavelength monitoring.

Plus, real-time high sensitivity scanning provides data you just can't get with diode array detectors. Programmability optimizes results automatically. And, to complete the picture, there's a full range of flow cells, from micro to macro.

Personal Chromatography Makes It Easy.

System Gold can be configured for any application. Easily. From methods development and quality control, to the analysis of amino acids, proteins, peptides and more.

Your System Gold LC could be up and running in just one hour. Installation is so easy you can do it yourself. And with electronic self-diagnostics, maintenance is a snap.

Go For The Gold.

Enough said. The best way to get to know the all-new System Gold series of Personal Chromatographs is to see them for yourself.

Just contact your local Beckman HPLC representative. In the U.S. call toll-free 800/742-2345. Or write Beckman Instruments, Inc., Altex Division, 2350 Camino Ramon, P.O. Box 5101, San Ramon, CA 94583. Offices in major cities worldwide.

Go for System Gold today.

Circle No. 204 on Readers' Service Card

BECKMAN

A SMITHKLINE BECKMAN COMPANY

Australia, Sydney (02) 816-5288. Austria, Vienna (0222) 32 25 57. Canada, Toronto (416) 673-9844. Denmark, Birkerød 2-81 92 11. France, Gagny (1) 381 93 00. Germany, München (089) 38871. Hong Kong, Aberdeen 5-539236. Italy, Milan 02/69911. Japan, Tokyo (03) 221-5931. Mexico, Mexico City 546-8109. Netherlands, Mijdrecht 02979-3651. Norway, Oslo 02-648040. Puerto Rico, Carolina (809) 762-3030. Singapore 5342508. Spain, Madrid (091) 723 16 66. Sweden, Bromma 08 98 53 20. Switzerland, Geneva (022) 35 64 80. Taiwan, Taipei (02) 700-2230. United Kingdom, High Wycombe (0494) 41181. IBM PC-AT, XT are trademarks of International Business Machines Corporation. NEC is a trademark of NEC Corporation.
©1987 Beckman Instruments, Inc. AX87-1008 A



INTRODUCING HPLC WITH A PERSONAL TOUCH.



DNA, FPLC.[®] PURE AND SIMPLE





FPLC improves plasmid purification, with high recovery and between 3 to 6-fold time savings. Without the high costs and concerns of a centrifuge, Superose[®] 6 prep grade purifies pBR322 in just two hours.

FPLC clearly resolves restriction fragments of up to 3000 base pairs with sample loads of up to 100 μ g. Hae III digest products of pBR322, including 123 and 124 base pairs fragments, have been separated to base line resolution on the ion exchanger Mono Q.[®]


FPLC excels at rapid, high resolution purification of synthetic oligonucleotides. Even large oligonucleotides of 100-mer or more can be purified with high recovery and greater than 98% purity by using Mono Q and ProRPC.[™]

To serve all your separation needs, high performance chromatography with FPLC now adds DNA to the already impressive list of proteins, peptides, monoclonal antibodies, and amino acids.


Request more information on DNA separated with FPLC...pure and simple.



Preparative purification of the plasmid pBR322.



Three-step purification of 100-mer oligonucleotide. ProRPC, Mono Q and ProRPC again.





Laboratory Separation Division
Piscataway, New Jersey 08854
Information: (800) 526-3618
In NJ: (201) 457-8000

Circle No. 190 on Readers' Service Card

PHARMACIA PRODUCTS WILL BE ON DISPLAY AT THE ASCB MEETING — BOOTHS 621, 623, 625

If safety evaluation is your responsibility, we're here to help you get the data you need.



Our Acute Toxicology Department specializes in short term toxicological studies, conducting as many as 800 each year.

We supply the data for such things as MSDS sheets, PMN's, FIFRA registration, product labeling, employee protection, OSHA compliance and DOT classifications. Our capabilities also include preliminary screening studies prior to initiation of comprehensive programs, and pre-market testing of a wide range of chemicals, pharmaceuticals and consumer products.

We routinely conduct acute toxicity, irritation and sensitization studies — and a full range of subchronic, neurotoxicity and phototoxicity studies.

And when it comes to people, we have more of the best. Credentials of our study directors, supervisory and technical staff include Society of Toxicology membership and certification by AALAS and the American Board of Toxicology.

So when a short term study is called for, contact us. We have the expertise to respond quickly, and provide accurate results on schedule.


Acute Toxicology Testing.
We do it right the first time.
On time.

Bio/dynamics, Inc. 
Results you can trust.

Mettlers Road, East Millstone, New Jersey 08873 • (201) 873-2550

TOXICOLOGY • INHALATION TOXICOLOGY • PATHOLOGY • METABOLIC AND ANALYTICAL CHEMISTRY • QUALITY ASSURANCE • FACILITIES MANAGEMENT

Circle No. 157 on Readers' Service Card



We don't
cut
corners.

***Baker's radius-cornered
one-piece interior
maximizes personnel
protection.***

Caution: inadequate containment of dangerous biological materials can be hazardous to your health.


The best way to reduce the risk? Baker Biological Safety Cabinets. With the best materials. The best workmanship. The best performance.

One-piece stainless steel and true metal rounded radiused corners make our interior work areas leak-tight and easy to clean. Other cabinets are *seamed* with metallic-colored caulking compound, which can shed, peel, and allow harmful contaminants and bacteria to build up.

At The Baker Company, we never compromise on personnel or product protection. Should you?

Call toll free for complete data. Or write to The Baker Company, P.O. Drawer E, Sanford, ME 04073.

1-800-992-2537



**the BAKER
company**



Cappel

Commitment...

For over forty years, Cappel has been dedicated to supplying high quality immunochemicals to the research community. Cappel continues to be the leading supplier of fluorochrome conjugates, enzyme conjugates, affinity purified antisera, antibody fractions, antibody fragments, immunocytochemical reagents, and animal sera.

The Cappel team total effort is directed towards supplying you with the products *you* need to perform *your* research. We routinely stock over 1,000 products to assure that your special requirements are met. Need something unique? Our Custom Service team will work with you to develop an antiserum production program to meet your laboratory's research reagent requirements.

Cappel continues the commitment to provide you the finest in immunoreagents for research use.

Call us today.

Organon Teknika—Cappel

One Technology Court
Malvern, PA 19355
Telephone: (215) 251-2000
(800) 523-7620
(800) 662-2440 in PA

Telex: 685-1293 CBI UW

Fax: (215) 251-2208

Outside the U.S., contact:

Organon Teknika N.V.

Veedijk 58-2300 Turnhout
Belgium
Telephone: 014/41.04.61
Telex: 71939 Otbel
Fax: 3214413893

ORGANON TEKNIKA

TOOLS FOR RESEARCH.

Labeled compounds from Merck Sharp & Dohme / Isotopes.



MSD ISOTOPES – research products that work, bringing results to researchers in biology, chemistry, physics, medicine and related fields.

Did you know that we now offer thousands of compounds labeled with deuterium, carbon 13, nitrogen 15 and other stable isotopes too? And did you know that even if you require a compound which is not available from us immediately 'off-the-shelf', we have the world's best facilities to custom synthesize it for you?

Fact is, that today more scientists around the world are getting more results than ever before with MSD ISOTOPES

'research tools' because of our unprecedented variety, quality and value.

So just call or write to us. Tell us what you need, (or even what you *think* you might need) and we will be happy to send you up to date information on the price and availability of your compound.

MSD ISOTOPES – Labeled compounds are our only business.

MSD

MSD ISOTOPES

DIVISION OF MERCK FROSST CANADA INC., Montreal, Canada

WEST COAST

P.O. Box 2951
Terminal Annex
Los Angeles, CA 90051
Los Angeles Area: (213) 723-9521
Outside State of CA: (800) 423-4977
State of CA: (800) 372-6454


EAST COAST & CENTRAL

4545 Oleatha Avenue
St. Louis, MO 63116
Outside state of MO: (800) 325-9034
State of Missouri: (314) 353-7000
Circle No. 257 on Readers' Service Card

CANADA

P.O. Box 899
Pointe Claire/Dorval, Quebec
Canada H9R 4P7
Telephone: (514) 697-2823

T H E D E C I S I O N




**Before you trust your life's work
to an ultra-low temperature freezer,
talk to someone who puts his trust in Harris,
and his reputation on the line.**

Talk to your independent Harris Sales Representative.

Because he operates his own company, it's his business to sell the best. He stakes his reputation on the quality and performance of each ultra-low temperature freezer he sells.

Which is why he sells Harris. A complete line of chest and upright freezers from -40°C to -150°C . A variety of models to suit every lab budget, every application. Each freezer built tough, backed by the best warranty in the industry and supported by more than 200 factory trained service technicians located throughout North America.

**Before you make your decision, talk to someone who deals from
experience and shares the risk, hands-on.
Your local Harris Sales Representative.**




Call us toll-free
and we'll put you in touch.

1-800-221-4201

Harris

Harris Manufacturing Company
P.O. Box 38, North Billerica, MA 01862
Phone 1-800-221-4201 Telex 577449

Circle No. 203 on Readers' Service Card



Your pad or ours?


If you perform calculations, the answer is obvious.

MathCAD 2.0.

It's everything you appreciate about working on a scratchpad—simple, free-form math—and more. More speed. More accuracy. More flexibility.

Just define your variables and enter your formulas anywhere on the screen. MathCAD formats your equations as they're typed. Instantly calculates the results. And displays them exactly as you're used to seeing them—in real math notation, as numbers, tables or graphs.

MathCAD is more than an equation solver. Like a scratchpad, it allows you to add



text anywhere to support your work, and see and record every step. You can try an unlimited number of what-ifs. And print your entire calculation as an integrated document that anyone can understand.

Plus, MathCAD is loaded with powerful

built-in features. In addition to the usual trigonometric and exponential functions, it includes built-in statistical functions, cubic splines, Fourier transforms, and more. It also handles complex numbers and unit conversions in a completely transparent way.

Yet, MathCAD is so easy to learn, you'll be using its full power an hour after you begin.

What more could you ask for? How about the exciting new features we've just added to MathCAD 2.0...

- Built-in equation solver
- Full matrix operations
- Two to four times increase in calculating speed
- Easier full-page text processing
- Auto-scaled plots
- Memory enhancements
- Additional printer and plotter support
- And more.

If you're tired of doing calculations by hand or writing and debugging programs, come on over to our pad. MathCAD. The Electronic Scratchpad.

Call for a detailed spec sheet and the name of a MathCAD dealer near you.

1-800-MathCAD (In MA: 617-577-1017).

Circle No. 213 on Readers' Service Card

Requires IBM PC® or compatible, 512KB RAM, graphics card.
IBM PC® International Business Machines Corporation.
MathCAD® MathSoft, Inc.

MathCAD®

MathSoft, Inc., One Kendall Sq., Cambridge, MA 02139

Automatic DNA sequence analysis

Introducing the new Gene-Master™ system, the first automated film reading and analysis system for DNA sequencing.

It's the only stand alone, totally integrated DNA sequence reader and analysis system. It features specially designed hardware and software for scanning autoradiograms in a matter of minutes! The scanner incorporates a state of the art CCD camera. The proprietary pattern recognition software locates the lanes, corrects for artifacts using expert rules, and assigns bases. All of this without operator intervention.

The Gene-Master system goes beyond autoreading. After the sequences are determined, the powerful Staden shotgun handler assembles the readings into a consensus. Then, additional Gene-Master software can analyze the sequences for protein coding regions and other biologically important signals.

Six reasons to automate with the Gene-Master system

- **Higher throughput.** Approximately 100,000 bases per day can be processed.
- **Superior accuracy and sensitivity.** The Gene-Master system combines the latest technology with superior chemistry. Accuracies of 99% have been achieved!
- **Original data retained.** Since the original autoradiogram is always available, you can refer to it or the stored image for verification.
- **Extremely low cost-per-base.** Because the Gene-Master system offers such high throughput, sequence reading costs less than a penny per base.



• **Use existing equipment and techniques.** Your existing power supplies, sequencing cells, and chemicals can be used with the Gene-Master system.

• **Comprehensive DNA and protein analysis software.** The Gene-Master system includes a selection of updatable software. Unique multi-tasking capability runs up to three analyses simultaneously, finding alignments, hairpin loops, and other signals in DNA. You can predict protein structure, hydroplots, and more. All results are displayed in high resolution graphics.

For those who prefer semi-automatic operation, Bio-Rad offers a digitizer based system which is upgradable to a fully automated system.

Sign up for a seminar in your area and find out how autoreading and analysis can give you the productivity and sensitivity edge in DNA sequencing. Call us directly at 800 4 BIORAD. For more information, circle the reader service number.

1. Elder, J.K., Green, D.K., and Southern, E.M. (1986). *Nuc. Acids Res.* 14:417-424.



Symposia on Automatic
Film Reading

BIO-RAD


**Chemical
Division**

1414 Harbour Way South
Richmond, CA 94804
(415) 232-7000
800-4-BIO-RAD

Also in Rockville Centre, NY; Hornsby, Australia; Vienna,
Austria; Mississauga, Canada; Watford, England; Munich,
Germany; Hong Kong; Milan, Italy; Tokyo, Japan; Utrecht,
The Netherlands; and Glattbrugg, Switzerland.

Circle No. 206 on Readers' Service Card

Answers.



Extend The Limits of Sequencing with . . .

SEQUENASE™

Reduces Gel Background
Gives Uniform Band Intensities
Eliminates Compressions Using dITP
Resolves 500, 600 or More Bases
Easy To Use

*Sequenase™ is a novel T7 DNA polymerase specifically engineered for DNA sequencing¹. Sequenase™ has been overproduced and purified. It is highly processive; it will polymerize thousands of nucleotides uninterrupted. It is not impeded by template secondary structure. It has no exonuclease activity, and efficiently incorporates nucleotide analogs (thio-dATP, dideoxynucleotides, dITP, etc.). Together, these properties make Sequenase™ ideal for determining the sequence of even "difficult" DNAs, and allow reading of sequence information well beyond the point where traditional sequencing enzymes fail.

Sequenase™ is used in combination with familiar M13 vectors and primers and either ³⁵S or ³²P labeled dATP. The products can be run on any standard sequencing gel apparatus with excellent results. However, it is recommended that very long or gradient gels be used for reading the maximum possible number of nucleotides. Specially prepared high purity reagents are included to insure optimal results.

The convenient Sequenase Kit includes:

- | | |
|--|-----------------------------|
| ■ Sufficient for 100 sequences | ■ Control DNA |
| ■ Buffer | ■ Stop solution |
| ■ Pre-mixed nucleotides including dITP | ■ Sequenase DNA Polymerase |
| ■ Universal Primer | ■ Complete protocol booklet |

Sequenase™ Kit, Product #70700; Enzyme only, Product #70722.

For a free copy of our Molecular Biology catalog, write USB at P.O. Box 22400, Cleveland, Ohio, 44122. Phone 800/321-9322. From Ohio or Canada, call collect to 216/765-5000. Telex: 980718. FAX: 216/464-5075.

*Patent Pending

1. Tabor, S. and Richardson, C.C., *Proc. Nat. Acad. Sci., USA* 84, No. 14, 4767-4771 (1987).

INTERNATIONAL DISTRIBUTORS:

CAMBRIDGE BIOSCIENCE
CAMBRIDGE, ENGLAND
223 316 855


RENNER GMBH
D-6701 DANNSTADT-SCHAUERNHELM 1,
GERMANY
0 62 31/44 47 oder 45 24

KEBO-AB
SPANGA, SWEDEN
46-8-760 00 20

TOUZART & MATIGNON
PARIS, FRANCE
46-80-85-21

LUCERNA CHEM AG
LUCERNE, SWITZERLAND
041-44 93 44

TOYOBO COMPANY, LTD.
OSAKA, JAPAN
06-348-3786



UNITED STATES BIOCHEMICAL CORPORATION

Circle No. 175 on Readers' Service Card

New International Journals from Academic Press

Brain, Behavior, and Immunity

Concerned with the interactions between the nervous system and the immune system at the molecular, cellular, and organismic levels

EDITOR-

IN-CHIEF: **Robert Ader**

ASSOCIATE

EDITORS: **Nicholas Cohen and David Felten**
*University of Rochester School of Medicine
and Dentistry, New York*

Volume 2 (1988), 4 issues ISSN 0889-1591
In the U.S.A. and Canada: \$100.00
All other countries: \$112.00

Genomics

International Journal of Gene Mapping and Nucleotide Sequencing Emphasizing Analyses of the Human and Other Complex Genomes

EDITORS-

IN-CHIEF: **Victor A. McKusick**

*The Johns Hopkins University, Baltimore,
Maryland*

Frank H. Ruddle

Yale University, New Haven, Connecticut

Volumes 2-3 (1988), 8 issues ISSN 0888-7543
In the U.S.A. and Canada: \$160.00
All other countries: \$176.00

Journal of Food Composition and Analysis

An official publication of INFOODS, the International Network of Food Data Systems, a UNU Organization

Devoted to all scientific aspects of the data on chemical composition of human foods

EDITOR: **Kent K. Stewart**

*Virginia Polytechnic Institute and State
University, Blacksburg, Virginia, U.S.A.*

Volume 1 (1987/1988), 4 issues ISSN 0889-1575
In the U.S.A. and Canada: \$70.00
All other countries: \$82.00

Journal of the Japanese and International Economies

Academic analyses of the Japanese economy and its interdependence with other national economies

EDITORIAL

BOARD:

Masahiko Aoki

University of Kyoto and Stanford University

Michihiro Ohyama, *Keio University*

Koichi Hamada, *Yale University*

Masahiro Okuno-Fujiwara

University of Tokyo

Volume 2 (1988), 4 issues ISSN 0889-1583
In the U.S.A. and Canada: \$100.00
All other countries: \$112.00

NEW IN 1988

Biomedical and Environmental Sciences

EDITORS-

IN-CHIEF:

Frederick Coulston

*White Sands Research Center, Albany,
New York*

Chen Chunming

*Chinese Academy of Preventive Medicine,
Beijing, China*

MANAGING

EDITOR:

Frank C. Lu, *Miami, Florida*

Volume 1 (1988), 4 issues
In the U.S.A. and Canada: \$100.00
All other countries: \$112.00

Journal of X-Ray Science and Technology

EDITOR:

Larry V. Knight

Brigham Young University, Provo, Utah

Volume 1 (1988), 2 issues
In the U.S.A. and Canada: \$50.00
All other countries: \$56.00

All prices are in U.S. dollars and are subject to change without notice.

For sample copies or more information, please write or call:




ACADEMIC PRESS, INC., Journal Promotion Department
1250 Sixth Avenue, San Diego, CA 92101, U.S.A. (619) 230-1840

Please visit the Academic Press booth at the American Society for Cell Biology meeting

S7185

Realize the potential
of your cell line.



Innovators in the many facets
of animal cell culture, Bio-Response advances
your product's commercial success.



FOR MORE INFORMATION ON OUR CONTRACT PRODUCTION SERVICES, CALL (415) 786-9744. BIO-RESPONSE INC., 1978 WINTON AVENUE, HAYWARD, CA 94545

Circle No. 167 on Readers' Service Card

The easiest and fastest way to synthesize peptides also happens to be the most affordable.

Meet the MilliGen 9000 Series.

Now you can realize all the advantages of the latest Fmoc synthesis methods coupled with the convenience of programmable automation. At a price that puts peptide synthesis within the reach of even the smaller facility.

Optimized Fmoc chemistry simplifies synthesis and eliminates HF. In contrast to earlier batchwise techniques which rely on the differential acid lability of protecting groups, the new continuous flow MilliGen instruments use a mild orthogonal protecting scheme. Preactivated amino acid esters simplify protocols and instrumentation. Deprotection with mild base solutions replaces potentially damaging repetitive acid treatments. And final side chain deprotection and peptide-resin cleavage is achieved with TFA, rather than hazardous HF. A choice of peptide linkage agents lets you easily generate free peptide acids, amides or fully protected peptides.

The faster route to purer peptides.

MilliGen's combination of straightforward, reliable Waters™ fluidics and innovative continuous flow methodology delivers high yields of high quality crude peptides. For standard protocols, total cycle times are less than an hour, and no double coupling is required. Post-synthesis procedures are fast, simple and safe. Syntheses are easily scaleable from 0.1 to 5.0 mmol, just by increasing reaction column diameter and pump flow rates.

Continuous flow operation minimizes reagent usage and greatly simplifies and speeds up production of peptide analogues. In addition, a built-in UV spectrometer provides a

real-time monitor of system performance. An integral IBM compatible computer works with MilliGen Express-Peptide™ software to take the guesswork out of reagent calculation, to control instrument operation and to give you complete documentation of your synthesis.


Select the MilliGen instrument that meets your needs and budget.

The price of the semi-automatic MilliGen 9020 is \$35,000; the complete, fully automated MilliGen 9050, \$55,000. All of our modular hardware and software can be easily and economically upgraded to meet your lab's changing needs.

To find out more, circle the Reader Service Card. Or call toll-free (800) 872-0071; in Massachusetts, (617) 275-5208.



MilliGen
Division of MILLIPORE



The new Eppendorf Micro Centrifuge.

*With 50% higher capacity,
variable speed,
quieter operation, and
quick-release rotor.*

Brand new and turning 18.

Higher capacity...plus.


The new 18-place Model 5415 Micro Centrifuge gives you important operating advantages—with unique Eppendorf quality.

Versatile in use.


Model 5415 has a variable-speed motor that reaches a maximum of 14,000 rpm with an RCF of 16,000 x g; a 30-minute timer; and a momentary button for short spins. It accepts 1.5 mL, 500 μ L, 400 μ L, and 250 μ L Eppendorf Microcentrifuge Tubes and blood collection microtubes, such as B-D Microtainer* Tubes.

New rotor design.

The enclosed rotor design reduces air turbulence for quieter operation. And the new quick-release feature lets you transport the rotor *with* tubes—especially convenient when the centrifuge is run in a cold room.



Enclosed rotor design reduces air turbulence and noise. Tubes are angled precisely at 45° to maximize pellet formation.



Quick-release feature allows the 18-position rotor to be easily transported even when loaded.

Microtainer Tubes is a registered trademark of Becton Dickinson and Company.

Safe and rugged.

The Eppendorf 5415 Micro Centrifuge is UL listed for safety. It's so rugged that an accidentally unbalanced load won't cause excessive vibration or motor damage.

For more information: call 800-645-3050; in New York, 516-334-7500. Or write Brinkmann Instruments, Inc., Cantiague Road, Westbury, NY 11590. (In Canada: 416-675-7911; 50 Galaxy Blvd., Rexdale, Ont. M9W 4Y5)

Specifications


Maximum speed	14,000 rpm
Maximum RCF	16,000 x g
Test-tube capacity	18
Time required for maximum speed	10 sec
Time required to stop	12 sec
Dimensions (L x W x H)	28 x 21 x 28.5 cm

For literature circle reader service number 165
For a demonstration circle reader service number 166

BRK-5458

eppendorf

Shaping the future. **Brinkmann**
INSTRUMENTS, INC.




THE IBI GOLD STANDARD

IBI is proud to introduce the first **BLOT QUALIFIED** line of **GENOMIC GRADE™** restriction enzymes specifically designed for **RFLP**, **gene mapping** and **OFAGE/FIGE** applications. In addition to our normal stringent quality control procedures, **each lot is assayed on high molecular weight GENOMIC DNA**. Samples prepared from both **peripheral blood and cultured cells** are digested to completion, transferred to a nylon membrane and subsequently probed. By incorporating into our Quality Control assay procedures the same substrate used in your own experiments, **GENOMIC DNA**, many commonly observed artifacts related to genomic analysis are eliminated. These unique functional assays guarantee:

■ ***STRONG UNIT EQUIVALENTS*** — The number of units required for complete digestion of DNA derived from both peripheral blood and cultured cells is reported for each lot of Genomic Grade enzymes. By taking into account due to substrate specificity (GENOMIC DNA vs. lambda), or the protocols employed in sample preparation (blood vs. cells), the likelihood of a **partial digestion in your hands is minimized**.

■ ***OPTIMAL SIGNAL-TO-NOISE RATIO*** — Minimum acceptable signal-to-noise ratio is determined by laser densitometry of blot autoradiograms. High levels of background or smeared lanes can obscure complex or faint patterns making data difficult or impossible to interpret. Each lot of Genomic Grade enzyme has passed this sensitive functional assay.




DNA was prepared from normal peripheral blood leukocytes, digested to completion with *Pvu* II (Genomic Grade) and transferred from a 1.5% agarose gel to a nylon filter membrane. A 1.2 Kb DR β cDNA clone contained in the *Pst* I site of pBR322 was used as the probe. Densitometry values are presented as ratios in the following format: filter background/lane background/band signal or, in this example, an average signal-to-noise ratio of 30:1.

Quality Biologicals from ...

Life on the Edge.

■ **NO DETECTABLE SPECIFIC ENDONUCLEASE CONTAMINATION** — The level of detection of extraneous fragments is increased **1,000-fold** under blotting conditions as compared to sensitivities achievable by visual examination of an ethidium bromide stained agarose gel (5 picograms as compared to 5 nanograms). These dramatically improved levels of sensitivity are a prerequisite if easily misinterpretable artifacts are to be avoided.

■ **NO DETECTABLE PLASMID CONTAMINATION** — Only intact-vector probes are used in IBI's Genomic Grade blotting assays to ensure that even trace amounts of vector DNA contaminants in the enzyme preparation are readily detected. This more closely resembles the conditions in your own lab where gel purified fragments invariably contain some vector material. The level of sensitivity exceeds 5 picograms plasmid DNA per lane, eliminating the possibility of false positives.



DNA was prepared from three different cultured β -lymphoblastoid cell lines, digested to completion with *Rsa* I (Genomic Grade) and transferred to a nylon filter membrane. A DR β clone contained in the plasmid pBR322 was then nick-translated and used as the probe. DR β was employed because of its highly polymorphic nature. As is evident in the autoradiogram, no signals are present in the molecular weight range >2.3 Kb, where one would expect to detect plasmid contamination.

Has your ability to interpret your southern blots become more of an art than a science? Protect the integrity of your precious samples, as well as your interpretations. **The highest level of confidence is assured with IBI Genomic Grade Enzymes. FUNCTIONAL TESTING** guarantees uncompromised enzyme performance and reliability in this specialized area of application.

GENOMIC GRADE™ restriction enzymes — ANOTHER IBI FIRST!

To find out more about this exciting new line of **GENOMIC GRADE** restriction enzymes please call 1-800-243-2555 today. International inquiries, call your local IBI affiliate or IBI at 203-786-5620. TELEX 643993. IBI CANADA, call TOLL-FREE 800-387-1525.

The Cutting Edge of Molecular Biology

Circle No. 238 on Readers' Service Card

IBI

Never before has a dual wavelength spectrophotometer done so much

Whether your samples are clear or turbid, the DW-2000™ UV/VIS spectrophotometer is the easiest, most precise way to perform simple or complex analyses. It combines the double beam mode of operation for routine samples with a dual wavelength capability for handling complex analyses.


It scans the entire UV/VIS range with automatic calibration and baseline correction.

Data is collected, stored in memory and displayed on the CRT with a speed that will startle you. Recall data from one of the eight storage locations, and perform a variety of post processing functions including cursor control and peak pick, scale expansion, difference or derivative spectra.

By combining fast and slow kinetics with wavelength scanning it's everything scientists have been asking for in an economical multi-user UV/VIS spectrophotometer.



for so many



In these days of tightening budgets, the DW-2000 is the *departmental* instrument of choice for easy operation and unparalleled results.

For example, in petrochemical laboratories, analyzing the aromatic content of complex hydrocarbons becomes easy *and* fast. Easy because the keyboard has a prompt display that leads you step by step through instrument set up. Fast because instrument parameters can be indefinitely stored, recalled, and modified at the push of a button.

These features provide real multi-user versatility and virtually eliminate training. All for the price of a single research-grade spectrophotometer.




so precisely

Scattering samples can be measured accurately since the PMT is end-on and only millimeters away from your sample. A patented zoom lens maintains focus automatically throughout the entire wavelength range. High stability, high intensity UV and visible light sources, combined with minimal reflective surfaces in the optical pathway result in impressive, high sensitivity and low noise and drift.

How impressive is it? You can pluck a leaf from a tree, slip it into the sample compartment, and press a few keys. You'll have a precise scan for the chlorophyll pigments.

Give us a call. We'll send you complete information on the DW-2000.

Never before has a spectrophotometer done so much, for so many, so precisely... and so economically.



SLM AMINCO®

SLM Instruments, Inc.
810 West Anthony Drive, Urbana, IL 61801 U.S.A.
(217) 384-7730/Telex 206079



DEMAND BETTER SEPARATION FROM YOUR GEL FILTRATION.



You must try Sephacryl® High Resolution (HR) from Pharmacia—a new range of gel filtration media, and the result of continuous development of its acclaimed predecessor Sephacryl Superfine, already the first choice of thousands.

Sephacryl S-200 HR, S-300 HR and S-400 HR together cover a wide separation range. Each gel has its own steep selectivity curve providing you with excellent resolution. The balanced particle size distribution (25–75 μm) exploits the improved kinetics of a smaller bead without generating high back pressure.




The rigid matrix arrives preswollen and is ready to use. It can be packed quickly and operates well under high flow rates at any scale, from analytical to industrial. Sephacryl HR is extremely stable to denaturing agents and detergents, and can be autoclaved and cleaned in place, without harming the matrix.

Sephacryl HR is compatible with your Standard Chromatography equipment and can, of course, be used to advantage with FPLC®. You won't find a better starting point if you are scaling up to pilot- or industrial scale.

Don't delay! Order your gels now. Further information available free on request.

ORDERING INFORMATION	CODE NO.	PACK SIZE
SEPHACRYL S-200 HR	17-0584-01	750 ML
SEPHACRYL S-300 HR	17-0599-01	750 ML
SEPHACRYL S-400 HR	17-0609-01	750 ML

Sephacryl HR can separate faster—at least two times faster—than Sephacryl Superfine, and with maintained resolution! The new gels can replace Sephacryl Superfine, with improved performance, in all applications.



Laboratory Separation Division
Piscataway, New Jersey 08854
Information: (800) 526-3618
In NJ: (201) 457-8000

Circle No. 189 on Readers' Service Card

Pharmacia products will be on display at the ASCB Meeting — Booths 621, 623, 625.

Linking Value to Performance.



See spot fly!

Our new "Flying Spot" scanner combines speed of scanning, high level quantitation, high resolution and wide dynamic range for TLC, gels and autoradiograms—all in one compact system.

It's versatile. For the researcher who's not quite prepared to purchase two separate, dedicated scanners, we have combined three photometric modes in one instrument: transmission, reflectance and fluorescence. (The new CS-9000 even has a laser option when higher sensitivity is required.) Whether you're scanning electrophoretic gels, TLCs, autoradiograms, paper, cellulose or films—you get accurate, reliable and reproducible results every time; time after time.

It's fast. Optics ingenuity and advanced computer technology provide a breakthrough in densitometry: High-speed scanning using a "flying spot". This patented system increases accuracy dramatically without com-

promising high speed scanning. And to keep up with the increased amount of datapoints the CS-9000 comes with a high speed, high resolution thermal head printer/plotter and a high resolution color video display.

It's accurate. The tuneable monochromator lets you choose the wavelength of best selectivity and sensitivity. The dual wavelength capability helps eliminate artifacts. And a host of signal processing options, background correction modes and the scan store/reprocess feature gives you quantitative accuracy and lets you optimize all parameters exactly as you need them.

And it's easy to use. There are enough automation features packed into the CS-9000 to make it ideal for high sample throughput and convenient enough to be used even by a novice operator. Yet it gives you all the manual control and programming options you want for more complex analyses.

So, before you buy just another scanner, look into the universal Shimadzu CS-9000 and "see our spot fly".

Shimadzu Scientific Instruments, Inc.
7102 Riverwood Drive, Columbia, MD 21046
(301) 381-1227

In Canada: Tekscience, Inc., (416) 844-1762



For information circle reader service number 208
For a demonstration circle reader service number 263

CHPS...



BIOSCAN Imaging Scanners provide rapid, quantitative radioisotope analysis without plate scraping or scintillation counting.

Imaging Scanners are the ultimate systems for analysis of TLC plates, gels, blots and paper strips . . . A powerful quantitative tool (yet surprisingly affordable) for research in:

- drug metabolism
- lipid biochemistry
- biotechnology
- enzymology
- radiochemistry
- toxicology

Call toll-free, 800-255-7226.

BIOSCAN
4590 MacArthur Blvd.
Washington, D.C. 20007

Have you ever wondered why some ideas fly and others don't?

Could it be that the right information at the right time makes a difference?

For instance, when it comes to finding information on industry trends, intelligence on setting standards, and applications of a material process, DIALOG® helps you locate sources quickly.

You can search for information in databases such as CA SEARCH with instant access to millions of items indexed by Chemical Abstracts Service.

Or in Inspec from the Institution of Electrical Engineers, covering electronics.

Or in Paperchem, containing technical literature and patent information on raw materials.

Or in Metadex, covering international literature on metallurgy and composites.

Or in Aerospace, with extensive coverage of advanced materials composites.

And in SCISEARCH, covering the world's

scientific and technical literature.

These are just a few examples of how DIALOG provides the most valuable databases in one easy to use service.

All you need is a personal computer or terminal and a telecommunications link. And you pay only for the time you use DIALOG. There's no minimum usage requirement.

Instant delivery of essential information will make a difference and get your ideas off the ground. In fact, your company may already be using DIALOG. Check to see.


To get acquainted, DIALOG offers \$100 worth of free time to first time users. Call 800-3-DIALOG. Or write DIALOG, 3460 Hillview Avenue, Palo Alto, CA 94304.

DIALOG® INFORMATION SERVICES, INC.
A SUBSIDIARY OF LOCKHEED CORP.

The world's largest online knowledgebank.

800-3-DIALOG

Ever wonder why some things never get off the ground?



Announcing a new journal that covers research in all aspects of microbial, plant and animal fouling

a call for papers

BIOFOULING

Editor-in-Chief:

Dr Len Evans,
Department of Plant Sciences,
The University, Leeds, LS2 9JT, UK

Regional Editors

North America:

Dr Paul Boyle,
Edgerton Research Laboratory,
New England Aquarium,
Central Wharf, Boston,
Massachusetts 02110-3309, USA

Dr Madilyn Fletcher,
C.O.M.B., University of Maryland,
600 E. Lombard Street,
Baltimore, Maryland 21202, USA

Europe:

Dr Maureen Callow,
Department of Plant Biology,
The University of Birmingham,
P.O. Box 363, Birmingham,
B15 2TT, UK

Australasia:

Professor Kevin Marshall,
School of Microbiology,
The University of
New South Wales,
P.O. Box 1, Kensington,
New South Wales 2033,
Australia

Editorial Board

Robert Baier (USA) • John
Barnacle (Australia) • Jacques
Berger (Canada) • Reg Bott
(UK) • James Bryers (USA) •
H. Busscher (The Netherlands)
• Keith Cooksey (USA) •
Robert Edyvean (UK) • Jonathan
French (USA) • Gill Geesey
(USA) • Kyle Hoagland (USA) •
Philip Marsh (UK) • Tom
McMeekin (Australia) • Ralph
Mitchell (USA) • Nancy Targett
(USA) • Ruth Turner (USA) •
Graham Walker (UK)

Aims and Scope

Biofouling is an international multidisciplinary journal on fouling and fouling organisms. It provides a forum for publishing pure and applied work of all kinds on microbial, plant or animal fouling. This may be the result of attachment and growth on natural or man-made surfaces in the aquatic (freshwater and marine) or aerial environments, and medical, industrial and agricultural biofouling are also included. The journal will contain a broad range of papers, particularly in the following areas:

- * *The effects of fouling*, such as roughness (ships), corrosion (biological, chemical, engineering significance); heat transfer (heat exchanger/cooling towers); hygiene (food, biomedical)
- * *Studies relevant to fouling*, such as ecology (colonization, development of microbial plant and animal communities); attachment studies; the physics and chemistry of bioadhesion
- * *Fouling control and prevention*, in areas such as antifouling coatings, 'non-stick' non-biocidal coatings and biological control including naturally-produced, non-toxic antifoulants; metals/alloys; biocides in addition to those in antifouling paints; anti-fungal/algal paints

A Call for Papers

The editors invite contributions to the journal from professionals and researchers working in biofouling, such as biologists and microbiologists, biomedics, chemists, engineers, metallurgists, naval architects and physicists. The journal will publish original research papers, short communications, review articles, book reviews and news items. Notes for contributors may be obtained either from the editors or directly from Harwood Academic Publishers. Manuscripts may be submitted to the Editor-in-Chief or appropriate regional editor.

Subscription Information

Volume 1, number 1 due Autumn 1987

4 issues per volume

Corporate subscription price per volume: \$286.00

University/academic library subscription price per volume: \$176.00

Individual subscription price per volume: \$88.00*

*The individual subscription price is available only to individuals who subscribe directly from the publisher and who pay by personal checks or credit cards.

**First issue due Autumn 1987
Send now for your free sample copy!**



Harwood Academic Publishers

P.O. Box 786 Cooper Station New York, N.Y. 10276, USA or
P.O. Box 197, London WC2E 9PX, UK

4800 RPM **4950** xg



INTRODUCING THE 218 THREE LITER BIOCONTAINED ROTOR FOR CENTRA-8/8R

- Twice the attainable force for significantly faster separation.
- Full three liter capacity . . . up to 148 RIA tubes per run.
- New Quick-Seal transparent cup domes for simple, reliable biohazard containment.
- Sealed transparent windshield lid provides added protection and permits visual inspection of samples.
- Unique mechanical "lift-off" feature enables easy rotor removal . . . entire rotor may be autoclaved without breaking lid seal.
- Accepts new high-strength 750 ml polypropylene bottles and Centrac Adapters for most tube sizes.



CENTRA-8R[®]

REFRIGERATED CENTRIFUGE

WE'VE JUST MADE THE WORLD'S BEST BENCH-TOP CENTRIFUGE EVEN BETTER!

Since its introduction, the IEC Centra-8R has fast become one of the world's most popular high performance refrigerated centrifuges. After all, it offers the capacity, power, versatility and accuracy of a floor model centrifuge with the convenience and economy of a bench-top unit.

Now we've *doubled* Centra-8R performance with the innovative 218 rotor system. So you get even more centrifuge productivity, safety and value than ever before.

Write today for complete information on the Centra-8R and its family of advanced centrifuge accessories.

ALSO AVAILABLE IN NON-REFRIGERATED
VERSION—CENTRA-8.



INTERNATIONAL EQUIPMENT COMPANY

A DIVISION OF DAMON

300 SECOND AVE., NEEDHAM, MA 02194 U.S.A. • TEL. (617) 449-0800 • FAX (617) 444-6743

DAMON/IEC (UK) LTD. DUNSTABLE, BEDFORDSHIRE, LU6 1BD, ENGLAND • TEL. (0582) 604669



Circle No. 238 on Readers' Service Card

Welcome to GIBCO SQUARE... the place to be.



New GIBCO SQUARE packaging can make life in your lab easier. After years of research and development, our fine line of cell culture media and related products is now available in new, convenient clear plastic bottles.

This new packaging offers some very positive advantages.

- **Easier to store.** Square shape allows you to get more in less space.
- **Enhances lab safety.** Flexible plastic eliminates breakage.
- **Easy to hold.** Light weight plastic material fits the hand better than rigid packaging.
- **Easier to pipette.** Wide neck allows better access to bottle content.
- **Convenient disposal.** Plastic bottles are completely disposable.

New GIBCO SQUARE plastic bottles do not affect product quality. Our extensive testing in a wide variety of cell culture systems has conclusively proven that the new PETG plastic bottles made for us by Nalge

are equivalent to the type 1 borosilicate glass bottles we currently use. The only difference is the added convenience for you.

If you want glass containers, we'll still offer them, but we think the advantages of GIBCO SQUARE make this new packaging a good choice for any laboratory.

Want more information? Easy enough. Just use this magazine's reader service card to make your request. If this publication doesn't have a card, call your GIBCO sales representative or write us at the address shown below. We'll send you all you need to see why GIBCO SQUARE is the place to be.

GIBCO[®]
LABORATORIES

GIBCO Laboratories
Life Technologies, Inc.
3175 Staley Road
Grand Island, NY 14072
Toll-free: (800) 828-6686
Tech-Line™: (800) 828-6613
(716) 773-0700

Visit us at ASCB, Booths #812, 814, 816.

Circle No. 207 on Readers' Service Card

GIBCO SQUARE plastic bottles are available only in the United States.

Synthetic peptides: pure and simple




Look what we've done to simplify peptide synthesis. For the first time we've brought together in one instrument the long awaited features of full automation, mild chemistry and continuous flow synthesis to offer you the most advanced peptide synthesiser ever developed – the LKB Biolynx 4170.

All the drawbacks inherent in conventional t-Boc peptide synthesisers have now been overcome by the Biolynx which incorporates Fmoc polyamide active ester chemistry, already recognised for its mild reaction conditions (no HF required), versatility in providing a choice of resin linking agents, and for the high purity of the peptide synthesised. Furthermore, the use of active esters eliminates the need to change solvents during synthesis and carry out 'in situ' preparation of symmetrical anhydrides.

The unique polyamide resin has also allowed us to use a continuous flow methodology to facilitate total automation of synthesis thereby enabling the Biolynx to operate completely unattended to produce up to three peptide analogues simultaneously. Moreover, the combination of computer control and constant on-line monitoring provides reassurance that the synthesis is proceeding satisfactorily.

Peptide synthesis from LKB: it's pure and simple. Contact your local LKB representative for further details of the above system or the Biolynx 4175, our alternative, lower cost model. And don't forget that the Biolynx system is fully backed up by high quality reagents, solvents and solid phase supports.

The Biolynx 4170 peptide synthesis system.



The microcomputer monitor keeps you constantly informed of instrument and synthesis status. A typical Biolynx fluidics display is shown here.

Analytical HPLC of crude acyl carrier protein decapeptide (65-74):
H₂N-Val-Gln-Ala-Ala-Ile-Asp-Tyr-Ile-Asn-Gly-OH. LKB TSK-ODS
120T (250 × 4.6 mm) column (5μ). Elution with 20% CH₃CN in
0.1% TFA at 1 mL/min.

LKB

LKB BIOCHROM LTD, SCIENCE PARK, CAMBRIDGE CB4 4FJ, ENGLAND

Athens Middle-East +30 (1) 894 73 96 · Beijing 89 06 21 · Brussels 02 242 46 60 · Copenhagen (01) 29 50 44 · Freiburg (0761) 49030 · Hongkong (852) 5-814 8421
London (01) 657 88 22 · Lucerne (041) 57 44 57 · Madras (044) 45 28 74 · Moscow (095) 255 6984 · Paris (01) 30 43 99 19 · Rome (06) 39 90 33 · Sao Paulo (011) 578 41 38
Stockholm (08) 764 75 50 · Tokyo (03) 406 62 61 · Turku (021) 678 111 · Vienna +43 (222) 92 16 07 · Washington (301) 963 3200 · Woerden 3480 77911

Over 60 qualified representatives throughout the world.



THE POWER OF ELECTROPHORESIS AT YOUR FINGERTIPS.

High-resolution electrophoresis, a powerful technique, is rarely used to its full potential due to time-consuming and often tedious methodology.

PhastSystem™, a technological breakthrough, transforms this powerful technique into a streamlined scientific tool.

To give you fast, reproducible, high-resolution results, PhastSystem integrates precision instrumentation, PhastGel™ media (see below) and micro-processor control of all separation and development conditions.

Try PhastSystem. You'll find IEF, SDS and native PAGE equally easy and reliable to perform in just 30 to 50 minutes. And with automated development at elevated temperatures, you'll find silver staining (60 min) as convenient to perform as coomassie staining (30 min).

PhastSystem, a new class of electrophoresis system for fast, reproducible, high-resolution results.

Actual gel size.



SDS-PAGE on PhastGel Gradient 10-15 with silver staining. Samples of meat and *E. Coli* proteins. FROM SAMPLE TO RESULT: 90 MINUTES.



2-D titration curve on PhastGel IEF 3-9 with coomassie staining. Sample of beef muscle proteins. FROM SAMPLE TO RESULT: 60 MINUTES.

Circle No. 244 on Readers' Service Card



Pharmacia LKB Biotechnology Inc.
800 Centennial Avenue
Piscataway, NJ 08854
Information: (800) 526-3618
In NJ: (201) 457-8000

WHEN THESE CUSTOMERS RAN THE NUMBERS ON MINISUPERCOMPUTERS, FPS CAME OUT ON TOP.

To learn why FPS has sold more minisuper-
computers than all other companies combined,
call 1-800-635-0938.




Floating Point Systems, Inc., P.O. Box 23489, Portland, OR 97223
Telex 4742018 FLOATPOIN BEAV.

Real performance.
The facts speak for themselves.

Circle No. 209 on Readers' Service Card

PRECAST MINI GRADIENT GELS



NEW! Mini Polyacrylamide Linear Gradient Electrophoresis Gels...compatible with Bio-Rad Mini Protean II™, Model 360, Hoefer Mighty Small™ SE200/SE250 and LKB Midget™ chambers.

SAVES UP TO 6 HOURS!

For complete details
write or phone
today...

JULÉ INC.
BIOTECHNOLOGY DIVISION

25 Science Park.
New Haven, CT 06511.
Call Collect: (203) 786-5490.

Circle No. 198 on Readers' Service Card


SCIENTIFIC/ENGINEERING GRAPHIC TOOLS for the IBM PC and compatibles

FORTRAN/Pascal tools: **GRAFATIC** (screen graphics) and **PLOTMATIC** (pen plotter driver)

These packages provide 2D and 3D plotting capabilities for programmers writing in a variety of FORTRAN/Pascal environments. We support MS, R-M, LAHEY FORTRAN and more. PLOTMATIC supports HP or Houston Instrument plotters. Font module available too!

Don't want to program? Just ask for **OMNILOT!** Menu-driven, fully documented integrated scientific graphics. Write or call for complete information and ordering instructions.


GRAFATIC—PLOTMATIC—OMNILOT [S] & [P]



Microcompatibles, 301 Prelude Drive, Silver Spring, MD 20901
(301) 593-0683

Circle No. 162 on Readers' Service Card

FREE SAMPLES



TRY THE EASY WAY TO RECORD TEMPERATURES FOR RAPID COMPUTER ANALYSIS.

The new Ryan TempMentor records up to 6,361 temperatures with 0.1° C resolution, then communicates to your computer for rapid analysis and charting on your printer. Records time/temperature changes in the lab, in the field or underwater for up to 530 days, completely unattended. Choose from 16 sampling rates. Download in field or office using an IBM-compatible computer or terminal.

We'll give you a **30-day free trial** use so you can take up to 6,361 free samples of critical temperatures in your work. Offer is limited to supply on hand. Write: Ryan Instruments, P.O. Box 599, Redmond, WA 98073. Call: Bill Evans at (206) 883-7926.



Ryan Instruments, Inc.
Ryan is a registered trademark of the
Pullman-Peabody Company

ABAQUS * ADINA * ADLPIPE * AMBER * AMPAC * ANSYS * AOS/
 MAGNETIC * APTEC IOC * ARC 2D * ARGUS * ASAS * ASKA *
 ASKAMESH * ASKAVIEW * BCSLIB * BEASY * CHARMm * COMIC *
 CSPI, CSPI/G * CSSL-IV * DIS * DISCOVER * DISPLAY II * DYNA3D
 * EASY 5 * ECLIPSE * EISI-EAL * EISPACK * ENDURE * EOS-PAK * EZBEA *
 * FE2000 * FIDAP * FIMESH * FIPOST * FIPREP * FL022 * FL052 * FL057 *
 FL059 * FLODYN * FMSLIB * GAMESS * GAUSSIAN 82 * GTSTRUDL *
 HARWELL SUBROUTINE LIBRARY * HCT * HSPICE * HYDRA * IMSL *
 ISOCROSS * LINPACK * LUSAS * MARC * MATHPACK * MENTAT * MOPAC
 * MSC/NASTRAN * NAG LIBRARY * NEKTON * NISAH * NISA/FLUID *
 NISAOPT * ODEPACK * OMNILOT * PAPEC * PCGPACK64 * PHOENICS *
 PIGS * PISCES * POLYDATA * POLYFLOW * POLYMESH * POLYPLOT *
 PORES * POST NEK * PPLLIB * PREFAS * PRENEK * QSPICE * SACS *
 SALE-3D * SAM * SAVFEM * SCORPIO * SEPS * SESAM * SILOS * SIMU-
 LATE-E * SINDA * SJE * SPARSPACK * SPECTRAL64 *
 SMPACK * SOLID * STO * STREAM * SUPREM-3 *
 * TPS10 * VESPA * VIP FAMILY *
 * VSAERO * VSPICE *

WHEN IT COMES TO MINISUPERCOMPUTER SUPPORT, THE FPS LIBRARY SPEAKS VOLUMES.

To learn what the largest library in minisuper-
 computing has in store for your application, call
 1-800-635-0938.



Floating Point Systems, Inc., P.O. Box 23489, Portland, OR 97223
 Telex 4742018 FLOATPOIN BEAV.

Real performance.
 The facts speak for themselves.

Circle No. 210 on Readers' Service Card

Variability and Management of Large Marine Ecosystems

Edited by Dr. Kenneth Sherman, *Director, Narragansett Laboratory, National Oceanic and Atmospheric Administration*, and Dr. Lewis M. Alexander, *Director, Center for Ocean Management Studies, University of Rhode Island*

Large marine ecosystems (LMEs) are being subjected to increasing stress from industrial and urban wastes, aerosol contaminants, and heavy exploitation of renewable resources. This book is a state-of-the-art review of effective means for measuring changes in populations and productivity, physical-chemical environments, and management options for LMEs. For the first time, this volume treats LMEs holistically as regional management units by bringing together the all too often fragmented efforts to optimize ocean resources. 319 pp., 1986.

\$31.85; AAAS members \$25.50 (include membership number from *Science*).

Order from: Westview Press, Dept. AAAS, 5500 Central Avenue, Boulder, CO 80301
Add \$2.50 postage and handling for the first book ordered; 75¢ for each additional book.

Published by Westview Press for AAAS

Tropical Rain Forests and the World Atmosphere

Edited by Ghilleen T. Prance, *Director, Institute of Economic Botany, New York Botanical Garden*


The ongoing destruction of tropical rain forests may have profound consequences for the global atmosphere. Increasing scientific knowledge finds that the role of the tropical rain forests in maintaining the equilibrium of the atmosphere may be far greater than previously believed. Based on a AAAS symposium, this volume reports on the urgent need for preserving the tropical rain forests and includes policy recommendations. 106 pp., 1986.

\$18.50; AAAS members \$14.80 (include membership number from *Science*).

Order from: Westview Press, Dept. AAAS, 5500 Central Avenue, Boulder, CO 80301
Add \$2.50 postage and handling for the first book ordered; 75 cents for each additional book.

Published by Westview Press for AAAS

6 NEW WAYS TO EXPRESS YOURSELF



Genes also available:

Interleukin 1 α
Interleukin 1 β
Interleukin 2
TNF α
TNF β
Interferon γ
Epidermal Growth Factor
Insulin-like Growth Factor 2

Genes available soon:

Interleukin 4
Interleukin 5
Interferon β ₂
M-CSF
TNF α (nat)
TNF β (nat)
Interleukin 1 β (nat)

BBL's range of genes is constantly expanding. Please contact us for further details of our new products and custom synthesis service.


British Bio-technology Limited

British Bio-technology Ltd, Brook House, Watlington Road, Cowley, Oxford OX4 5LY. Telephone (0865) 718817. Fax (0865) 717598. Telex 838083 BIOTEC G.

Don't buy any microplate reader until you've seen its software.


Watch all 96 kinetic plots appear as the reactions proceed.

Adjust your OD plot limit (y-axis): shorten it for fast reactions, extend it for slower ones. SOFTmax™ will replot the data and recalculate rates using only the initial linear portion of the reaction.




TEMPLATE SET-UP works like a chemical spreadsheet.

Assign wells to group types (blank, standard, control, unknown, empty). Run your assay or overlay the template on existing data and recalculate. Use it to mask outliers and set up multiple standard curves on the same microplate.



For quantitation, try any of 5 standard curve-fitting algorithms.

Linear. Quadratic. Semi-log. Log/logit. 4-parameter logistic. SOFTmax can also format your data in ASCII for use with Lotus 1-2-3® or your own compatible data processing program.



Flexible data displays.

Use Grey Scale for cytoproliferation studies, Ranged data for hybridoma screening, Threshold display for serology. All display options may be used with kinetic and endpoint data.



Moving from data to results requires more than instrumentation. The new SOFTmax™ controls the Vmax™ Kinetic Microplate Reader, collects and analyzes endpoint and kinetic data, and delivers many useful options for computing, plotting, displaying, printing and managing microplate assay data. All in one simple, menu-driven package.

Evaluate SOFTmax on your IBM® or IBM® compatible microcomputer. And find out how our powerful microplate assay software can make a difference in managing data *and* results. Contact Molecular Devices today for a free SOFTmax demonstration diskette.

Vmax and SOFTmax are trademarks of Molecular Devices Corporation. IBM is a registered trademark of International Business Machines Corporation. Lotus 1-2-3® is a registered trademark of Lotus Development Corporation.


Molecular Devices

Molecular Devices Corporation
3180 Porter Drive
Palo Alto, CA 94304
(415) 493-0166
(800) 635-5577
(800) 635-5588 (within California)

For \$19⁹⁵, RS/1[®] will give you a world of understanding.

Thousands of engineers and scientists are using RS/1[®] software to get a better understanding of their technical data. And now you can find out why.

For the first time, you can see why the RS/1 package has become an industry standard in data analysis software. All on your own PC.



RS/1 software displays results in tabular and graphical forms.

Just send for our new introductory trial version of RS/1 software. For \$19.95, you get a special demonstration version of RS/1 software with data storage limitations. With this trial package,

you can perform all sorts of data analysis in a fraction of the time it would take with conventional methods.

RS/1 Data Analysis System

- Data management: easy entry, import, storage, and retrieval of data
- A full range of data analysis tools
- Presentation quality graphics: graphs, histograms, 3Ds and more
- Easy-to-use statistics: t-tests, 1- and 2-way Anova, non-parametric, and other standard tests
- Spreadsheets and models
- Curve-fitting: linear, polynomial, multiple stepwise, and nonlinear - with goodness-of-fit statistics
- A powerful programming language for customizing applications

Experience easy-to-use statistical and curve-fitting capabilities. And presentation quality graphics. Plots, graphs, histograms. This trial package does it all.

Before long, you'll discover how RS/1 software can help you improve laboratory research, cut product development costs, and increase manufacturing yields. And that's just the beginning. There's no telling what you can do with such a fast, powerful data analysis tool.

And the \$19.95 you pay for the trial version can be applied to the purchase of the complete RS/1 software package for the IBM PC AT and XT.*


Call 1-617-499-8200 for a credit card purchase, or fill out the coupon and send it in.

And find out how much of your world you can understand.

Circle No. 260 on Readers' Service Card

BBN Software Products

10 Fawcett Street Cambridge, MA 02238

A Subsidiary of Bolt Beranek and Newman Inc. 

*Or full equivalent with 10MB hard disk and 512KB memory.

RS/1 is a registered trademark of BBN Software Products Corporation.

Enclosed find a check for \$19.95 for the introductory version of RS/1 for the IBM PC. SCI 1087

AMEX # _____ Exp. Date _____

Signature _____

Name _____


Title _____

Company _____

Address _____

City _____ State _____ Zip _____

BBN Software Products, 10 Fawcett Street, Cambridge, MA 02238



Membrane vs. polystyrene based assays?

S&S asked proteins to decide.

We at S&S had a head start in developing an efficient and economic alternative to standard ELISA testing.

Because we know a lot about the behavior of proteins.

We know their properties. What they bind to. And what they don't bind to.

It comes from the experience of thousands of dot-blots and "westerns" done with S&S products.

So when we modified the Minifold I system for protein assays, we already had the "opinions" of trillions of proteins.


And this is what they told us.

HIGHER SIGNAL - LESS NOISE.

Polystyrene 96-well plates are excellent binders of proteins. Unfortunately, they are also good at binding just about any other organic particle with which they happen to come in contact. The result - higher background noise.

With the S&S Minifold I system, modified for incubation, you can choose the substrate that has the highest binding capacity for the protein you are studying. Without inhibiting biological activity.

And with the Minifold system, noise and false positives are greatly reduced because of the low non-specific binding characteristics of S&S papers and membranes.



Plus an assay run with the Minifold I system requires less processing and incubation time than standard ELISA systems.

Best of all, the Minifold system is priced thousands less than standard polystyrene-based ELISA systems, which means that the system is well within the reach of anyone who wishes to add protein assay capabilities to their lab.

FIND OUT ABOUT MINIFOLD I.

Call for a paper comparing the features of both polystyrene-based and filter-based assays. Or just to talk about protein binding in general. We always welcome your questions and comments.

That's because in addition to finding out what a protein might tell us, we're interested in your opinion, too.

MORE ACCURATE AND REPEATABLE RESULTS.


You can't count on the binding characteristics of plastic supports. Binding varies from maker to maker. Plate to plate. Protein to protein. Ergo, results vary just as widely.

But with Minifold filter-based assays, known quantities can be bound. And binding characteristics vary little from batch to batch.

And results can be quantified by colorimetric and isotopic techniques, fluorography, scintillation counting or densitometry.

BETTER ASSAYS, FOR THOUSANDS LESS.

Unlike other systems, the Minifold I features both filtration and incubation in the same unit.



The Minifold I consists of (1) the 96-well sample plate, (2) the filter support plate and (3) the vacuum plenum. The accessory incubation plate (4) replaces the filter support during incubation.

An EC910 Transmission Densitometer from the EC Apparatus Corporation was used for the dot-blot reading.

Schleicher & Schuell

Keene, NH 03431 800-245-4024 • 603-352-3810

See us at ASCB, Booth #638, 640.

Circle No. 205 on Readers' Service Card

Introducing **The Fluoro IV™** Spectrofluorometer

Performs to your high standards
Circle No. 127 on Readers' Service Card

Reinforcing our commitment to listening to our customers, we have combined measuring power, flexibility and simplified operation into a single instrument. The Fluoro IV...developed to meet the demands of researchers, day in and day out.



The Fluoro IV features include:

- 4-Position Thermostatted Cell Holder
- Printer/Plotter for Both Tabular and Graphical Data Presentation
- Excitation, Emission, and Synchroscans
- Time-Based Studies
- HPLC Adapter
- RS232...and Much More.

For more information, call us toll-free at **800-GILFORD**.


CIBA-CORNING

Ciba Corning Diagnostics Corp.
Gilford Systems



Environmental Control From Percival

Have control of the environment with a biological, insect rearing, dew, or plant growth chamber from Percival. Whether you need a table model or a giant walk-in, Percival has the knowledge, experience, and quality product line to meet your specific needs. Or, Percival will prepare a recommendation to fit your requirements. Write today for more information and the complete Percival catalog or call collect 515/432-6501.



Percival

Percival Manufacturing Company
P.O. Box 249 Dept. S-10 Boone, IA 50036

...The name to remember for versatile Biological Incubators, Insect Rearing Chambers, Dew Chambers, and Plant Growth Chambers.


Circle No. 3 on Readers' Service Card

A New Slant on Statistics

\$390. That's the cost of joining the many others who have already discovered STATA. Of course, quantity discounts, leases, and discounts for qualified educational institutions are available. And there's a 30-day money-back guarantee.

\$20. That's the cost of our Demonstration System. Call 1-800-STATAPC, or, from California, Alaska, and Hawaii, call 1-213-470-4341. Place an order or just ask questions.

- ☒ ANOVA and ANCOVA
- ☒ Chi-square tests
- ☒ Cumulative functions
- ☒ Cross-tabulation
- ☒ Covariance & Correlation matrices
- ☒ F-tests & t-tests
- ☒ Influence statistics
- ☒ Logit & probit
- ☒ Multiple-comparison tests
- ☒ Non-parametric tests
- ☒ Predicted and residual values
- ☒ Random numbers
- ☒ Regression
- ☒ Univariate statistics
- ☒ Weighted & unweighted data
- And, more.



STATA

Computing Resource Center
10801 National Boulevard
Los Angeles, California 90064
1-800-STATAPC
(In California, Alaska, and Hawaii:
1-213-470-4341)
Telex: 263-199-STATAPC

1-800-STATAPC

Circle No. 216 on Readers' Service Card

Micromanipulators, Microscope Stages and Precision Laboratory Tools


C A T A L O G

Manual, DC or stepper motor micromanipulators and microscope stages

Custom-made precision slides and positioners

Micromanipulators and bases in-stock for immediate delivery

Send today for this catalog featuring hundreds of general and special purpose micromanipulators.



The single source for all your micromanipulator needs.

STOELTING


1350 S. Kostner Ave. Chicago, Illinois 60623
(312)522-4500

Circle No. 38 on Readers' Service Card

Announcing a Conference on

DIRECTIONS IN

BIOSCIENCE



APRIL 11-15, 1988

Newark, Delaware, and Philadelphia, Pennsylvania

Speakers

Kari Alitalo	Klaus Mosbach
Margareta Baltscheffsky	Peter Nowell
Baruch Blumberg	Jorma Palo
Kari Cantell	Karl Paul
John Eveleigh	Per Petterson
Kjell Fuxe	Jan Ponten
Paul Greengard	Darwin Prockop
Goran Holm	Garry Rechnitz
Leroy Hood	Eero Saksela
Leevi Kaariainen	Jouni Uitto
Kari Kivirikko	Torsten Wiesel
George Klein	Hans Wigzell
J.M. Lindstrom	Morton Wikstrom


This international event will be convened in commemoration of the 350th anniversary of Swedish and Finnish settlement in America and in recognition of the significant developments in human bioscience that have been made by researchers from the three countries. The speakers will focus, however, not on past achievements but on the directions bioscience will take in the next decade. Following are the session topics:

- Human Molecular Genetics
- Regulation of Immune Response
- Biotechnology: Tools and Applications
- Regulation of Normal and Abnormal Cell Growth
- Frontiers in Neuroscience
- Molecular Bioenergetics

The conference is sponsored by the Academy of Finland, the Fox Chase Cancer Center, the Swedish National Committee for New Sweden 88, the Thomas Jefferson University, the University of Delaware, and the University of Lund and the Bioscience 88 Conference.

For registration materials, send name, address, telephone number, and organization to the following address:

DIRECTIONS IN BIOSCIENCE
Office of Research & Patents
210 Hullihen Hall
University of Delaware
Newark, DE 19716 U.S.A.
(302) 451-2136



Microelectrodes For Measuring
 pH, O₂, CO₂, Na⁺, K⁺, Ca⁺⁺,
 NH₃, ORP

Send for *free* information!

MICROELECTRODES, INC.
 OAK HILL PARK, LONDONDERRY, NEW HAMPSHIRE 03053
 (603) 668-0692

Circle No. 230 on Readers' Service Card

PICO-INJECTOR

PILI-100

**for all your injection needs...
 extracellular-intracellular**

- Reproducible delivery • Digital displays • Balance/clear pressures
- Vacuum fill/hold • Optional foot switches • Injection count

The Medical Systems Corporation PICO-INJECTOR is a highly versatile single channel pressure unit that delivers volume from microliters to femtoliters using micropipettes. Its reproducible delivery allows a beginning or experienced user to master all aspects of microinjection quickly.




THE RIGHT TECHNOLOGY
Medical Systems Corp.

Call One Plaza Road, Greenvale, NY 11548 • (516) 621-9190
 800-654-5406 TELEX: 645031 MED SYSTEMS • FAX: 516-621-8503

Circle No. 104 on Readers' Service Card

SECOND ANNOUNCEMENT



ENZON'S 2ND CONFERENCE IN A SERIES "Innovations in Protein Therapeutics: From Research to the Clinic"

December 8-11, 1987
Walt Disney World Resort
Lake Buena Vista, Florida USA

The field of peptide and protein therapy has expanded tremendously over the past few years. Both heterologous and homologous proteins are presently being considered for treatment of a variety of diseases. This conference is intended to address improvements in the pharmacodynamics, pharmacokinetics and immunology of proteins through the following methodologies:

Stabilization
 Modification: Polyethylene glycol (PEG), albumin
 Encapsulation
 Controlled release

Presentations will either review the current state of each technology, or discuss the specific applications from the viewpoint of basic research, preclinical or clinical studies.

For further information regarding this conference, please call or write: M. L. Nucci, ENZON, Inc., 300 C Corporate Court, South Plainfield, NJ 07080 USA. (201) 668-1800. Early registration ends November 13, 1987.

Circle No. 251 on Readers' Service Card

BURN YOUR REFERENCE CARDS!

REF-11TM

Computerizes your REFERENCES
 and prepares your BIBLIOGRAPHIES

- ☐ Maintains a data base of references
- ☐ Searches for any combination of authors, years of publication, reference title, publication title, keywords or abstract
- ☐ Formats bibliographies exactly as you want them
- ☐ Reads your paper, inserts citations into the paper, and prepares a bibliography of the references cited (optional)
- ☐ Downloads references from MedLine data bases such as NLM, BRS and DIALOG (optional)

IBM PC/XT/AT, MS-DOS, CP/M 80 ... **\$195⁰⁰**

RT-11, TSX-Plus, RSX-11, P/OS **\$250⁰⁰**

VAX/VMS (native mode) **\$350⁰⁰**



ANY **\$15⁰⁰**
 MANUAL

ANY **\$20⁰⁰**
 MANUAL
 & DEMO

322 Prospect Ave., Hartford, CT 06106
 (203) 247-8500

Connecticut residents add 7 1/2 % sales tax.

Circle No. 16 on Readers' Service Card

 Polaroid

It took us four years
to develop this slide.

Now you can do it
in 4 minutes.

PolaBlue Instant Slide

**Introducing PolaBlue.
The first and only 35mm
instant white-on-blue
slide film.**

It took Polaroid to develop this tremendous breakthrough. White-on-blue 35mm slides that you can now make instantly. Unique new PolaBlue 35mm instant film from Polaroid provides sharp, bright, high contrast slides of the highest quality for presentation or documentation in just minutes.

You'll also find new PolaBlue is far more cost-effective than traditional methods of making blue slides—less than \$1.25 per slide.* So with PolaBlue and the Polaroid 35mm Instant Slide System it's easy to make white-on-blue slides in-house, in minutes, within budget.

Pick up the phone now for technical assistance, more information or the Polaroid Professional Film Center nearest you. Call toll-free 800-343-5000.

For more information, mail this coupon to Polaroid Corporation, Dept. 676, P.O. Box 5011, Clifton, NJ 07015.

Please send me information about:

- ☐ New PolaBlue 35mm Instant Slides
☐ The complete Polaroid 35mm Instant Slide System

SC 10/16

Name _____

Company _____

Address _____

City _____

State _____ Zip _____

Telephone _____

 **Polaroid**




When it comes to high-quality peptides, Bachem wrote the book.

And the book, we're pleased to report, contains more information and is easier to use than ever.

As in previous years, the Bachem Catalog has the most up-to-date line of Bioactive Peptides, Substrates and Synthetic Intermediates available.

But the 1987-88 catalog has even more. We've added molecular weights for every product. We've added an easy-to-use alphabetical index. And to assist you in researching new products, we've added literature references on new peptides.

In addition, the new Bachem Catalog features hundreds of new products, including: *Magainins*, *Hypercalcemia of Malignancy Factor (1-40)*, *ω -Conotoxin GVIA*, and specific antagonists for *Vasopressin*, *Bradykinin* and *Bombesin*.



For a free copy of the new Bachem 1987-88 Catalog, call us today. You'll see that finding high-quality peptides is now easier than ever.

b BACHEM INC.
FINE CHEMICALS

BACHEM INC.
3132 Kashiwa Street, Torrance, CA 90505
Telephone (213) 539-4171
TWX 910-347-7320 BACHEM TRNC
Telefax (213) 530-1571

BACHEM (UK) LTD.
69 High Street, Saffron Walden
Essex CB10 1AA, England
Telephone (0799) 26465
TLX 817205 BACHEM G Telefax (0799) 26351

PHIDEA S.R.L., in Italy
Via Ariosto, 34, 20145 Milano, Italia
Tel. (02) 4988405-434061
Telex 321315 PHIDEA I
Telefax (02) 4697387

Circle No. 194 on Readers' Service Card

YOU GET A LOT OF KODAK WITH EVERY CATALOG.

NEW PACKETS, NEW KITS FOR ELECTROPHORESIS

Now you can order new convenience from the fullest electrophoresis line on the market. Full enough to warrant a separate section in the newest Kodak catalog that includes high quality reagents, photographic products and analytical instrumentation.



New BufferEZE™ Packets For Electrophoresis

Packet convenience simplifies handling, because you've got better things to do than mix. BufferEZE Formula 1 TRIS/Glycine is the first in a full line of pre-mixed buffers. A single packet makes 4 litres. **Call or circle No. 143.**


New Gel Duplicating Kits For Film Or Paper

Quality reproductions. No camera needed. Contact photographic process takes just 15 minutes. You can even enhance faint protein bands. EDF film for silver and nickel stained gels. EDP paper for gels stained with blue dye. Distinctly Kodak.

Pre-packaged starter kits contain developing trays, box of electrophoresis film or paper, photochemical kit, electronic timer, clips and stirrers. We also offer the film, paper and photochemicals separately. **Call or circle No. 144.**




NEW AND FREE



You'll find the usual *and* the unusual in the new KODAK Laboratory and Research Products Catalog No. 53. All the specialties you expect from Kodak in six expanded sections—General Organics, Electrophoresis, Chromatography, Dyes and Stains. LabLEADER™ Analytical Standards and Biological Probes. And we're one supplier that's ready to produce and deliver in about any quantity you want—gram, kilo or more. Lower pricing for volume. Call for a quote.

EXPANDED LINE OF MERRIFIELD PEPTIDE RESINS

Chloromethylated Polystyrene 1% and 2% Divinylbenzene Copolymer Beads are the new high quality resins for peptide synthesis. *Manufactured and characterized by Kodak.* Available in gram to kilo quantities. **Call or circle No. 145.**



Call 1-800-225-5352

Toll free from anywhere in the U.S. (New York included!) to order or request catalog.

Call 1-800-225-5572

Toll free from anywhere in the U.S. (New York included!) to talk with a technical specialist.

Or circle No. 146.

to get your copy of the catalog anywhere in the world. To order products outside the U.S., check our dealer list in the catalog or call 1-716-458-4014.



LABORATORY AND
RESEARCH PRODUCTS DIVISION
EASTMAN KODAK COMPANY
ROCHESTER, NY 14650

OLYMPUS®

The Image of Quality

Introducing the IMT-2.

The new inverted tissue culture microscopy system with extended routine and research capabilities.

The quality of any microscope system relies on the interrelationship of optical performance, mechanical design and the component support system that increases its versatility as your needs grow. That is the Olympus IMT-2 advantage!


Whether used for laboratory routine or research study, the IMT-2 offers greater capability, its low-profile fixed stage and focusing

optics permitting use with a large variety of vessels from ordinary slide to Petri dish and flask. A wide choice of components offers brightfield, phase contrast, differential interference contrast, Nomarski, epi-fluorescence and polarizing microscopy, as well as parfocal multiple-format photomicrography and multiple unit viewing. Combined with built-in magnification changer, 3-way light path selector and superior long barrel LB optics, IMT-2 represents a significantly advanced new generation in inverted microscope design ideally suited to cell manipulation, monoclonal antibody and tissue culture research.


See for yourself by contacting Olympus Corporation, Precision Instrument Division, 4 Nevada Drive, Lake Success, New York 11042-1179. Phone toll-free 1 (800) 446-5967 for demonstration.

In Canada: W. Carsen Co., Ltd., Ontario.

OLYMPUS®



For information circle reader service number 29
For demonstration circle reader service number 30



Cut the sequencing story short while relying on your current know-how with our automatic sequencer

Genoscope™ helps you gain time in sequencing thanks to:

automatic reading of the band pattern represented on the screen in real time. Your intervention is only required to remove ambiguities. The computer will assist you in data analysis.

well mastered techniques used: Maxam and Gilbert as well as Sanger.

Genoscope's built-in detectors. You no longer need long gels and gradients.

shortened run time due to:

- 16 cm long gels,
- continuously optimized migration speed of the bands.

increased reliability due to:

- automatic control by computer,
- elimination of gel handling.

Moreover, GENOSCOPE™ offers you **increased safety** by reducing the handling and the quantity of labelled products used.



For further information please contact,

GENOSFIT

CP 239 -CH-1212 Grand-Lancy 1/Genève, Switzerland – Tel. (022) 71 44 44 – Telex 423 082 gfit ch.

Circle No. 43 on Readers' Service Card

What the 210,625 best supplied scientists are reading.

The brand new VWR Scientific Catalog, complete with 50,000 items to meet your every laboratory need. From the basics such as chemicals, pH meters and laboratory furniture to sophisticated instrumentation and supplies for cell culture and hybridoma work, electrophoresis, chromatography and ultrafiltration. Specialty products for molecular biology, immunology and protein chemistry...

To make sure these products reach your laboratory when you need them, VWR maintains 27 offices across the country, and large inventories at local warehouses. This means that most of your orders are filled immediately from stock.

We also want our customers to be the first to know about new products, special promotions, seminars and important new reference books. Six times a year, our *VWR Laboratory Showcase* arrives at laboratories,


packed with information about the latest developments—in laboratory safety, computerized procurement and inventory management... A recent issue, for example, features an expanded section on our growing line of products used in biotechnology.

Become the 210,626th best supplied scientist. Send in your reader service card and ask for a FREE subscription to *The VWR Laboratory Showcase*.



VWR.
Today's
Biotechnology
Supplier.

P.O. Box 7900
San Francisco, CA 94120
800-634-5423 619-259-1015 (In California)



VWR Scientific
A VWR COMPANY

Circle No. 159 on Readers' Service Card

The care to meet to know, to understand.

Ares-Serono Symposia is an independent foundation, created in 1971, to promote scientific research in all disciplines which contribute towards improving human health. This aim is pursued by means of congresses, courses, seminars and specialized studies.

Some of the international meetings scheduled for 1988 are:

5th European Workshop on the Molecular and Cellular Endocrinology of the Testis
Brighton, UK / April 13-16
Scientific Organization: B.A. Cooke (UK)

Endocrinology Under 35
Florence, Italy / May 23-25
Scientific Organization: M. Maggi (I) and C.J. Johnston (USA)

Platelets and Vascular Occlusion
Rome, Italy / June 1-3
Scientific Organization: G.A. FitzGerald (USA) and C. Patrono (I)

Advances in Biotechnology of Membrane Ion Transport
L'Aquila, Italy / September 19-20
Scientific Organization: R. Verria (I), P.L. Jørgensen (DK) and R.P. Garay (F)

Third Conference on Differentiation Therapy
Villasimius (CA), Sardinia/September 6-9
Scientific Organization: G.B. Rossi (I), F. Takaku (J) and S. Waxman (USA)

The Adrenal and Hypertension: From Cloning to Clinic
Tokyo, Japan / July 25-26
Scientific Organization: E.G. Biglieri (USA), J. Funder (AUS), F. Mantero (I) and R. Takeda (J)

No form of advertising or promotional activity by any pharmaceutical company is allowed at Ares-Serono Symposia Congresses.



ARES
SERONO
SYMPOSIA

Do you know all the advantages that the Ares-Serono Symposia Congress Card offers?
Please fill in this form and mail to Ares-Serono Symposia for information.

CONGRESS
CARD
1988



ARES
SERONO
SYMPOSIA

swissair  Official Carrier

I would like to receive information about:
☐ Ares-Serono Symposia Congress Card
☐ All the above Congresses
☐ In particular the Congress on

☐ Ares-Serono Symposia Publications

Name

Address

Institution

Please send to Ares-Serono Symposia,
Via Ravenna 8 - 00161 Rome - Italy

sintesi informazione - roma

GET SERIOUS

NEW MENUING
INTERFACE AND
RELIABILITY PROCEDURE



...ABOUT ANALYZING YOUR DATA.

You might be spreading your spreadsheet a little too thin. Or maybe you're starting from scratch. But if you're serious about data analysis, you're ready for SPSS/PC+™—a full software family that brings you eight high-powered ways to complete any data analysis task.

Enter it. SPSS Data Entry II™ is a fully integrated data entry, cleaning and editing tool.

Analyze it. The SPSS/PC+ base package provides a powerful array of statistical and reporting procedures.

Picture it. SPSS/PC+ Graph-in-the-Box™ featuring New England Software's Graph-in-the-Box™ offers full color "snapshot" graphics.

Examine it. SPSS/PC+ Advanced Statistics™ lets you get more serious with your data.

Predict it. SPSS/PC+ Trends™—our latest option—is the complete time series analysis/forecasting tool.

Table it. SPSS/PC+ Tables™ produces presentation-ready tables instantly.

Chart it. SPSS/PC+ Graphics™ featuring Microsoft® Chart creates show-stopping graphs and charts.

Map it. SPSS/PC+ Mapping™ featuring MAP-MASTER™ creates maps where vast amounts of data can be summarized and presented in one, simple picture.

SPSS/PC+ products are being put to productive use by serious fact finders in business, government and education. For countless purposes such as market research. Wage and salary studies. Survey analysis. And quality control. Plus each product is superbly documented and supported by SPSS Inc., a leader in statistical software for nearly 20 years. While specially tailored customer support is available through the VALUE PLUS™ plan. And SPSS now offers a SPSS/PC+ version for Novell LANs.

So if you're serious about data analysis, step up to SPSS/PC+. For details, contact our Marketing Department.

CALL 1/312/329-3319

SPSS inc.

SPSS Inc. • 444 North Michigan Avenue, Suite 3000 • Chicago, Illinois 60611


In Europe: SPSS Europe B.V. • P.O. Box 115 • 4200 AC Gorinchem, The Netherlands • Telephone: +31183036711 • TWX: 21019

SPSS/PC+ runs on IBM PC/XT/AT's with hard disk. Contact SPSS Inc. for compatible microcomputers. SPSS/PC+, SPSS Data Entry II, SPSS/PC+ Graph-in-the-Box, SPSS/PC+ Advanced Statistics, SPSS/PC+ Trends, SPSS/PC+ Tables, SPSS/PC+ Graphics and SPSS/PC+ Mapping are trademarks of SPSS Inc. VALUE PLUS is a trademark of SPSS Inc. Chart and Microsoft are trademarks of Microsoft Corporation. MAP-MASTER is a trademark of Ashton-Tate. Graph-in-the-Box is a trademark of New England Software, Inc.


© 1987, SPSS Inc.

Circle No. 249 on Readers' Service Card


BUY THREE




6, 12, 24 wells. 96 well TC & ELISA



15, 50, 250 mL



35, 60, 60 (w/grid), 100 mm



25, 75, 150 cm²


GET ONE FREE

Buy three cases of two or more of these quality Corning Orange plastic products.

Get a Corning Orange T shirt free.

Take your pick from the full array of Corning multiple well plates, centrifuge tubes, and tissue culture dishes or flasks.

Mix and match two or more of



these products for a case total of three. We'll send you a Corning Orange T shirt absolutely free!

See your Corning or distributor rep for full product details. Ask for our latest biotechnology catalog.

Send coupon with proof of purchase for your free T shirt from Corning. Offer good through 30 November, 1987.



SEND IT TODAY

Here's my proof of purchase for a minimum of three cases of these Corning Orange plastic products. Send me one T shirt free.

Size: ☐ Small ☐ Medium ☐ Large ☐ X Large
Offer good through 30 November, 1987. USA, Canada only.

Mail coupon and proof of purchase to Corning Orange T Shirt Offer, P.O. Box 10484, Rochester, NY 14610.

Name _____
Title _____
Company _____
Address _____
City _____ State _____ Zip _____
Phone (____) _____

SE-10-87-CB

Circle No. 183 on Readers' Service Card

Science Products
a **CORNING** Laboratory Sciences Company

The Most Trusted Tools Of Science

CORNING



Your NMR work deserves our RF amplifiers.

Conservatively-rated power amplifiers, with the noise-blanking capability that pulsed NMR demands, have been a specialty of ours for well over a decade. Whether your needs for clean rf power are at the 200- to 500-watt level (as supplied by our Model 200L shown here) or up in the kilowatt range, we have the pulse power systems to ensure your peace of mind.

During pulse operation (at duty cycles up to 25%), the 200L can deliver up to 500 watts over a bandwidth of 1-200 MHz; yet when blanked with a +5V signal it reduces noise 30 dB in less than 5 microseconds. We know how important that noise-free environment is to the integrity of your results.

If you're upgrading an existing system or moving into high-power spectrometry for solid-material experiments, we suggest you work for a few moments with an AR amplifier. Enjoy the instant frequency response without need for tuning or bandswitching; the total immunity to any degree or phase of load mismatch; the assurance that nowhere within the bandwidth will the output power be less than the rated minimum. (When we say minimum, we mean *minimum*.)

Call us to discuss your present and expected applications. Or write for our NMR Application Note 0013 and the informative booklet "Your guide to broadband power amplifiers."



160 School House Road, Souderton, PA 18964-9990 USA
Phone 215-723-8181 • TWX 510-661-6094

7816

Course Announcement

FUNDAMENTAL NEUROBIOLOGICAL APPROACHES TO DISEASE: Neurodegenerative Diseases

January 17-22, 1988

The Neurosciences Institute will sponsor a course emphasizing the disease process as an object for biological research. The purpose is to introduce scientists to the basic physiology and pathology of a selected group of human neurodegenerative diseases and to the opportunities they provide for solving fundamental biological problems.

This intensive course will focus on Alzheimer's, Parkinson's, and motor neuron diseases, reviewing both clinical phenomena and pertinent areas of basic neuroscience. The course format includes lectures, clinical case presentations, and small group discussions. Instructors include David Anderson, Anders Björklund, Benjamin Bunney, John Growden, Mark Gurney, James Gusella, Zach Hall, Thomas Jessell, Joseph Martin, Donald Price, Peter Spencer, Stephen Waxman, Mark Willard, and Michael Zigmond.

Applications are invited from scientists with advanced training and experience in neuroscience or related biological sciences. Participants' expenses will be reimbursed by the Institute.

For application forms or further information, please write to Dr. W. Einar Gall, Research Director, The Neurosciences Institute, 1230 York Avenue, New York, New York 10021, or telephone (212) 570-8975.



The Neurosciences Institute
of the Neurosciences Research Program

PRECISION® SHIPPING IN 72 HOURS OR YOU SAVE \$100

New Precision® Incubators. Stocked To Deliver.

We guarantee your new Precision® incubator will be shipped within 72 hours of the time your distributor gets the order. If your distributor is temporarily out of stock, call Precision and we'll ship for your distributor within 72 hours. If we fail, you'll save \$100 on your next purchase of Precision products.

We can make this most timely guarantee because we've just built a huge volume of newly updated Precision incubators, and our distributors are stocking up so they can serve you quickly.

You won't have to wait long to enjoy the quality performance new Precision incubators deliver.

It's Right From Precision Scientific

We've built more quality constant temperature equipment for the lab than anybody else in the business. Made in America, for nearly 70 years. Fast, responsive service is just a phone call away.

Talk to your distributor rep. Call us at 800-621-8820. Or send coupon for more.

The Precision Scientific 72 hour shipping guarantee is good throughout all of 1987.


New digital model 4L

Get $\pm 0.5^{\circ}\text{C}$ uniformity at 37°C , measured at 9 places throughout the chamber.

Platinum sensor rests high in the chamber for a most accurate reading.

Tempered inner glass doors let you check samples without disturbing uniformity.

Low watt density heater brings you a broad surface area for steadier heat that's easier to control.




New microprocessor control carries a full five year warranty. It continuously monitors setpoint against actual $^{\circ}\text{C}$ to minimize overshoot. Controller cycles heat in split second bursts near setpoint to minimize fluctuation.

Just touch the pad to enter your temperature to 0.1°C . No knobs, no guesswork. Bright LED constantly displays $^{\circ}\text{C}$ in real time, never an average.

Circle No. 215 on Readers' Service Card

Your Performance Demands
PRECISION®

Visit us at ASCB Nov. 16-20, Booth #'s 629-631.



Make It Quick


Get our new Precision incubator literature and details on our new 72 hour shipping guarantee.

Mail coupon on your letterhead with name, title, and address to Precision Scientific, 3737 West Cortland Street, Chicago, Illinois 60647.

S101687

MEMBER
SAMA
SCIENTIFIC APPARATUS
MAKERS ASSOCIATION

Precision Scientific



```

DNA Homology Search
FILE: HURH2BPS.DNA 3121 214 bp 500 NOR 1 - 214 PRINT ON
NO TARGET FILE DEFINITION INIT OPT
1 XENHIS2A.DNA 298 312
66.3% identity in 178 bp overlap

      10      20      30      40      50      60
HURH2BPS.DNA  GGTAGCGACGGCGGTTTTCGGCGCCTTTCGGATTGCCAGCGAGGAGTTTCTCTCGTGAC
                : : : : :
XENHIS2A.DNA  AATATCAGTAAGTAGGGGAGTGCAGCTTCAGTCTACACATCTTCTGATTGGTTGAT
      80      90      100      110      120      130

      70      80      90      100
HURH2BPS.DNA  TACTATCGCTGT---CATGTCTGGTGGTGGCAGCGAGGAGGCA-----GGC
: : : : : : : : : : : ^ : : : : : v : : : : :
XENHIS2A.DNA  TTGTAGCAGCTAATCATGTCTGGAGAGGCGAAGCAGCGCGCAGACTCCGCTAGGCG
      140      150      160      170      180      190

      110      120      130      140      150      160
HURH2BPS.DNA  CAGATCGCGCTCGTCCGCGCTGGCTTCAGTTCGGGTAGGAGGAGTGCATCGCTTGT
: : : : : : : : : : : : : : : : : : : : : :
XENHIS2A.DNA  AAGAGCTCGCTCATCTCGAGCGCGGCTGCAGTTCCGAGTCCGCGGTGTTCCCGGCTCTT
      200      210      220      230      240      250
Command (U: up/D: down/P: print/RETURNing: Exit)?
TARGET FILE XENHIS2A.DNA 3121 748 FOUND 1 WAIT

```

Above: A reconstruction of Ichthyostega, the "missing link" between fishes and amphibia. Left: An example of a MacroGene homology search with DNASIS software and a CD-ROM disk. In the LKB MacroGene Workstation, there are no "missing links".

The short-legged fish on the facing page is Ichthyostega, the oldest known four-footed animal. Until fossils of this meter-long creature were found, he was one of evolution's missing links. Too bad there's none of his DNA to analyze! On LKB's new MacroGene Workstation, his evolutionary relationships could be evaluated quickly with the sequence information on the unique CD-ROM Laser Reference Disk.

Not only does the LKB MacroGene Workstation let you study intergenic relationships easily — it fully equips you to achieve optimum results at every step in DNA sequencing. From start to finish.

You can cast flawless ultrathin and wedge-shaped gels, as thin as 0.1 mm, in seconds. Maximize both resolution and the number of readable bases by drying your gels down to a very thin film. Carry out dideoxy reactions in a single dish. Load radioactive samples into 0.2 mm gel slots, with no risk of breaking a tip.

You can also read a longer DNA sequence rapidly and transfer it directly into your personal computer for analysis with user-friendly DNASIS* software.


DNASIS has all the functions needed to manipulate and evaluate large amounts of data in little time at all. It enables maximum matching, auto-connection of shotgun sequence fragments, control element search, automatic generation of restriction maps from gel data, reformatting of external files, and scores of other routines. And DNASIS alone offers you instantaneous access to the complete GenBank and NBRF databases, both on one laser disk.

MacroGene is the only DNA sequencing system you can buy which has no missing links. See for yourself! Ask your nearest LKB representative today for a MacroGene brochure and a DNASIS demonstration disk.

* DNASIS is a trademark of Hitachi Software Engineering Systems, Ltd. GenBank is a registered trademark of NIH.

** MacroMould is based on an invention of Dr. W. Ansorge. Rights are held by The European Molecular Biology Laboratory, Heidelberg, FRG.

There are no missing links in a MacroGene Workstation



For sequencing reactions: the 60-well LKB MicroSample Plate.




*For casting and moulding sequencing gels: the LKB MacroMould™ System.***



For loading: the LKB Sequencing Syringe, with unbreakable, exchangeable needle.



For heat control: the LKB Thermostatic Plate or Isothermal Plate.



For fast data input: the LKB MacroRead™ Digitizer and Light Table.



For analysis and comparison. DNASIS software and the compact laser disk.

LKB

Circle No. 172 on Readers' Service Card


Pharmacia LKB Biotechnology AB, 751 82 Uppsala, Sweden, Tel. 018-16 30 00, Telex 760 70, Telefax 018-14 38 20.

LKB products: Box 305, S-161 26 Bromma, Sweden, Tel. 08-799 80 00, Telex 10492, Telefax 08-98 63 64.

Australia (02)888 36 22 • Austria (0222)69 16 01 0 • Austria-Eastern European Countries +43 (222)92 16 07 • Belgium (02)242 4660
Canada (0514)684 88 81 • Denmark (02)26 52 00 • Finland (021)67 81 11 • France (01)64 46 36 36 • Greece-Middle East +30(1)894 73 96
Great Britain (0908)66 11 01 • Holland (034)80 77 911 Hong Kong-Far East (852)5-814 84 21 • India (0812)57 37 88 • Italy (02)253 913 21 • Japan (03)444-4811
Norway (02)54 90 95 • People's Republic of China (Beijing) 89 06 21 • Sovjet Union (095)255-6984 • Spain (03)230 05 01 • Sweden (08)799 80 00
Switzerland (01)821 18 16 • United States (0201)457 80 00 • West Germany (0761)490 30. Over 60 qualified representatives throughout the world.

Your Personal Copy of the New ICN Biomedicals Combined Reference Catalog Is Waiting.

**Ask for your
free copy
today.**



The ICN Biomedicals Combined Reference Catalog contains over 16,000 products in the following areas:
Monoclonal Antibodies, Blood Proteins, Film Badge and TLD Services, Enzymes, Human IgG Subclass Kits, Growth Factors and other Tissue Culture Reagents, Laboratory Accessories, NMR Solvents, Electrophoresis Reagents, Ultra Pures, Aluminas, Silicas, Animal Research Diets & Components, Pipetting Instruments, Gamma Counters, Immunodiagnostic Reagents, Automated RIA, Immunochemicals, Iodinated Products, Standards/Sources, Nuclides, Liquid Scintillation Cocktails, ^{35}S -, ^{32}P -, ^{14}C -, and ^3H Labeled Compounds.

ICN Biomedicals, Inc.

(800) 854-0530

P.O. Box 19536, Irvine, CA 92713



Circle No. 71 on Readers' Service Card

NIOSH ANNOUNCES RESIDENT & POSTDOCTORAL ASSOCIATESHIP AWARDS



FOR Ph.D & M.D. SCIENTISTS
The National Institute for Occupational Safety and Health announces outstanding opportunities for research in:

- APPLIED PSYCHOLOGY
- BIOACOUSTICS
- ERGONOMICS
- HEALTH PHYSICS
- EPIDEMIOLOGY

The National Research Council in cooperation with the National Institute for Occupational Safety and Health, Division of Behavioral and Biomedical Research in Cincinnati, Ohio will offer competitive research awards in 1988 for independent scientific research. Recipients will be free to work on problems of their own choice as guest investigators at NIOSH. This program is open to experienced senior scientists as well as to recent PhD and MD graduates. Stipends begin at \$27,150 for recent graduates and are appropriately higher for seniors. Relocation reimbursement and funds for professional travel on tenure are available, and a health insurance program is offered.

For further information and application materials please write to:

Associateship Programs (GF1, Room 424) NI2
NATIONAL RESEARCH COUNCIL
2101 Constitution Avenue, N.W.
Washington, D.C. 20418

Deadlines for application are
January 15, April 15 and August 15.


ENGINEERING FOUNDATION AND GENENCOR, INC.

are pleased to announce
the
**1987 Award in
Enzyme Engineering
To
Ephraim
Katchalski-Katzir**

**1985 – Klaus Mosbach
1983 – Ichiro Chibata**

ENGINEERING FOUNDATION, New York


deFonbrune MICROMANIPULATOR



- Exclusive pneumatic joystick control
- X, Y and Z axis movement
- Essentially **zero drift**

Circle No. 21 on Readers' Service Card


deFonbrune MICROFORGE



- 6 movement controls
- Accomplish contact fusion/melting, remote fusion, micro-forging, fracturing, micro-glass blowing, distant/contact stretching
- Forge/fabricate with glass, metals, alloys, organic crystals, polymers, other fusible materials

Circle No. 22 on Readers' Service Card

VIBRATOME® Series 1000



- Cuts fresh or fixed tissue evenly from 20 microns up
- Sections in fluid bath with vibrating blade
- Simple, precise adjustments

Circle No. 23 on Readers' Service Card


Detailed literature available.



**TECHNICAL
PRODUCTS
INTERNATIONAL, INC.**

13795 RIDER TRAIL, SUITE 104, ST. LOUIS, MO 63045
(314) 291-0198

PC3D™ three dimensional reconstruction software



Reconstruction courtesy of E. Effmann, M.D. and B. Smith, Ph.D.

Analyze and sort out complex structures (e.g. intertwined neurons). Components may be removed or replaced at will, allowing computer "dissection" of complex objects. Up to 400 polygons per section, 1000 serial sections per image. Applications in electron and optical microscopy, macroscopic sectioning, CAT scans, positron emission tomography, NMR imaging, autoradiographs, and more.

Reconstruct multi-colored three dimensional computer models of any object from serial section data. Display the reconstructed image, rotating and viewing it from any angle or orientation. Automatically calculate volumes of elements. Output 3-D images to plotter.

Free brochure **800-374-1888**
In CA **415-924-8640**
FAX: 415-924-2850 / Telex: 4931977
In Europe: R.J.A. Handers, Germany
Ph #02101/666268

Jandel
SCIENTIFIC

"Microcomputer Tools
for the Scientist"

65 Koch Road ■ Corte Madera, CA 94925

Circle No. 219 on Readers' Service Card



Forum '87: Students and Science Learning

20-21 November 1987
Hyatt Regency Crystal City

(Adjacent to National Airport, Washington, DC)

The **National Forum for School Science** provides an annual opportunity for teachers, scientists, science educators, and other individuals to analyze and discuss key issues, and to focus attention on the most enduring concerns and most promising new directions in science education.

Forum '87 Program

OPENING ADDRESS: What Do We Want Students to Learn?

PANEL: Perspectives on Science Learners — The speakers in this session will present some of the different ways of thinking about learners: as constructors of knowledge, as processors of information, and as products of social and cultural influences. How these different perspectives affect science education policy and practice will be the topic of the ensuing discussion.

LUNCHEON ADDRESS: Diversity Among Learners

PANEL: Information Needed to Develop Practice and Policy — Following a brief analysis by educators and policymakers, Forum participants will break into groups to discuss the types of information they need about students and learn-

ing. The focus will be on identifying research (such as demographic data, science aptitude, attitudes toward science, etc.) that will enable them to make informed policy and practice decisions.

PANEL: Environments and Strategies for Effective Science Learning — The environments and strategies for science learning can be based on our conceptions of science, or on our knowledge about how people learn. This panel will explore the differences and similarities between these designs, and ask if there are optimal systems. Particular attention will be paid to the practical considerations imposed by school structures.

SUMMARY DISCUSSION: Developing Practice and Policy That Adapt to Science Learning Goals

Register today by completing and returning the form on the next page. For more information, contact the AAAS Office of Science and Technology Education, 1333 H Street, NW, Washington, DC 20005; (202) 326-6620.

American Association for the Advancement of Science

Advance Registration Form

NS3

AAAS Forum '87: Students and Science Learning 20 – 21 November 1987 ♦ Hyatt Regency Crystal City, Arlington, VA (Adjacent to National Airport, Washington, DC)

Please Print or Type

Name _____
(last) (first and initial)

Affiliation _____

Mailing Address _____
(street and number)

(city) (state & zip) (telephone number)

☐ Check enclosed or charge to my ☐ VISA or ☐ MasterCard
(no other credit cards accepted)

Card No. _____ Expiration Date _____

Cardholder's signature _____

☐ Check here if you need special services due to a handicap. We will contact you before the meeting.

Registration Fees

\$125 Full (meals included) \$ _____

\$90 Partial (without meals) _____

\$ 45 Student (without meals) _____

Separate Meal Tickets

\$ 20 Lunch, Friday (20 Nov.) _____

\$ 20 Lunch, Saturday
(21 Nov.) _____

TOTAL AMOUNT: \$ _____

All registrants will receive a program, registration badge, and *This Year in School Science 1987: Students and Science Learning*. Registration is limited; preregistration form and payment **must be received by 30 October** to ensure receipt of advance materials. Registrations received after 30 October will be held at the AAAS registration desk in the Hyatt Regency. **Refund Policy:** If cancellations are received by 30 October, advance registration and meal ticket fees will be refunded after the Forum. **No refunds will be made on cancellations received after 30 October.**

Mail top half (registration form) to: AAAS, Forum '87 Registration, Room 830
1333 H Street, NW, Washington, DC 20005

Hotel Reservation ♦ Hyatt Regency Crystal City

AAAS Forum '87 ♦ 20 – 21 November 1987

Reservations received after 30 October cannot be guaranteed.

Send confirmation to:

Name _____ Street _____

City _____ State _____ Zip _____ Telephone No. _____

Other occupants of room: Name _____ Name _____

Special housing needs due to handicap: _____

Room: ☐ Single (\$75)* ☐ Double (\$85)* *Add 4.5% VA sales tax and 5% occupancy tax.

Arrival: Date _____ Time _____ **Departure:** Date _____ Time _____

Please list definite arrival and departure dates and times. Check-in time is 3:00 p.m.; check-out time is 12 noon. Reservations will be held until 6 p.m. only unless guaranteed by one of the credit cards below or one night's deposit is received. Failure to cancel 24 hrs. prior will result in one night's charges.


☐ American Express ☐ Diner's Club Number _____ Expiration Date _____

☐ Carte Blanche ☐ MasterCard Cardholder's signature _____

☐ Deposit (check) for \$ _____ enclosed. Make check payable to Hyatt Regency Crystal City.

Mail bottom half (reservation form) to: Reservations, Hyatt Regency Crystal City
2799 Jefferson Davis Hwy., Arlington, VA 22202

The right chromatograph is going to cost you less



\$10,965

Binary/ternary gradient programmer, programmable pump, variable detector, recorder, valve, column.

\$16,640

Dual programmable pumps, variable detector, valve and column. Includes ChemResearch™ data management/gradient software and color monitor option. Now with IBM System 2/Model 30 and printer.

Isco chromatographs offer top value. An Isco single-pump ternary or dual-pump binary gradient system complete with IBM computer can give you greater performance and capability than the two top-selling systems that cost much more. You'll get great value with any Isco HPLC system: analytical, microbore, or prep; stainless or inert.

With top performance. An Isco system gives you a precision pump that can be programmed to change flow rate based on either time or volume delivered. A

sensitive, stable 190-to-750 nm detector with flow cells down to 0.12 μ l. And a complete, powerful data management system that is actually easier to use than many simple integrators found in other chromatographs.

And top reliability. We've made over 50,000 detectors, pumps, and related instruments to earn our reputation for dependability. We back it up with the best system warranty offered anywhere—three years parts and labor.


Find out how much value Isco gives you by sending for a comparison of the leading HPLC systems. We'll also include our



80-page general catalog and 48-page HPLC catalog. Call toll free (800) 228-4250 (except NE, AK) or write: Isco, Inc., P.O. Box 5347, Lincoln, NE 68505.

Circle No. 182 on Readers' Service Card

HELLMA®



Your one stop for Cells and D₂ Lamps.

A full line of quality cells:

- Stock and custom
- Dye laser

Plane Optics

UV light sources for all major instruments:

- D₂ and Hg.

CUV-O-STIR

- Mini-magnet cuvette stirrer
- Variable speed
- Multiple cell capability.

HELLMA®

Box 544, Borough Hall Station, Jamaica, NY 11424
(718) 544-9534 or (718) 544-9166.

Circle No. 176 on Readers' Service Card

TMC MICRO-g® LAB TABLES

Higher performance...more options, more features

The new TMC Micro-g Series 63-500 lab tables with advanced Gimbal Piston® isolators achieve highest attenuation of both vertical and horizontal vibration. Features include new leveling feet, new 2" thick, high-stiffness, highly damped stainless steel top with rounded seamless edges and corners.

Options offer added utility and convenience: armrests, shelving, choice of top surfaces, casters, equipment support bars, Faraday cages. Sizes from 25" x 36" to 36" x 60". Write or call for new product bulletin No. MG30-86 and price list today.

Technical Manufacturing Corporation
15 Centennial Drive • Peabody, MA 01960
Telephone: 617-532-6330 • Telex: 951408



Micro-g® Series 63-500
vibration isolation tables

TMC-25

Circle No. 252 on Readers' Service Card

Try A New Approach to Life Sciences—Search STN International®!

Search STN International online to get effective answers to your life science questions. Stay current on new developments in fast-changing fields such as biotechnology, pharmaceuticals, forensic science, and immunology.

More than one-third of our CAS ONLINE® database contains life science information. Try BIOSIS Previews®, and our newest file, BIOCAS, the link between BIOSIS Previews and CAS ONLINE. All searchable with only one command language on STN!

Please fill out and clip this coupon today.

YES: Please rush me my FREE STN information packet.

Name _____

Organization _____

Address _____


Mail to:

STN International
Marketing, Dept. 33886
2540 Olentangy River Road
P. O. Box 02228
Columbus, Ohio 43202 USA



Circle No. 7 on Readers' Service Card

Consistent swirling action,
accurate results,
versatile
operation...



BURRELL WRIST-ACTION® SHAKER

As long as necessary, at selected speeds, the Burrell Shaker duplicates hand mixing's swirling action, assuring faster, more complete mixing. The swirling motion is consistent at every speed, so you can replicate exact operations.

Burrell's Shaker Build-Up® System lets you add separatory funnel supports, longer side-arms, platforms, and clamps. The Wrist-Action® Shaker can hold from 4 to 24 Erlenmeyer flasks, up to eight 250 ml flasks, and separatory funnels as large as 2,000 ml!

Write or call for our brochure on the only Wrist-Action® Shaker on the market.

BURRELL


BURRELL CORPORATION

2223 Fifth Avenue, Pittsburgh, PA 15219
Telephone: (412) 471-2527

Circle No. 154 on Readers' Service Card

DNASTAR PULSE

The FIELD INVERSION GEL
ELECTROPHORESIS CONTROLLER



AVAILABLE ONLY FROM DNASTAR

Under license from Field Inversion Technologies Corp.
St. Louis, MO (G. Carle, M. Olson patent pending)


DNASTAR

565-569 Chiswick High Rd.
London W4 3AY ENGLAND
(UK) 01-994-0619

1801 University Ave.
Madison WI 53705 USA
(608) 233-5525

Circle No. 233 on Readers' Service Card

COOL UNDER PRESSURE



SEM COLD STAGE

MMR Technologies proudly introduces the SEM Cold Stage System, a unique, SEM-attachable cold stage that allows you to study mounted samples at precisely controlled temperatures – **without** the use of liquid nitrogen.

The system includes MMR's K-20 Temperature Controller with IEEE-488 (GPIA) and RS-232 computer interfaces.


Other features include:

- Built-in microminiature Joule-Thomson refrigerator (80K-400K)
- Radiation shielded Cold Stage with 1cm² grounded cold pad
- Couples with standard SEM specimen mount; offers full rotation and tilt of stage
- 5 minute cooldown/2 watt capacity now available on ALL MMR refrigerators

MMR

MMR Technologies, Inc.
1400 Sterlin Road, #A-5
Mountain View, California 94043-1312, U.S.A.
Telephone (415) 962-9620 Telex 184817 Fax (415) 962-9647
England: AG Electro-Optics, Tarporley
Japan: ABE Trading Company, Ltd., Osaka
India: V.S. Scientific Instruments AG, New Delhi

Ask about our Hall System, Low Temperature Micro Probe for DLTS and X-Ray Diffractometer.



(weak anion exchange)

Analytical Conditions:

Column: 10 mm x 250mm, 15- μ m BAKERBOND WP-PEI (Cat. #7278-0)
Mobile Phase: A = 25mM Tris, pH 7.5
Gradient: B = 3 M NH₄OAc, pH 5.8
Flow Rate: 0% B to 100% B in 1.5 hours
Pressure: 8 ml/min
Detection: 50 psi
Sample: UV at 280 nm; 2 AUFS
Peaks: 48 mg total protein

1. Cytochrome c (Horse heart type VI)	8. α_1 -acid glycoprotein (Human)
2. α -Chymotrypsinogen	9. β -Lactoglobulin B (Bovine milk)
3. Carbonic anhydrase (Bovine erythrocyte)	10. β -Lactoglobulin A (Bovine milk)
4. Impurity (?)	11. Pepsin (Repurified isozyme)
5. Impurity (?)	12. Pepsin (Repurified isozyme)
6. Conalbumin (Bovine milk)	13. Pepsin (Repurified isozyme)
7. α -lactalbumin (Bovine milk)	14. Calmodulin (Bovine brain)

Turn each of your separations into a masterpiece...use BAKERBOND[®].

The BAKERBOND Universal Family of Chromatography Media offers different particle sizes of wide-pore silica using the same surface chemistry...to facilitate scale-up, obtain high recoveries, and increase throughput.


For technical assistance or ordering information on the entire BAKERBOND family of silica-based bonded phases for protein analysis and purification, call toll free 1-800-JTBAKER (1-800-582-2537); in NJ 201-859-2151; FAX 201-859-2865 Telex 299514 BAKR UR.




J.T. Baker Chemical Co.
 222 Red School Lane
 Phillipsburg, NJ 08865, USA

J.T. Baker Chemicals B.V., Rijsterborgherweg 20, 7400 AA Deventer, Holland, P.O. Box 1, Tel (5700) 87500. Telex 49072.
 *BAKERBOND and BAKERBOND Wide-Pore (WP) are trademarks of J.T. Baker Chemical Co.
 © 1987 J.T. Baker Chemical Co. All rights reserved.
 +D.R. Nau, BioChromatography (submitted)

Circle No. 239 on Readers' Service Card



Racemic propranolol separation on Chiralcel[®] OD



Racemic malic acid separation on Chiralpak[®] WE(-)

Introducing 3 new chiral columns from Daicel Chemical Industries

- CHIRALCEL OD for β -adrenergic blocking agents
- CHIRALPAK WE(-) for amino acid, malic acid separations
- CHIRALCEL CA-1 for a wide variety of aromatic compounds.

BAKERBOND Chiral Phase[®] products can be used to separate a wide variety of enantiomers.

O Type
 Comprised of optically active polymethacrylate (CHIRALPAK) or cellulose derivatives (CHIRALCEL). These highly resolving chiral phases offer unique separations of enantiomers with either non-polar or polar mobile phases.

CHIRALPAK WH, WM, and WE(-)
 Comprised of ligand exchangers. Provides resolution of many unmodified racemic α -amino acids.

PIRKLE-TYPE COLUMNS
 Comprised of DNBPG and DNBLeu. Provides the largest number of documented separations with non-derivatized and derivatized compounds.

For technical assistance or ordering information on the entire line of BAKERBOND Chiral Phase Products, call or write J.T. Baker Inc. 1-800-JTBAKER (1-800-582-2537); in NJ 201-859-2151; FAX 201-859-2865, Telex 299514 BAKR UR.




J.T. Baker Inc.
 222 Red School Lane
 Phillipsburg, NJ 08865, USA


J.T. Baker Chemicals B.V., Rijsterborgherweg 20, 7400 AA Deventer, Holland, P.O. Box 1, Tel (5700) 87500. Telex 49072. *BAKERBOND Chiral Phase is a trademark of J.T. Baker. Chiralcel[®] and Chiralpak[®] are registered trademarks of Daicel Chemical Industries. © 1987 J.T. Baker Inc. All rights reserved.

Circle No. 255 on Readers' Service Card


From the beginning,
 We have been
 your source for:



DNA Transilluminators



DNA Photodocumentation Systems



DNA Analysis Systems

Call 1-800-DNA-FOTO



Write or call today for a free copy of our 1987/88 catalog. In Wisconsin or outside the US, call 414-786-9550.

FOTODYNE
 INCORPORATED

16700 W. Victor Road
 New Berlin, WI 53151-4131
 Telex 260127

Circle No. 168 on Readers' Service Card

Designed to Freeze.
Monitored For Safety.
Priced To Sell.

NEW!
-135°C Bio-Freezer®

Bio-Freezers. Zero to -135°C From Forma

Now, there are even more reasons to invest in Forma ultra-low temperature technology.

Our new -135°C cryogenic preservation system provides a safe, economical alternative to liquid nitrogen for long term storage with better cell viability. Improved operating efficiencies and inventory configurations mean lower comparative storage costs per vial.

Plus, Forma's family of ultra-low temperature chest and upright freezers still give you an extra measure of performance, stability and reliability you can trust... all in a complete line for the broadest range of applications.

Operating Temperature	Cabinet Configuration	
	Chest	Upright
-135°C	X	
-120°C	X	
-85°C	X	X
-75°C	X	X
-40°C	X	X

Call or write for the details about Forma cryopreservation technology and refrigeration products... tested with patience and guaranteed to work longer, harder and quieter.

Free Offer!

A handcrafted silk rose for your time and interest in the energy efficient upright and chest freezers. Call us today, toll-free.



Call Us!



Forma Scientific
DIVISION OF MALLINCKRODT, INC.

BOX 649 MARIETTA, OH 45750 TELEX 29-8205

Toll-Free USA 800-848-3080

In Ohio 800-848-3080

In Canada 800-848-3080

Circle No. 5 on Readers' Service Card