

AMERICAN
ASSOCIATION FOR THE
ADVANCEMENT OF
SCIENCE

SCIENCE

5 MAY 1989

\$3.50

VOL. 244 ■ PAGES 505-620



The World's Most Effective Separation Tool

It's the Beckman TL-100 tabletop ultracentrifuge for rapid, efficient runs.

The proof is in the applications literature:

- **Isolate both DNA and RNA in one 3-hour CsCl gradient run in the TLS-55 rotor**
P. L. Iversen et al., *BioTechniques* 5, 521, 523-524 (1987)
- **Prepare rabbit reticulocyte fractions with four times the protein synthesis activity of conventional preparations in 20 minutes using the TLA-100.2.**
J. W. B. Hershey, Application Note DS-711. Beckman Instruments, Spinco Division, 1987.
- **Isolate DNA in 4-5 hours on CsCl-ethidium bromide gradients run in the TLA-100.2 rotor.**
D. P. Weeks et al., *Anal. Biochem.* 152, 376-385 (1986)
- **Pellet RNA through a CsCl cushion in 2.5 hours using the TLS-55 rotor**
M. G. Peterson and J. F. B. Mercer, *Eur. J. Biochem.* 160, 579-585 (1986)
- **Separate lipoprotein fractions by a single 25-minute spin in the TLV-100 rotor**
J. P. Segrest et al., *Beckman TL-100 News* 1:2, 1-2 (Dec. 1986)
- **Prepare a protein-free plasma fraction for drug binding studies by a 1-hour run in the TLA-100.2 rotor**
T. Zysset, *Ther. Drug Monit.* 9, 433-438 (1987)
- **Clarify solubilized membrane proteins in 15 minutes in the TLA-100.2 rotor**
P. H. Lin et al., *Biochemistry* 26, 731-736 (1987)
- **Separate lipoprotein fractions using single or serial 2.5-hour runs in the TLA-100.2 rotor**
H. K. Naito, Applications Note DS-693. Beckman Instruments, Inc., Spinco Division, 1986.
- **Estimate sedimentation coefficient of a protein from a 4-hour sucrose gradient run using the TLS-55 rotor**
D. L. Gard and M. W. Kirschner, *J. Cell Biol.* 105, 2203-2215 (1987)


Small wonder the TL-100 has become the fastest selling ultracentrifuge in the world!

For more information, including our new annotated applications bibliography with over 100 literature references, contact Beckman Instruments, Inc., Spinco Division, 1050 Page Mill Road, Palo Alto, CA 94304. (800) 742-2345. Offices worldwide.

BECKMAN

A SMITHKLINE BECKMAN COMPANY















Circle No. 95 on Readers' Service Card



NEW

TLA-45 fixed angle rotor holds 12 1.5-mL conical-tip tubes at a 45° angle and spins them at 45,000 rpm/125,000 g. Refrigerated microcentrifuge capability for the TL-100!

New! From OMEGA

CN1200 Series  Universal Digital Controller From \$250	Model HH110  Digital Thermometer \$85	Model CN310  Temperature Controller \$199	Model CL-307  Thermocouple Calibrator \$495	CN350 Series  1/16 DIN Temperature Controllers From \$99	CN9000 Series  Miniature Temperature Controller From \$199
Thermocouple  Wire Teflon® Coated Priced by the Spool	Model HH-82  Handheld Digital Thermometer \$169	TX90 Series  Miniature Transmitters From \$125	DP10 Series  Digital Panel Meter From \$425	Model 115  Digital Thermometer \$155	5-Dot Labels  Non-Reversible Temp. Labels \$19/10-pack
Model DP462  Temperature Indicator \$328	OS60 Series  Fixed Infrared Temperature Sensor From \$545	88000 Series  Surface Probes From \$110	SA1 Series  Self-Adhesive Thermocouples \$60/5-pack	CN380 Series  PID Auto-Tune Controllers From \$299	 All These Products and Many More! For FREE Handbook Call (203) 359-RUSH.

Widest Selection of Temperature Measurement and Control Products Available

OMEGA Engineering provides customers the world over with the technology and innovation they need to master their temperature measurement and control requirements. OMEGA Handbooks and Encyclopedias have long been the trusted source for quality sensors and instrumentation for Industry, Research, and Education.

Our knowledgeable sales representatives and product engineers support OMEGA's unyielding commitment to complete customer satisfaction. Wherever you are, ordering from OMEGA is easy

with toll-free sales and technical assistance during business hours coast-to-coast in the U.S. and Canada, 24-hour FAX ordering, and a full staff of international specialists.

For your custom design requirements, let our applications engineering staff show you how we can custom design the right product to your exact specifications.

Call Sales Toll-Free at
1-800-TC-OMEGA
 IN CT (203) 359-1660

Name _____
 Address _____
 City _____ State _____ Zip _____

Please send me the following FREE OMEGA Handbooks:

- ☐ Temperature ☐ pH
☐ Pressure ☐ Data Acquisition
☐ Flow ☐ Electric Heaters

Mail To
OMEGA ENGINEERING
 One Omega Drive, P.O. Box 4047
 Stamford, CT 06907
 (203) 359-1660

© COPYRIGHT 1989, OMEGA ENGINEERING, INC. ALL RIGHTS RESERVED, PRINTED IN USA

511 This Week in *Science*

Editorial

513 Drunk Driving and Statistical Morality

Letters

515 Planet Exploration: L. FRIEDMAN ■ Science Education: G. M. CARROW ■ Alar and Apples: T. H. JUKES; L. ROBERTS

News & Comment

517 Clean Air? Don't Hold Your Breath ■ The Political Battle Over Clean Air
520 NIH Reopens Baltimore Inquiry
521 Watson Floats a Plan to Carve Up the Genome
522 Utah Looks to Congress for Cold Fusion Cash ■ Cold Water from Caltech
524 NAS Elects New Members
Peoples Selected Head of Fermilab
Nobelists Back Animal Research

Research News

525 Space Science on the Rebound? ■ The Unveiling of Venus
527 Biologists Disagree Over Bold Signature of Nature
528 The Circle Can Be Squared!
529 Another Movement in the Dance of the Plates

Articles

541 Rewarding Performance That Is Hard to Measure: The Private Nonprofit Sector: B. A. WEISBROD
546 Studies of Inositol Phospholipid-Specific Phospholipase C: S. G. RHEE, P.-G. SUH, S.-H. RYU, S. Y. LEE

Research Articles

551 Interleukin-2 Receptor β Chain Gene: Generation of Three Receptor Forms by Cloned Human α and β Chain cDNA's: M. HATAKEYAMA, M. TSUDO, S. MINAMOTO, T. KONO, T. DOI, T. MIYATA, M. MIYASAKA, T. TANIGUCHI

Reports

557 Microearthquake Imaging of the Parkfield Asperity: P. E. MALIN, S. N. BLAKESLEE, M. G. ALVAREZ, A. J. MARTIN
559 Long-Range Electronic Perturbations Caused by Defects Using Scanning Tunneling Microscopy: H. A. MIZES AND J. S. FOSTER

■ **SCIENCE** is published weekly on Friday, except the last week in December, and with an extra issue in February by the American Association for the Advancement of Science, 1333 H Street, NW, Washington, DC 20005. Second-class postage (publication No. 484460) paid at Washington, DC, and at an additional entry. Now combined with **The Scientific Monthly**® Copyright © 1989 by the American Association for the Advancement of Science. The title **SCIENCE** is a registered trademark of the AAAS. Domestic individual membership and subscription (51 issues): \$75. Domestic institutional subscription (51 issues): \$120. Foreign postage extra: Canada \$46, other (surface mail) \$46, air mail via Amsterdam \$85. First class, airmail, school-year, and student rates on request. **Single copy sales:** Current issue, \$3.50; back issues, \$5.00; Biotechnology issue, \$6.00 (for postage and handling, add per copy \$0.50 U.S., \$1.00 all foreign); Guide to Biotechnology Products and Instruments, \$18 (for postage and handling add per copy \$1.00 U.S., \$1.50 Canada, \$2.00 other foreign). Bulk rates on request. **Change of address:** allow 6 weeks, giving old and new addresses and 11-digit account number. Authorization to photocopy material for internal or personal use under circumstances not falling within the fair use provisions of the Copyright Act is granted by AAAS to libraries and other users registered with the Copyright Clearance Center (CCC) Transactional Reporting Service, provided that the base fee of \$1 per copy plus \$0.10 per page is paid directly to CCC, 21 Congress Street, Salem, Massachusetts 01970. The identification code for *Science* is 0036-8075/89 \$1 + .10. **Postmaster:** Send Form 3579 to *Science*, P.O. Box 1722, Riverton, NJ 08077. *Science* is indexed in the *Reader's Guide to Periodical Literature* and in several specialized indexes.

■ The American Association for the Advancement of Science was founded in 1848 and incorporated in 1874. Its objects are to further the work of scientists, to facilitate cooperation among them, to foster scientific freedom and responsibility, to improve the effectiveness of science in the promotion of human welfare, and to increase public understanding and appreciation of the importance and promise of the methods of science in human progress.



COVER Adult mudpuppy (*Necturus maculosus*). This large aquatic salamander has been used frequently in physiological studies because cells throughout its body are exceptionally large. The eyes of this animal are used to show that glial cells in the retina normally buffer light-evoked changes in intraretinal potassium concentration. See page 578. [Photograph by W. B. Pavlik, Department of Psychology, University of Georgia].

- 562 Detection of C₅ in the Circumstellar Shell of IRC+10216: P. F. BERNATH, K. H. HINKLE, J. J. KEADY
- 564 Diode-Laser Absorption Spectroscopy of Supersonic Carbon Cluster Beams: The ν_3 Spectrum of C₅: J. R. HEATH, A. L. COOKSY, M. H. W. GRUEBELE, C. A. SCHMUTTENMAER, R. J. SAYKALLY
- 566 AP1/*jun* Function Is Differentially Induced in Promotion-Sensitive and Resistant JB6 Cells: L. R. BERNSTEIN AND N. H. COLBURN
- 569 Selective Amplification and Cloning of Four New Members of the G Protein-Coupled Receptor Family: F. LIBERT, M. PARMENTIER, A. LEFORT, C. DINSART, J. VAN SANDE, C. MAENHAUT, M.-J. SIMONS, J. E. DUMONT, G. VASSART
- 572 Murine MHC Polymorphism and T Cell Specificities: S. ROY, M. T. SCHERER, T. J. BRINER, J. A. SMITH, M. L. GEFTER
- 575 Interferon- α But Not AZT Suppresses HIV Expression in Chronically Infected Cell Lines: G. POLI, J. M. ORENSTEIN, A. KINTER, T. M. FOLKS, A. S. FAUCI
- 578 Spatial Buffering of Light-Evoked Potassium Increases by Retinal Müller (Glial) Cells: C. J. KARWOSKI, H.-K. LU, E. A. NEWMAN
- 580 Diffusible Factors Essential for Epidermal Cell Redifferentiation in *Catharanthus roseus*: B. A. SIEGEL AND J. A. VERBEKE
- 582 Commitment of Mouse Fibroblasts to Adipocyte Differentiation by DNA Transfection: S. CHEN, L. C. TEICHER, D. KAZIM, R. E. POLLACK, L. S. WISE
- 585 Synaptic Connections in Vitro: Modulation of Number and Efficacy by Electrical Activity: P. G. NELSON, C. YU, R. D. FIELDS, E. A. NEALE

AAAS Meetings

- 588 1990 Annual Meeting: Call for Poster Papers

Book Reviews

- 590 Proceedings of a Workshop on Optical Surveys for Quasars, *reviewed by* D. C. KOO ■ Seabed Pockmarks and Seepages, K. A. KVENVOLDEN ■ The Bioinorganic Chemistry of Nickel, M. K. JOHNSON ■ Some Other Books of Interest ■ Books Received

Products & Materials

- 593 Image Synthesis, Graphics Generation Software ■ Cell Culture Starter Kit ■ Transfer Solution for RNA Blotting ■ Disposable Filter for Media ■ Calculator Replaces Statistical Tables ■ Optical Micropositioner ■ Literature

Board of Directors

Walter E. Massey
*Retiring President,
Chairman*

Richard C. Atkinson
President

Donald N. Langenberg
President-elect

Mary Ellen Avery
Francisco J. Ayala
Floyd E. Bloom
Mary E. Clutter
Eugene H. Cota-Robles
Joseph G. Gavin, Jr.
John H. Gibbons
Beatrix A. Hamburg
William T. Golden
Treasurer

Richard S. Nicholson
Executive Officer

Editorial Board

Elizabeth E. Bailey
David Baltimore
William F. Brinkman
E. Margaret Burbidge
Philip E. Converse
Joseph L. Goldstein
Mary L. Good
F. Clark Howell
James D. Idol, Jr.
Leon Knopoff
Oliver E. Nelson
Helen M. Ranney
David M. Raup
Howard A. Schneiderman
Larry L. Smarr
Robert M. Solow
James D. Watson


Board of Reviewing Editors

John Abelson
Qais Al-Awqati
Don L. Anderson
Stephen J. Benkovic
Floyd E. Bloom
Henry R. Bourne
James J. Bull
Kathryn Calame
Charles R. Cantor
Ralph J. Ciccerone
John M. Coffin
Robert Dorfman
Bruce F. Eldridge
Paul T. Englund
Fredric S. Fay
Theodore H. Geballe

Roger I. M. Glass
Stephen P. Goff
Robert B. Goldberg
Corey S. Goodman
Jack Gorski
Stephen J. Gould
Richard M. Held
Gloria Heppner
Eric F. Johnson
Konrad B. Krauskopf
Charles S. Levings III
Richard Losick
Karl L. Magleby
Philippa Marrack
Joseph B. Martin
John C. McGiff
Mortimer Mishkin
Carl O. Pabo

Yeshayau Pocker
Michael I. Posner
Dennis A. Powers
Russell Ross
James E. Rothman
Erkki Ruoslahti
Ronald H. Schwartz
Vernon L. Smith
Robert T. N. Tjian
Virginia Trimble
Emil R. Unanue
Geerat J. Vermeij
Bert Vogelstein
Harold Weintraub
Irving L. Weissman
George M. Whitesides
Owen N. Witte
William B. Wood

Choose One. Or Choose Them All.



The IntelliGenetics Team Gives Protein Chemists and Molecular Biologists a Choice of Software *and* Hardware.

On computers ranging from PCs to Macs to mainframes, IntelliGenetics software products operate independently or linked together.

The IntelliGenetics Suite™

Running on VAX™ and MicroVAX II™ computers and Sun® workstations, the IntelliGenetics Suite offers multiuser access to programs and sequence databases. The speed and performance of these computers make the IntelliGenetics Suite the obvious choice for database searches, sequence comparisons, and assembling sequencing projects.

PC/GENE™

High resolution plots and tables display the results of over sixty nucleic acid and protein analysis programs. PC/GENE now features a powerful new sequence assembly program. Particularly strong in protein functionality, PC/GENE runs on IBM XTs, ATs, PS/2s, and compatibles. The latest releases of the EMBL and SWISS-PROT data banks are available on CD-ROM for use with PC/GENE.

GeneWorks™

To be introduced this year, GeneWorks is sequence analysis software for the Macintosh®. It offers a comprehensive range of functions—from restriction mapping, to structure and function prediction, to data base searching on CD-ROM. Taking full advantage of the Macintosh interface and graphics, GeneWorks gives scientists a flexible tool for analysis and publication.

Linked Systems

By linking IntelliGenetics microcomputer programs to the IntelliGenetics Suite on a VAX or Sun, you can distribute computing tasks and still take advantage of the larger machine's power. Together, these systems provide the most complete set of molecular biology programs available.

Solve your research problems by choosing software from the IntelliGenetics team. For more information about our products, please call or write.

IntelliGenetics Suite and GeneWorks are trademarks of IntelliGenetics, Inc. PC/GENE is a trademark of IntelliGenetics, Inc./Genofit, S.A. VAX and MicroVAX II are trademarks of the Digital Equipment Corporation. Sun is a registered trademark of Sun Microsystems, Inc. Macintosh is a registered trademark of Apple Computer, Inc.



IntelliGenetics, Inc.

700 East El Camino Real
Mountain View, CA 94040
415-962-7354

Circle No. 110 on Readers' Service Card

This Week in SCIENCE

C₅ in the stars

IN the envelope surrounding the star IRC+10216 there are molecules that consist of linear chains of five carbon atoms—C₅ (page 562). This star is the prototype of a mass-losing carbon star: such a star is a cooling red giant that, late in stellar evolution, expels material from its envelope. The infrared spectrum of IRC+10216 was recorded with ground-based instruments—telescope and spectrometer—at Kitt Peak National Observatory. A high resolution vibration-rotation spectrum attributable to C₅ was recorded and is in good agreement with laboratory measurements (see below). The data suggest that between each pair of atoms the bond lengths are approximately equal (1.283 angstroms). Bernath *et al.* speculate that the formation of C₅ and carbon chains of other lengths depends directly or indirectly on photochemistry: such chains may be generated from carbonaceous materials that have fragmented as a result of the absorption of photons and they may also be the building blocks of such materials. C₅ is a factor of 11 less abundant than C₃, which had earlier been detected in this star. The carbon chains appear to form in the cooler outer part of the circumstellar shell unlike other carbon compounds—CO, C₂H₂, HCN, and CS—that are present throughout the envelope.

C₅ in the laboratory

CARBON chains, rings, and clusters are found not only in association with stars but also as components of soot. Thus an understanding of pure carbon compounds is of interest both for astrophysical chemistry and for combustion chemistry. The infrared spectroscopic signature of C₅ was recorded in the laboratory by Heath *et al.* (page 564), confirming the identification of C₅ recorded from the star IRC+10216 (see above). The major difficulty of studying carbon clusters—generating large chains—was overcome through the use of supersonically cooled carbon; through laser vaporiza-

tion of graphite, C₅ chains were produced and their infrared rotation-vibration absorption spectra were recorded at high resolution. The chains appear to be linear and to be symmetric around a center; the data support the theoretical values of carbon-carbon bond lengths at approximately 1.283 angstroms. Carbon chains of other lengths should be amenable to generation and study by these procedures.

Combination drug therapy for AIDS

INFECTIONS with HIV-1 have acute and chronic phases. Thus Poli *et al.* recommend that a comprehensive clinical approach to controlling the AIDS virus should include attacks on the virus in both phases (page 575). They have studied how two drugs, AZT and IFN- α , affect HIV-1 replication in chronically infected cells in tissue culture. Two cell lines were used: one is a monocytic line (monocytes are known reservoirs for HIV) and the other is a T cell line (T cells are both virus reservoirs and virus “factories” during acute infections). Both cell lines normally produce small amounts of virus: the monocytes make defective noninfectious particles and the T cells produce infectious virions. However, both cell lines could be induced to produce large amounts of virus. IFN- α blocked the induced release of the viruses from each of the chronically infected lines; AZT did not. AZT, in contrast, prevented acute HIV infections in susceptible cells. Thus, if a proper balance of these drugs could be struck, it might be possible to block both the de novo infection of cells with HIV-1 and the release and spread of the virus from chronically infected cells.

Plant cell fates

MOST of the epidermal cells of plants are committed to an epidermal future fairly early in ontogeny. Exceptions to this rule, however, are the epidermal cells that form carpels in some plants; carpels are the

structures that enclose ovules of flowering plants. Siegel and Verbeke studied how carpel epidermal cells of Madagascar periwinkles redifferentiate to parenchymal cells (page 580); the two cell types are different in shape, in the density of their cytoplasm (epidermal cells have dense cytoplasm whereas parenchymal cells have many vacuoles), and in their planes of division. Through cell manipulations that exposed carpel epidermal cells to diffusible factors, the epidermal cells could be made to redifferentiate to parenchymal cells, whereas cells that had been removed from such exposure were stopped from redifferentiating. The diffusible factors could be trapped in agar barriers from which it should be possible to elute and characterize them.

Making neuronal connections

HOW nerve cells link up and fire during development affects how the electric circuitry of the nervous system will later operate (page 585). Nelson *et al.* show that early electric stimulation increases the strength and number of synapses (electrical connections) between neurons. In a three-chambered dish, sensory neurons in “peripheral” chambers sent axons beneath the barriers to link up with spinal cord neurons in the “central” chamber; in this way, nervous system circuitry was partially reconstructed in vitro. After the physical connections had been established, axons of the sensory neurons were electrically stimulated (or not) for several days. Stimulation promoted the attachment of axons to the central neurons; it also enhanced the amplitude of the excitatory potential in the central cells. Competition between stimulated and unstimulated sensory neurons that were linked to the same central cell apparently diminished the strength of the connection on the non-stimulated side. Thus, the developing nervous system is a pragmatic system: links that are actually made determine what the final functional architecture of the circuitry will be.

Here's why Gilson may be the only source you need for liquid handling, LC and HPLC instruments

Whether you're considering a new instrument or putting together a complete HPLC system, start with *Instruments for Liquid Chromatography and Sample Preparation*. You're likely to find everything you need in this new guide to Gilson instruments and systems. It presents an overview of the entire line—our well-known fraction collectors and Minipuls® peristaltic pumps, plus HPLC systems, detectors and automated sample processors. Here are a few examples.


A fraction collector with 98% problem-free performance

Most people associate Gilson with reliable fraction collectors. Many FC 80 and FC 100 micro-fractionators have been in use more than 15 years. But reliability doesn't stop with these instruments. Our newest collector, the FC 203, has a 98% problem-free performance record. A stationary rack system is one feature that assures reliable performance. Other features simplify operation. A special HELP key gives on-line instruction for time, drop or peak mode operation. And the FC 203 is just one of five Gilson fraction collectors.


A versatile detector with 254/280 nm built-in

Among the five Gilson detectors, one of the more versatile is the 112 UV/VIS detector. It has built-in 254 nm and 280 nm wavelengths to meet the most common detection requirements. And with optional lamp/filter assemblies, it can detect at 214 nm, 229 nm, and in the visible range of 340-640 nm. A choice of flow cells is offered for prep and analytical sample loads. Other Gilson detectors include variable wavelength UV, electrochemical, refractive index and fluorometric detectors.


Multi-mode fraction collector




UV/VIS detector



Automatic solid phase extraction system



Gradient auto-analytical HPLC system



A new system that automates solid phase extraction

The ASPEC system is one of the newest Gilson products. It automates sample cleanup with extraction columns...and then automatically injects the prepared sample into an on-line HPLC system. Other Gilson systems and instruments for sample prep and injection range from simple sample changers to large capacity auto-sampling injectors for up to 540 samples.


Eight complete HPLC systems to meet specific needs

Besides individual instruments, we also offer complete HPLC systems. Eight complete series are offered. These are based on gradient or isocratic elution, analytical or prep sample loads, automatic or manual injection. System controllers and software provide both control and data analysis. And peripheral components such as fraction collectors and detectors can be easily added. In fact, the modular system design facilitates upgrading at any time.

Request your free guide today

Instruments for Liquid Chromatography and Sample Preparation provides an overview of these and other Gilson instruments. It

includes general specifications to help you determine the suitability of an instrument for your application. For your free copy, call toll-free 1-800-445-7667 (in Wisconsin, 608-836-1551). Or mark the reader service number below.



GILSON

Gilson Medical Electronics, Inc., Box 27, 3000 W. Beltline Hwy., Middleton, WI 53562 USA, Tel: 608-836-1551, TLX: 26-5478, FAX: 608-831-4451
 Gilson Medical Electronics (France) S.A., 72 rue Gambetta, B.P. No. 45, 95400 Villiers-le-Bel, France, Tel: (1) 34.29.50.00, TLX: 606682, FAX: (1) 39.94.51.83

Circle No. 94 on Readers' Service Card

American Association for the Advancement of Science

Science serves its readers as a forum for the presentation and discussion of important issues related to the advancement of science, including the presentation of minority or conflicting points of view, rather than by publishing only material on which a consensus has been reached. Accordingly, all articles published in *Science*—including editorials, news and comment, and book reviews—are signed and reflect the individual views of the authors and not official points of view adopted by the AAAS or the institutions with which the authors are affiliated.

Publisher: Richard S. Nicholson

Editor: Daniel E. Koshland, Jr.

News Editor: Ellis Rubinstein

Managing Editor: Patricia A. Morgan

Deputy Editors: Philip H. Abelson (*Engineering and Applied Sciences*); John I. Brauman (*Physical Sciences*)

EDITORIAL STAFF

Assistant Managing Editor: Nancy J. Hartnagel

Senior Editor: Eleanor Butz

Associate Editors: Keith W. Brocklehurst, Martha Coleman,

R. Brooks Hanson, Barbara Jasny, Katrina L. Kelner, Edith

Meyers, Linda J. Miller, Phillip D. Szurmi, David F. Voss

Letters Editor: Christine Gilbert

Book Reviews: Katherine Livingston, *editor*; Susan Millius

This Week in Science: Ruth Levy Guyer

Contributing Editor: Lawrence I. Grossman

Chief Production Editor: Ellen E. Murphy

Editing Department: Lois Schmitt, *head*; Mary McDaniel,

Patricia L. Moe, Barbara E. Patterson

Copy Desk: Joi S. Granger, Jane Hurd, MaryBeth Shartle,

Beverly Shields

Production Manager: Karen Schools Colson

Assistant Production Manager: James Landry

Art Director: Yolanda M. Rook

Graphics and Production: Holly Bishop, Catherine S. Siskos

Covers Editor: Grayce Finger

Manuscript Systems Analyst: William Carter

NEWS STAFF

Correspondent-at-Large: Barbara J. Culliton

Deputy News Editors: Roger Lewin, Colin Norman

News and Comment/Research News: William Booth, Greg-

ory Byrne, Mark H. Crawford, Constance Holden, Richard A.

Kerr, Elliot Marshall, Jean L. Marx, Robert Pool, Leslie Rob-

erts, Marjorie Sun, M. Mitchell Waldrop

European Correspondent: David Dickson

Contributing Writer: John Walsh

BUSINESS STAFF

Circulation Director: John G. Colson

Fulfillment Manager: Ann Ragland

Business Staff Manager: Deborah Rivera-Wienhold

Classified Advertising Supervisor: Karen Morgenstern

Guide to Biotechnology Products and Instruments:

Sauna S. Roberts

ADVERTISING REPRESENTATIVES

Director: Earl J. Scherago

Traffic Manager: Donna Rivera

Traffic Manager (Recruitment): Gwen Canter

Advertising Sales Manager: Richard L. Charles

Marketing Manager: Herbert L. Burklund

Employment Sales Manager: Edward C. Keller

Sales: New York, NY 10036: J. Kevin Henebry, 1515 Broad-

way (212-730-1050); Scotch Plains, NJ 07076: C. Richard

Callis, 12 Unami Lane (201-889-4873); Chicago, IL 60914:

Jack Ryan, 525 W. Higgins Rd. (312-885-8675); San Jose, CA

95112: Bob Brindley, 310 S. 16th St. (408-998-4690); Dorset,

VT 05251: Fred W. Dieffenbach, Kent Hill Rd. (802-867-5581);

Damascus, MD 20872: Rick Sommer, 11318 Kings Valley Dr.

(301-972-9270); U.K., Europe: Nick Jones, +44(0647)52918;

Telex 42513; FAX (0647) 52053.

Information for contributors appears on page XI of the 31 March 1989 issue. Editorial correspondence, including requests for permission to reprint and reprint orders, should be sent to 1333 H Street, NW, Washington, DC 20005. Telephone: 202-326-6500. **Advertising correspondence** should be sent to Tenth Floor, 1515 Broadway, New York, NY 10036. Telephone 212-730-1050 or WU Telex 968082 SCHERAGO, or FAX 212-382-3725.

Drunk Driving and Statistical Morality

The picture pulls at the heartstrings: a grieving family at the side of a daughter in a coma. It appeared last week in the newspapers, as one more victim of a drunk driver was counted. She joins a growing group, including many killed and others seriously injured. The best estimates are that 50% of motor vehicle accidents, which kill more people each year than were killed during all of the Vietnam War, are caused by drunk drivers. In too many such cases the driver is a repeater whose license was not revoked for previous offenses, despite a law on the books, by a kindly but muddleheaded judge who did not understand statistics.

What may be called "statistical morality" is the precept that a given course of action that may cause some harm to individuals now will result in greater benefit to more individuals in the future. Tests of a new drug or vaccine on a limited number of volunteers, half of whom are given placebos, is one example. A tough law on drunk driving is another. It is inevitable that such a law will cause some hardship on offenders who are subjected to an automatic sentence, whether it be prison, suspension of a driver's license, or a large fine. Careful consideration of the punishment to make the penalty commensurate with the crime must be made at the time that the law is designed, but once the law is enacted, every judge who suspends a sentence should be aware that he or she erodes the deterrence effect and ensures that some future victim will pay for the mercy shown.

Mercy to the lawbreaker in the dock is frequently justified by the argument, "We can't do anything more to help the victim, so let's not be too mean to the criminal." In fact, the overwhelming purpose of punishment is to deter the next act, and not, as claimed by some who wish to confuse the issue and whip up sentiment for the criminal, an act of vengeance. The murderer, the drunk driver, the rapist—all appear in court after committing an act that cannot be undone. Society correctly says that, if convicted, they must pay the penalty to diminish the chances of the future murder, driving incident, or rape. Deterrence depends on the sureness of punishment. Unfortunately, there are far too many sitting on judicial benches who do not understand this and who rationalize that "in this case, the defendant will suffer undue hardship if the penalty is enforced." The argument of being kind to the violator today is often rationalized on the basis that one does not know the name and telephone number of the victim who will be injured and therefore it is not our responsibility. Such thinking is truly immoral because it allows a sanctimoniousness in the present to obscure the harsh reality that a future victim will certainly become a statistic. To some extent the future injury will be caused by the failure to implement a carefully thought-out law.

In some Scandinavian countries, penalties for drunk driving are so severe that revelers select in advance a driver who abstains from alcohol during an evening of merriment. In the United States there are tough laws on the books in many states, but these laws are so often overstepped by judges that the impression is created that they are easy to evade. That leads to exactly the situation that we have today; far fewer people here are deterred from drinking and driving than is the case in other countries.

Statistical morality implies that a painstakingly structured procedure will minimize the total victim count in our society. It requires that the initial plan be well considered but, once implemented, there should be no exceptions based on the false assumption that small deviations today are less serious than the consequences to the unnamed but statistically certain victims of the future. Statistical morality applies to tests for drugs and vaccines—for instance, the administration of a vaccine like the whooping cough. It is known that when this vaccine is given to 1 million children about 100 will have serious side reactions, including death. When public outcries led to a cancellation of the public requirement for vaccination in England, the number of deaths from whooping cough climbed so alarmingly that the vaccine requirement was reinstated. Statistical morality would say, "Require the test—more people will be saved." Those who believe that "we are only responsible for what we cause" would have said, "No. We are not responsible for the deaths caused by nature's scourges, only for those caused by the vaccine we required." Morality in the new age should state that the errors of omission are as important as the errors of commission. The more we know about risks, whether from drunk drivers or vaccines for whooping cough, the more we will require a statistical morality that ensures the least harm to the most people.

—DANIEL E. KOSHLAND, JR.

Three ways to give yourself more time for research:

Account 7-89918				
Sub	Sub-account Name	Budget	Expenses	Balance
1000	LAB SUPPLIES	40,000.00	7,331.25	32,668.75
2000	EQUIPMENT	15,000.00	5,865.00	9,135.00
3000	FACILITY SALARIES	35,000.00	9,775.00	25,225.00
4000	RESEARCH ASSOCIATE	15,000.00	14,662.50	337.50
5000	FRINGE BENEFITS	45,000.00	4,887.50	2,112.50
9000	OVERHEAD		6,353.75	38,646.25
Totals (Direct Costs).....		112,000.00	42,521.25	69,478.75
Totals (Indirect Costs)....		45,000.00	6,353.75	38,646.25

Keep track of your budgets with Grant Manager.

A no-nonsense accounting and ordering system for investigators and departments.


Type orders into your computer, and let Grant Manager print them on your university's order form. To make future ordering even faster, the program records your orders in vendor catalogs. Reorder with the press of a key!

Meanwhile, Grant Manager is tracking your account balances so that you know what's left to spend at any time -- not just once a month.

For IBM PC and compatibles, Macintosh 512KE, Plus, SE, II, or IIX. Requires 512K RAM and hard disk. \$425. Demo - \$10.

Project salary costs with Personnel Manager.

Working out the charging of people's salaries to your grants can be quite a puzzle. Salaries go up. People come and go. Grants have different time frames. You want to spend your grants, but not overspend.



Personnel Manager helps you solve this puzzle. It keeps a salary database, and enables you to play "what-if" so that you can quickly find the optimal allocation of people to grants over time. Works with Grant Manager.

For IBM PC and compatibles. Requires 512K RAM. \$425. Demo - \$10.

Paleo - Word	
The fossil in question was generally regarded to be the 220 million year-old remains of <i>Tyrannosaurus</i> . Several years ago, however, this hypothesis was overturned by the discovery of yet another species, <i>Neotyrannosaurus</i> , which is now the undisputed progeny (Argus, 1984 #1; Turbouse, 1987 #4) of the species at hand. While the true nature of this large mammal is questionable, the <i>Paleo Library</i> from central <i>EndNote</i> Showing 1143 references out of 1143	
Where the of	Billoski 1987 Triceratops extinction linked
mystery is	Billoski 1988 A stochastic model of the pr
that other r	Bowstein 1984 Micro-mammal prototypes in
Some think d	Jones 1976 The Theory of Foraging
Schwartz, 19	1977 Statistical confidence in pale
currently in	1978 Statistical confidence in pale
that the end	1979 Tooling ecosystems: role in
[Whiting, 19	Schwartz 1976 Tooling ecosystems: role in


Add bibliographic power to your word processor with EndNote.

Still typing bibliographies on your Macintosh? EndNote will do them for you.

You can store up to 32,000 references in EndNote's database. There's plenty of room for information. For example: up to 32,000 characters of abstracts for *each* reference! Searching is simple and extremely fast.

EndNote works with Microsoft Word, WriteNow, WordPerfect and MacWrite. Use the desk accessory to access the database while you write. Copy citations like [Smith, 1986] from your database and paste them into your text. When you finish writing, select a bibliographic style, and EndNote does the rest: It scans the paper, finds the citations you pasted, and adds a perfectly formatted bibliography to the end of your paper. It even replaces the in-text citations with numbers or (author, date).

For Macintosh 512KE, Plus, SE, II or IIX. \$129. Demo - \$5.




Niles & Associates
2200 Powell, Suite 765
Emeryville, CA 94608
(415) 655-6666

Trademarks: IBM: IBM Corp., Macintosh: Apple Computer Inc., Microsoft: Microsoft Corporation, WriteNow: licensed to T/Maker, Corp. MacWrite: Claris Corp. WordPerfect: WordPerfect Corp.

For Grant Manager or Personnel circle no. 8 on readers' service card.

For EndNote circle no. 9 on readers' service card.

Immobilize Your Resources



The new HiPAC™ Aldehyde Activated Kit makes it easy to create custom affinity supports.


The kit contains three Luer Tip Quick (LTQ™) Columns and an activated HPLC column—all packed with our high-performance HiPAC Aldehyde silica.

To determine optimal immobilization conditions, prepare an LTQ column and inject your protein or peptide ligand. Once the ligand has been immobilized to the silica, it only takes about 5 minutes to purify your sample.

When you're ready, transfer the immobilization protocol directly to the HiPAC analytical or prep-scale HPLC column.

The HiPAC Aldehyde Kit—it's the fastest way to get your separations moving.

For complete information, call toll-free.



ChromatoChem, Inc.
2837 Fort Missoula Road
Missoula, MT 59801
Call toll-free: 800/426-7227

Circle No. 149 on Readers' Service Card

BURN YOUR REFERENCE CARDS!

REF-11™

Computerizes your REFERENCES
and prepares your BIBLIOGRAPHIES

- ☐ Maintains a data base of references
- ☐ Searches for any combination of authors, years of publication, reference title, publication title, keywords or abstract
- ☐ Formats bibliographies exactly as you want them
- ☐ Reads your paper, inserts citations into the paper, and prepares a bibliography of the references cited (optional)
- ☐ Downloads references from MedLine data bases such as NLM, BRS and DIALOG (optional)

IBM PC/XT/AT, MS-DOS, CP/M 80 ... **\$195⁰⁰**

RT-11, TSX-Plus, RSX-11, P/OS **\$250⁰⁰**

VAX/VMS (native mode) **\$350⁰⁰**



ANY MANUAL **\$15⁰⁰**

ANY MANUAL & DEMO **\$20⁰⁰**

322 Prospect Ave., Hartford, CT 06106
(203) 247-8500


Connecticut residents add 7 1/2 % sales tax.

Circle No. 29 on Readers' Service Card

JENA

BIRTHPLACE OF
MODERN OPTICS

Fine, Infinity-corrected Instruments From JENA





JENAVAL

Transmitted Light Research

- Large wide field optical components with a built-in magnification changer for chrome-free super wide field index (25)
- Adjustable phase medium of contrast as needed
- DIK differential interference contrast, darkfield, polarizing accessories
- Photographic equipment utilizes large (25) field index

JENAVAL INTERPHAKO

- High measuring accuracy with Shearing, half-shade plate, interference fringe, and photomeric measuring techniques
- Also equipped for brightfield, positive & negative phase contrast, interference contrast, and interphako
- Photographic capabilities



Surgery Microscope

- Crystal-clear, infinity-corrected wide field optics with 3x-30x magnification
- Coaxial halogen illumination
- Monocular & binocular co-observation
- C-mount, TV, and 35mm adapters
- Counterbalanced floor stand

For more information on these and other JENA microscopes, call **1-800-444-7952**



SEILER

INSTRUMENT & MANUFACTURING CO. 170 E. Kirkham Ave. • St. Louis, Missouri 63119 • (314) 968-2282
See us at Booth #714 at FASEB

Dealerships Available

Circle No. 128 on Readers' Service Card

“We’re
used to
being
quoted.”

When scientists gather to talk science, they very often end up talking SCIENCE.

With a circulation almost 500% larger than the second-ranked journal, SCIENCE gets into more conversations because it gets into more hands, reaching more doctoral scientists in R&D, more NIH Grant recipients, more science policy-makers than any other publication in the field.

Just as the elite readership of SCIENCE responds to its editorial content, it responds to its advertising as well. If you’re involved in scientific products or recruitment advertising, the best place for that advertising is SCIENCE. You’ll enjoy being quoted. And acted upon.

• THE GLOBAL WEEKLY OF RESEARCH
SCIENCE

ADVERTISING REPRESENTATIVES
SCHERAGO ASSOCIATES, INC.
1515 Broadway, New York, N.Y. 10036 • (212) 730-1050

gangliosides

PURITY

**A wide range of gangliosides
of highest purity.**

For neurology and oncology research.

In large

quantities as well.

From BioCarb Chemicals.

*Please, circle the Readers' Service Card for more
information about gangliosides of highest purity.*

 **BioCarb
Chemicals**

BioCarb Chemicals AB S-223 70 Lund, Sweden
Phone + 46-46 18 25 23 Telefax + 46-46 15 13 95 Teletex 830 52 40

Interleukin Too.

We're introducing not only a new Interleukin-2 RIA, but also new Interleukin-1 α and Interleukin-1 β RIAs. They're all highly specific and highly sensitive. And only Advanced Magnetix has all three.


All are overnight assays that save washing and incubation time. All deliver the high performance and reliability that make Advanced Magnetix a leader in immunoassay.

We also offer a new highly specific Mouse IL-1 α antiserum for use in your own assays. Besides cytokines, we have a large selection of eicosanoids, phosphoinositides, and cyclic nucleotide assays, too. Satisfaction is guaranteed. Send for our complete catalog of research products.




Advanced Magnetix, Inc.
61 Mooney Street
Cambridge, MA 02138 U.S.A.
1-800-RIA-KITS 1-800-343-1346
617-497-2070 FAX 617-497-6927

For research only. Not for use in human clinical diagnosis.



Circle No. 72 on Readers' Service Card



Your life line.

BIOSIS Connection®

Keeping your finger on the pulse of biological and biomedical research information is as easy as pressing a button. The **BIOSIS Connection** is your low-cost direct line to the world's largest abstracting and indexing service for the life sciences, and you can be part of it!

This exciting and invaluable life sciences online end-user system connects you to all research and professional information needs, with current awareness databases covering journal

literature, patents, books, meetings, theses, AIDS research, upcoming events and employment openings.

All it takes to extend your life line is a modem-equipped personal computer or terminal and a phone line. And what's more, you need no special search skills to use the **BIOSIS Connection**!


Put your fingers on the life line... Subscribe now to the **BIOSIS Connection**!

BIOSIS Connection is a trademark of BIOSIS.



For more information, a Sign-Up Kit, or a free demonstration disk, contact BIOSIS, 2100 Arch Street, Dept. S589LL, Philadelphia, PA 19103-1399 USA. Telephone (215) 587-4800 worldwide; toll free 1-800-523-4806 (USA except PA); Telex 831739; Fax (215) 587-2016.

Circle No. 56 on Readers' Service Card



Precision control peristaltic pump

Tris™ is the peristaltic pump that's different. You get three channels, precision control with 2% reproductibility, robust metal construction, electronic intercommunication with your other instruments, a three-year warranty, and a price lower than you're asked to pay for someone else's single channel plastic model. The snap-out tubing cassettes give you uniform channel-to-channel flow and long tubing life—all without manual adjustment.

To order or get more information, phone toll free (800)228-4250. Or write Isco, Inc., P.O. Box 5347, Lincoln, NE 68505, USA.



Circle No. 75 on Readers' Service Card

The State-of-the-Art PCR Laboratory Manual!

PCR PROTOCOLS

A Guide to Methods and Applications

edited by

**Michael A. Innis, David H. Gelfand,
John J. Sninsky, and Thomas J. White**

Designed for use at the laboratory bench, this is the most comprehensive manual on PCR available today. Over 50 chapters provide precise instructions on procedures, with advice on primer design. All of the techniques described, from amplification and direct sequencing of genomic DNA through cDNA cloning and quantitation of mRNA, are tested, current, and supplemented with helpful notes and illustrations.

You'll also learn how to:

- Optimize novel applications
- Avoid many cumbersome molecular biological techniques
- Set up the laboratory to avoid contamination.

Affordable, well written, and self-contained—all biological scientists using PCR will find this spiralbound volume indispensable.

CONTENTS (Section Headings): Part I: Basic Methodology. Part II: Research Applications. Part III: Genetics and Evolution. Part IV: Diagnostics and Forensics. Part V: Instrumentation and Supplies.

Spiralbound: \$39.95 (tentative)

ISBN: 0-12-372181-4

Casebound: \$69.50 (tentative)

ISBN: 0-12-372180-6

July 1989, c. 320 pages



ACADEMIC PRESS

Harcourt Brace Jovanovich, Publishers

Book Marketing Department #36059

1250 Sixth Avenue, San Diego, CA 92101


CALL TOLL FREE 1-800-321-5068

Prices subject to change without notice. MC/ES #36059.
© 1989 by Academic Press, Inc. All Rights Reserved.

Circle No. 74 on Readers' Service Card

PURIFY MORE DNA

Purification of a
5'-DMT 36mer



PRP-1 Polymeric Reversed Phase Columns are ideal for purification of protected oligomers. They provide:

- **Superior Recovery-** Greater than 95% is recovered
- **High Capacity-** Load 10 mg or more of sample on a 15 cm column
- **Longer Life-** The polymeric packing easily outlasts silica based columns
- **pH Stability-** Purify oligomers with strong secondary structure at pH 12.7

For more information on PRP-1 HPLC columns for synthetic DNA purification, contact us today.

HAMILTON

Hamilton Company
P.O. Box 10030
Reno, Nevada 89520
Telephone: (702) 786-7077
TOLL FREE: (800) 648-5950


In Europe, circle no. 107
In the U.S.A., circle no. 108

SigmaScan™

measures areas, lengths
coordinates, angles,
slopes, and more —

IBM PC and Compatibles

NEW
VERSION 3.9



Quickly and accurately
digitize, measure and analyze
photomicrographs, strip charts,
X-rays, maps, and more — using
your own PC.

- Areas, lengths, angles, and slopes
- X,Y digitizing
- User-defined units
- Standard ASCII data output
(use in Lotus, dBase, etc.)
- Keyboard macros
- User-defined data transforms

Save hundreds of hours annually over
manual measurement techniques.
Automate complex analyses. Comes
complete with software, choice of
digitizing tablet, money-back guarantee
and full year hardware warranty.

Sigma-Scan™ software is also
available separately.

Free brochure **800-874-1888**
In CA **415-924-8640**

FAX: 415-924-2850/Telex: 4931977
In Europe: R.J.A. Handels GmbH, Germany
Phone: 2101/666268 FAX: 2101/64321

Jandel
SCIENTIFIC

"Microcomputer Tools
for the Scientist"

■ 65 Koch Road ■ Corte Madera, CA 94925

Circle No. 140 on Readers' Service Card

Astronomy & Astrophysics

Edited by
Morton S. Roberts

Selected from articles published in *Science* between 1982–1984, these 24 papers are divided into four parts: Solar System, Structure and Content of the Galaxy, Galaxies and Cosmology, and Instrumentation. The volume reveals a broad, coherent, and contemporary picture of our astronomical universe. Includes two Nobel Prize lectures. Fully indexed and illustrated. 1985; 384 pages.

Hardcover \$29.95
Softcover \$17.95
(AAAS members deduct 20%)

Order from AAAS Marketing,
Dept. Q, 1333 H Street, NW,
Washington, DC 20005.
Add \$1.50 postage and handling.

**American Association for the
Advancement of Science**

THE SECOND CAROLINA CONFERENCE ON PROTEIN ENGINEERING

The 1989 Conference
will be held September 8
through 11 at the Pine
Needles Resort,
Southern Pines
North Carolina



**THE
PROGRAM
IN
MOLECULAR
BIOLOGY &
BIOTECHNOLOGY**

KEYNOTE ADDRESS: BRIAN MATTHEWS (*Univ. of Oregon*)

SESSION TOPICS AND CHAIRS:

STRUCTURE: CRYSTALLOGRAPHY (*Charles Carter, Univ. of North Carolina*)

STRUCTURE: NMR (*Gary Pielak, Univ. of North Carolina*)

MOLECULAR DYNAMICS (*Jan Hermans, Univ. of North Carolina*)

MODELING AND VISUALIZATION (*Mike Carson, Univ. of Alabama*)

BACK TO ENGINEERING (*Joseph Bonaventura, Duke Marine Laboratories*)

POSTER SESSION (*Michael Cory, Burroughs Wellcome*)

ROUND TABLE DISCUSSION:

SUPERCOMPUTING TECHNOLOGY (*Steve Wallach, Convex Computer Corporation, co-sponsor of the conference*)

The total cost for room, board and registration is \$400. The number of participants is limited. Apply by **July 1** by writing a brief letter describing your research interests and their relevance to this meeting. For further information contact Jo Ann McPherson, The Program in Molecular Biology & Biotechnology, CB* 7100, 130 MBBRL, The University of North Carolina at Chapel Hill, Chapel Hill, NC 27599-7100.

There Are No Substitutes

We are Electron Microscopy Sciences, Fort Washington

You know us for our product quality, our reasonable prices, and our prompt service. There are no substitutes for these advantages. Now for the first time ever, Electron Microscopy Sciences is offering a complete line of chemicals, supplies and accessories for your EM, LM, SEM, Biological and Materials Science needs. Call or write for our catalog today.




Electron Microscopy Sciences

321 Morris Road, Box 251
Fort Washington, PA 19034
Toll-free: 1-800-523-5874
In PA: (215) 646-1566
Fax: (215) 646-8931
Telex: 510-661-3280

Circle No. 58 on Readers' Service Card

DIALOG. THINK OF IT AS A CATALYST.

CALL TODAY FOR
YOUR FREE SAMPLE



Dialog puts ideas in motion. It's the force that changes the way you look at problems and solutions.

We're the largest online knowledgebank in the world. Our databases have the chemical property, toxicity, patent and business information you need.

We have chemical substance, substructure and identification information online. Dialog also offers databases covering market data, products and companies.

But what really sets Dialog apart is how thoroughly we cover so many different industries from chemical to phar-

maceutical to petrochemical. We can make a big difference when a project demands more than chemical information. Since it's all on one service, you'll have what you need, quickly and easily.

With Dialog you get more than in-depth information. You get quality service and training.

We're reasonable, too. All you pay for is the time you spend using Dialog. There is no minimum usage requirement.

So when you need a catalyst to get ideas going, or to push them further, turn to Dialog. And watch how your mind reacts.

Circle No. 63 on Readers' Service Card

See us at The American Society of Microbiology tradeshow—Booth #s 135, 136.


Ask your librarian or information center if your company has DIALOG. If not call 800-3-DIALOG (800-334-2564) today. Or write to us at 3460 Hillview Ave., Palo Alto, CA 94304. 415-858-3785. Fax: 415-858-7069.

DIALOG® INFORMATION SERVICES, INC.
A Knight-Ridder Company

The world's largest online knowledgebank.

800-3-DIALOG


Dialog Information Service Inc., 1988. All rights reserved.
Dialog is a service mark of Dialog Services Inc. Reg. U.S.
Patent and Trademark office.



**SCIENCE has 153,192*
good reasons for being
number one.**

The pre-eminence of SCIENCE as the definitive life science journal starts with our circulation, a full 435% greater than our nearest rival! But that's only one reason for our supremacy. SCIENCE reaches more doctoral scientists, more NIH/Hughes Institute grant recipients, more U.S. science policy leaders than any other publication. For working scientists who must keep in touch with the life sciences universe, there is no choice but to read SCIENCE. For businesses with products or services to sell them, there is no choice but to advertise in SCIENCE.

*ABC Audit, 6/88 †BPA Audit, 12/87








**Number 2 has 35,192†
reasons for coming
in second.**

• THE GLOBAL WEEKLY OF RESEARCH
SCIENCE

ADVERTISING REPRESENTATIVES
SCHERAGO ASSOCIATES, INC.
1515 Broadway, New York, N.Y. 10036 • (212) 730-1050

from the pages of

SCIENCE



Neuroscience

Edited by Philip H. Abelson, Eleanore Butz, and Solomon H. Snyder

A collection of 27 *Science* articles on neuroscience research—ranging from genetic engineering to clinical therapy. Provides an integrative treatment of brain anatomy, physiology, and chemistry, and addresses fundamental questions concerning nervous system functioning. 1985; 454 pp.; hardcover \$29.95 (\$23.95 for AAAS members); softcover \$14.95 (\$11.95 for AAAS members); AAAS Publication #84-13.

Biotechnology: The Renewable Frontier

Edited by Daniel E. Koshland, Jr.

Covers the latest and most important topics at the forefront of biological research in 26 papers from *Science*. Illustrates the extraordinary cross-disciplinary aspects of modern biology and its tremendous impact on the future. Topics include virology, immunology, developmental biology, cancer, hormones, sensory phenomena, and more. 1986; 384 pp.; hardcover \$29.95 (\$23.95 for AAAS members); softcover \$17.95 (\$14.35 for AAAS members); AAAS Publication #85-26.

AIDS: Papers from Science, 1982-1985

Edited by Ruth Kulstad

Presents some of the most frequently cited research papers and news reports on AIDS, published in *Science* between 1982 and 1985. Chronicles how far AIDS research has come and provides an indication of the

direction in which it might go. An overview of AIDS research is provided by Myron Essex, chairman of the Department of Cancer Biology, Harvard University School of Public Health. 1986; 654 pp.; hardcover \$32.95 (\$26.35 for AAAS members); softcover \$19.95 (\$15.95 for AAAS members); AAAS Publication #85-23.

Frontiers in the Chemical Sciences

Edited by William Spindel and Robert M. Simon

Describes recent progress in all areas of chemistry, including new work in ultrafast reactions, chemical catalysis, organic synthesis, enzyme mechanisms, theoretical chemistry, analysis, and the chemistry of materials. Compiled from 38 *Science* articles. 1986; 624 pp.; hardcover \$29.95 (\$23.95 for AAAS members); softcover \$17.95 (\$14.35 for AAAS members); AAAS Publication #85-10.

Astronomy & Astrophysics


Edited by Morton S. Roberts

From our own solar system to the pulsars at the very edge of the cosmos, this collection of 24 articles from *Science* reveals a broad, coherent, and contemporary picture of our astronomical universe. Topics include galaxies and cosmology, the solar system, instrumentation, and the structure and content of the galaxy. 1985; 384 pp.; hardcover \$29.95 (\$23.95 for AAAS members); softcover \$17.95 (\$14.35 for AAAS members); AAAS Publication #85-4.

Order from: AAAS Marketing Dept. SM, 1333 H Street, NW, Washington, DC 20005. Please specify the appropriate AAAS Publication #(s), add \$2.00 postage & handling per order, and allow 4-6 weeks for delivery. VISA and MasterCard accepted. For shipments to CA, add applicable sales tax.

Published by the American Association for the Advancement of Science

Is protein chemistry an art or a science?



The art of protein chemistry is the experimental design. The science is the correct interpretation of results. Just as the protein chemist integrates art with science, Applied Biosystems integrates chemistry and instrument design with application expertise to become partners with you in developing solutions to your problems.

The complexity of the protein molecule poses a unique analytical challenge. Applied Biosystems' broad product portfolio—sequencers, synthesizers, amino acid analyzers,

HPLC, high performance electrophoresis, capillary electrophoresis—represents a few of the techniques needed for structural elucidation. Leading-edge instrumentation must be combined with clear insight and intuition to be successful.

Breakthroughs in protein chemistry derive from intuition and systematic analyses. Applied Biosystems recognizes this. We invite you, our partners in research, to share your viewpoint. Contact us today.



U.S.A. 850 Lincoln Centre Drive, Foster City, CA 94404. Tel: (415) 570-6667.
(800) 874-9868, in California (800) 831-3582. Telex: 470052 APBIO UI. Fax: (415) 572-2743.
Mississauga, Canada. Tel: (416) 821-8183. Fax: (416) 821-8246.
Warrington, U.K. Tel: 0925-825650. Telex: 629611 APBIO G. Fax: 0925-828196.
Weiterstadt, West Germany. Tel: 06151-87940. Telex: 4197318 Z ABI D. Fax: 06151-84899.
Paris, France. Tel: (1) 48 63 24 44. Telex: 230458 ABIF. Fax: (1) 48 63 22 82.
Milan, Italy. Tel: 02-835 4920. Fax: 02-832 1655.
Rotterdam, The Netherlands. Tel: (0) 10 452 35 11. Telex: 26249. Fax: (0) 10 452 91 57.
Burwood, Australia. Tel: (03) 288-7777. Fax: (03) 887-1469.
Tokyo, Japan. Tel: (03) 699-0700. Fax: (03) 699-0733.

For information circle reader service number 77
For sales call circle reader service number 78

CALL FOR POSTER PAPERS

AAAS ANNUAL MEETING

15-20 February 1990
New Orleans

Share your research with your colleagues across the scientific community at the world's most prestigious general science forum.

Gain visibility in your field by sharing your research with over 5,000 colleagues throughout the scientific community. This is an ideal opportunity, particularly for those developing their careers (such as entry-level scientists and graduate students), to get professional exposure at the world's most prestigious general science forum. Posters are informal, visually oriented presentations that uniquely allow for one-on-one discussions.

If you are a AAAS member or fellow, or are sponsored by one, submit an abstract of your research paper (as specified

below) by 9 November 1989. Appropriate topics span the physical, life, behavioral, and social sciences. Accepted abstracts will be published and distributed to all Meeting registrants. Each contributor will be assigned, by general discipline, to a poster session, and will be provided with a bulletin board on which to display graphics (charts, graphs, photos) and large, easy-to-read text for 1½ hours.

REGISTRATION: Poster presenters must be registered for the Meeting

Instructions

ABSTRACTS

We do not typeset abstracts. The copy you submit is reduced by one-third in size and printed. Adherence to the standard format (illustrated at right) is necessary to produce a readable publication.

- Abstracts must be no larger than a 5" square.
- Use only a typewriter or letter-quality printer on plain white paper (dot matrix print does not reproduce well).
- Indent, space, underline, and capitalize as in the example.
- Do not double-space the body of the text.
- Use reproducible black ink for all hand-lettering.
- Do not box abstract or paste it onto another piece of paper.

SUBMISSION

Above the 5" square, type the name of the broad discipline that encompasses the subject matter and provide 3 index words to describe the area within that discipline. Below and to the left of the square, type the full name, address, and phone number of the person to be contacted regarding status and scheduling. Below and to the right, type the name, affiliation, and complete membership number (appearing above the name on SCIENCE mailing labels) of the member or fellow endorsing the paper with his/her signature.

Send the original plus one copy of the abstract no later than 9 November 1989 to

**Contributed Papers, AAAS Meetings Office
1333 H Street, NW, Washington, DC 20005**

ABSTRACT DEADLINE: 9 November 1989

EXAMPLE:

5 inches (12.7cm)

Indent Five Spaces and Type Title in Upper and Lower Case Letters and Underline. AUTHOR'S NAME IN UPPER CASE (Institution Name in Upper and Lower Case), SECOND AUTHOR (Institution).*

(skip one line). Type abstract. The full width of the column of typed material should be 5 inches (12.7cm) and must not extend beyond that. The total length of the material, from top of title to bottom of footnotes must not exceed 5 inches (12.7cm). Abstracts which exceed these parameters will be returned. All special symbols and signs which must be hand lettered (e.g., π) should be rendered in reproducible black ink as clearly and carefully as possible. The entire submission should be of camera-ready quality so that it can be photographed, turned into a plate, and printed. The printed abstract will be about 2/3 the size of the typed version. Avoid paragraphing as this wastes space. However, you may use your allotted space to neatly letter in equations and diagrams as you deem necessary.

as indicated in this example.

(skip one line). Type footnotes.

AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE

Announcing... **HUMAN GENOME I**

**October 2 – 4, 1989
Town & Country Hotel • San Diego, CA**

Human Genome I, the conference on the largest biological project ever contemplated, is set for October 2-4, 1989 in San Diego, CA. It is an international conference on the status and future of research on the human genome.

Co-chaired By

Charles R. Cantor, Ph.D.
Director, Human Genome Center
Lawrence Berkeley Laboratory

Daniel E. Koshland Jr., Ph.D.
Editor of **SCIENCE**

PARTIAL LIST OF SPEAKERS FOR HUMAN GENOME I

Sydney Brenner
*Director, MRC Molecular Genetics Unit
Cambridge, England*

Eric Lander
*Fellow of the Whitehead Institute for Biomedical Research
Cambridge, Massachusetts*

Robert Moyzis
*Genetics Group Leader
Los Alamos National Laboratory*

Charles Cantor
*Professor of Molecular Biology
University of California, Berkeley
Director, Human Genome Center
Lawrence Berkeley Laboratory*

James Watson
*Director, Cold Spring Harbor Laboratory
Associate Director, NIH Human Genome Project*

Victor McKusick
*University Professor of Medical Genetics
Johns Hopkins University School of Medicine
President, International Human Genome Organization (HUGO)*

Francis Collins
*Associate Investigator, Howard Hughes Medical Institute
Chief, Division of Medical Genetics
University of Michigan Medical Center*

Tasuku Honjo
*Professor of Medical Chemistry
Kyoto University Faculty of Medicine*

Jean Dausset
*Professor, College de France
President, Human Polymorphism Study Center (CEPH)*

Hans Zachau
*Professor, Institute for Physiological Chemistry
University of Munich*

Ronald Davis
*Professor of Biochemistry
Stanford University School of Medicine*

Peter Pearson
*Chairman, Department of Human Genetics
Sylvius Laboratories, Leiden*

Allan Wilson
*Professor of Biochemistry
University of California, Berkeley*

Thomas Caskey
*Professor and Director, Institute for Molecular Genetics
Investigator, Howard Hughes Medical Institute
Baylor College of Medicine*

Russell Doolittle
*Professor of Chemistry and Biology
Center for Molecular Genetics
University of California, San Diego*

Cassandra Smith
*Associate Professor of Molecular Biology
University of California, Berkeley*

Peter Dervan
*Professor of Biochemistry
California Institute of Technology*

David Cox
*Associate Professor of Psychiatry
University of California, San Francisco*

George Church
*Assistant Professor of Genetics
Harvard Medical School*

Raymond White
*Investigator, Howard Hughes Medical Institute
University of Utah*

Norman Arnheim
*Professor of Molecular Biology
Ahmanson Center for Biological Research
University of Southern California*

Renato Dulbecco
*Distinguished Research Professor and Acting President
Salk Institute*

Michio Oishi
*Professor of Applied Microbiology
University of Tokyo*

PARTIAL LIST OF SESSION TOPICS

- Organization of the Human Genome Project
- State of our Current Knowledge
- Advances in Technology: New Methods for Mapping and Sequencing
- Progress in Interesting Regions of the Human Genome
- The Diversity of the Human Genome in Different Populations

Sponsored By
SCIENCE MAGAZINE

**A Publication of the American Association for the Advancement of Science
and**

The Human Genome Organization (HUGO)

The conference will be accompanied by an exhibition open to commercial firms offering or planning to offer appropriate equipment. The magnitude of the project and its requirements for the latest in scientific laboratory equipment, supplies and services will make the exhibit a significant contribution both to scientists and the companies that serve the scientific community.

SCHERAGO ASSOCIATES, INC.
1515 Broadway, New York, NY 10036
Tel: (212) 730-1050 • Fax: (212) 382-3725