

A microscopic image of cells, possibly from a developing embryo, showing various stages of cell division and differentiation. The cells are stained with a variety of colors, including bright red, green, blue, and yellow, creating a vibrant, almost abstract pattern. The background is a light, textured surface, and the overall composition is dense and detailed.

AMERICAN
ASSOCIATION FOR THE
ADVANCEMENT OF
SCIENCE

Science

12 FEBRUARY 1993
VOL. 259 • PAGES 865-1092

\$6.00

EVOLUTION OF ATMOSPHERES



**STELLAR
BREAKTHROUGH**

**INTRODUCING SPEED AND UNIFORMITY
IN TEMPERATURE CYCLING -
THE ROBOCYCLER™ 40 FROM STRATAGENE.**

- Reduce your cycling times by 30%
- Well-to-well temperature uniformity of $\pm 0.1^{\circ}\text{C}$

The RoboCycler™ 40 Temperature Cycler takes the most innovative approach ever to thermal cycling. While conventional systems ramp the temperature of a single block, the four-block RoboCycler 40 Temperature cycler uses a robotic arm to quickly move samples from one block to another. The results? A dramatic 30% reduction in cycle times. And unmatched well-to-well temperature uniformity of $\pm 0.1^{\circ}\text{C}$.

The RoboCycler 40 Temperature Cycler completely avoids the problem of temperature nonuniformity. Single-block cycler blocks develop hot spots because of the need to quickly ramp temperature. In contrast, temperatures of the RoboCycler 40 temperature cycler-well blocks remain constant throughout a cycling program. And upon completion, the samples are safely stored in the cold block.

The system's user-friendly interface simplifies programming. The hot blocks can be programmed from 25° to 99°C , the cold block from 6° to 20°C . Program information and run status are displayed on the LCD or printed to an optional serial printer.

The RoboCycler 40 Temperature Cycler...The answer when speed and consistency are critical.

Please contact Stratagene for the distributor nearest you.



STRATAGENE®

USA:
Corporate Headquarters
Telephone: 800-424-5444
Fax: 619-535-0034

Germany:
Stratagene GmbH
Telephone: (06221) 40 06 34
Telefax: (06221) 40 06 39

United Kingdom:
Stratagene Ltd.
Telephone: (0223) 42 09 55
Telefax: (0223) 42 02 34

France:
Stratagene France
Telephone: (0590) 72 36
Telefax: (1) 44 28 19 00

Switzerland:
Stratagene GmbH
Telephone: 01-3641106
Telefax: 01-3657707

STATISTICA/w

Windows,
DOS, Macintosh

STATISTICA/w™ (for Windows) Complete Statistical System with thousands of on-screen customizable, presentation-quality graphs fully integrated with all procedures ■ Complete Windows 3.1 support, DDE, OLE, TT-fonts, multiple toolbars, right mouse button support ■ Unlimited numbers of data-, results-, and graph-windows ■ Inter-window integration: data, results, and graphs can be treated as objects and converted into one another in a number of ways ■ The largest selection of statistics and graphs in a single system; comprehensive implementations of: Exploratory techniques; multi-way tables with banners (presentation-quality reports); nonparametrics; distribution fitting; multiple regression; general nonlinear estimation; stepwise logit/probit; general ANCOVA/MANCOVA; stepwise discriminant analysis; log-linear analysis; factor analysis; cluster analysis; multidimensional scaling; canonical correlation; item analysis/reliability; survival analysis; time series modeling; forecasting; lags analysis; quality control; process analysis; experimental design (with Taguchi); and much more ■ Manuals with comprehensive introductions to each procedure and examples ■ Hypertext-based Stats Advisor expert system ■ Extensive data management facilities (spreadsheet with long formulas, block operations, advanced Clipboard support, DDE hot links, relational merge, data verification, powerful programming language) ■ Batch command language and macros also supported, "turn-key system" options ■ All output displayed in Scrollsheets™ (dynamic, customizable, presentation-quality tables with toolbars, pop-up windows, and instant 2D, 3D and multiple graphs) ■ Extremely large analysis designs (e.g., correlation matrices up to 32,000x32,000, unlimited ANOVA designs) ■ Megafiler Manager with up to 32,000 variables (8 Mb) per record ■ Unlimited size of files; extended ("quadruple") precision; unmatched speed ■ Exchanges data and graphs with other applications via DDE or an extensive selection of file import/export facilities ■ Hundreds of types of graphs, including categorized multiple 2D and 3D graphs, matrix plots, icons, and unique multivariate (e.g., 4D) graphs ■ Facilities to custom design new graphs and add them permanently to menu ■ On-screen graph customization with advanced drawing tools, interactive stretching and resizing of complex objects, interactive embedding of graphs and artwork, special effects, icons, maps, multi-graphics management, page layout control for slides and printouts; unmatched speed of graph redraw ■ Interactive rotation, perspective and cross-sections of all 3D and 4D graphs ■ Extensive selection of tools for graphical exploration of data: fitting, smoothing, overlaying, spectral planes, projections, layered compressions, marked subsets ■ Price \$995.

Quick STATISTICA/w™ (for Windows) A comprehensive selection of basic statistics and the full graphics capabilities of STATISTICA/w ■ Price \$495.

STATISTICA/dos™ (for DOS) A STATISTICA/w-compatible data analysis system ■ Price \$795.

Quick STATISTICA/dos™ (for DOS) A subset of STATISTICA/dos statistics and graphics ■ Price \$295.

Domestic sh/h \$10 per product; 14-day money back guarantee.

Circle No. 22 on Readers' Service Card



StatSoft™

2325 E. 13th St. • Tulsa, OK 74104 • (918) 583-4149
Fax: (918) 583-4376

Overseas Offices: Statsoft of Europe (Hamburg, FRG), ph: 040/4200347, fax: 040/4911310; StatSoft UK (London, UK), ph: 0462/482822, fax: 0462/482855; StatSoft Pacific (Melbourne, Australia), ph: (03) 663 6580, fax: (03) 663 6117; StatSoft Canada-CCO (Ontario), ph: 416-849-0737, fax: 416-849-0918. Available From: CORPORATE SOFTWARE and other Authorized Representatives Worldwide: Holland: MAB Julsing, 071-230410; France: Version US (1) 40590913; Sweden: AkademiData 018-240035; Hungary: Dagent Kft 186-5782; Belgium: TEXMA 10 61 16 28; South Africa: Osiris 12 663-4500; Japan (Macintosh): Three's Company, Inc., 03-3770-7600; Japan (DOS and Windows): Design Technologies, Inc., 03-3667-1110.

StatSoft, STATISTICA/w, Quick STATISTICA/w, STATISTICA/Mac, Quick STATISTICA/Mac, STATISTICA/DOS, Quick STATISTICA/DOS, and Scrollsheet are trademarks of StatSoft, Inc.




STATISTICA/Mac™ (for Macintosh) A STATISTICA/w-compatible, comprehensive data analysis and graphics system designed for the Macintosh ■ Large selection of statistical methods fully integrated with presentation-quality graphics (incl. EDA, multiplots, a wide selection of interactively rotatable 3D graphs; MacDraw-style tools) ■ Unlimited size of files ■ Full support for System 7, incl. "Publish and Subscribe" ■ Price \$595.

Quick STATISTICA/Mac™ (for Macintosh) A subset of STATISTICA/Mac a comprehensive selection of basic statistics and the full graphics capabilities of STATISTICA/Mac ■ Price \$295.



T H E B E C K



Number 1 in a Series

CENTRIFUGATION

*Why do 7 out of 10 choose
L/XL Optima™ Series Ultracentrifuges?*

Dr. John Noti, Ph.D. of the Guthrie Foundation in Sayre, Pennsylvania, USA, puts it succinctly: *"The XL-90 is in a class by itself. The short run times we're achieving give us more time to perform actual experiments. Programming is simple and fast. This instrument is very reliable and service is, as always, excellent. A remarkable machine!"*

RUN TIMES SHORTENED BY AS MUCH AS 90%!

Beckman's ESP™ Efficient Sedimentation Program lets you use high performance rotors to full advantage without danger of precipitating gradient-forming salts. ESP designs protocols for optimum efficiency by simulating run conditions at the rotor's maximum speed. The software determines how the salts will redistribute during the run and automatically adjusts speed to avoid precipitation. By applying significantly higher forces at the early stages of the separation process, dramatic reductions in run times are achieved.

To simplify the set-up of subsequent runs, optimized protocols designed by ESP can be stored for future recall.

ENVIRONMENTALLY SENSITIVE - THE ONLY ULTRACENTRIFUGE WITH NO CFC'S

As the leading supplier of ultracentrifuges, Beckman's concern about the environmental impact of these instruments led to the development of the L and XL Optima Series. These instruments

are designed with no chlorofluorocarbons. Instead they have an air-cooled drive and solid-state thermoelectric cooling and heating modules.


Power consumption is reduced to 3100 BTU/hr with the L and XL — less than a third that of other ultracentrifuges. And the benchtop Optima TLX generates a low 2000 BTU/hr.

THE WORLD'S FIRST AND ONLY IMBALANCE-TOLERANT DRIVE

The driving force behind Optima Series instruments is Beckman's proven, vacuum-encased induction drive, now brought to a new level of stress tolerance. This drive is so resilient it eliminates the need for tedious sample weighing. Tubes can be conveniently balanced by eye to within $\pm 10\%$ or 5 mL, whichever is greater.

MORE USEFUL INNOVATIONS FROM BECKMAN

Beckman pushes the limits of centrifuge technology in a commitment to increase the efficiency of this essential laboratory process. Still the most widely used separation technique, continual Beckman innovation allows its forces to be applied in new ways.




ESP-designed protocol with VTI 90 Rotor reduced run times from 36 hours to just 2.5 hours.

For more information contact: Beckman Instruments, Inc., 2500 Harbor Blvd., Fullerton, CA, 92634-3100. (800) 742-2345. Service and distribution offices worldwide.

BECKMAN

Circle No. 18 on Readers' Service Card



NEWS & COMMENT

- Hughes' Tough Stand on Industry Ties** 884
Conflict Confusion: Five Views on Equity
Federal Conflict Rules Nearing Completion

- A Shot in the Arm for TB Research** 886

- NASA Stakes Its Reputation on Fix for Hubble Telescope** 887

- Scripps-Sandoz Deal Comes Under Fire** 889

- AIDS Research: Shalala Backs Reorganization** 889

RESEARCH NEWS

- Will Future Computers Be All Wet?** 890
Chemists Work Their Way From Pores to Riches

- Unexpected Intelligence Turns Up in a Cellular Gel** 893
Polymer Gels Get Smarter

- Mathematicians Gather to Play the Numbers Game** 894

- Cell Communication Failure Leads to Immune Disorder** 896
Tyrosine Kinase Defect Also Causes Immunodeficiency

SPECIAL SECTION

Evolution of Atmospheres

NEWS

- Is the Geological Past a Key to the (Near) Future?** 906

- Searching for Clues to Ancient Carbon Dioxide** 908

- An 'Outrageous Hypothesis' for Mars: Episodic Oceans** 910

- Earth Scientists Look NASA's Gift Horse in the Mouth** 912

- The Sad Saga of Small Satellites** 913

ARTICLES

- Atmospheric Evolution of the Terrestrial Planets** 915
D. M. Hunten

- Earth's Early Atmosphere** 920
J. F. Kasting

- The Ice Record of Greenhouse Gases** 926
D. Raynaud, J. Jouzel, J. M. Barnola, J. Chappellaz, R. J. Delmas, C. Lorius

- The Global Carbon Dioxide Budget** 934
E. T. Sundquist

DEPARTMENTS

- THIS WEEK IN SCIENCE** 873

- EDITORIAL** 875
Evolution of Atmospheres

- LETTERS** 877
Asteroid or Volcano: Have the Volcanists Been Heard?: D. M. McLean; D. E. Koshland, Jr. • Basic Research and Society: R. W. Briehl; R. G. Sachs • Correction: A Li-Fraumeni Syndrome p53 Mutation: D. Malkin and S. H. Friend

- SCIENCESCOPE** 883

- RANDOM SAMPLES** 898

- BOOK REVIEWS** 995
Revolutions in Mathematics, reviewed by J. McCleary
• *Containing the Atom*, J. M. Jasper • Vignettes: Book Promotion • Books Received

- PRODUCTS & MATERIALS** 998

AAAS Board of Directors

Leon M. Lederman
Retiring President,
Chairman
F. Sherwood Rowland
President
Eloise E. Clark
President-elect
Mary Ellen Avery
Francisco J. Ayala
Robert A. Frosch

Florence P. Haseltine
Alan Schriesheim
Jean'ne M. Shreeve
Chang-Lin Tien
Warren M. Washington

William T. Golden
Treasurer
Richard S. Nicholson
Executive Officer

John Abelson
Frederick W. Alt
Don L. Anderson
Stephen J. Benkovic
David E. Bloom
Floyd E. Bloom
Henry R. Bourne
James J. Bull
Kathryn Calame
C. Thomas Caskey
Dennis W. Choi

John M. Coffin
Bruce F. Eldridge
Paul T. Englund
Richard G. Fairbanks
Douglas T. Fearon
Harry A. Fozzard
Victor R. Fuchs
Theodore H. Geballe
Margaret J. Geller
John C. Gerhart
Roger I. M. Glass

Stephen P. Goff
Corey S. Goodman
Stephen J. Gould
Ira Herskowitz
Eric F. Johnson
Stephen M. Kosslyn
Michael LaBarbera
Charles S. Levings III
Harvey F. Lodish
Richard Losick
Anthony R. Means


Mortimer Mishkin
Roger A. Nicoll
William H. Orme-Johnson III
Stuart L. Pimm
Yeshayau Pocker
Dennis A. Powers
Ralph S. Quatrano
V. Ramanathan
Douglas C. Rees
Erkki Ruoslahti
Ronald H. Schwartz

Terrence J. Sejnowski
Thomas A. Steitz
Richard F. Thompson
Robert T. N. Tjian
Emil R. Unanue
Geerat J. Vermeij
Bert Vogelstein
Harold Weintraub
Zena Werb
George M. Whitesides
Owen N. Witte
Keith Yamamoto

COVER

Thin section of Antarctic ice (from a depth of 87 meters) observed under polarized light; the variety of colors reflects different orientations of the ice crystals. Air bubbles (black, about 1 millimeter across) provide a tiny record of ancient atmospheres that can be preserved in

the ice for more than 100,000 years. See page 926. Other aspects of atmospheric evolution are discussed in the special section beginning on page 905. [Photograph: Michel Creseveur]



PERSPECTIVES

- The Pluses of Subtraction 942
R. M. Myers

- RNA Polymerase Marching Backward 944
G. A. Kassavetis and E. P. Geiduschek

RESEARCH ARTICLE

- Cloning the Differences Between Two Complex Genomes 946
N. Lisitsyn, N. Lisitsyn, M. Wigler

REPORTS

- Existence of a Flat Phase in Red Cell Membrane Skeletons 952
C. F. Schmidt, K. Svoboda, N. Lei, I. B. Petsche, L. E. Berman, C. R. Safinya, G. S. Grest

- Photoinduced Polymerization of Solid C_{60} Films 955
A. M. Rao, P. Zhou, K.-A. Wang, G. T. Hager, J. M. Holden, Y. Wang, W.-T. Lee, X.-X. Bi, P. C. Eklund, D. S. Cornett, M. A. Duncan, I. J. Amster

- Electrochemical Processing of Conducting Polymer Fibers 957
S. Li, C. W. Macosko, H. S. White

- Colicin E1 Binding to Membranes: Time-Resolved Studies of Spin-Labeled Mutants 960
Y.-K. Shin, C. Levinthal, F. Levinthal, W. L. Hubbell

- The Secretory Granule Matrix: A Fast-Acting Smart Polymer 963
C. Nanavati and J. M. Fernandez

- Structure of the Thiamine- and Flavin-Dependent Enzyme Pyruvate Oxidase 965
Y. A. Muller and G. E. Schulz

- Deletion of IRF-1, Mapping to Chromosome 5q31.1, in Human Leukemia and Preleukemic Myelodysplasia 968

C. L. Willman, C. E. Sever, M. G. Pallavicini, H. Harada, N. Tanaka, M. L. Slovak, H. Yamamoto, K. Harada, T. C. Meeker, A. F. List, T. Taniguchi

- Anti-Oncogenic and Oncogenic Potentials of Interferon Regulatory Factors-1 and -2 971
H. Harada, M. Kitagawa, N. Tanaka, H. Yamamoto, K. Harada, M. Ishihara, T. Taniguchi

- "Infectious" Transplantation Tolerance 974
S. Qin, S. P. Cobbold, H. Pope, J. Elliott, D. Kioussis, J. Davies, H. Waldmann

- Carboxyl Methylation of Ras-Related Proteins During Signal Transduction in Neutrophils 977
M. R. Philips, M. H. Pillinger, R. Staud, C. Volker, M. G. Rosenfeld, G. Weissmann, J. B. Stock


- Altered Prevalence of Gating Modes in Neurotransmitter Inhibition of N-Type Calcium Channels 980
A. H. Delcour and R. W. Tsien

- Shutdown of Class Switch Recombination by Deletion of a Switch Region Control Element 984
S. Jung, K. Rajewsky, A. Radbruch

- An Adenovirus Vector for Gene Transfer into Neurons and Glia in the Brain 988
G. Le Gal La Salle, J. J. Robert, S. Berrard, V. Ridoux, L. D. Stratford-Perricaudet, M. Perricaudet, J. Mallet

- CD40 Ligand Gene Defects Responsible for X-Linked Hyper-IgM Syndrome 990
R. C. Allen, R. J. Armitage, M. E. Conley, H. Rosenblatt, N. A. Jenkins, N. G. Copeland, M. A. Bedell, S. Edelhoff, C. M. Disteche, D. K. Simoneaux, W. C. Fanslow, J. Belmont, M. K. Spriggs

960
Radical insertion





■ Indicates accompanying feature

■ **SCIENCE** (ISSN 0036-8075) is published weekly on Friday, except the last week in December, by the American Association for the Advancement of Science, 1333 H Street, NW, Washington, DC 20005. Second-class postage (publication No. 484460) paid at Washington, DC, and additional mailing offices. Copyright © 1993 by the American Association for the Advancement of Science. The title **SCIENCE** is a registered trademark of the AAAS. Domestic individual membership and subscription (51 issues): \$87 (\$47 allocated to subscription). Domestic institutional subscription (51 issues): \$205. Foreign postage extra: Mexico, Caribbean (surface mail) \$50; other countries (air assist delivery) \$95. First class, airmail, student and emeritus rates on request. Canadian rates with GST available upon request, GST #1254 88122. **Change of address:** allow 6 weeks, giving old and new addresses and 11-digit account number. **Postmaster:** Send change of address to *Science*, P.O. Box 2033, Marion, OH 43305-2033. **Single copy sales:** \$6.00 per issue prepaid includes surface postage; Guide to Biotechnology Products and Instruments, \$20.

Bulk rates on request. **Authorization to photocopy** material for internal or personal use under circumstances not falling within the fair use provisions of the Copyright Act is granted by AAAS to libraries and other users registered with the Copyright Clearance Center (CCC) Transactional Reporting Service, provided that the base fee of \$1 per copy plus \$0.10 per page is paid directly to CCC, 27 Congress Street, Salem, MA 01970. The identification code for *Science* is 0036-8075/83 \$1 + .10. *Science* is indexed in the *Reader's Guide to Periodical Literature* and in several specialized indexes.

■ The American Association for the Advancement of Science was founded in 1848 and incorporated in 1874. Its objectives are to further the work of scientists, to facilitate cooperation among them, to foster scientific freedom and responsibility, to improve the effectiveness of science in the promotion of human welfare, to advance education in science, and to increase public understanding and appreciation of the importance and promise of the methods of science in human progress.

³³P-nucleotides. Safety, sensitivity and resolution all in one.



Get high sensitivity and low background in *in situ* hybridizations.

[α -³³P]dATP (Cat. No. NEG-312H) was used to nick translate cDNA restriction fragments of various tubulin mRNAs. The resulting probes were hybridized to 10 μ m sections of rat dorsal root ganglia and frozen brain sections as described by Oblinger and Pickert (*DuPont Biotech Update*, October 1992) and subsequently visualized with NTB2 emulsion.

A Rat dorsal root ganglia were hybridized to a ³³P-labeled β -tubulin cDNA probe. Only the neurons (stained yellow) show labeling. The glia cells (stained purple) display no apparent expression of β -tubulin mRNA. Exposure time: 1 day.

B A ³³P-labeled β -tubulin cDNA probe was used to detect presence of β -tubulin mRNA in neurons from rat facial motor nucleus brain stem (heavily labeled large cells). The smaller, purple stained glia cells do not show any apparent labeling. Exposure time: 2 days. Data courtesy Dr. Monica Oblinger, Chicago Medical School.

New ³³P-nucleotides have β -emissions that are five-fold weaker than ³²P, so they can be handled routinely with no elaborate shielding. And with higher emission energies than ³⁵S, they require shorter exposure times. You get faster sequencing autoradiographs and maintain equivalent band resolution. ³³P-nucleotides are ideal for emerging technologies such as DNA cycle sequencing, single-stranded conformational polymorphism and *in situ* hybridization.

And like all NEN[®] nucleotides, new ³³P-nucleotides offer consistent results. Thanks to two DuPont advances – our patented radiochemical stabilizer,

Tricine, and our MicroSpec specific activity analysis. In fact, DuPont is the only manufacturer to accurately measure the specific activity of its nucleotides using microbore HPLC.

DuPont offers a choice of radiolabeled nucleotides as well as nonradiometric fluorescent and biotin-based labels, so you can choose the nucleotide best suited to your needs.

For additional information by fax, 24 hours a day, 7 days a week, call DuPont FaxBack[®] at 1-800-666-6527 (or 302-892-0616) and request document #2002.

United States 1-800-551-2121 • Canada 1-800-387-8391 • Australia +61 (08) 257149 • Belgium +32 (02) 724 2717 • Denmark +45 31506610 • France +33 (01) 69825450 • Germany +49 (06172) 87-2600 • Italy +39-2-25-302-1 • Japan +81 (03) 3224-8763 • Latin America/Asia Pacific Fax +1 (617) 426-2464 • Netherlands +31 (073) 206550 • Sweden +46 (08) 7503700 • Switzerland +41 (01) 8410330 • United Kingdom +44 (0438) 734865



Flat membranes

Two-dimensional polymer membranes exhibit a variety of phase transitions. A flat phase has been predicted in numerical simulations to have enhanced stiffness because it is roughened by thermal fluctuations, much like the way corrugations strengthen sheet metal. Schmidt *et al.* (p. 952) used small-angle x-ray and light scattering to study the flat phase in a system that models the classic tethered membrane—the spectrin skeletons of red blood cells. The calculated roughness of the skeletons from computer simulations agreed well with the experimental values.

Bucky polymers

Individual buckyballs may be joined in a polymer network by irradiating thin films of C_{60} with ultraviolet light. Rao *et al.* (p. 955) used several probes, including mass spectroscopy, Raman scattering, and scanning electron microscopy, to observe an irreversible transformation of the films into a different solid phase. Vacuum-deposited films of C_{60} were exposed to light from a mercury arc lamp or an argon-ion laser. Changes in vibrational modes were consistent with photopolymerization of the films. The C_{60} molecules may be connected by photochemically created $=C=C=$ bridges.

Electrochemical polymer processing

Conducting polymers, such as polyaniline and polypyrrole, have potential applications that include use in lightweight batteries and in display devices. Most of these materials are difficult to process in that they are insoluble and tend to decom-

Balancing act in tumor suppression

Interferons (IFNs) are proteins that can inhibit cell growth. The expression of IFNs is controlled by two mutually antagonistic transcription factors, interferon regulatory factors-1 and -2 (IRF-1 and IRF-2). Harada *et al.* (p. 971) show that a balanced expression of these two factors is critical to the maintenance of normal cell growth. When IRF-2 is overexpressed in NIH 3T3 cells, the cells become transformed and exhibit increased tumorigenicity in mice; simultaneous overexpression of IRF-1, however, causes the cells to revert to a more normal phenotype. The anti-oncogenic activity of IRF-1 is supported by the cytogenetic studies of Willman *et al.* (p. 968), who show that the gene for IRF-1 maps to human chromosome 5q31.1, a region that is frequently deleted in leukemia and preleukemic myelodysplasia.

pose rather than melt. Li *et al.* (p. 957) show that fibers of poly(3-methylthiophene) can be grown electrochemically by using a capillary tube to flow a solution containing the monomer past the electrode. By changing the capillary geometry, fibers of different shapes could be formed.

Moving into the membrane

Steps in the adsorption and unfolding of a protein at a membrane surface were followed by time-resolved electron paramagnetic resonance experiments. Shin *et al.* (p. 960) selectively spin-labeled the 20-kilodalton channel-forming fragment of a bacterial cytotoxin, colicin E1, with nitroxide compounds. The fragment likely consists of eight amphipathic α helices that surround two very hydrophobic α helices. The signal from spin-labeled residues in solvent-exposed surface loops of colicin E1 was attenuated within a few seconds after mixing with phospholipid vesicles, reflecting the loss of rotational mobility after adsorption. In contrast, the signal from a residue within one of the buried hydrophobic helices had a time constant greater than

100 seconds, reflecting the rate-limiting insertion of the helix into the membrane.

Cofactor binding in pyruvate oxidation

Pyruvate plays a central role in both aerobic and anaerobic metabolism. Muller and Schulz (p. 965) describe the structure of the tetrameric enzyme pyruvate oxidase from *Lactobacillus*, which oxidatively decarboxylates pyruvate to produce the energy-storage metabolite acetylphosphate. One of the cofactors, thiamine pyrophosphate (TPP), is bound by a metal ion and a $\beta\alpha\beta$ unit. The spatial relation between the TPP cofactor and the flavin adenine dinucleotide (FAD) cofactor suggests that oxidation proceeds through a two-step transfer of two electrons.

Ras relations

Posttranslational modifications of Ras-related proteins are required to target these proteins to the plasma membrane where they function in signal transduction. Philips *et al.* (p. 977) find that the carboxyl methylation of Ras-related low molecular weight GTP-binding pro-

teins is stimulated after the activation of neutrophils with the chemoattractant FMLP. Blockage of neutrophil response to FMLP by inhibitors of prenyl cysteine carboxyl methylation suggests that receptor-mediated signal transduction is linked to carboxyl methylation of Ras-related proteins.

Regulating switching

During B cell development, the antigen-binding variable part of the immunoglobulin gene can bind to different constant regions, a process known as class switching. Two reports provide insights into factors controlling switching. Allen *et al.* show that (p. 990) defects in a T cell protein block switching in B cells and induce an inherited immunodeficiency disease (see also news report by Marx, p. 896). Jung *et al.* (p. 984) used mutant mice to identify a DNA sequence that regulates a particular switch combination.

Gene transfer in the brain

Viral vectors can be used to transfer genes into differentiated cells that no longer undergo cell division, but in the case of cells in the nervous system, expression from transferred genes often declines within a few weeks. Le Gal La Salle *et al.* (p. 988) show that a replication-defective adenovirus vector can be used to introduce a gene encoding β -galactosidase into the brain. Injection of the vector into the hippocampus and substantia nigra of the rat-infected neurons, microglial cells, and astrocytes without apparent cytotoxic effects. Neuronal expression of β -galactosidase was stable for 2 months.

With Chemiluminescent
Sample Detection

Major Breakthrough!

Phosphor Image Analysis for Molecular Biology



New Model GS-250 Molecular Imager™ System. The Model GS-250 Molecular Imager system represents a breakthrough in rapid imaging of radioactivity and chemiluminescence. Now you can image your gels, dot blots, slot blots, or any radioactive or chemiluminescent sample quickly, without the need for x-ray film, intensifiers, or low temperature film storage and exposure.

The Model GS-250 Molecular Imager system offers the following advantages over traditional autoradiography:

- Reduced Exposure Times
- Increased Sensitivity
- More Accurate Quantitation
- Erasable and Reusable Storage Phosphor Screens
- Lower Cost per Image

(Microsoft, and Windows, are registered trade marks of Microsoft Corporation)

How it Works:

Simply expose your radioactive or chemiluminescent sample to the proprietary storage phosphor screens. (The screens are sensitive to not only ^{32}P , but also ^{35}S , ^{14}C , ^{125}I , X-rays, AMPPD, and others.) In minutes, you will have generated a latent image which can be quickly scanned by the patented fiber optic laser scanner. Displaying, analyzing, and outputting the resulting data from your 1-D or 2-D protein and DNA gels or blots is fast and simple.

High Performance Software: Great technology is wasted without great software! We offer sophisticated software packages for both Windows® compatible and Macintosh® computers; you choose the graphical user interface you prefer. Bio-Rad software will let you quickly quantify your DNA, RNA, or protein samples; perform molecular weight analysis; plot 1-D density profiles or histograms; and perform a myriad of image viewing and processing manipulations.

For more information call 1-800-4BIORAD.

(Macintosh is a registered trademark of Apple Computers, Inc.)

BIO-RAD

**Bio-Rad
Laboratories, Inc.**

Life Science
Group

U.S. (800) 4BIORAD • California Ph. (510) 741-1000 • New York Ph. (516) 756-2575 • Australia Ph. 02-805-5000 • Austria Ph. 0222-877 89 01 • Belgium Ph. 091-85 55 11 • Canada Ph. (416) 624-0713 • China Ph. 2563146 • France Ph. 01-49 60 68 34 • Germany Ph. 089-318 84 0 • Italy Ph. 02-21609 1 • Japan Ph. 03-3534-7515 • Hong Kong Ph. 7893300 • The Netherlands Ph. 08385-40666 • New Zealand Ph. 09-443 3099 • Scandinavia Ph. 46 (0) 8 590-73489 • Spain Ph. (91) 661 70 85 • Switzerland Ph. 01-810 16 77 • United Kingdom Ph. 0800 181134

Circle No. 9 on Readers' Service Card

LIGHT BASED DETECTION OF BIOLOGICAL MOLECULES

ECLTM

Light Years Ahead in Western Blotting

Western blotting



Greater sensitivity than any
other technology -
results in minutes,
often seconds!

The Proof

For Western blotting applications, ECLTM systems are simply the best! - but don't just take our word for it. After just two years, ECL has become the most referenced light-based system, proving it works.

No other system is as sensitive or as fast as ECL.

- Orders of magnitude more sensitive
- Results in minutes, often seconds
- Simple multiple reprobing
- Hard copy results

Let ECL Western blotting systems help you get ahead.

Extensive bibliography -
the proof that it works



Amersham International plc
Amersham England HP7 9NA
Tel: 0494 544000

For further information
contact your local office.

Australia
Sydney (02) 888-2288
Belgium
Gent (091) 41-52-80
Canada
Oakville (416) 847-1166
Denmark
Birkerød (45) 82-02-22
Far East
Hong Kong (852) 802-1288

France
Les Ulis (1) 69-1828-00
Germany
Braunschweig (053) 07 2 06-0
Iberica
Madrid (91) 304-42-00
Italy
Milan (2) 5097740
Japan
Tokyo (03) 38-16-6161

Nederland
Hertogenbosch (073) 41-85-25
Norway
Gjettum (02) 54-63-18
Sweden
Solna (08) 734-0800
UK Sales
Aylesbury (0296) 395222
USA Arlington Heights
IL (800) 323-9750


ECL, ECL+device and Amersham are trademarks of Amersham International plc

Circle No. 2 on Readers' Service Card

© Amersham International plc 1992. All rights reserved.

FREE PEPTIDE*

from the world's price leader in peptide products for the research community
...Chiron Mimotopes U.S.



PEPTIDES \$150

Quality custom peptides at a fraction of the cost of conventional suppliers. Full range of quantity and purification levels. Individual or multiple peptide synthesis, peptide-conjugates, anti-peptide sera, multi-synthesis kits.

Also offering contract testing, epitope mapping and other services.
And now two offices to serve you.

What better reason to Celebrate ... Coast to Coast!

Call or write today

Order Toll Free:

800-733-7025

West Coast

800-633-8161

East Coast

* Free peptide offer limited to one free peptide per customer (≤ 15 -mer in length, 1mg \geq 70% purity) on next purchase of three or more peptides at our regular low price. Expires March 31, 1993.



CHIRON MIMOTOPES U.S.

TODAY'S LEADER IN PEPTIDE PRODUCTS AND SERVICES.


West Coast: 4560 Horton Street, Emeryville, CA 94608. Tel: (800) 733-7025 or (510) 601 3316. Fax: (510) 601-3300.

East Coast: 172 Mine Lake Ct., Suite 100, Raleigh, NC 27615. Tel: (800) 633-8161 or (919) 848 4323. Fax: (919) 870-9352.

Project CM 0555.1.93

Circle No. 21 on Readers' Service Card

A WORLD OF MICROMANIPULATION AND MICROINJECTION.



CHOICE. FLEXIBILITY. EXPERIENCE.

Carl Zeiss offers the most complete selection of state-of-the-art manipulation and cell injection workstations available. Whether you have a Zeiss inverted, upright or stereo microscope, we have the right fit for you.

- Each system is optimized for ergonomic, high performance operation.
- Equipment from many manufacturers is available including Narishige and Eppendorf.
- Complete training and technical support are provided.

Why compromise your research and use anything less? Zeiss has a micromanipulation/microinjection workstation to fit every application. And every budget.

For further details, contact your Zeiss representative. Or call 800-233-2343; fax 914-681-7446.

Carl Zeiss, Inc.

Microscope Division
One Zeiss Drive
Thornwood, New York 10594
800-233-2343 • Fax 914-681-7446



Brightfield
Microscopy

Darkfield
Microscopy

Phase
Contrast
Microscopy

Fluorescence
Microscopy

Differential
Interference
Contrast
Microscopy

Polarizing
Microscopy

Video
Microscopy

Confocal
Microscopy

Microscope
Photometry

SEE
WHAT
YOU'RE
MISSING

Talent runs in our family.

NEW Fluorescence Systems from SLM-AMINCO.



- **Photometry**
- **Microphotometry**
- **Fluorescence Lifetime Measurements**
- **Quantitative Ratio Fluorescence Imaging**

Our complete line of research fluorescence measurement systems covers the widest range of life science applications. Wavelength and time drive scanning, T-Optics™ polarization measurements, intracellular studies using fluorescent probes, and remote fiber optic spectroscopy are just a few examples of the many techniques these systems can employ.

With a number of mechanical and optical design enhancements, these systems offer improved stray light rejection and greater sensitivity. The new mouse-driven multitasking software combines exceptional power and flexibility with incredible ease of use. New software features include real time 3D acquisitions, quantitation, and both preset and custom routines for intracellular probe measurements.


We're proud of this family, and we'd like to tell you more. Call or write to us today for new literature.

SLM-AMINCO®

A Milton Roy Company

810 West Anthony Drive, Urbana, IL 61801 USA
(217) 384-7730 (800) 637-7689 FAX: (217) 384-7744

SUTTER INSTRUMENT COMPANY



PURVEYORS OF FINE GLASS TUBING

Offering a wide variety of sizes in quartz, borosilicate, and aluminosilicate compositions for micropipette fabrication. Available with or without internal filament.

Quality, Consistency, and Availability you can count on—from the leader in micropipette technology.



Call us for a complete listing.

Sutter Instrument Company

40 Leveroni Court, Novato, CA 94949 Ph: (415) 883-0128 Fax: (415) 883-0572

Circle No. 3 on Readers' Service Card

CUSTOM DNA SYNTHESIS

PURE & SIMPLE
(and now more economical than ever)

☐ PRICE \$ **2.50**
+ per base plus setup

☐ SUPERB TECHNICAL SUPPORT
+
☐ IMPECCABLE QUALITY
+
☐ WORLD'S FASTEST SERVICE

MIDLAND

THE UNDISPUTED #1
CUSTOM DNA SYNTHESIS SERVICE

THE MIDLAND CERTIFIED REAGENT COMPANY
3112-A WEST CUTHBERT AVENUE
MIDLAND, TEXAS, 79701

PHONE 1-800-247-8766 FAX 1-915-694-2387

Circle No. 7 on Readers' Service Card

ATTOFLUOR™ REAL-TIME DIGITAL FLUORESCENCE ANALYZER



FASTER. SIMPLER. MORE SOPHISTICATED.

State-of-the-art computerized workstation ideally suited for kinetic measurements of calcium, pH, cyclic AMP and other ions.

- Zeiss AXIOVERT for exceptional fluorescence transmission.
- Highly sensitive, compact intensified CCD camera.
- Vibration-free excitation filter changer for up to four ratios per second.
- Powerful 486 class computer.

ATTOFLUOR provides effortless image acquisition, processing, archiving and data file manipulation. It's offered as an integrated turnkey system, or as individual components.

For immediate detailed information and pricing, contact Roxanne Moosher at 914-681-7653, fax 914 681-7446. Or call 800-233-2343.

Brightfield
Microscopy

Darkfield
Microscopy

Phase
Contrast
Microscopy

Fluorescence
Microscopy

Differential
Interference
Contrast
Microscopy

Polarizing
Microscopy

Video
Microscopy

Confocal
Microscopy

Microscope
Photometry


Carl Zeiss, Inc.

Microscope Division
One Zeiss Drive
Thornwood, New York 10594
800-233-2343 • Fax 914-681-7446



SEE
WHAT
YOU'RE
MISSING

Circle No. 25 on Readers' Service Card



**The system:
Superior performance
fluorescence microscopy**

**The objective:
Nikon CFN
Plan Apo 60X
1.4 N.A.**

Exceptional throughput, resolution and sensitivity clearly demonstrate Nikon's superiority in fluorescence microscopy. Investigate over an extended spectrum with diffraction limited microscopy, employing very low fluorochrome doses.

High UV transmission and large numerical apertures are the rule throughout the wide range of Nikon CF® fluors and CFN plan apochromats. The 60X CFN plan apo, for instance, provides the brightest, highest resolution images possible over the broadest spectrum.


Call today to schedule a demonstration:

1-800-52-NIKON, ext. 133.

See us at the Biophysical Society Meeting,
Booth #'s 805 & 807.


Nikon

AT THE HEART OF
THE SYSTEM YOU NEED



Academic Press

Visit Booth 207 at the 1993 Biophysical Society Meeting



Guide to Electroporation and Electrofusion

edited by
Donald C. Chang, Bruce M. Chassy,
James A. Saunders, and Arthur E. Sowers

"Excellent instructions...Anyone who wants to use electroporation or electrofusion in anything but a routine fashion will find this book very useful."

—ANALYTICAL BIOCHEMISTRY

Comb Bound: \$45.00/ISBN: 0-12-168041-X
Casebound: \$89.00/ISBN: 0-12-168040-1
1992, 581 pp.

Physical Forces and the Mammalian Cell

edited by
John A. Frangos

This book examines the physical forces—fluid shear, stretch, and gravity—that play a role in the physiology of tissues and cellular functions.

February 1993, 400 pp., \$105.00/ISBN: 0-12-265330-0

Encyclopedia of Microbiology

Editor-in-Chief
Joshua Lederberg



Hailed by leading scientists around the world as "excellent" and "outstanding," the *Encyclopedia of Microbiology* will be a standard text for many years. The clearly written, four-volume set describes current advances in environmental microbiology, genetic and molecular microbiology, immunology, and many other areas in rich detail.

Four-Volume Set
October 1992, 2,518 pp., \$695.00/ISBN: 0-12-226890-3

Optical Microscopy Emerging Methods and Applications

edited by
Brian Herman and John J. Lemasters

Optical Microscopy: Emerging Methods and Applications covers recent technical advances and new approaches to monitoring and altering cell physiology, observing multiple cellular activities, examining membrane



cytoarchitecture, confocal imaging of live cell function, observation of intact organ physiology, and lifetime imaging, plus automated clinical imaging cytometry.

December 1992, 456 pp. \$79.95/ISBN: 0-12-342060-1

Bioenergetics 2

David G. Nicholls and Stuart J. Ferguson

CONTENTS: Chemiosmotic Energy Transduction Ion Transport Across Energy-Conserving Membranes. Quantitative Bioenergetics: The Measurement of Driving Forces. The Chemiosmotic Proton Circuit. Respiratory Chains. Photosynthetic Generators of Protonmotive Force. The ATP Synthase. Secondary Transport.

Paperback: \$34.50
August 1992, 272 pp./ISBN: 0-12-518124-8

Nonisotopic DNA Probe Techniques

edited by
Larry J. Kricka

Nonisotopic DNA Probe Techniques presents recently developed alternatives to traditional detection methods that are inexpensive, safe, and ecologically responsible. This book contains chapters on colorimetric, bioluminescent, chemiluminescent, fluorescent, and time-resolved fluorescent detection methods.

Comb Bound: \$39.95/ISBN: 0-12-426296-1
Casebound: \$69.95/ISBN: 0-12-426295-3
1992, 358 pp.


Capillary Electrophoresis

Theory and Practice

edited by
Paul D. Grossman and Joel C. Colburn

This well-illustrated book is designed to help you make the most of your CE separations. It includes comprehensive information on: electroosmosis, separation efficiency; Joule heating, detection methods; in-depth discussion of the separation principles and capabilities of the major modes of CE; and detailed treatment of the application of CE to a wide range of molecules.

1992, 352 pp., \$69.95/ISBN: 0-12-304250-X



Immobilized Affinity Ligand Techniques

Greg T. Hermanson, A. Krishna Mallia, and Paul K. Smith

Comb Bound: \$39.95
September 1992, 454 pp./ISBN: 0-12-342330-9

Volumes in the METHODS IN ENZYMOLOGY Series

Volume 207

Ion Channels

edited by
Bernardo Rudy and Linda E. Iverson

This book describes modern electrophysiological techniques for recording and analysis of ion channels in various systems including cells in culture, *Xenopus* oocytes, brain slices, cellular organelles, and microorganisms.

1992, 917 pp., \$99.00/ISBN: 0-12-182108-0

Volumes 202 and 203

Molecular Design and Modeling: Concepts and Applications

edited by
John J. Langone

Part A Proteins, Peptides, and Enzymes

1991, 824 pp., \$95.00/ISBN: 0-12-182103-X

Part B Antibodies and Antigens, Nucleic Acids, Polysaccharides, and Drugs

1991, 764 pp., \$95.00/ISBN: 0-12-182104-8

Volume 216

Recombinant DNA, Part G

edited by
Ray Wu

CONTENTS (Section Headings): Isolation, Synthesis, and Detection of DNA and RNA. Enzymes and Methods for Cleaving and Manipulating DNA. Reporter Genes. Vectors for Cloning Genes. Author Index. Subject Index.

November 1992, 736 pp., \$90.00/ISBN: 0-12-182117-X

Volume 210

Numerical Computer Methods

edited by
Ludwig Brand and Michael L. Johnson

1992, 718 pp., \$95.00/ISBN: 0-12-182111-0

Order from your local bookseller or directly from

ACADEMIC PRESS

HBJ Order Fulfillment Dept. DM17915
6277 Sea Harbor Drive, Orlando, FL 32887

CALL TOLL FREE

1-800-321-5068

FAX 1-800-336-7377

Prices subject to change without notice. ©1993 Academic Press, Inc.
All Rights Reserved JE/JW/AB 23023

WORKING WITH CONGRESS

A PRACTICAL
GUIDE FOR
SCIENTISTS AND
ENGINEERS

William G. Wells, Jr.

This new guide will help scientists and engineers gain an understanding of how Congress operates and how to communicate effectively with its members and staff.

William G. Wells, Jr., formerly a senior official in Congress and the Executive Office of the President, is presently a Professor of Management Science at George Washington University.

Dr. Wells offers a practical and authoritative orientation to the workings of Congress. He includes advice and step-by-step suggestions on everything from writing your Representative to how to present formal testimony.

A concise, easy-to-read manual that is must reading—especially during this transition period!

AVAILABLE NOW FOR
ONLY \$12.95
(AAAS MEMBERS, \$10.36)


Mail order to AAAS Books,
Box 753, Department A72,
Waldorf, MD 20604

Add \$4 shipping per order. If you prefer, order by phone (VISA/MasterCard only) 301-645-5643 (9am-4pm ET) and ask for AAAS or Fax 301-843-0159.

Circle No. 38 on Readers' Service Card



*The
choice
is yours.*



cDNA SYNTHESIS

AMV
Reverse Transcriptase

The RiboClone® cDNA Synthesis and Cloning Systems simplify the steps involved in synthesizing double-stranded cDNA utilizing RNA from any source. First strand synthesis is driven by AMV Reverse Transcriptase using either random primers or an oligo (dT)-containing primer, followed directly by second strand replacement synthesis using RNase H, and DNA polymerase I. Features:

- Choice of four primers to fit your cloning strategy
 - Xba I Primer-Adaptor (Cat. #C1001) For directional cloning
 - Oligo (dT)₁₅ Primer (Cat. #C2100)
 - Not I Primer-Adaptor (Cat. #C3000) For directional cloning
 - Random Primers (Cat. #C1450) For poly (A-) RNA
- Sufficient reagents to process up to 40µg RNA
- AMV Reverse Transcriptase is the choice for RNAs with high secondary structure.
- A comprehensive *Protocols and Applications Guide* and technical services department are available for your assistance.
- All enzymes available separately for your convenience.



Promega

Promega Corporation

2800 Woods Hollow Road
Madison, WI 53711-5399 USA
Telephone 608-274-4330
Fax 608-273-6967
Telex 62057092
In USA

800-356-9526
Fax 800-356-1970


Copyright 1993, Promega Corporation.
All Rights Reserved.

A complete line of cloning systems and a packaging system are available including Adaptor Ligation Systems, Lambda Vectors and Arms ready for cloning and Packagene® Lambda DNA Packaging System.

Please call: Austl. 02-565-1100, 008-22-5123 (long distance toll free); Fr. 05-48-79-99 (toll free); Neth. 071-217121; Switz. 01-830-70-37; U.K. 0800-378994

Circle No. 33 on Readers' Service Card

SAVE YOUR COPIES OF SCIENCE



These custom-made, imprinted cases and binders are ideal for protecting your valuable *Science* copies from damage. Each binder or case holds one volume of *Science*, or 13 weekly issues – order four binders or cases to hold a complete year of issues. Constructed from reinforced board and covered with durable, leather-like red material and stamped in gold, the cases are V-notched for easy access; binders have a special spring mechanism to hold individual rods which easily snap in.

Quantity	Cases	Binders
One	\$ 7.95	\$ 9.95
Two	\$14.95	\$18.95
Four	\$27.95	\$35.95

SCIENCE

Jesse Jones Industries, Dept. SCE
499 East Erie Avenue
Philadelphia, PA 19134

Enclosed is \$_____ for _____
Cases; _____ Binders. Add \$1 per case/
binder for postage & handling. Outside
USA \$2.50 per case/binder (US funds
only). PA residents add 7% sales tax.

Print

Name _____

Address _____

No P.O. Box Numbers Please

City _____


State/Zip _____

CHARGE ORDERS (Minimum \$15): Am
Ex, Visa, MC, DC accepted. Send card
name, #, Exp. date.

CALL TOLL FREE 7 days, 24 hours
1-800-825-6690

— SATISFACTION GUARANTEED —

*The
choice
is yours.*



cDNA SYNTHESIS

M-MLV

Reverse Transcriptase

Promega introduces four new additions to the cDNA Synthesis and Cloning family using M-MLV Reverse Transcriptase, RNase H Minus for first strand synthesis. Available with choice of primer to suit your cloning strategy, these systems complement the original RiboClone® AMV Reverse Transcriptase Systems.

The systems are designed for the efficient synthesis of double-stranded cDNA using RNA from any source.

All systems are rigorously quality control tested and provided with a complete technical manual to take you effortlessly through the procedure.

M-MLV Reverse Transcriptase
RNase H Minus or
AMV Reverse Transcriptase -
the choice is yours!

■ Four primers available to fit your cloning scheme-

Xba I Primer-Adaptor
(Cat. #C1620) For directional cloning

Oligo (dT)₁₅ Primer
(Cat. #C1630)

Not I Primer-Adaptor
(Cat. #C1660) For directional cloning

Random Primers (Cat. #C1640)
For poly A(-) RNA

■ M-MLV Reverse Transcriptase RNase H Minus is the choice for longer length cDNA's

■ Sufficient reagents to process up to 40 µg RNA

■ All enzymes available separately for your convenience



Promega

Promega Corporation

2800 Woods Hollow Road
Madison, WI 53711-5399 USA
Telephone 608-274-4330
Fax 608-273-6967
Telex 62057092
In USA

800-356-9526
Fax 800-356-1970


Copyright 1993, Promega Corporation.
All Rights Reserved.

A complete line of cloning systems and a packaging system are available including Adaptor Ligation Systems, Linker Ligation Systems, Lambda Vectors and Arms ready for cloning and Packagene® Lambda DNA Packaging System.

Please call: Austl. 02-565-1100, 008-22-5123 (long distance toll free); Fr. 05-48-79-99 (toll free); Neth. 071-217121; Switz. 01-830-70-37; U.K. 0800-378994

Circle No. 34 on Readers' Service Card

NOW PUBLISHED BY
APPLETON & LANGE



Principles of Neural Science

Third Edition

Eric R. Kandel, MD
James H. Schwartz, MD, PhD
Thomas M. Jessell, PhD

*Center for Neurobiology and Behavior,
College of Physicians & Surgeons of
Columbia University; The Howard Hughes
Medical Institute, New York, New York*

with 26 renowned neuroscience contributors

The preeminent neuroscience textbook, the third edition of **Principles of Neural Science**, continues its tradition as the standard textbook in the field. This state-of-the-science resource offers a definitive, comprehensive view of neuroscience, featuring clear writing, instructive illustrations and boxes highlighting key issues for easy referral. It's the perfect reference for every practitioner and researcher.

The extraordinary advances that are revolutionizing neuroscience are featured, including

- Newly identified mechanisms important in the development of the nervous system including the generation of cell lineages, target recognition, and synapse formation
- Insight into the genetic basis of diseases such as muscular dystrophy, retinoblastoma, and neurofibromatosis

1991, 1,135 pp., 850 illus., case, ISBN 0-8385-8034-3, A8034-9, \$69.50

To order this or any Appleton & Lange book, visit your local bookstore or call toll-free
1-800-423-1359 (in CT 838-4400).



APPLETON & LANGE
P.O. Box 5630
25 Van Zant Street
Norwalk, CT 06856
Simon & Schuster Business and Professional Group


1/93

Circle No. 17 on Readers' Service Card

Optimize Your PCR* And Get The Job Done Right.

The PCR* Optimizer™ from Invitrogen lets you determine optimal conditions for your PCR* reaction. This results in higher success rates and increased reliability. With PCR* Optimizer™, you can get the job done right the first time.

The success of your PCR* reactions cannot be predetermined by primer/template specificity alone. Other conditions such as Mg⁺⁺ concentration, pH and the presence or absence of DMSO also effect the successful outcome of PCR*.



Before Optimization After Optimization

By identifying the optimal Mg⁺⁺ concentrations, pH and DMSO levels you can;

- Optimize any Primer/Template Combination.
- Generate the maximum yield and purity of PCR* products.

Each PCR* Optimizer™ kit contains sufficient reagents for evaluating 100 PCR* reactions as well as detailed instructions and buffer concentrations, so you can get the job done right the first time.

Make certain your next PCR* reaction is a success, call...

1-800-955-6288 EXT-728

and ask for the Optimizer™, Cat. No K1220-01.



3985 • B Sorrento Valley Blvd.,
San Diego, CA 92121

Phone: (619) 597-6200 • FAX: (619) 597-6201

BRITISH BIO-TECHNOLOGY
PRODUCTS LTD

UK - Tel: (0865) 781045 Fax: (0235) 533420

France - Numéro Vert 05 90 72 49

Sweden - 020-Linjen 020 793149

Norway - Ring Grønt Nummer 050 11033

Denmark - Grønt Nummer 80 01 85 92

Belgium - Numéro Vert/Groen Nummer 78 11 04 68

Italy

Tel: 39-238103171

Fax: 39-2381014651



Germany - Tel: 62-21303907

Fax: 62-21303511

Switzerland - Tel: 155-5044



Japan

Tel: 81-356841622

Fax: 81-356841633




FUKAKOSHI CO., LTD.

BDH INC., CANADA - TEL: 800-565-7933 • BIO-TRADE, AUSTRIA - TEL: 43-1-8891819 • FINNZYMES, FINLAND - TEL: 35804208077 • MEDOS, AUSTRALIA - TEL: 61-38089077 • SANBIO BV, NETHERLANDS - TEL: 31-413251115 • TAL RON, ISRAEL - TEL: 972-8-472563 • TDI, SPAIN - TEL: 34-14091251 • ECOGEN, SPAIN - TEL: 34-934560607 • UNITED RESEARCH/GOODMAN BIOTECHNOLOGIES, INDIA - TEL: 59 1107

* PCR is covered by U.S. Patent No. 4,683,202 issued to Cetus Corporation.

Circle No. 10 on Readers' Service Card

SlideWrite™ Plus makes your data look terrific.



Try a free, full-featured trial copy of SlideWrite Plus. You can even use your own data! No risk or obligation!

Scientists and technical users! Now you can try the only presentation graphics program designed specially for *you*! SlideWrite Plus, for DOS or Windows™, gives you technical graphing, word charting, science figure packs, scalable fonts, freehand drawing, and more! Supports all popular printers, plotters, and film recorders.

New DOS Version Includes: Advanced curve fitting. 100 linear and 16 parametric nonlinear functions. **Statistics and data analysis.** t-Test module, histogram, and function evaluator. **Table editor.** Makes it easy to draw tables of text and numbers.

Call toll free **1-800-795-4754** or FAX 619-931-9313

Include name, company, address, and business phone. Specify DOS or Windows, disk size.

Advanced Graphics Software, Inc.

5825 Avenida Encinas, Suite 105, Carlsbad, CA 92008

In Europe: BIS Bedrijfs Informatie Systemen, 2980 GB Ridderkerk, The Netherlands.
Ph. +31-1804-30000. Fax: +31-1804-12377

Our free trial disk will prove it!



For DOS circle reader service number 29
For Windows circle reader service number 30

New from AAAS Press!

Two books that are shaking the world
of science...*with laughter!*

**Chalk Up Another One
The Best of Sidney Harris**
by Sidney Harris

This "best-of-the-best" collection includes many of his classics. Harris proves once again that science is very funny, indeed!



Big Science
by Nick Downes

Anything can happen when cartoonist Nick Downes looks to *Big Science* for explanations. Don't miss out on the weird and wonderful world of Nick Downes.

Available in bookstores or by mail for only \$10.95 each
or to AAAS members for \$8.75 each.

Mail to: AAAS Books, P.O. Box 753
Department A69, Waldorf, MD 20604

Add \$4 shipping per order. If you prefer, order by phone (VISA/MasterCard only)
(301) 645-5643 (9am-4pm ET) and ask for AAAS or Fax (301) 843-0159.

American Association for the Advancement of Science

Circle No. 37 on Readers' Service Card

The John Ugelstad Conference II



MIT, Cambridge, MA
19-20 April 1993

Magnetic Separation in Molecular and Cellular Biology

Topics will include:

• Subtraction & Hybridization Strategies • Genome Sequencing • Diagnostic Sequencing • Forensic Analysis • Microbiology • HIV Infection • Subcellular Separation • Gene Therapy • Future Technologies

The Massachusetts Institute of Technology is honored to host the second John Ugelstad Conference on April 19-20, 1993. The focus of this meeting will be magnetic separation techniques in Molecular and Cellular Biology. Distinguished speakers from around the world will present results on the latest magnetic separation techniques applied to virology, gene therapy, bacteriology, clinical immunology and the human genome project. This two day meeting will promote the continued advancement of magnetic separation technology in clinical and basic research.

Abstract Presentations:

Researchers are encouraged to submit poster abstracts. Deadline for submission of abstracts is March 1, 1993.

Special Events:

• The John Ugelstad Award • John Ugelstad Review Lecture • Conference Banquet

Conference Speakers:

Anthony J. Sinskey
MIT

Robert B. Weiss
University of Utah

Annemarie Poustka
German Cancer Res. Ctr.

Richard A. Gibbs
Baylor College of Medicine

Katsumi Isono
Kobe University, Japan

Lloyd M. Smith
Univ. of Wisconsin-Madison

Mathias Uhlen
Royal Institute, Sweden

S. Mizutani
*National Children's Med.
Res. Ctr., Japan*

Kathryn Howell
University of Colorado

Charles R. Cantor
Boston University

Jean-Claude Chermann
*National Institute of Med. Res.,
France*

James J. Murtagh
Emory University

Orjan Olsvik
CDC, Atlanta

Edward F. Rossomondo
University of Connecticut

Bruce Trapnell
Genetic Therapy, Inc.

Kevin M. Sullivan
Forensic Science Service, U.K.

Ignatio Rodriguez
National Institutes of Health

The John Ugelstad Conference II Registration Form

19-20 April 1993, MIT, Cambridge, MA

Registration Fees: \$200 U.S. before March 15, 1993 (\$300 U.S. Late Registration Fee)

To register, mail this form with payment to: The Secretariat of The John Ugelstad Conference II, MIT Conference Services Office, Room 7-111, Massachusetts Institute of Technology, Cambridge, MA 02139 Tel: 617-253-1700 Fax: 617-253-7002

Name: _____ Title: _____

Organization: _____

Address: _____

City/State/Zip: _____ Country: _____

Telephone: _____ ☐ Check enclosed (Please make payable to The Secretariat of the John Ugelstad Conference II)

☐ Please charge my credit card

☐ VISA ☐ MasterCard Account #: _____

Signature: _____ Expiration Date: _____

Cancellation Policy: All cancellations must be received in writing. A refund, minus a \$25 processing fee, will be made for cancellations received prior to April 1, 1993. No refund will be made for cancellations received after this date. Substitutions from the same institution can be made at any time.

Sponsor's Liability: In the event of a cancellation of this conference, the sponsor's liability is limited to the return of the registration fee.

