


AMERICAN  
ASSOCIATION FOR THE  
ADVANCEMENT OF  
SCIENCE

# SCIENCE

26 MARCH 1993  
VOL. 259 • PAGES 1797-1992

\$6.00







106 Charles Lindbergh Blvd. Uniondale, NY 11553

# The SOURCE™ for PhosphoSHure™ SH2/SH3 Discovery Tools

Biologically engineered to detect and isolate activated phosphotyrosyl-proteins

## New! SH2/SH3 Products

For the First Time Ever!

Obtain specific information on the SH2/SH3 binding characteristics of proteins involved in signal transduction pathways.

Available as agarose and biotinylated conjugates and in two sampler kits!

In addition to these 6 new products, we are The Source™ for antibodies to many signal transduction proteins, including *src*, *trk*, and *abl*, as well as antibodies to phosphotyrosine, phosphatases and growth factor receptors including EGF-R, PDGF-R, and *erb* B2(*neu*).

• Identify • Capture • Purify

### ORDERING INFORMATION

PRODUCT	CAT. #	FORM	WORKS FOR
<b>Agarose Coupled Reagents:</b>			
PhosphoBind <i>src</i> SH2	SH01A	1mg/0.5 ml	IP, Purification
PhosphoBind PLCγ-N SH2	SH02A	1 mg/0.5 ml	IP, Purification
PhosphoBind p85-C SH2	SH03A	1 mg/0.5 ml	IP, Purification
PhosphoBind Sampler Kit	SHK01	0.1 ml x 3 Agarose Reagents	
<b>Biotin Coupled Reagents:</b>			
PhosphoBlot <i>src</i> SH2	SH01B	100 µg	Blotting
PhosphoBlot PLCγ-N SH2	SH02B	100 µg	Blotting
PhosphoBlot p85-C SH2	SH03B	100 µg	Blotting
PhosphoBlot Sampler Kit	SHK02	20 µg x 3 Biotin Reagents	

Biotinylated products are supplied in 0.5 ml of 0.05 M sodium phosphate buffer, pH 7.4, 0.2% gelatin, 0.1% azide. Agarose products are supplied as 0.5 ml of packed beads in sodium phosphate buffer, pH 7.4, 0.1% azide, 10% glycerol, 5 mM DTT

## Oncogene Science:

The SOURCE™ for over 1,000 Antibodies, Probes, ELISA's, and the Reagents to Make Them Work!

**For technical information CALL 800-662-2616**

Circle 16 for a Complete Catalog

Circle 99 for SH2 Product Information

Circle No. 100 on Readers' Service Card



# Extremely Sensitive Immunoassays for the Quantitation of Cytokines.

## QUANTIKINE<sup>TM</sup>



### Quantikine Kit

*The Tool To Measure Cytokines:*

IL-1 $\alpha$	TNF $\alpha$
IL-1 $\beta$	TNF $\beta$
IL-2	G-CSF
IL-3	GM-CSF
IL-4	New: LIF
IL-6	IL-1ra
IL-7	Erythropoietin*
IL-8	Transferrin Receptor

The R&D Systems' Quantikine series of immunoassays is a complete set of reagents designed to accurately measure the concentration of cytokines in fluids. The high sensitivity and wide standard curve range offered by each Quantikine Kit enables precise measurement of cytokines at both high and low concentrations in less than 4½ hours.

**FOR RESEARCH USE ONLY. NOT FOR USE IN  
DIAGNOSTIC PROCEDURES.**


\*CLINIGEN ERYTHROPOIETIN KIT – FOR IN VITRO DIAGNOSTIC USE.

In Europe contact:

**British Biotechnology, Ltd.**  
4-10 The Quadrant, Barton Lane  
Abingdon, Oxon OX14 3YS  
Telephone: +44 (0865) 781045  
Fax: +44 (0235) 533420

In Japan contact:

**Funakoshi Co., Ltd.**  
9-7, Hongo 2-Chome  
Bunkyo-ku, Tokyo 113  
Telephone: +81 (03) 56841622  
Fax: +81 (03) 56841633



# R&D SYSTEMS


## 1-800-343-7475

Circle No. 45 on Readers' Service Card

**R&D Systems**  
614 McKinley Place N.E.  
Minneapolis, MN 55413  
In Minnesota: (612) 379-2956  
Fax: (612) 379-6580



# DNA THERMAL CYCLER 480 SYSTEM





# Advancing Gene Expression Research Through RNA PCR.

Innovations in PCR amplification are leading the way with a powerful tool for gene expression research. Today the DNA

Thermal Cyclers 480 System makes rapid amplification and quantitation of mRNA routine.

By integrating state-of-the-art instrumentation, consumables and RNA PCR reagents, the DNA Thermal Cyclers 480 System ensures

efficient amplification of all types of mRNA sequences, including complex secondary structures. The exceptional temperature uniformity of the DNA Thermal Cyclers 480, coupled with rapid heat transfer facilitated by Thin-Walled GeneAmp™ Reaction Tubes, provides the accuracy and precision required for advanced PCR quantitation.

Above all, the DNA Thermal Cyclers 480 System for RNA PCR comes with our PCR Performance Guarantee. A commitment backed by the comprehensive experience of Perkin-Elmer's worldwide technical support organization.

**Free technical note and bibliography.** To receive free literature on the DNA Thermal Cyclers 480 and RNA PCR, contact your local Perkin-Elmer sales representative or call 1-800-762-4000. To order, call 1-800-762-4002. Outside the U.S., contact your Perkin-Elmer sales representative.

## PERKIN ELMER


Europe Vaterstetten, Germany Tel: 49-8106-381-115 Fax: 49-8106-6697  
Canada Montreal, Canada Tel: 514-737-7575 Fax: 514-737-9726  
Far East Melbourne, Australia Tel: 61-3-560-4566 Fax: 61-3-560-3231  
Latin America Mexico City, Mexico Tel: 52-5-651-7077 Fax: 52-5-593-6223

Perkin-Elmer PCR reagents are developed and manufactured by Roche Molecular Systems, Inc., Branchburg, New Jersey, U.S.A.




Perkin-Elmer is a registered trademark of The Perkin-Elmer Corporation. GeneAmp is a registered trademark/trademark of Hoffmann-La Roche Inc. and F. Hoffmann-La Roche A.G. The GeneAmp PCR process is covered by U.S. patents owned by Hoffmann-La Roche Inc. and F. Hoffmann-La Roche A.G.

Circle No. 47 on Readers' Service Card



DETECTION OF MARINE SPECIES METALLOTHIONEIN MRNAs (114 bp) WITH INTERNAL STANDARDS (122 bp) USING THE GENEAMP® RNA PCR KIT AND THE DNA THERMAL CYCLER 480 WITH THIN-WALLED GENEAMP REACTION TUBES. COURTESY OF DR. JOSEPH F. CRIVELLO, UNIVERSITY OF CONNECTICUT.





1816

Washington focuses  
on technology

## NEWS & COMMENT

**R&D Policy That Emphasizes the 'D'** 1816  
Technology Boosting: A Checkered History  
NIST: Measuring Up to a New Task

**NIH Adds an Extra Layer of Review for Sensitive Grants** 1820  
Eyeing a Project's Ethics

**AIDS Vaccines: MicroGeneSys Withdraws From Trial** 1821

**Engineering Academy Elects New Members** 1822

## RESEARCH NEWS

**Plugging a Cosmic Information Link** 1824

**Aging Twins Offer Clues to Late-Onset Diseases** 1826  
What the Twin Studies Might Tell Us

**Enzyme May Blunt Cocaine's Action** 1828

**The Search for Liver Stem Cells Picks Up** 1829

## POLICY FORUM

**Integrated Assessment of Climate Change** 1813  
H. Dowlatabadi and M. G. Morgan

## PERSPECTIVE

**Protein Prenylation: A Mediator of Protein-Protein Interactions** 1865  
C. J. Marshall

## ARTICLES

**Dynamics of Elliptical Galaxies** 1867  
D. Merritt

**What Do We Learn from Neutrinos?** 1872  
J. Steinberger

## RESEARCH ARTICLE

**Regional Ground-Water Mixing and the Origin of Saline Fluids: Midcontinent, United States** 1877  
M. Musgrove and J. L. Banner

## DEPARTMENTS

**THIS WEEK IN SCIENCE** 1805

**EDITORIAL** 1807  
The Great Overcoat Scare

**LETTERS** 1809  
AIDS Virus History: J. C. Gluckman • Japanese Support of U.S. Research: T. S. B. Yen • Thoroughly Modern Reptiles: R. T. Moore • All in the Family: P. Mitchell • Kuwait Oil Fires: Correction: P. V. Hobbs and L. F. Radke

**SCIENCESCOPE** 1815

**RANDOM SAMPLES** 1823

Administration Drops Indirect Cost Cuts • Europe First With New X-Ray Source • A New Look at Racehorse Genetics

**INSIDE AAAS** 1922

**BOOK REVIEWS** 1923  
*Heisenberg's War*, reviewed by J. Bernstein • *The Major Gods of Ancient Yucatan*, J. Kowalski • *Linguistic Diversity in Space and Time*, W. H. Baxter • *Pi in the Sky*, J. Borwein and P. Borwein • Vignettes: Obscure Vices • Books Received

## AAAS Board of Directors

F. Sherwood Rowland  
*Retiring President,  
Chairman*  
Eloise E. Clark  
*President*  
Francisco J. Ayala  
*President-elect*

Robert A. Frosch  
Florence P. Haseltine  
William A. Lester, Jr.

Alan Schriesheim  
Jean'ne M. Shreeve  
Chang-Lin Tien  
Warren M. Washington  
Nancy S. Wexler

William T. Golden  
*Treasurer*  
Richard S. Nicholson  
*Executive Officer*

John Abelson  
Frederick W. Alt  
Don L. Anderson  
Stephen J. Benkovic  
David E. Bloom  
Floyd E. Bloom  
Henry R. Bourne  
James J. Bull  
Kathryn Calame  
C. Thomas Caskey  
Dennis W. Choi

John M. Coffin  
Bruce F. Eldridge  
Paul T. Englund  
Richard G. Fairbanks  
Douglas T. Fearon  
Harry A. Fozzard  
Victor R. Fuchs  
Theodore H. Geballe  
Margaret J. Geller  
John C. Gerhart  
Roger I. M. Glass

Stephen P. Goff  
Corey S. Goodman  
Stephen J. Gould  
Ira Herskowitz  
Eric F. Johnson  
Stephen M. Kosslyn  
Michael LaBarbera  
Charles S. Levings III  
Harvey F. Lodish  
Richard Losick  
Anthony R. Means

Mortimer Mishkin  
Roger A. Nicoll  
William H. Orme-Johnson III  
Stuart L. Pimm  
Yeshayau Pocker  
Dennis A. Powers  
Ralph S. Quatrano  
V. Ramanathan  
Douglas C. Rees  
Erkki Ruoslahti  
Ronald H. Schwartz

Terrence J. Sejnowski  
Thomas A. Steitz  
Richard F. Thompson  
Robert T. N. Tjian  
Emil R. Unanue  
Geerat J. Vermeij  
Bert Vogelstein  
Harold Weintraub  
Zena Werb  
George M. Whitesides  
Owen N. Witte  
Keith Yamamoto


## Board of Reviewing Editors



## COVER

Transparent surface models of the helical nucleoprotein filaments formed by the yeast RAD51 protein (top) and the bacterial RecA protein (bottom). Within RAD51 filaments, DNA (red) is extended and untwisted in a manner similar to its extension and untwist-

ing in RecA filaments. The homology of these filaments suggests that the structures for recombination in eukaryotes are similar to those in prokaryotes. See page 1896. [Computer graphic: Edward H. Egelman]



## REPORTS

### The Mechanical Response of Gold Substrates Passivated by Self-Assembling Monolayer Films

1883

R. C. Thomas, J. E. Houston, T. A. Michalske, R. M. Crooks

### C<sub>60</sub>H<sub>2</sub>: Synthesis of the Simplest C<sub>60</sub> Hydrocarbon Derivative

1885

C. C. Henderson and P. A. Cahill

### Fabrication and Properties of Free-Standing C<sub>60</sub> Membranes

1887

C. B. Eom, A. F. Hebard, L. E. Trimble, G. K. Celler, R. C. Haddon

### Outgassed Water on Mars: Constraints from Melt Inclusions in SNC Meteorites

1890

H. Y. McSween, Jr., and R. P. Harvey

### Structural Relationship of Bacterial RecA Proteins to Recombination Proteins from Bacteriophage T4 and Yeast

1892

R. M. Story, D. K. Bishop, N. Kleckner, T. A. Steitz

### Similarity of the Yeast RAD51 Filament to the Bacterial RecA Filament

1896

T. Ogawa, X. Yu, A. Shinohara, E. H. Egelman

### Antibody-Catalyzed Degradation of Cocaine

1899

D. W. Landry, K. Zhao, G. X.-Q. Yang, M. Glickman, T. M. Georgiadis

### Translocation of TCR $\alpha$ Chains into the Lumen of the Endoplasmic Reticulum and Their Degradation

1901

J. Shin, S. Lee, J. L. Strominger

### Germ Line Transmission of a Yeast Artificial Chromosome Spanning the Murine $\alpha_1$ (I) Collagen Locus

1904

W. M. Strauss, J. Dausman, C. Beard, C. Johnson, J. B. Lawrence, R. Jaenisch

### Cyclin-Dependent Regulation of G<sub>1</sub> in Mammalian Fibroblasts

1908

M. Ohtsubo and J. M. Roberts

### NF- $\kappa$ B Controls Expression of Inhibitor I $\kappa$ B $\alpha$ : Evidence for an Inducible Autoregulatory Pathway

1912

S.-C. Sun, P. A. Ganchi, D. W. Ballard, W. C. Greene

### Activity-Dependent Regulation of Conductances in Model Neurons


1915

G. LeMasson, E. Marder, L. F. Abbott

### An Essential Heparin-Binding Domain in the Fibroblast Growth Factor Receptor Kinase

1918

M. Kan, F. Wang, J. Xu, J. W. Crabb, J. Hou, W. L. McKeehan



1867  
Star dance

1904  
Expression of YAC DNA transmitted through the germ line to generate transgenic mice


■ Indicates accompanying feature

■ **SCIENCE** (ISSN 0036-8075) is published weekly on Friday, except the last week in December, by the American Association for the Advancement of Science, 1333 H Street, NW, Washington, DC 20005. Second-class postage (publication No. 484460) paid at Washington, DC, and additional mailing offices. Copyright © 1993 by the American Association for the Advancement of Science. The title **SCIENCE** is a registered trademark of the AAAS. Domestic individual membership and subscription (51 issues): \$87 (\$47 allocated to subscription). Domestic institutional subscription (51 issues): \$205. Foreign postage extra: Mexico, Caribbean (surface mail) \$50; other countries (air assist delivery) \$95. First class, airmail, student and emeritus rates on request. Canadian rates with GST available upon request, GST #1254 88122. **Change of address:** allow 6 weeks, giving old and new addresses and 11-digit account number. **Postmaster:** Send change of address to *Science*, P.O. Box 2033, Marion, OH 43305-2033. **Single copy sales:** \$6.00 per issue prepaid includes surface postage; Guide to Biotechnology Products and Instruments, \$20.

Bulk rates on request. **Authorization to photocopy** material for internal or personal use under circumstances not falling within the fair use provisions of the Copyright Act is granted by AAAS to libraries and other users registered with the Copyright Clearance Center (CCC) Transactional Reporting Service, provided that the base fee of \$1 per copy plus \$0.10 per page is paid directly to CCC, 27 Congress Street, Salem, MA 01970. The identification code for *Science* is 0036-8075/93 \$1 + .10. *Science* is indexed in the *Reader's Guide to Periodical Literature* and in several specialized indexes.

■ The American Association for the Advancement of Science was founded in 1848 and incorporated in 1874. Its objectives are to further the work of scientists, to facilitate cooperation among them, to foster scientific freedom and responsibility, to improve the effectiveness of science in the promotion of human welfare, to advance education in science, and to increase public understanding and appreciation of the importance and promise of the methods of science in human progress.





# Ultrapure Plasmid DNA for Transfections


## QIAGEN Plasmid DNA Purification means:

- multiple plasmid preps in 2 hours.
- plasmid yields up to 2.5 mg.
- no CsCl/EtBr gradient ultracentrifugation.
- no phenol/chloroform extractions.

With gravity flow QIAGEN-tips, large quantities of ultrapure plasmid DNA can be obtained in just 2 hours. DNA is selectively purified on QIAGEN-tips containing our unique anion exchange resin. RNA, protein and polysaccharides are efficiently removed to yield pure supercoiled, plasmid DNA.

DNA purified with QIAGEN-tips performs better than DNA prepared by CsCl density gradient centrifugation. The DNA is ideally suited for demanding applications such as transfection of eukaryotic cell lines (Figure 1) and automated sequencing.

QIAGEN Plasmid Kits come ready-to-use with gravity flow QIAGEN-tips, RNase A, all reagents and buffers and a detailed handbook.



**QIAGEN purified DNA for transfection:** NIH 3T3 cells were transfected with 1, 2 and 3 µg of the plasmid pRSVlacZ according to the standard calcium phosphate method. Plasmid DNAs used for transfection were prepared with QIAGEN Plasmid Mega Kit (blue bar) or by CsCl gradient centrifugation (red bar). The transfection efficiencies were determined by measurement of β-galactosidase activity. Samples 1 and 2: 1 µg DNA; samples 3 and 4: 2 µg DNA; samples 5 and 6: 3 µg DNA. Data kindly provided by Frauke Ehlert, Institute for Molecular Biology and Tumor Research, Marburg, Germany.

## QIAGEN Plasmid Kits

QIAGEN Plasmid Midi Kit	Cat. No. 12143
QIAGEN Plasmid Maxi Kit	Cat. No. 12162
QIAGEN Plasmid Mega Kit	Cat. No. 12181

**DIAGEN GmbH** Max-Volmer-Straße 4, 4010 Hilden, Germany, Phone (0)2103-892-0, Orders (0)2103-892-230, Fax (0)2103-892-222, Technical Service (0)2103-892-240

**QIAGEN Inc.** 9259 Eton Avenue, Chatsworth, CA 91311, USA, Orders 800-426-8157, Fax 818-718-2056, Technical Service 800-DNA-PREP (800-362-7737)

**Distributors:** AUSTRALIA: PHOENIX Scientific Industries Ltd. (03) 544 8022, AUSTRIA: BIO-TRADE (222) 889 18 19, BENELUX: Westburg B. V. (33) 95 00 94, DENMARK: KEBO Lab A/S (44) 88 72 00, FINLAND: KEBO OY (90) 435 420 03, FRANCE: Coger (1) 45 32 35 17, GREECE: Bio Analytica (01) 644 0547, HONG KONG/CHINA: Diagnostics Co., Ltd. (892) 542 0566, ISRAEL: BIO-LAB Ltd. (02) 524 447, ITALY: Genenco (M-Medical srl) (055) 5001871, JAPAN: Funakoshi Co., Ltd. (3) 5684-1622, KOREA: LRS Laboratories, Inc. 924 8697, MALAYSIA: Filarite Company (03) 755 8582, NORWAY: KEBO Lab AS (02) 30 17 18, PORTUGAL: Izasa Portugal (01) 757 0740, RSA: Whitehead Scientific Supplies (021) 981-1560, SINGAPORE: TIBS TRADING PTE Ltd. 292 9783, SPAIN: Izasa S.A. (93) 401 01 01, SWEDEN: KEBO Lab AB (08) 621 34 00, SWITZERLAND: KONTRON Instruments AG (01) 733-5-733, TAIWAN: Formo Industrial Co., Ltd. (02) 736 7125, UK: Hybaid Ltd. (081) 977 3266



Circle No. 3 on Readers' Service Card



## Hot galaxies

Elliptical galaxies are elliptical not because they are rotating, but because the stars in them have distinctly anisotropic velocity distributions. As stellar systems, they are "hot," their shapes maintained against their own gravity by random motions characterized by three different velocity dispersions along their principal axes. Merritt (p. 1867) explains how *N*-body simulations have been used to calculate stellar orbits in a triaxial gravitational potential and show how those orbits can be populated so as to produce a mass distribution that leads to the presumed potential. Observational studies allied with these models have allowed astronomers to test for the presence of dark matter and massive black holes in elliptical galaxies.

## Free-standing buckyfilms

Solid fullerenes are held together by van der Waals forces. Eom *et al.* (p. 1887) have found that these forces, usually considered among the weaker interactions, are strong enough to create free-standing films of fullerene molecules. The films were fabricated by micromachining a silicon wafer after thermal deposition of a C<sub>60</sub> layer. The layer was uniformly etched down to 2000 angstroms in a CF<sub>4</sub> plasma. Measurements of the mechanical properties of the films suggest possible applications, including sensors and filters.

## Tracing water in Mars

The amount of water originally contained in Mars and its degassing history are uncertain. McSween and Harvey (p. 1890) searched for clues of the amount

## Tracing ground water in the crust

Three ground-water systems of different ages and origins mix on a regional scale in the midcontinental United States. Musgrove and Banner (p. 1877) use chemical analyses to trace the origins of the waters and to characterize the mixing. Interestingly, two of the fluids are saline but acquired their salinity in distinct ways. Ground water originating in the Colorado Front Range and flowing to eastern Kansas becomes saline by mixing with brines leaching from halite-bearing rocks above the primary aquifers. Near northern Oklahoma, these waters mix with brines that may have originated from trapped Paleozoic seawater. Finally, these waters mix with dilute meteoric waters flowing west from the Ozark Dome in southern Missouri.

of water in Mars' interior in the only direct samples of Martian rocks available—the so-called SNC meteorites, which are believed to be igneous rocks crystallized from a Martian magma and then ejected into space by a large asteroid impact. Analysis of melt inclusions frozen in crystals suggest that the water content of progenitor magmas was about 1.4 percent by weight. Degassing of this amount from all of the igneous rocks observed on Mars could yield a global ocean about 200 meters deep, which yields a volume that is still somewhat less than that required to account for the large flood structures on Mars.

## RecA relatives

In the presence of adenosine triphosphate, RecA polymerizes on DNA in prokaryotes to form a helical nucleoprotein filament in which double-stranded DNA is stretched and unwound. Within the filament, RecA catalyzes recognition, pairing, and strand-exchange between two homologous DNA molecules. Ogawa *et al.* (p. 1896) show that a yeast protein, RAD51, which has 30 percent sequence homology with RecA, can polymerize on DNA and form nucleoprotein filaments

like those formed by RecA. The similarity of these types of filaments suggests that RAD51 may function like RecA. Story *et al.* (p. 1892) discuss the structural similarity of recombination proteins from yeast and bacteriophage to RecA.

## Getting large genes into stem cells

Introducing large segments of DNA into the germ line cells of mice can be useful for functional studies of large genes and for identifying genes within the segment through their ability to complement lost genetic functions in mutant strains. Strauss *et al.* (p. 1904) used yeast artificial chromosomes to introduce a purified 150-kilobase DNA segment comprising the gene for  $\alpha_1$ (I) collagen into the germ line of a mouse. Modifications of existing technology resulted in significantly higher transfection efficiencies for transfer of lipid-bound DNA into embryonic stem (ES) cells. These cells were introduced into host-mouse blastocytes to generate founder chimeras for transgenic mouse strains. Fluorescence *in situ* hybridization and ribonuclease protection analysis confirmed that the intact transgene was present at one chromoso-

mal location in each ES clone and was expressed at levels similar to that of the endogenous collagen gene.

## Neuron model

Neurons manage to produce stable electrical signals despite changes in cell shape, replacement of ion channels, and extracellular conditions. LeMasson *et al.* (p. 1915) have developed a mathematical model for neuronal conductances in which the activity of the cell, as represented by the concentration of intracellular calcium, helps control the maximum conductances of ion currents. For example, a perturbation in extracellular potassium concentration, which induces continuous firing of potassium channels, increases intracellular calcium and restores burst firing of the channels by slowly decreasing inward currents and increasing outward currents. Such processes can also allow the neuron to maintain its pattern of electrical activity within a network. The model also allows the neuron to change its state in response to presynaptic stimulation.

## Autoregulated factors

The transcription factor NF- $\kappa$ B is involved in gene activation in the immune, acute phase, and inflammatory responses. NF- $\kappa$ B is maintained in the cytoplasm in an inactive form when coupled to an inhibitory protein, I $\kappa$ B. When T cells are activated by stimuli such as mitogens, NF- $\kappa$ B is released from I $\kappa$ B and is translocated to the nucleus. Sun *et al.* (p. 1912) have identified an autoregulatory pathway in which NF- $\kappa$ B controls the expression of the I $\kappa$ B gene.



**Not I  
Price  
Reduction!**  
Recombinant Not I Now Available.  
60% Lower Cost Per Unit.  
2.5x More Units Per Package.

# Cloned enzymes put the future in your hands.

*Cloned enzymes are the best tools available today for your research tomorrow.*

At New England Biolabs, we are helping you shape the future of molecular biology by offering the largest selection of cloned enzymes available today. You know why you clone: that's why we clone. And with our 12 years of experience in the cloning and overexpression of restriction/modification systems, we have introduced not only substantial cuts in reagent costs but also unsurpassed improvements in both product quality and purity.

GET 2.5x  
MORE **Not I**  
FOR THE  
SAME LOW  
PRICE, A  
60% LOWER  
COST PER UNIT.



And, New England Biolab's cloned enzymes provide consistency you can count on. In fact, we guarantee it. So, no matter what your application: basic or applied research, paternity testing, DNA typing, diagnostics, or even glyco-biology...cloned enzymes from New England Biolabs offer improved consistency, outstanding purity and unmatched affordability.

*Call 1-800-NEB-LABS and  
put the future in your  
laboratory today.*

- New England Biolabs Inc. 32 Tozer Road, Beverly, MA 01915 USA  
1-800-NEB-LABS (US and MA) Tel. (508) 927-5054 Fax (508) 921-1350
- New England Biolabs Ltd., Canada Tel. (800) 387-1095 (416) 672-3370 Fax (416) 672-3414






LIGHT BASED DETECTION OF BIOLOGICAL MOLECULES



## Light Years Ahead in Western Blotting

### Western blotting



Greater sensitivity than any  
other technology -  
results in minutes,  
often seconds!

### The Proof

For Western blotting applications, ECL™ systems are simply the best! - but don't just take our word for it. After just two years, ECL has become the most referenced light-based system, proving it works.

No other system is as sensitive or as fast as ECL.

- Orders of magnitude more sensitive
- Results in minutes, often seconds
- Simple multiple reprobings
- Hard copy results

Let ECL Western blotting systems help you get ahead.

Extensive bibliography -  
the proof that it works



Amersham International plc  
Amersham England HP7 9NA  
Tel: 0494 544000

For further information  
contact your local office.

**Australia**  
Sydney (02) 888-2288  
**Belgium**  
Gent (091) 41-52-80  
**Canada**  
Oakville (416) 847-1166  
**Denmark**  
Birkerød (45) 82-02-22  
**Far East**  
Hong Kong (852) 802-1288

**France**  
Les Ulis (1) 69-1828-00  
**Germany**  
Braunschweig (053) 07 2 06-0  
**Iberica**  
Madrid (91) 304-42-00  
**Italy**  
Milan (2) 5097740  
**Japan**  
Tokyo (03) 38-16-6161

**Nederland**  
Hertogenbosch (073) 41-85-25  
**Norway**  
Gjetsum (02) 54-63-18  
**Sweden**  
Solna (08) 734-0800  
**UK Sales**  
Aylesbury (0296) 395222  
**USA Arlington Heights**  
IL (800) 323-9750

ECL, ECL+device and Amersham are trademarks of Amersham International plc

Circle No. 2 on Readers' Service Card



**NEW!**  
Rat & Mouse Cell  
Adhesion Antibodies

CMP-Sialic  
Acid Synthetase



Endoglycoceramidase



**NEW!**  
Endoglycosidase F<sub>2</sub>



Endoglycosidase H



**COMING SOON!**  
 $\alpha$ 2,3 Sialyltransferase



O-Glycanase®



**NEW!**  
 $\alpha$ 2,6 Sialyltransferase



Neuraminidase



**NEW!**  
BSA-Free  
O-Glycanase®



Recombinant  
N-Glycanase®



Natural  
N-Glycanase®



## Carbohydrate problems? Here's food for thought.

Complex carbohydrates come in all shapes and sizes. That's why Genzyme offers an ever-expanding menu of enzymes to suit every taste.

Now, you can sample these four new varieties:

- BSA-free O-Glycanase®
- Endoglycosidase F<sub>2</sub>
- $\alpha$ 2,3 Sialyltransferase
- $\alpha$ 2,6 Sialyltransferase

*To find out more about our complete line of enzymes for glycoconjugate research — and for your FREE copy of our new Deglycosylation Selection Guide — just call 1-800-332-1042 or contact your local Genzyme distributor.*

**genzyme**

GlycAnalysis™ Systems


**Genzyme Corporation**, One Kendall Square, Cambridge, MA 02139 1-800-332-1042 FAX 617-252-7600

N-Glycanase® and O-Glycanase® are registered trademarks of Genzyme Corporation. ©1992, Genzyme Corporation

Circle No. 62 on Readers' Service Card



PAIN AND DISTRESS ARE ON THEIR WAY OUT.



## TiterMax<sup>®</sup> is in.

Now there is an ANIMAL-SAFE, EFFECTIVE ADJUVANT.

All the performance of the traditional adjuvant with MINIMAL TOXIC EFFECTS.

### TiterMax:

- produces reliably high titer levels
- maintains high titer levels over time
- can be boosted safely when needed

For more TiterMax information complete the coupon or call 800-345-2987.

### Vaxcel<sup>™</sup>, Inc.

150 Technology Parkway  
Technology Park/Atlanta  
Norcross, GA 30092

Vaxcel, Inc. is a subsidiary of  
CytRx<sup>®</sup> Corporation.

Vaxcel<sup>™</sup> and TiterMax<sup>®</sup> are trademarks of  
CytRx Corporation.

**CAUTION:** TiterMax is a new adjuvant for  
investigational use only in laboratory research  
animals or for tests in vitro. Not for use in  
humans.

©1993 CytRx Corporation. All Rights Reserved.  
Printed in USA.

PLEASE SEND MORE TiterMax<sup>®</sup> INFORMATION.

NAME..... TITLE.....

COMPANY/  
INSTITUTION.....

STREET.....

CITY/STATE/ZIP.....

TELEPHONE (.....)..... EXT.....

RETURN TO: **Vaxcel<sup>™</sup>, Inc.**

150 Technology Parkway  
Technology Park/Atlanta  
Norcross, GA 30092

Attn: Technical Services Group

Circle No. 52 on Readers' Service Card



Scientists reluctant to incorporate nonradioactive labeling and detection into their research techniques have historically accused nonradioactive products of lacking sufficient levels of sensitivity. However, recent technological advancements — chemiluminescent detection, to be exact — have allowed select products to meet or exceed the sensitivity levels achievable with  $^{32}\text{P}$  labeling and detection.

### Southern blotting: single-copy gene detection in 30 minutes


The accepted wait for accurate single-copy gene detection using  $^{32}\text{P}$  labeled probes has always been one to three days. Now identical results can be achieved in 30 minutes with the Genius System and its chemiluminescent substrate.

In addition, the probes are stable, allowing researchers to reuse them

employed to prepare radioactive probes, and multiple filters can be processed simultaneously.

### Dot blots: guaranteed results in under one hour


With exposure times reduced from overnight to 30-60 minutes, critical answers are relayed quickly, helping research progress as soon as possible.



Comparison of slot blot hybridizations using Genius-labeled oligonucleotide probe (left panel, 60 minute film exposure) and  $^{32}\text{P}$ -labeled probe (right panel, overnight film exposure)

## New information identifies ideal techniques for nonradioactive labeling and detection.


Those products, part of the Genius™ System with Lumi-Phos™ 530 from Boehringer Mannheim, guarantee detection of 0.1 pg of target nucleic acid in less than one hour. In addition, they allow detection of single-copy genes in human DNA in less than 30 minutes when used in Southern blotting techniques.



The Genius products rely on the very same procedures already used for radioactive labeling, yet provide results not just hours earlier, but **days** earlier.

Results can be recorded on X-ray film and, because the chemiluminescence lasts more than 48 hours, multiple exposures are possible. Those films can be easily scanned and quantified by densitometry.

multiple times without relabeling. Furthermore, Southern blots developed with the Genius System can be repeatedly stripped and reprobed using the same procedures as the older radioactive method.




Human genomic Southern blot (RFLP) showing chemiluminescent detection of single-copy genes following a 2.5 minute X-ray film exposure

### Library screening: positive colonies visible in one hour

Using probes labeled with the Genius System, researchers are detecting positive colonies and plaques within one hour versus the overnight exposures required with radioactive techniques.

Colorimetric  
detection of  
recombinant  
lambda gt 10 phage



Because the Genius System can use DNA, RNA or oligonucleotide probes, many research needs can be met with a single system. Probe labeling utilizes the same techniques currently

Performance is guaranteed as the Genius System products are function tested in dot blot (also referred to as slot blot) assays.

The unique digoxigenin-labeled probes are stable for at least one year and can be reused multiple times, offering savings in terms of the time and money normally spent relabeling decayed radioactive probes. By using the same probe, experimental variables are reduced and reproducibility between experiments is provided.

**0.1 pg DNA  
Guaranteed**

### Technique-specific information and bibliography available upon request

To receive additional information about the Genius System products, their applications and specifications, or to receive a bibliography of publication references, call Boehringer Mannheim at **800-428-5433**.

If you would like to speak directly with a specialist about your individual research needs, call **800-262-4911**.



**Genius™  
System**

**BOEHRINGER  
MANNHEIM  
CORPORATION**






# Functional Features

## of the Forma 3250 Cell Culture Incubator




### **Optimum Culturing Chamber**

Chamber design provides optimum culturing environment. Ultra-flat electropolished shelves are square for convenience and perforated for uniform air movement. A stainless steel humidity pan is provided standard to prevent desiccation and facilitate cleaning. Unit interior is comprised of three components which are easily removable without the use of tools.



### **Cleanability**

A crevice-free, stainless steel interior with mirrored finish features 100 percent 1/2" radius coved corners. No harsh finish to harbor contaminants. Construction of the chamber promotes aseptic culturing conditions.




### **Precision Control**

A door-mounted control panel features a microcomputer with vacuum fluorescent alpha/numeric display. Programming is easy with the touch key pad. Self-prompting messages allow you to control and monitor the critical performance parameters of the incubator. Proportional and integral (PI) algorithms combine with fuzzy logic characteristics to provide a stable culturing environment. A patented interactive dual temperature probe system provides precise balance of chamber air for ultimate temperature control and uniformity.



### **Functional Design**

The 3250 water jacketed incubator gives you the flexibility to address any lab space requirements. Set on a countertop or stack them for maximum utilization of space. The incubator door swing can be reversed to accommodate lab layout or for side-by-side operation.




### **Interface Protection**

Remote monitoring of temperature and CO<sub>2</sub> can be achieved through remote alarm contacts, millivolt output signal and optional RS232 interface. Standard electrical outlet is provided for connection to accessory items, such as temperature recorder, gas guard, etc.



### **Product Compliance**

U.L. Listed and C.S.A. Certified, every 3250 incubator is manufactured to conform with the intense guidelines of these testing facilities. Copies of U.L. and C.S.A. documentation are available from Forma for compliance verification.



**Forma Scientific, Inc.**

Box 649, Marietta, OH 45750 1-614-373-4763

Fax: 1-614-373-6770 Telex: 29-8205


**Toll Free 1-800-848-3080**

An ISO 9001 Company

Circle No. 69 on Readers' Service Card



# SEQUENASE<sup>®</sup> IMAGES<sup>™</sup> NON-ISOTOPIC SEQUENCING



## New Sequenase<sup>®</sup> Images<sup>™</sup> Non-Isotopic DNA Sequencing System

For fast, non-radioactive detection and reliable results.

- **Features Sequenase Version 2.0 T7 DNA Polymerase**
- **A simple blotting procedure using Sequenase Images membrane**
- **Hard-copy results on standard X-ray film**
- **Reactions to results in a single 8-hour day**
- **Compatible with Mn Buffer<sup>®</sup>, dITP, Sequence Extending Mix and Glycerol-Tolerant Gel Buffer<sup>†</sup>**

The Sequenase Images System (Product No. 71350) includes the Sequenase Images Non-Isotopic DNA Sequencing Kit (Product No. 71352) and the Sequenase Images Detection Kit (Product No. 71355)

Features of the Images Product Line include:

### Gene Images<sup>™</sup> Detection System

Product No. 75250

#### Gene Images<sup>™</sup> Kit

Product No. 75275

Single copy gene detection

Detection of Northern blots

#### Random Primed Images<sup>™</sup> Biotin Labeling Kit

Product No. 76150

Highly efficient biotin labeling of probes with Exonuclease-free Klenow

### Colony Images<sup>™</sup> Kit

Product No. 75350

Rapid screening of recombinant DNA libraries

### Protein Images<sup>™</sup> Kit

Product No. 71900


Sensitive detection of Western blots

Contact United States Biochemical Corporation,  
P.O. Box 22400, Cleveland, Ohio 44122.

Phone: 800-321-9322. Fax: 800-535-0898.

International: 216-765-5000.

Telex: 980718. Fax: 216-464-5075.



\*Sequenase is a registered trademark of United States Biochemical Corporation. This reagent (kit) contains materials covered by or suitable for use under one or more U.S. Patent Nos.: 4,795,699; 4,946,786; 4,942,130; 4,962,020; 4,994,372, 5,145,776 and 5,173,411. Patents pending in the U.S. and other countries.

\*\*Mn Buffer - Purchase includes a nonexclusive sublicense solely for use with Sequenase<sup>®</sup> T7 DNA Polymerase. No other license is granted expressly, impliedly, or by estoppel.

†Glycerol Tolerant Gel Buffer - Pat. pending.



# USB<sup>™</sup>

**United States Biochemical**


USB INTERNATIONAL OFFICES: Cambridge, England • Bad Homburg, Germany INTERNATIONAL DISTRIBUTORS: Alpha Scientific - Thessaloniki, Greece • BioTrade - Austria • Cambridge BioScience - Cambridge, England • Getter (2000) Ltd. - Tel-Aviv, Israel • Integro B.V. - Holland • Kabo - Spanga, Sweden • Ballerup, Denmark • Espoo, Finland • Oslo, Norway • Labex Corp. - Calcutta/New Delhi, India • Lucerna-Chem AG - Lucerne, Switzerland • Quimigranel-SA - Spain/Portugal • Societa Italiana Chimici - Roma, Italy • Sopachem - Belux • Touzart et Matignon - Paris, France • Toyobo Company, Ltd. - Osaka, Japan • Trace Scientific - Melbourne, Australia

Circle No. 32 on Readers' Service Card



**The system:  
Analytical microscopy workstation**

**The instrument:  
Nikon Microphot-SA**



The advanced Microphot-SA research microscope is designed to be the heart of your analytical workstation, today and tomorrow.

Capable of being run from a host computer, it is uniquely adapted to complex multi-imaging requirements. For instance, perform simultaneous real-time image processing and microspectrophotometry. Or put world famous Nikon CF® optics to work for both high resolution video microscopy and photomicrography, in all techniques including DIC and fluorescence.

To arrange for a demonstration, call **1-800-52-NIKON, ext. 144.**

See us at the FASEB Meeting,  
Booth No. 1619

Circle No. 36 on Readers' Service Card

**Nikon**

AT THE HEART OF  
THE SYSTEM YOU NEED



# The Film Critics' Choice



**Automatic scanning  
of up to 10 films**

## New! FilmMate™ DNA Film Scanning System.

Those with a critical eye for DNA sequencing films will appreciate the

FilmMate DNA Film Scanning System. • Its high-capacity image

processing capability eases your workload. Load up to 10 films and it auto-

matically scans, digitizes, and stores the images for analysis. • The FilmMate scanner

uses 12-bit processing to capture autoradiogram images for the highest resolution film

imaging. Variable density settings enable you to optimize resolution for a broad

array of films. • Automatic base calling and ambiguity marking with the new



**Accurate BaseWorks  
sequencing software**

BaseWorks™ software give you accurate control of sequence

identification. You can also take advantage of powerful sequence comparison

and alignment capabilities, as well as versatile report formatting. And it's



**Highest resolution  
film imaging**

Macintosh®-based and System 7™ savvy, so networking and

multi-tasking are incredibly easy. • IntelliGenetics/Betagen has over a decade of

field-proven experience providing products, services, support, and training for

DNA sequencing. For a complete demonstration or for more information, contact your

local representative at 1-800-876-9994.



**Industry-leading support**



**IntelliGenetics**



**Betagen**

Please visit us at  
Experimental Biology '93,  
Booth 812

IntelliGenetics/Betagen, Mountain View, CA: 415.962.7300 fax 415.962.7302 IntelliGenetics, Europe: +32.3.219.5352 fax +32.3.219.5354  
FilmMate and BaseWorks are trademarks of IntelliGenetics, Inc. Macintosh is a registered trademark and System 7 is a trademark of Apple Computer, Inc.

Circle No. 64 on Readers' Service Card





## There Are Places You Expect To Find Water Hazards. Your Lab Shouldn't Be One Of Them.

Water hazards on a golf course can cost you a couple of strokes. In a lab, however, water hazards in the form of ions, pyrogens and organics can cost you months of research work.

That's why it's critical for you to have a water purification system that removes contamination instead of one that contributes to it. If you don't, you could end up like one of these facilities:

1. A lab doing protein purification with HPLC found that endotoxins in its water system contaminated proteins.
2. Organics in another facility's water system kept mammalian cells from growing in defined medium.
3. A semiconductor plant testing for copper had its chromatographic results ruined by contamination in their water system.
4. In yet another lab, phosphates



being used to clean the water system got into the water used for polyacrylamide gels, rendering them useless.

In each of these cases, when the lab switched to a Milli-Q® Plus Water Purification System, all their problems disappeared.

That's because the Milli-Q Plus system delivers the highest quality water available: 1.5 liters per minute at 18 megohm and < 10ppb TOC (with RO pretreatment).

To ensure you're always getting high quality water, the system tells you when the filter module needs changing.

And unlike other systems, changing our module is easy. Instead of separate bowls, our QPAK™ purification module has an all-in-one design that takes hardly any time to

change out and prevents you from ever having to come in contact with wetted surfaces.

Also, because we don't use separate bowls, the system has virtually no hold-up volume. This results in faster flush-ups, along with reduced extractables and bacterial growth.

To further prevent water degradation, the system even recirculates water automatically every hour for approximately five minutes.

So beware of less expensive water purification systems that claim to do as good a job as the Milli-Q Plus system: They may prove hazardous to your work.

For additional information about the Milli-Q Plus system, including a brochure, technical bulletins and data sheet, call 1-800-225-1380.



**MILLIPORE**

Circle No. 46 on Readers' Service Card

©1992 Millipore Corporation



# NATURE POSES SOME ENORMOUS CHALLENGES

One of today's greatest challenges to life science researchers is the full characterization of the carbohydrate moieties of glycoproteins.

The fundamental procedures required to characterize protein glycosylation – glycan release and profiling – have now been automated by Oxford GlycoSystems.

To learn more about how the automation of these procedures can help your research, call Oxford GlycoSystems today for free technical and applications literature.



Tel: +44 (0)235 553066 Fax: +44 (0)235 554701  
Tollfree from:

North America and Canada: 1-800 722 2597

UK: 0800 212061; Germany: 0130 81 37 48


France: 0590 86 08; Switzerland: 155 2739

Circle No. 23 on Readers' Service Card





# Bio-Rad: The Evolution of Gene Transfer



At Bio-Rad we don't have a gene transfer instrument—we have a gene transfer program. One that Bio-Rad has pioneered, researched, and provided to the life science community. Plant and animal cells, bacteria, other cells: Bio-Rad not only gives you the tools to efficiently transform them, but also supports you throughout your research with unparalleled sup-

port, free telephone consultation, and unlimited access to published protocols and yet unpublished data. We have been dedicated to gene transfer since 1986, and we will continue to further your success with high quality electroporation and particle delivery products.

Call 1-800-4BIORAD today to discuss your particular gene transfer needs.

\* Tang, D., DeVit, M. and Johnston, S.A., *Nature*, 356, 152-154 (1992). Biolistic is a registered trademark of the DuPont Company.

**BIO-RAD**

**Bio-Rad  
Laboratories, Inc.**

**Life Science  
Group**

U.S. (800) 4BIORAD • California Ph. (510) 741-1000 • New York Ph. (516) 756-2575 • Australia Ph. 02-805-5000 • Austria Ph. 0222-877 89 01 • Belgium Ph. 091-85 55 11 • Canada Ph. (416) 624-0713 • China Ph. 2563146 • France Ph. 01-49 60 68 34 • Germany Ph. 089-318 84 0 • Italy Ph. 02-21609 1 • Japan Ph. 03-3534-7515 • Hong Kong Ph. 7893300 • The Netherlands Ph. 08385-40666 • New Zealand Ph. 09-443 3099 • Scandinavia Ph. 46 (0) 8 590-73489 • Spain Ph. (91) 661 70 85 • Switzerland Ph. 01-810 16 77 • United Kingdom Ph. 0800 181134

Circle No. 41 on Readers' Service Card



# TO CREATE YOUR FINEST WORK YOU NEED THE RIGHT TOOLS

Life science researchers are at the frontier of one of today's greatest challenges – a deeper understanding of the biological role of carbohydrates. To meet that challenge, investigators need the right tools to probe and characterize oligosaccharides to the finest level of detail possible – tools which include instrumentation, reagents, kits and analytical services.

Oxford GlycoSystems is dedicated to providing the highest quality glycobiology research tools plus technical literature, applications information and the latest research news. To help you create your finest work and ensure that you have the right Tools for Glycobiology, call us today for free technical information.



**Oxford GlycoSystems**

Tel: +44 (0)235 553066 Fax: +44 (0)235 554701  
Tollfree from:  
North America and Canada: 1-800 722 2597  
UK: 0800 212061; Germany: 0130 81 37 48  
France: 0590 86 08; Switzerland: 155 2739

**Circle No. 24 on Readers' Service Card**





From Jandel Scientific, the makers of **Java**<sup>®</sup> and **SigmaPlot**<sup>®</sup>


**NEW!** **Mocha**<sup>™</sup> Image Analysis Software for Windows<sup>™</sup>

# Simply the most powerful software for automating your research.

Did you know your desktop PC has the power to analyze your color and gray scale images? With MOCHA<sup>™</sup> image analysis software for Windows<sup>™</sup> you can measure and quantify up to 60,000 objects. Find intensities, areas, angles, lengths, distances and more. Then analyze and plot your results – all on your own PC! Whether you're in biological, industrial, or other research, MOCHA will greatly enhance the accuracy and efficiency of your analysis.

## Powerful Measurement Capabilities

Count fibers, size particles, quantify gels, measure neurons, and more with MOCHA's



Analyze color and black and white images.

broad selection of point, line, and area measurements. Four overlay planes provide unparalleled versatility. Collect image data, draw major/minor axes, and number objects, for example, all in one step! Just specify your measurements and let MOCHA do the rest!

## Multiple Image Formats

MOCHA accepts a range of image formats so you can work with images from video cameras, scanners, or even from disk. With features such as pseudo-clearfield equalization, multiple image splicing, floodfill measurements, and a host of enhancement filters, you can process and analyze your images swiftly, with or without access to the original image source!

## Fully Integrated Analysis

MOCHA's exclusive, SigmaPlot<sup>®</sup> compatible data worksheet, built-in plot functions, and convenient annotation tools relieve you from using cumbersome business spreadsheets and drawing programs. MOCHA has the versatility to process images, analyze the data, then graph and print results all within one package! No other image analysis software offers this level of integration.


## Configure To Your Needs

MOCHA can analyze saved or scanned color

and b/w images on a standard PC (386 or better) using a single monitor, with or without a frame grabber. If you work with live images, MOCHA supports a variety of popular color and monochrome frame grabbers in

both single and dual monitor configurations. This flexibility means you can use MOCHA

for desktop image analysis (no additional hardware required), or as a workstation for capture and analysis of live video images.



Intensities, areas, counts and more!

## Call For Demonstration


See a demonstration and discover how MOCHA image analysis software can automate *your* research. Jandel Scientific backs MOCHA with a network of authorized dealers, free technical support, and a 90-day money back guarantee.

For a FREE brochure and name of the dealer nearest you call toll-free (800) 874-1888 (within the U.S.) or (415) 453-6700.

**Jandel**  
SCIENTIFIC  
"Microcomputer Tools for the Scientist"

2591 Kerner Blvd.  
San Rafael, CA 94901  
(800) 874-1888  
(415) 453-6700  
Fax: 415-453-7769


In Europe:  
Schimmelbuschstr. 25  
4006 Erkrath, Germany  
Phone: +2104/36098  
Fax: +2104/33110




Integrated data worksheet.




# One Fluorimeter Fits All.




*DNA Quantitation*




*Gel Scanning*




*Macromolecular Anisotropy Measurements*




*Intracellular Ion Determinations*




*Enzyme Kinetic Analysis*



*Microplate Assays*



*Stopped-Flow Kinetic Analysis*



*Membrane Fluidity Measurements*

## FLUORESCENCE

Enzymes or membranes. Ions or DNA. Perkin-Elmer's new Model LS-50B Luminescence Spectrometer lets you look at cells in many different ways. This unique, computer-controlled spectrometer measures fluorescence, bioluminescence, chemiluminescence, and phosphorescence, faster, more accurately and with less photo bleaching.

No other luminescence spectrometer gives you more options to meet all of your needs. You can measure multiple fluorescent probes — FURA-2, BCECF, INDO-1, SNARF®-1, SNAFL®-2 — in milliseconds with the Fast Filter Accessory. Monitor temperature, stir the sample, and

mark when a reagent is added with the biokinetic attachment. Read 96-well microplates or scan TLC plates or gels with the Microplate Reader. Analyze membrane fluidity by measuring changes in anisotropy. Measure the fluorescence of living cells right on the coverslip. And more.

With Perkin-Elmer's worldwide service, technical support and the LS-50B, you've got the perfect fit for any

size bioresearch lab. To order or for more information on the Model LS-50B Luminescence Spectrometer, contact your local Perkin-Elmer office. For product literature in the U.S., call 1-800-762-4000.




*The Model LS-50B Luminescence Spectrometer is PC-controlled with Fluorescence Data Manager Software.*

**PERKIN ELMER**

The Perkin-Elmer Corporation, Norwalk, CT 06859-0012 U.S.A. SNARF and SNAFL are registered trademarks of Molecular Probes, Inc.

For information circle reader service number 38  
For sales call circle reader service number 39





## INTRODUCING THE FIRST INVERTED CONFOCAL LASER SCAN MICROSCOPE DESIGNED BY EAR.

"More lasers." "Additional detectors." "Increased flexibility." We heard your requests and sent them straight to our engineers. This is their response—the new LSM 410 Inverted.

The LSM 410 features previously unheard of benefits. Simultaneous laser scanning and detection of up to three different fluorescence signals at up to 1024 x 1024 pixels with 24-bit resolution. Time series software for high-speed ratio imaging for kinetic analyses of spatial changes in cell metabolism. Add 3-D projection with animation and it's far and away the most advanced system available.

Multiple lasers and detectors are just part of the story. With Windows™-based operating software, the LSM 410 makes it easy to accommodate multiple users and applications. The price makes it hard to pass up.

The new LSM 410 Inverted Confocal Laser Scan Microscope. Isn't it about time you got what you've been asking for? Call 800-233-2343.

Multiple  
Fluorescence

3-Color  
Overlay (Live)

Time Series  
Kinetic Analysis

Multi-User  
Operation

3-D Projection  
With Animation

Macroprogramming

Carl Zeiss, Inc.

Microscope Division  
One Zeiss Drive  
Thornwood, New York 10594  
800-233-2343 • Fax 914-681-7446



Circle No. 77 on Readers' Service Card


Windows is a trademark of Microsoft Corporation.

SEE  
WHAT  
YOU'RE  
MISSING



# **THE CAMFOLIO APPROACH TO UNDERSTANDING THE COMPLEX PHENOMENA OF LEUKOCYTE- ENDOTHELIAL INTERACTIONS.**

Begin with the most comprehensive selection of unique monoclonal antibodies to human cell adhesion molecules. Portfolios of complementary molecules, formulated for multiple assays, give you the flexibility to explore multiple pathways of interest. Clone-specific references and documented functionality take the guesswork out of experimental design. By letting you define the variables, CAMFolio monoclonal antibodies give you more control over your research. For information on CAMFolio monoclonal antibodies relevant to leukocyte-endothelial interactions or for a complete catalog of CAMFolio monoclonal antibodies, contact Becton Dickinson today.



**ELAM-1 (E-SELECTIN)**  
**GMP-140 (P-SELECTIN)**  
**ICAM-1 (CD54)**  
**INTEGRIN  $\alpha$ 4 (VLA-4)**  
**LECAM-1 (L-SELECTIN)**  
**LFA-1 $\alpha$  (CD11A)**  
**LFA-1 $\beta$  (CD18)**  
**MAC-1 (CD11B)**  
**SIALYL-LE<sup>x</sup>**  
**VCAM-1**

**OUR EXPANDING PORTFOLIO  
OF MONOCLONAL ANTIBODIES  
NOW INCLUDES  $\beta$ 1,  $\alpha$ 3 AND  
EIGHT OTHER INTEGRINS.**

Becton Dickinson  
Advanced Cellular Biology  
2350 Qume Drive  
San Jose, CA 95131-1807  
Ordering (800) 223-8226  
Customer Support (800) 952-3222  
Fax (408) 954-2009

Becton Dickinson Canada, Inc.  
(416) 822-4820  
Fax (416) 855-1243

Becton Dickinson  
European HQ  
Erembodegem-Aalst, Belgium  
(32) 53-720211  
Fax (32) 53-720450

Nippon Becton Dickinson  
Company, Ltd.  
Tokyo, Japan  
(81) 3-3403-9991  
Fax (81) 3-3403-5008

Becton Dickinson Worldwide Inc.  
Singapore  
(65) 861-0633  
Fax (65) 860-1590

CAMFolio is a trademark of  
Becton Dickinson and Company.

**BECTON  
DICKINSON**

Circle No. 43 on Readers' Service Card

For research use only. Not for use in  
diagnostic or therapeutic procedures.




NEW EXTRA CLEAN  
TIP RACK



# The right choice for DNA sequencing

Comfortable pipetting, fine nanolitre control, confidence in reproducibility... Just a few of the reasons why the new Pipetman Models P2/P10 are fast becoming the world's preferred pipettes for DNA sequencing and other biotechnology applications.

Used with Gilson's new miniature C10 tips, Models P2 and P10 offer maximum security, even at 0.1  $\mu\text{L}$ !



Pipetman<sup>®</sup> micro pipettes  
in volumes down to  
0.1  $\mu\text{L}$



 **GILSON<sup>®</sup>**



*Ares-Serono Symposia is an independent foundation, created in 1971, to promote scientific research in all disciplines which contribute towards improving human health. This aim is pursued by means of congresses, courses, seminars and scientific publications.*

pietrasante ray subiali - rome

# THE CARE TO MEET, TO KNOW, TO UNDERSTAND

## SOME OF THE INTERNATIONAL MEETINGS SCHEDULED FOR 1993 ARE:

*"Experimental and Clinical  
Precancerous Lesions: Approaches to  
Cancer Prevention and Early  
Diagnosis"*

P. Marks (USA) and R. Weil (CH)  
**Montreux (CH), March 29-31**

*"Molecular Diagnosis and Monitoring  
of Leukaemia and Lymphoma"*

F. Grignani (I)  
**Perugia (I), April 15-17**

*"Molecular Basis of Inflammation"*

J. Navarro (USA)  
**Heidelberg (D), April 21-23**

*"Metabolism in the Female Life Cycle"*

M.P. Diamond and F. Naftolin (USA)  
**Taormina (I), May 17-18**

*"Recent Advances on Monoclonal  
Gammopathies and Related  
Malignancies"*

B. Barlogie (USA) and F. Dammacco (I)  
**Evian (F), June 3-5**

*"Inhibin and Inhibin-Related Proteins"*

H.G. Burger (AUS)  
**Siena (I), June 17-18**

*"Cell and Molecular Biology of the  
Testis"*

M.L. Dufau (USA) and A. Isidori (I)  
**Majorca (E), September 13-14**

*"GTPase-Controlled Molecular  
Machines"*

D. Corda, S. Garattini and A. Luini (I)  
**S. Maria Imbaro (I), Sept. 22-25**

*"Developmental Endocrinology"*

M.L. Aubert and P.C. Sizonenko (CH)  
**Geneva (CH), Sept. 30 - Oct. 2**

*"The Challenge of Biotechnology: from  
Laboratory Diagnosis to Clinical  
Therapy"*

S.A. Aaronson (USA) and R. Verna (I)  
**Rome (I), October 11-12**

*"Molecular Basis of Endocrine  
Diseases"*

C. Pavía (E)  
**Barcelona (E), November 18-19**

## **G**ordon Research Conferences

*Ares-Serono Symposia is pleased to  
announce its support of the Gordon  
Research Conferences in Italy for 1993.*

*"Biodegradable Polymers"*

**San Miniato (I), May 2-7**

*"Biological Structure and Gene  
Expression"*

**Volterra (I), May 2-7**

*"Organic Superconductors"*

**Il' Ciocco (I), May 9-14**

No form of advertising or promotional activity by any pharmaceutical company is allowed at Ares-Serono Symposia Congresses.



ARES  
SERONO  
SYMPOSIA

Circle No. 13 on Readers' Service Card

I would like to receive information about:

☐ All the Congresses

☐ In particular the Congress on .....

Name .....

Address .....


Institution .....

Please mail or fax to Ares-Serono Symposia,  
Via Ravenna 8 - Rome Italy (fax 6/44291324)



# *Now You Have a Choice!*

**Choose a model from Meridian's ACAS 570 family  
to fit your fluorescence applications...  
and your budget.**



*3-D reconstruction of soybean cell labeled with Hexyl Rhodamine. (Courtesy of Dr. Mel Schindler, Michigan State University, East Lansing, MI.)*

The premier fluorescence instrument for biological applications—the Meridian **ACAS 570** Interactive Laser Cytometer—is now available in a family of five separate models, each configured for a different set of applications. No matter what your quantitative fluorescence work involves, Meridian's new **ACAS 570** family has a model for you:

- *The ACAS 570C* for UV confocal microscopy.
- *The ACAS 570K* for ratio kinetics (imaging and quantitation).
- *The ACAS 570I* for fluorescence imaging and quantitation.
- *The ACAS 570T* for toxicology screening (e.g., "FRAP" or Fluorescence Redistribution After Photobleaching, and cell-cell communication).
- *The ACAS 570S* for adherent cell sorting.

The new **ACAS 570** family makes it easy for you to choose an instrument that fits your needs exactly. Features can be mixed *and* matched, and UV laser excitation can be added to any configuration. Call us today to find the **ACAS 570** for your needs.

***Meridian's ACAS 570 family—  
A model for every fluorescence  
application.***

**MERIDIAN**

Instruments, Inc.

2310 Science Parkway  
Okemos, Michigan 48864  
Phone: 800-247-8084 • 517-349-7200  
Fax: 517-349-5967

**Meridian Instruments Europe, Inc.**  
Industriepark—West 75  
B-9100 St. Niklaas BELGIUM  
Phone: 32-3-7801760 • Fax: 32-3-7781727



**Meridian Instruments Far East K.K.**  
The Sonehara Building (2nd Floor)  
2-23-3 Higashi Nihonbashi  
Chuo-ku, Tokyo 103 JAPAN  
Phone: 03-5820-3315 • Fax: 03-5820-3316

**Distributors** located in Australia, Austria, Benelux, France, Germany, Hong Kong, Israel, Italy, Korea, New Zealand, Portugal, Spain, Sweden, Switzerland, Taiwan, and the United Kingdom.

Circle No. 37 on Readers' Service Card

©1993 Meridian Instruments, Inc.



HUMAN

MONOCLONAL

CLONE	BINDS TO AMINO ACID RESIDUES NO.	NEUTRALIZING
<i>Anti-gp41 anti-gp160 monoclonal antibodies</i>		
IAM 41-3D6	env 604-617	No
IAM 41-4D4	env 590-600	No
IAM 41-4B3	env 580-604	No
IAM 41-25C2	(a)	No
IAM 41-2F5	env 662-667	Yes
IAM 41-3H12	env 693-856	No
IAM 41-5F3	(a)	No
<i>Anti-gp120 anti-gp160 monoclonal antibodies</i>		
IAM 120-2G12	(b)	Yes
IAM 120-3D12	(b)	Yes
IAM 120-1B1	(b) (blocks CD4 binding)	Yes
IAM 120-2G6	(b)	No
<i>Anti-p24 monoclonal antibodies</i>		
IAM 24-37G12	structural epitope	No
IAM 24-3A6	p 24 122-149	No
IAM 24-1D7	unknown	No

(a) Binding site at residues 539-542, 556-566, or 683-693 of gp160, or structural epitope.  
(b) Binds to carbohydrate moieties anchored between residues 001-201, 216-254, 285-302, 328-414, or 434-475 of gp160.  
All IAM anti-gp41 and anti-gp120 antibodies also bind to gp160.

HIV-1

Targeted to the



scourge of mankind

ANTIBODIES


\*Research Use Only

Each of our monoclonal HIV-1 antibodies comes in sterile, lyophilized 0.1mg portions with a technical data sheet, priced at \$300 - \$370 per vial.

FOR ORDERS AND INFORMATION CALL:  
VTS, 1-800-323-2437

Circle No. 82 on Readers' Service Card





## Reach for the Peak of GMP Peptide Synthesis

### Reaching New Heights in GMP Capability

For 20 years Bachem California has set the standard for peptide purity and customer service, bringing you products of unequalled quality for your research and clinical applications. No one can match our experience or credentials. Bachem is the world's leading producer of GMP peptides for clinical applications with over 30 Drug Master Files with the FDA and other regulatory agencies throughout the world. Bachem California is the only company worldwide that has been inspected and approved by both the U.S. FDA and England's MCA for GMP production of pharmaceutical grade peptides.

### Bachem's Proven Track Record Speaks for Itself

- The rigorous standards we have developed for every GMP-grade product are your guarantee of quality and consistency.
- The unmatched experience we have gained in supporting hundreds of projects involving GMP-grade peptides means you can count on Bachem for the documentation and regulatory support your project demands.
- Our unique capability in large-scale synthesis and purification is your assurance of a reliable product supply now, and in the future, as your requirements grow.

**Bachem California - Exceeding Your Expectations for GMP-Grade Peptide Production**



**Bachem California**  
3132 Kashiwa Street • Torrance, CA 90505  
Tel: (310) 539-4171  
Telefax: (310) 530-1571  
**Circle No. 84 on Readers' Service Card**


**Bachem (UK), Ltd.**  
A subsidiary of Bachem California  
69 High Street • Saffron Walden  
Essex CB10 1AA, England  
Tel: (0799) 526465  
Telefax: (0799) 526351

**Saxon Biochemicals, GmbH**  
A subsidiary of Bachem California  
Feodor-Lynen-Strasse 14  
3000 Hannover 61, Germany  
Tel: (0511) 530080 • Telefax: (0511) 5300888



# FETALCLONE®

A Family of High Performance FBS Alternatives



## Cell Performance At Its Best for About Half the Cost of FBS

Decreasing the amount of fetal bovine serum (FBS) used in cell culture is a significant factor in lowering costs. HyClone's new FETALCLONE® products are designed to do just that — deliver both cost-effectiveness and high performance.

FETALCLONE® products are manufactured and processed with the same care as other HyClone sera using our "closed-system" technology and "true pool" methods in large lot sizes. These steps result in consistent performance.

FETALCLONE® is available in two options: FETALCLONE®I, optimized for the growth of hybridomas and other related mammalian cell lines; and FETALCLONE®II, optimized for growing CHO-K1 and other epithelial cells. Both products demonstrate many of the same growth characteristics inherent in FBS, and with comparable IgG levels.



Low cost, high performance, and stable supply make FETALCLONE® THE ideal FBS alternative.

Contact us for FREE samples for evaluation.

### Start Switching Your Cells to FETALCLONE® and Watch Your Savings Grow!

For example, assume current usage of FBS is at 50 liters per year, at \$320/L\*, total FBS cost = \$16,000.

For every liter of FetalClone used in place of FBS, there will be savings of about \$160\*\*. Consider the following table:

Current FBS Requirements	FBS Cost (Price @ \$320/L*)	% Switch to FETALCLONE®	Total Cost of Using FETALCLONE® & FBS	COST SAVINGS
50 Liters	\$16,000			
		10%	\$15,200	\$ 800
		30%	\$13,600	\$ 2,400
		50%	\$12,000	\$ 4,000
		70%	\$10,400	\$ 5,600
		90%	\$ 8,800	\$ 7,200

\*Hypothetical pricing based on current FBS market prices only. Not to be interpreted as actual price.

\*\* Reflects average current U.S. FETALCLONE®.




"Innovative Products for Cell Culture"

1725 South HyClone Road, Logan, Utah 84321 • Phone: 1-800-HYCLONE (492-5663) • FAX: 1-800-533-9450

Circle No. 44 on Readers' Service Card



## EASILY PURIFY PROTEINS WITH PRP-X500 COLUMNS



PRP-X500 SUPER-  
FICIALLY POROUS  
ANION EXCHANGE  
COLUMNS COMBINE  
THE BEST OF  
POROUS AND NON-  
POROUS PACKINGS

**HIGH CAPACITY**  
0.2mg PER RUN FOR  
EACH SAMPLE  
COMPONENT

**GOOD RECOVERY**  
PROTEIN RECOVERY IS  
THE SAME AS NON-  
POROUS SUPPORTS

**FAST SEPARATIONS**  
ANALYZE FOUR  
PROTEINS BY ANION  
EXCHANGE IN 3  
MINUTES

FOR MORE INFOR-  
MATION ABOUT PRP-  
X500 PROTEIN ANALYSIS  
COLUMNS CALL (800)  
648-5950 OR CIRCLE THE  
READER SERVICE CARD.

See us at Experimental Biology  
Meeting, Booth 928-930.

**HAMILTON**  
THE MEASURE OF EXCELLENCE

HAMILTON COMPANY  
P.O. BOX 10030  
RENO, NEVADA 89520-0012  
800-648-5950 OR 702-786-7077

# PFANSTIEHL

**We've added something new:  
Thioglycosides.**

Thioglycosides are a convenient source of the glycosyl donors needed for complex carbohydrate synthesis. They are highly crystalline and stable to a wide range of blocking/deblocking conditions. Thioglycosides can easily be converted to glycosyl halides for use under Koenigs-Knorr conditions or activated directly with a variety of reagents. Stereochemical control of the newly-formed glycosidic bond can be achieved by the appropriate group on C-2. Data sheet available.

**PFANSTIEHL LABORATORIES, INC.**

*The source for carbohydrate chemistry*

1219 Glen Rock Avenue/Waukegan, IL 60085-0439  
Tel.: 1-708/623-0370/Toll Free: 1-800/383-0126  
FAX: 708/623-9173  
68-W



## MAKING THINGS WORK

Circle No. 27 on Readers' Service Card

## <sup>14</sup>C MEVALONOLACTONES <sup>3</sup>H ISOPRENOIDS 20-30% OFF\*

CAT #	PRODUCT DESCRIPTION	QTY	AMSM	DPNT	ARC
ART-158	Dolichol, [1- <sup>3</sup> H]	250μCi	\$----	\$----	\$399
ARC-620	Farnesyl-PP, [1,2- <sup>14</sup> C]	5μCi			599
ARCD-101	Farnesyl-PP, "cold"	100μg			109
ART-324	Farnesyl-PP, [1- <sup>3</sup> H]	250μCi	610	650	439
ARCD-102	Geranylgeranyl-PP, "cold"	100μg			109
ART-348	Geranylgeranyl-PP, [1- <sup>3</sup> H]	250μCi	910	630	469
ARC-541	Isopentenyl-PP, [1- <sup>14</sup> C]	50μCi	905	----	679
ART-377	Isopentenyl-PP, [1- <sup>3</sup> H]	50μCi			439
ARC-702	Mevalonic acid-5-P, [2- <sup>14</sup> C]	10μCi			599
ART-406	Mevalonic acid-5-P, [5- <sup>3</sup> H]	50μCi			999
ARC-522	Mevalonic acid, RS- [2- <sup>14</sup> C]	50μCi	260	310	199
ART-334	Mevalonic acid, RS- [5- <sup>3</sup> H]	1mCi	----	950	599
ART-427	Mevalonic acid, R- [5- <sup>3</sup> H]	1mCi	----	950	769
ARC-452	Mevalonolactone, RS- [2- <sup>14</sup> C]	50μCi	240	305	189
ART-315	Mevalonolactone, RS- [5- <sup>3</sup> H]	5mCi	----	1380	999
ARC-540	Squalene, [ <sup>14</sup> C]	1μCi			299
ART-355	Squalene, [ <sup>3</sup> H]	50μCi	----	1150	859

\*Amersham (AMSM) &

DuPont's (DPNT) List Prices

Complete Satisfaction or 100%

Money Back Guarantee

**We do <sup>14</sup>C & <sup>3</sup>H Custom Synthesis**

**AMERICAN RADIOLABELED CHEMICALS, INC.**

11624 Bowling Green Drive, St. Louis, MO 63146

Phone: 1-314-991-4545 FAX: 1-314-991-4692


Toll Free (US & Canada): 1-800-331-6661



Circle No. 93 on Readers' Service Card



# Visit Booth 307 at Experimental Biology '93.




*The New Standard of Excellence.*

## Academic Press Dictionary of Science and Technology

The **Academic Press Dictionary of Science and Technology** offers you, in one volume, the largest collection of scientific terms available in any single reference source of its kind.

1992, 2,240 pp., \$115.00/ISBN: 0-12-200400-0



## Pulsed Field Gel Electrophoresis

*A Practical Guide*

Bruce Birren and Eric Lai

This book provides everything a researcher or advanced student needs to know in order to understand and carry out pulsed field gel experiments.

**Comb Bound: \$39.95**

April 1993, c. 244 pp./ISBN: 0-12-101290-5

## Genome Research in Molecular Medicine and Virology

*edited by*

Kenneth W. Adolph

This volume emphasizes:

- Genome research—DNA sequences and gene structure of chromosomes
- Molecular medicine—gene therapy, cancer research, oncogenes, hereditary disorders, normal DNA, and model systems
- Virology—retroviruses, herpesviruses, and poliovirus

February 1993, 355 pp., \$75.00/ISBN: 0-12-044300-7

## Molecular Biology and Pathology

*A Guidebook for Quality Control*

*edited by*

Daniel H. Farkas


**Comb Bound: \$49.95**, February 1993, 299 pp./ISBN: 0-12-249100-9

## The Leucocyte Antigen FactsBook

A. Neil Barclay, Marian L. Birkeland, Marion H. Brown, Albertus D. Beyers, Simon J. Davis, Chamorro Somoza, and Alan F. Williams

A Volume in the FACTSBOOK Series

September 1992, 448 pp., \$42.00/ISBN: 0-12-178180-1



## Gene Expression in Neural Tissues

*edited by*


P. Michael Conn

A Volume in the METHODS IN NEUROSCIENCES Series

**Comb Bound: \$49.95**/ISBN: 0-12-185268-7

Casebound: \$95.00/ISBN: 0-12-185267-9

1992, 492 pp.



## Immunoassay Automation

*A Practical Guide*

*edited by*

Daniel W. Chan

**Immunoassay Automation: A Practical Guide** describes automation of immunoassay from the practical viewpoint of the clinical laboratory.

**Comb Bound: \$49.95**

1992, 367 pp./ISBN: 0-12-167632-3

## Immunochemical Techniques Laboratory Manual

John Goers

This manual, designed for courses and workshops in immunochemical techniques, contains student-tested protocols for affinity chromatography, Western blotting, ion exchange chromatography, immunostaining of cells, ELISAs, and more than 30 other methods.

**Comb Bound: \$29.95**, January 1993, 216 pp./ISBN: 0-12-287048-4

## Veterinary Virology

Second Edition

Frank J. Fenner, E. Paul Gibbs, Frederick A. Murphy, Rudolf Rott, Michael J. Studdert, and David O. White

February 1993, 666 pp., \$75.00/ISBN: 0-12-253056-X

## Astrocytes

*Pharmacology and Function*

*edited by*

Sean Murphy

February 1993, 458 pp., \$99.00/ISBN: 0-12-511370-6

## Reconstitution of Intracellular Transport

*edited by*

James E. Rothman

A Volume in the METHODS IN ENZYMOLOGY Series

This volume includes sections on reconstitution in cell-free extracts, reconstitution using semi-intact and perforated cells, and identification of transport intermediates.

November 1992, 512 pp., \$75.00/ISBN: 0-12-182120-X

## RNA Methodologies

*A Laboratory Guide for Isolation and Characterization*

Robert E. Farrell, Jr.

**Comb Bound: \$49.95**, February 1993, 317 pp./ISBN: 0-12-249700-7



Order from your local bookseller or directly from

**ACADEMIC PRESS**

HBJ Order Fulfillment Dept. DM17915  
6277 Sea Harbor Drive, Orlando, FL 32887

Prices subject to change without notice. ©1993 by Academic Press, Inc. All Rights Reserved. JE/AB/DC — 26033


CALL TOLL FREE

**1-800-321-5068**

FAX 1-800-336-7377



# Capillary Electrophoresis



## It's Time You Took A Second Look At CE Sensitivity

New levels of sensitivity and reproducibility have been achieved by Applied Biosystems' Capillary Electrophoresis System (Model 270A-HT). We've now eclipsed our own industry-leading standard by more than 30-fold, ensuring results you can't attain from any other system.

Whether you need high-efficiency separations for research or quality control, the automated 270A-HT provides the flexibility to handle a range of chemistries and applications. And now we've developed a new,


proprietary optical technology especially for our CE systems—increasing sensitivity by more than an order of magnitude. You can achieve this unprecedented sensitivity without modifying applications or protocols you currently use on our systems. Additional sensitivity-enhancing strategies, such as stacking techniques, can increase performance even further.

With the 270A-HT you can measure reaction rates—instead of evaporation rates—because our specially designed


sample vials and thermostatted sample handling system eliminate evaporation for precise, run-to-run reproducibility.

There are more reasons to take another look at Applied Biosystems capillary electrophoresis systems. In addition to providing the highest performance CE instrumentation, we're committed to keeping our customers on the leading edge with a full selection of specialty application kits, responsive service and expert technical assistance. For free product information phone Applied Biosystems today at: **Australia** (03) 808-7777, **Benelux** (0) 3465-74868, **Canada** (800) 668-6913, **France** (1) 48 63 24 44, **Germany** (0) 6150/101-0, **Italy** (02) 8912 6011, **Japan** (03) 3699-0700, **U.K.** (0925) 825650, **U.S.** (800) 345-5ABI.

©1993 Applied Biosystems, Inc.



*Antihypertensive drug formulation analysis with impurity sensitivity capability of <0.01%.*



*10 pmol/mL digest of  $\beta$ -lactoglobulin.  
Top graph: High Sensitivity option.  
Bottom graph: Standard Instrument.*


**Applied Biosystems**

*A Division of Perkin-Elmer Corporation*


Circle No. 40 on Readers' Service Card



# Spin less and spend less with CHROMA SPIN<sup>TM</sup>.




You've been there—  
waiting for samples  
to finish. You should  
be getting more done.



But instead, you're looking at your watch. **Obtain purified DNA or RNA in just 12 minutes with CHROMA SPIN Columns.** Not only do they work faster, they are more versatile. CHROMA SPIN Columns are based on the principle of gel filtration, and come in five matrix pore sizes and three different

buffers. They can thus be used for buffer exchange, purification, and size selection—a 3-fold size differential is all you need for good separation. Which means you can use CHROMA SPIN to purify DNA, RNA, or proteins quickly and efficiently—in half the time, and without the hassle of syringes or ethanol precipitation. **Why wait?**

**Think how much faster you'll finish your experiments.** You'll save preparation time and get results that much sooner. Don't wait any longer.



**Stop going round-and-round.**  
Call 1-800-662-CLON or your local distributor  
for our free booklet.

Come visit us at Experimental Biology Booth No. 901, 903, 1000 and 1002

4030 Fabian Way, Palo Alto, California 94303 USA • 800/662-CLON • 415/424-8222 • FAX: 415/424-1064 • 800/424-1350

Circle No. 71 on Readers' Service Card

**CLONTECH**



## CALL FOR ABSTRACTS

# SEVENTH ANNUAL NORTH AMERICAN CYSTIC FIBROSIS CONFERENCE

Loews Anatole Hotel  
Dallas, Texas  
October 13-16, 1993

*A multidisciplinary  
conference for  
researchers, physicians  
and allied health  
professionals on  
the latest advances in  
CF research and care.*

### PROGRAM CHAIRMEN

Bonnie W. Ramsey, M.D.  
Jeffrey A. Whitsett, M.D.

## ▲ PLENARY SESSIONS

*The CF Gene: Old Questions—  
New Insights*  
Francis S. Collins, M.D., Ph.D.

*Cystic Fibrosis: Electrolyte  
Transport Revisited*  
Michael J. Welsh, M.D.

*Gene Therapy for CF: A Glimpse  
Into the Future*  
Richard C. Boucher, M.D.  
Ronald G. Crystal, M.D.  
Jeffrey A. Whitsett, M.D.  
James M. Wilson, M.D., Ph.D.

▲ **SYMPOSIA:** Rites of Passage—Protein Folding and Transit • Vector Systems • Management of Infections in CF • Controversies in Physical Therapy • The Biliary Tree & Liver Disease in CF • Ages & Stages: The Developmental Perspective • Neutrophil-Mediated Oxidant Injury to the Lung & Anti-oxidant Therapy • Etiology of Energy Metabolism in CF • New Perspectives on Reproductive Issues in CF • Pharmacology of ATP-sensitive Ion Channels • Model Systems • Submucosal Gland Biology • Advances in Clinical Therapies • Macro & Micronutrients • Topology and Domains • Host Vector Interaction in Gene Therapy • On Being the CF Family • Late-Breaking Science • Topics in Transplantation

▲ **SHORT COURSES:** Human Subjects Research: Design & Grantsmanship • Modern-Day Molecular Genetics for the Clinician • Pulmonary Physiology • The Microbiological Basis for Anti-microbial Therapy • Technological Advances in CF Home Care • A Comprehensive Educational Program for CF Patients/Caregivers

Conference proceedings will be published as a supplement to *Pediatric Pulmonology*. **ABSTRACT DEADLINE: May 20, 1993.**


For abstract forms and registration information, contact:



6931 Arlington Road  
Bethesda, MD 20814  
Attn: Caroline McPherson  
(301) 951-4422 Fax (301) 951-6378

Circle No. 87 on Readers' Service Card

NEW




## PROPHET

A National Computing Resource for Life Science Research

1	2	3
C1	C2	C3
1	12	a
2	24	b
3	36	c

Table



Graph

$$\sum x^2$$

$$R^2$$

Analysis

cgu ggc ugc

Arg Gly Cys

Sequence

NC(=O)CNC(=O)C

Molecule

### Making Life Easier for the Life Scientist

The new PROPHET software package gives Life Scientists tools for storing, analyzing, and presenting *Data Tables*, *Graphs*, and *Statistical Analyses*. PROPHET's point-and-click graphical user interface lets you choose commands from pull down menus, enter information through dialog boxes tailored for each task, and invoke a variety of analysis tools. *Molecular Modeling* and *Sequence Analysis* capabilities, currently available with a command interface, will soon be supported under the graphical interface. PROPHET's new and easy-to-use interface lets you concentrate on life science – not computer science.

For only \$510 you can get PROPHET in your lab or office.

**The PROPHET Group**  
**BBN Systems and Technologies**  
 10 Moulton Street Cambridge, MA 02138  
 (617) 873-2669 prophet-info@bbn.com

The National Center for Research Resources, NIH, sponsors the continuing improvement of PROPHET.

Circle No. 83 on Readers' Service Card



# New Grant?

**...Then get the latest about  
what's new in physiological  
research instrumentation**

**New** automated syringe pumps  
**New** gas anesthesia  
**New** micromanipulators  
**New** stereotaxic instrumentation  
**New** MacLab™ data acquisition

**Call (708) 860-9700** or send today for your free copy of Stoelting's 80 page, **1993** catalog for physiological and pharmacological research.



**STOELTING**

Stoelting Co., 620 Wheat Lane  
Wood Dale, IL 60191 USA  
**(708) 860-9700**

Circle No. 20 on Readers' Service Card



SEE US AT EXPERIMENTAL BIOLOGY, BOOTH #1403

# How Far Does Your Stats Package Go?



"If you're looking for the most fully functional statistical package, look no further than StatView 4.0." — MacUser

StatView® 4.0 is the only software package that offers all the tools that scientists and researchers need to analyze data and present results — in a single application. StatView 4.0 integrates data management, statistical analysis, graphing, drawing, and color presentation into one powerful, elegant program. From data entry to presentation, StatView 4.0 dramatically saves you time and effort.

## NOT JUST ANOTHER STATS PACKAGE

With StatView 4.0, Abacus Concepts introduces many unique features not found in traditional data analysis software, including the ability to:

- access data, analyses, and results via floating browsers
- update results automatically when data changes
- generate new data columns using other variables and/or over 150 built-in functions
- specify and save simple or complex criteria to restrict analyses to data subsets
- customize and annotate graphs and tables with a full set of

drawing and presentation tools

Best of all, you can save your work as templates and apply them to any other data. Because templates remember everything you do — analysis and graph specifications, formatting, text, drawings — you save time and never have to repeat your work.

StatView 4.0 was developed so that new statistical and graphing features may be easily added as modular extensions. Look for our first module — Quality Control Tools — this summer.


These features (and many more), combined with seamless integration and unparalleled ease of use, make StatView much more than just another statistics package.

## CRITICS AND USERS AGREE: STATVIEW IS THE BEST


The editors of *MacUser* awarded StatView the 1992 Eddy for Best New Science/Engineering Product. Our users have voted StatView five *Macworld* World Class Awards. And it's the most popular statistics software on the Macintosh. Before you try the rest, try the best.

## FIND OUT MORE

StatView 4.0 has a list price of \$595 (U.S.). Academic discounts are available. To learn more, please call Abacus Concepts at 1-800-666-STAT.



"Abacus Concepts scores with new enhancements that push this analysis program beyond the competition." — MacWeek



"StatView offers a nice mix of statistical features and operating convenience, coupled with outstanding control over printed output." — Macworld

**ABACUS**  
**CONCEPTS**

1984 Bonita Avenue, Berkeley, CA 94704-1038  
Tel: 510-540-1949; Fax: 510-540-0260

## Upgrade Customers!


If you own an older version of StatView, you should upgrade to StatView 4.0 at our great low price. Call for details today!

© 1993 Abacus Concepts, Inc. StatView is a registered trademark of Abacus Concepts. All other trademarks are property of their respective holders. Overseas offices: France: (+33) 76-41-84-30; Switzerland: (+41) 32-23-54-65; Belgium: (+32) 25-34-45-35; UK & Ireland: (+44) 0 865-784800; Japan: (+81) 03-590-2311; Sweden: (+46) 08 - 753 22 23; Germany: (+49) 030 31 59 13-0; Australia & New Zealand: (+61) 03 326 5255.

Circle No. 91 on Readers' Service Card



# The Plate.



## MicroBeta™ serves all of your micro-plate label assays right through to results print-out

In applications like, *cytotoxicity, coated plate immunoassays, homogeneous immunoassays/SPA, labelled cell assays, membrane receptor assays, Cerenkov counting, multi-label assays and work with coated beads*, the use of microtitration plates has become universal. Accessories and sample handling equipment are designed for this format. Now, at last, the counting end is secured too. With MicroBeta, we can offer you the best available counting technology, moulded to suit the way you work - on standard microtitration plates.



### MicroBeta™ - another 1st from Wallac

- ◆ Allows fast multidetector counting of 16 plates, any combination of 24- or 96-well, per loading
- ◆ Free plate counting - does not require special, expensive opaque plates
- ◆ Saves space - it's only a quarter the size of a conventional LSC
- ◆ Assures all the reliability of the liquid scintillation counting technique
- ◆ Is supported by everything you need in the way of accessories and applications info
- ◆ Comes in a version to suit your exact needs

MicroBeta is a trademark of Wallac Oy, Finland


Wallac Inc., 9238 Gaither Road, Gaithersburg,  
MD 20877, Telephone: 301-963-3200,  
800-638-6692  
Telecopier: 301-963-7780

**wallac**

Circle No. 33 on Readers' Service Card

WE MAKE BETTER SCIENCE POSSIBLE.

## The Company.




1st rack-based counters  
- RackGamma (LKB-Wallac, 1975), RackBeta™ (1977)



1st ultra low-level liquid scintillation counter  
- Quantulus™ (LKB-Wallac, 1984)



1st time-resolved fluorometry immunodiagnostic system  
- DELFIA® (Wallac, 1985)




1st LSC specifically for harvested samples  
- Betaplate (Wallac, 1986)



1st LSC for microtitration plates  
- MicroBeta (Wallac, 1990)





# One Hour Protocol

including DNA elution/purification from agarose gels and ligation reaction, achieves subcloning or size-fractionated deletion cloning within one hour.

■ By combining SUPREC™-01, SUPREC™-02 and the DNA Ligation Kit, the ONE HOUR PROTOCOL facilitates DNA recovery followed by cloning with speed and simplicity. Furthermore, the amount of DNA required for this protocol could be at the level comparable to what would generally be used for analytical electrophoresis (0.5 ~ 2 µg for subcloning and 5 µg for size fractionated deletion cloning).

## Products Required

### 1. SUPREC™-01 (Code No.9040)

Microcentrifuge cartridge for DNA recovery from agarose gel.

### 2. SUPREC™-02 (Code No.9041)

Microcentrifuge cartridge for sample concentration and/or for removal of small DNA fragments.

### 3. DNA Ligation Kit (Code No.6021)

Two-component system containing all the reagents necessary to perform DNA ligation in 30 minutes with considerably high transformation frequency.



Please contact your local representative for details and pricing.

## One Hour Protocol

### Step 1

Cut 1~2 µg of DNA with restriction enzyme and separate by agarose gel electrophoresis.

Cut out the gel containing the DNA band of interest.

Treat 5 µg of plasmid DNA with Kilo-Sequence Deletion Kit (Code No.6030) and separate by electrophoresis.

Cut out the gel containing the DNA fraction of required size.

### Step 2

Put the agarose gel slice into the inner cartridge of SUPREC™-01 and centrifuge at 10,000 rpm for 5 min.

### Step 3

Transfer the filtrate into the inner cartridge of SUPREC™-02 and centrifuge at 6,000 rpm for 3 min. Add 15 µl of TE and collect DNA retained on the filter membrane by pipetting.

### Step 4

Add 80 µl of A solution and 15 µl of B solution of the Ligation Kit into the 10 µl of recovered DNA solution. Allow to stand at room temperature (16~25°C) for 15 min.

### Step 5

Add the 10 µl of Ligation solution to 100 µl of competent cells and leave in the ice bath for 15 min. After heat shock (at 42°C for 30 sec), add 890 µl of SOC culture medium and cultivate with shaking at 37°C for 15 min.

### Step 6

Concentrate the culture medium into 200 µl and plate 10 µl and all the rest respectively.



**TAKARA SHUZO CO., LTD.**  
Shijo-Takakura, Shimogyo-ku,  
Kyoto 600, Japan  
Phone: +81 75-241-5177  
Fax: +81 75-241-5208

**TAKARA BIOCHEMICAL INC.**  
719 Allston Way, Berkeley  
California 94710, USA  
Phone: +510-649-9895  
Fax: +510-649-8933  
Toll Free: 1-(800)544-9899

**Distributors:**  
**Europe**  
**PBS-Organics**  
Parc de l'Innovation BP 209,  
67405 Illkirch Cedex France  
Tel: 88 67 08 30  
Fax: 88 67 38 61

**ITC Biotechnology, GmbH**  
Postfach 103026, 6900 Heidelberg, Germany  
Tel: (06221)303907  
Fax: (06221)303511

**British Bio-technology Products Ltd.**  
4-10, The Quadrant, Barton Lane, Abingdon  
OX14 3YS, UK  
Tel: (0235)529449  
Fax: (0235)533420

**Taiwan**  
**Cheng Chin Trading Co., Ltd.**  
No. 75 Sec. 1, Hankow St. Taipei Taiwan  
Tel: (02)331-3111 Fax: (02)382-2197

**Canada**  
**BIO/CAN SCIENTIFIC**  
2170 Dunwin Drive, Unit 5  
Mississauga, Ontario CANADA L5L 1C7  
Tel: 416-828-2455  
Toll Free: 1-800-387-8125  
Fax: 416-828-9422





"We're making headway on those buckyballs."    *"This could be a real coup for us."*    "I wonder if anyone's hot on our tail."



*"Doesn't look like it."*    "Not even the French lab?"    *"No way. We're months ahead of them."*



"Really? You have connections in Paris?"    *"No, I have a much better source."*

# THE ANSWER IS IN DIALOG.

When research is moving fast, the best way to keep up with developments—and with competitors—is with Dialog®, the world's most comprehensive online database service.

Dialog takes you beyond the basic research facts of other services to deliver late-breaking news and information. It's the one source you

can count on for everything you need to know, now. And Dialog's superior searching capabilities enable you to find more thorough and precise answers in less time, making it your best value in an online service.

**1-800-3-DIALOG**

Outside U.S., 415-858-3785. Fax 415-858-7069

You can even alert Dialog to notify you automatically of any new developments in your field.

Talk to your company information center about Dialog, or call for free information on Dialog services for the chemistry field.

©1992 Dialog Information Services, Inc., 3460 Hillview Avenue, Palo Alto, California 94304. All rights reserved. DIALOG is a servicemark of Dialog Information Services, Inc., Registered U.S. Patent and Trademark Office.

**DIALOG INFORMATION SERVICES, INC.**


A KNIGHT-RIDDER COMPANY

Circle No. 18 on Readers' Service Card



# VISIT PERGAMON AT EXPERIMENTAL BIOLOGY '93, BOOTHS 319-320!

*Free Sample Copies Available Upon Request!*



## CELL TRANSPLANTATION

Official Journal of the Cell Transplant Society

**Editor-in-Chief:** Paul R. Sanberg, Division of Neurosurgery, Department of Surgery, University of South Florida College of Medicine, Tampa, FL 33612


### Section Editors:

**Neuroscience:** Olle Lindvall; **Islets and Other Endocrines:** David W. Scharp; **Tissue Engineering and Bioartificial Organs:** Patrick Aebischer; **Hematopoietic System:** Robert Peter Gale; **Bone, Cartilage and Muscle:** A. Hari Reddi; **Blood Vessels, Skin and Other Tissues:** Bruce E. Jarrell; **Methods and New Technologies:** Camillo Ricordi

**Cell Transplantation** publishes original research and review articles on the subject of cell transplantation and its application to human diseases. Reports deal with a wide range of topics including physiological, medical, preclinical, tissue engineering and device-oriented aspects of transplantation of nervous system, endocrine, growth factor-secreting, bone marrow, epithelial, endothelial and genetically engineered cells, among others. Reports of basic clinical studies, relevant technological advances, immunology research, patents and ethical and regulatory considerations of cell transplants are all welcome.

### Cell Transplantation features:

- **Original Contributions:** Peer-reviewed, high-quality research investigations that represent new and significant contributions to science.
- **Review Articles:** Reviews of major areas in cellular transplantation. These may be of any length and are peer-reviewed.
- **Brief Communications:** Timely and brief, peer-reviewed articles.
- **Letters to the Editor:** Readers' comments on journal articles and other matters of interest to transplant researchers.
- **Announcements and News:** Notices of upcoming meetings, conferences, seminars, and other events of interest to those in the field.



## CALL FOR PAPERS!

## CELLULAR AND MOLECULAR BIOLOGY RESEARCH

(Incorporating Cellular and Molecular Biology)

**Editor-in-Chief:** M.A.Q. Siddiqui, State University of New York Health Science Center at Brooklyn, New York

**Cellular and Molecular Biology Research (Incorporating Cellular and Molecular Biology)** publishes original articles emphasizing the wide range of topics concerned with cellular and molecular mechanisms involved in cell growth, development and differentiation. Among those areas of coverage are *in vitro* and *in vivo* studies of structure and function of biologically important macromolecules in eukaryotic organisms, gene regulatory mechanisms, transcription, structure and function of regulatory factors, carcinogenesis, oncogenes and tumor suppressor genes in growth and differentiation, signal transduction mechanisms, structure and function of cell organelles, control of growth and development, cell-cell interaction, embryogenesis, and structure and function of viruses and bacteriophages. Papers with emphasis on recombinant DNA technology, transgenes, and genetic manipulation of cells and organs are of special interest.

## CALL FOR PAPERS

The journal considers full-length, original articles, short communications, summary articles, techniques papers, and commentaries.

The Editorial Board of **CMBR** is committed to objective peer review and rapid publication. Manuscripts are published in the shortest time possible commensurate with high scientific quality. Please complete the form below to receive complete guidelines for electronic submission and manuscript preparation.

Please send me: ☐ Free sample copy ☐ Instructions for Contributors

☐ **Cell Transplantation** ISSN: 0963-6897 Vol. 2, 1993 6 issues per annum

☐ **Cellular and Molecular Biology Research** (Incorporating Cellular and Molecular Biology) ISSN: 0968-8773 Vol. 39, 1993 8 issues per annum

Name \_\_\_\_\_

Address \_\_\_\_\_

Institution/Department \_\_\_\_\_

City/State/Zip \_\_\_\_\_

Return to: **PERGAMON PRESS**

US: Marketing Department, 660 White Plains Road, Tarrytown, NY 10591-5153, USA. Fax: 914-333-2444

UK: Marketing Department, Headington Hill Hall, Oxford, OX3 0BW, England. Fax: 011-44-865-743952


A member of the Elsevier Science Publishing Group

Circle No. 9 on Readers' Service Card




# If you don't like the numbers you get from a phosphor plate imager,

Phosphor Plate Imager



AMBIS Imaging System



## just rerun the sample.

Which of their numbers should you believe? Run 1? Run 2? Run 3? None of the above? If reproducible data is important to you, an AMBIS Quantitative Imaging System is what you need.

Our imager directly measures radioactivity in counts per minute (cpm). So the results are statistically valid and can be compared easily between



experiments and from lab to lab.

Phosphor plate imaging uses an indirect approach, much like x-ray film, so it provides arbitrary and highly variable units. To prove it, run the same sample a few times and watch the numbers change! Think what that could do to your research reports.

So if you want your results to count, use the only imaging system that counts. AMBIS Direct Imaging.

Our automated sample changer and easy-to-use software make AMBIS the productivity choice as well.

Don't settle for data that won't hold still. Call AMBIS today, at 1-800-882-6247.



AMBIS, Inc. 3939 Ruffin Road, San Diego, CA 92123 • (800) 882-6247 • (619) 571-0113

\* Identical samples, 1 hour exposures. Data on file at AMBIS, Inc.


©1993 AMBIS, Inc.

Circle No. 75 on Readers' Service Card



# It's Key To Serum-Free !

*Rapid Attachment  
and Reduced  
Lag Time Will  
Speed Adaptation  
To Serum-Free  
Media!*



## PROnECTIN<sup>®</sup>F

**PROnECTIN<sup>®</sup>F can:**

*Eliminate the need for  
fetuin, serum fractions  
or blood-derived  
attachment factors.*

*Allow very low total  
protein or protein-free  
formulations to be  
successfully used.*

*Simplify the formulation  
of custom media*

### APPLICATIONS:

- Routine Adherent Cell Culture
- Recombinant Protein and Vaccine Production
- Virology
- In Vitro Diagnostics
- In Vitro Toxicology
- In Situ Hybridization
- DNA Transfection

If you have been struggling to get your anchorage-dependent cell lines adapted to serum-free media... If you have considered going serum-free but dread the laborious serum weaning protocols... If you're looking for better performance out of your current media... **PROnECTIN<sup>®</sup>F** is the innovative answer you've been waiting for.

**PROnECTIN<sup>®</sup>F** is a genetically engineered protein polymer that presents multiple RGD receptor-specific binding sites on coated cell culture substrates which allows the immediate attachment and spreading of numerous cell types—the critical step in the successful culture of adherent cell lines.

For application reports, additional information or to place an order, please call:

**800-755-0407**




## Protein Polymer

TECHNOLOGIES

10655 Sorrento Valley Road, San Diego, CA 92121  
(619) 558-6064 Fax (619) 558-6477

Circle No. 90 on Readers' Service Card



Vero cells were passed directly from DMEM with 10% FBS to serum-free media in 75 cm<sup>2</sup> tissue culture flasks, both uncoated and precoated with PROnECTIN<sup>™</sup>F. Subsequent passages were in the two serum-free media as noted. Initial seedings were 5 x 10<sup>6</sup> cells per ml. After 2 and 4 days, cells were released and counted with a hemacytometer.

For Research Use Only. Patents Pending.  
©1992 Protein Polymer Technologies, Inc.  
RGD Sequence Licensed from Telios Pharmaceuticals, Inc.



# THE EMBO JOURNAL

**Your chance to save over 80% on a  
personal subscription to The EMBO Journal**



*A subscription to the EMBO Journal in 1993 will include at least 5200 pages, in  
12 monthly issues. Plus an extra bonus issue in December.*

- **the world's most intensively cited monthly journal of original molecular biological research (ISI, October 1992)**
- **world-wide contributions and editorial board ensure a truly international perspective**
- **covers current research into all areas of molecular and cellular biology from gene structure and virology to molecular genetics and immunology**
- **rapid review system to ensure swift publication of important information**

## **Now available to personal subscribers at a special discount**

Why wait for your turn on the departmental list? A personal subscription offers you the chance to read each issue as it comes out – and to have all the back numbers there when you need them. *What is more, first time personal subscribers receive an extra discount.*

**Institutional rate**  
\$625

**Personal subscriber rate**  
\$149 (save \$476)

**First time personal subscriber rate**  
\$110 (save \$515)

<b>SPECIAL PERSONAL SUBSCRIPTION ORDER FORM – SAVE OVER 80%</b>		SC1/93/1
<p>I would like my own subscription to The EMBO Journal. Please airmail 1993 issues to the following address:</p> <p>Name <input style="width: 80%;" type="text"/></p> <p>Institution <input style="width: 80%;" type="text"/></p> <p>Address <input style="width: 80%;" type="text"/></p> <p><input type="checkbox"/> I enclose my check payable to Oxford University Press</p> <p><input type="checkbox"/> Please charge my Visa/Mastercard/American Express/Diners Card No. <input style="width: 80%;" type="text"/></p> <p>Expiry date <input style="width: 20%;" type="text"/></p> <p>Signature <input style="width: 80%;" type="text"/></p>	<p><input type="checkbox"/> I am a first time subscriber (please check)</p> <p>Please return this form to: Oxford University Press Inc., 2001 Evans Road, Cary, North Carolina 27513, USA. or Oxford University Press, Walton Street, Oxford OX2 6DP, UK.</p> <p> For a faster service, you can call us, in the US, on (919) 677 0977 or fax a copy of this order form to us on (919) 677 1714. Outside the USA, call us in the UK on +44 865 56767 or fax us on +44 865 267782.</p>	
<div style="display: flex; justify-content: space-between; align-items: center;"> <div style="background-color: black; color: white; padding: 5px 20px; font-weight: bold;">OXFORD JOURNALS</div> </div>		

Circle No. 60 on Readers' Service Card



With Chemiluminescent  
Sample Detection

# Major Breakthrough!

## Phosphor Image Analysis for Molecular Biology



**New Model GS-250 Molecular Imager™ System.** The Model GS-250 Molecular Imager system represents a breakthrough in rapid imaging of radioactivity and chemiluminescence. Now you can image your gels, dot blots, slot blots, or any radioactive or chemiluminescent sample quickly, without the need for x-ray film, intensifiers, or low temperature film storage and exposure.

The Model GS-250 Molecular Imager system offers the following advantages over traditional autoradiography:

- Reduced Exposure Times
- Increased Sensitivity
- More Accurate Quantitation
- Erasable and Reusable Storage Phosphor Screens
- Lower Cost per Image

(Microsoft, and Windows, are registered trade marks of Microsoft Corporation)

### How it Works:

Simply expose your radioactive or chemiluminescent sample to the proprietary storage phosphor screens. (The screens are sensitive to not only  $^{32}\text{P}$ , but also  $^{35}\text{S}$ ,  $^{14}\text{C}$ ,  $^{125}\text{I}$ , X-rays, AMPPD, and others.) In minutes, you will have generated a latent image which can be quickly scanned by the patented fiber optic laser scanner. Displaying, analyzing, and outputting the resulting data from your 1-D or 2-D protein and DNA gels or blots is fast and simple.

**High Performance Software:** Great technology is wasted without great software! We offer sophisticated software packages for both Windows® compatible and Macintosh® computers; you choose the graphical user interface you prefer. Bio-Rad software will let you quickly quantify your DNA, RNA, or protein samples; perform molecular weight analysis; plot 1-D density profiles or histograms; and perform a myriad of image viewing and processing manipulations.

For more information call 1-800-4BIORAD.

(Macintosh is a registered trademark of Apple Computers, Inc.)

**BIO-RAD**


**Bio-Rad  
Laboratories, Inc.**

**Life Science  
Group**

U.S. (800) 4BIORAD • California Ph. (510) 741-1000 • New York Ph. (516) 756-2575 • Australia Ph. 02-805-5000 • Austria Ph. 0222-877 89 01 • Belgium Ph. 091-85 55 11 • Canada Ph. (416) 624-0713 • China Ph. 2563146 • France Ph. 01-49 60 68 34 • Germany Ph. 089-318 84 0 • Italy Ph. 02-21609 1 • Japan Ph. 03-3534-7515 • Hong Kong Ph. 7893300 • The Netherlands Ph. 08385-40666 • New Zealand Ph. 09-443 3099 • Scandinavia Ph. 46 (0) 8 590-73489 • Spain Ph. (91) 661 70 85 • Switzerland Ph. 01-810 16 77 • United Kingdom Ph. 0800 181134

Circle No. 35 on Readers' Service Card





## **Look to NEN<sup>®</sup> Research Products. The leader in quality, safety and support.**

The surest way to keep pace with complex life science research is to work with a partner who can help you meet the challenge. That partner is DuPont NEN Research Products.

When you needed safer, more convenient radiochemical packaging, we set a new standard with the NENSURE<sup>™</sup> system.

We responded to environmental concerns by starting the first recycling program for plastic foam packaging and lead shielding for radioactive materials. Currently available in the

U.S., this program will soon be extended to other countries.

And, of course, we offer you a choice between radiolabeled products or nonradiometric labels. All backed by DuPont's customer service and technical support.

For more details, call your DuPont representative. Or for information on our newest products by fax, 24 hours a day, seven days a week, call DuPont FaxBack<sup>®</sup> at 1-800-666-6527 (or +1 302 892-0616) and request document #1107.

United States 1-800-551-2121 • Canada 1-800-387-8391 • Australia +61 (008) 257149 • Belgium +32 (02) 724 2717 • Denmark +45 31506610  
France +33 (01) 69825450 • Germany +49 (06172) 87-2600 • Italy +39 (055) 247-8044 • Japan +81 (03) 3224-8763 • Latin America/Asia Pacific Fax  
+1 (617) 426-2464 • Netherlands +31 (073) 206550 • Sweden +46 (08) 7503700 • Switzerland +41 (01) 8410330 • United Kingdom +44 (0438) 734865

