


AMERICAN
ASSOCIATION FOR THE
ADVANCEMENT OF
SCIENCE

SCIENCE

22 APRIL 1994
VOL. 264 • PAGES 477-628

\$6.00



ADVICE TO THE STATISTICALLY WEARY.

Just call 1-800-4-JANDEL to order SigmaStat today!
-SigmaStat Advisor

Dear SigmaStat Advisor,
Help! I'm a scientist, NOT a statistician. I need a foolproof way to pick the best statistical tests for my research data.
Signed,

Looking for Answers
P.S. I'm running on Windows.

Dear Looking,
Look no further! The Jandel SigmaStat Advisor can help. Simply type your research goals into new SigmaStat for Windows and the "Advisor" feature automatically suggests the best statistical test for your data. SigmaStat also checks your assumptions and evaluates the validity of your results. It can even suggest the ideal sample size. New SigmaStat for Windows, puts your own expert statistician right on your desktop!
-SigmaStat Advisor

Dear SigmaStat Advisor,
Help! I tripped and dropped a test tube! What do I do now that my data is unbalanced?
Signed,
Statistically Unbalanced

Dear S.U.,
You need Jandel SigmaStat! Other statistical packages either use erroneous substitution methods or make you throw out valuable data points. SigmaStat for Windows uses a sophisticated general linear model to give you the most accurate results possible from your "messy" data.
-SigmaStat Advisor


Dear SigmaStat Advisor,
How do I get a copy of new SigmaStat for Windows? (I heard from a trusted colleague that SigmaStat is the fastest, easiest way to do statistical analysis on my research data.)
Signed,
Searching for SigmaStat

Dear SigmaStat Advisor,

Try New Jandel SigmaStat™ for Windows! New SigmaStat for Windows comes complete with your own expert statistician right in the box! The SigmaStat Advisor appears with a click of your mouse to make your statistical worries go away. What's the ideal sample size? Which statistical test is best? Are your assumptions and results valid? SigmaStat has the answers.

The SigmaStat Advisor Automatically Picks the Best Test for You. SigmaStat automatically points you to the appropriate tests. You simply type in your goals, and SigmaStat's Advisor automatically suggests the best statistical test for your data!

SigmaStat Gives You Better Results from "Messy Data." Other statistical software packages throw out valuable data. Not SigmaStat for Windows. SigmaStat uses a sophisticated general linear model to fill data gaps, so you get more complete results and better conclusions.




SigmaStat takes full advantage of Windows environment.

SigmaStat Checks Your Assumptions. SigmaStat automatically checks that your data fit the assumptions underlying the statistical methods you're using. If there's a problem, SigmaStat will automatically suggest a more appropriate test!

SigmaStat Lets You Print Detailed Reports and Stunning Graphs. SigmaStat automatically creates a detailed report of your test results—complete with its own analysis! Plus you can automatically print publication-quality graphs of your SigmaStat results or cut and paste to your favorite Windows application.

**TO ORDER SIGMASTAT TODAY
CALL 1-800-4-JANDEL**

Competitive and customer upgrades available. Quantity, academic, and student pricing available.



Jandel
SCIENTIFIC SOFTWARE

2591 Kerner Blvd., San Rafael, CA 94901 Phone: (800) 4-JANDEL, (415) 453-6700 Fax: (415) 453-7769, Internet: sales@jandel.com


In Europe: Jandel Scientific, Schimmelbuschstrasse 25, 40699 Erkrath, Germany 02104/36098 02104/33110 (FAX)
For international dealers call: Australia 2 958 2688, Denmark 45-42150544, Canada 519 767 1061, France 059-03755, Japan 3 3590 2311, Switzerland 41-617121616, Taiwan 2 785 3202, UK 0800-894982

© 1994, Jandel Corp. All companies and/or product names are trademarks of their respective companies.

Circle No. 40 on Readers' Service Card




Cytokine Immunoassays



Sensitive...Precise...Accurate.

Standard curve for IL-6 in serum diluent using the 2 hour shaker option.



Cytokines exert their effects at extremely low concentrations in the presence of many other bioactive proteins. Quantikine immunoassays provide precise determination of cytokines in very complex biological samples.

- Sensitivity in the picogram range
- Tested for cross reactivity with all available cytokines
- Assay time - 2.0 - 4.5 hours
- Precoated 96 well strip plates
- Prediluted HRPO-conjugated second antibody
- Diluents for Serum/Plasma and Cell Culture Media
- Low inter-and intra-assay CV's

To obtain a catalog, detailed product information or to place an order:

In North America Contact:

R&D Systems
614 McKinley Place N.E.
Minneapolis, MN 55413
Telephone: 800-343-7475
Fax: (612) 379-6580

Freephone:

Belgium: 078 11 04 68
Denmark: 80 01 85 92

In Europe Contact:

R&D Systems Europe
4-10 The Quadrant, Barton Lane
Abingdon, Oxon OX14 3YS
Telephone: +44 (0235) 531074
Fax: +44 (0235) 533420

France: 05 90 72 49
Netherlands: 060 225 607

In Japan Contact:

Funakoshi Co., Ltd.
9-7, Hongo 2-Chome
Bunkyo-ku, Tokyo 113
Telephone: +81 (03) 56841622
Fax: +81 (03) 56841633

Norway: 800 11 033
Sweden: 020 79 31 49


R&D
SYSTEMS

1-800-343-7475

Circle No. 36 on Readers' Service Card

GENEAMP PCR INSTRUMENT SYSTEMS





There's No Limit To What You Can Achieve With Our Thermal Cyclers.

In research and product development, being first makes all the difference. GeneAmp® PCR Instrument Systems get you started quickly and easily. They are available with pretested reagents and protocols. Our standard of quality—the highest in the industry—translates into consistently superior performance and unmatched reliability.

In fact, our track record in reliability has allowed us to reduce prices on all GeneAmp PCR Instrument Systems and offer an extended warranty.

What about the future? GeneAmp PCR Instrument Systems are the tools for continuing innovations in laboratories around the world, including our own. This means new methodologies can be applied to your research as soon as they are available. Our new Applied Biosystems Division integrates our expertise and experience in PCR with a reputation for innovation, leadership and support in nucleic acid synthesis, genetic analysis and protein research. Together, we offer you the most comprehensive range of systems, technologies and support for your research.

Invest in the PCR instrument system that won't limit your achievements. For information about our new reduced prices and extended warranty, call 1-800-345-5224 in the U.S., or contact your local Perkin-Elmer sales representative.

PERKIN ELMER

Europe Weiterstadt, Germany Tel: 49-6150-101-0 Fax: 49-6150-101-101
Canada Mississauga, Canada Tel: 800-668-6913 Fax: 905-821-8246
Japan Tokyo, Japan Tel: 81-4-7380-8500 Fax: 81-4-7380-8505
Latin America Mexico City, Mexico Tel: 52-5-651-7077 Fax: 52-5-593-6223
Australia Melbourne, Australia Tel: 61-3-212-8585 Fax: 61-3-212-8502

Perkin-Elmer PCR reagents are developed and manufactured by Roche Molecular Systems, Inc., Branchburg, New Jersey, U.S.A.




Perkin-Elmer is a registered trademark of The Perkin-Elmer Corporation.
GeneAmp is a registered trademark of Roche Molecular Systems, Inc.
The GeneAmp PCR process is covered by U.S. patents owned by Hoffmann-La Roche Inc. and Hoffmann-La Roche Ltd.

Circle No. 28 on Readers' Service Card



496
Bruce Alberts
reshapes the NAS



558
Blue light
phototropism

NEWS & COMMENT

- Dr. Alberts Comes to Washington 496
Making Public Schools the Place to Be
- NIH Tightens Clinical Trials Monitoring 499
- Pesticides and Breast Cancer: No Link? 499
- Cancer Prevention: Beta-Carotene: Helpful or Harmful? 500
- Astronomy in India: Big Science in a Developing Country 501

RESEARCH NEWS

- Vaccines Get a New Twist 503
- Pulsing Star Confirms More Planets in the Universe 506
Other Planets, Other Searches
- Gene Transfer to Spark a Failing Heart 507
- A New Laser Promises to Put an End to Band Gap Slavery 508

- Structure Meets Function at Materials Gathering 510

PERSPECTIVES

- Sulfate Aerosols and Polar Stratospheric Cloud Formation 527
M. A. Tolbert
- The Prion Connection: Now in Yeast? 528
C. Weissmann
- Structural Clues to Prion Replication 530
F. E. Cohen, K.-M. Pan, Z. Huang, M. Baldwin, R. J. Fletterick, S. B. Prusiner

ARTICLE

- Diet and Health: What Should We Eat? 532
W. C. Willett

RESEARCH ARTICLE

- Confirmation of Earth-Mass Planets Orbiting the Millisecond Pulsar PSR B1257+12 538
A. Wolszczan

DEPARTMENTS

- THIS WEEK IN SCIENCE** 485
- EDITORIAL** 487
Chemical Ecology
- LETTERS** 489
Livermore Lab Head: C. B. Tarter • Viral Recombination in Transgenic Plants: M. Mellon and J. Rissler; G. Bruening and B. W. Falk • Adolescent Health Study: D. Alexander; Eds. • Proposed NSF Violence Center: J. A. Mercy • Element Naming: E. K. Hulet • Inhibition of Rev-Mediated HIV-1 Expression: Retraction: M. Campbell, B. K. Felber, G. Nasioulas, E. Afonina, G. N. Pavlakis
- SCIENCESCOPE** 495
- RANDOM SAMPLES** 512
Running Afoul of German Biotech Regs • Horsing With High Tech, etc.
- BOOK REVIEWS** 601
Biological Relationships Between Africa and South America, reviewed by S. D. Webb • *Contrast Sensitivity*, L. O. Harvey Jr. • *The Great Ape Project*, R. H. Tuttle • *Translational Regulation of Gene Expression 2*, D. J. Goss • Vignettes • Books Received
- PRODUCTS & MATERIALS** 609

Board of Reviewing Editors

Frederick W. Alt
Don L. Anderson
Michael Ashburner
Stephen J. Benkovic
David E. Bloom
Floyd E. Bloom
Piet Borst
Henry R. Bourne
Michael S. Brown
James J. Bull

Kathryn Calame
C. Thomas Caskey
Dennis W. Choi
John M. Coffin
Paul J. Crutzen
Robert Desimone
Bruce F. Eldridge
Paul T. Englund
Richard G. Fairbanks
Douglas T. Fearon

Harry A. Fozzard
Klaus Friedrich
Theodore H. Geballe
John C. Gerhart
Roger I. M. Glass
Stephen P. Goff
Peter N. Goodfellow
Corey S. Goodman
Ira Herskowitz
Eric F. Johnson

Stephen M. Kosslyn
Michael LaBarbera
Nicole Le Douarin
Charles S. Levings III
Alexander Levitzki
Harvey F. Lodish
Richard Losick
Diane Mathis
Anthony R. Means
Shigetada Nakanishi

Roger A. Nicoll
Stuart L. Pimm
Yeshayau Pocker
Dennis A. Powers
Ralph S. Quatrano
V. Ramanathan
Douglas C. Rees
T. M. Rice
David C. Rubie
Erkki Ruoslahti

Gottfried Schatz
Jozef Schell
Ronald H. Schwartz
Terrence J. Sejnowski
Ellen Solomon
Thomas A. Steitz
Michael P. Stryker
Robert T. N. Tjian
Emil R. Unanue
Geerat J. Vermeij

Bert Vogelstein
Harold Weintraub
Arthur Weiss
Zena Werb
George M. Whitesides
Owen N. Witte
William A. Wulf
Keith Yamamoto

In the quantum cascade laser, alternating semiconductor layers (red and gray) create quantum wells for electrons. As electrons (light blue streaks) cascade from well to well (left to right), they jump from a higher to a lower energy level, emitting light (yellow band).

The straight blue lines in the well regions represent the energy levels, and the bell-shaped curves, the probability distributions for occupation of that level. See page 553 and the news story on page 508. [Illustration: Keith D. Drake and Frank A. Antalec, AT&T Bell Labs]



REPORTS

Antarctic Total Ozone in 1958 543
P. A. Newman

Strontium Isotopic Composition of Mid-Cretaceous Seawater 546
B. L. Ingram, R. Coccioni, A. Montanari, F. M. Richter

Origin of Saturn's E Ring: Self-Sustained, Naturally 550
D. P. Hamilton and J. A. Burns

Quantum Cascade Laser 553
J. Faist, F. Capasso, D. L. Sivco, C. Sirtori, A. L. Hutchinson, A. Y. Cho

A Mechanism of Lithium Storage in Disordered Carbons 556
K. Sato, M. Noguchi, A. Demachi, N. Oki, M. Endo

A Putative Role of the Xanthophyll, Zeaxanthin, in Blue Light Photoreception of Corn Coleoptiles 558
M. A. Quiñones and E. Zeiger

The Role of Interleukin-6 in Mucosal IgA Antibody Responses in Vivo 561
A. J. Ramsay, A. J. Husband, I. A. Ramshaw, S. Bao, K. I. Mattheai, G. Koehler, M. Kopf

Loudness-Coding Mechanisms Inferred from Electric Stimulation of the Human Auditory System 564
F.-G. Zeng and R. V. Shannon

[URE3] as an Altered URE2 Protein: Evidence for a Prion Analog in *Saccharomyces cerevisiae* 566
R. B. Wickner

Requirement of Vascular Integrin $\alpha_v\beta_3$ for Angiogenesis 569
P. C. Brooks, R. A. F. Clark, D. A. Cheresh

***Bacillus subtilis* SpoIIIE Protein Required for DNA Segregation During Asymmetric Cell Division** 572
L. J. Wu and J. Errington

A Neurochemically Distinct Third Channel in the Macaque Dorsal Lateral Geniculate Nucleus 575
S. H. C. Hendry and T. Yoshioka

Gating Current Noise Produced by Elementary Transitions in *Shaker* Potassium Channels 578
D. Sigg, E. Stefani, F. Bezanilla

Enhanced Myocardial Function in Transgenic Mice Overexpressing the β_2 -Adrenergic Receptor 582
C. A. Milano, L. F. Allen, H. A. Rockman, P. C. Dolber, T. R. McMinn, K. R. Chien, T. D. Johnson, R. A. Bond, R. J. Lefkowitz


TECHNICAL COMMENTS

Muscle Biochemistry and a Genetic Study of Myotonic Dystrophy 587
A. D. Roses; Y.-H. Fu, J. A. Pearlman, A. Pizzuti, R. G. Fenwick Jr., D. L. Friedman, M. B. Perryman, S. Richards, R. A. Gibbs, T. Ashizawa, G. Scarlato, C. T. Caskey

A Truncated Erythropoietin Receptor and Cell Death: A Reanalysis 588
Y. Nakamura and H. Nakauchi

546

Fishing for ancient seawater strontium



AAAS Board of Directors

Eloise E. Clark
Retiring President,
Chairman
Francisco J. Ayala
President
Rita R. Colwell
President-elect

William A. Lester Jr.
Simon A. Levin
Anna C. Roosevelt


Alan Schriesheim
Jean'ne M. Shreeve
Chang-Lin Tien
Warren M. Washington
Nancy S. Wexler

William T. Golden
Treasurer
Richard S. Nicholson
Executive Officer

Indicates accompanying feature

■ **SCIENCE** (ISSN 0036-8075) is published weekly on Friday, except the last week in December, by the American Association for the Advancement of Science, 1333 H Street, NW, Washington, DC 20005. Second-class postage (publication No. 484460) paid at Washington, DC, and additional mailing offices. Copyright © 1994 by the American Association for the Advancement of Science. The title **SCIENCE** is a registered trademark of the AAAS. Domestic individual membership and subscription (51 issues): \$92 (\$50 allocated to subscription). Domestic institutional subscription (51 issues): \$215. Foreign postage extra: Mexico, Caribbean (surface mail) \$50; other countries (air assist delivery) \$95. First class, airmail, student and emeritus rates on request. Canadian rates with GST available upon request, GST #125488122. Change of address: allow 6 weeks, giving old and new addresses and 11-digit

account number. Postmaster: Send change of address to *Science*, P.O. Box 2033, Marion, OH 43305-2033. Single copy sales: \$6.00 per issue prepaid includes surface postage; Guide to Biotechnology Products and Instruments, \$20. Bulk rates on request. Authorization to photocopy material for internal or personal use under circumstances not falling within the fair use provisions of the Copyright Act is granted by AAAS to libraries and other users registered with the Copyright Clearance Center (CCC) Transactional Reporting Service, provided that the base fee of \$1 per copy plus \$0.10 per page is paid directly to CCC, 27 Congress Street, Salem, MA 01970. The identification code for *Science* is 0036-8075/93 \$1 + .10. *Science* is indexed in the *Reader's Guide to Periodical Literature* and in several specialized indexes.




Supercoiled DNA for Transfections

Plasmid DNA Purification using QIAGEN means:

- multiple, ultrapure plasmid preps in 2 hours.
- plasmid yields up to 10 mg.
- no CsCl/EtBr gradient ultracentrifugation.
- no phenol/CHCl₃ extractions.

With the gravity flow QIAGEN-tips, up to 10 mg of ultrapure, supercoiled plasmid DNA can be obtained. Plasmid DNA is selectively purified on QIAGEN-tips which contain our unique anion exchange resin. RNA, protein and polysaccharides are efficiently removed. Ultrapure QIAGEN DNA performs better than or equal to double banded CsCl purified DNA.



Comparison of transfection efficiencies. NIH 3T3 cells were transfected with the plasmid pRSVlacZ which was prepared with QIAGEN Plasmid Mega Kit (blue bar) or by CsCl (gray bar). Data from Ehler et. al. (1993) *BioTechniques* **14**: 546

QIAGEN purified DNA is not exposed to toxic reagents which can potentially inhibit cell growth. QIAGEN-tips are therefore ideally suited for isolating DNA for such demanding applications as transfections, microinjection and gene therapy research.

QIAGEN Plasmid Kits are delivered ready-to-use and include gravity flow QIAGEN-tips, RNase A, all necessary reagents and buffers, and a detailed handbook.

To receive additional information about the QIAGEN Plasmid Kits, or to have a product specialist contact you, please call QIAGEN GmbH, QIAGEN Inc. or contact your local distributor.

Title photo: digitalized reproduction of EM photo of supercoiled DNA prepared with QIAGEN resin.

QIAGEN GmbH Max-Volmer-Straße 4, 40724 Hilden, Germany, Orders (0)2103-892-230, Fax (0)2103-892-222, Technical Service (0)2103-892-240
QIAGEN Inc. 9600 De Soto Avenue, Chatsworth, CA 91311 USA, Orders 800-426-8157, Fax 818-718-2056, Technical Service 800-DNA-PREP (800-362-7737)

DISTRIBUTORS: **AUSTRALIA:** PHOENIX Scientific Pty. Ltd. (03) 868 6888 **AUSTRIA:** BIO-TRADE (0222) 889 18 19 **BENELUX:** Westburg B.V. NL: (033) 95 00 94, B: (07811) 9815 (toll free) **DENMARK:** KEBO Lab A/S (44) 88 72 00 **FINLAND:** KEBO OY (90) 804 4900 **FRANCE:** Cogem (1) 45 32 35 17 **GREECE:** Bio-Analytica (01) 64 62 748 **HONG KONG/CHINA:** Diagnostics Co., Ltd. (892) 542 0566 **INDIA:** Genetix: (11)-550 6463 **ISRAEL:** BIO-LAB Laboratories Ltd. (02) 524 447 **ITALY:** Genenco (M-Medical srl) (055) 5001871 **JAPAN:** Funakoshi Co., Ltd. (3) 5684-1622 **KOREA:** LRS Laboratories, Inc. 924 8697 **NORWAY:** KEBO Lab AS (02) 30 17 18 **PORTUGAL:** Izasa Portugal (01) 758 07 40 **RSA:** Whitehead Scientific Supplies (021) 981-1560 **SINGAPORE:** TIBS TRADING PTE LTD (Biotech Branch) 292 9783 **SPAIN:** Izasa S. A. (3) 401 01 01 **SWEDEN:** KEBO Lab AB (08) 621 34 00 **SWITZERLAND:** KONTRON Instruments AG (01) 733-5-733 **TAIWAN:** Formo Industrial Co., Ltd. (02) 736 7125 **UK:** Hybaid Ltd. (081) 977 3266

Circle No. 1 on Readers' Service Card



Telltale teeth

The strontium isotopic composition of seawater reflects the interaction of several processes, including continental weathering and hydrothermal circulation through ocean crust. The Cretaceous was a time of rapid formation of oceanic crust, but finding fossils that faithfully record seawater Sr compositions despite burial and diagenesis has been difficult. Ingram *et al.* (p. 546) show that fossil fish teeth provide such a record. The overall Sr isotopic variation in the fish teeth accounts for the volume of large Cretaceous oceanic plateaus but is inconsistent with high crustal production rates that have been proposed.

■

Cascade laser


In conventional semiconductor lasers, such as those found in compact disc players, the emitted light is produced as a result of electron-hole recombination across a band gap. Faist *et al.* (p. 553; see cover and news story by Taubes, p. 508) report a quantum well laser that operates on a different principle: the laser action occurs between quantum-confined conduction-band states in a semiconductor heterostructure. The lasers are grown by molecular beam epitaxy, so the dimensions of the quantum structures can be tailored to yield different lasing wavelengths. Lasing occurs at about 4.2 micrometers, a wavelength of interest to molecular spectroscopists and one that is difficult to obtain by other means.

■

Extra lithium

Atomic intercalation of metals into carbon solids, such as lithium intercalation into graphite and alkali metal insertion into

fullerenes, is well known. Sato *et al.* (p. 556) use nuclear magnetic resonance and electron microscopy to show that in a disordered carbonaceous material, which was made by pyrolyzing polyphenylene, lithium could be introduced electrochemically not only in an atomic form but also as discrete Li_2 dimers. The identification of this covalent



site for lithium storage could lead to higher energy densities for batteries; the measured discharge capacity is about double that found for C_6Li electrodes.

■

Defense secretions

A first line of defense against the entry of pathogens through the mucosal linings of the body is the secretion of immunoglobulin A (IgA) antibodies from antibody-containing cells. Ramsay *et al.* (p. 561) show that interleukin-6 (IL-6) is critical for the development of such immunity. Disrupting the IL-6 gene in mice greatly diminished IgA responses in the gut to viral challenges. Normal responses could be restored in the lungs of mice by infection with a recombinant vaccinia virus that carried the IL-6 gene.

Distant pulsar planets

Two years ago, timing variations were found in the radio signals from the isolated 6-millisecond pulsar PSR B1257+12 that suggested the existence of at least two orbiting planets. That discovery has now been confirmed. Wolszczan (p. 538; see news story by Travis, p. 506) explains how a total of 3 years of data show additional systematic variations in the pulse arrival times due to orbital changes caused by the gravitational interaction of the two Earth-mass planets. There is also a third inner planet which has a mass about that of the moon.

Sensing blue light

Blue light affects plants in a variety of ways. Some of the plant's responses to blue light can include developing new chloroplasts, opening stomata, and phototropism, or growth toward the source of the light. The receptor for the various blue light responses has, however, remained elusive. Quiñones and Zeiger (p. 558) manipulated the amount of the carotenoid zeaxanthin in maize coleoptiles and showed that the amount of zeaxanthin correlates with the phototropic response of coleoptiles to blue light. As zeaxanthin is present in very small concentrations, it may be a candidate for a blue light receptor.

■

Unevenly divided

When the bacterium *Bacillus subtilis* is starved, an asymmetric division is initiated, and transcription of certain genes is confined to one or the other of the daughter cells. Wu and Errington (p. 572) analyzed this process and found that the SpoIIIE protein is apparently required for movement of the chromosome into the smaller spore cell, a process that remains incomplete in SpoIIIE mutants. It may be the incomplete segregation of the chromosome in SpoIIIE mutants that leads to some of the normal sporulation—specific genes being transcribed, and others not, depending upon their location on the chromosome.

Shouts and whispers

That the intensity of visual and auditory events varies over many orders of magnitude does not prevent humans from accurately sensing these stimuli. The relation between sensation and intensity is a power function, which reduces the dynamic range to a physiologically acceptable window. Zeng and Shannon (p. 564) used auditory prostheses that had been implanted in the cochlea or brainstem of deaf patients to examine where such transformations might occur between the periphery—the cochlea—and the central nervous system. They suggest that high-frequency stimuli undergo a logarithmic transformation in the cochlea, whereas low-frequency sounds are transmitted linearly until a logarithmic transformation occurs in the cochlear nucleus of the central nervous system.

■


Prying into prions

Prions are infectious proteins that cause diseases of the central nervous system—scrapie in sheep and Creutzfeldt-Jakob disease in humans. How prions infect cells without the involvement of nucleic acids has been puzzling, but in a Perspective, Cohen *et al.* (p. 530) discuss how an abnormal prion protein may act directly on normal host proteins to change their conformation to an abnormal, disease-causing one. In a related report, Wickner (p. 566; see accompanying Perspective by Weissmann, p. 528) discusses how yeast get into the prion business. An obscure yeast mutant, URE, seems to owe its unusual phenotype to a protein, Ure2p, with many of the same properties as prion proteins. This yeast "prion" could prove more amenable to genetic and biochemical analysis.

Most Likely to Succeed.

For non-radioactive nucleic acid detection, the NEBlot™ Phototope™ Kit offers both increased sensitivity and speed.

The NEBlot Phototope Kit rapidly detects nucleic acids in standard Southern and Northern blotting. Biotin is incorporated by use of a unique biotinylated random octamer primer and by incorporation of a biotinylated nucleotide. Resulting probes are then hybridized and detected by chemiluminescence with NEB's Phototope Detection kit.



NEBlot Phototope Kit Catalog # 7550


Random primer biotin labeling for chemiluminescent Northern and Southern blotting and plaque hybridizations. Kit contains reagents for 25 reactions.

Phototope Detection Kit Catalog# 7006, 7030, 7060


Reagents for chemiluminescent detection. Includes streptavidin, biotinylated alkaline phosphatase and Lumigen-PPD®. Available in 3 sizes.

Non-radioactive Phototope detection offers:

- ▶ **Sensitivity:** Equal or greater than radioactive detection.
- ▶ **Speed:** 50 minute process produces films identical to radioactive detection.
- ▶ **Convenience:** No fixing or drying gels required.
- ▶ **Stability:** Biotinylated probes and reagents are stable for extended periods.
- ▶ **Simple Reprobing/Multiple Exposures**



X-ray film detection of Benton-Davis plaque lift using NEBlot Phototope Kit. Exposure time: 10 minutes.



Southern Blot using NEBlot Phototope Kit. Exposure time: 10 minutes.

To experience successful, rapid DNA detection every time, call New England Biolabs: **1-800-NEB-LABS.**

And for those who like it hot...

NEB offers the NEBlot Kit (Cat # 1550-25, 1550-50) which produces high specific activity labeled DNA probes ideal for use in hybridization and blotting procedures.



- **New England Biolabs Inc.**
32 Tozer Road, Beverly, MA 01915 USA
1-800-NEB-LABS
Tel. (508) 927-5054 Fax (508) 921-1350
- **New England Biolabs Ltd., Canada**
Tel. (800) 387-1095 (905) 672-3370
Fax (905) 672-3414
- **New England Biolabs GmbH,
Federal Republic of Germany**
Tel. (06196) 3031 Fax (06196) 83639
- **New England Biolabs Ltd., UK**
Tel. (0800) 31 84 86 (0462) 420616
Fax (0462) 421057

DISTRIBUTORS: Australia (075) 94-0299; Belgium 078 119815; Brazil (011) 66-3565; Denmark (31) 56 20 00; Finland (90) 420-8077; France (1) 34 60 24 24; Greece (01) 5226547; India (542) 311473; Israel (03) 535-1205; Italy (02) 38103171; Japan (03) 3272-0671; Korea (02) 5560311; Mexico (5) 519-3463; The Netherlands (033) 95 00 94; New Zealand (09) 418-3039; Norway 22 22 04 11; PR of China (1) 2544399; Portugal (01) 858 15 64; Spain (01) 594 0806; Sweden (08) 734 8300; Switzerland (061) 481 47 13; Taiwan (02) 8964437

High Definition Real-Time 3D Microscope

INTRODUCING THE EDGE R400:
EXPERIENCE MICROSCOPY'S TRUE DIMENSIONS



The R400 system consists of the Microscope (pictured with Nikon camera, trinocular head, and objectives), an outboard Control Box (containing lamps, filters, and control circuitry), and an outboard power supply (not pictured). Four fiber optic cables deliver "cool" light.

Thanks to its patented* four-beam Multiple Oblique Illumination system, the Edge R400 is the only research light microscope that offers

DIRECTLY OBSERVABLE REAL-TIME 3D

images in true color up to the highest magnification available.

But that's not all. The R400 also delivers significantly better resolution and sharpness, higher contrast, and greater penetrating power than conventional axial-light microscopes.

If two-dimensional images permit recognition, then

3D IMAGES PROMOTE UNDERSTANDING.

Edge believes that, given the choice, researchers would prefer seeing their samples in three dimensions - especially since the R400's true-color 3D images are superior to what can be obtained from comparable 2D microscopes. But when a 2D image is desired, the R400 will present it with unequalled clarity.

So if you'd like to


EXPERIENCE MICROSCOPY'S TRUE DIMENSIONS

please call Edge at (310) 391-3050.

MULTIPLE OBLIQUE ILLUMINATION
REAL-TIME 3D
3D PHOTOMICROGRAPHY
INCREASED RESOLUTION
HIGHER CONTRAST
GREATER PENETRATION



ULTRA-FINE MANUAL FOCUS
MOTORIZED AUXILIARY FOCUS
Z-AXIS FOCUS MEASUREMENT
FIBER OPTIC LIGHTING
PRECISION COMPONENTS
DESIGNED AND BUILT IN THE USA



EDGE SCIENTIFIC INSTRUMENT CORPORATION, 1630 17TH STREET, SANTA MONICA, CA 90404
TEL (310) 391-3050 FAX (310) 391-8715

Circle No. 63 on Readers' Service Card

*ADDITIONAL U.S. AND FOREIGN
PATENTS PENDING



MADE IN
THE USA

THE ICE AGE IS OVER.

New EasyTides™ ambient shipped radiolabeled nucleotides from DuPont NEN.



- **Easy to use.** Shipped in NENSURE™ vials without dry ice. No need to thaw, just store in the refrigerator at 4°C.
- **Easy to see.** Colored with an inert dye for easier pipetting.
- **Easy on the environment.** Less bulky packaging, and no Styrofoam® to recycle.

EasyTides perform like the products you use now. Available in ³²P, ³³P, and ³⁵S-nucleotides. Try EasyTides the next time you order. For information or to order, call 1-800-551-2121. Or call FaxBack at 1-800-666-6527 and request document 9001.

United States 1-800-551-2121 • Canada 1-800-387-8391 • Australia +61 (008) 257149 • Belgium +32 (02) 7242717
Denmark +45 31506610 • France (01) 6982 5450 • Germany +49-6172-872600 • Italy (02) 25302 481/483
Japan +81 (03) 3224-8763 • Latin America/Asia Pacific Fax +1 (617) 426-2464 • Netherlands +31 (073) 206550
Sweden +46 (08) 7503700 • Switzerland +41 (01) 8410330 • United Kingdom +44 (0438) 734865

Manufactured under an ISO 9001 quality system registered by UL® and approved by BSI.
Styrofoam® is a registered trademark of the Dow Chemical Company.



DuPont NEN®



To Dispose Or Not To Dispose... That is the Question.

Now there is an easy answer to this sample preparation dilemma: Packard's New MultiPROBE™.

Carry-Over Elimination... Your Way.

Packard's MultiPROBE offers you a choice: washable Teflon®-coated sampling probes or low cost disposable tips. The MultiPROBE can be configured with conductive liquid level sensing tips or standard tips. Carry-over requirements change? Only a MultiPROBE can be upgraded to meet your changing needs.

High Throughput Liquid Handling.

Automatic variable spacing between multiple sampling tips guarantees high throughput and flexibility for different methodologies regardless of the combination of test tubes, microplates or vials used.

Easy-To-Use.

A graphical user interface with pull-down menus provides versatile, uncomplicated programming for all of today's clinical and life science applications.

On-Site Support.

Add the reputation of Packard's worldwide service and support network and the choice for automating your liquid handling procedures becomes obvious:

The Packard MultiPROBE Family of Automated Liquid Handling Systems... Choices Without Compromise.



PACKARD

A Canberra Company


Packard Instrument Company, 800 Research Parkway, Meriden, CT 06450 U.S.A.
Tel: 203-238-2351 Toll Free: 1-800-323-1891 TX: 643251 FAX: 203-639-2172



Packard International Offices:

Australia, Mt Waverley 61-3-543-4266; **Austria**, Vienna 43-1-302504-0; **Belgium**, Brussels 32-2-4668210; **Canada**, Ontario 1-800-387-9559; **Denmark**, Greve 45-42909023; **France**, Rungis (33) 1 46 86 27 75; **Germany**, Frankfurt (49-69) 663010; **Italy**, Milano 39-2-33910796/7/8; **Japan**, Tokyo 81-3-3866-5850; **Netherlands**, Groningen 31-50-413360; **Tilburg** (013) 423900; **Russia**, Moscow, 7-095-238-7335; **Switzerland**, Zurich (01) 481 69 44; **United Kingdom**, Pangbourne, Berks (44) 0734 844981.

Circle No. 73 on Readers' Service Card



- 152°C ultra-low temperature freezers
- 135°/- 86°C ultra-low temperature freezers
- 40°/- 35°/- 30°C medical freezers

Blood bank refrigerators
Pharmaceutical refrigerators
Cooling rooms
Freezers/refrigerators with totally reliable internal temperatures.

Multi-gas incubators
CO₂ incubators
Cooled incubators
Heated incubators
Plant growth chambers
Incubators for advanced research and testing.



Automatic tablet counting and packing machines
Based on microcomputer advances, these compact units enable accurate tablet dispensing.

Laboratory ware jet washers
Sanyo's new washing technology for scientific use is based on jet water action and microcomputer control.

**Electric ovens
Dry heat sterilizers
Autoclaves**
Equipped with electronic control systems, these models limit temperature fluctuations.

**Bio-clean rooms
Bio-clean benches
Clean benches
HEPA-filter units**
Sanyo's know-how extends to bio and industrial clean rooms.

-152°C



MDF1155
An ultra-low - 152°C temperature freezer with special refrigerant and refrigerating oil for safe, low temperature operation.



MCO-17A
Microprocessor-controlled air-jacketed CO₂ incubator with accurate temperature control and alarm functions. Easy to install and maintain.



MLR-350
Versatile plant growth chamber with programmable temperature and lighting functions and large capacity plus a slim, space saving design.



MPR-411F
Exclusively tailored refrigeration/freezing systems provide a stable preservation environment for reagents and pharmaceuticals.

A wide range of advanced Sanyo technology in support of medicine and research.

Using state-of-the-art electronics technology and control systems, Sanyo makes it possible to perform a wide variety of environmental and medical experiments that are vital for biotechnology research. Our advanced technology used in compressors and ultra-low temperature freezers enables

long-term preservation of cells, bacteria, bone marrow, and other materials essential to the advancement of science. Thanks to the use of CFC-free technology, Sanyo offers an array of outstanding products that are safe for the environment.



■ SANYO Electric Co., Ltd.

■ SANYO Electric Trading Co., Ltd.

5-15, Hiyoshi-cho 2-chome, Moriguchi-City, Osaka, 570 Japan Tel: (06)-992-1521
Fax: (06)-992-1565 or (06)-992-2874 Telex: SANYO J63363, J63353





























■ SANYO Sales & Supply (USA) Corporation

900 N. Arlington Heights Rd. - Suite 300 Itasca, IL 60143, U.S.A. Tel: 708-775-0404 Fax: 708-775-0044

TOLL FREE - 1-800-858-8442

Circle No. 10 on Readers' Service Card

Our Competitors' Enzymes Can't Pass The Quality Testing Performed On GIBCO BRL Restriction Endonucleases[†]

Assay Performed	GIBCO BRL	Company A	Company B	Company C
DNA Endonuclease Assay	All Pass	 	All Pass	
DNA Exonuclease Assay	All Pass	All Pass	 	
DNA Ligation/Recut Assay	All Pass	 	 	All Pass
Nicking Assay	All Pass		 	 
Unit Assay	All Pass	    	  	   
 = One Failure <p>Summary: Seven out of ten restriction endonucleases from each competitor failed one or more of our quality tests.</p>				

[†]Enzymes tested in the manufacturers' recommended buffers: *Bam*H I, *Bgl* I, *Cla* I, *Eco*R I, *Hind* III, *Kpn* I, *Not* I, *Pst* I, *Sal* I, and *Sst* I (*Sac* I)

In a recent quality comparison between Life Technologies and three leading manufacturers, seven out of ten restriction endonucleases from each competitor failed one or more of our quality tests.

Life Technologies offers the highest quality restriction endonucleases on the market because of our unparalleled system of rigorous and complete quality testing. Up to eight different assays are performed at four separate critical stages of production: during the purification process, on the bulk product in storage buffer, on the diluted product, and on the final product you receive.

For reliable, reproducible performance, GIBCO BRL restriction endonucleases beat the competition hands down. Call now for our free restriction endonuclease brochure or for a summary of the comparative data.

To Order/TECH-LINESM: (800) 828-6686

LIFE  TECHNOLOGIESTM

Corporate Headquarters
8717 Grovemont Circle
P.O. Box 6009
Gaithersburg, MD 20884-9980 U.S.A.
(301) 840-8000


U.S.A. Orders
P.O. Box 68
Grand Island, NY 14072-0068 U.S.A.
Facsimile: (800) 331-2286
To Order/TECH-LINE: (800) 828-6686
Industrial Customers: (800) 874-4226


U.S.A. Office for Latin America
To Order/TECH-LINE: (301) 840-4027
Facsimile: (301) 258-8238
Telex: 64210 BRL GARG UW

Canada Orders
Burlington, Ontario L7P 1A1 Canada
To Order/TECH-LINE: (416) 335-2255

*Taq DNA Polymerase is sold under licensing arrangements with Roche Molecular Systems and The Perkin-Elmer Corporation. Purchase of this product is accompanied by a license to use it in the Polymerase Chain Reaction (PCR) process in conjunction with an authorized thermal cycler. TECH-LINESM and the Life Technologies logo are marks of Life Technologies, Inc. For research use only. Not intended for human or animal diagnostic or therapeutic use. 94-001

Circle No. 12 on Readers' Service Card





**The system:
Superior performance
fluorescence microscopy**

**The objective:
Nikon CFN
Plan Apo 60X
1.4 N.A.**

Exceptional throughput, resolution and sensitivity clearly demonstrate Nikon's superiority in fluorescence microscopy. Investigate over an extended spectrum with diffraction limited microscopy, employing very low fluorochrome doses.

High UV transmission and large numerical apertures are the rule throughout the wide range of Nikon CF® fluors and CFN plan apochromats. The 60X CFN plan apo, for instance, provides the brightest, highest resolution images possible over the broadest spectrum.

Call today to schedule a demonstration:
1-800-52-NIKON, ext. 133.

**Fluorochrome
of your choice:
DAPI, AMCA,
FITC, TRITC,
CY5, TR ...**

Nikon

AT THE HEART OF
THE SYSTEM YOU NEED

NEW Control System

UNICORN™ and BioPilot®
With UNICORN and BioPilot, scale up is smooth and predictable, maximizing your return on investment in system control and validation.




A smoother path from laboratory bench to production floor

Smooth transfer from one scale of operation to another is vital when bringing products to market. With UNICORN and BioPilot that transfer is both swift and efficient.


UNICORN is the new control system for BioPilot, the chromatography system for process development, scale up and small scale production. UNICORN complies fully with GLP/GMP requirements and includes validation support.

UNICORN is issued with a Validation Certificate and Report.



UNICORN Control System:

- Developed according to current FDA guidelines for software validation
- Simultaneously supervises up to four liquid handling units of any scale and complexity with control at unit level
- Displays the actual liquid flow path and the readings of the process monitors on a unique dynamic flow diagram
- Keeps the same documentation and user interface during scale up to routine production on BioProcess™ System



UNICORN control system is scaleable. It supervises BioPilot and BioProcess systems keeping documentation and user interface constant.

Ask for new information



Pharmacia Biotech. For information in US: 1-800-526-3593.

In Canada: 1-800-463-5800.

Visit us at Experimental Biology—Booth #1221
Circle No. 34 on Readers' Service Card


Secure and Stable

redivue ^{32}P nucleotides:
now shipping in the
redivial without
dry ice or styrofoam

Ambient shipment of ^{32}P nucleotides demands a secure vial design and a stable formulation. Amersham has launched both with its new Redivial™ and the already successful Redivue™ ^{32}P nucleotides.


The novel vial design uses a splash guard to contain the radioactive solution at the bottom of the "V" insert. This prevents contamination of the vial cap during ambient shipment, when the solution is not frozen. In addition, Redivue ^{32}P nucleotides are stable for days, even at summer shipping temperatures.

redivue - stable at 42°C, 37°C



[^{32}P]dCTP was stored at +42°C for 2 days and then at +37°C for 3 weeks. The control was stored at -20°C. Radiochemical purity was determined by HPLC and biological performance by DNA labeling using the random prime protocol.

redivial



Redivue ^{32}P nucleotides shipped in the Redivial provides the benefits of:-

- ▶ Ready-to-use radiolabel that can be stored in the fridge
- ▶ Ambient shipment without dry ice or styrofoam
- ▶ Identical performance with enhanced security and stability

Shouldn't you buy your nucleotides from the company that thinks through the issues and provides innovative solutions?

FREE BROCHURE



Amersham International plc
Little Chalfont England HP7 9NA
Tel: 0494 544000
Fax: 0494 542929

For further information
contact your local office.


Amersham LIFE SCIENCE

Amersham Corporation
2636 South Clearbrook Drive
Arlington Heights
Illinois 60005
Tel: (800) 323-9750
Fax: (708) 437-1640

Amersham Canada Ltd
1166 South Service Road West
Oakville
Ontario L6L5T7
Tel: (905) 847-1166
Fax: (905) 847-7790

SCIENCE

and HUGO Present



HUMAN GENOME

1 9 9 4



THE GENES AND BEYOND

HORIZONS IN MEDICINE,
BASIC RESEARCH, AND TECHNOLOGY

OCTOBER 2 - 5, 1994
WASHINGTON, D.C.

Mark your calendar now! For the fifth year, Science magazine and The Human Genome Organisation (HUGO) are pleased to sponsor the most comprehensive update on the genome project. The meeting will span issues in mapping, development, signal transduction, genome analysis, medical ethics, and policy.

The genome project is having a revolutionary effect on medical diagnosis, biotechnology, and investigations in cell and molecular biology.

- What are the newest advances in technology?
- How can the genome project be applied to understanding fundamental cellular processes?
- What are the latest developments in identifying disease genes?
- What are the social and policy issues facing us in the future?

For answers to these questions and a complete update on the genome project, don't miss this meeting!

Send me Human Genome 1994 info today!

I'm interested in: ☐ Conference ☐ Exhibition

Name

Title

Organization

Address

City

State Zip

Phone


Fax

Mail or fax to: Global Trade Productions, Inc., 5203 Leesburg Pike,
Suite 1313, Falls Church, VA 22041. Fax (703) 671-7695.

Call (703) 671-1400 today or send
in the coupon for more information.

Unlocking the Secrets
of Programmed Cell Death


Only New ApopTag™ Gives You This Insight Into Apoptosis *In Situ*.



Spontaneous spermatogonial apoptosis in rat testis¹, ApopTag™ Peroxidase Kit and methyl green (X400)



Developmental cell deletion in remodeling of mouse forelimb bud (14-day embryo), ApopTag™ Peroxidase Kit and methyl green (X200)²



Apoptosis of human peripheral blood lymphocytes *in vitro*, ApopTag™ Fluorescein Kit (yellow) and propidium iodide (red) (X400)³

In the most dynamic branches of the life sciences, apoptosis is now recognized as a process of mainstream significance.

At Oncor, we immediately noted the emerging importance of this biological phenomenon and anticipated your need to explore it. Our foresight enables us to put the key to apoptosis investigation in your hands *today*.

Here's why you'll find our new ApopTag™ *In Situ* Apoptosis Detection Kit far superior to other methods:

- ApopTag™ gives you quick, easy and permanent staining of apoptotic cells in the context of the surrounding tissue.
- Able to detect individual apoptotic cells, ApopTag™ provides more information than DNA ladder analysis.
- ApopTag™ is in a complete, simplified, dropper bottle format that saves time and effort.

Call 1-800-77-ONCOR and ask for ApopTag™. Get the first in a series of kits from the company that gives you easy access to apoptosis research.

For Research Use Only. Not for use in diagnostic procedures. ApopTag™ is a trademark of Oncor, Inc.

¹ DJ Allan, BV Harmon & SA Roberts, (1992) *Cell Proliferation*; 25:241-250.

² JFR Kerr, J Searle, BV Harmon & CJ Bishop, in: Potten, CS (ed) (1987) *Perspectives in mammalian cell death*. Oxford U. Press, pp. 93-128. Z Zakeri, D Quaglino, T Latham & R Lockshin, (1993) *FASEB Journal*; 7:470-478; and manuscripts submitted.

³ X Li, W James, F Traganos & Z Darzynkiewicz, (1993) manuscript submitted.



oncor®

209 Perry Parkway
Gaithersburg,
MD 20877

1-800-77-ONCOR
301-963-3500
Fax: 301-926-6129

GENOSYS

CUSTOM OLIGONUCLEOTIDE AND PEPTIDE SYNTHESIS

Arbitrary 10-mer primer kits for fingerprinting and genetic mapping

Each kit includes ten arbitrary 10-mers, selected on the basis of a 50%, 60%, 70% or 80% G+C content, your choice. The arbitrary sequences avoid the presence of complementary ends and palindromes. Additional kits contain 10-mers with common restriction sites, so that amplified products may be digested and cloned for further study. All primers are supplied dry with a guaranteed yield of 0.5 O.D. units. Prices start at \$65 per kit.



With every Genosys oligo, the proof is in the package. Your Quality Assurance Certificate includes a digitized PAGE image, measured yield in O.D.s and µg, and valuable information about your oligo, such as melting temperature, molecular weight and %AT.

**Custom Oligos
\$2.50 Per Base,
Plus \$20. Setup**

Turn down the heat with non-radioactively labeled oligos.

If it's too hot in the lab, Genosys can help cool things down — with non-isotopically labeled custom oligos.

The process is simple. We synthesize your oligo, attach the label of your choice, and deliver your purified product in 7–10 days.

For chemiluminescent or colorimetric detection, any custom oligo can be labeled with biotin or digoxigenin (licensed from Boehringer Mannheim).

For fluorescent detection, choose from fluorescein, the ABI dyes, Texas

red, or rhodamine. You always receive at least 1.5 O.D.s of labeled product.

Standard Sequencing Primer Kits

All of the most commonly used sequencing primers are also available labeled with either biotin, fluorescein, or all four ABI dyes. These off-the-shelf primers are shipped within 24 hours. And they're supplied purified, dried, and ready-to-use — enough for over 100 sequencing reactions.

For more information, fax back or return the coupon below. Or call 1-800-234-5362.

NOW YOU CAN ORDER WHOLE GENES AS EASILY AS PROBES AND PRIMERS

We'll work with you to design a coding sequence to suit your application. We'll clone it into the vector of your choice, verify the sequence, and deliver a guaranteed product, complete with restriction maps and QC autoradiographs.

ORDER NOW AND GET A FREE STANDARD PRIMER OR PROBE!

Return or fax this coupon for details, or call

1-800-234-5362

- ☐ Free standard primer or probe ☐ Custom oligo synthesis ☐ Non-isotopic labeled oligos
☐ Custom gene synthesis ☐ Arbitrary 10-mer primers ☐ Non-isotopic sequencing primer kits

Name _____

Company/Institution _____

Address _____

City _____ State _____ Zip _____

Phone _____ Fax _____

GENOSYS

Genosys Biotechnologies, Inc.

1442 Lake Front Circle, Suite 185, The Woodlands, TX 77380-3600
Phone: (713) 363-3693 Fax: (713) 363-2212

Copyright © 1994 Geosys Biotechnologies, Inc.

Circle No. 8 on Readers' Service Card

Molecular Biology Certified Media

Premixed, Premeasured Capsules

Save Time, Enhance Media Consistency.

BIO 101 designs Yeast and Bacteriological growth media specifically for recombinant DNA applications. Unlike microbiology growth media manufacturers, we focus on plasmids, phagemids, and phages during growth.

We study amplification, fidelity, effect of temperature and media components, and the optimal culturing conditions of the plasmid and the host.

Media is capsulated for convenience, longer shelf-life and guaranteed performance. No measuring, no mess, and no pH adjustment. Simply add them to water and autoclave.

Come On In, Molly.
Let's Autoclave!

We now have optimized powdered formulations for *Pichia pastoris* and Dual Hybrid Systems.

For more information, call BIO 101 or your nearest distributor.



"Check with BIO 101...
We have your time in mind!"

1070 Joshua Way, Vista, CA 92083
619-598-7299 • Fax: 619-598-0116

Internet Mail:
BIO101@BIOTECHNET.COM

800-424-6101


PARTIAL LIST OF CAPSULATED PRODUCTS:

Circlegrow®
L-Broth
L-Broth Agar
L-Top Agar
LB-Medium
LB-Agar Medium
LB-Top Agar
Tryptone Broth
Tryptone Agar
Tryptone Top Agar
NZYCM Broth
NZYCM Agar
NZYCM Top Agar
NZYM Broth
NZYM Agar
NZYM Top Agar
NZM Broth
NZM Agar
NZM Top Agar
YT-Broth
YT-Broth Agar
YT-Top Agar
2xYT-Broth
2xYT-Broth Agar
2xYT-Top Agar
SOB
SOB Agar
"Super Broth"
Casamino Acids
Yeast Extract-Y
YPD (YEPD) Broth
YPD (YEPD) Agar
Mock YEPD Broth
Mock YEPD Agar
YEP Broth
YEP Agar
YE
YE Agar
YES
YES Agar
ME
ME Agar

Media for
P. pastoris
Dual Hybrid Systems
Now available
Call 800-424-6101


International Distributors:

• **Australia** • Integrated Sciences Tel: (2) 436-2611, Fax: (2) 436-4442 • **Austria** • Bio-Trade Tel: (1) 889 18 19, Fax: (1) 889 18 19 20 • **Benelux** • Westburg B.V. Tel: (33) 950094, Fax: (33) 951222 • **Canada** • Bio/Can Scientific Tel: (905) 828-2455, Fax: (905) 828-9422 • **France** • Ozyme Tel: (1) 34 60 24 24, Fax: (1) 34 60 92 12 • **Germany** • Dianova G.m.b.H. Tel: (40) 32 30 74, Fax: (40) 32 21 90 • **Israel** • Tamar, Ltd. Tel: (2) 520 279, Fax: (2) 527 318 • **Italy** • Stebbio s.r.l. Tel: (51) 63 43 340, Fax: (51) 63 43 332 Dasit S.p.A. Tel: (2) 9027381, Fax: (2) 90361209 • **Japan** • Funakoshi, Co. Ltd. Tel: (3) 5684-1622, Fax: (3) 5684-1633 • **Korea** • Bio-Medical & Science Tel: (2) 569 6902, Fax: (2) 553 9670 • **Mexico** • Bioelec Tel: (525) 341 7764, Fax: (525) 341 7403 • **Singapore** • TIBS Trading PTE Ltd. Tel: (65) 2929783, Fax: (65) 2932517 • **Scandinavia** • Kebo Lab Tel: (8) 621 34 00, Fax: (8) 621 35 68 • **Spain** • Ramon Cornet SA Tel: 237 5562, Fax: 415 3466 • **Switzerland** • Lucerna Chem AG Tel: (41) 36 96 36, Fax: (41) 36 96 56 • **Taiwan** • Cold Spring Biotechnology Co. Ltd. Tel: (2) 694 0066, Fax: (2) 694 3204 • **United Kingdom** • Stratech Scientific, Ltd. Tel: (582) 481884, Fax: (582) 481895



Smile! Renaissance™ non-rad DNA labeling kits give you reproducible results, not high backgrounds.

Are you repeating experiments just to reduce backgrounds? Then look into Renaissance™ non-radioactive DNA labeling and detection products from DuPont NEN®. And get low backgrounds and reproducible results the first time, and every time.



- Sensitive HRP-luminol systems for colony plaque lifts and Northern, Southern, and Western blots.
- Results in minutes.
- Guaranteed one-year shelf life.
- Backed by full protocol and comprehensive technical service.
- Now available: Random Primer Fluorescein dUTP Labeling Kit and Oligonucleotide 3' End Labeling kit (Fluorescein ddUTP).


DuPont gives you a choice of radioactive and now Renaissance non-radioactive labeling and detection products. For information, or orders call 1-800-551-2121. For information by fax, call 1-800-666-6527 and request number 9002.

All DuPont products are manufactured under an ISO 9001 quality system registered by UL and approved by BSI.

United States 1-800-551-2121 • Canada 1-800-387-8391 • Australia +61 (008) 257149 • Belgium +32 (02) 7242717
Denmark +45 31506610 • France (01) 6982 5450 • Germany +49-6172-872600 • Italy (02) 25302 481/483
Japan +81 (03) 5421-1354 • Latin America/Asia Pacific Fax +1 (617) 426-2464 • Netherlands +31 (073) 206550
Sweden +46 (08) 7503700 • Switzerland +41 (01) 8410330 • United Kingdom +44 (0438) 724037



DuPont NEN®



Our UNIX Solutions Are Winning Solutions.

"We selected Digital's Alpha AXP™ systems as our initial development platform for our next generation of computer-aided emergency dispatching decision support software. The Alpha AXP systems provide BAPSS with the best system price/performance. This, combined with a 64-bit hardware architecture, fast graphics performance and an industry-leading standards-compliant UNIX® operating system, which resulted in our decision to employ Digital's technology."

*David Riffelmacher, President
Bell Atlantic Public Sector Systems*

"The system installation went very smoothly – no glitches. Overall, we would say that Digital has done an outstanding job in their 64-bit OSF/1® implementation."

*Jeffrey S. Jalbert, President
JCC Consulting, Inc.*

"Digital's unified UNIX, DEC OSF/1, is emerging as a strong contender in the commercial UNIX market. Oracle™ 7 for Alpha AXP DEC OSF/1 in its very first quarter resulted in revenue significantly higher than Oracle expects for a new product line. Our results demonstrate the high quality of Oracle 7 and the competitiveness of DEC OSF/1 in the commercial UNIX market."

*Bonnie Crater, Director, Product Management
Midrange Products Division, Oracle Corporation*

"Over the past three years, Digital has invested a considerable amount of money and development effort in UNIX. And today, the work they have done has put them in a leadership position within the UNIX Operating System Market."

*Judith Hurwitz, President
Hurwitz Consulting Group, Inc.*

"We have ported Sybase's Gain Momentum® product to DEC OSF/1 V1.3 and run extensive testing. As a result, we are confident that DEC OSF/1 is of such high quality and performance that it will prove an excellent addition to the world of commercial, mission-critical UNIX computing."

*Noel Leca, Vice President,
Business Unit Manager for the Gain Multimedia
Products Group, Sybase, Inc.*

For more information call our network of Authorized Resellers or Digital's Alpha AXP Fast FAX at 1-800-842-7027 in the U.S. and Canada.

Putting Technology To Work

digital

© Digital Equipment Corporation 1994. DIGITAL, the DIGITAL logo, DEC, the AXP logo and AXP are trademarks of Digital Equipment Corporation. UNIX is a registered trademark of UNIX System Laboratories, Inc., a wholly owned subsidiary of Novell, Inc. OSF and OSF/1 are registered trademarks and Open Software Foundation is a trademark of the Open Software Foundation. Oracle is a trademark of the Oracle Corporation. Gain Momentum is a registered trademark of Sybase, Inc.

Circle No. 68 on Readers' Service Card

Half-Day DNA. And That's Just The Half Of It.



Synthesizing primers, probes, linkers, or gene fragments? Cut your time — and your reagent use — in half with the Millipore Expedite™ Nucleic Acid Synthesis System.

For instance, a pair of 20-mer PCR primers can be synthesized in about an hour and worked up in about 15 minutes, so you can start using your new DNA before noon.

How's that possible? Because no matter how many columns you're running — 1, 2, 3, or 4 — the total cycle time is under 4 minutes. The Expedite system's exclusive alternating phase synthesis allows columns to act independently, so each is always synthesizing.

Expedite reagents speed things along further. Oligonucleotides made with Expedite brand

chemistry can be deprotected in only 15 minutes at 55°C or 2 hours at room temperature.

But speed is just half the story. The Expedite system's patented microfluidics plate provides for extremely low internal volume. Reagent delivery to the column is cut by as much as 50% over other systems. So, you get more couplings per chemistry changeout; more DNA synthesized per dollar.


Now you've got the whole picture. Just give us a call and we'll send you more information, expeditiously. In the U.S. and Canada 1-800-872-0071.

In Japan (03) 3471-8191 and Europe (44) 0923 816375.




MILLIPORE

If you're working with glycoproteins or carbohydrates, the FACE® Imager should be working for you.



PHOTOGRAPH OF FACE GEL SHOWING VISIBLE TRACKING DYES



PHOTOGRAPH OF COOMASSIE STAINED PROTEIN GEL

IMAGE OF FACE GEL SHOWING FLUORESCENT OLIGOSACCHARIDES RELEASED FROM GLYCOPROTEINS






IMAGE OF COOMASSIE STAINED PROTEIN GEL




Say "hello" to the fastest, most reliable, and cost-effective way available to image all your minigels. The FACE (Fluorophore-Assisted Carbohydrate Electrophoresis) Imager allows you to analyze your oligosaccharides and monosaccharides on specially formulated pre-cast gels without radiolabelling, staining, or exhaustive sample preparation.

And, say "good-bye" to guesswork and lab books filled with cracked protein gels. The FACE Imager can analyze your Coomassie Blue and silver-stained protein gels, allowing you to determine relative quantities and molecular weights in seconds. You can conveniently use the FACE Imager for ethidium bromide DNA minigels.

It's easy to put the FACE Imager to work for you. Just call Glyko at 415 382 6653, fax 415 382 7889, or call 800 33 GLYKO (334 5956) toll-free.



Bachem California Offers the Absolute Purest Foundation for Your Research



The quality of the materials you use is the cornerstone upon which your research is built. Inferior or sub-standard raw materials add unexpected variables to your work, costing you valuable time and money. At Bachem California, we understand the vital importance of providing you with products of consistent quality and reliability. That's why for over 23 years we have gained a well-deserved reputation for peptides, growth factors, and protected amino acids that set the standard in the industry. Our commitment to quality is your guarantee of reliability.

Only Bachem California combines uncompromising commitment to quality with the most advanced peptide technology and 23 years' experience in the peptide field. This enables us to support every aspect of your project from basic research to clinical applications and commercial development.

No one can match our credentials.

- Complete, State-of-the-Art Selection of Peptides and Related Products
- Proven Track Record in GMP Production and Regulatory Support
- Unmatched Scale-Up Capabilities
- Over 30 Drug Master Files with the FDA
- Technology to Succeed on the Most Complex Sequences
- Success in Supporting Hundreds of IND and NDA Applications
- Confidential, Quality Custom Synthesis

Bachem California means uncompromising quality, purity, reliability, and service—the ultimate foundation for your research.

Call us at 310-539-4171 today to find out how we can support your work.



Bachem (UK), Ltd.
A subsidiary of Bachem California
69 High Street, Saffron Walden
Essex CB10 1AA, England
Tel: (0799) 526465 • Telefax: (0799) 526351

Bachem California
3132 Kashiwa Street
Torrance, CA 90505
Tel: (310) 539-4171
Telefax: (310) 530-1571

Saxon Biochemicals GmbH
A Subsidiary of Bachem California
Feodor-Lynen-Strasse 14
D-30625 Hannover 61, Germany
Tel: (0511) 546080 • Telefax: (0511) 5460888

Circle No. 65 on Readers' Service Card

THE NEW BIOMEK® 2000 AUTOMATES LIFE SCIENCE... YOUR WAY

Automate virtually any liquid-handling procedure *the way you want* with the Biomek® 2000 Laboratory Automation Workstation from Beckman. Choose the labware, tools and accessories you require to perform your research with the powerful and easy-to-use Windows®-based BioWorks™ operating system.

EXTENSIVE APPLICATIONS

Years of reliability, versatility and innovation have made Biomek the benchmark system for BioRobotics in Pharmaceutical, Human Genome and Bioresearch laboratories worldwide. This new generation Workstation provides the productivity and throughput you need today for:

- *Immunology* — ELISA and RIA
- *DNA* — Sequencing, Mapping and PCR** Processes
- *Drug Discovery* — Receptor Binding and Cytotoxicity Assays
- *Biochemistry* — Total Protein and Enzymatic Assays
- *Molecular Diagnostics* — HLA Typing and Prenatal Screening
- *Quality Control* — Endotoxin Testing and Liquid Extractions
- *Tissue Culture* — Cell-based and Sterile Techniques
- *Clinical Research* — HIV and Hepatitis Screening


INTELLIGENT PERFORMANCE

Intuitive icon-based BioWorks software intelligently automates your protocols in minutes. Advanced workstation features and fast motion translate operating system commands into efficient performance. The Biomek easily integrates other software programs and devices into a complete automation solution that includes temperature control, detection and quantitation.

For greatest flexibility in automating repetitive tasks, up to 12 user-defined positions are available on the large Biomek worksurface. The unique interchangeable liquid-handling tools deliver accurate volumes from 1.0 µL to 150 mL with exceptional precision. And the modular design allows you to add options as you need them, such as the high-capacity Side Loader system for extended periods of unattended operation.

Find out how you can Automate Life Science...Your Way by contacting your local Beckman office for *free* technical information about the new Biomek 2000 BioRobotics System.

* Windows is a trademark of Microsoft Corporation.
** PCR is covered by patents owned by Hoffmann-La Roche Inc.



BECKMAN

Worldwide Offices: Africa, Middle East, Eastern Europe (Switzerland) (22) 994 07 07. Australia (61) 02 816-5288. Austria (2243) 85656-0. Canada (800) 387-6799. China (861) 5051241-2. France (33) 1 43 01 70 00. Germany (49) 89-38871. Hong Kong (852) 814 7431. Italy (39) 2-953921. Japan 3-3221-5831. Mexico 525 575 5200, 525 575 3511. Netherlands 02979-85651. Poland 408822, 408833. Singapore (65) 339 3633. South Africa (27) 11-805-2014/5. Spain (1) 358-0051. Sweden (8) 98-5320. Switzerland (22) 994 07 07. Taiwan (886) 02 378-3456. U.K. (0494) 441181. U.S.A. 1-800-742-2345.

©1994 Beckman Instruments, Inc.

Circle No. 50 on Readers' Service Card




With a couple of cards, the Genesys 5 spectrophotometer can do a lot more tricks.

You don't have to be a performer to know high-performance. The new Milton Roy (high-performance) SPECTRONIC GENESYS 5 SPECTROPHOTOMETER with optional programmable SoftCards combines the artistry of quality and accuracy with time-saving and space-saving innovations. It comes with a 5 nm slit width • a graphic LCD screen • a built-in automatic 8-position cell holder • RS-232-C interface • optional built-in printer • linear and non-linear standard curves • scanning • multiple wavelength tests. **For under \$6,000...this could be a tough act to follow.**



For our free video, call 1-800-654-9955

 **MILTON ROY**



The ^{32}P Kinase Assay is about to dye.

Introducing a Non-Radioactive PKA assay that's sensitive to 0.02 units of enzyme activity—the Pierce SpinZyme™ Non-Radioactive PKA Assay.

Sensitive. The Pierce SpinZyme™ Non-Radioactive PKA Assay measures as little as 0.02 units of PKA activity.

Non-Radioactive. No more tedious wash procedures . . . no more radioactive waste disposal . . . no more waiting for the scintillation counter to see your results. The SpinZyme™ Non-Radioactive PKA Assay gives you the sensitivity you need to measure PKA activity . . . without the hassles and hazards associated with ^{32}P .

Specific. Unlike the ^{32}P assay, the SpinZyme™ Non-Radioactive PKA Assay

detects only the phosphorylation of a dye-labeled peptide substrate. There's no interference from the phosphorylation of endogenous protein in physiological samples.

Fast. In as little as 20 minutes, you can perform a complete spin-column separation—including all binding and elution steps. The Pierce SpinZyme™ separation unit incorporates a unique affinity membrane, which quickly and efficiently separates the phosphorylated dye-labeled peptide from the non-phosphorylated peptide.

Wash and elution fractions can be measured spectrophotometrically or fluorometrically. Results can be measured both qualitatively and quantitatively.

 **PIERCE**


a Perstorp Biotec Company



3747 N. Meridian Rd.
P.O. Box 117
Rockford, IL 61105 U.S.A.
Tel: 800.874.3723 or 815.968.0747

Pierce Europe B.V.
P.O. Box 1512
3260 BA Oud Beijerland
The Netherlands
Tel: 31.1860.19277

Circle No. 60 on Readers' Service Card



I

Leukosorb™. A New Filtration Medium.

Leukosorb is a new filtration medium from Pall which is designed to immobilize white cells from whole blood or packed red cells. The process eliminates all tedious centrifugation steps, and can be completed in less than one hour.

I-Blood Volumes As Little As 50 μ l.

A wide range of blood volume can easily be processed from as little as 50 μ l to as much as 50ml.



II

The blood sample is placed in a syringe. The white cells are adsorbed by the Leukosorb filter as the blood passes through the fibrous matrix (Figure I).

III-Ready For DNA Analysis.

The process is now completed. The resulting DNA

lysates can be used directly for hybridization or for PCR amplification (Figure III). Total time from start to finish is less than one hour. For a free sample pack, contact Pall at one of the locations listed below.

Extract DNA In Under An Hour Without A Centrifuge.



After adsorption of the white cells, the filter medium is washed in a hypotonic saline solution. After the wash, no red cells remain on the Leukosorb medium.

II-Extract DNA.

Now, the Leukosorb filter disks are placed in a microfuge tube with proteinase enzyme to free the DNA from the immobilized white blood cells.



III


USA
(516)
484-5400

United Kingdom
(0705) 75-3545

Japan
(03) 3495-8343



Pall BioSupport Division
East Hills, New York 11548



For the ^{32}P Kinase Assay that won't dye...

Try the fastest, easiest radioactive assay available—the Pierce SpinZyme™ PKA Assay Kit, featuring our exclusive Basic Separation Units.

Fast. In less than 5 minutes, the SpinZyme™ Basic Separation Units will separate excess [γ - ^{32}P]ATP from the phosphorylated peptide. These exclusive Pierce units feature a phosphocellulose membrane in a spin column format, which eliminates the tedious wash steps of the standard procedure.

Low Background. Researchers who have tested these separation units report significantly lower background levels compared with standard assay methods using phosphocellulose squares.

Convenient Handling, Easy Disposal. The sample buckets containing the phosphocellulose membrane and bound phosphorylated peptide

are easily transferred to scintillation vials for counting. Radioactive liquid waste is contained in microcentrifuge tubes for easy disposal.

Reduced Radioactive Waste. Only 1 ml of radioactive liquid waste is generated per sample.

The SpinZyme™ PKA Assay Kit (Product #29530) contains sufficient materials to perform 100 assays including: Reaction Buffer, Activator Solution, Kemptide (PKA Substrate), PKA Dilution Buffer, 100 SpinZyme™ Basic Separation Units and 100 Additional Receptacles.

The SpinZyme™ Basic Separation Units (Product #29520) includes 100 SpinZyme™

Basic Separation Units and 100 Additional Receptacles. Suitable for use in kinase assay protocols using peptide substrates that will bind to phosphocellulose.

 **PIERCE**

a Perstorp Biotec Company

3747 N. Meridian Rd.
P.O. Box 117
Rockford, IL 61105 U.S.A.
Tel: 800.874.3723 or 815.968.0747

Pierce Europe B.V.
P.O. Box 1512
3260 BA Oud Beijerland
The Netherlands
Tel: 31.1860.19277



PERSONAL BIBLIOGRAPHIC DATABASES...

These cost more:

ENDNOTE • DMS 4 CITE • PRO-CITE
REFERENCE MANAGER • REF-11

This does more:

PAPYRUS™

Version 7!

- Manages up to 2 million reference citations. Stores up to 16,000 characters and 100 keywords per reference.
- Dozens of predefined output formats, plus the ability to easily design your own.
- 100% compatible with WordPerfect*, Microsoft Word*, Ami Pro*, WordStar, XyWrite, Signature, ChiWriter, TeX. *including Windows™ versions
- Can also be used with virtually all other word processors.
- Fast, powerful search capabilities.
- Able to import references from national databases, CD-ROM files, monthly diskette services, or almost any other database or text file.
- Allows an unlimited number of Notecards for each reference.
- Powerful new user interface.
- Fully compatible with Windows™.

for IBM-PC and compatibles
also available for VAX-VMS
Macintosh version under development

Research
SOFTWARE DESIGN

2718 SW Kelly Street, Suite 181
Portland, OR 97201
(503) 796-1368 FAX: 503-241-4260

Complete System \$99


Full money-back guarantee
on purchase of Complete System.

Demo System \$25

Demo price credited toward subsequent Complete System purchase.

Outside North America, add \$20 shipping charge - U.S. funds, on a U.S. bank.

Circle No. 62 on Readers' Service Card



New Grant?

...Then get the latest about what's new in physiological research instrumentation

- New** automated syringe pumps
- New** gas anesthesia
- New** micromanipulators
- New** stereotaxic instrumentation
- New** MacLab™ data acquisition

Call (708) 860-9700 or send today for your free copy of Stoelting's 80 page, 1993 catalog for physiological and pharmacological research.



Stoelting Co., 620 Wheat Lane
Wood Dale, IL 60191 USA
(708) 860-9700

Circle No. 4 on Readers' Service Card

FREE REFERENCE POSTER ON CELL ADHESION

Make your next experiment in inflammation or tumor metastasis a success.

This new teaching and reference tool includes a current listing of cell adhesion molecules, ligands and molecular characteristics. You will also find a bibliography as well as a synopsis of practical considerations for optimizing your experiments. Contact Becton Dickinson for your free copy today. 1-800-952-3222 (U.S.)

Becton Dickinson
Cellular Imaging Systems
Ordering (800) 223-8226
Customer Support (800) 952-3222
Fax (408) 954-2156

Becton Dickinson Canada, Inc.
(905) 822-4820
Fax (905) 855-1243

Becton Dickinson European HQ
Erembodegem-Aalst, Belgium
(32) 53-720211
Fax (32) 53-720450

Nippon Becton Dickinson Co., Ltd.
Tokyo, Japan
(81) 3-3403-9991
Fax (81) 3-3403-5008

Becton Dickinson Worldwide Inc.
Singapore
(65) 861-0633
Fax (65) 860-1590

CAMFolio is a trademark of
Becton Dickinson and Company.
For research use only. Not for use in
diagnostic or therapeutic procedures.

**BECTON
DICKINSON**

Circle No. 25 on Readers' Service Card

CELL ADHESION RESEARCH: EXPERIMENTAL CONSIDERATIONS

CAMFOLIO
TOOLS FOR EXPLORING CELLULAR INTERACTIONS AT THE MOLECULAR LEVEL

Cell Type	Adhesion Molecule	Characteristics	Experimental Considerations	References
Epithelial	E-cadherin	120 kDa, transmembrane, Ca ²⁺ -dependent	Use anti-E-cadherin antibody for Western blotting	Greenhalgh et al., 1991
Endothelial	ICAM-1	80 kDa, transmembrane, Ca ²⁺ -independent	Use anti-ICAM-1 antibody for immunoprecipitation	Springer et al., 1990
Leukocyte	L-selectin	140 kDa, transmembrane, Ca ²⁺ -dependent	Use anti-L-selectin antibody for flow cytometry	Springer et al., 1990
Neuronal	Neural cell adhesion molecule (NCAM)	140 kDa, transmembrane, Ca ²⁺ -independent	Use anti-NCAM antibody for immunoprecipitation	Springer et al., 1990
Epithelial	Integrin α5β1	150 kDa, transmembrane, Ca ²⁺ -independent	Use anti-α5β1 antibody for immunoprecipitation	Springer et al., 1990
Endothelial	Integrin αvβ3	150 kDa, transmembrane, Ca ²⁺ -independent	Use anti-αvβ3 antibody for immunoprecipitation	Springer et al., 1990
Leukocyte	Integrin αLβ2	150 kDa, transmembrane, Ca ²⁺ -independent	Use anti-αLβ2 antibody for immunoprecipitation	Springer et al., 1990
Neuronal	Integrin α3β1	150 kDa, transmembrane, Ca ²⁺ -independent	Use anti-α3β1 antibody for immunoprecipitation	Springer et al., 1990

Commonly Used Techniques

- Flow Cytometry:** Measure cell surface expression of adhesion molecules.
- Immunoprecipitation:** Isolate adhesion molecules from cell lysates.
- Western Blotting:** Detect specific adhesion molecules in cell lysates.
- Immunofluorescence:** Visualize adhesion molecules on cell surfaces.
- Cell Adhesion Assays:** Measure cell-cell and cell-matrix interactions.

Cell Adhesion Assays:

- Adhesion to Immobilized Ligands:** Measure cell attachment to ligands immobilized on a solid support.
- Adhesion to Immobilized Cells:** Measure cell attachment to cells immobilized on a solid support.
- Adhesion to Immobilized Extracellular Matrix (ECM):** Measure cell attachment to ECM components immobilized on a solid support.

POWER PCR

from

TaKaRa

TaKaRa PCR *in vitro* Mutagenesis Kit.....

Provides a Simple and Easy Way to Introduce Site-Specific Mutations into DNA

Cloned in either pUC18/118, pUC19/119, or their Derivatives.

* TaKaRa is the first Licensee possessing the PCR Licence in these various fields; Research, Forensics and Human Identity, Environmental Testing, Plant Diagnostics, Animal Diagnostics, Animal Identity Testing and Positive Trait Breeding, and Food Testing.

Code No. R016; for 20 Preparations

TaKaRa PCR *in vitro* Mutagenesis Kit


Features:

An Assortment of Six Primers for Mutagenesis: Provides a wide selection of pre-made mutagenesis primers for introducing sequence alterations into DNA cloned at any position in a vector, as well as an operational convenience that reduces user preparation of mutant induction primers from two to one for each mutagenesis.

Generalized Mutagenesis via PCR: Facilitates the introduction of mutations to such an extent that multiple mutations can be easily and rapidly achieved, which may enhance systematic studies on the structure-function relationship of enzymes, proteins and so on.

Automatic Mutant Selection: Eliminates *in vivo* selection for a mutant strand through transformation. In this procedure only the mutant strand is cut for recloning, and clearly distinguished from the non-mutant strand. Mutagenic primers are designed to eliminate a restriction enzyme site by causing a single-base change in one of the following sites; *EcoR* I, *SacI*, *BamH* I, *Xba* I, *Sal* I, *Acc* I, *Hinc* II, *Sph* I, or *Hind* III.

A Control Set Included in the Kit: Composed of Control Template, Control Primer, and Restriction Enzyme *EcoR* I, verifies and assures kit performance.



POWER PCR products from TaKaRa are sold under licensing arrangements with Roche Molecular Systems and F. Hoffman -La Roche Ltd. and The Perkin-Elmer Corporation. Purchase of these products is accompanied by a license to use them in the Polymerase Chain Reaction (PCR) process in conjunction with an Authorized Thermal Cycler.



TAKARA SHUZO CO., LTD.

Biomedical Group
Otsu, Shiga, Japan
Tel: +81 775-43-7247
Fax: +81 775-43-9254

DISTRIBUTORS for More Information and Orders:

PanVera Corporation (USA and South America)
Tel: 608-233-5050 Fax: 608-233-3007
Toll Free: 1-(800) 791-1400

R&D Systems Europe Ltd. (UK)
Tel: (0235) 529449 Fax: (0235) 533420

ITC Biotechnology GmbH (Germany)
Tel: (06221) 303907 Fax: (06221) 303511

CIS bio international (France)

Tel: (1) 69 85 73 90 Fax: (1) 69 85 73 65

CIS Diagnostici S.p.A. (Italy)


Tel: (0161) 912292 Fax: (0161) 912308

Cheng Chin Trading Co., Ltd. (Taiwan)

Tel: (02) 331-3111 Fax: (02) 382-2197

Protech Technology Enterprise Co., Ltd. (Taiwan)

Tel: (02) 381-0844 Fax: (02) 311-8524



MOLECULAR RESEARCH CENTER, INC.
Introduces a splendid isolation for life science and clinical research

TRI REAGENT™

*The most capable reagent in molecular biology
simultaneously isolates total RNA, poly A⁺ RNA, DNA and Proteins*.*

The simultaneous isolation of RNA, DNA and proteins from the same sample by TRI REAGENT™ is the newest version of the single-step liquid-phase separation method (P. Chomczynski, *Biotechniques*, Vol. 15 No. 3 1993). This highly reliable method performs well with small and large quantities of tissues or cultured cells, and allows simultaneous processing of a large number of samples of human, animal, plant, yeast, bacterial and viral origin. For the simultaneous isolation of RNA, DNA and proteins from liquid samples such as whole blood, serum and CSF we offer TRI REAGENT™ LS.

- **Total RNA** is isolated in 1 hour and is ready for Northern (See Figure) and dot blotting, RT-PCR (Hoffmann LaRoche Corp. patent), RNase protection assay and molecular cloning. 1ml of the reagent isolates 50-200 µg of total RNA.
- **Poly A⁺ RNA** is isolated in 1.5 hours by application of the aqueous phase containing total RNA on the oligo dT-cellulose column.
- **High molecular weight DNA** is isolated in 3 hours and is ready for restriction analysis, Southern (See Figure) and dot blotting, PCR and molecular cloning.
- **Proteins** are isolated in 3 hours and are ready for Western blotting (See Figure).

Introductory discount of 20% for first time orders (U.S. only): 50ml TRI REAGENT™ (cat. no. TR-118-IS), 50 ml TRI REAGENT™ LS (cat. no. TS-120-IS), 5 reusable oligo dT-cellulose columns (cat. no. OT-125). We guarantee that TRI REAGENT™ will outperform RNA, poly A⁺ RNA, and DNA isolation methods currently used in your laboratory. For more information or to place an order call 1-800-462-9868, 513-841-0900, or contact us by Fax 513-841-0080.

* Patent pending

North America

MRC

MOLECULAR RESEARCH CENTER, INC.
5645 Montgomery Rd.
Cincinnati, OH 45212

United Kingdom

MOLECULAR RESEARCH CENTRE OXFORD
Oxford, U.K.
PHONE: 44-865-351-511
FAX: 44-865-351-511

Italy

BIO-OPTICA MILANO SpA
Milano, Italy
PHONE: 39-2-26-40-274
FAX: 39-2-21-53-000

Israel

TAL-RON
Rehovot, Israel
PHONE: 972-8-472563
FAX: 972-8-471156


Taiwan

PAN ASIA BIOMEDICAL TECHNOLOGY, INC.
Taipei, Taiwan
PHONE: 886-2-7418169
FAX: 886-2-7764372

Northern blot of rat
growth hormone
(rGH) mRNA (1 kb)



Southern blot of EcoR1
rGH gene fragment
(11.4 kb)




Western blot of rGH
(22 kD)



Circle No. 18 on Readers' Service Card

DNA SEQUENCING



**Sequencing strategies
custom-tailored for your
projects and budget.**

- Priced from 99¢/base to \$3.20/bp
- Manual sequencing (Sanger Method)
- Rapid turn-around time
- Sequences guaranteed
- Complete confidentiality assured

Services to support your Research Projects:

- Sub-Cloning
- In-vitro Mutagenesis
- Phage & Plasmid DNA Production
- Probe Labeling
- Oligonucleotide Synthesis




Lofstrand Labs Limited

"Highest quality sequencing at the lowest price"


Circle No. 45 on Readers' Service Card

Call
1-800-5410-DNA
to discuss how
we can work as an extension
of your research laboratory.
(Outside U.S. call:
+1-301-330-0111)

FINDING THE PERFECT PROTEIN PARTNER IS EASIER THAN YOU THINK.



If you're looking to positively identify interacting proteins, *and* immediately access their cloned genes, search no further. With MATCHMAKER™ (#K1605-I), the only commercially-available kit employing the innovative two-hybrid system, you'll dramatically accelerate the process of discovery by tracing protein interaction pathways *in vivo*. This direct, sensitive assay is simple, yet extraordinarily powerful—CLONTECH also offers premade MATCHMAKER cDNA and genomic libraries which are cloned into GAL4 activation domain vectors. So, if you're exploring cell cycle control, cell differentiation, signal transduction, or any complex



Results of a MATCHMAKER Two-Hybrid assay for interaction of the murine p53 protein with SV40 large T-antigen. The two candidate proteins were cloned into a DNA-binding domain vector and an activation domain vector, respectively. Transformants were assayed for β -galactosidase expression. Rows 1-3: negative controls. Row 4: test interaction.


protein interaction, MATCHMAKER will help you find exact protein partners faster and easier than you've ever thought possible. Just call 1-800-662-CLON or contact your local distributor.

DON'T EXPERIMENT WITH ANYTHING ELSE. **CLONTECH**

4030 Fabian Way, Palo Alto, California 94303 USA • Fax: 800/424-1350 415-424-1064 • Phone: 800/662-2566 415/424-8222
©1994 CLONTECH Laboratories, Inc. MATCHMAKER is a trademark of CLONTECH Laboratories, Inc.

Circle No. 37 on Readers' Service Card

ANNOUNCING MAJOR DEVELOPMENTS IN 1994!



CBP COMPARATIVE BIOCHEMISTRY AND PHYSIOLOGY, Parts A, B, & C

An International Journal

An Official Journal of the:

Japanese Society for Comparative Physiology & Biochemistry
European Society for Comparative Physiology & Biochemistry
Canadian Society for Zoology (CPB Section)

Editors-in-Chief:

Peter W. Hochachka, University of British Columbia, Canada
and Thomas P. Mommsen, University of Victoria, Canada

In order to better serve its readership, **CBP** is undergoing a major transformation. **CBP** now features new editors, new associate editors, new aims and scope, redefined sections and new society affiliations.

NEW AIMS AND SCOPE

- To publish original articles emphasizing comparative aspects of biology.
- Adaptation and evolution as organizing principles are encouraged, as well as comparative and mammalian studies.


REDEFINED SECTIONS:

A: Physiology — scope has been expanded from the classical molecular, cellular & organ-specific physiology to include *ecological physiology*.

B: Biochemistry and Molecular Biology — the expression of functional mechanisms and evolutionary principles at the molecular level.

C: Pharmacology, Toxicology and Endocrinology — concerned with drug action at different levels of cellular organization including *endocrinology* and message transduction. A molecular approach to these fields is encouraged but not a prerequisite.

CALL FOR PAPERS!



CELLULAR AND MOLECULAR BIOLOGY RESEARCH

Editor-in-Chief: M.A.Q. Siddiqui, State University of
New York Health Science Center at Brooklyn, New York

Cellular and Molecular Biology Research publishes original articles emphasizing the wide range of topics concerned with cellular and molecular mechanisms involved in cell growth, development and differentiation. Among those areas of coverage are *in vitro* and *in vivo* studies of structure and function of biologically important macromolecules in eukaryotic organisms, gene regulatory mechanisms, transcription, structure and function of regulatory factors, carcinogenesis, oncogenes and tumor suppressor genes in growth and differentiation, signal transduction mechanisms, structure and function of cell organelles, control of growth and development, cell-cell interaction, embryogenesis, and structure and function of viruses and bacteriophages. Papers with emphasis on recombinant DNA technology, transgenes, and genetic manipulation of cells and organs are of special interest.

CALL FOR PAPERS

The journal considers full-length, original articles, short communications, summary articles, techniques papers, and commentaries.

The Editorial Board of **CMBR** is committed to objective peerreview and rapid publication. Manuscripts are published in the shortest time possible commensurate with high scientific quality. Please complete the form below to receive complete guidelines for electronic submission and manuscript preparation.

Free Sample Copies Available Upon Request!

Please send me: () Sample Copy () Instructions for Contributors

Name

Institution/Department

Address

City/State/Zip

A41595

CBP COMPARATIVE BIOCHEMISTRY AND PHYSIOLOGY

() Part A () Part B () Part C

ISSN: 2090-0000 Vol. 107-109, 1994 33 issues per annum

() CELLULAR AND MOLECULAR BIOLOGY RESEARCH

ISSN: 0968-8773 Vol. 40, 1994 8 issues per annum



Return to: **PERGAMON** An imprint of Elsevier Science

US: 660 White Plains Road, Tarrytown, NY 10591-5153, USA
Fax: 914-333-2444

UK: The Boulevard, Langford Lane, Kidlington, Oxford OX5 1GB
Fax: 011-44-865-843952