

AMERICAN  
ASSOCIATION FOR THE  
ADVANCEMENT OF  
SCIENCE


# SCIENCE

2 SEPTEMBER 1994  
VOL. 265 • PAGES 1333-1496

\$6.00

VACCINES

FRONTIERS IN  
MEDICINE





# hNT Neurons

## STRATAGENE



## Terminally Differentiated Postmitotic Human Neurons\*

**Internet Mail**  
tech\_services@stratagene.com

USA

**Corporate Headquarters**  
(800) 424-5444  
(619) 535-5400  
Fax: (619) 535-0045

GERMANY

**Stratagene GmbH**  
Telephone: (06221) 400634  
Telefax: (06221) 400639

UNITED KINGDOM

**Stratagene Ltd.**  
Telephone: (0223) 420955  
Telefax: (0223) 420234

SWITZERLAND

**Stratagene GmbH**  
Telephone: (01) 3641106  
Telefax: (01) 3657707

- elaborate processes typical of mature postmitotic human neurons
- express many neuronal markers<sup>1</sup>
- express functional NMDA and non-NMDA glutamate receptor channels<sup>2</sup>

1. Pleasure, S.J., and Lee, V.M.Y. (1993) *J. Neurosci. Res.* 35:585-602

2. Pleasure, S.J., et al. (1992) *J. Neuroscience* 12: 1802-1815

- express and secrete amyloid precursor peptides<sup>3</sup>
- can be transplanted into rodent brains<sup>4</sup>
- are permissive for HIV infection<sup>5</sup>
- can be transfected

3. Yankin, D.P., et al. (1993) *Proc. Natl. Acad. Sci. USA* 90: 2174-2178

4. Werkin, R., et al. (1993) *Proc. Natl. Acad. Sci. USA* 90: 9513-9517

5. Hlinka, G., et al. (1991) *J. Virol.* 65: 2732-2735

Call Stratagene at 1-800-424-5444  
for information about placing standing orders and volume discounts

hNT Neuron Kit: Catalog #204100  
NT2 Precursor Cells: Catalog #204101  
Circle No. 34 on Readers' Service Card

\*U.S. Patent No. 5,175,103





# In LSC, It's Experience That Counts



## Building on Forty Years of Innovations, the New Packard Tri-Carbs Are Ready For the 21st Century

In 1954, Packard introduced the world's first liquid scintillation counter (LSC): the Tri-Carb 314. In both name and model number, it heralded a new era of biological research: sensitive, accurate, and convenient measurement of  $^3\text{H}$  and  $^{14}\text{C}$ . For 40 years, Packard LSC innovations have consistently set the benchmarks against which other LSC's are compared.

Now, there's a whole new generation of Tri-Carbs: the 2100, 2300, and 2700 series LSC's. Different by design, today's Tri-Carbs are 386 or 486 computer controlled analyzers. These new Tri-Carbs reduce background by a factor of ten, cut cocktail use in half, slash waste disposal costs by up to 80%, and increase throughput.

These instruments even test their own performance automatically.

Today's Tri-Carbs are designed for the future with multi-user, single key counting, sample analysis, data processing with your choice of software, and final answer output. Direct networking lets you automatically send your results where you want them. And, with the available OS/2® operating system and Packard's Windows-based ChroniCal™ database software, you can easily customize your own complete assay quality control system.

There's even more, so call today and ask for information about the new Tri-Carb 2100, 2300, and 2700 LSC's.



Packard Instrument Company, 800 Research Parkway, Meriden, CT 06450 U.S.A.  
Tel: 203-238-2351 Toll Free: 1-800-323-1891 TX: 643251 FAX: 203-639-2172




**Packard International Offices:**

**Australia**, Mt Waverley 61-3-543-4266; **Austria**, Vienna 43-1-302504-0; **Belgium**, Brussels 32-2-4668210; **Canada**, Ontario 1-800-387-9559; **Denmark**, Greve 45-42909023; **France**, Rungis (33) 1 46.86.27.75; **Germany**, Frankfurt (49-69) 663010; **Italy**, Milano 39-2-33910796/7/8; **Japan**, Tokyo 81-3-3866-5850; **Netherlands**, Groningen 31-50-413360; **Tilburg** (013) 423900; **Russia**, Moscow, 7-095-238-7335; **Switzerland**, Zurich (01) 481 69 44; **United Kingdom**, Pangbourne, Berks (44) 0734 844981.

Circle No. 21 on Readers' Service Card



I  
like  
to draw.





# Structure searching on the STN network from CAS. Be creative.

We know that when it comes to pharmaceutical research, a structure represents a substance in the way you best understand.

We also know that when armed with our STN Express® software at your personal computer, you can easily search databases on STN®. And when you view structure graphics, they look as they would in a book. Or a journal.

That is why we made structures and structure searching key advantages of using STN for pharmaceutical research.

So why not give it a try? All you have to do is draw a chemical structure. Search it, and the substance will be instantly identified by a CAS Registry Number® from our extensive CAS REGISTRY<sup>SM</sup> database.

And with that number, you'll be able to find all kinds of related information about the identified substance from our many STN files. Everything from R&D to drug launch. You'll be able to locate references to prior research and numeric property data. Even molecular connection tables for modeling. Plus, STN includes an enormous amount of online biosequence data and patent information.


STN is also quite possibly the largest storehouse of reaction information in the world. And it's just as easy to find as substance information.

Of course, you'll also be able to search by text in databases like Chemical Abstracts, International Pharmaceutical Abstracts, the Merck Index Online, EMBASE, and MEDLINE. Files from PJB Publications...including Pharmaprojects. IMSWorld R&D Focus files. And PROMT and Pharmaceutical News Index.

For more about structure searching and pharmaceutical information on the STN network from CAS, call **1-800-933-4350** or **614-447-3731**. Or fax **614-447-3751**.

In this world, there is some information you absolutely have to know.

We have it. And we'd like to share it with you.



A Division Of The American Chemical Society®


Circle No. 17 on Readers' Service Card





1360

Motley crowd  
of planets



1432

Connecting California  
earthquakes

## NEWS & COMMENT

- Forensic DNA Goes to Court With O.J. 1352  
German Lab Points the Way in Hunt for 'Hot' Plutonium 1355  
Scopes and Squirrels Return to Court 1356  
Bromley Memoirs Take Off the Gloves 1357  
Buckyballs and Kin Are Still on a Roll 1357

## RESEARCH NEWS


- The Solar System's New Diversity 1360  
Micro-Quasars Found in Our Galaxy 1362  
Making Molecular Filters More Reactive 1363  
Link to Hereditary Melanoma Brightens Mood for p16 Gene 1364  
Alzheimer's: Could There Be a Zinc Link? 1365

## RESEARCH ARTICLE

- Structures of Active Conformations of  $G_{101}$  and the Mechanism of GTP Hydrolysis 1405  
D. E. Coleman, A. M. Berghuis, E. Lee, M. E. Linder, A. G. Gilman, S. R. Sprang

## REPORTS

- Molecular Recognition by Self-Assembled Monolayers of Cavitand Receptors 1413  
K. D. Schierbaum, T. Weiss, E. U. Thoden van Velzen, J. F. J. Engbersen, D. N. Reinhoudt, W. Göpel



## FRONTIERS IN MEDICINE: VACCINES

### NEWS REPORTS

- Bumps on the Vaccine Road 1371  
Are Researchers Racing Toward Success, Or Crawling? 1373  
U.S. National Program Is Going Nowhere Fast 1375  
Childrens' Vaccine Initiative Stumbles 1376

### POLICY FORUM

- The United States Needs a National Vaccine Authority 1378  
B. R. Bloom

### PERSPECTIVES

- Malaria Vaccines: Multiple Targets 1381  
R. S. Nussenzweig and C. A. Long  
Vaccines for Varicella-Zoster Virus and Cytomegalovirus: Recent Progress 1383  
S. A. Plotkin

- Nanocatalysis by the Tip of a Scanning Tunneling Microscope Operating Inside a Reactor Cell 1415  
B. J. McIntyre, M. Salmeron, G. A. Somorjai

- A Redox Fuel Cell That Operates with Methane as Fuel at 120°C 1418  
S. H. Bergens, C. B. Gorman, G. T. R. Palmore, G. M. Whitesides

## DEPARTMENTS

- THIS WEEK IN SCIENCE 1341  
EDITORIAL 1343  
LETTERS 1345  
Embryo Research Guidelines: H. Varmus • Primates and New Viruses: J. S. Allan • The Sobering D<sub>2</sub> Story: K. Blum and E. P. Noble; C. Holden • Machine-Written Summary: Missed Opportunity?: N. R. Olsen • Correction: Ferricyanide and Pattern Formation: K.-J. Lee, W. D. McCormick, Q. Quyang, H. L. Swinney

- SCIENCESCOPE 1351  
RANDOM SAMPLES 1358  
BOOK REVIEWS 1468  
Richard Owen, reviewed by C. Limoges • *Evolution of Infectious Diseases*, J. J. Bull and B. R. Levin • *Remote Geochemical Analysis*, B. M. Jakosky • Vignettes • Books Received  
PRODUCTS & MATERIALS 1472

## Board of Reviewing Editors

Frederick W. Alt  
Don L. Anderson  
Michael Ashburner  
Stephen J. Benkovic  
David E. Bloom  
Floyd E. Bloom  
Piet Borst  
Henry R. Bourne  
Michael S. Brown  
James J. Bull

Kathryn Calame  
C. Thomas Caskey  
Dennis W. Choi  
John M. Coffin  
Paul J. Crutzen  
James E. Dahlberg  
Robert Desimone  
Bruce F. Eldridge  
Paul T. Englund  
Richard G. Fairbanks

Douglas T. Fearon  
Harry A. Fozzard  
Klaus Friedrich  
Theodore H. Geballe  
John C. Gerhart  
Roger I. M. Glass  
Stephen P. Goff  
Peter N. Goodfellow  
Corey S. Goodman  
Ira Herskowitz

Eric F. Johnson  
Stephen M. Kosslyn  
Michael LaBarbera  
Nicole Le Douarin  
Charles S. Levinger III  
Alexander Levitzki  
Harvey F. Lodish  
Richard Losick  
Diane Mathis  
Anthony R. Means

Shigetada Nakanishi  
Roger A. Nicoll  
Stuart L. Pimm  
Yeshayau Pocker  
Dennis A. Powers  
Ralph S. Quatrano  
V. Ramanathan  
Douglas C. Rees  
T. M. Rice  
David C. Rubie

Erkki Ruoslahti  
Gottfried Schatz  
Jozef Schell  
Ronald H. Schwartz  
Terrence J. Sejnowski  
Ellen Solomon  
Thomas A. Steitz  
Michael P. Stryker  
Robert T. N. Tjian  
Emil R. Unanue


Geerat J. Vermeij  
Bert Vogelstein  
Harold Weintraub  
Arthur Weiss  
Zena Werb  
George M. Whitesides  
Owen N. Witte  
William A. Wulf



## COVER

This *Frontiers in Medicine* issue addresses the science and politics of vaccines. The hypodermic syringe, a classic tool of vaccination, may soon be out-of-date. A global effort is under way to produce a "supervaccine" that would confer, in a single oral dose, lifetime immunity against

many of the major infectious diseases. A special section begins on page 1371 and related Reports on pages 1448 and 1451. [Images: needle, Peter Steiner/The Stock Market; globe, Tom Van Sant/Geosphere Project/The Stock Market. Illustration: E. Carroll]



1405

A  $G\alpha$  protein transition state analog complex:  $GDP\cdot AIF_4^-$

**Pneumococcal Disease: Prospects for a New Generation of Vaccines** 1385  
G. R. Siber

**Cholera Vaccines: Fighting an Ancient Scourge** 1387  
J. J. Mekalanos and J. C. Sadoff

**Rotavirus Vaccines: Success by Reassortment?** 1389  
R. I. Glass, J. Gentsch, J. C. Smith

**Measles Vaccine: Do We Need New Vaccines or New Programs?** 1391  
S. L. Katz and B. G. Gellin

**Prospects for a Respiratory Syncytial Virus Vaccine** 1393  
C. B. Hall

## ARTICLES

**Lymphocyte Life-Span and Memory** 1395  
J. Sprent and D. F. Tough

**Vaccine Technologies: View to the Future** 1401  
N. Regina Rabinovich, P. McInnes, D. L. Klein, B. F. Hall

**A Spectroscopic Measurement of the Coronal Density of Procyon** 1420  
J. H. M. M. Schmitt, B. M. Haisch, J. J. Drake

**Morphology and Ionization of the Interstellar Cloud Surrounding the Solar System** 1423  
P. C. Frisch

**Chemical Sequence Control of  $\beta$ -Sheet Assembly in Macromolecular Crystals of Periodic Polypeptides** 1427  
M. T. Krejchi, E. D. T. Atkins, A. J. Waddon, M. J. Fournier, T. L. Mason, D. A. Tirrell

**Stress Triggering of the 1994  $M = 6.7$  Northridge, California, Earthquake by Its Predecessors** 1432  
R. S. Stein, G. C. P. King, J. Lin

**Seismic Evidence for a Lower-Crustal Detachment Beneath San Francisco Bay, California** 1436  
T. M. Brocher, J. McCarthy, P. E. Hart, W. S. Holbrook, K. P. Furlong, T. V. McEvilly, J. A. Hole, S. L. Klemperer

**Redistribution of Intracellular  $Ca^{2+}$  Stores During Phagocytosis in Human Neutrophils** 1439  
O. Stendahl, K.-H. Krause, J. Krischer, P. Jerström, J.-M. Theler, R. A. Clark, J.-L. Carpentier, D. P. Lew

**Ku80: Product of the XRCC5 Gene and Its Role in DNA Repair and V(D)J Recombination** 1442  
G. E. Taccioli, T. M. Gottlieb, T. Blunt, A. Priestley, J. Demengeot, R. Mizuta, A. R. Lehmann, F. W. Alt, S. P. Jackson, P. A. Jeggo

**Thrombocytopenia in *c-mpl*-Deficient Mice** 1445  
A. L. Gurney, K. Carver-Moore, F. J. de Sauvage, M. W. Moore

**Engineering Poliovirus as a Vaccine Vector for the Expression of Diverse Antigens** 1448  
R. Andino, D. Silvera, S. D. Suggett, P. L. Achacoso, C. J. Miller, D. Baltimore, M. B. Feinberg

**Prophylactic Vaccines, Risk Behavior Change, and the Probability of Eradicating HIV in San Francisco** 1451  
S. M. Blower and A. R. McLean

**Degradation of  $G\alpha$  by the N-End Rule Pathway** 1454  
K. Madura and A. Varshavsky

**Ras-Dependent Growth Factor Regulation of MEK Kinase in PC12 Cells** 1458  
C. A. Lange-Carter and G. L. Johnson

**Genetics of a Pheromonal Difference Contributing to Reproductive Isolation in *Drosophila*** 1461  
J. A. Coyne, A. P. Crittenden, K. Mah

**Rapid Induction of Alzheimer  $A\beta$  Amyloid Formation by Zinc** 1464  
A. I. Bush, W. H. Pettingell, G. Multhaup, M. d. Paradis, J.-P. Vonsattel, J. F. Gusella, K. Beyreuther, C. L. Masters, R. E. Tanzi

## AAAS Board of Directors

Elaine E. Clark  
Retiring President,  
Chairman  
Francisco J. Ayala  
President  
Rita R. Colwell  
President-elect

William A. Lester Jr.  
Simon A. Levin  
Anna C. Roosevelt

Alan Schriesheim  
Jean'ne M. Shreeve  
Chang-Lin Tien  
Warren M. Washington  
Nancy S. Wexler

William T. Golden  
Treasurer  
Richard S. Nicholson  
Executive Officer

■ **SCIENCE** (ISSN 0036-8075) is published weekly on Friday, except the last week in December, by the American Association for the Advancement of Science, 1333 H Street, NW, Washington, DC 20005. Second-class postage (publication No. 484460) paid at Washington, DC, and additional mailing offices. Copyright © 1994 by the American Association for the Advancement of Science. The title **SCIENCE** is a registered trademark of the AAAS. Domestic individual membership and subscription (51 issues): \$92 (\$50 allocated to subscription). Domestic institutional subscription (51 issues): \$215. Foreign postage extra: Mexico, Caribbean (surface mail) \$50; other countries (air assist delivery) \$95. First class, airmail, student and emeritus rates on request. Canadian rates with GST available upon request, GST #1254 88122. Printed in the U.S.A.

## ■ Indicates accompanying feature

**Change of address:** allow 6 weeks, giving old and new addresses and 11-digit account number. **Postmaster:** Send change of address to *Science*, P.O. Box 2033, Marion, OH 43305-2033. **Single copy sales:** \$6.00 per issue prepaid includes surface postage; Guide to Biotechnology Products and Instruments, \$20. Bulk rates on request. **Authorization to photocopy** material for internal or personal use under circumstances not falling within the fair use provisions of the Copyright Act is granted by AAAS to libraries and other users registered with the Copyright Clearance Center (CCC) Transactional Reporting Service, provided that the base fee of \$1 per copy plus \$0.10 per page is paid directly to CCC, 27 Congress Street, Salem, MA 01970. The identification code for *Science* is 0036-8075/93 \$1 + .10. *Science* is indexed in the *Reader's Guide to Periodical Literature* and in several specialized indexes.



# ANTIPEPTIDE ANTIBODIES

Can you be  
more specific?

Definitely.

Improved antisera can advance your research. Immunizing with peptides gives you higher titers and more specific antibodies. Now you can see even better results—with multiple antigenic peptides, or MAPs.

Antisera produced with MAPs contain specific antibodies with a monoclonal nature. Yet they're much easier to get. You can easily make exactly the immunogen you need with the Synergy™ Peptide Synthesizer and our special MAP resin.



*Create exactly the immunogens you need—overnight—with Synergy and our MAP resin.*

*Galactosyl receptors in sperm labeled with antibodies prepared using Synergy and MAPs.*

With 4 or 8 peptide arms branching from a non-immunogenic core, MAPs maximize the mass and weight percent of antigen, dramatically improving immunogenic response.

And MAPs eliminate the need for a carrier protein, so you save time and effort.

When we designed the Synergy system, we optimized automated peptide synthesis to match the needs—and budgets—of bioresearchers. After a 10-minute set up, an overnight run, and simple post-synthesis work-up, you're ready to immunize.

The Synergy system and the MAP resin are backed by the leaders in peptide synthesis. With our integrated resources, we provide the most comprehensive and advanced range of bioresearch systems, technologies and support in the world.

Find out how researchers like you are accelerating their research with Synergy. **In the U.S., phone Perkin-Elmer's Applied Biosystems Division at 1-800-345-5224.** Outside the U.S., contact your local Perkin-Elmer representative. And discover how easy it is to be more specific.

**PERKIN ELMER**

**Europe** Weiterstadt, Germany Tel: 49-6150-101-0 Fax: 49-6150-101-101

**Canada** Mississauga, Canada Tel: 800-668-6913 Fax: 905-821-8246

**Japan** Tokyo, Japan Tel: 81-473-80-8500 Fax: 81-473-80-8505

**Latin America** Mexico City, Mexico Tel: 52-5-593-7077 Fax: 52-5-593-6223

**Australia** Melbourne, Australia Tel: 61-3-212-8585 Fax: 61-3-212-8502

Photo courtesy of Dr. A. Kierszenbaum, City University of New York Medical School.

Synergy is a trademark of Perkin-Elmer's Applied Biosystems Division.

Perkin-Elmer is a registered trademark of The Perkin-Elmer Corporation.

Circle No. 22 on Readers' Service Card



## In their cups

Many organic molecules that adsorb onto a clean metal surface are held in place by relatively weak intermolecular forces. Schierbaum *et al.* (p. 1413) have used self-assembly to adsorb organic molecules that, through molecular recognition, can selectively adsorb a small nonpolar organic molecule,  $C_2Cl_4$ . The cup-like resorcin[4]arene monolayer adsorbed on a gold surface increased the desorption temperature for adsorbed  $C_2Cl_4$  molecules. Such self-assembled layers could have applications as highly selective sensors.

## Catalyst tips

Numerous scanning tunneling microscope (STM) studies have been reported of surfaces that can act as catalysts. McIntyre *et al.* (p. 1415) have used a platinum-rhodium STM tip itself as a catalyst to hydrogenate carbonaceous layers on a platinum surface. When short voltage pulses were applied to these tips operating in predominantly or pure hydrogen environments, hydrocarbon clusters near the tip were rehydrogenated.

## Fuel for thought

Practical fuel cells are powered by hydrogen, which is not a primary energy source. For example, alkanes such as methane, which is abundant in natural gas, are steam-reformed to generate hydrogen; direct oxidation of alkanes in fuel cells is slow. Bergens *et al.* (p. 1418) have taken advantage of reaction in which platinum black adsorbed on graphite oxidizes methane to carbon dioxide while converting iron (III) to iron (II) in concentrated sulfuric acid solutions. They coupled

## G protein structures and GTP hydrolysis


A crucial step in signal transduction by G proteins is the hydrolysis of guanosine triphosphate (GTP) to guanosine diphosphate (GDP). However, hydrolysis of GTP by  $G_{\alpha}$  subunits such as  $G_{i\alpha 1}$  is actually quite slow, usually only a few hundred times faster than the background rate, and unraveling the mechanism has proven difficult. Coleman *et al.* (p. 1405) present crystal structures of the  $G_{i\alpha 1}$  subunit complexed with  $GTP\gamma S$  and with  $GDP\cdot AlF_4^-$ . These studies show that two highly conserved residues, Gln<sup>204</sup> and Arg<sup>178</sup>, stabilize the transition state—the former helps orient the nucleophilic water molecule, while the latter helps stabilize negative charge that builds up on the  $\gamma$  phosphate oxygen atoms.

this reaction with a vanadium reaction that reduces oxygen to create a working fuel cell based on methane. This system operates at relatively low temperature (120°C) compared to hydrogen cells (700° to 900°C).

rium, a state which at the low densities involved must take million of years to attain.

## Our place in space

The solar system resides in a cloud of interstellar gas that appears to be the quiescent remnant of a supernova explosion of uncertain age. The ionized gas in the cloud imposes absorption lines on spectra of nearby stars, and Frisch (p. 1423) used the wavelength and strength of



lines in several stellar spectra to reconstruct the geometry of the cloud, the sun's position within it, and its ionization state. Our system appears to be skating along just inside the cloud, suggesting that we entered it no more than a few thousand years. Comparison of lines due to different elements indicates that the gas is in ionization equilib-

## Underlying causes

Two reports examine possible temporal and spatial connections between earthquakes and faults in California. Stein *et al.* (p. 1432) modeled how stress was transferred due to large earthquakes in Los Angeles. Their results suggest that quakes in 1933 and 1952 slightly increased stresses on the fault that ruptured in the 1971 San Fernando quake, and that this quake increased stresses on the fault that ruptured in the 1994 Northridge quake. Brocher *et al.* (p. 1436) examined the deep structure of the San Andreas, Hayward, and other faults near San Francisco. Seismic data revealed a prominent subhorizontal reflector at a depth of about 15 kilometers that may be a detachment surface that truncates the shallow faults.

## Role in repair

The Ku protein, a heterodimer that was originally identified as an autoantigen in sera, is known to bind to the free ends of double-stranded DNA. Taccioli *et al.* (p. 1442) show that its 80-

kilodalton subunit, Ku80, is the product of the human DNA repair gene XRCC5. They found that introducing the gene that encodes Ku80 complemented the *xrs-6* mutation in Chinese hamster ovary cells that leads to defects in repair of DNA double-strand breaks and V(D)J recombination. Study of the role of Ku protein in these processes may also shed light on the action of DNA-dependent protein kinase, in which Ku serves as a DNA-binding component.

## In the blood

Recently the cytokine that regulates platelet production, thrombopoietin (TPO), was purified and cloned. Gurney *et al.* (p. 1445) have investigated the specificity of TPO for its receptor, c-Mpl, in mice that have had this gene deleted. These mice were normal except for a marked decrease in platelet numbers and high concentrations of circulating TPO. Such results are important in considering the possible side effects of administering TPO to immunodeficient patients.

## Who's who

A major unsolved mystery in evolution is the genetic basis of speciation. One significant cause of speciation may be mating discrimination. Coyne *et al.* (p. 1461) show that one factor in sexual isolation for several *Drosophila* species is differences in the hydrocarbons of the female cuticle. When females of different fly species are crowded together, the cuticular hydrocarbons can be transferred between females of different species and can subsequently alter male sexual behavior toward females of another species.




# OLIGO 5.0 PCR SOFTWARE. WITH 22 NEW FEATURES, WHO KNOWS WHAT YOU MIGHT COME UP WITH?

It was huge. Some called it colossal.

With the release in 1989 of Oligo Primer Analysis Software, researchers using PCR and DNA sequencing were getting a powerful new tool. Soon it became the standard by which all others were measured.

Now, NBI has ushered in a worthy successor: Oligo 5.0. And all of a sudden, everything else seems prehistoric.




So powerful, yet so easy. Just click on any oligo for instant analysis.

## BETTER PCR PRIMERS, BY DESIGN.

If you're big on PCR, here's something you can sink your teeth into. Because with Oligo 5.0, you have the ultimate program for selecting *all* optimal PCR primer pairs from any DNA/RNA sequence file. In seconds.

Its easy-to-use search function selects primers which are not only unique, highly specific, and cross-compatible, they're also free of dimer, hairpin, and sequence repeats. Oligo 5.0 even searches for sequencing primers and hybridization probes. So you can design the best possible primers by simply lifting a finger.



All optimal PCR primer pairs are listed with complete analyses.

What's more, it calculates optimal annealing temperature, PCR product composition, various  $T_m$  calculations, and extinction coefficients.

But perhaps most impressive is what we've added: Like LCR analysis, a GC clamp, primer  $T_m$  matching, highly accurate false priming checks, and nested primer selection.

## THE EVOLUTION OF OLIGO: HIGHER INTELLIGENCE.

Oligo 5.0 – for Windows or Mac – works so efficiently, so

intelligently, and so intuitively, it could render all lesser programs extinct. Here are some other notable advancements:

With a new oligo database and order feature, you can download primer selections and user data to an order form. Just print and send it to your synthesis provider.

In addition, hybridization time can now be calculated for any DNA segment. You also gain back-translation capabilities, which no other program offers. Plus a host of other eye-popping enhancements.

## SEEING IS BELIEVING.

To fully appreciate the magnitude of Oligo 5.0, ask us for a free demo disk

(PC or Mac). In North America and elsewhere, call (800) 747-4362 or (612) 550-2012. Or fax (612) 550-9625. In Europe, you can call us at 47 22 20 01 37 or fax 47 22 20 01 89. In Japan, call 0775-43-7231 or fax 0775-43-9254. You may also contact us by E-mail: [nbi@biotechnet.com](mailto:nbi@biotechnet.com).

With software this powerful, imagine the possibilities. Some could even be gigantic.



**Oligo**

Primer Analysis Software from  
National Biosciences, Inc.



# STATISTICA/w

Windows,  
DOS, Macintosh

**STATISTICA/w™** (for Windows) Complete Statistical System with thousands of on-screen customizable, presentation-quality graphs fully integrated with all procedures ■ Complete Windows 3.1 support, DDE, OLE, TT-fonts, multiple toolbars, right mouse button support ■ Unlimited numbers of data-, results-, and graph-windows ■ Inter-window integration: data, results, and graphs can be treated as objects and converted into one another in a number of ways ■ The largest selection of statistics and graphs in a single system; comprehensive implementations of: Exploratory techniques; multi-way tables with banners (presentation-quality reports); nonparametrics; distribution fitting; multiple regression; general nonlinear estimation; stepwise logit/probit; general ANCOVA/MANCOVA; stepwise discriminant analysis; log-linear analysis; factor analysis; cluster analysis; multidimensional scaling; canonical correlation; item analysis/reliability; survival analysis; time series modeling; forecasting; lags analysis; quality control; process analysis; experimental design (with Taguchi); and much more ■ Manuals with comprehensive introductions to each procedure and examples ■ Hypertext-based Stats Advisor expert system ■ Extensive data management facilities (spreadsheet with long formulas, block operations, advanced Clipboard support, DDE hot links, relational merge, data verification, powerful programming language) ■ Batch command language and macros also supported, "turn-key system" options ■ All output displayed in Scrollsheets™ (dynamic, customizable, presentation-quality tables with toolbars, pop-up windows, and instant 2D, 3D and multiple graphs) ■ Extremely large analysis designs (e.g., correlation matrices up to 32,000x32,000, unlimited ANOVA designs) ■ Megafiler Manager with up to 32,000 variables (8 Mb) per record ■ Unlimited size of files; extended ("quadruple") precision; unmatched speed ■ Exchanges data and graphs with other applications via DDE or an extensive selection of file import/export facilities ■ Hundreds of types of graphs, including categorized multiple 2D and 3D graphs, matrix plots, icons, and unique multivariate (e.g., 4D) graphs ■ Facilities to custom design new graphs and add them permanently to menu ■ On-screen graph customization with advanced drawing tools, interactive stretching and resizing of complex objects, interactive embedding of graphs and artwork, special effects, icons, maps, multi-graphics management, page layout control for slides and printouts; unmatched speed of graph redraw ■ Interactive rotation, perspective and cross-sections of all 3D and 4D graphs ■ Extensive selection of tools for graphical exploration of data: fitting, smoothing, overlaying, spectral planes, projections, layered compressions, marked subsets ■ Price **\$995**.

**Quick STATISTICA/w™** (for Windows) A comprehensive selection of basic statistics and the full graphics capabilities of STATISTICA/w ■ Price **\$495**.

**STATISTICA/dos™** (for DOS) A STATISTICA/w-compatible data analysis system ■ Price **\$795**.

**Quick STATISTICA/dos™** (for DOS) A subset of STATISTICA/dos statistics and graphics ■ Price **\$295**.

Domestic sh/h \$10 per product; 14-day money back guarantee.

Circle No. 5 on Readers' Service Card




**StatSoft™**

2325 E. 13th St. • Tulsa, OK 74104 • (918) 583-4149  
Fax: (918) 583-4376

Overseas Offices: StatSoft of Europe (Hamburg, FRG), ph: 040/4200347, fax: 040/4911310; StatSoft UK (London, UK), ph: 0462/482822, fax: 0462/482855; StatSoft Pacific (Melbourne, Australia), ph: (03) 663 6580, fax: (03) 663 6117; StatSoft France (Paris), ph: (1) 45 66 97 00, fax: (1) 45 66 06 51; Available from other Authorized Representatives worldwide: Sweden: AkademiData Scientific AB ph: 018-210035, fax: 018-210039; Finland: Statcon Oy ph: 24-334678, fax: 24-333867; Belgium: Texma Newtech 10 61 16 28; South Africa: Osiris Technical Systems 12-663-4500; Japan (Macintosh): Three's Company, Inc., 03-3770-7600; Japan (Windows): Design Technologies, Inc., 03-3667-1110.

StatSoft, STATISTICA/Mac, Quick STATISTICA/Mac, STATISTICA/w, Quick STATISTICA/w, and Scrollsheet are trademarks of StatSoft, Inc.; Macintosh, Mac IIx, Excel and MacDraw are trademarks of their respective companies.



**STATISTICA/Mac™** (for Macintosh) A STATISTICA/w-compatible, comprehensive data analysis and graphics system designed for the Macintosh ■ Large selection of statistical methods fully integrated with presentation-quality graphics (incl. EDA, multiplots, a wide selection of interactively rotatable 3D graphs; MacDraw-style tools) ■ Unlimited size of files ■ Full support for System 7, incl. "Publish and Subscribe" ■ Price **\$695**.

**Quick STATISTICA/Mac™** (for Macintosh) A subset of STATISTICA/Mac: a comprehensive selection of basic statistics and the full graphics capabilities of STATISTICA/Mac ■ Price **\$395**.





**FOURTH INTERNATIONAL  
SYMPOSIUM ON IMMUNOTOXINS  
KINGSTON RADISSON RESORT  
MYRTLE BEACH, SOUTH CAROLINA**

**JUNE 8-11, 1995**

**SESSIONS/WORKSHOPS  
and CHAIRPEOPLE**

1. **Toxin Structure/Function:** S. Olsnes, J. Robertus, D. FitzGerald, J. Collier
2. **Single Chain Fv's Structure/Function:** S. Carroll, B. Bird
3. **Cytokine Structure/Function:** T. Ciardelli
4. **Toxin Modification:** W. Blatter, E. Vitetta, I. Pastan, J. Murphy, F. Stirpe
5. **Novel Ligands:** P. Thorpe, F. Preijers, R. Wiley
6. **Novel Toxins:** S. Ramakrishnan, S. Lepia, A. Tonevitsky
7. **Animal Studies:** D. Vallera, D. Bigner, C. Siegall, R. King
8. **Vascular Leak Syndrome:** M. Willingham, J. Schnitzer
9. **Tumor Penetration:** C. Sung, M. Rosenblum
10. **Immunogenicity:** M. Tepper
11. **Cell Resistance:** R. Bast, R. Kreitman
12. **Clinical Studies/Leukemia and Lymphoma:** E. Sausville, S. Esseltine, F. Foss, B. Fallini, A. Engert, F. Uckun, J. Kersey, D. Flavell
13. **Clinical Studies/Solid Tumors:** R. Youle, L. Pai, J. Nichols, T. Lynch
14. **Clinical Studies/Autoimmune Disease:** T. Woodworth, J. Nichols, P. Scannon
15. **Biological and Toxin Weapons in a Changing World:** R. Manchee, T. Hudson

For Information Contact:

**Arthur E. Frankel, MD**  
Hollings Cancer Center, Room 306  
Medical University of South Carolina  
171 Ashley Avenue  
Charleston, South Carolina 29425  
(803) 792-1450

# PFANSTIEHL

## The Start Of Something Big: Ribonolactone

D-Ribono-1, 4-lactone and its 2,3-O-isopropylidene derivative have been the starting points for the synthesis of a wide variety of natural products, C-glycosides and carbocycles. We now offer both of these important chiral starting materials, in addition to our wide selection of other sugar lactones. All are available in quantities to meet your needs; research through scale-up. Write or call for a technical brochure and price list and see how Ribonolactone can be the start of something big for you.

### PFANSTIEHL LABORATORIES, INC.

*The source for carbohydrate chemistry*

1219 Glen Rock Avenue/Waukegan, IL 60085-0439  
Tel.: 1-708/623-0370/Toll Free: 1-800/383-0126  
FAX: 708/623-9173  
77-W



# MAKING THINGS WORK

Circle No. 15 on Readers' Service Card

## With B-I-T-S®, the Choice is Yours.

### A Customized Literature Search Service From BIOSIS.

**W**ith an increasing volume of journals, authors and topics to track, you need to use your time and money wisely. Be sure you choose **B-I-T-S** (*BIOSIS Information Transfer System*), the only life science information service focused to match your needs.

#### **B-I-T-S: The Smart Choice.**

Most information services force you to pay for references from standard profiles, but not **B-I-T-S**. To ensure that **B-I-T-S** references reflect **your** research focus, we create a personalized search strategy from a specially designed Profile Development Form. The result – targeted information that makes **B-I-T-S** your most cost-effective choice.

**B-I-T-S** supplies references from over 7,000 international life science publications. Comprehensive coverage, combined with a customized search strategy, is our formula for providing you with the information you require.

#### **Put B-I-T-S to Work with Our Special Offer.**

Call today for our \$150 three-month introductory offer or a no-obligation estimate for an annual subscription:


**1-800-523-4806** (USA and Canada)  
**(215) 587-4847** (Worldwide)



**BIOSIS®**

*Information for Today's Decisions and Discoveries*

BIOSIS is a registered trademark of Biological Abstracts, Inc. S994BT




Circle No. 11 on Readers' Service Card



# Sigma.

Supporting  
the Future  
of Research  
With a  
Tradition  
of Service.



What distinguishes a chemical company: First and foremost - service; which is Sigma's approach to supporting the future of research. Should technical questions arise, a Sigma specialist is available to assist you with the use or application of our 33,000-plus products. In addition - quality; every one of Sigma's products is thoroughly assayed in state-of-the-art analytical laboratories using the most demanding analytical criteria.



Another distinguishing factor is expertise; as in Sigma's family of experts in life science research. Sigma's chemists are constantly at work on a wide spectrum of specialty chemicals for targeted areas of research, adding over 2,000 new products each year. And they're continually refining and improving the products that established Sigma's reputation over 50 years ago.

Service, quality and expertise distinguish Sigma as the leading choice of laboratories worldwide. Weigh all the advantages and consider Sigma as your single source for research biochemicals and organic compounds.

**SIGMA** 

Where Science and Service Come Together.

Call collect: 314-771-5750,  
800-325-3010,  
or contact your local Sigma office.

AUSTRIA AUSTRALIA BELGIUM BRAZIL CZECH REPUBLIC FRANCE GERMANY HUNGARY INDIA ITALY  
MEXICO NETHERLANDS JAPAN SOUTH KOREA SPAIN SWITZERLAND UNITED KINGDOM UNITED STATES

A Sigma-Aldrich Company

Circle No. 10 on Readers' Service Card





Designed to meet the growing demands for integrated systems in life science research, <sup>TM</sup>Vistra DNA Systems brings together the strengths of Molecular Dynamics' instrumentation capability and Amersham's expertise in innovative reagent kit chemistry.



# Sequencing Automation from Cells to Gels

The <sup>TM</sup>Vistra DNA Labstation 625 is the only system that provides 'walkaway' automation in DNA sequencing. Template preparation and sequencing reactions are completely integrated on the Labstation, with flexible protocols available for all major sequencing strategies.

## NO CENTRIFUGATION WITH FMP

Central to the design of the Labstation is FMP<sup>TM</sup>, a novel separation method based on magnetic purification that removes the need for centrifugation. A complete range of reagent kits ensure optimum performance with either M13 or plasmid DNA templates.



## INTEGRATION WITH FLUORESCENT SEQUENCERS


Designed for use with four color single lane fluorescent sequencers, the Labstation offers unparalleled automation, consistency and accuracy in DNA sequencing.

For complete automation whatever your sequencing strategy

Call now: **1-800-9-VISTRA** (in the U.S.)

In France, call (33) 1-6086-6513. In Germany (49) 2151-83870. In the U.K. and other European countries (44) 732-762565. Molecular Dynamics Inc. • 928 East Arques Avenue • Sunnyvale, California 94086-4520 USA


VISTRA IS AN ALLIANCE OF



© 1994, Molecular Dynamics. All rights reserved. Product specifications subject to change without notice. Amersham, FMP and Vistra DNA are trade marks of Amersham International plc. Molecular Dynamics is a trademark of Molecular Dynamics Corp.

Circle No. 19 on Readers' Service Card

Ad 9440





# supernova



At NOVABIOCHEM we are especially proud to celebrate our 10th Anniversary with the NEW 1994/95 Catalog and Peptide Synthesis Handbook, featuring...

- Over 120 innovative NEW products
- Detailed technical specifications
- Highly competitive prices
- Clear, easy to use layout
- Unique 'Synthesis Notes' a complete review of peptide methodology

Our aim has been to produce a catalog that will become an invaluable resource for every lab making or using peptides.

To receive your NOVABIOCHEM Catalog and Peptide Synthesis Handbook please call your local distributor below.

**Australia**  
Calbiochem-  
Novabiochem Pty Ltd.  
Tel: (02) 318 0322  
Fax: (02) 319 2440

**Germany**  
Calbiochem-  
Novabiochem GmbH  
Tel: (06196) 63 955  
Fax: (06196) 62 361

**Japan**  
Calbiochem-  
Novabiochem Japan Ltd.  
Tel: (03) 5443 0281  
Fax: (03) 5443 0271

**Switzerland**  
Calbiochem-  
Novabiochem AG  
Tel: (062) 69 19 22  
Fax: (062) 69 19 92

**United Kingdom**  
Calbiochem-  
Novabiochem (UK) Ltd.  
Tel: (0800) 622935  
Fax: (0602) 430951

**USA & Canada**  
Calbiochem-  
Novabiochem Corp.  
Tel: (800) 228-9622  
Fax: (800) 432-9622

**Available locally worldwide**  
 CALBIOCHEM  
NOVABIOCHEM

Austria (01) 889 18 19 Belgium (010) 41 24 55 Denmark 44 94 88 22 Finland (0) 70 80 900 France (1) 46.26.78.70 Greece +41-62-69 19 22 India (011) 559 6820  
Israel (09) 90 45 14 Italy (02) 5521 3005 Japan (03) 3279 6363 Netherlands (076) 79 57 95 Norway (069) 261 777 Portugal (01) 458 2209 Spain (91) 551 5403 Sweden (08) 742 02 50

Circle No. 18 on Readers' Service Card



## FLAMING / BROWN MICROPIPETTE PULLER



**IF YOUR RESEARCH DEPENDS ON PIPETTES  
... DEPEND ON THE P-87 PIPETTE PULLER**

**MICROPIPETTES                      INJECTION NEEDLES  
   PATCH PIPETTES**

Microprocessor controlled and fully programmable, the P-87 is an easy to use, highly reproducible puller.

With storage for multiple pipette programs, the P-87 is ideal as a general purpose lab puller, or for the investigator who utilizes a variety of techniques.

Our technical support and quality control standards are the highest in the industry.



**Sutter Instrument Company**

40 Leveroni Court, Novato, CA 94949 Ph: (415) 883-0128 Fax: (415) 883-0572


Circle No. 14 on Readers' Service Card

# ISCAN®

## EYE MOVEMENT MONITORING SYSTEMS

ISCAN eye tracking products set the pace for performance and affordability in the research marketplace. Applications include fixation monitoring, vision research, advertising research, psycho-visual experiments, reading studies, toxicology, human factors evaluation and drug or alcohol response measurement.

Non-invasive, real-time video image processors automatically measure pupil position & size



and compute subject's point of gaze. Benchtop, head mounted & high speed versions as well as computer data analysis software are available. Complete Eye Movement Monitoring Laboratory or component systems may be provided to suit any research application.

**ISCAN, Inc.**

125 Cambridgepark Drive  
Cambridge, MA 02140, U.S.A.

Tel: (617) 868-5353 Fax: (617) 868-9231


Circle No. 7 on Readers' Service Card

## Another first from SCIENCE... The 1995 SCIENCE Calendar

From the first experiments with electricity to the genetic revolution, scientists around the globe have looked to the pages of SCIENCE for leading-edge research and scientific news. Now SCIENCE announces another publishing event: **The SCIENCE 1995 Calendar**. This full-color calendar features:

- 12 full-color mini-poster SCIENCE covers
- Oversized 9"x 12" format
- Large date blocks with major holidays
- An important moment in the history of science highlighted each month

Whether you're planning scientific research or researching the perfect gift for your favorite scientist, you'll want to have your own copy of the **1995 SCIENCE Calendar**!



### It's Easy to Order!

I want to show my support for SCIENCE. Please send me \_\_\_\_\_ copies of the 1995 SCIENCE Calendar, for \$12.50 (AAAS Member price \$11.50). I have enclosed a check for the amount of the calendar (plus shipping charges for non-US orders). I understand my order is refundable if I am not completely satisfied.

### Shipping Address (please print)

Name \_\_\_\_\_

(last name)

(first name)

Address \_\_\_\_\_

City/State/Postal Code \_\_\_\_\_

Please allow 4-6 weeks for delivery. Prices and shipping charges subject to change without notice. Calendars returned must be in salable condition.

Wash., DC residents please add sales tax.

**SCIENCE**

### Payment Information

☐ US dollar check    ☐ VISA    ☐ Master Card

Important: All payments must be made in US dollars. \$25 minimum order for all credit card orders. Make check payable to AAAS Science Publications, Inc.

Credit card number \_\_\_\_\_

Expiration date \_\_\_\_\_

Signature \_\_\_\_\_

Send orders to: Attn: Corine Harris, AAAS Science Publications, Inc.  
1333 H St., NW, Washington, DC 20005  
(202) 326-6527 • FAX (202) 682-0816

Subtotal \_\_\_\_\_

DC Sales Tax \_\_\_\_\_

Outside the US: \_\_\_\_\_

Add \$3.50 per Calendar air delivery \_\_\_\_\_

Total Enclosed \_\_\_\_\_



# If these are your questions:

■ Is my protein glycosylated?

■ Are the carbohydrates N-linked or O-linked?

■ What types of monosaccharides are present on my glycoprotein?

■ How consistent is the glycosylation between my production lots?

■ How do various culture conditions affect my protein glycosylation?

## This is your **\$795** answer:

### INTRODUCING THREE GLYKO FACE® CARBOHYDRATE ANALYSIS STARTER KITS

■ N-Linked Oligosaccharide Profiling ■ O-Linked Oligosaccharide Profiling ■ Monosaccharide Analysis

Now you can choose one of these easy-to-use Starter Kits for only \$795. Each is guaranteed, and everything you'll need to start getting answers to your carbohydrate questions is included.\*

Additional supplies are as close as your telephone. Just ask for details.




To order your Starter Kit for \$795, phone 1 415 382-6653 or fax 1 415 382 7889. In U.S. only, phone 1 800 33 GLYKO (334 5956). Sorry, one kit per customer only at this special price.

\*Power supply, recirculating chiller, and long-wave UV light box not included. ©1994 Glyko, Inc. FACE is a registered trademark of Glyko, Inc. 81 Digital Drive, Novato, CA 94949

Circle No. 31 on Readers' Service Card

**GLYKO**





# **Biomolecular Interaction Analysis - the shortest route between structure and function**

*Comparison of a structural and  
functional binding domain in human  
growth hormone using real-time BIA*

*B. C. Cunningham and J. A. Wells*

*(J. Mol. Biol. (1993) 234, 554-563)*

From a structural point of view, the complex between human growth hormone (hGH) and the extracellular domain of its receptor (hGHbp) is one of the best characterized hormone-receptor complexes. One molecule of hGH binds sequentially to two molecules of hGHbp, involving separate sites on the hormone molecule and resulting in dimerization of the receptor.

The work of Brian C Cunningham and James A Wells of Genentech Inc, USA, used Biomolecular Interaction Analysis (BIA) to measure the progress of the macromolecular interactions in real-time, allowing the evaluation of the kinetics as well as the affinity of the interactions. By this method it was possible to determine the effect of replacing each of the 30 contact residues in the hGH site 1 structural binding domain with alanine.

The results indicated that only one quarter of the residues account for the major part of the binding energy. Thus the functional binding domain is considerably smaller than the structural binding domain. The results are potentially valuable in the design of hormone analogues for therapeutic purposes, identifying critical residues for the binding interaction and implying that it might be possible to design smaller hormone mimics.

Establishment of a clearly defined experimental situation for the real time BIA studies was aided by the specific immobilization chemistry on the sensor chip. By immobilizing the receptor through a single cysteine residue introduced at a chosen position, the properties of the surface-bound interaction could be closely controlled. This work clearly demonstrates the value of real-time BIA in analysing functional aspects of macromolecular interaction.

With BIAcore™ and BIAcore™, BIA technology is available in a range of automation and sample processing capabilities at a price level to fit every budget.

If you want to see more of what's going on in your biomolecular interactions contact Pharmacia Biosensor today.

See more of what's going on



**Biosensor**

Pharmacia Biosensor

USA 1-800-BIA-2599, France +33-1-3064 3400  
Germany +49-7614 9030, UK +44-727-814 075  
Japan 03-3492-9499, Australia toll free 008 252 265  
Head Office: Sweden +46-18-165376

Circle No. 27 on Readers' Service Card