



AMERICAN  
ASSOCIATION FOR THE  
ADVANCEMENT OF  
SCIENCE

# SCIENCE

27 SEPTEMBER 1996  
VOL. 273 • PAGES 1765-1968

\$7.00





# The Future of Microplate Assays Is Here — 384-well Radioisotopic and Luminescence Counting!

## TopCount HTS Revolutionizes Microplate Counting — Again

- Measure radioisotopic and luminescence assays in both 96- and 384-well microplates.
- Process 12 luminescence samples simultaneously for throughputs of 10,000 samples per hour, or more.
- Stack up to 40 plates or 15,000 samples for unattended counting.
- Interface with robotic systems for complete automation.

## Meeting the Demands of Microplate Counting for Today and for the Future

Assays have been scaled down from test tubes to 24-well microplates and the 96-well format. Now, the next step in assay miniaturization is here with the TopCount HTS. Continue today's assays in 96-well plates while converting assays to the higher density 384-well plates — the screening format of the future. With the TopCount HTS you can do both without risk or compromise.

## Less Reagents, Less Handling, More Samples

The Packard 384-well plate accommodates assay volumes of up to 80  $\mu$ L. This allows practical miniaturization of assays such as scintillation proximity, in-plate binding, liquid scintillation, Cerenkov and luminescence.

## TopCount Harnesses the Power of Luminescence

Packard has pioneered the field of high throughput microplate luminescence. Now up to 12 wells can be measured simultaneously in 96- or 384-well plates for throughputs never before possible. Ultra-sensitive photon counting is combined with innovative reagents, such as LucLite with a five hour half-life, for measuring the luciferase reporter gene. Greater capacity, greater throughput.

To complement sample throughput, TopCount now has an extended 40 plate stacker. Load up to 40 plates with the TopCount ZipLoad stacker, carry them safely to the TopCount, load and count. It's that easy. If your experiments demand accuracy, versatility or high throughput, TopCount has the answer for your laboratory.

Contact Packard at 800-323-1891 (U.S. only) or 203-238-2351 for more details. Or, visit our Web Site at <http://www.packardinst.com>.



Circle No. 21 on Readers' Service Card

Packard Instrument Company, 800 Research Parkway, Meriden, CT 06450 U.S.A.  
Tel: 203-238-2351 Toll Free: 1-800-323-1891 FAX: 203-639-2172  
Web Site: <http://www.packardinst.com> Email: [webmaster@packardinst.com](mailto:webmaster@packardinst.com)




### Packard International Offices:

**Australia**, Mt Waverley 61-3-9543-4266; **Austria**, Vienna 43-1-2702504; **Belgium**, Brussels 32-2-4668210; **Canada**, Ontario 1-800-387-9559; **Central Europe**, Schwadorf, Aus. 43 456 2230 015; **Denmark**, Greve 45-42909023; **France**, Rungis (33) 1 46.86.27.75; **Germany**, Dreieich (49) 6103 385-151; **Italy**, Milano 39-2-33910796/7/8; **Japan**, Tokyo 81-3-3866-5850; **Netherlands**, Groningen 31-50-5413360; Tilburg (013) 5423900; **Russia**, Moscow, 7-095-238-7335; **Switzerland**, Zurich (01) 481 69 44; **United Kingdom**, Pangbourne, Berks (44) 01734 844981.



# Assays for Biological Redox Systems



## Generation of reactive oxygen species and the defense mechanisms against damage by reactive oxygen

### Colorimetric Technology

- ▶ Glutathione
- ▶ Lipid Peroxidation
- ▶ Superoxide Dismutase

### Immunoassay Technology

- ▶ Myeloperoxidase
- ▶ Lactoferrin
- ▶ Plasma Selenium Glutathione Peroxidases

R&D Systems also produces and markets over 600 cell biology reagents, including cytokines and adhesion molecules, corresponding monoclonal and polyclonal antibodies (unlabeled and labeled), human and mouse cytokine ELISAs, as well as genes and probes for cytokines.

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES.**

**U.S.A. and Canada**  
**R&D Systems, Inc.**  
 614 McKinley Place NE  
 Minneapolis, MN 55413, USA.  
 Tel: 612 379-2956  
 Fax: 612 379-6580

**Europe**  
**R&D Systems Europe Ltd.**  
 4-10 The Quadrant, Barton Lane  
 Abingdon, OX14 3YS, UK.  
 Tel: +44 (0)1235 531074  
 Fax: +44 (0)1235 533420

**Germany**  
**R&D Systems GmbH**  
 Borsigstrasse 7  
 65205 Wiesbaden, Germany  
 Tel: +49 06122 90980  
 Fax: +49 06122 909819

**International Distributors – Australia:** (61) 02-521-2177. **Austria:** (43) 02 292 35 27. **Chile:** (56) 2-671-9369.  
**Israel:** (972) 9 764 8787. **Italy:** (39) 02 25 75 377. **Korea:** (850) 82-2-569-0781. **Mexico:** (52) 5-612-0085.  
**Spain & Portugal:** (34) 01 448 84 86 or 03 456 97 06. **Taiwan:** (866) 2-368-3600.

**International Freefone Numbers – Belgique/België:** 078 11 04 68. **Danmark:** 80 01 85 92. **Deutschland:** 013011 0169.  
**France:** 05 90 72 49. **Nederland:** 060-225607. **Norge:** 800 1103. **Sverige:** 02079 31 49. **Switzerland:** 155 2482.

Circle No. 31 on Readers' Service Card


**R&D**  
**SYSTEMS**  
**1-800-343-7475**

[www.rndsystems.com](http://www.rndsystems.com)



# BioDoc II™ The ultimate video documentation system

**SPECIAL  
KEYBOARD  
FOR EASY  
HANDLING**



**SUPERIOR SENSITIVITY** by electronic control of the 'on chip' integrating CCD varying exposure time, background, brilliance, gain and edge enhancement:  $4 \times 10^{-6}$  lux and 63 dB signal-to-noise ratio.

**UNSURPASSED USER FRIENDLINESS** by completely new design, easy-to-cleanse keyboard with easy-to-read and self explaining touch keys.

**UNIQUE** specially developed auto-exposure function - six preset protocols and ten user defined protocols guaranteeing extremely fast handling of sensitive specimens.

Data storage on either floppy, hard disc or the network in both BMP or TIFF format.

Freeze function allows for additional improvement of the digitised images in direct or inverted presentation.

Printing of live, frozen or opened images on high resolution digital thermoprinter.

High resolution VGA colour display and execution of controlling software under Windows.

Compatible with the Biometra analysis softwares BioDocAnalyse and ScanPack 2.0.

## Biometra®

Circle No. 14 on Readers' Service Card

Germany	Biometra GmbH	Tel: 0551/50 68 60	Fax: 0551/50 68 666	e-mail: <a href="mailto:info@biometra.co.uk">info@biometra.co.uk</a>
UK	Biometra Ltd	Tel: 01622 678872	Fax: 01622 752774	e-mail: <a href="mailto:sales@biometra.co.uk">sales@biometra.co.uk</a>
US	Biometra Inc	Tel: 1-800 932 7250	Fax: 813 282 1936	Internet: <a href="mailto:biometra@gate.net">biometra@gate.net</a>

# We're proving that positive thinking can cure disease.

The first step in accomplishing anything is the belief that it can be done. We believe there are cures for today's "incurable" diseases and we're working fast to find them. At Pharmacia & Upjohn we share this belief with more than 30,000 colleagues worldwide.


Along the way, we're finding better treatments for diseases like cancer, AIDS, diabetes, Parkinson's, growth disorders and arthritis.

Believing there is a cure is the first step in finding it. By sharing that belief you too can join in the fight. That's the power of positive thinking. That's the power of Pharmacia & Upjohn.



**Pharmacia  
& Upjohn**





**1795**

Microbes star in biofilms



**1850**

How old is Dominican  
amber?

## NEWS & COMMENT

Whose Genome Is It, Anyway? 1788

France: Research Is Spared the Budget Ax 1790

Germany: No Respite for International Labs 1790

Computer Networks: Explosive Growth Helps Japan Take Its Place on Global Team  
New Networks Aren't Good Enough 1791

A Plea to Protect Threatened Collections 1792

Congressional Scorecard Sparks Furor 1793

Undergraduate Education: In Boston, a Revolutionary Experience 1794

## RESEARCH NEWS

Biofilms Invade Microbiology 1795  
Linking Microbiology and Engineering 1796

AIDS Research: Receptor Mutations Help Slow Disease Progression  
Quick Work Draws Scientific Praise, Colleagues' Complaints 1797  
1798

Hubble Spectrograph Takes a Hard Look in the Ultraviolet 1799

Plastics May Add New Colors to Lasers' Light Show 1800

How the Songbird Makes His Song 1801

Politics Meets Technology in Next Collider 1802

## PERSPECTIVES

Flux Lattice Melting 1811  
D. J. Bishop

Is Neural Noise Just a Nuisance? 1812  
D. Ferster

A Parallel Spliceosome 1813  
T. W. Nilsen

## ARTICLES

Noble Gases and Earth's Accretion 1814  
C. L. Harper Jr. and S. B. Jacobsen

Copper Sources, Metal Production, and Metals Trade in Late Postclassic Mesoamerica 1819  
D. Hosler and A. Macfarlane

## DEPARTMENTS

THIS WEEK IN SCIENCE 1773

EDITORIAL 1779  
Cooperation, Competition, and Science Policy  
E. Bloch

LETTERS 1781  
Interleukin-1 Research: T. P. Hopp • On Sociological Biophobia: A. Abbott; J. H. Gagnon • Morality Play: L. Z. Benet; W. R. Debler; J. R. Thomen; E. B. Hook; C. Eckert; D. S. Zinberg

SCIENCESCOPE 1787

RANDOM SAMPLES 1805  
Health of Native Americans • Rabies Death for

Poe? • Micromanaging Atoms

BOOK REVIEWS 1806  
*Frontiers of Illusion*, reviewed by D. H. Guston • *The Evolution and Extinction of the Dinosaurs and Dinosaur Extinction and the End of an Era*, D. A. Russell • *The Great Dinosaur Extinction Controversy* • Vignettes • Books Received

AAAS NEWS & NOTES 1880

PRODUCTS & MATERIALS 1882

QUARTERLY INDEX 1884

## AAAS Board of Directors

Rita R. Colwell  
Retiring President, Chairman  
Jane Lubchenco  
President  
Mildred S. Dresselhaus  
President-elect

Sheila Jasanoff  
William A. Lester Jr.  
Simon A. Levin  
Marcia C. Linn  
Michael J. Novacek  
Anna C. Roosevelt  
Jean E. Taylor  
Nancy S. Wexler

William T. Golden  
Treasurer  
Richard S. Nicholson  
Executive Officer

SCIENCE (ISSN 0036-8075) is published weekly on Friday, except the last week in December, by the American Association for the Advancement of Science, 1200 New York Avenue, NW, Washington, DC 20005. Periodicals Mail postage (publication No. 484460) paid at Washington, DC, and additional mailing offices. Copyright © 1996 by the American Association for the Advancement of Science. The title SCIENCE is a registered trademark of the AAAS. Domestic individual membership and subscription (51 issues): \$102 (\$55 allocated to subscription). Domestic institutional subscription (51 issues): \$250. Foreign postage extra: Mexico, Caribbean (surface mail) \$55; other countries (air assist delivery) \$90. First class, airmail, student, and emeritus rates on request. Canadian rates with GST available upon request, GST #1254 88122. Printed in the U.S.A.



## COVER

Methane clathrate hydrate, an icy compound capable of sustaining its own combustion, displays intriguing formation kinetics as well as unusual plastic flow properties. Although naturally occurring methane hydrate

may harbor the largest untapped reservoir of natural gas on Earth, the material properties of this ephemeral compound are not yet well understood. See page 1843. [Photo: John Pinkston and Laura Stern]



## RESEARCH ARTICLE

- Highly Diverged U4 and U6 Small Nuclear RNAs Required for Splicing Rare AT-AC Introns 1824  
W.-Y. Tarn and J. A. Steitz

## REPORTS

- Semiconducting Polymers: A New Class of Solid-State Laser Materials 1833  
F. Hide, M. A. Díaz-García, B. J. Schwartz, M. R. Andersson, Q. Pei, A. J. Heeger

- Nanorod-Superconductor Composites: A Pathway to Materials with High Critical Current Densities 1836  
P. Yang and C. M. Lieber

- Methane Hydrate and Free Gas on the Blake Ridge from Vertical Seismic Profiling 1840  
W. S. Holbrook, H. Hoskins, W. T. Wood, R. A. Stephen, D. Lizarralde, Leg 164 Science Party

- Peculiarities of Methane Clathrate Hydrate Formation and Solid-State Deformation, Including Possible Superheating of Water Ice 1843  
L. A. Stern, S. H. Kirby, W. B. Durham

- A Trinuclear Intermediate in the Copper-Mediated Reduction of O<sub>2</sub>: Four Electrons from Three Coppers 1848  
A. P. Cole, D. E. Root, P. Mukherjee, E. I. Solomon, T. D. P. Stack

- Age and Paleogeographical Origin of Dominican Amber 1850  
M. A. Iturralde-Vinent and R. D. E. MacPhee

- Initiation of Runaway Cell Death in an *Arabidopsis* Mutant by Extracellular Superoxide 1853  
T. Jabs, R. A. Dietrich, J. L. Dangl

- Genetic Restriction of HIV-1 Infection and Progression to AIDS by a Deletion Allele of the *CKR5* Structural Gene 1856  
M. Dean, M. Carrington, C. Winkler, G. A. Huttley, M. W. Smith, R. Allikmets, J. J. Goedert, S. P. Buchbinder, E. Vittinghoff, E. Gomperts, S. Donfield, D. Vlahov, R. Kaslow, A. Saah, C. Rinaldo, R. Detels, Hemophilia Growth and Development Study, Multicenter AIDS Cohort Study, Multicenter Hemophilia Cohort Study, San Francisco City Cohort, ALIVE Study, S. J. O'Brien

- CD40 Ligand-Dependent T Cell Activation: Requirement of B7-CD28 Signaling Through CD40 1862  
Y. Yang and J. M. Wilson

- Requirement for CD40 Ligand in Costimulation Induction, T Cell Activation, and Experimental Allergic Encephalomyelitis 1864  
I. S. Grewal, H. G. Foellmer, K. D. Grewal, J. Xu, F. Hardardottir, J. L. Baron, C. A. Janeway Jr., R. A. Flavell

- Dynamics of Ongoing Activity: Explanation of the Large Variability in Evoked Cortical Responses 1868  
A. Arieli, A. Sterkin, A. Grinvald, A. Aertsen

- Temporal Hierarchical Control of Singing in Birds 1871  
A. C. Yu and D. Margoliash

- Conformational States of the Nuclear Pore Complex Induced by Depletion of Nuclear Ca<sup>2+</sup> Stores 1875  
C. Perez-Terzic, J. Pyle, M. Jaconi, L. Stehno-Bittel, D. E. Clapham


## TECHNICAL COMMENTS

- Sm-Nd Isotopic Data and Earth's Evolution 1878  
S. Moorbath and M. J. Whitehouse; S. A. Bowring and T. Housh



1875

Nuclear pore complex



## ■ Indicates accompanying feature

**Change of address:** allow 4 weeks, giving old and new addresses and 8-digit account number.  
**Postmaster:** Send change of address to *Science*, P.O. Box 1811, Danbury, CT 06813-1811.  
**Single copy sales:** \$7.00 per issue prepaid includes surface postage; bulk rates on request.  
**Authorization to photocopy:** material for internal or personal use under circumstances not falling within the fair use provisions of the Copyright Act is granted by AAAS to libraries and other users registered with the Copyright Clearance Center (CCC) Transactional Reporting Service, provided that \$4.00 per article is paid directly to CCC, 222 Rosewood Drive, Danvers, MA 01923. The identification code for *Science* is 0036-8075/93 \$4.00. *Science* is indexed in the *Reader's Guide to Periodical Literature* and in several specialized indexes.

## On the Web

You don't have to wait anymore!

Grab any scientific paper on the day it hits the mail:

<http://www.sciencemag.org/science/home/faq-text.shtml>



# How to Find the Perfect Transfection Tool.

If you're looking for a transfection tool that provides maximum transfection efficiency, look no further. The PerFect Lipid™ Transfection Kit allows you to test a panel of eight unique lipids so you can quickly and easily find the one that yields the highest efficiency in your cell line.

## Maximize Efficiency

The efficiency of lipid-mediated transfection varies greatly depending on the cell line and lipid used (Figure 1).




Figure 1: Percent blue cells produced by transfection of CHO cells with a *lacZ* control vector using PerFect Lipids™ (pFex 1 - pFex 8).

That's why we developed the PerFect Lipid™ Transfection Kit. Now finding the right lipid for your cell line is easy. Simply transfect the supplied  $\beta$ -galactosidase reporter plasmid with each of the eight PerFect Lipids™ (pFex™1-8) and assay for the percent of blue cells (Figure 2).




Figure 2: Transfected CHO cells expressing  $\beta$ -galactosidase, stained with X-gal.

## Minimize Cost

Once you've found the perfect lipid you can continue to achieve optimum transfection efficiencies. Each PerFect Lipid™ is available in individual pack sizes that are convenient and economical to use.

## Tested for Results

To make it even easier for you to choose the right lipid, we've tested the PerFect Lipids™ in a variety of common cell lines. To see for yourself, check out our web site (<http://www.invitrogen.com>) or call our Technical Services Department.

For perfect transfections with every cell line, there's no better choice than the PerFect Lipid™ Transfection Kit and PerFect Lipids™. For more information or to place your order, call today.



Telephone (619) 597-6200  
Fax (619) 597-6201  
<http://www.invitrogen.com>  
1-800-955-6288

Circle No. 19 on Readers' Service Card

## EUROPEAN HEADQUARTERS:

**Invitrogen BV**  
De Schelp 12, 9351 NV Leek  
The Netherlands  
Tel: +31 (0) 594 515 175  
Fax: +31 (0) 594 515 312  
Email: [tech\\_service@invitrogen.nest.nl](mailto:tech_service@invitrogen.nest.nl)

## Distributors:

Australia 03 9562 6888	Korea 822 569 6902
Hong Kong 886 2 381 0844	Singapore 64 482 3888
Israel 02 524 447	or 29 29 783
Italy 02 33 91 00 72	Spain 03 450 2601
Japan 03 5684 1616	Taiwan 886 2 381 0844

From all other countries, please contact our European headquarters at +31 (0) 594 515 175  
PerFect Lipids™ are manufactured by Genta, Inc. and licensed to Invitrogen Corporation.


## TOLL FREE TELEPHONE NUMBERS:

Belgium 0800 111 73	The Netherlands 06 022 8848
Denmark 800 188 67	Norway 800 113 70
Finland 990 31 800 5345	Sweden 020 795 369
France 19 31 800 5345	Switzerland 155 1966
Germany 0130 8100 43	UK 0800 96 61 93



## Copper(III) chemistry

Several enzymes that can reduce  $O_2$  to  $H_2O$  contain a trinuclear copper core. Cole *et al.* (p. 1848) have synthesized a small trinuclear copper com-



plex that breaks the O—O bond and does so with this unusual 3:1 Cu: $O_2$  stoichiometry. The four-electron oxidation produces a bridging copper atom in its unusual trivalent state.

## Deliberate defects

High-temperature superconductors (HTSCs) can carry high current densities only if the crystalline grains are well aligned and if flux lines are pinned by defects with the material. Often such defects are introduced after the material has been grown, but Yang and Lieber (p. 1836) report a method for incorporating nanometer-scale magnesium oxide columns in the bismuth HTSCs. High critical current densities comparable to those seen for the best polycrystalline thick films were observed in materials where the rods were oriented perpendicular to the copper oxide planes.

## Methane hydrates

Abundant methane is stored as a solid hydrate on continental shelves and also may also be present in icy moons. On Earth, the phase transition between the hydrate and free methane plus water is a prominent seismic

## Receptor mutations and HIV disease progression

Chemokines receptors such as CKR5 have recently been identified that act as co-receptors with CD4 for the entry of HIV-1 (human immunodeficiency virus—type 1) into human cells, and defects in these receptors could lead to resistance against viral infection and disease progression. Dean *et al.* (p. 1856; see the news story by Cohen, p. 1797) mapped the CKR5 structural gene to human chromosome 3 and identified a 32-base pair deletion that occurs in 10% of the Caucasian population. In examining 1955 individuals at high risk for HIV-1 infection, they found that 17 individuals homozygous for this allele were HIV-1 negative, and that disease progression is slower in individuals heterozygous for this allele than in individuals who lack this allele.

reflector. Estimation of the abundance of trapped methane requires knowledge of the meaning of the seismic signal and the phase relations. Experiments by Stern *et al.* (p. 1843; see cover) that examined the phase relations of methane clathrate hydrate and suggest that superheating of the solid may be possible. Holbrook *et al.* (p. 1840) compared drill core records and seismic images through the methane hydrate layer off of South Carolina and conclude that less methane hydrate may be there than previously thought.

## Formative stages

A simple two-stage model of Earth's formation may explain how the Earth can presently have an oxidized mantle that is not in equilibrium with its reduced iron-rich core. Harper and Jacobsen (p. 1814) used noble gas isotopic compositions of basalts and primitive meteorites to develop an accretion model for the Earth. They suggest that the mantle and core were initially equilibrated under the reducing conditions of a massive hydrogen-rich proto-atmosphere that blanketed the proto-Earth's magma ocean and

allowed the incorporation of the light noble gases. During and after the removal of the proto-atmosphere by the sun, the heavier noble gases and other components continued to accrete and oxidize the mantle.

## Amber range of age

Amber from the Dominican Republic and related sites in the Caribbean has provided key specimens for studying the evolution of many trapped plant and animal species and may be the only viable source for extraction of very ancient DNA, but interpretations have been hampered because the age of the amber has been poorly known; estimates range back to more than 65 million years ago. Iturralde-Vinent and MacPhee (p. 1850) reconstructed the stratigraphy of the region and conclude that all of the key amber-bearing rocks were deposited in a single sedimentary basin 15 to 20 million years ago.

## Overkill in plants

Plant cells can resist pathogens by inducing cell death locally around a site of infection, but in certain mutants cell death is

unchecked and can engulf an entire leaf. Jabs *et al.* (p. 1853) studied such mutants in *Arabidopsis* and found that superoxide, a reactive oxygen intermediate, alone could initiate this process. Signals that normally would suppress superoxide production appear to be blocked, which then leads to amplification of the superoxide signal.

## Effects of CD40L

Two reports illustrate the important role of the ligand to CD40 (CD40L) in stimulating the immune response. Grewal *et al.* (p. 1864) show that induction of an autoimmune disease, experimental allergic encephalomyelitis, could not be induced by antigen injection if they had no CD40L. Yang and Wilson (p. 1862) show that CD40L-deficient mice also did not have the expected adverse reaction to a gene therapy vector. These two different approaches show that CD40L works through induction of B7 on antigen-presenting cells, which increases the costimulation signal that the T cells need to increase their response.

## Nuclear pore control


Nuclear pores punctuate the nuclear envelope and contain a channel through which molecules enter and leave the nucleus. When the nuclear envelope is depleted of calcium, nuclear import is inhibited. Perez-Terzic *et al.* (p. 1875) have visualized conformational changes in the structure of the nuclear pore itself which explain why, when calcium is depleted from the nuclear envelope, the diameter of the channel through the pore shrinks.



SEE

WHAT YOU'VE BEEN MISSING.

CAPTURE IT, ENHANCE IT,  
SHARE IT.



Pixera Professional System/John Bagley Photography



**Pixera Studio:** Compose your picture with a through-the-lens motion viewfinder. Manipulate and enhance the acquired image and regions of interest with a wide variety of built-in high-performance tools, filters and effects. Compress it and store it in JPEG format or convert it into TIFF, PICT, BMP, or DIB.



**Pixera Album:** Store, catalog and manage an unlimited number of images, sound files, even video clips. Organize, sort, query and filter stored objects using a powerful database tool set. Capture, edit or enhance new images with one-button access to the full capability of Pixera Studio.

Sometimes, it's not what you see but what you don't that's critical. You need the best image you can get your eyes on. And no digital camera can give you better image quality at a lower price than a Pixera Professional Visual Communication System. It provides the resolution and image quality of systems costing many times more. Yet you get everything you need at one low price -- including the first digital camera system with million-pixel resolution and a full suite of application software.

The Pixera Professional captures two- or three-dimensional images and can connect to any optical microscope with a standard C-mount interface. Choose from PCI or PCMCIA versions for Windows 95 or Macintosh.

For more information, call toll-free 1-888-4PIXERA or visit our web page at [www.pixera.com](http://www.pixera.com).



**pixera**  
visual SYSTEM communication

© 1996 Pixera Corporation, 10211 Bubb Road, Cupertino, CA 95014  
408/342-0110, [sales@pixera.com](mailto:sales@pixera.com)  
All trademarks are the property of their respective holders.

Circle No. 26 on Readers' Service Card



# STRATAGENE

## XL2-Blue

### ULTRACOMPETENT CELLS


DO YOU WANT TO MAKE LARGE REPRESENTATIVE PLASMID LIBRARIES?

DO YOU EVER HAVE TROUBLE GETTING THE CLONES YOU WANT?

DO YOU EVER HUNT FOR RARE MESSAGES?

DO YOU DO PCR CLONING, BLUNT-END LIGATIONS OR MULTIPLE FRAGMENT LIGATIONS?


THEN YOU NEED XL2-BLUE ULTRACOMPETENT CELLS!



The Most Competent *E. Coli* in the World

*XL2-Blue Ultracompetent Cells  
are 5 times more efficient than  
DH5 MAX cells, so you can  
get 5 times more colonies!*

*When others are MAXed out,  
XL2-Blue cells are still transforming.*



Circle No. 30 on Readers' Service Card

USA  
Corporate Headquarters  
(800) 424-5444  
(619) 535-5400  
Fax: (619) 535-0045

GERMANY  
Stratagene GmbH  
Telephone: (06221) 400634  
Telefax: (06221) 400639

UNITED KINGDOM  
Stratagene Ltd.  
Telephone: (01223) 420955  
Telefax: (01223) 420234

SWITZERLAND  
Stratagene GmbH  
Telephone: (01) 3641106  
Telefax: (01) 3657707

Internet Mail  
techservices@stratagene.com




**STRATAGENE**



# At Jandel, you supply the data, we provide the solution.

Whether you're creating technical graphs for publication or searching for the perfect equation to describe empirical data, Jandel offers time-saving software tools to help you do better research, accurately and more easily. That's why over 150,000 scientists and engineers count on Jandel products for their research solutions.



## GRAPHING SigmaPlot®

Award-winning SigmaPlot is used by more scientists and engineers than any other technical graphing program to visualize, analyze and publish their data. Choose from a complete set of 2D and 3D graph types and then customize every detail to create the exact graph you need. Analyze your data with SigmaPlot's extensive mathematical transforms and summary statistics. SigmaPlot automatically fits and graphs polynomial regressions up to the 10th order and 95% or 99% confidence and prediction intervals. Microsoft Office and OLE 2 compatibility means SigmaPlot is even more versatile – you can even use an Excel worksheet inside SigmaPlot! From data import to final output SigmaPlot gives you graphing flexibility every step of the way.





## DATA ANALYSIS

### SigmaStat®, TableCurve™ & PeakFit™

Analyzing and validating your results is a crucial part of the research process. SigmaStat offers powerful, comprehensive statistical procedures and an "expert system" design that features automatic assumption checking, messy data handling, reports with an explanation of your results, presentation-quality graphs and an Advisor Wizard that recommends the correct test to run. SigmaStat gives you results you can trust.

TableCurve 2D and 3D combine powerful data fitting capabilities with the ability to find the ideal equation to describe empirical data. Extremely fast and easy to use, TableCurve 2D and 3D automatically locate best fit candidate equations for your X,Y and X,Y,Z data from a built-in set of functions. Or enter your own functions to be fitted and ranked along with the built-in equations. PeakFit is also available to separate and fit up to 100 overlapping, difficult to resolve peaks on your spectroscopy, chromatography and electrophoresis data.



## IMAGE ANALYSIS

### SigmaScan™, SigmaScan Pro & SigmaGel™

SigmaScan and SigmaScan Pro are two powerful image measurement and analysis software tools that allow you to analyze your digital images quickly and accurately. Just scan in your image, enhance, calibrate or annotate your image, measure it using any of 44 different measurement features, analyze your data and output the information. SigmaGel is also available to perform quantitative gel and spot analysis on your own PC.

Jandel – The right solution for Scientists and Engineers. Since 1982, scientists and engineers worldwide have been using Jandel products to complete their work.

**CALL 1-800-4-JANDEL (1-800-452-6335)  
AND PUT JANDEL PRODUCTS TO WORK FOR YOU TODAY!**

OR VISIT OUR WEBSITE AT [HTTP://WWW.JANDEL.COM](http://www.jandel.com)

Jandel Scientific Software 2591 Kerner Blvd. San Rafael, CA 94901 Phone: 1-800-4-JANDEL, (415) 453-6700 FAX (415) 453-7769  
e-mail: [sales@jandel.com](mailto:sales@jandel.com) Internet: <http://www.jandel.com> Jandel Scientific Software has dealers and distributors worldwide. Call, fax or email Jandel U.S. for a representative in your country.

All company and/or product names are trademarks of their respective companies. ©1996 Jandel Corporation.

Circle No. 33 on Readers' Service Card

**Jandel**  
SCIENTIFIC SOFTWARE



#### MEMBERSHIP/CIRCULATION

**Director:** Michael Spinella  
**Deputy Director:** Marlene Zendell  
**Member Services:** Donald Crowl, *Manager*; Mary Curry, *Supervisor*; Pat Butler, Laurie Baker, *Representatives*  
**Marketing:** Dee Valencia, *Manager*; Jane Pennington, *Europe Manager*; Hilary Baar, *Assistant Manager*  
**Research:** Renuka Chander, *Manager*  
**Business and Finance:** Robert Smariga, *Manager*; Nina Araujo de Kobes, Kimberly Parker, *Coordinators*

#### ADVERTISING AND FINANCE

**Associate Publisher:** Beth Rosner  
**Advertising Sales Manager:** Susan A. Meredith  
**Recruitment Advertising Manager:** Janis Crowley  
**Business Manager:** Deborah Rivera-Wienhold  
**Finance:** Randy Yi, *Senior Analyst*; Connie Dong, *Financial Analyst*  
**Marketing:** John Meyers, *Manager*; Allison Pritchard, *Associate*  
**Electronic Media:** David Gillikin, *Manager*; Pamela Sams, *Production Associate*  
**Product Advertising:** Carol Maddox, *Traffic Manager*; Natalie Britt, *Sales Associate*  
**Product Advertising Sales:** **East Coast/E. Canada:** Richard Teeling, 201-904-9774, FAX 201-904-9701 • **Midwest/Southeast:** Elizabeth Mosko, 312-665-1150, FAX 312-665-2129 • **West Coast/W. Canada:** Neil Boylan, 415-673-9265, FAX 415-673-9267 • **UK, Scandinavia, France, Italy, Belgium, Netherlands:** Andrew Davies, (44) 1-457-838-519, FAX (44) 1-457-838-898 • **Germany/Switzerland/Austria:** Tracey Peers, (44) 1-260-297-530, FAX (44) 1-260-271-022 • **Japan:** Mashy Yoshikawa, (81) 3-3235-5961, FAX (81) 3-3235-5852  
**Recruitment Advertising:** Terri Seiter Azie, *Sales and Production Operations Manager*; Celeste Miller, *Sales Supervisor*; Eric Banks, Bethany Ritchey, *Sales*; Debbie Cummings, *European Sales Manager*; Wendy Green, *Production Associate*; Nicole Robinson, *Advertising Assistant*  
**Recruitment Advertising Sales:** **US:** Janis Crowley, 202-326-6532, FAX 202-289-6742 • **Europe:** Debbie Cummings, (44) 1223-302067, FAX (44) 1223-576208 • **Australia/New Zealand:** Keith Sandell, (61) 02-922-2977, FAX (61) 02-922-1100 • **Japan:** Mashy Yoshikawa, (81) 3-3235-5961, FAX (81) 3-3235-5852  
**Assistant to Associate Publisher:** Nancy Hicks  
**Permissions:** Lincoln Richman  
**Exhibits Coordinator:** Arlene Ennis  
 Send materials to *Science* Advertising, 1200 New York Avenue, NW, Washington, DC 20005.

#### SCIENCE EDITORIAL BOARD

Charles J. Arntzen  
 David Baltimore  
 J. Michael Bishop  
 William F. Brinkman  
 E. Margaret Burbidge  
 Pierre-Gilles de Gennes  
 Joseph L. Goldstein  
 Mary L. Good  
 Harry B. Gray

John J. Hopfield  
 F. Clark Howell  
 Paul A. Marks  
 Yasutomi Nishizuka  
 Helen M. Ranney  
 Bengt Samuelsson  
 Robert M. Solow  
 Edward C. Stone  
 James D. Watson  
 Richard N. Zare

■ Published by the **American Association for the Advancement of Science (AAAS)**, *Science* serves its readers as a forum for the presentation and discussion of important issues related to the advancement of science, including the presentation of minority or conflicting points of view, rather than by publishing only material on which a consensus has been reached. Accordingly, all articles published in *Science*—including editorials, news and comment, and book reviews—are signed and reflect the individual views of the authors and not official points of view adopted by the AAAS or the institutions with which the authors are affiliated.

■ The American Association for the Advancement of Science was founded in 1848 and incorporated in 1874. Its objectives are to further the work of scientists, to facilitate cooperation among them, to foster scientific freedom and responsibility, to improve the effectiveness of science in the promotion of human welfare, to advance education in science, and to increase public understanding and appreciation of the importance and promise of the methods of science in human progress.

## INFORMATION RESOURCES

#### SUBSCRIPTION SERVICES

For change of address, missing issues, new orders and renewals, and payment questions, please contact AAAS at Danbury, CT: 800-731-4939 or Washington, DC: 202-326-6417, FAX 202-842-1065. Mailing addresses: AAAS, P.O. Box 1811, Danbury, CT 06813 or AAAS Member Services, 1200 New York Avenue, NW, Washington, DC 20005 • **Other AAAS Programs:** 202-326-6400

#### REPRINTS & PERMISSION

**Reprints:** Ordering/Billing, 800-407-9191; Corrections, 202-326-6501 • **Permissions:** 202-326-7074, FAX 202-682-0816

#### INTERNET ADDRESSES

science\_editors@aaas.org (for general editorial queries); science\_news@aaas.org (for news queries); science\_letters@aaas.org (for letters to the editor); science\_reviews@aaas.org (for returning manuscript reviews); science@science-int.co.uk (for the Europe Office); membership@aaas.org (for member services); science\_classifieds@aaas.org (for submitting classified advertisements); science\_advertising@aaas.org (for product advertising)

#### INFORMATION FOR CONTRIBUTORS

See pages 93–94 of the 5 January 1996 issue or

access <http://www.sciencemag.org/science/home/con-info.shtml>.

#### EDITORIAL & NEWS CONTACTS

##### North America

**Address:** 1200 New York Avenue, NW, Washington, DC 20005  
**Editorial:** 202-326-6501, FAX 202-289-7562  
**News:** 202-326-6500, FAX 202-371-9227 • **Bureaus:** **Berkeley, CA:** 510-841-1154, FAX 510-841-6339, **San Diego, CA:** 619-942-3252, FAX 619-942-4979, **Chicago, IL:** 312-360-1227, FAX 312-360-0537, **Boston, MA:** 617-566-7137, FAX 617-734-8088

##### Europe

**Headquarters:** 14 George IV Street, Cambridge, UK CB2 1HH; (44) 1223-302067, FAX (44) 1223-302068  
**Paris Correspondent:** (33) 1-49-29-09-01, FAX (33) 1-49-29-09-00

##### Asia

**Japan Office:** Carl Kay, Esaka-cho 5-chome 11-10, Suita-shi, Osaka 564 Japan; (81) 6-368-1925, FAX (81) 6-368-6905; science@magical.egg.or.jp • **News Bureau:** (81) 3-3335-9925, FAX (81) 3-3335-4898 • **China Office:** Hao Xin, (86)10-6255-9478; haoxin@sun.ihep.ac.cn

#### BOARD OF REVIEWING EDITORS

Frederick W. Alt  
*Children's Hospital, Boston*  
 Don L. Anderson  
*California Institute of Technology*  
 Michael Ashburner  
*Univ. of Cambridge*  
 Stephen J. Benkovic  
*Pennsylvania State Univ.*  
 Alan Bernstein  
*Mount Sinai Hospital, Toronto*  
 David E. Bloom  
*Harvard Institute for International Development*  
 Piet Borst  
*The Netherlands Cancer Institute*  
 Henry R. Bourne  
*Univ. of California, San Francisco*  
 Michael S. Brown  
*Univ. of Texas Southwestern Medical Center*  
 James J. Bull  
*Univ. of Texas at Austin*  
 Kathryn Calame  
*Columbia Univ. College of Physicians & Surgeons*  
 Dennis W. Choi  
*Washington Univ. School of Medicine, St. Louis*  
 David Clapham  
*Mayo Foundation*  
 Adrienne E. Clarke  
*Univ. of Melbourne, Parkville*  
 John M. Coffin  
*Tufts Univ. School of Medicine*  
 F. Fleming Crim  
*Univ. of Wisconsin, Madison*  
 Paul J. Crutzen  
*Max-Planck-Institut für Chemie*  
 James E. Dahlberg  
*Univ. of Wisconsin Medical School, Madison*  
 Robert Desimone  
*National Institute of Mental Health, NIH*  
 Paul T. Englund  
*Johns Hopkins Univ. School of Medicine*

G. Ertl  
*Max-Planck-Gesellschaft, Boston*  
 Richard G. Fairbanks  
*Lamont-Doherty Earth Observatory*  
 Robert E. Fay  
*U.S. Bureau of the Census*  
 Douglas T. Fearon  
*Univ. of Cambridge*  
 Harry A. Fozzard  
*The Univ. of Chicago*  
 Roger I. M. Glass  
*Centers for Disease Control*  
 Stephen P. Goff  
*Columbia Univ. College of Physicians & Surgeons*  
 Peter N. Goodfellow  
*SmithKline Beecham, UK*  
 Peter Gruss  
*Max Planck Institute of Biophysical Chemistry*  
 Philip C. Hanawalt  
*Stanford Univ.*  
 Paul Harvey  
*Univ. of Oxford*  
 M. P. Hassell  
*Imperial College at Silwood Park*  
 Nobutaka Hirokawa  
*Univ. of Tokyo*  
 Tomas Hökfelt  
*Karolinska Institutet*  
 Tasuku Honjo  
*Kyoto Univ.*  
 Susan D. Iversen  
*Univ. of Oxford*  
 Eric F. Johnson  
*The Scripps Research Institute*  
 Elliott Kieff  
*Harvard Univ.*  
 Stephen M. Kosslyn  
*Harvard Univ.*  
 Michael LaBarbera  
*The Univ. of Chicago*  
 Nicole Le Douarin  
*Institut d'Embryologie Cellulaire et Moléculaire du CNRS*  
 Charles S. Levings III  
*North Carolina State Univ., Raleigh*  
 Harvey F. Lodish  
*Whitehead Institute for Biomedical Research*

Richard Losick  
*Harvard Univ.*  
 Reinhard Lührmann  
*Institut für Molekularbiologie und Tumorforschung der Philipps-Universität*  
 Ruth Lynden-Bell  
*Queen's Univ., Belfast*  
 Seth Marder  
*California Institute of Technology*  
 Diane Mathis  
*Institut de Chimie Biologique, Strasbourg*  
 Anthony R. Means  
*Duke Univ. Medical Center*  
 Shigetada Nakanishi  
*Kyoto Univ.*  
 Kim Nasmyth  
*Research Institute of Molecular Pathology, Vienna*  
 Roger A. Nicoll  
*Univ. of California, San Francisco*  
 Staffan Normark  
*Swedish Institute for Infectious Disease Control*  
 Stuart L. Pimm  
*The Univ. of Tennessee, Knoxville*  
 Yeshayau Pocker  
*Univ. of Washington, Seattle*  
 Ralph S. Quatrano  
*Univ. of North Carolina, Chapel Hill*  
 Martin Raff  
*Univ. College London*  
 V. Ramanathan  
*Scripps Institute of Oceanography*  
 Douglas C. Rees  
*California Institute of Technology*  
 T. M. Rice  
*ETH-Hönggerberg, Zürich*  
 David C. Rubie  
*Universität Bayreuth*  
 Erkki Ruoslahti  
*The Burnham Institute, CA*  
 Gottfried Schatz  
*Biozentrum, Basel*  
 Jozef Schell  
*Max-Planck-Institut für Zuchtungsforschung*

Ronald H. Schwartz  
*National Institute of Allergy and Infectious Diseases, NIH*  
 Terrence J. Sejnowski  
*Salk Institute*  
 Christopher R. Somerville  
*Carnegie Institute of Washington*  
 Thomas A. Steitz  
*Yale Univ.*  
 Michael P. Stryker  
*Univ. of California, San Francisco*  
 Tomoyuki Takahashi  
*Univ. of Tokyo*  
 Masatoshi Takeichi  
*Kyoto Univ.*  
 Keiji Tanaka  
*RIKEN Institute*  
 David Tilman  
*Univ. of Minnesota, St. Paul*  
 Robert T. N. Tjian  
*Univ. of California, Berkeley*  
 Yoshinori Tokura  
*Univ. of Tokyo*  
 Emil R. Unanue  
*Washington Univ. School of Medicine, St. Louis*  
 Geerat J. Vermeij  
*Univ. of California, Davis*  
 Bert Vogelstein  
*Johns Hopkins Oncology Center*  
 Arthur Weiss  
*Univ. of California, San Francisco*  
 Zena Werb  
*Univ. of California, San Francisco*  
 George M. Whitesides  
*Harvard Univ.*  
 Owen N. Witte  
*Univ. of California, Los Angeles*



# Filtration Problems?



## The Solution Is Clear — Millipore® Filters From Techware



When it comes to filtration problems, there's no clearer solution than Millipore® Filtration Products from Sigma-Aldrich Techware.

One phone call is all it takes to gain immediate access to Millipore's comprehensive family of ready-to-use, disposable filter devices...as well as over 10,000 other brand-name laboratory equipment and supply products.


In addition to this wide selection of top quality products, Sigma-Aldrich Techware backs each purchase with

unequaled service and support, including:

- Same day shipping
- Comprehensive technical support
- Competitive pricing

So for solutions that are perfectly clear, rely on Millipore filtration products from Sigma-Aldrich Techware — your single, most convenient source for brand name laboratory supplies and equipment.

Circle the number below for your free Filtration Product Guide.



Call collect: 314-771-5750, 800-325-3010, or contact your local Sigma office.

AUSTRALIA AUSTRIA BELGIUM BRAZIL CZECH REPUBLIC FRANCE GERMANY HUNGARY INDIA ITALY JAPAN  
KOREA MEXICO NETHERLANDS POLAND SPAIN SWEDEN SWITZERLAND UNITED KINGDOM UNITED STATES

Millipore is a registered trademark of Millipore Corporation.

A Sigma-Aldrich Company

Circle No. 1 on Readers' Service Card




# LIPOFECTAMINE™ Reagent



## Highest Transfection Efficiencies

GIBCO BRL LIPOFECTAMINE Reagent provides previously unattainable transfection efficiencies—up to 90%  $\beta$ -Galactosidase-expressing cells. And as figure 1 shows, LIPOFECTAMINE Reagent is far superior to other transfection reagents.



**FIGURE 1.** Relative transfection efficiencies of LIPOFECTAMINE Reagent and two cationic lipid reagents from other suppliers. The amount of reagent and DNA were optimized according to manufacturers' recommended conditions. Peak activities are shown normalized relative to LIPOFECTAMINE Reagent.

Cell Lines and Cell Types			
NIH/3T3	293	COS-7	Fibroblasts
P19	HeLa	CHO-K1	Keratinocytes
PC 12	D5	BHK-21	Vero

LIPOFECTAMINE Reagent provides unparalleled transfection efficiencies for a wide variety of cell lines, including many that have traditionally been hard to transfect.

- Up to 90% positive cells.
- Superior results in both stable or transient transfections.
- Effective in the absence or presence of serum.
- Preformulated to be easy to use.
- Long shelf life at 4°C.
- Now available in 4-ml pack size.

Product	Cat. No.	Size
LIPOFECTAMINE Reagent	18324-012	1 ml
	18324-020	4 ml <b>New Size</b>

To Order/TECH-LINE™: **(800) 828-6686**  
 U.S. Industrial Orders/TECH-LINE: **(800) 874-4226**  
 Internet: <http://www.lifetech.com>



Producer of GIBCO BRL Products

**U.S.A. Orders**  
 To Order/TECH-LINE: (800) 828-6686  
 Fax: (800) 331-2286

**U.S. Industrial Orders**  
 To Order/TECH-LINE: (800) 874-4226  
 Fax: (800) 352-1468

**Latin America Orders**  
 To Order/TECH-LINE: (301) 840-4027  
 Fax: (301) 258-8238

**Canada Orders**  
 To Order: (800) 263-6236  
 TECH-LINE: (800) 757-8257  
 Fax: (800) 387-1007

**Internet**  
[info@lifetech.com](mailto:info@lifetech.com)  
<http://www.lifetech.com>



# Remodeling the extracellular matrix? There's only one place to find just the right enzyme.

COLLAGEN FIBRILS

LAMININ: CELL ATTACHMENT  
GROWTH PROMOTOR

COLLAGEN MOLECULE

HYALURONAN

TYPE VIII  
COLLAGEN

**Sigma has the specific enzyme you need for ECM research.**

Exploring the tissue architecture of the ECM is like any other project. The more appropriate and precise the tool for each application and activity, the more successful the result.

That's why you'll turn to Sigma. And select exactly what you need from a group of enzymes with an extensive range of specific activities.

Our comprehensive collection of enzymes for ECM research includes high purity collagenases,

leukocyte elastases, disases, hyaluronidases, tissue plasminogen activator and cathepsin G.

For guidance, contact our Technical Service experts. They'll rush literature to you by fax and mail to help you choose just the right enzymes to dismantle, rebuild and remodel the ECM.

**SIGMA®**

Call 1-800-325-5832,  
314-771-5750 or e-mail  
sigma-techserv@sial.com.

Sigma Chemical Co. P.O. Box 14508 St. Louis, MO 63178 USA 314-771-5750  
<http://www.sigma.sial.com>

Circle No. 23 on Readers' Service Card





**Get**

**Get Access to the Arrayed  
I.M.A.G.E. Consortium (LLNL)<sup>1</sup>  
Human and Mouse cDNA  
Libraries from Genome Systems.**

**Human and Mouse  
cDNAs**

**Now!**

For over four years,  
Genome Systems, Inc.  
has been supplying

researchers with high quality genomic clones and

unequaled technical support. Our commitment to you has resulted in our being  
chosen as a distribution site for the I.M.A.G.E. cDNA clones.

**How to Get a cDNA Clone** that matches  
your sequence using the internet world  
wide web.

**Access to Clones:** Over 300,000 5' and 3'  
sequences have been deposited in dbEST.

**Homology Search:** Perform a homology search  
against the dbEST database. Use the Blast  
program on the ncbi world wide web server. Type,  
<http://www.ncbi.nlm.nih.gov/dbEST/index.html>  
and search the dbEST database with your  
sequence or with a particular text string.

**Clone ID:** Provide Genome Systems  
with the dbEST database clone ID plus  
the genbank number and we will ship  
out the clone for next day delivery.

**Cost:** \$17.50/clone. Individual clones will  
be shipped to you 24 hours after adequate  
growth. \$14.00/384 well plate for not-for-profit  
organizations, \$19.00/384 well plate for  
commercial entities.

**Hybridization Screening:** Screen the library by  
hybridization with a DNA fragment. Genome  
Systems supplies high density gridded filters.  
Each set of filters contains the entire library with  
each clone double spotted on each filter that  
allows the investigator to easily identify the clone  
address. If you prefer, let Genome Systems  
perform the screen and ship you the clones.

**Break Down Your Genomic Clones:**


Genome Systems will help you isolate the  
cDNA clones whose sequences are  
contained within your genomic clone.

**Cost:** \$1,995/screen. Or you can purchase  
our double spotted high density gridded  
filters. **Cost:** \$2,350/filter set

**Characteristics of the Library**

**Tissue Source:** The I.M.A.G.E. cDNA clones  
are derived from ovary, retina, pineal gland, fetal  
spleen, liver, fetal liver spleen, breast, placenta,  
infant brain, adult brain and lung.

**Related Services:** P1 or PAC/BAC™ library  
screening services.



**GenomeSystemsInc™**

8620 Pennell Drive  
St. Louis, Missouri 63114, USA

(800) 430-0030 or, (314) 692-0033

Facsimile: (314) 692-0044

France: Appel gratuit, 0590-2104

Germany: Rufen sie uns an zum ortstarif, 0130-81-9081

UK: Call us free on, 0800-89-3733

email address: [Genome@MO.net](mailto:Genome@MO.net)

World Wide Web: <http://www.genomesystems.com>

1. "The I.M.A.G.E. Consortium: An Integrated Molecular Analysis of  
Genomes and their Expression", Lennon, G.G., Auffray, C.,  
Polymeropoulos, M., and Soares, M. B. [1995] Genomics.

Circle No. 16 on Readers' Service Card





<http://www.sciencemag.org>

# The world of SCIENCE... On-line every week!



- *SCIENCE Table of Contents*
- *The SCIENCE Editorial*
- *This Week in SCIENCE*
- *Beyond the Printed Page:  
On-line interactive projects & articles*
- *Abstracts of all research papers*
- *Synopsis of news stories*
- *SCIENCE Professional Network:  
On-line classified advertising*
- *The SCIENCE Electronic Marketplace:  
On-line product information*

**SCIENCE**  
COVERS THE WORLD



# Best Choice in DIGITAL IMAGING



- Document
- Archive
- Analyze

AlphaEase™ software makes image acquisition, image enhancement and image analysis a breeze.

**Easy**

**Fast**

Instant prints and quantitative analysis within **1-3 minutes!**



## Applications

- Fluorescence Imaging
- Visible Light Imaging
- Chemiluminescence – *low light imaging*
- Microscopy Imaging

**NEW**


For additional information and demonstration in your lab call

**1-800-795-5556**



**Alpha  
Innotech  
Corporation**

## Versatile



Circle No. 37 on Readers' Service Card