



AMERICAN
ASSOCIATION FOR THE
ADVANCEMENT OF
SCIENCE

SCIENCE


22 NOVEMBER 1996

\$7.00

VOL. 274 • PAGES 1269-1432



At Last!



A Robotic Liquid Handler That's Mastered the Art of Change.

New MultiPROBE® VersaTip™ technology automatically adapts pipetting tips to your assay procedures. Disposable tips, fixed washable tips, or both? Microliter or milliliter volumes? Liquid level sensing of ionic or nonionic solutions? Different disposable tip sizes and types? With MultiPROBE, you get it all in one system.

The new VersaTip automatically changes from fixed, washable tips to a variety of disposable tip sizes or types — within the same protocol and without any user-intervention. Now you can minimize consumable costs by using washable tips whenever possible, and at the same time optimize performance by automatically switching from small volume to larger volume disposable tips. In other words, MultiPROBE with VersaTip adapts to your application, instead of your having to adapt the application to it.

Drug Discovery Work

For drug discovery work, VersaTip will sense small sample volumes in microplates, and handle both ionic and nonionic solutions such as DMSO. With four- and eight-tip MultiPROBE systems, you can stop wasting precious compounds for all solubilization, distribution, and screening applications.


Molecular Biology

For molecular biology applications, VersaTip can automatically pick up micro tips to handle very small volumes, switch to larger

tips when required, or use washable fixed tips whenever possible. Ultra-sensitive liquid level sensing tips further eliminate any risk for DNA cross-contamination.

Clinical Testing

For clinical testing, VersaTip will automatically adapt to the carry-over requirements of your assays. Four-tip systems with Varispan



can switch from disposable to washable tips, and transfer samples from any size test tube or vial, to microplates or other labware — all without user intervention. Combine this versatility with a large deck capacity and a bar code reader that can read

256 tubes in 90 seconds, and you have the most efficient robotic liquid handler your lab can own.

Now, there's no need to compromise. Demand the MultiPROBE with VersaTip — the first robotic liquid handler that adapts to your assays.



PACKARD

A Canberra Company

Packard Instrument Company, 800 Research Parkway, Meriden, CT 06450 U.S.A.

Tel: 203-238-2351 Toll Free: 1-800-323-1891 FAX: 203-639-2172

Web Site: <http://www.packardinst.com> Email: webmaster@packardinst.com



Packard International Offices:

Australia, Mt Waverley 61-3-9543-4266; **Austria**, Vienna 43-1-2702504; **Belgium**, Brussels 32-2-4668210; **Canada**, Ontario 1-800-387-9559; **Central Europe**, Schwadorf, Aus. 43 456 2230 015; **Denmark**, Greve 45-42909023; **France**, Rungis (33) 1 46.86.27.75; **Germany**, Dreieich (49) 6103 385-151; **Italy**, Milano 39-2-33910796/7/8; **Japan**, Tokyo 81-3-3866-5850; **Netherlands**, Groningen 31-50-5413360; Tilburg (013) 5423900; **Russia**, Moscow, 7-095-238-7335; **Switzerland**, Zurich (01) 481 69 44; **United Kingdom**, Pangbourne, Berks (44) 01734 844981.

Circle No. 36 on Readers' Service Card

ANG	IL-10 sR ♦
AR	IL-11
BTC	IL-12
C10	IL-12 p40 ♦
CNTF	IL-13
β-ECGF	IL-15 ♦
EGF	IP-10
ENA-78	JE
Epo	KC ♦
FGF acidic	LAP
FGF basic	LIF
FGF-4	LIFR ♦
FGF-5	M-CSF
FGF-6	MCP-1
FGF-7	MCP-2 ♦
FGF-9 ♦	MCP-3 ♦
G-CSF	MIF ♦
GDNF ♦	MIP-1α
GM-CSF	MIP-1β
gp130 ♦	MIP-2 ♦
GROα	MK
GROβ	β-NGF
GROγ	NT-3 ♦
HB-EGF	NT-4 ♦
HRG-α ♦	OSM
I-309	PD-ECGF
IFN-γ	PDGF-AA
IGF-I	PDGF-AB
IGF-II	PDGF-BB
IL-1α	PI/GF
IL-1β	PTN
IL-1ra	RANTES
IL-1sRI ♦	SCF
IL-1sRII ♦	SLPI
IL-2	TGF-α
IL-2 sRα	TGF-β sRII ♦
IL-2 sRβ	TGF-β sRIII ♦
IL-2 sRγ ♦	TGF-β1
IL-3	TGF-β1.2
IL-3 Rα ♦	TGF-β2
IL-4	TGF-β3
IL-4 sRα ♦	TGF-β5
IL-5	TNF-α
IL-5 sRα ♦	TNF-β
IL-6	sTNF RI
IL-6 sR	sTNF RII
IL-7	Tpo ♦
IL-8	VEGF ₁₂₁ ♦
IL-9	VEGF ₁₆₅ ♦
IL-9R ♦	VEGF/PI/GF ♦
IL-10	♦ denotes NEW

R&D
SYSTEMS
1-800-343-7475

Searching for Cytokines and Antibodies?

R	I	N	T	E	R	F	E	R	O	N	V	G	M	C	S	F	X	S	C	T	B
&	N	C	H	E	M	O	K	I	N	E	S	X	S	M	H	Q	L	T	1	U	L
D	T	G	R	O	A	L	P	H	A	I	3	0	9	I	U	P	I	E	0	M	E
S	E	E	O	4	H	E	R	E	G	U	L	I	N	D	I	J	F	M	P	O	D
Y	R	S	M	F	X	I	N	T	E	R	L	E	U	K	I	N	8	C	L	R	J
S	L	A	B	G	P	G	R	O	G	A	M	M	A	I	W	D	N	E	E	N	V
T	E	R	O	F	A	R	6	F	G	F	9	X	N	N	P	O	O	L	I	E	X
E	U	E	P	5	H	O	C	Y	T	O	K	I	N	E	S	M	F	L	O	C	P
M	K	F	O	L	P	B	N	H	U	I	L	2	Y	3	N	G	E	F	T	R	2
S	I	I	I	I	L	E	T	A	F	U	Q	7	L	I	F	H	J	A	R	O	S
5	N	C	E	V	A	T	F	F	G	E	V	I	T	A	T	I	G	C	O	S	E
Q	1	U	T	K	F	A	J	E	H	Q	P	A	C	O	C	L	C	T	P	I	I
H	X	L	I	J	G	N	R	K	D	1	T	I	D	Y	I	1	S	O	H	S	D
R	G	D	N	F	T	I	L	1	0	S	D	N	M	W	S	5	F	R	I	F	O
N	E	R	Y	T	H	R	O	P	O	I	E	T	I	N	A	P	3	T	N	A	B
A	4	Y	6	P	Z	W	H	C	C	A	Z	E	G	F	B	D	N	F	K	C	I
W	L	L	M	B	R	A	N	T	E	S	M	I	F	V	F	G	K	Q	Z	T	T
I	I	A	N	G	I	O	G	E	N	I	N	H	B	E	G	F	N	T	4	O	N
L	P	L	A	C	E	N	T	A	G	R	O	W	T	H	F	A	C	T	O	R	A
5	B	E	T	A	C	E	L	L	U	L	I	N	D	T	G	F	B	E	T	A	Z

There are 60 cytokines hidden in this puzzle

The first place to look ...

R&D Systems offers the most extensive line of cytokines and soluble receptors as well as antibodies to cytokines and cytokine receptors.

Cytokines & Soluble Receptors

- Over 130 available – *human and mouse*
- Full biological activity
- >97% pure

Antibodies

- Over 150 available – *anti-human and anti-mouse*
- Polyclonal and monoclonal
- Highly active and specific

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES.

U.S.A. and Canada
R&D Systems, Inc.
614 McKinley Place NE
Minneapolis, MN 55413, USA.
Tel: 612 379-2956
Fax: 612 379-6580

Europe
R&D Systems Europe Ltd.
4-10 The Quadrant, Barton Lane
Abingdon, OX14 3YS, UK.
Tel: +44 (0)1235 531074
Fax: +44 (0)1235 533420


Germany
R&D Systems GmbH
Borsigstrasse 7
65205 Wiesbaden, Germany
Tel: +49 06122 90980
Fax: +49 06122 909819

Japan
Funakoshi Co., Ltd.
9-7, Hongo 2-Chome
Bunkyo-ku, Tokyo 113, Japan
Tel: +81 (03) 56841622
Fax: +81 (03) 56841633

International Distributors – Australia: (62) 008-25-1437. **Austria:** (43) 02 292 35 27. **Chile:** (56) 2-671-9369.

Israel: (972) 02 9230048. **Italy:** (39) 02 25 75 377. **Korea:** (850) 82-2-569-0781. **Mexico:** (52) 5-612-0085.


Spain & Portugal: (34) 01 448 84 86 or 03 456 97 06. **Taiwan:** (866) 2-368-3600.

 www.rndsystems.com


International Freefone Numbers – Belgique/België: 078 11 04 68. **Danmark:** 80 01 85 92. **Deutschland:** 013011 0169.

France: 05 90 72 49. **Nederland:** 060-225607. **Norge:** 800 1103. **Sverige:** 02079 31 49. **Switzerland:** 155 2482.

Circle No. 43 on Readers' Service Card



**AVANTI® MEANS
HIGH
ACCELERATION
CENTRIFUGATION**



**AVANTI® MEANS
HIGH
PRODUCTIVITY
CENTRIFUGATION**



BETTER QUALITY, INCREASED OPERATOR SAFETY

Accelerate up to $110,000 \times g$ in half the time of conventional high speed drives. Designed and tested as a system, Avanti J systems meet the world's highest standards for quality and operator safety. Even at maximum g -force, high-speed centrifugation is a breeze, certified by CE marking, UL and CSA approval.

We have two High Performance Centrifuge systems to choose from — the Avanti J-25 Series and our fastest model, Avanti J-30I. Both are engineering breakthroughs, and both offer longer warranties, diagnostic notification of instrument performance, and are designed for easy servicing.

Whichever Avanti J system you choose, you can be assured of a safe, quality drive with unmatched accel/decel rates, in a high-performance system that's driven to productivity.

BECKMAN

Worldwide Offices: Africa, Middle East, Eastern Europe (Switzerland) (41) 22 994 07 07. Australia (61) 2 844 6000. Austria (43) 222 7292164. Canada (800) 387-6799. China (861) 5051241-2. France (33) 1 43 01 70 00. Germany (49) 89-358700. Hong Kong, PRC (852) 814 7431. Italy (39) 2-953921. Japan 3-5352-2825. Mexico (525) 559-1635. Netherlands (0297-230630. Poland (48) 22 408822, 408833. Singapore (65) 339 3633. South Africa (27) 11-805-2014/5. Spain (1) 358-0051. Sweden (8) 445 11 20. Switzerland (41) 22 994 07 07. Taiwan (886) 02 378-3456. U.K. (01494) 441181. U.S.A. 1-800-742-2345.

© 1996 Beckman Instruments, Inc.

Circle No. 56 on Readers' Service Card



EASIER SET-UPS, FASTER SEPARATIONS

Break away from the pack. At up to $110,000 \times g$, Avanti J systems offer a range of separations broader than that of conventional high speeds. Plus, Avanti J systems are as easy to use as they are powerful. Every system is optimized by analog or digital display of set parameters and one-touch functions that accelerate your throughput while cutting run time.

We have two High Performance Centrifuge systems to choose from — the Avanti J-25 Series and our fastest model, Avanti J-30I. Both are engineering breakthroughs, and both offer longer warranties, diagnostic notification of instrument performance, and are designed for easy servicing.

Whichever Avanti J system you choose, you can be assured of a fast, easy way to produce more separations in less time.

BECKMAN

Worldwide Offices: Africa, Middle East, Eastern Europe (Switzerland) (41) 22 994 07 07. Australia (61) 2 844 6000. Austria (43) 222 7292164. Canada (800) 387-6799. China (861) 5051241-2. France (33) 1 43 01 70 00. Germany (49) 89-358700. Hong Kong, PRC (852) 814 7431. Italy (39) 2-953921. Japan 3-5352-2825. Mexico (525) 559-1635. Netherlands (0297-230630. Poland (48) 22 408822, 408833. Singapore (65) 339 3633. South Africa (27) 11-805-2014/5. Spain (1) 358-0051. Sweden (8) 445 11 20. Switzerland (41) 22 994 07 07. Taiwan (886) 02 378-3456. U.K. (01494) 441181. U.S.A. 1-800-742-2345.


© 1996 Beckman Instruments, Inc.

Circle No. 57 on Readers' Service Card



1292

A research
agency run by
researchers



1338

Cool galaxy clusters

NEWS & COMMENT

- Scientists With Clout **1292**
Generous Funding Wins a Seat at the
Genome Top Table **1293**
- Cluster Mission to Rise From the Ashes **1295**
- School Achievement: Asia and Europe Top
in World, But Reasons Are Hard to Find **1296**
- Mars Loss Could Sink Planetary Probes **1297**
- EU Stops Fiddling While Cows Burn **1297**
- A Rare Glimpse of an Early Human Face **1298**

RESEARCH NEWS

- The Case of the Missing Migrants **1299**
A Brewing Controversy Over Bird-Friendly
Coffee **1300**
- Homing In On a Prostate Cancer Gene **1301**
- Illusion Reveals Pain Locus in Brain **1301**
- Tumor Cells Fight Back to Beat
Immune System **1302**
- New Way to Read the Record
Suggests Abrupt Extinction **1303**
- Hedgehog's Patterning Call Is Patched
Through, Smoothly **1304**
- Quasar Pairs: A Redshift Puzzle? **1305**

- PER Protein in Silkworms Marches to
Different Drummer **1306**

PERSPECTIVES

- Apomixis: The Asexual Revolution **1322**
J.-P. V. Calzada, C. F. Crane, D. M. Stelly
- Hats Off to the Tricorn Protease **1323**
C. Schneider and F. U. Hartl
- Heavy Ozone—A Difficult Puzzle to
Solve **1324**
D. Krankowsky and K. Mauersberger
- The 110% Solution **1325**
D. Voss

ARTICLE

- Microwave Spectroscopy at the
Dissociation Limit **1327**
A. Carrington

REPORTS

- The Coulomb Blockade in Coupled
Quantum Dots **1332**
C. Livermore, C. H. Crouch, R. M. Westervelt, K. L.
Campman, A. C. Gossard
- Diffuse Extreme-Ultraviolet Emission **1335**
from the Coma Cluster: Evidence for Rapidly
Cooling Gases at Submegakelvin Temperatures
R. Lieu, J. P. D. Mittaz, S. Bowyer, J. O. Breen, F. J.
Lockman, E. M. Murphy, C.-y. Hwang

DEPARTMENTS

- THIS WEEK IN SCIENCE** **1277**
- EDITORIAL** **1283**
Materials Research and Applications
- LETTERS** **1285**
Activation of Estrogen Receptors: M. Joffe; J. A.
McLachlan • "Quantum Voodoo": E. Geislinger •
Flights from Reason?: W. P. Handwerker • Arsenic
Poisoning in West Bengal: D. P. Saha; K. S.
Subramanian • In Defense of Nannobacteria: R. L.
Folk • Cold Neutron Production: J. D. Jorgensen, J.
M. Carpenter, G. Aepli • "Quetzal" Coatings:
- S.H. Li, W. Messier
- SCIENCESCOPE** **1291**
- RANDOM SAMPLES** **1307**
- BOOK REVIEWS** **1320**
Krakatau, reviewed by T. J. Case • *Supernovae and
Nucleosynthesis*, G. J. Mathews • *Browsings*
- GORDON RESEARCH CONFERENCE** **1396**
- PRODUCTS & MATERIALS** **1403**

AAAS Board of Directors

Rita R. Colwell
Retiring President, *Chairman*
Jane Lubchenco
President
Mildred S. Dresselhaus
President-elect

Sheila Jasanoff
William A. Lester Jr.
Simon A. Levin
Marcia C. Linn
Michael J. Novacek
Anna C. Roosevelt
Jean E. Taylor
Nancy S. Wexler

William T. Golden
Treasurer
Richard S. Nicholson
Executive Officer

■ SCIENCE (ISSN 0036-8075) is published weekly on Friday, except the last week in December, by the American Association for the Advancement of Science, 1200 New York Avenue, NW, Washington, DC 20005. Periodicals Mail postage (publication No. 484460) paid at Washington, DC, and additional mailing offices. Copyright © 1996 by the American Association for the Advancement of Science. The title SCIENCE is a registered trademark of the AAAS. Domestic individual membership and subscription (51 issues): \$102 (\$55 allocated to subscription). Domestic institutional subscription (51 issues): \$250. Foreign postage extra: Mexico, Caribbean (surface mail) \$55; other countries (air assist delivery) \$90. First class, airmail, student, and emeritus rates on request. Canadian rates with GST available upon request, GST #1254 88122. Printed in the U.S.A.

COVER

Cretaceous-Tertiary (K-T) boundary near the town of Zumaya, Spain. The boundary occurs immediately below the ledge (lower center) formed by the erosion-resistant Tertiary strata. The boulders in the foreground lie in the pre-K-T gap. Statistical analysis of the fossil

distributions in this and other Bay of Biscay sections reveals a complex of extinction patterns in the latest Cretaceous. See page 1360 and the News story on page 1303. [Photo: Peter D. Ward, University of Washington]



Extreme-Ultraviolet Flux from the Virgo Cluster: Further Evidence for a 500,000-Kelvin Component 1338

S. Bowyer, M. Lampton, R. Lieu

Stratospheric Mean Ages and Transport Rates from Observations of Carbon Dioxide and Nitrous Oxide 1340

K. A. Boering, S. C. Wofsy, B. C. Daube, H. R. Schneider, M. Loewenstein, J. R. Podolske, T. J. Conway

An Explanation for Symmetry-Induced Isotopic Fractionation in Ozone 1344

G. I. Gellene

Oceanic Carbon Dioxide Uptake in a Model of Century-Scale Global Warming 1346

J. L. Sarmiento and C. Le Quééré

Red-Emitting Semiconductor Quantum Dot Lasers 1350

S. Fafard, K. Hinzer, S. Raymond, M. Dion, J. McCaffrey, Y. Feng, S. Charbonneau

Chemiluminescence in the Agglomeration of Metal Clusters 1353

L. König, I. Rabin, W. Schulze, G. Ertl

Glacial to Interglacial Fluctuations in Productivity in the Equatorial Pacific as Indicated by Marine Barite 1355

A. Paytan, M. Kastner, F. P. Chavez

Stability of Perovskite (MgSiO₃) in the Earth's Mantle 1357

S. K. Saxena, L. S. Dubrovinsky, P. Lazor, Y. Cerenius, P. Häggkvist, M. Hanfland, J. Hu

Sudden and Gradual Molluscan Extinctions in the Latest Cretaceous of Western European Tethys 1360

C. R. Marshall and P. D. Ward

Melanoma Cell Expression of Fas(Apo-1/CD95) Ligand: Implications for Tumor Immune Escape 1363

M. Hahne, D. Rimoldi, M. Schröter, P. Romero, M. Schreier, L. E. French, P. Schneider, T. Bornand, A. Fontana, D. Lienard, J.-C. Cerottini, J. Tschopp

Structure of the A Site of *Escherichia coli* 16S Ribosomal RNA Complexed with an Aminoglycoside Antibiotic 1367

D. Fourmy, M. I. Recht, S. C. Blanchard, J. D. Puglisi

Major Susceptibility Locus for Prostate Cancer on Chromosome 1 Suggested by a Genome-Wide Search 1371

J. R. Smith, D. Freije, J. D. Carpten, H. Grönberg, J. Xu, S. D. Isaacs, M. J. Brownstein, G. S. Bova, H. Guo, P. Bujnovsky, D. R. Nusskern, J.-E. Damber, A. Bergh, M. Emanuelsson, O. P. Kallioniemi, J. Walker-Daniels, J. E. Bailey-Wilson, T. H. Beaty, D. A. Meyers, P. C. Walsh, F. S. Collins, J. M. Trent, W. B. Isaacs

RAC Regulation of Actin Polymerization and Proliferation by a Pathway Distinct from Jun Kinase 1374

T. Joneson, M. McDonough, D. Bar-Sagi, L. Van Aelst

Uncoupling of Obesity from Insulin Resistance Through a Targeted Mutation in aP2, the Adipocyte Fatty Acid Binding Protein 1377

G. S. Hotamisligil, R. S. Johnson, R. J. Distel, R. Ellis, V. E. Papaioannou, B. M. Spiegelman

Liver Failure and Defective Hepatocyte Regeneration in Interleukin-6-Deficient Mice 1379

D. E. Cressman, L. E. Greenbaum, R. A. DeAngelis, G. Ciliberto, E. E. Furth, V. Poli, R. Taub

Neuroprotection by Aspirin and Sodium Salicylate Through Blockade of NF- κ B Activation 1383

M. Grilli, M. Pizzi, M. Memo, P. F. Spano

Tricorn Protease—The Core of a Modular Proteolytic System 1385

T. Tamura, N. Tamura, Z. Cejka, R. Hegerl, F. Lottspeich, W. Baumeister

Control of *C. elegans* Larval Development by Neuronal Expression of a TGF- β Homolog 1389

P. Ren, C.-S. Lim, R. Johnsen, P. S. Albert, D. Pilgrim, D. L. Riddle

TECHNICAL COMMENTS

Resistance to *Leishmania major* in Mice 1392


P. Demant, M. Lipoldova, M. Svobodova; *Response:* M. Güler, K. Murphy, J. Gorham

Role of β -Chemokines in Suppressing HIV Replication 1393

C. E. Mackewicz, E. Barker, J. A. Levy; *Response:* F. Cocchi, A. L. DeVico, A. Garzino-Demo, P. Lusso, R. C. Gallo

1323 & 1385

Tricorn protease



Indicates accompanying feature

Change of address: allow 4 weeks, giving old and new addresses and 8-digit account number. **Postmaster:** Send change of address to *Science*, P.O. Box 1811, Danbury, CT 06813-1811. **Single copy sales:** \$7.00 per issue prepaid includes surface postage; bulk rates on request. **Authorization to photocopy** material for internal or personal use under circumstances not falling within the fair use provisions of the Copyright Act is granted by AAAS to libraries and other users registered with the Copyright Clearance Center (CCC) Transactional Reporting Service, provided that \$4.00 per article is paid directly to CCC, 222 Rosewood Drive, Danvers, MA 01923. The identification code for *Science* is 0036-8075/83 \$4.00. *Science* is indexed in the *Reader's Guide to Periodical Literature* and in several specialized indexes.

On the Web


SCIENCE NOW

Get your daily news fix from *Science's* news team: <http://www.sciencenow.org/>

Introducing
AmpliTaq Gold™

PERFORMANCE
GUARANTEE

Where There's Gold,
You'll Find
Performance.



Co-amplification of 8 human STR loci. L1, 2: male control DNA; L3, 4: female control DNA; L5: AmpliTaq Gold negative control; L6, 7: male control DNA; L8, 9: female control DNA; L10: negative control.



Amplification of HIV-1 Control DNA. L2: 0 copies, AmpliTaq DNA Polymerase, No Hot Start; L3: 10 copies, AmpliTaq DNA Polymerase, No Hot Start; L4: 10 copies, AmpliTaq DNA Polymerase, manual Hot Start; L5: 10 copies, AmpliTaq Gold.

For PCR performance with higher yield, better specificity and more reliable results, discover AmpliTaq Gold™.

This new version of AmpliTaq® DNA Polymerase provides the specificity of Hot Start PCR, without all the extra steps. In most cases, you can substitute AmpliTaq Gold directly in existing PCR amplification protocols—without re-optimization.

You'll find AmpliTaq Gold saves time and money with dramatically lower drop-out rates, improved specificity, and easier multiplexing.

It also gives you consistently better PCR results. Because AmpliTaq Gold remains inactive until heated, conditions that lead to primer-dimer formation and mispriming are eliminated.

And of course, you have the continued assurance of knowing that AmpliTaq Gold

is backed by PE Applied Biosystems' exclusive PCR Performance Guarantee.

So discover AmpliTaq Gold, and discover high performance PCR. To request information, call 1-800-327-3002. Outside the U.S. and Canada, contact your local PE Applied Biosystems representative. On the Internet, visit our home page at <http://www.amplitaqgold.com>, or e-mail pebio@perkin-elmer.com.

PE Applied Biosystems

Europe Langen, Germany Tel: 49 6103 708 301 Fax: 49 6103 708 310
Japan Tokyo, Japan Tel: (0473) 80-8381 Fax: (0473) 80-8505
Latin America Mexico City, Mexico Tel: 52-5-651-7077 Fax: 52-5-593-6223
Australia Melbourne, Australia Tel: (03) 9212-8585 Fax: (03) 9212-8502

Perkin-Elmer PCR reagents are developed and manufactured by Roche Molecular Systems, Inc., Branchburg, New Jersey, U.S.A.




AmpliTaq Gold is a trademark and AmpliTaq is a registered trademark of Roche Molecular Systems, Inc. The PCR process is covered by patents owned by Hoffmann-La Roche, Inc. and F.Hoffmann-La Roche Ltd. PE Applied Biosystems is a trademark and Perkin-Elmer is a registered trademark of The Perkin-Elmer Corporation. The Perkin-Elmer Corporation is ISO 9001 certified.

Circle No. 41 on Readers' Service Card

Artificial molecules

A quantum dot is a semiconductor structure in which electrons are confined within a small volume and have discrete energy levels resembling those



of atoms. Livermore *et al.* (p. 1332) have constructed coupled quantum dots in which electrons can tunnel between dots, thus creating "artificial molecules." The conductances observed in both the weak and strong tunneling limits agree with predictions from many-body theory.

Glowing clusters

Chemical reactions may be accompanied by emission of visible light, or chemiluminescence, if the energy created by the reaction is stored initially in an excited state that later decays. König *et al.* (p. 1353) report that formation of metal clusters can be accompanied by chemiluminescence. During cluster agglomeration in a noble gas matrix, formation of unstable intermediates is proposed to lead to emission of excited fragments, which decay while emitting visible light.

Fresh air?

Transport of air from the troposphere into and out of the stratosphere, and its residence time in the stratosphere, can determine the rates at which ozone-destroying compounds reach the ozone layer and the

Insulin resistance, diabetes, and obesity

Obesity can lead to insulin resistance and diabetes, and the two appeared inseparable in animal models. Hotamisligil *et al.* (p. 1377) now report that in mice lacking the gene encoding AP2, the fatty acid-binding protein from adipocytes, dietary obesity fails to cause insulin resistance or diabetes. Somehow, AP2 must be critical for the metabolic pathway that leads from obesity to insulin resistance. The results provide a focus possibly intervening in the process that causes abnormal glucose homeostasis and symptoms of diabetes as a consequence of obesity.

effects of aircraft emissions. However, details of the global transport of air are difficult to obtain. Boering *et al.* (p. 1340) measured various gases on board the NASA ER-2 aircraft at tropospheric and stratospheric altitudes between 1992 and 1996 and showed that air enters the stratosphere continuously throughout the year and is distributed rapidly. The measurement allowed the determination of the mean age of stratospheric air, which is related to the residence time of pollutants in the stratosphere.

Not so stable

Earth's lower mantle has generally been thought to be composed mostly of perovskite [(Mg,Fe)SiO₃], and many geophysical models of the lower mantle are based on this assumption. Some earlier work had indicated, however, that perovskite containing some iron might not be stable throughout the pressure range of the lower mantle. Saxena *et al.* (p. 1357) performed synchrotron x-ray studies on high-pressure, high-temperature samples of the end-member MgSiO₃. The results suggest that the perovskite broke down to MgO (periclase) and SiO₂ (stishovite) at pressures of about 60 gigapascals.

Aspirin and glutamate

The neurotransmitter glutamate can actually be toxic to neurons if its levels are elevated for too long—as may occur during stroke. Grilli *et al.* (p. 1383) describe how aspirin, at doses already in frequent use for the treatment of arthritis, can help to protect rat neurons in primary culture and in hippocampal slices from glutamate-induced neurotoxicity.

Tumor evasion

Activated T cells are normally eliminated after an immune response is completed by expression of the Fas ligand (FasL, also called Apo-1 or CD95 ligand), and immune-privileged sites in the body such as the eye also express FasL. Hahne *et al.* (p. 1363; see the news story by Williams, p. 1302) show that unlike normal skin cells, malignant melanoma cells express FasL to avoid the immune response. Injection of mouse melanoma cells expressing FasL led to rapid tumorigenesis in normal mice but not in mice deficient in Fas.

Restoring the liver

Adult liver cells (hepatocytes) can replicate rapidly, thus allowing the liver to recover from

toxic-induced damage and to regenerate in a few days after surgical procedures that remove most of its mass. Cressman *et al.* (p. 1379) have shown that interleukin-6 (IL-6) is a critical cytokine in this recovery. Mice lacking the gene for IL-6 were unable to regenerate liver tissue unless IL-6 was administered exogenously after tissue removal. The necessity of IL-6 is an important consideration in strategies that produce decreases in cytokine activity in order to control liver damage, such as those used in treating cirrhosis.


Hats off

There has been a recent explosion in knowledge about how the cell uses selective proteolysis to control a variety of functions. Tamura *et al.* (p. 1385; see the Perspective by Schneider, p. 1323) describe the discovery of a protease complex likely to form part of a multicatalytic complex. The protease subunits have a distinctive structure in electron micrographs and resemble a three-cornered (tricorn) hat.

Stress signal

When the nematode *Caenorhabditis elegans* finds its environment inhospitable, it develops into a dauer larva adapted for survival in adverse conditions. Ren *et al.* (p. 1389) show that response to the pheromone that induces the dauer phase includes alterations in transcription of a growth factor related to transforming growth factor- β . The transcriptional modulation occurs within certain chemosensory neurons of the larva.

Leave No Stone Unturned.



You know the answers you need are out there. Lost under masses of data. To uncover them, leverage the world's most complete collection of scientific information right from your desktop with Knight-Ridder Information.

Search everything from technical and patent literature to industry information....online with DIALOG®, on CD-ROM, on the Web, or use our document delivery service. Whether you're searching individually or organization-wide, you'll have the freedom to transform information into insight, the impulse of science. You don't need to turn the world upside-down.

Find out more with Knight-Ridder Information –
more easily, more quickly, more productively.


Call toll-free in the US: 800-334-2564.

In Canada and California: 415-254-8800. Europe e-mail: enquiries@dm.krinfo.ch

Or visit us at <http://www.krinfo.com>.



**KNIGHT-RIDDER
INFORMATION**



INFORMATION SOLUTIONS ONLINE, ON CD-ROM, ON THE WEB, ON PAPER

©1996 Knight-Ridder Information, Inc. 2440 W. El Camino Real, Mountain View CA 94040. All rights reserved. DIALOG is a service mark of Knight-Ridder Information, Inc. Registered U.S. Patent and Trademark Office. scitod

Circle No. 44 on Readers' Service Card

High-Efficiency Lambda Cloning


Lambda ZAP® Vectors for Easy Library Screening and Amplification

Stratagene's Lambda ZAP® vectors* combine the high efficiency of lambda cloning—for easier library screening and amplification—with the convenience of a plasmid system.^{1,2} In addition, the insert size bias inherent in plasmid libraries is not found with libraries constructed in lambda phage. Stratagene offers several derivatives of the Lambda ZAP vector, each tailored to meet your specific cloning needs. The Lambda ZAP II vector has six unique cloning sites that accommodate inserts from 0-10 kb in length, and recombinants can be screened with either DNA or antibody probes. *In vivo* excision of the pBluescript® plasmid allows for rapid characterization of inserts in a plasmid system without time-consuming phage preparations or sub-cloning steps. In addition, entire libraries can be excised for screening and analysis.

The Uni-ZAP XR® vector, Lambda ZAP II vector digested for unidirectional cloning, ensures that all clones are in the proper orientation for protein expression. The ZAP Express™ vector allows unidirectional cloning, both eukaryotic and prokaryotic expression, and increased cloning capacity up to 12 kb. For studying signal transduction, cell growth and differentiation, gene expression, secretion and metabolism, Stratagene provides a complete line of HybriZAP® two-hybrid system products for generating cDNA or genomic libraries.

REFERENCES

1. Short, J.M., Fernandez, J.M., Sorge, J.A., and Huse, W.D. (1988) *Nucl. Acids Res.* 16: 7583-7600.



Lambda ZAP® II Vector Excision

Individual lambda plaques or an entire lambda library is allowed to infect cells that are co-infected with filamentous helper phage. Inside the cell, trans-acting proteins from the helper phage recognize two separate domains (initiator and terminator) positioned within the Lambda ZAP vector arms. Both of these signals are recognized by the helper phage gene III protein and a new DNA strand is synthesized, displacing the existing strand. The displaced strand is packaged as a filamentous phage by the helper phage proteins, and secreted from the cell. pBluescript plasmids are recovered by infecting an F' strain and growing in the presence of ampicillin.

Construct Directional cDNA Libraries with Stratagene's cDNA Synthesis Kits**

Stratagene's cDNA Synthesis Kit is the only cDNA synthesis kit quality controlled to produce a library with 2×10^6 primary clones. The method of choice for construction of directional cDNA libraries is the cDNA Synthesis Kit from Stratagene. This kit is designed to make directional cDNA libraries in your choice of innovative Lambda ZAP® cloning vectors. The kit uses 5-methyl-dCTP during first-strand synthesis, eliminating the need for site-specific methylases. All of Stratagene's cDNA synthesis kits are provided with *Pfu* DNA polymerase, instead of Klenow polymerase, to create blunt-ended cDNA before adaptor ligation. Studies show that using *Pfu* DNA polymerase to end polish cDNA creates more efficient adapter ligation, resulting in more primary clones. Fragments that have been cloned into any Lambda ZAP vector can be quickly excised to generate subclones or entire libraries in versatile phagemid vectors. This means no more time-consuming subcloning experiments.

For construction of high-quality directional cDNA libraries, Stratagene provides complete kits that include your choice of powerful Lambda ZAP vector, the cDNA Synthesis Kit and Gigapack® III Gold packaging extract. All components are also available separately.

*U.S. Patent Nos. 5,128,256 and 5,286,636, and European Patent No. 286200B1

**Patent Pending

Don't Settle for Less!

The Highest Efficiency Packaging Extract

Gigapack® III packaging extracts offer the highest packaging efficiencies available *plus* a convenient, single-tube format and simplified protocol. Gigapack III Gold extract is the only packaging extract that guarantees efficiencies of 2×10^9 pfu/μg DNA. Gigapack III Gold extract lacks all known *E. coli* restriction systems, allowing packaging of methylated DNA at efficiencies 15-fold higher than those obtained with restriction-positive extracts.

Since 1984, Gigapack extracts have offered the same consistent high performance. Stratagene's Gigapack III Gold extracts are the best extracts for high-efficiency construction of cDNA and genomic libraries. Don't trust your valuable libraries to anything less than Gigapack III Gold packaging extracts.

ZAP-cDNA® Gigapack® III Gold Cloning Kit	#200450
ZAP Express™ cDNA Gigapack® III Gold Cloning Kit	#200451
HybriZAP® cDNA Gigapack® III Gold Cloning Kit	#235612



UNITED STATES
Stratagene Headquarters
800-424-5400
INTERNET MAIL:
techservices@stratagene.com

AUSTRALIA: 1800-252-204
AUSTRIA: 1 5332666
BRAZIL: 11 530 7833
CANADA: 905 713 1201
DENMARK: 8 6101055

FRANCE: 1 34 60 2424
GERMANY: 0130-840911
HONG KONG: 5785839
ISRAEL: 3 7521687
ITALY: 2580 134 50

JAPAN: (Toyoko) 3 3660 4819
KOREA: 2 5560311
MALAYSIA: 37031888
THE NETHERLANDS: 33 4950094
NEW ZEALAND: 9 443 5867

NORWAY: 22715090
PORTUGAL: 1-4581641
ROC: 2 694 0066 or 2 720 2215
SINGAPORE: 2730898
SPAIN: 1 729 03 33


SWEDEN: 8 680 0845
SWITZERLAND: 061-6930540
THAILAND: 23080611
UNITED KINGDOM: 0800 585370
OTHER COUNTRIES CALL STRATAGENE U.S.A.



TAKE THE GUESSWORK OUT OF PROTEIN PURIFICATION

ProScale — a revolutionary way to achieve rapid, predictable results for protein purification. By integrating unique chemistry, intelligent software and hardware, ProScale gives you a logical and efficient path to optimize purification. From research to process development, all the way to production, ProScale offers you a no-fail formula for success with your protein.

Plus, ProScale is the only system that comes with proven, scalable HyperD[®] media. Matched with HyperD, our BioSys[™] protein purification workstations are designed to deliver more purified protein, in less time, while retaining maximum activity.



*Sample: Cytochrom c, BSA, Transferrin.
Separations on 10, 20 and 35 µm
Q HyperD particles showing
the same selectivity.*

Our intelligent ProScale software takes much of the guesswork out of method development and optimization, thereby saving



you time and effort. By running simulation experiments within the software, much of your trial-and-error work is eliminated.

Take the guesswork out of protein purification and put ProScale in your lab. For additional information, just call your local Beckman representative or visit our Web site at <http://www.beckman.com>.

* HyperD is a trademark of BioSeptra, Inc.

BECKMAN

Worldwide Offices: Africa, Middle East, Eastern Europe (Switzerland) (41) 22 994 07 07. Australia (61) 02 816 5288. Austria (43) 222 7292164. Canada (800) 387-6799. China (861) 5051241-2. France (33) 1 43 01 70 00. Germany (49) 89-358700. Hong Kong (852) 814 7431. Italy (39) 2-953921. Japan 3-5352-2825. Mexico (525) 5-59-16-35. Netherlands 02979-85651. Poland (48) 22 408822, 408833. Singapore (65) 339 3633. South Africa (27) 11-805-2014/5. Spain (1) 358-0051. Sweden (8) 98-5320. Switzerland (41) 22 994 07 07. Taiwan (886) 02 378-3456. U.K. (01494) 441181. U.S.A. 1-800-742-2345.

©1996 Beckman Instruments, Inc.

Circle No. 22 on Readers' Service Card

MEMBERSHIP/CIRCULATION

Director: Michael Spinella
Deputy Director: Marlene Zendell
Member Services: Donald Crowl, *Manager*; Mary Curry, *Supervisor*; Pat Butler, Laurie Baker, Yvette Carter *Representatives*
Marketing: Dee Valencia, *Manager*; Jane Pennington, *Europe Manager*; Hilary Baar, *Assistant Manager*; Lauri Sirois, *Coordinator*
Research: Renuka Chander, *Manager*
Business and Finance: Robert Smariga, *Manager*; Nina Araujo de Kobes, Kimberly Parker, *Coordinators*
Computer Specialist: Charles Munson

ADVERTISING AND FINANCE

Associate Publisher: Beth Rosner
Advertising Sales Manager: Susan A. Meredith
Recruitment Advertising Manager: Janis Crowley
Business Manager: Deborah Rivera-Wienhold
Finance: Randy Yi, *Senior Analyst*; Connie Dang, *Financial Analyst*
Marketing: John Meyers, *Manager*; Allison Pritchard, *Associate*
Electronic Media: David Gilliken, *Manager*; Mark Croatti, Crystal Young, *Production Associates*
Product Advertising: Carol Maddox, *Traffic Manager*; Natalie Britt, *Sales Associate*
Product Advertising Sales: **East Coast/E. Canada:** Richard Teeling, 201-904-9774, FAX 201-904-9701 • **Midwest/Southeast:** Elizabeth Mosko, 773-665-1150, FAX 773-665-2129 • **West Coast/W. Canada:** Neil Boylan, 415-673-9265, FAX 415-673-9267 • **UK, Scandinavia, France, Italy, Belgium, Netherlands:** Andrew Davies, (44) 1-457-838-519, FAX (44) 1-457-838-898 • **Germany/Switzerland/Austria:** Tracey Peers, (44) 1-260-297-530, FAX (44) 1-260-271-022 • **Japan:** Masay Yoshikawa, (81) 3-3235-5961, FAX (81) 3-3235-5852
Recruitment Advertising: Terri Seiter Azie, *Sales and Production Operations Manager*; Celeste Miller, *Sales Supervisor*; Eric Banks, Troy Benitez, Bren Peters-Minnis, Bethany Ritchey, *Sales*; Debbie Cummings, *European Sales Manager*; Wendy Green, *Production Associate*; Nicole Robinson, *Advertising Assistant*
Recruitment Advertising Sales: **US:** Janis Crowley, 202-326-6532, FAX 202-289-6742 • **Europe:** Debbie Cummings, (44) 1223-302067, FAX (44) 1223-576208 • **Australia/New Zealand:** Keith Sandell, (61) 02-922-2977, FAX (61) 02-922-1100 • **Japan:** Masay Yoshikawa, (81) 3-3235-5961, FAX (81) 3-3235-5852
Assistant to Associate Publisher: Nancy Hicks
Permissions: Lincoln Richman
Exhibits Coordinator: Arlene Ennis
Send materials to Science Advertising, 1200 New York Avenue, NW, Washington, DC 20005.

SCIENCE EDITORIAL BOARD

Charles J. Arntzen
David Baltimore
J. Michael Bishop
William F. Brinkman
E. Margaret Burbidge
Pierre-Gilles de Gennes
Joseph L. Goldstein
Mary L. Good
Harry B. Gray
John J. Hopfield
F. Clark Howell
Paul A. Marks
Yasutomi Nishizuka
Helen M. Ranney
Bengt Samuelsson
Robert M. Solow
Edward C. Stone
James D. Watson
Richard N. Zare

■ Published by the American Association for the Advancement of Science (AAAS), *Science* serves its readers as a forum for the presentation and discussion of important issues related to the advancement of science, including the presentation of minority or conflicting points of view, rather than by publishing only material on which a consensus has been reached. Accordingly, all articles published in *Science*—including editorials, news and comment, and book reviews—are signed and reflect the individual views of the authors and not official points of view adopted by the AAAS or the institutions with which the authors are affiliated.

■ The American Association for the Advancement of Science was founded in 1848 and incorporated in 1874. Its objectives are to further the work of scientists, to facilitate cooperation among them, to foster scientific freedom and responsibility, to improve the effectiveness of science in the promotion of human welfare, to advance education in science, and to increase public understanding and appreciation of the importance and promise of the methods of science in human progress.

INFORMATION RESOURCES

SUBSCRIPTION SERVICES

For change of address, missing issues, new orders and renewals, and payment questions, please contact AAAS at Danbury, CT: 800-731-4939 or Washington, DC: 202-326-6417, FAX 202-842-1065. Mailing addresses: AAAS, P.O. Box 1811, Danbury, CT 06813 or AAAS Member Services, 1200 New York Avenue, NW, Washington, DC 20005 • **Other AAAS Programs:** 202-326-6400

REPRINTS & PERMISSION

Reprints: Ordering/Billing, 800-407-9191; Corrections, 202-326-6501 • **Permissions:** 202-326-7074, FAX 202-682-0816

INTERNET ADDRESSES

science_editors@aaas.org (for general editorial queries); science_news@aaas.org (for news queries); science_letters@aaas.org (for letters to the editor); science_reviews@aaas.org (for returning manuscript reviews); science@science-int.co.uk (for the Europe Office); membership@aaas.org (for member services); science_classifieds@aaas.org (for submitting classified advertisements); science_advertising@aaas.org (for product advertising)

INFORMATION FOR CONTRIBUTORS

See pages 93-94 of the 5 January 1996 issue or

access <http://www.sciencemag.org/science/home/con-info.shtml>.

EDITORIAL & NEWS CONTACTS

North America

Address: 1200 New York Avenue, NW, Washington, DC 20005
Editorial: 202-326-6501, FAX 202-289-7562
News: 202-326-6500, FAX 202-371-9227 • **Bureaus:** **Berkeley, CA:** 510-841-1154, FAX 510-841-6339, **San Diego, CA:** 619-942-3252, FAX 619-942-4979, **Chicago, IL:** 312-360-1227, FAX 312-360-0537, **Boston, MA:** 617-566-7137, FAX 617-734-8088

Europe

Headquarters: 14 George IV Street, Cambridge, UK CB2 1HH; (44) 1223-302067, FAX (44) 1223-302068
Paris Correspondent: (33) 1-49-29-09-01, FAX (33) 1-49-29-09-00

Asia

Japan Office: Carl Kay, Esaka-cho 5-chome 11-10, Suita-shi, Osaka 564 Japan; (81) 6-368-1925, FAX (81) 6-368-6905; science@magical.egg.or.jp
News Bureau: (81) 3-3335-9925, FAX (81) 3-3335-4898 • **China Office:** Hao Xin, (86)10-6255-9478; science@public3.bta.net.cn

BOARD OF REVIEWING EDITORS

Frederick W. Alt <i>Children's Hospital, Boston</i>	G. Ertl <i>Max-Planck-Gesellschaft</i>	Richard Losick <i>Harvard Univ.</i>	Gottfried Schatz <i>Biozentrum, Basel</i>
Don L. Anderson <i>California Institute of Technology</i>	Richard G. Fairbanks <i>Lamont-Doherty Earth Observatory</i>	Reinhard Lührmann <i>Institut für Molekularbiologie und Tumorforschung der Philipps-Universität</i>	Jozef Schell <i>Max-Planck-Institut für Zuchtungsforschung</i>
Michael Ashburner <i>Univ. of Cambridge</i>	Robert E. Fay <i>U.S. Bureau of the Census</i>	Ruth Lynden-Bell <i>Queen's Univ., Belfast</i>	Ronald H. Schwartz <i>National Institute of Allergy and Infectious Diseases, NIH</i>
Stephen J. Benkovic <i>Pennsylvania State Univ.</i>	Douglas T. Fearon <i>Univ. of Cambridge</i>	Seth Marder <i>California Institute of Technology</i>	Terrence J. Sejnowski <i>Salk Institute</i>
Alan Bernstein <i>Mount Sinai Hospital, Toronto</i>	Harry A. Fozzard <i>The Univ. of Chicago</i>	Diane Mathis <i>Institut de Chimie Biologique, Strasbourg</i>	Christopher R. Somerville <i>Carnegie Institute of Washington</i>
Seth Blair <i>Univ. of Wisconsin, Madison</i>	Roger I. M. Glass <i>Centers for Disease Control</i>	Susan K. McConnell <i>Stanford Univ.</i>	Thomas A. Steitz <i>Yale Univ.</i>
David E. Bloom <i>Harvard Institute for International Development</i>	Stephen P. Goff <i>Columbia Univ. College of Physicians & Surgeons</i>	Anthony R. Means <i>Duke Univ. Medical Center</i>	Michael P. Stryker <i>Univ. of California, San Francisco</i>
Piet Borst <i>The Netherlands Cancer Institute</i>	Peter N. Goodfellow <i>SmithKline Beecham, UK</i>	Stanley Meisel <i>Univ. of California, Davis</i>	Cliff Tabin <i>Harvard Medical School</i>
Henry R. Bourne <i>Univ. of California, San Francisco</i>	Peter Gruss <i>Max Planck Institute of Biophysical Chemistry</i>	Shigetada Nakanishi <i>Kyoto Univ.</i>	John Jen Tai <i>Academia Sinica, Taiwan</i>
Michael S. Brown <i>Univ. of Texas Southwestern Medical Center</i>	Philip C. Hanawalt <i>Stanford Univ.</i>	Kim Nasmyth <i>Research Institute of Molecular Pathology, Vienna</i>	Tomoyuki Takahashi <i>Univ. of Tokyo</i>
James J. Bull <i>Univ. of Texas at Austin</i>	Paul Harvey <i>Univ. of Oxford</i>	Roger A. Nicoll <i>Univ. of California, San Francisco</i>	Masatoshi Takeichi <i>Kyoto Univ.</i>
Kathryn Calame <i>Columbia Univ. College of Physicians & Surgeons</i>	M. P. Hassell <i>Imperial College at Silwood Park</i>	Staffan Normark <i>Swedish Institute for Infectious Disease Control</i>	Keiji Tanaka <i>RIKEN Institute</i>
Dennis W. Choi <i>Washington Univ. School of Medicine, St. Louis</i>	Nobutaka Hirokawa <i>Univ. of Tokyo</i>	Bert W. O'Malley <i>Baylor College of Medicine</i>	David Tilman <i>Univ. of Minnesota, St. Paul</i>
David Clapham <i>Mayo Foundation</i>	Tomas Hökfelt <i>Karolinska Institutet</i>	Stuart L. Pimm <i>The Univ. of Tennessee, Knoxville</i>	Robert T. N. Tjian <i>Univ. of California, Berkeley</i>
Adrienne E. Clarke <i>Univ. of Melbourne, Parkville</i>	Tasuku Honjo <i>Kyoto Univ.</i>	Yeshayau Pocker <i>Univ. of Washington, Seattle</i>	Yoshinori Tokura <i>Univ. of Tokyo</i>
John M. Coffin <i>Tufts Univ. School of Medicine</i>	Susan D. Iversen <i>Univ. of Oxford</i>	Ralph S. Quatrano <i>Univ. of North Carolina, Chapel Hill</i>	Emil R. Unanue <i>Washington Univ. School of Medicine, St. Louis</i>
F. Fleming Crim <i>Univ. of Wisconsin, Madison</i>	Eric F. Johnson <i>The Scripps Research Institute</i>	Martin Raff <i>Univ. College London</i>	Derek van der Kooy <i>Univ. of Toronto</i>
Paul J. Cruzten <i>Max-Planck-Institut für Chemie</i>	Elliott Kieff <i>Harvard Univ.</i>	V. Ramanathan <i>Scripps Institute of Oceanography</i>	Geerat J. Vermeij <i>Univ. of California, Davis</i>
James E. Dahlberg <i>Univ. of Wisconsin Medical School, Madison</i>	Judith Kimble <i>Univ. of Wisconsin, Madison</i>	Douglas C. Rees <i>California Institute of Technology</i>	Bert Vogelstein <i>Johns Hopkins Oncology Center</i>
Robert Desimone <i>National Institute of Mental Health, NIH</i>	Stephen M. Kosslyn <i>Harvard Univ.</i>	T. M. Rice <i>ETH-Hönggerberg, Zürich</i>	Arthur Weiss <i>Univ. of California, San Francisco</i>
Paul T. Englund <i>Johns Hopkins Univ. School of Medicine</i>	Michael LaBarbera <i>The Univ. of Chicago</i>	David C. Rubie <i>Universität Bayreuth</i>	Zena Werb <i>Univ. of California, San Francisco</i>
	Nicole Le Douarin <i>Institut d'Embryologie Cellulaire et Moléculaire du CNRS</i>	Erkki Ruoslahti <i>The Burnham Institute, CA</i>	George M. Whitesides <i>Harvard Univ.</i>
	Charles S. Levings III <i>North Carolina State Univ., Raleigh</i>		Owen N. Witte <i>Univ. of California, Los Angeles</i>
	Harvey F. Lodish <i>Whitehead Institute for Biomedical Research</i>		

Picture Perfect

BioMax Film
Maximum Performance For Biological Imaging



KODAK Products and Techware — the perfect match for scientific imaging.



When it comes to your biological imaging and analysis product needs, there's no better choice than KODAK products from Sigma-Aldrich Techware.

A single phone call is all it takes to gain immediate access to many KODAK Scientific Imaging Systems products...as well as over 10,000 other brand name laboratory equipment and supply products.

In addition to a wide selection of top quality laboratory products, Sigma-Aldrich Techware backs each

purchase with unequalled service and support, including:

- Same day shipping
- Comprehensive technical support
- Competitive pricing



If your research calls for biological imaging and analysis products that are picture perfect, choose KODAK products from Sigma-Aldrich Techware—your single, most convenient source for brand name laboratory equipment and supplies.

SIGMA-ALDRICH
Techware
FROM BENCHTOP TO BOOKSHELF

SIGMA[®] 


Where Science and Service Come Together.

Call collect: 314-771-5750,
Toll Free: 800-325-3010,
or contact your local Sigma office.
World Wide Web: <http://www.sigma.sial.com>


AUSTRALIA AUSTRIA BELGIUM BRAZIL CZECH REPUBLIC FRANCE GERMANY HUNGARY INDIA ISRAEL ITALY
JAPAN KOREA MEXICO THE NETHERLANDS POLAND SPAIN SWEDEN SWITZERLAND UNITED KINGDOM UNITED STATES

A Sigma-Aldrich Company
Circle No. 2 on Readers' Service Card


Biotinylated dextran injected iontophoretically into the auditory cortex of rat brain to label somata and neuronal profiles. Tissue treated with KPL's Streptavidin-Peroxidase, followed with TrueBlue[®] Peroxidase Substrate. Counterstained with Contrast RED. Photographed using a 3-dimensional light microscope.




Microplate wells coated with 1.5 µg/ml of rabbit IgG excepting row A, then blocked. Alkaline phosphatase-labeled goat anti-rabbit IgG (H+L) added to row A at a concentration of 50 ng/ml, with a dilution series of conjugate ending in row H at a concentration of 2.5 ng/ml. Incubated for 20 minutes with KPL's BluePhos[®] Substrate (left side) and pNPP Phosphatase Substrate (right side).



Why use radioactivity when you can get results like this without it?



Southern Blot of EcoR1 digested pE1.4 detected using KPL's Biotin Detection Kit. Lane 1: 10 ng pE1.4. Lanes 2-6: 5 ng, 1 ng, 500 pg, 100 pg, 50 pg. Lanes 7-11: 10 pg, 5 pg, 1 pg, 0.5 pg, 0.1 pg.



Western Blots of HIV-1 viral lysate electrophoreted onto nitrocellulose membrane. Immunodetection with HIV-positive human serum followed by enzyme-labeled antibody. Immunodetected with appropriate KPL substrate.

Detect, measure and characterize with KPL's high-performance, non-radioactive kits and reagents.


You don't need another lecture on the usage and disposal problems associated with radioactivity. But you do need to know about viable, reliable alternatives. KPL offers over 600 products for applications in ELISA; Western, Southern and Northern blotting; and immunohistochemistry.

- Nucleic Acid Labeling and Detection Systems
- Substrates for Microwell ELISA, Immunoblotting and Immunohistochemistry
- Western Blot Systems
- HistoMark[®] Biotin Streptavidin Kits
- Secondary Antibodies to Human and Animal Immunoglobulins
- BacTrace[®] Antibodies to Selected Bacteria

Because we make our own products, we can provide customers with immediate access to in-depth expertise via our Technical Service line and FastFacts[™]. And we're happy to discuss custom and bulk orders.

Granted, sometimes you have no choice but to use radioactive methods. However, most times you do. That's when you should count on KPL.

For a catalog and more information, call 1-800-638-3167.



2 Cessna Court, Gaithersburg, MD 20879-4174 U.S.A. • Phone: 800-638-3167 or 301-948-7755 • Fax: 301-948-0169

TrueBlue[™], BluePhos[™], HistoMark[®], BacTrace[®] and FastFacts[™] are trademarks of KPL, Inc.

Circle No. 29 on Readers' Service Card

GO GERMLINE™



a la carte

MOUSE KNOCKOUT SERVICES

Genome Systems, Inc. introduces a complete line of mouse knockout services to fit each investigator's unique requirements. Whether you need our complete knockout services or just a "taste", Genome Systems provides the answer.

Isolation of 129/SvJ genomic clones from our BAC library: All knockout experiments start with isolation of your gene. Genome systems is the **only** company that makes the libraries it screens. \$995/screen. catalog # BAC-4921

Engineering of the targeting vector: We prefer to start with our own 129/SvJ DNA. However, if you have already isolated your gene, you can supply us with DNA from another 129 strain. We will subclone fragments, extensively map your gene and do all sequencing necessary to make a good targeting construct for homologous recombination. \$10,000/experiment. catalog # MK-2110

Electroporation and antibiotic selection of ES cells: We will take your targeting vector, electroporate our "Go Germline" RW-4 ES cells and plate them on our primary MEF feeder cells. We do all tissue culture and antibiotic selection. You receive a minimum of 150 ES clones, in duplicate, for Southern analysis and eventual blastocyst injection or morula aggregation. \$3,500/150 clone minimum. catalog # MK-2120

Southern analysis of antibiotic selected ES clones: We take ES clones and perform all Southern analysis so we can determine which clone contains the knocked out allele. \$7,650/150 clones. catalog # MK-2130; or you provide the targeting construct \$3,825/150 clones. catalog # MK2130A

Karyotyping: We will examine 20 metaphases and send back a numerical count and structural analysis of the chromosomes. \$750 per clone. catalog # MK-2150

Morula aggregation: We use the latest technology to perform Morula aggregation to create highly chimeric mice. U.S. patent (#5,449,620) is licensed exclusively to Genome Systems. Send us your targeted ES cells and we will send you at least 3 highly chimeric male mice so you can obtain founder mice. \$7,000/experiment catalog # MK-2140

GenomeSystemsInc™

8620 Pennell Drive
St. Louis, Missouri 63114, USA
(800) 430-0030 or, (314) 692-0033
Facsimile: (314) 692-0044

Circle No. 61 on Readers' Service Card

France: Appel gratuit, 0590-2104


Germany: Rufen sie uns an zum ortstarif, 0130-81-9081

UK: Call us free on, 0800-89-3733

email address: Genome@MO.net

World Wide Web: <http://www.genomesystems.com>






**OUR NEW
cDNA SUBTRACTION KIT
HAS ONE IMPORTANT
FEATURE...
IT WORKS.**

Subtracted

Unsubtracted

PCR-Select™ cDNA Subtraction Kit

Above photo: Southern blot analysis showing enrichment of IL-2 receptor, known to be activated in human Jurkat T-cells by PHA/PMA treatment. Lane 1: unsubtracted cDNA, treated cells. Lane 2: unsubtracted cDNA, untreated cells. Lane 3: subtracted cDNA, treated cells. A dramatic reduction of the abundant housekeeping gene, G3PDH, was also seen (data not shown).



Northern blot analysis of two selected representative cDNA clones from a testis-specific library. Testis cDNA was subtracted against a mixture of cDNAs from 10 different tissues. The subtracted testis-specific cDNA was cloned into a plasmid vector, and 10 randomly selected clones were used to probe Human Multiple Tissue Northern Blots containing 2 µg of poly A⁺ RNA from the indicated tissues. All 10 clones hybridized only to testis RNA. The exposure times needed to generate signal ranged from 5 hr (Panel A) to 7 days (Panel B), indicating that abundant and relatively rare cDNAs were obtained. Lanes 1-16: heart (1), brain (2), placenta (3), lung (4), liver (5), skeletal muscle (6), kidney (7), pancreas (8), spleen (9), thymus (10), prostate (11), testis (12), ovary (13), small intestine (14), colon (15), peripheral blood leukocyte (16).

Introducing the **PCR-Select™ cDNA Subtraction Kit** (#K1804-1), a revolutionary new method for finding differentially expressed genes. Actually, it has several important features. First, you selectively amplify only your genes of interest—in fact, **no physical subtraction** is needed. Further, this unique method allows you to **equalize and subtract in one procedure**, to dramatically increase the probability of obtaining cDNAs that correspond to rare transcripts—with **greater than 1,000-fold enrichment**. Thus, you can directly clone the subtracted cDNAs or use them as hybridization probes to screen cDNA or genomic libraries. And, it only requires 0.5–2.0 µg of poly A⁺ RNA. But, best of all, it works. For more information or to order, call 1-800-662-CLON or contact your local distributor.

Circle No. 49 on Readers' Service Card

DON'T EXPERIMENT WITH ANYTHING ELSE.

CLONTECH

1020 East Meadow Circle, Palo Alto, California 94303 USA

FAX: 800/424-1350 415/424-1064 • TEL: 800/662-CLON 415/424-8222


e-mail: tech@clontech.com • orders@clontech.com

©1996 CLONTECH Laboratories, Inc.

PCR is covered by patents owned by Hoffmann-La Roche and F. Hoffmann-La Roche, A.G.

(AD61036)

New, Fast, Sensitive and Safe Labeling Technology



Introducing BioLite™: A luminescent labeling and detection assay kit for the quantification of cells, particles, and microorganisms.

BioLite¹, a long-lived “glow” type signal assay kit, allows you to carry out high throughput cell adhesion, chemotactic or infection studies in the microplate format. With BioLite, you simply pre-label your cells, particles, or microorganisms of interest using the label provided. Perform your assay; wash any non-adherent, non-bound, or non-invasive labeled cells or microorganisms free; add the BioLite detection reagent; and measure the produced signal.

BioLite labeling is rapid, permanent, cell-type independent, and does not require long incubation times. The assay simplicity and high detection sensitivity with BioLite enables the use of microplate technology for easy handling and high throughput analysis. And, with a half-life of several hours, BioLite allows you to prepare multiple microplates at the same time for compound screening assays or to measure plates at multiple time points.

¹Patent pending

Combined with the TopCount™ microplate scintillation and luminescence counter, BioLite assays can be performed in either the 96- or 384-well microplate format for batch processing of cellular or microbiological assays. TopCount enables you to analyze 12 samples simultaneously. It also provides temperature control and stackers to hold up to 40 microplates.



Switch on BioLite for:

- Rapid, non-isotopic labeling and detection
- Long-lived luminescent signal (half-life of several hours)
- High sensitivity
- A label that binds equally well to all cell types
- A label that does not interfere with cell adhesion kinetics or membrane receptor function



PACKARD
A Canberra Company


Packard Instrument Company, 800 Research Parkway, Meriden, CT 06450 U.S.A.
Tel: 203-238-2351 Toll Free: 1-800-323-1891 FAX: 203-639-2172
Web Site: <http://www.packardinst.com> Email: webmaster@packardinst.com



Packard International Offices:

Australia, Mt Waverley 61-3-9543-4266; **Austria**, Vienna 43-1-2702504; **Belgium**, Brussels 32-2-4668210; **Canada**, Ontario 1-800-387-9559; **Central Europe**, Schwadorf, Aus. 43 456 2230 015; **Denmark**, Greve 45-42909023; **France**, Rungis (33) 1 46.86.27.75; **Germany**, Dreieich (49) 6103 385-151; **Italy**, Milano 39-2-33910796/7/8; **Japan**, Tokyo 81-3-3866-5850; **Netherlands**, Groningen 31-50-5413360; Tilburg (013) 5423900; **Russia**, Moscow, 7-095-259-9632; **Switzerland**, Zurich (01) 481 69 44; **United Kingdom**, Pangbourne, Berks (44) 01734 844981.

Circle No. 45 on Readers' Service Card



Looking into protein folding?

Begin with our purified, performance-tested GroEL and GroES chaperonins and antibodies.

Now you can get everything you need to study polypeptide folding and assembly from one convenient source: Sigma.

Our comprehensive selection of heat shock proteins includes chaperonin 60 (GroEL), chaperonin 10 (GroES), and a 1:1 mixture of the two. Purified to >95% as determined by SDS-PAGE, they're also assayed for folding activity and ATPase activity.

To help you detect and study these proteins, we offer antibodies to GroEL and GroES.

Mono-specific, they do not cross-react with other chaperonins. They're tested in both immunoblotting and immunoprecipitation. And Anti-GroEL is also available labeled with peroxidase for direct immunoblotting.

Because we develop these products together and test them as a system, you can be sure of their performance, each and every time.

Look into protein folding with reliable, reasonably-priced Sigma reagents. And explore with confidence.

SIGMA®

To learn more, call Technical Service, 1-800-262-9141 or your local Sigma office.

Circle No. 35 on Readers' Service Card


CyDye fluorescent reagents

Direct labelling for the brightest results


If you want the ease of direct non-radioactive labelling coupled with high sensitivity and bright, multicolour results, you need look no further than the extensive range of CyDye™ fluorescent probes from Amersham™.

Based on the cyanine fluor, the seven CyDye fluors all offer intense colours with a narrow emission spectrum, allowing easy discrimination between multiple probes on the same sample. Over 180 CyDye products are available - either as individual reactive NHS-esters or as part of optimized kits. All are designed to meet the needs of specific applications including chromosome painting, DNA labelling, protein labelling and immunocytochemistry.

**Make the right choice - the bright choice -
CyDye from Amersham.**



Three colour confocal laser scanning microscope analysis of islet from neonatal rat pancreas. FITC anti-insulin (green), Cy™3 anti-somatostatin (red), and Cy5 anti-glucagon (blue).



Whole insect embryos directly labelled with Cy3 dye (red) and Cy5 dye (blue).
Image courtesy of Dr. T.C. Brelje,
University of Minnesota Medical School.



Circle No. 31 on Readers' Service Card

Amersham International plc
Amersham Place, Little Chalfont,
Buckinghamshire, HP7 9NA England
UK 0800 515 313
Europe (01494) 544 000
USA (800) 323 9750
Japan (03) 38 16 1091

For further information
contact your local office

Location: <http://www.amersham.co.uk/life>

See our worldwide web site
for news, information and
local office contact details.



CyDye

Multicolour fluorescent
labelling reagents and
directly labelled biologicals

All goods and services are sold subject to the terms and conditions of sale of the company within the Amersham Group which supplies them. A copy of these terms and conditions is available on request. Amersham, Cy and CyDye are trademarks of Amersham International plc.
© Amersham International plc 1996. All rights reserved.

MCYD

Amersham LIFE SCIENCE

We're proving that positive thinking can cure disease.

The first step in accomplishing anything is the belief that it can be done. We believe there are cures for today's "incurable" diseases and we're working fast to find them. At Pharmacia & Upjohn we share this belief with more than 30,000 colleagues worldwide.

Along the way, we're finding better treatments for diseases like cancer, AIDS, diabetes, Parkinson's, growth disorders and arthritis.

Believing there is a cure is the first step in finding it. By sharing that belief you too can join in the fight. That's the power of positive thinking. That's the power of Pharmacia & Upjohn.



**Pharmacia
& Upjohn**

Get the Message...

...RNA Purification
Products from Promega.



PolyATtract® System 1000

- Isolation of high quality mRNA directly from cells or tissue
- Protocol complete in 45 minutes or less
- Start with as little as 10mg of tissue or 10^5 cultured cells
- Advanced magnetics-based isolation: no organic solutions

PolyATtract® Series 9600™ mRNA Isolation System

- Simultaneous isolation of mRNA from up to 96 samples
- From tissue or cells to cDNA in 3 hours

PolyATtract® System I-IV


- "Get the Message" from total RNA
- Magnetics-based isolation of mRNA from 100µg to 5mg of total RNA

RNAgents® Total RNA Isolation System


- The best value in total RNA isolation
- Demonstrated performance in RT-PCR and Northern blotting
- Protocol complete in 90 minutes with as little as 5mg starting material

PolyATtract and RNAgents are trademarks of Promega Corporation and are registered with the U.S. Patent and Trademark Office. Series 9600 is a trademark of Promega Corporation.

Please contact Promega or your nearest Branch Office or Distributor to receive your RNA Purification products folder.



RT-PCR of mRNA purified from Mouse Liver tissue using the PolyATtract® System 1000.



Northern blot of mRNA prepared using the PolyATtract® System 1000.



For more information about the RNAgents®, PolyATtract® Systems and other Promega products, visit our Home Page on the Internet at <http://www.promega.com>



Promega Corporation • Toll Free in USA 800-356-9526 • 2800 Woods Hollow Road • Madison, Wisconsin • 53711-5399 • Toll Free Fax in USA 800-356-1970 • Phone 608-274-4330 • Fax 608-277-2516
Internet <http://www.promega.com> • In Europe <http://www.euro.promega.com> • Promega products may also be purchased in the USA from **Fisher Scientific** • Toll Free 800-766-7000 • Fax 800-926-1166
Branch Offices: Australia 1 800 225 123 • France 0800 48 79 99 • Japan 03 3669 7981 • Netherlands (+31) (0) 71 5324244 • Switzerland 01 830 70 37 • United Kingdom 0800 378994 • **Joint Ventures:** China 21 6483 5136 • 379 393 5607
International Distributors: Argentina 1 384 7962 • Austria 1 891 460 • Belgium 0800 18098 • Brazil 031 275 1240 • Canada 800 267 7424 • Chile 2 334 0253 • China 10 6256 3159 • Colombia 1 310 8579 • Czech Republic 02 20610151
Denmark 44 94 88 22 • Egypt 2 245 1785 • Finland 0 350 9250 • France 1 45 33 67 17 • Germany 06221 50 24 70 • Greece 01 643 6138 • Hong Kong 2751 9488 • Hungary 01 251 0344 • Iceland see Sweden • India 11 684 6565 • Indonesia 21 723 1018 • Ireland 01 8426644 • Israel 8 9406 530 • Italy 055 5001871 • Japan 3 3270 0536 • Korea 2 3461 5911 • Luxembourg 31 71 324244 • Malaysia 3 718 3600 • Mexico 5 519 3463 • New Zealand 9 570 3200 • Norway 035 541999 • Poland 58 31 2807 • Portugal 1 458 1641 • Russia 095 135 4206 • Singapore 775 7284 • Slovak Republic 95 632 4729 • South Africa 11 476 9457 • Spain 93 404 52 14 • Sweden 0346 83050 • Taiwan 02 733 6877 • Thailand 2 294 7777 • Turkey 216 385 8321 • Uruguay see Argentina • Venezuela 2 265 0891

©1996 Promega Corporation. All Rights Reserved. Prices and specifications subject to change without prior notice. Rev.092796

Circle No. 42 on Readers' Service Card

The World's Best Science Journal Just Got Better


INTRODUCING

Full Text SCIENCE ONLINE

• IMMEDIACY

• INTERACTIVE

• SEARCHABILITY



SEE FOR YOURSELF

<http://www.sciencemag.org>

THE CREDIT CARD YOU'LL CARRY INTO THE NEXT CENTURY



AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE MBNA® *PLATINUM PLUS* CREDIT CARD

IT'S LIKE NO CREDIT CARD YOU CURRENTLY CARRY.

- A credit line up to \$100,000, No Annual Fee, and a low introductory 5.9% Annual Percentage Rate (APR) for cash advance checks and balance transfers†
- Toll-free MBNA *Platinum Plus* service 24 hours a day
- *Platinum Passage*—a 24-hour toll-free travel service that guarantees the lowest available published airfare at the time of booking
- Express delivery for card replacement at no additional cost. Free additional card for family and friends.
- Free Lost Card Registry
- Purchase protection against theft or damage
- \$1,000,000 Common Carrier Travel Accident Insurance at no additional cost*
- Free Year-End Summary of Charges
- Emergency cash and airline tickets, up to your available credit line, with free express delivery
- Credit line increase decisions in 15 minutes or less

Get the new standard in credit cards.

CALL TOLL-FREE 1-800-457-3714

(Please mention priority code FVHX when calling)

MBNA®
PLATINUM PLUS
t h e n e w s t a n d a r d



† The Annual Percentage Rate (APR) for purchases and ATM and Bank cash advances is 15.65%, which may vary. The current promotional APR offer for cash advance checks and balance transfers is 5.9% through your first five statement closing dates, commencing the month after your account is opened. When your minimum monthly payment is not received by the close of the first complete billing cycle following its Payment Due Date, or when the promotional offer expires, whichever occurs first, your APR for both new and outstanding cash advance balances (consisting of cash advance checks and balance transfer transactions) will be calculated using the Variable-Rate Information disclosures accompanying your card. The current indexed APR for cash advance checks and balance transfers is 15.65%, which may vary. Transaction fee for Bank and ATM cash advances: 2% of each cash advance (minimum \$2). Transaction fee for credit card cash advance checks: 1% of each cash advance (minimum \$2, maximum \$10). Transaction fee for the purchase of wire transfers, money orders, bets, lottery tickets, and casino gaming chips: 2% of each such purchase (minimum \$2). Cash advances and balance transfers cannot be used to pay off or pay down any MBNA account. We may allocate your monthly payments to your promotional APR balance(s) before your nonpromotional APR balance(s).

* Certain restrictions apply to these and other benefits described in the Portfolio of Services sent soon after your account is opened. The information about the costs of the card described in this advertisement is accurate as of 8/96. The information may have changed after that date. To find out what may have changed, call MBNA at 1-800-457-3714. TTY users, call 1-800-833-6262.

Platinum Passage travel services are provided to MBNA *Platinum Plus* Customers by, and are the responsibility of, an independently owned and operated travel agency. Visa is a federally registered service mark of Visa U.S.A. Inc., used pursuant to license. MBNA is a federally registered service mark of MBNA America Bank, N.A.
©1996 MBNA America Bank, N.A.
AAAS(ITTf)9/96

Circle No. 27 on Readers' Service Card

New Telomerase PCR ELISA Offers Simplified, Nonradioactive TRAP Assay for Measuring Telomerase, A Potential Marker for Cancer Research

Boehringer Mannheim is now offering a Telomerase PCR ELISA for the highly sensitive, nonradioactive detection of telomerase activity in extracts from cell cultures and tissue samples.

Telomerase as an important parameter in cancer research

Telomeres, the specialized DNA/protein structures at the end of eukaryotic chromosomes, contain tandemly repeated DNA sequences that are believed to protect genomic DNA from degradation and deleterious recombination events. During normal somatic cell proliferation, telomeric ends are progressively shortened with each replication cycle, which may play a role in limiting the proliferative capacity of normal cells. Germline cells, many tumor cells, and "immortalized" cell lines are believed to circumvent this telomere shortening using telomerase, a ribonucleoprotein that adds new repeats to the ends of chromosomes. Telomerase activity has recently been identified in many cancers (e.g., prostate cancers [1], advanced-stage breast cancers [2], neuroblastomas [3], and

primary lung cancer tissues [4]) that have been confirmed by other methods (e.g., histochemical staining). Thus, telomerase reactivation may allow cells to escape from the proliferative limitations of cellular senescence and could be further investigated as a potential marker for the development of malignant tumor cells.

Telomerase PCR ELISA improves upon previous TRAP assays

Telomerase activity is most frequently detected by the Telomeric Repeat Amplification Protocol (TRAP) of Kim *et al.* (5), in which the telomerase-reaction product is amplified by PCR. However, the conventional TRAP assay achieves full sensitivity only when performed with a hazardous radioactive label, and visualization of results requires time-consuming gel electrophoresis and autoradiography. The new Telomerase PCR ELISA[†] combines a one-step/one-tube TRAP assay with nonradioactive detection in a highly sensitive photometric ELISA (Figure 1).




Figure 1. Detection of telomerase activity with the Telomerase PCR ELISA.

- Step 1. Telomerase, if present, adds multiple 6-nucleotide telomeric repeats to a biotinylated synthetic primer.
- Step 2. The telomerase reaction product is amplified by PCR, using a biotinylated primer.
- Step 3. After denaturation, the PCR product hybridizes to a digoxigenin-labeled probe specific for the telomeric repeat.
- Step 4. The DNA hybrid binds to a streptavidin-coated microtiter plate, and anti-digoxigenin-peroxidase and TMB substrate generate a colored product measurable with a microplate reader.

Note: If desired, the TRAP reaction product from Step 2 can also be detected by the traditional gel electrophoresis method.




Figure 2. Correlation of results obtained with the Telomerase PCR ELISA and conventional, radioisotopic TRAP assays. Samples from known carcinomas, normal specimens (negative control), benign prostatic hyperplasia (BPH) specimens, and immortalized cell lines were tested with the Telomerase PCR ELISA and conventional, radioisotopic TRAP assays. In all sample types, the methods were able to identify the same number of samples featuring telomerase activity.

Data provided by M. Müller and R. Heicappell (6) and by H. J. Sommerfeld.

Easy-to-use ELISA delivers results in less time

The Telomerase PCR ELISA delivers results within 6 hours, eliminating the need for laborious, time-consuming gel electrophoresis and autoradiography techniques. Its ready-to-use TRAP reaction mix (telomerase substrate, amplification primers, nucleotides, *Taq* DNA polymerase, reaction buffer) eliminates the need to prepare multiple solutions and minimizes the risk of assay failure caused by contamination. Up to 96 TRAP reactions can be simultaneously analyzed with an ELISA plate reader.

Sensitive results correspond closely with those of radioactive TRAP assays

Besides avoiding the use of hazardous radioisotopes, the Telomerase PCR ELISA produces sensitive results comparable to those of the radioisotopic TRAP assay (Figure 2). The kit's optimized detection probe and hybridization conditions maximize both specificity and sensitivity.

Additionally, optimized primer sequences eliminate the need for "hot start" PCR while avoiding amplification artifacts (e.g., primer dimers).

The Telomerase PCR ELISA is currently available

The Telomerase PCR ELISA (96 tests; Cat. No. 1 854 666) is now available from Boehringer Mannheim Biochemicals representatives. Additional information can also be found at <http://biochem.boehringer-mannheim.com>.

References:

1. Sommerfeld, H. J. *et al.* (1996) *Cancer Research* **56**:218-222.
2. Hiyama, E. *et al.* (1996) *J. National Cancer Institute* **88**:116-122.
3. Hiyama, E. *et al.* (1995) *Nature Medicine* **1**:249-255.
4. Hiyama, K. *et al.* (1995) *J. National Cancer Institute* **87**:895-902.
5. Kim, N. W. *et al.* (1994) *Science* **266**:2011-2015.
6. Müller, M. *et al.* (1996) *Int. J. Oncology* **9**: in press.

[†]Licensed from Geron Corporation. Patents pending.

*Purchase of this product is accompanied by a limited license to use it in the Polymerase Chain Reaction (PCR) process in conjunction with a thermal cycler whose use in the automated performance of the PCR process is covered by the up-front license fee, either by payment to Perkin-Elmer or as purchased, i.e., an authorized thermal cycler.

Helping biomedical research become medical practice.

**BOEHRINGER
MANNHEIM**



Join the Partnership



MAXWELL MACKENZIE

"I have played a diverse set of roles in the scientific and engineering community, and believe that we are at our best when we are linking our work to the goals of society. The Center will enhance our ability to do that, and I encourage all of you to be a part of this effort."



H. Guyford Stever
Engineer
Former Director of the
National Science Foundation

For the latest on the Center for Science and Engineering, contact the AAAS Development Office at NEWBUILDING@AAAS.ORG or (202) 326-6636. Visit us on the web at WWW.AAAS.ORG.

Yes, I want to support the mission of AAAS and the Center for Science and Engineering

Honor a colleague, mentor, or family member with a gift ...

- \$1,000 will inscribe your name on the Center's Wall of Honor.
- \$5,000 also places your name on a handsome engraved plaque on the back of an auditorium chair.

Yes, I/We would like to support the **Center for Science and Engineering** with

☐ a gift of \$ _____ or

☐ a pledge of \$ _____

Enclosed is the first installment of \$ _____, with balance to be contributed:

☐ annually ☐ semi-annually ☐ quarterly beginning _____

I make this gift

☐ in my name

☐ in honor/memory of: _____

☐ and I will determine if my employer or my spouse's employer will match this contribution to AAAS.

Signature

Printed Name

Address

City, State, ZIP

Phone

Please return to:
AAAS Development Office
Room 730
1200 New York Avenue, NW
Washington, DC 20005




AMERICAN ASSOCIATION FOR THE
ADVANCEMENT OF SCIENCE

SC1196

CAMBRIDGE - Leaders in Scientific Scholarship

Physical Sciences



The MATHEMATICA® Book

Third Edition

Stephen Wolfram

With the release of the third edition of *Mathematica* software in 1996, Cambridge presents the manual that explains its every feature. The *Mathematica* Book contains a complete description of how to take advantage of the software's ability to solve everything from algebra and calculus problems to matrices and transformations.

1996 1400 pp. 58888-X Paperback \$44.95

A History of Modern Planetary Physics

Volume 1: Nebulous Earth:

The Origin of the Solar System and the Core of the Earth from Laplace to Jeffreys

Stephen G. Brush

Nebulous Earth follows the development of Laplace's Nebular Hypothesis, its connection with ideas about the interior of the Earth, and its role in the establishment of the "evolutionary" worldview that dominated science in the latter part of the nineteenth century.

1996 324 pp.
44171-4 Hardback \$54.95

Volume 2: Transmuted Past:

The Age of the Earth and the Evolution of the Elements from Lyell to Patterson

Stephen G. Brush

"An excellent work that covers a limited topic thoroughly."

—Science Books & Films

Transmuted Past follows the development of theories of stellar evolution and nucleosynthesis in the twentieth century and describes radiometric methods for estimating the age of the Earth.

1996 144 pp.
55213-3 Hardback \$44.95

Volume 3: Fruitful Encounters:

The Origin of the Solar System and of the Moon from Chamberlin to Apollo

Stephen G. Brush

Fruitful Encounters follows the eventual refutation of the encounter theory and the subsequent revival of a modernized Nebular Hypothesis. Brush also discusses the role of findings from the Apollo space program, especially the analysis of lunar samples, culminating in the establishment of the "giant impact" theory of the Moon's origin in the 1980s.

1996 366 pp.
55214-1 Hardback \$54.95

Available in bookstores or from

Chaotic Dynamics

An Introduction

Second Edition

**Gregory L. Baker and
Jerry P. Gollub**

Widely praised for its clarity and style, the previous edition of this text was the first to provide a quantitative introduction to chaos and nonlinear dynamics at the undergraduate level. This edition includes additional material on the analysis and characterization of chaotic data and applications of chaos.

1996 270 pp.
47106-0 Hardback \$70.00
47685-2 Paperback \$24.95

The Tectonic Evolution of Asia

An Yin and Mark Harrison, Editors

This volume is a collection of 21 contributions on the tectonic evolution of Asia. The book is divided into five parts: geodynamic models of the Cenozoic deformation in Asia, seismotectonics, geological evolution of the Himalaya-Karakorum Ranges, tectonics of the Cenozoic Indo-Asia collision, and Mesozoic-Paleozoic assembly of Asia.

World and Regional Geology Series 8

1996 678 pp.
48049-3 Hardback \$200.00

Biological Sciences

Designing Conservation Projects

Julian Caldecott

The author, a field conservationist who has worked with NGOs, governments, and donor agencies, provides a first-hand account of conservation in action across tropical Asia, Africa and Latin America. The book provides an analysis of where modern conservation has come from, where it is going, and what to do next.

1996 321 pp.
47328-4 Hardback \$64.95

The Trophic Cascade in Lakes

**Stephen R. Carpenter and
James F. Kitchell, Editors**

"I found the book to be clearly written, and appropriate for its target audience, researchers in ecology and resource management...The breadth of approaches brought to bear in the study is extremely impressive...I recommend this book as a fine case history of a pluralistic approach to an important ecological issue."

—Tim Wootton, Ecology

Cambridge Studies in Ecology

1996 399 pp.
43145-X Hardback \$85.00
56684-3 Paperback \$29.95

Taking Animals Seriously

Mental Life and Moral Status

David DeGrazia

"...David DeGrazia is one of the best and brightest of the current generation of philosophers to tackle this question, as he clearly demonstrates in this careful and lucid philosophical analysis of the issues and arguments. All those who are concerned about these issues would benefit greatly from reading this book."

—Andrew N. Rowan, Tufts University

School of Veterinary Medicine

1996 312 pp.
56140-X Hardback \$59.95
56760-2 Paperback \$18.95

Limb Regeneration

Panagiotis A. Tsonis

This is the first book that describes and analyzes the mechanisms of limb regeneration and patterning by incorporating the information obtained from older experiments with the recent advances in molecular biology.

Developmental and Cell Biology Series 31

1996 253 pp.
44149-8 Hardback \$74.95

Telling Lives in Science: Essays On Scientific Biography

**Michael Shortland and
Richard Yeo, Editors**

This collection of original essays explores for the first time the nature and development of scientific biography and its importance in forming our ideas about what scientists do, how science works, and why scientific biography remains popular. It is written by historians of science and science biographers in a scholarly but accessible style.

1996 309 pp.
43323-1 Hardback \$75.00

**CAMBRIDGE
UNIVERSITY PRESS**

40 West 20th Street, New York, NY 10011-4211 Call toll-free 800-872-7423.


The Edinburgh Building, Shaftesbury Road, Cambridge CB2 2RU, U.K.

Web site: <http://www.cup.org> MasterCard/VISA accepted. Prices subject to change.

Circle No. 1 on Readers' Service Card

NEW ENGLAND BIOLABS INTRODUCES...

ONE STEP NATIVE PROTEIN PURIFICATION



Purification of Maltose-Binding Protein (MBP) with >98% Purity in one step.

Std: NEB Broad Range Protein Marker
Lane 1: uninduced cell extract. Lane 2: induced cell extract showing expressed fusion protein (97 kDa). Lane 3: Flow through. Lane 4 and 5: MBP fractions eluted after inducing cleavage overnight at 4°C.

SO FAST IT'S SURREAL

Introducing **IMPACT™**, a new kit for obtaining pure native protein after a single affinity step. Developed from basic research at New England Biolabs on protein splicing, **IMPACT™** uses a new class of proteins called "inteins" modified to undergo directed cleavage at 4°C when induced by thiol-reagents such as DTT. Simply clone your target protein into the provided expression vector to create a *target protein-intein-chitin binding domain (CBD)* fusion. Bind the fusion protein to a chitin bead column, wash, induce on-column intein cleavage by adding DTT and then elute your native target protein at > 98% purity.

Advantages of the IMPACT™ system

- One step affinity purification
- Procedure carried out at 4°C
- Yields protein with native sequence
- No polyhistidine tags
- No proteolytic cleavage of fusion protein
- Requires only chitin beads
- Can label target protein C-terminus

Circle No. 51 on Readers' Service Card

For more information, call 1-800-NEB-LABS or visit the NEB web site: <http://www.neb.com>

THE NEXT STEP

- New England Biolabs Inc. 32 Tozer Road, Beverly, MA 01915 USA 1-800-NEB-LABS Tel. (508) 927-5054 Fax (508) 921-1350 email: info@neb.com
- New England Biolabs Ltd., Canada Tel. (800) 387-1095 (905) 672-3370 Fax (905) 672-3414 email: info@ca.neb.com
- New England Biolabs GmbH, Federal Republic of Germany Tel. (0130) 83 30 31 (06196) 3031 Fax (06196) 83639 email: info@de.neb.com
- New England Biolabs (UK) Ltd. Tel. (0800) 31 84 86 (01462) 420616 Fax (01462) 421057 email: info@uk.neb.com

DISTRIBUTORS: Australia (07) 5594-0299; Belgium (0800) 1 9815; Brazil (011) 66-3565; Denmark (31) 56 20 00; Finland (9) 584-121; France (1) 34 60 24 24; Greece (01) 5226547; Hong Kong 2649-9988; India (542) 311473; Israel (03) 5351205; Italy (02) 381951; Japan (03) 5820-9408; Korea (02) 556-0311; Mexico (5) 519-3463; Netherlands (033) 495 00 94; New Zealand (09) 418-3039; Norway 22 22 04 11; Singapore 4457927; Sweden (08) 30 60 10; Switzerland (061) 481 47 13; Taiwan (02) 8802913

NEW ENGLAND
BioLabs Inc.