

Science

19 March 1999

Vol. 283 No. 5409
Pages 1805-1968 \$8




AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE

NEW
ENGLAND
BIOLABS

finally

the end of primer walking,
random subcloning
and nested deletions

Overview of the
GPS-1 Genome Priming System



Transprimer is inserted randomly into the target DNA, generating a set of insertion products with only one Transprimer per molecule.

choose **GPS™-1 Genome Priming System** and complete your sequencing project in no time.

New England Biolabs introduces GPS-1, a fast, *in vitro* Tn7 transposon-based method for generating a population of DNA sequencing templates containing primer-binding sites (Transprimer element) at randomly interspersed positions. Simply supply your DNA of interest (0.08 µg per 20 µl reaction mix), follow a short protocol, transform into *E. coli*, miniprep colonies and then sequence.

- Transprimer is inserted randomly into the target DNA: no "hot spots" for insertion
- High efficiency: 10^1 - 10^3 insertion products recovered per reaction; entire protocol takes less than 90 minutes
- Unique rare-cutting restriction enzyme sites within Transprimer facilitate insert mapping
- DNA sequence can be obtained from both strands of target DNA using supplied primers
- Faster than primer walking, random subcloning and nested deletion methods

GPS-1 Genome Priming System

#7100

10 reactions

Includes: pGPS1 (Kan^r) and pGPS2 (Cam^r) Transprimer donor plasmids, TnsABC⁺ Transposase, control target plasmid, sequencing primers for the two ends of the Transprimer, and buffers.

ALSO AVAILABLE: GPS-M Mutagenesis System

#7101

10 reactions

Customize and insert a transposable element into target DNA at random locations

For further information on GPS-1, visit www.neb.com or call 1-800-NEB-LABS

- New England Biolabs Inc. 32 Tozer Road, Beverly, MA 01915 USA 1-800-NEB-LABS Tel. (978) 927-5054 Fax (978) 921-1350 email: info@neb.com
- New England Biolabs Ltd., Canada Tel. (800) 387-1095 (905) 672-3370 Fax (905) 672-3414 email: info@ca.neb.com
- New England Biolabs GmbH, Federal Republic of Germany Tel. 0800/BIOLABS (06196) 3031 Fax 0800/BIOLABX (06196) 83639 email: info@de.neb.com
- New England Biolabs (UK) Ltd. Tel. (0800) 31 84 86 (01462) 420616 Fax (01462) 421057 email: info@uk.neb.com


DISTRIBUTORS: Australia (07) 5594-0299; Belgium (0800) 1 9815; Brazil (011) 366-3565; Denmark (39) 56 20 00; Finland (9) 584-121; France (1) 34 60 24 24; Greece (01) 5226547; Hong Kong 2649-9988; India (542) 311473; Israel (03) 9229013; Italy (02) 381951; Japan (03) 5820-9408; Korea (02) 556-0311; Netherlands (033) 495 00 94; New Zealand 0800 807809; Norway 22 09 15 00; Singapore 4457927; Spain (03) 902 20 30 90; Sweden (08) 30 60 10; Switzerland (061) 486 80 80; Taiwan (02) 28802913

Circle No. 45 on Readers' Service Card

NEW ENGLAND
BioLabs Inc.

Adhesion Molecules and Cytokines

plus antibodies, ELISAs, and Primer Pairs™



R&D Systems' 700+ cytokine-related reagents are versatile tools that can be used in many applications. All are rigorously tested to ensure the highest quality, activity and reproducibility. Cytokines and adhesion molecules are usually supplied with BSA as a carrier protein. Antibodies are tested for many applications, including Neutralization, Western blot, and ELISA. Many of our antibodies have also been shown to perform in immunohistochemistry applications. PCR* Primer Pairs are a new addition to our range of reagents. Each pair has been carefully tested to identify optimum amplification conditions.

ADHESION MOLECULES	Human ELISA	Recombinant human protein	Anti-human polyclonal antibody	Anti-human monoclonal antibody	PCR Primer Pairs
sICAM-1	✓	✓	✓	✓	✓
sVCAM-1	✓	✓	✓	✓	✓
sE-Selectin	✓	✓	✓	✓	✓
sL-Selectin	✓	✓	✓	✓	✓
sP-Selectin	✓	✓	✓	✓	✓
sCD31 (PECAM-1)	✓	✓	✓	✓	✓

CYTOKINES	Human ELISA	Recombinant human protein	Anti-human polyclonal antibody	Anti-human monoclonal antibody	PCR Primer Pairs
IFN-γ	✓	✓	✓	✓	✓
IL-1α	✓	✓	✓	✓	✓
IL-1β	✓	✓	✓	✓	✓
IL-4	✓	✓	✓	✓	✓
IL-8	✓	✓	✓	✓	✓
MCP-1	✓	✓	✓	✓	✓
MIP-1α	✓	✓	✓	✓	✓
MIP-1β	✓	✓	✓	✓	✓
TNF-α	✓	✓	✓	✓	✓

*PCR is covered by US Patent Nos 4683195 and 4683202 assigned to Hoffmann-LaRoche

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES.

North America
R&D Systems, Inc.
614 McKinley Place NE
Minneapolis, MN 55413, USA.
Tel: 612 379-2956
Fax: 612 379-6580
info@rndsystems.com

Europe
R&D Systems Europe Ltd.
4-10 The Quadrant, Barton Lane
Abingdon, OX14 3YS, UK.
Tel: +44 (0)1235 551100
Fax: +44 (0)1235 533420
info@rndsystems.co.uk

Germany
R&D Systems GmbH
Borsigstrasse 7
65205 Wiesbaden, Germany
Tel: +49 (0)6122 90980
Fax: +49 (0)6122 909819
info@rndsystems.co.uk

Japan
Funakoshi, Co., Ltd.
9-7, 2-Chome
Hongo, Bunkyo-ku
Tokyo 113, Japan
Tel: 81-3 5684-1622
Fax: 81-3 5684-1633

International Distributors – Argentina: (54-1) 942-3654. Australia: (61-2) 9521-2177. Austria: 01 292 35 27.
Chile: (56-2) 671-9369. Greece: 031 322 525. Hong Kong: (852) 2649-9988. Israel: (972-9) 764-8787. Italy: 39 2 25 75 377.
Korea: (82-02) 569-0781. Mexico: (52-5) 612-0085. New Zealand: (64-9) 377-3336. Spain: 34 1 535 39 60. Poland: 48 22 720 44 54.
Portugal: 01 352 87 74. South Africa: 021 981 1560. Taiwan: (866-02) 368-3600. Venezuela: (58-2) 239-07-80.

Europe Free Phone – Belgique/België: 0800 10 468. Danmark: 80 01 85 92. France: 0800 90 72 49. Nederland: 060 225607.
Norge: 800 11033. Sverige: 020 79 31 49. Switzerland: 0800 55 2482.

Circle No. 46 on Readers' Service Card

R&D
SYSTEMS


1-800-343-7475

 www.rndsystems.com

Science

www.sciencemag.org

COVER Molecular data indicate that acoel flatworms are not an order of the phylum Platyhelminthes, as previously thought, but are members of the earliest divergent Bilateria (hence an early evolutionary branch of triploblastic Metazoa). Acoels are characteristic in their lack of anatomical features, and few are as attractive as these, as yet unclassified, specimens (~6 mm long) found on coral at Ngerdwais, a coastal lagoon reef in Palau. [Image: Coral Reef Research Foundation, Palau]



1835

Buried treasures

NEWS

NEWS OF THE WEEK

- 1822 **HUMAN GENOME: Academic Sequencers Challenge Celera in a Sprint to the Finish**
- ▼1823 **EVOLUTION: From a Flatworm, New Clues on Animal Origins**
1919
- ▼1825 **PLANT BIOLOGY: Key Molecular Signals Identified in Plants**
1911
- ▼1826 **GEODYNAMICS: A Lava Lamp Model for the Deep Earth**
1881
1885
1888
- 1827 **ASTRONOMY: NASA Plans Earlier Hubble Rescue**
- 1827 **EUROPEAN UNION: Cresson Resigns in Wake of Fraud Report**
- 1828 **HUMAN EVOLUTION: Genetic Study Shakes Up Out of Africa Theory**

- 1829 **RESEARCH ETHICS: Sweden Considers More Oversight of Research**
- 1829 **AIDS VACCINES: Nabel to Head NIH Vaccine Research Center**
- 1831 **CANCER RESEARCH: A Surprising Partner for Angiostatin**

NEWS FOCUS

- 1832 **NEUROPSYCHOLOGY: Nurture Helps Mold Able Minds**
- 1835 **ARCHAEOLOGY: Restorers Reveal 28,000-Year-Old Artworks**
- 1836 **ASTRONOMY: Submillimeter Astronomy Reaches New Heights**
- 1837 **MICROBIOLOGY: Forging a Link Between Biofilms and Disease**


RESEARCH

RESEARCH ARTICLE

- 1877 **The Origin of Chaos in the Outer Solar System** N. Murray and M. Holman

REPORTS

- ▼1881 **Compositional Stratification in the Deep Mantle** L. H. Kellogg, B. H. Hager, R. D. van der Hilst
1826
1885
1888
- ▼1885 **Compositional Heterogeneity in the Bottom 1000 Kilometers of Earth's Mantle: Toward a Hybrid Convection Model** R. D. van der Hilst and H. Kárason
1826
1881
1888
- ▼1888 **Dipping Low-Velocity Layer in the Mid-Lower Mantle: Evidence for Geochemical Heterogeneity** S. Kaneshima and G. Helffrich
1826
1881
1885
- 1892 **Chemical Transformations in Individual Ultrasmall Biomimetic Containers** D. T. Chiu, C. F. Wilson, F. Ryttsén, A. Strömberg, C. Farre, A. Karlsson, S. Nordholm, A. Gaggar, B. P. Modi, A. Moscho, R. A. Garza-López, O. Orwar, R. N. Zare
- 1895 **Nonequilibrium Self-Assembly of Long Chains of Polar Molecules in Superfluid Helium** K. Nauta and R. E. Miller



1905

Pumping electrons in mesoscopic systems

- 1897 **Epitaxial Cubic Gadolinium Oxide as a Dielectric for Gallium Arsenide Passivation** M. Hong, J. Kwo, A. R. Kortan, J. P. Mannaerts, A. M. Sergeant
- 1900 **Degradation Mechanism of Small Molecule-Based Organic Light-Emitting Devices** H. Aziz, Z. D. Popovic, N.-X. Hu, A.-M. Hor, G. Xu

DEPARTMENTS

NETWATCH
1811

THIS WEEK IN SCIENCE
1813

SCIENCESCOPE
1825

RANDOM SAMPLES
1841

CONTACT SCIENCE
1847

NEW PRODUCTS
1939



AMERICAN
ASSOCIATION FOR THE
ADVANCEMENT OF
SCIENCE

SCIENCE (ISSN 0036-8075) is published weekly on Friday, except the last week in December, by the American Association for the Advancement of Science, 1200 New York Avenue, NW, Washington, DC 20005. Periodicals Mail postage (publication No. 484460) paid at Washington, DC, and additional mailing offices. Copyright © 1999 by the American Association for the Advancement of Science. The title SCIENCE is a registered trademark of the AAAS. Domestic individual membership and subscription (51 issues): \$110 (\$62 allocated to subscription). Domestic institutional subscription (51 issues): \$325; Foreign postage extra: Mexico, Caribbean (surface mail) \$55; other countries (air assist delivery) \$90. First class, airmail, student, and emeritus rates on request. Canadian rates with GST available upon request, GST #1254 88122. Publications Mail Agreement Number 1069624. Printed in the U.S.A.

SCIENCE'S COMPASS

EDITORIAL

- 1847 **The Universality of Science** I. A. Lerch

LETTERS

- 1849 **Science Over Politics** R. P. Lanza *et al.* **Basic Research in China** S. Zhang. **Free Electrons?** J. M. Williams. **Monkey Numeration** P. L. Stocklin. **Response** E. M. Brannon and H. S. Terrace. **Corrections and Clarifications**

POLICY FORUM


- 1855 **CONSERVATION: Sustaining Natural and Human Capital: Villagers and Scientists** W. M. Getz, L. Fortmann, D. Cumming, J. du Toit, J. Hilty, R. Martin, M. Murphree, N. Owen-Smith, A. M. Starfield, M. I. Westphal
- 1856 **CONSERVATION: Capitalizing on Nature: Protected Area Management** A. Inamdar, H. de Jode, K. Lindsay, S. Cobb

BOOKS ETAL.

- 1858 **ECOLOGY: Untangling Ecological Complexity** *The Macroscopic Perspective* B. A. Maurer, reviewed by I. Hanski
- 1858 **Browsings**
- 1859 **EXHIBITIONS: Epidemic! The World of Infectious Disease**

1861

Creation of the moon



PERSPECTIVES

- ▼ 1860 **NEUROBIOLOGY: Dissecting Dendrite Dynamics** S. J. Smith
- 1861 **ORIGIN OF EARTH AND MOON: Colliding Theories** A. N. Halliday and M. J. Drake
- 1863 **ORIGIN OF EARTH AND MOON: A Couple of Uncertain Age** F. A. Podosek
- ▼ 1864 **CONDENSED MATTER PHYSICS: Pumping Electrons** B. L. Altshuler and L. I. Glazman

TECH.SIGHT

- 1867 **DNA SEQUENCING: Sequencing the Genome, Fast** J. C. Mullikin and A. A. McMurray
- 1869 **SOFTWARE: SYSTAT 8.0**, reviewed by A. Rawitch
- 1871 **TechSightings**

SCIENCE ONLINE
www.scienceonline.org

SCIENCE
THE JOURNAL
www.sciencemag.org

SCIENCE NOW
DAILY NEWS SERVICE
www.sciencenow.org

NEXT WAVE
WEEKLY CAREER UPDATES
www.nextwave.org

GRANTSNET
RESEARCH FUNDING DATABASE
www.grantsnet.org

NEUROAIDS
EXPERIMENTAL WEB SITE
www.sciencemag.org/NAIDS

- 1903 **Phase-Separated Composite Films for Liquid Crystal Displays** V. Vorflusev and S. Kumar

- ▼ 1905 **An Adiabatic Quantum Electron Pump** M. Switkes, C. M. Marcus, K. Campman, A. C. Gossard

- 1908 **Structural Maturation of Neural Pathways in Children and Adolescents: In Vivo Study** T. Paus, A. Zijdenbos, K. Worsley, D. L. Collins, J. Blumenthal, J. N. Giedd, J. L. Rapoport, A. C. Evans

- ▼ 1911 **Signaling of Cell Fate Decisions by CLAVATA3 in Arabidopsis Shoot Meristems** J. C. Fletcher, U. Brand, M. P. Running, R. Simon, E. M. Meyerowitz

- 1914 **Crystal Structure of Human ZAG, a Fat-Depleting Factor Related to MHC Molecules** L. M. Sánchez, A. J. Chirino, P. J. Bjorkman

- ▼ 1919 **Acoel Flatworms: Earliest Extant Bilaterian Metazoans, Not Members of Platyhelminthes** I. Ruiz-Trillo, M. Riutort, D. T. J. Littlewood, E. A. Herniou, J. Baguña

- ▼ 1923 **Rapid Dendritic Morphogenesis in CA1 Hippocampal Dendrites Induced by Synaptic Activity** M. Maletic-Savatic, R. Malinow, K. Svoboda

- 1927 **Variability in Spike Trains During Constant and Dynamic Stimulation** A.-K. Warzecha and M. Egelhaaf

- 1931 **The Role of Local Actin Instability in Axon Formation** F. Bradke and C. G. Dotti

- 1935 **Reversible Conversion of Monomeric Human Prion Protein Between Native and Fibrillogenic Conformations** G. S. Jackson, L. L. P. Hosszu, A. Power, A. F. Hill, J. Kenney, H. Saibil, C. J. Craven, J. P. Waltho, A. R. Clarke, J. Collinge

TECHNICAL COMMENTS

- North American Carbon Sink** E. A. Holland and S. Brown; C. S. Potter and S. A. Klooster. **Response** S. Fan, M. Gloor, J. Mahlman, S. Pacala, J. Sarmiento, T. Takahashi, P. Tans

www.sciencemag.org/cgi/content/full/283/5409/1815a



1911

Signals for developmental fate

Change of address: allow 4 weeks, giving old and new addresses and 8-digit account number. Postmaster: Send change of address to *Science*, P.O. Box 1811, Danbury, CT 06813-1811. Single copy sales: \$8.00 per issue prepaid includes surface postage; bulk rates on request. Authorization to photocopy material for internal or personal use under circumstances not falling within the fair use provisions of the Copyright Act is granted by AAAS to libraries and other users registered with the Copyright Clearance Center (CCC) Transactional Reporting Service, provided that \$4.00 per article is paid directly to CCC, 222 Rosewood Drive, Danvers, MA 01923. The identification code for *Science* is 0036-8075/83 \$4.00. *Science* is indexed in the *Reader's Guide to Periodical Literature* and in several specialized indexes.

That's just what you get with our pre-packed HiTrap™ columns. Results. Time after time, giving you peace of mind and time to concentrate on other important aspects of your work.

New HiTrap™ columns for hassle-free purification

HiTrap frees you from the needless details of many purification procedures because we've done all the work. Pre-packed columns and easy-to-follow instructions ensure fast start-up and method optimisation. There are even recommendations for sample handling.

In addition to applications for His-tagged fusion proteins, monoclonal and polyclonal IgG and DNA-binding proteins, the HiTrap range now includes immobilised lectins, immobilised streptavidin and media for purification of IgM and IgY.

Solvent delivery can be as simple as a syringe or as sophisticated as an ÄKTA™ design system. No other range offers such a choice of application from affinity to desalting.

Give us a call: in Europe +44 (0)1494 544550; in the US 1-800 526 3593; in Japan +81 3 5331 9336 and from the rest of the world +44 (0)1494 544100.

Or visit our web site: www.apbiotech.com/hitrap
— it'll be an eye-opener.



Blind you with science?
I'd rather dazzle
you with results.


Circle No. 52 on Readers' Service Card

HiTrap™ and ÄKTA™ are trademarks of Amersham Pharmacia Biotech Limited or its subsidiaries. Amersham is a trademark of Nycomed Amersham plc. Pharmacia and Drop Design are trademarks of Pharmacia & Upjohn Inc. All goods and services are sold subject to the terms and conditions of sale of the company within the Amersham Pharmacia Biotech group which supplies them. A copy of these terms and conditions is available on request. (c) Amersham Pharmacia Biotech UK Limited 1998. All rights reserved. Amersham Pharmacia Biotech UK Limited, Amersham Place, Little Chalfont, Buckinghamshire, England HP7 250. Amersham Pharmacia Biotech AB, SE-751 84 Uppsala, Sweden. Amersham Pharmacia Biotech Inc., 800 Centennial Avenue, PO Box 1327, Piscataway, NJ 08855, USA.



amersham pharmacia biotech


QIAGEN® tools for maximum transfection efficiencies



Advanced transfection reagents and ultrapure plasmid DNA

QIAGEN is the only company offering both the latest advances in transfection technology and the highest purity plasmid preparations for reproducibly high transfection efficiencies.

Lay the foundation for success with the proven plasmid DNA quality offered by QIAGEN®, QIAfilter™, and EndoFree™ Plasmid Kits.



Each bar represents the average efficiency from 4 replicates of 2 different plasmid preparations each. Plasmid: pCMVβ; Host strain: DH5α™

Build up your transfection results with SuperFect™ and Effectene™ Transfection Reagents for outstanding efficiencies, high reproducibility, and low cytotoxicity.

The QIAGEN combination of plasmid purification kits and innovative transfection reagents provide the best transfection results possible.

Let the best tools work for you! Call QIAGEN today at 800-426-8157 ext. 230 for your free Transfection Reagent Selector Kit*.

* The Transfection Reagent Selector Kit, cat. no. 301399, contains 0.3 ml SuperFect Reagent and 0.3 ml Effectene Reagent. Offer valid until May 31, 1999 while supplies last. Offer valid in the following countries: Australia, Canada, France, Germany, Japan, Switzerland, UK, and USA. Offer limited to 1 free kit per customer. Void where prohibited by institutional policy.

Trademarks: QIAGEN®, QIAfilter™, EndoFree™, Effectene™, SuperFect™ (QIAGEN). DH5α™ (Life Technologies, Inc.) © 1999 QIAGEN, all rights reserved.

<http://www.qiagen.com>

Circle No. 57 on Readers' Service Card

Germany:
QIAGEN GmbH
Tel. 02103-892-240
Fax 02103-892-255

USA:
QIAGEN Inc.
Tel. 800-426-8157
Fax 800-718-2056

Australia:
QIAGEN Pty Ltd
Tel. 03-9489-3666
Fax 03-9489-3888

Canada:
QIAGEN Inc.
Tel. 800-572-9613
Fax 905-501-0373

France:
QIAGEN S. A.
Tel. 01-60-920-930
Fax 01-60-920-925

Japan:
QIAGEN K.K.
Tel. 03-5805-7261
Fax 03-5805-7263

Switzerland:
QIAGEN AG
Tel. 061-319-3031
Fax 061-319-3033

UK:
QIAGEN Ltd.
Tel. 01293-422-999
Fax 01293-422-922

DISTRIBUTORS: Austria/Hungary/Slovenia Austria (01) 889 18 19 Belgium/Luxemburg 0800-1-9815 Brazil (11) 543 1455 or 0800 55 1321 Central & South America USA (305) 828-3818 China (852)2896-6263 Czech Republic (02)4447 1239 Denmark 43 86 87 88 Finland (09)-804 551 Greece (01)-643 61 38 India (011)-542 1714 Israel 02 584 1111 Italy (055) 500 1871 Korea (02) 924-86 97 Malaysia (03)-7312099 Mexico USA (661) 294-7940 The Netherlands (033)-4950094 New Zealand (09)9806700 or 0800933966 Norway 22 90 00 00 Poland (071) 51 41 06 or 51 54 20 Portugal (1)-751 6000 Singapore 445 7927 Slovak Republic (07) 54 01 336 South Africa (021) 615166 Spain (93) 902.20.30.90 Sweden (08) 621 34 00 Taiwan (02) 2880 2913 Thailand (02) 412-5672 In other countries contact: QIAGEN GmbH, Germany



PLSTR02993US

A SINKING SLAB OF CRUST?

Global tomographic models have revealed velocity anomalies in the mantle that may represent remnants of subducted slabs. Kaneshima and Helffrich (p. 1888) used short-period seismic waves to image a small-scale, low-velocity layer in the lower mantle beneath the Mariana and Izu-Bonin subduction zones in the Western Pacific. The authors suggest that this small, previously undetected dipping layer [about 500 kilometer (km) by 300 km in extent but only 8 km thick] is an undeformed fragment of oceanic crust that was subducted about 170 million years ago along the paleo-Indonesian subduction zone. This interpretation implies that mantle mixing was too weak to destroy this small, isolated slab remnant.

CHAOS IN THE OUTER SOLAR SYSTEM


Previous numerical simulations have showed that orbits could be chaotic for bodies in the outer solar system, but the mechanism for producing chaos was not determined and it was also possible that the simulated chaos was a numerical artifact. Murray and Holman (p. 1877) derived an analytical solution to the problem and showed in numerical simulations that chaos in the outer solar system can be caused by interactions between the orbits of Jupiter, Saturn, and Uranus or of Saturn, Uranus, and Neptune. The simulated chaotic zones are very narrow, so there is about a 20% chance that the jovian planets are currently in a chaotic zone.

MOLECULAR CHAINS

Dipolar molecules might be expected to align and form chains if dipole-dipole interactions dominate. However, other forces, such as van der Waals forces, would favor the formation of more compact structures. Nauta and Miller (p. 1895) show that when the strongly dipolar hydrogen cyanide molecules self-assemble inside cold, superfluid liquid helium droplets, linear chains of up to eight molecules form. In the gas phase, such clusters are known to form more compact cyclical structures. It is unlikely that the weak interactions with the superfluid helium "solvent" change the relative stabilities of the structures significantly. Rather, the formation of structures that do not represent the global energy minimum is attributed to their assembly mechanism, which is strongly influenced by the dipole-dipole interactions.

A DEEPER MANTLE BOUNDARY?

Geochemical data have implied that Earth's mantle is layered, with the boundary between an upper and a lower convective system marked by a change in seismic velocities at a depth of 660 kilometers. However, recent geophysical observations seem to show that subducted ocean crust penetrated this boundary and extended to much greater depths (see the news story by Kerr). Van der Hilst and Kárason (p. 1885) show that recent tomographic models imply that the seismic structure of the mantle changes at a much greater depth of 1600 kilometers. Kellogg *et al.* (p. 1881) present results from a series of models to argue that this depth may mark the convective boundary in the mantle.



AN ELECTRON HEART

A bias voltage is generally used to create electron currents in devices. Switkes *et al.* (p. 1905; see the Perspective by Altshuler and Glazman) abandon this usual method of bias-driven electron flow and report on a mechanism in which the electrons are pumped between two electron reservoirs by means of a periodic asymmetric deformation of an intervening quantum dot. In this system, the flow of electrons is not quantized but depends on the magnitude and frequency of the deformation.

AN INSULATOR FOR GALLIUM ARSENIDE

Gallium arsenide (GaAs) has several properties that make it superior for optoelectronics compared to silicon (Si), including higher charge carrier mobilities and a wider, direct band gap that allows high-temperature operation. However, GaAs-based devices are less commonly used than their Si counterparts because of the notorious difficulties in fabricating defect-free insulating oxide barriers on GaAs. Hong *et al.* (p. 1897) describe the successful growth by molecular beam epitaxy of a thin (~15 angstroms) gadolinium-based oxide on GaAs that has a low density of interface states, exhibits low leakage currents across the interface, and is stable against surface stress.

LESS ASSEMBLY REQUIRED

The familiar liquid-crystal displays of laptop computers are often the costliest component because of the numerous steps associated with their fabrication. Vorflusev and Kumar (p. 1903) show that some steps could be eliminated by preparing liquid-crystal films through phase separation in a polymer-liquid-

crystal mixture induced by ultraviolet illumination. In particular, the alignment of the liquid crystal can be achieved with a single substrate, whereas two are usually required. In addition, the resulting cells switched their display state faster than conventional cells and exhibited a gradation in their optical response.

FORMING YOUNG MINDS

The formation of appropriate connections between neurons in the brain, for example, in the visual cortex, is influenced by a blend of nature (molecular specifications of neighbors) and nurture (incoming visual input after birth), but what of the connections between areas of the brain? Paus *et al.* (p. 1908) present structural magnetic resonance imaging data from a sample of 111 children, ages 4 to 17, and describe age-related changes in interareal fiber tracts, which were assessed as increases in white matter density. Maturation of the fronto-temporal pathway was significantly greater in the left hemisphere than in the right, a result consistent with the maturation of brain areas involved in language.

PEPTIDES PATTERN FLOWERS

The pattern of petals in a flower, and even the pattern of flowers in a cluster, is determined by the development of the shoot meristems. In the plant *Arabidopsis*, the *CLAVATA1* gene, which encodes a receptor-like molecule, and the *CLAVATA3* gene are both expressed during development of the meristem but in different layers of cells. Fletcher *et al.* (p. 1911; see the news story by Barinaga) show that the *CLAVATA3* gene encodes a

CONTINUED ON PAGE 1815

Time in a Bottle. Expression-Ready Human Genes.

Over 100
Human Kinases

Chances are, your research project is time-sensitive. Cloning and expressing the human genes you need will take weeks, and there are no guarantees. GeneStorm™ expression-ready human clones give you a faster, better way.

GeneStorm™ Clones are recombinant human genes cloned into a state-of-the-art mammalian expression vector. Each vector has practical features including a polyhistidine sequence for convenient purification and an epitope tag for simplified detection. Once cloned, screened for orientation, and transfected into mammalian cells, gene expression is meticulously verified

for each of the GeneStorm™ clones. All this is done so that GeneStorm™ clones are immediately ready to use in your project. GeneStorm™ clones are ideal for dozens of research applications, from screening chemical libraries to antigen generation. And there are hundreds of GeneStorm™ human clones available from functional groups you need like kinases, transcription factors, and receptors, to name a few.

To learn more about these indispensable time-saving research tools, contact Invitrogen or visit our continually-updated GeneStorm™ web site (www.invitrogen.com/genestorm).



Circle No. 31 on Readers' Service Card

European Headquarters:

Invitrogen BV
De Schelp 12, 9351 NV Leek
The Netherlands
Tel: +31 (0) 594 515 175
Fax: +31 (0) 594 515 312
email: tech_service@invitrogen.nl

Toll Free Phone Numbers:

Belgium 0800 111 73
Denmark 800 188 67
Finland 990 31 800 5345
France 00 31 800 5345
Germany 0130 8100 43
The Netherlands 0800 022 88 48
Norway 800 113 70
Sweden 020 795 369
Switzerland 0800 551 966
United Kingdom 0800 96 61 93

Distributors:

Australia 1 800 882 555
China 010 6255 3477
Hungary 01 280 3728
India 91 80 839 1453
Israel 02 652 2102
Italy 02 38 19 51
Japan 03 5684 1616
Korea 82 2 569 6902

Malaysia 03 432 1357

Poland 058 41 42 26
Portugal 01 453 7085
Singapore 65 779 1919
Slovak Republic 07 3707 368
Spain 03 450 2601
Taiwan 080 231 530
Thailand 246 7243

From all other countries, contact
our European headquarters at
+31 (0) 594 515 175.

United States Headquarters:

Invitrogen
1600 Faraday Avenue
Carlsbad, California 92008
Tel: 1-800-955-6288
Fax: 760-603-7201
Email: tech_service@invitrogen.com
<http://www.invitrogen.com>

THIS WEEK IN SCIENCE

CONTINUED FROM PAGE 1813

small protein that may serve as a messenger between cells, perhaps signaling to the *CLAVATA1* protein found in neighboring cells. Thus, a peptide signal may be critical to pattern formation of tissues derived from the meristem.

MHC MIMIC

How similar to a major histocompatibility complex (MHC) class I protein does a molecule have to be in order to look like one? Sánchez *et al.* (p. 1914) analyzed the crystal structure of ZAG, a protein that stimulates lipolysis in fat cells. Although only 30 to 40% identical in sequence, and unable to bind the class I light chain β_2 -microglobulin or peptides, ZAG retains the same general tertiary structure, down to a pocket available for ligand binding. This finding illustrates the versatility of the basic structure and raises questions about the function and evolution of the ancestral protein.

BECOMING AN ANIMAL

The evolution of animals from protozoans involved several radical changes, including the formation of multicellular structures, the shift from radial to bilateral symmetry, and the development of digestive cavities. Ruiz-Trillo *et al.* (p. 1919; see the cover and the news story by Pennisi) have used molecular analyses to clarify the position of certain curiously primitive worms in this family tree. Previously thought to be secondarily simplified flatworms of the Platyhelminthes order, worms of the Acoela group are now seen to represent the earliest step in evolution of bilaterally symmetrical organisms.

SWITCHING PRION STATES

The human prion protein, implicated in neurodegenerative disorders, is found in a nonpathological soluble conformation rich in α -helix and in a pathological, fibrillar conformation rich in β sheet. Jackson *et al.* (p. 1935) have identified conditions in vitro that allow soluble human prion protein to switch between an α -helix rich conformation and a β sheet rich-conformation. The retention of the "pathological" conformation in solution may be important in understanding both the pathology and the biochemistry of this important protein.

DENDRITE AND NEURITE DYNAMICS

As a neuron develops, one of its neurites elongates to become the axon, and the others become dendrites. Bradke and Dotti (p. 1931) found that neurites that had growth cones with low levels of assembled actin were most likely to become axons, and that by locally breaking down the actin cytoskeleton in a dendrite, could be induced to form an axon. Subsequent electrical activity can change the morphology and topology of the remaining dendrites. Maletic-Savatic *et al.* (p. 1923; see the Perspective by Smith) present evidence from two-photon scanning imaging, which shows that, after brief periods of high-frequency synaptic stimulation, the dendrites of postsynaptic neurons start to grow new filopodia. This sprouting requires stimulation at high frequencies and can be antagonized by blocking *N*-methyl-D-aspartate receptors, which links this response to activity-dependent neuronal development and functional plasticity.

TECHNICAL COMMENT SUMMARIES

North American Carbon Sink

The full text of these comments can be seen at www.sciencemag.org/cgi/content/full/283/5409/1815a

S. Fan *et al.* (Reports, 16 Oct., p. 442) used "atmospheric and oceanic carbon dioxide data and models" to estimate "the spatial distribution of terrestrial carbon dioxide uptake." They concluded that North America "is the best constrained continent," with a large "carbon sink."


E. A. Holland and S. Brown comment that "direct estimates of forest C uptake" suggest that the report may have overestimated the magnitude of the sink. C. S. Potter and S. A. Klooster also studied the question using "newly derived terrestrial C fluxes" predicted from "forward modeling" simulations and arrived at "merely a transient C sink pattern in North America."

Fan *et al.* respond with an "intercomparison" between recent studies of C uptake. They state that, unlike their report, the comments "estimate terrestrial C uptake by invoking specific causes. The difference between [the various] estimates is the subject of intense research and may be resolved by a growing network of atmospheric and ecological observations."


3 Reasons Life Scientists choose SigmaPlot.

More than 100,000 of your colleagues have discovered how SigmaPlot can help them analyze and present their data with clear, compelling graphs.

Create exact graphs for publication




Analyze your enzyme kinetics data or any raw data with ease




Curve-fit any of more than 100 carefully selected built-in equations or your own using the Regression Wizard.

Automate repetitive reports of your experiments



Record macros to graph and analyze your data and create customized reports.



Try it FREE!

www.spss.com/software/science/SigmaPlot

(800) 345-4740


SPSS Science
233 S. Wacker Dr., 11th Floor
Chicago, IL 60606-6307
+1 262.551.3000; +800.621.1393
Fax: +1 262.551.3668; +800.841.0064
Distributors worldwide:
www.spss.com/order

SPSS Science Software GmbH
Schimmelbuschstrasse 25
40699 Erkrath, Germany
+49 (0) 2104.9540
Fax: +49 (0) 2104.95410
France: +0800.90.37.55

SIGMAPLOT
Exact Graphs for Exact Science

YSP5AD-0399-S

A4740



**“Now I CAN
BUY MY
REAGENTS ONLINE.
AND SAVE
MORE TIME
FOR MY
RESEARCH.”**


chemdex.com
Search less. Discover more.

Location: <http://www.chemdex.com>

About Chemdex | Chemdex News | User Account

Product Search | Your Favorites | Track Orders | Support | Your Profile | Current Orders

28	29	30	31
----	----	----	----



Announcing chemdex.com

The Faster, Easier Way to Purchase Biological and Chemical Reagents.

Spend less time searching for reagents and more time pursuing your research. When you purchase your reagents through chemdex.com, you get:

- The world's largest source of biological and chemical reagents for life science research
- Ordering with the click of a mouse—in a secure environment
- The suppliers you know—and trust
- A powerful, precision search engine to quickly find what you need

Purchase your reagents faster. Easier. ONLINE.
At www.chemdex.com.



chemdex.com

**Search less.
Discover more.**

Chemdex Corporation • 3950 Fabian Way, Palo Alto, CA 94303 • 650-813-0300
© 1999 Chemdex Corporation. Chemdex is a registered trademark of Chemdex Corporation.

EDITOR-IN-CHIEF
Floyd E. Bloom

EDITOR
Ellis Rubinstein

MANAGING EDITOR
Monica M. Bradford

EDITORIAL

DEPUTY MANAGING EDITORS: Richard B. Gallagher (Biological Sciences), R. Brooks Hanson (Physical Sciences), Katrina L. Kelner (Compass); SUPERVISORY SENIOR EDITORS Linda J. Miller, Phillip D. Szurorn; SENIOR EDITORS Gilbert J. Chin, Pamela J. Hines, Barbara Jasny, Paula A. Kiberstis, L. Bryan Ray; ASSOCIATE EDITORS Lisa D. Chong, Beverly A. Purnell, Linda R. Rowan; EDITORIAL SUPPORT Candace Gallery, Carolyn Kyle, Elise Laffman, Patricia M. Moore, Anita Wynn; ADMINISTRATIVE SUPPORT Sylvia Kihara

SCIENCE'S COMPASS

SENIOR EDITORS/PERSPECTIVES Orla Smith, Julia Uppenbrink; ASSOCIATE BOOK REVIEW EDITOR Sherman J. Suter; CONTRIBUTING EDITORS Kevin Ahern, Richard Peters, Robert Sikorski, David F. Voss;

PUBLISHER

Richard S. Nicholson

ASSOCIATE PUBLISHER
Beth Rosner

MEMBERSHIP/CIRCULATION DIRECTOR
Michael Spinella

MEMBERSHIP/CIRCULATION

SUBSCRIPTION SERVICES: 202-326-6417; 800-731-4939

DEPUTY DIRECTOR Marlene Zendell; MEMBER SERVICES: MANAGER Michael Lung; SUPERVISOR Mary Curry; REPRESENTATIVES Laurie Baker, Pat Butler, Christine Ford, Mari Pope, Jantell Smith; MARKETING: MANAGER Scott Oser; COORDINATOR Lauri Sirois; EUROPE MANAGER Jane Pennington; ASSOCIATE Ruth Jackson; RESEARCH: MANAGER Renuka Chander; BUSINESS AND FINANCE: MANAGER Dwight Theall; ASSISTANT Susan Maxim; COMPUTER SPECIALIST Charles Munson

FINANCE AND ADVERTISING

BUSINESS AND FINANCE: BUSINESS MANAGER Deborah Rivera-Wienhold; SENIOR ANALYST Randy Yi; FINANCIAL ANALYST Lisa Donovan RIGHTS AND PERMISSIONS: ASSOCIATE Lincoln Richman; ASSISTANT Emilie David MARKETING: DIRECTOR John Meyers; ASSOCIATE Allison Pritchard ELECTRONIC MEDIA: MANAGER David Gillikin; COMPUTER SPECIALIST Wendy Green; PRODUCTION ASSOCIATES Mark Croatti, Ellen McGuire

PUBLICATION ASSISTANTS Brent Gendelman, Jeffrey Hearn; INFORMATION SPECIALIST Janet Kegg; LETTERS AND TECHNICAL COMMENTS: EDITOR Christine Gilbert; ASSOCIATE EDITOR Steven S. Lapham; PUBLICATION ASSISTANT Charlene King

NEWS

NEWS EDITOR Colin Norman; FEATURES EDITOR Tim Appenzeller; DEPUTY NEWS EDITORS Elizabeth Culotta (contributing editor), Jean Marx, Jeffrey Mervis, Richard Stone; NEWS WRITERS Martin Enserink, Michael Haggmann (intern), Constance Holden, Jocelyn Kaiser, Richard A. Kerr, David Kestenbaum, Andrew Lawler, David Malakoff, Eliot Marshall, Elizabeth Pennisi, Robert F. Service, Gretchen Vogel BUREAUS: BERKELEY, CA Marcia Barinaga (contributing correspondent); SAN DIEGO, CA Jon Cohen; CHICAGO, IL James Glanz; COPY EDITORS Linda B. Felaco, Daniel T. Helgerman; CONTRIBUTING CORRESPONDENTS Barry A. Cipra, Ann Gibbons, Charles C. Mann, Anne Simon Moffat, Virginia Morell, Evelyn Strauss, Gary Taubes, Ingrid Wickelgren; ADMINISTRATIVE SUPPORT Scherraine Mack, Fannie Groom

EDITING AND PROOFREADING

DIRECTOR Dawn McCoy; SUPERVISOR Cara Tate; SENIOR COPY EDITORS Cay Butler, Harry Jach, Barbara Ordway, Christine M. Pearce; COPY EDITORS Jeffrey E. Cook, Etta Kavanagh, Jason Llewellyn, Joshua Marcy; COPY DESK Joi S. Granger, Monique Martineau, Ellen E. Murphy, Beverly Shields; ASSISTANT Kathy Libal

PRODUCTION

DIRECTOR James Landry; MANAGER Wendy K. Shank; ASSISTANT PRODUCTION MANAGER, INTERNET Elizabeth A. Harman; ASSOCIATES Vicki J. Jorgensen, Tara L. Kelly, Jessica K. Moshell, Rebecca Thomas

PRODUCT ADVERTISING

NATIONAL SALES MANAGER NORTHEAST AND E. CANADA Richard Teeling: 973-904-9774, FAX 973-904-9701 • MIDWEST/SOUTHEAST Elizabeth Mosko: 773-665-1150, FAX 773-665-2129 • WEST COAST/W. CANADA Neil Boylan: 415-673-9265, FAX 415-673-9267 • MID ATLANTIC AND U.S. INSIDE SALES Christopher Breslin: 410-273-1007, FAX 410-273-1591 • INTERNET SALES Laura Tellez 202-326-6599, FAX 202-682-0816 UK/SCANDINAVIA/France/ITALY/BELGIUM/NETHERLANDS Andrew Davies: (44) 7-071-226-216, FAX (44) 7-071-226-233 • GERMANY/SWITZERLAND/AUSTRIA Tracey Peers: (44) 1-260-297-530, FAX (44) 1-260-271-022 JAPAN Masah Yoshikawa: (81) 3-3235-5961, FAX (81) 3-3235-5852 • TRAFFIC MANAGER Carol Maddox; SALES ASSOCIATES Sheila Myers, Sandra Walls; ADMINISTRATIVE SUPPORT Jessica Tierney

RECRUITMENT ADVERTISING

SALES AND PRODUCTION OPERATIONS MANAGER Terri Seiter Azie U.S.: SALES MANAGER Gabrielle Boguslawski: 718-491-1607, FAX 202-289-6742; SALES SUPERVISOR Daryl Anderson; SALES REPRESENTATIVES Troy Benitez, Beth Dwyer, Bren Peters-Minnis, Kristin West-apher; ASSISTANTS Erika Bryant, Kathleen Clark; Christina Geiger PRODUCTION: SENIOR ASSOCIATE Jennifer Rankin; ASSOCIATES Chriss Borkowski, Elizabeth Lenox COPY EDITOR/PROOFREADER Chris Filiatreau UK/EUROPE: SALES MANAGER Debbie Cummings; SALES EXECUTIVE Sabine Lenuid; ASSISTANT Elisabeth Py: (44) 1-223-326-5000, FAX (44) 1-223-326-532 AUSTRALIA/NEW ZEALAND: Keith Sandell: (61) 02-9922-2977, FAX (61) 02-9922-1100 JAPAN: Masah Yoshikawa: (81) 3-3235-5961, FAX (81) 3-3235-5852

DEPUTY EDITORS

Philip H. Abelson (*Engineering and Applied Sciences*); John I. Brauman (*Physical Sciences*); Thomas R. Cech (*Biological Sciences*)

BOARD OF REVIEWING EDITORS

Frederick W. Alt <i>Children's Hospital, Boston</i>	David Clapham <i>Children's Hospital, Boston</i>	Peter Gruss <i>Max Planck Institute of Biophysical Chemistry</i>	Norman L. Letvin <i>Beth Israel Hospital, Boston</i>	Bert W. O'Malley <i>Baylor College of Medicine</i>	National Taiwan University
Don L. Anderson <i>California Institute of Technology</i>	Adrienne E. Clarke <i>University of Melbourne, Parkville</i>	Philip C. Hanawalt <i>Stanford University</i>	Harvey F. Lodish <i>Whitehead Institute for Biomedical Research</i>	Roy R. Parker <i>University of Arizona, Tucson</i>	Tomoyuki Takahashi <i>University of Tokyo</i>
Michael Ashburner <i>University of Cambridge</i>	F. Fleming Crim <i>University of Wisconsin, Madison</i>	Paul Harvey <i>University of Oxford</i>	Richard Losick <i>Harvard University</i>	Stuart L. Pimm <i>The Univ. of Tennessee, Knoxville</i>	Masatoshi Takeichi <i>Kyoto University</i>
Frank S. Bates <i>Univ. of Minnesota, Minneapolis</i>	Paul J. Crutzen <i>Max-Planck-Institut für Chemie</i>	M. P. Hassell <i>Imperial College at Silwood Park</i>	Seth Marder <i>University of Arizona</i>	Yeshayahu Pocker <i>Univ. of Washington, Seattle</i>	Keiji Tanaka <i>RIKEN Institute</i>
Stephen J. Benkovic <i>Pennsylvania State University</i>	James E. Dahlberg <i>University of Wisconsin Medical School, Madison</i>	Nobutaka Hirokawa <i>University of Tokyo</i>	Diane Mathis <i>Institut de Chimie Biologique, Strasbourg</i>	Martin Raff <i>University College London</i>	David Tilman <i>Univ. of Minnesota, St. Paul</i>
Alan Bernstein <i>Mount Sinai Hospital, Toronto</i>	Robert Desimone <i>National Institute of Mental Health, NIH</i>	Tasuku Honjo <i>Kyoto University</i>	Susan K. McConnell <i>Stanford University</i>	Douglas C. Rees <i>California Institute of Technology</i>	Robert T. N. Tjian <i>Univ. of California, Berkeley</i>
Michael J. Bevan <i>University of Washington, Seattle</i>	Hans Eklund <i>Swedish Univ. of Agricultural Sciences, Uppsala</i>	Susan D. Iversen <i>University of Oxford</i>	Anthony R. Means <i>Duke University Medical Center</i>	T. M. Rice <i>ETH-Hönggerberg, Zürich</i>	Yoshinori Tokura <i>University of Tokyo</i>
Seth Blair <i>University of Wisconsin, Madison</i>	Paul T. Englund <i>Johns Hopkins University School of Medicine</i>	Eric F. Johnson <i>The Scripps Research Institute</i>	Stanley Meisel <i>University of California, Davis</i>	David C. Rubie <i>Universität Bayreuth</i>	Derek van der Kooy <i>University of Toronto</i>
David E. Bloom <i>Harvard Institute for International Development</i>	G. Ertl <i>Max-Planck-Gesellschaft</i>	Hans Kende <i>Michigan State University</i>	Douglas A. Melton <i>Harvard University</i>	Erkki Ruoslahti <i>The Burnham Institute, CA</i>	Geerat J. Vermeij <i>University of California, Davis</i>
Piet Borst <i>The Netherlands Cancer Institute</i>	Richard G. Fairbanks <i>Lamont-Doherty Earth Observatory</i>	Jeffrey T. Kiehl <i>National Center for Atmospheric Research, Boulder</i>	Andrew Murray <i>Univ. of California, San Francisco</i>	Gottfried Schatz <i>Biozentrum, Basel</i>	Bert Vogelstein <i>Johns Hopkins Oncology Center</i>
Henry R. Bourne <i>Univ. of California, San Francisco</i>	Douglas T. Fearon <i>University of Cambridge</i>	Judith Kimble <i>University of Wisconsin, Madison</i>	Elizabeth G. Nabel <i>The Univ. of Michigan Medical Center</i>	Jozef Schell <i>Max-Planck-Institut für Zuchtungsforchung</i>	Gerhard Wegner <i>Max-Planck-Institut für Polymerforschung</i>
James J. Bull <i>University of Texas at Austin</i>	Harry A. Fozzard <i>The University of Chicago</i>	Stephen M. Kosslyn <i>Harvard University</i>	Shigetada Nakanishi <i>Kyoto University</i>	Ronald H. Schwartz <i>National Institute of Allergy and Infectious Diseases, NIH</i>	Arthur Weiss <i>Univ. of California, San Francisco</i>
Kathryn Calame <i>Columbia Univ. College of Physicians & Surgeons</i>	Roger I. M. Glass <i>Centers for Disease Control</i>	Michael LaBarbera <i>The University of Chicago</i>	Kim Nasmyth <i>Research Institute of Molecular Pathology, Vienna</i>	Terrence J. Sejnowski <i>The Salk Institute</i>	Zena Werb <i>Univ. of California, San Francisco</i>
Dennis W. Choi <i>Washington Univ. School of Medicine, St. Louis</i>	Peter N. Goodfellow <i>SmithKline Beecham, UK</i>	Antonio Lanzavecchia <i>Basel Institute for Immunology</i>	Roger A. Nicoll <i>Univ. of California, San Francisco</i>	Christopher R. Somerville <i>Carnegie Institute of Washington, Stanford, CA</i>	George M. Whitesides <i>Harvard University</i>
Joanne Chory <i>The Salk Institute</i>	Jack F. Greenblatt <i>University of Toronto</i>	Nicole Le Douarin <i>Institut d'Embryologie Cellulaire et Moléculaire du CNRS</i>	Staffan Normark <i>Swedish Institute for Infectious Disease Control</i>	Michael P. Stryker <i>Univ. of California, San Francisco</i>	Ian A. Wilson <i>The Scripps Research Institute</i>
			Kiyotaka Okada <i>Kyoto University</i>	Cliff Tabin <i>Harvard Medical School</i>	Alan P. Wolfe <i>National Institute of Child Health and Human Development, NIH</i>
				John Jen Tai	Martin Zatz <i>National Institute of Mental Health, NIH</i>

ART

DESIGN DIRECTOR Amy Decker Henry; ART DIRECTOR C. Faber Smith; ASSOCIATE ART DIRECTOR Elizabeth Carroll; SCIENTIFIC ILLUSTRATOR Katharine Sutliff; GRAPHICS ASSOCIATES Holly Bishop, Preston Morrighan, Darcel Pugh, Patricia M. Riehn; PHOTO RESEARCHER Leslie Blizard

SCIENCE INTERNATIONAL EUROPE

EDITORIAL: OFFICE HEAD Richard B. Gallagher; ASSOCIATE EDITORS Stella M. Hurtley, Ian S. Osborne, Peter Stern; EDITORIAL ASSOCIATE Belinda Holden; SCIENCE'S COMPASS: SENIOR EDITOR/PERSPECTIVES Julia Uppenbrink; NEWS: EDITOR Daniel Clery; CONTRIBUTING CORRESPONDENTS Michael Balter (Paris) Robert Koenig (Bern); SCIENCE'S NEXT WAVE, UK EDITOR: John MacFarlane; ADMINISTRATIVE SUPPORT Janet Mumford, Liz Ellis

ASIA

JAPAN NEWS BUREAU: Dennis Normile (contributing correspondent); CHINA REPRESENTATIVE Hao Xin; Pallava Bagla (contributing correspondent, India)

SCIENCE NOW: www.sciencenow.org

EDITOR Erik Stokstad

SCIENCE'S NEXT WAVE: www.nextwave.org

MANAGING EDITOR Wendy Yee; SENIOR EDITOR Nicole Ruediger WRITER Melissa Merti; CANADA EDITOR Charles Boulakia; ASSISTANT Suzanne Moore

AAAS BOARD OF DIRECTORS

RETIRING PRESIDENT, CHAIR M. R. C. Greenwood
PRESIDENT Stephen Jay Gould
PRESIDENT-ELECT Mary Lowe Good
TREASURER William T. Golden
EXECUTIVE OFFICER Richard S. Nicholson

Lewis M. Branscomb; Robert D. Goldman; Alice S. Huang; Sheila Jasanoff; Sally Gregory Kohlstedt; Marcia C. Linn; Neena B. Schwartz; David E. Shaw

Published by the American Association for the Advancement of Science (AAAS), *Science* serves its readers as a forum for the presentation and discussion of important issues related to the advancement of science, including the presentation of minority or conflicting points of view, rather than by publishing only material on which a consensus has been reached. Accordingly, all articles published in *Science*—including editorials, news and comment, and book reviews—are signed and reflect the individual views of the authors and not official points of view adopted by the AAAS or the institutions with which the authors are affiliated.

The American Association for the Advancement of Science was founded in 1848 and incorporated in 1874. Its objectives are to further the work of scientists, to facilitate cooperation among them, to foster scientific freedom and responsibility, to improve the effectiveness of science in the promotion of human welfare, to advance education in science, and to increase public understanding and appreciation of the importance and promise of the methods of science in human progress.

DHC
technology

Better Results
Easy Optimization
Higher Efficiencies
Double-Your-Satisfaction Guarantee


Plate: Image of NIH 3T3 transfected with green fluorescent protein (GFP).
Background: CHO-K1 transfected with β -Gal plasmid.

Finally, a transfection reagent that makes you smile, every time!

Introducing

GenePORTER™ Transfection Reagent*

GenePORTER™ reagent, which incorporates direct hydrophilic conjugation (DHC) technology*, is the latest innovation from the lab of Dr. Philip Felgner, inventor of lipofection. We at Gene Therapy Systems are so confident as to the performance of our product, we stand behind it with a double-your-satisfaction guarantee† on cell types listed. In a variety of cell lines, GenePORTER reagent consistently delivers higher levels of transfection than any other commercially available product. Without compromising transfection efficiency, this robust reagent performs in a wide range of conditions, including different ratios of plasmid and reagent. GenePORTER reagent is easy to use and does not require enhancers or special handling of cells, saving time, cost, and reagents. Order GenePORTER reagent today. The results will make you smile.



Transfected Cell Types†

HeLa S3	BHK-21
293	CHO-K1
Jurkat	CV-1
NIH 3T3	COS-1
B16-F0	COS-7
PC-12	

* Patents pending.

† If you do not get better transfection results following our protocol, kit purchase price will be refunded or double-the-kit price will be applied to your GTS, Inc. account toward future purchases. Offer valid in USA only and ends 3/31/99. Guarantee is limited to cell types in this ad. The purchased GenePORTER reagent must be returned using the GTS, Inc. return policy. This offer is limited to one T201007 or T201015 kit per customer.

GenePORTER™
Transfection
Reagent

75 reactions • T201007
150 reactions • T201015
750 reactions • T201075



To Order: 888-428-0558

Fax: 619-587-1499

10190 Telesis Court
San Diego, CA 92121, USA

For the complete product list
check out the GTS Website @


<http://www.genetherapysystems.com>

Circle No. 39 on Readers' Service Card

GTS International Distributors

Australia • ASTRAL +61 2 9540 2055 Austria • BIO-TRADE +43 1 889 18 19
Denmark, Finland, Norway, Sweden • KEBO Lab +46 8 621 35 07 France •
OZYME +1 30 85 92 92 Germany • Biozym Diagnostik GmbH +49 5152-9020
Italy • Duotech +39 02 3310 6630 Japan, Rep. of China, Korea • Funakoshi +81
3 5684 1622 Switzerland • Axon Lab AG +41 56 484 80 80 Taiwan • PROTECH,
Ltd. 886-2-23810844 United Kingdom • Lifescreeen LTD. +44 0 1923 241515

Nucleic Acid Labeling Made Easy



Single-step, nonradioactive labeling just got better!

Mirus
Label IT® Nucleic Acid Labeling Kits are now
available with $Cy3^{TM}$ and $Cy5^{TM}$ dyes.

Versatile: Label your nucleic acid with Rhodamine, Fluorescein, Digoxin, Biotin, AND NOW $Cy3$ and $Cy5$ dyes using just a single reagent.

Fast: Purified, labeled DNA or RNA is obtained in less than 2 hours.

Simple: Just mix the *Label IT* Reagent with your DNA or RNA.

Scalable: Label virtually any concentration of DNA or RNA without destroying your sample, and get quantitative recoveries of labeled product.

Efficient: The *Label IT* Reagents can add at least one label per 10 base pairs of DNA.

From Now On, *Label IT*!

Worldwide distributors:
Austria genexpress
Belgium Sopactem n.v.
Canada: BioCan Scientific (905) 828-2455
France: Takara Biomedical Europe S.A. 33-1-41-47-01-14
Germany: Mobitec GmbH 49-551-371062
Greece: Bega Ltd. 30-1-3300403




India Zonal, Ltd.
Israel Takara Shuzo Co., Ltd.
Japan JUMO Supply AG
Switzerland Cambridge Bioscience
United Kingdom Fisher Scientific
United States
Cyt3 and Cyt5 are trademarks of Amersham Pharmacia Biotech.
Label IT is a registered trademark of Mirus Corporation.

www.panvera.com
© 1999 Panvera Corporation. All rights reserved.

800.791.1400
In the USA and Canada

Panvera Corporation • 545 Science Dr., Madison WI 53711
Phone: 608.233.9450 • Fax: 608.233.3007
E-Mail: info@panvera.com




The New **Mini-Prep 24** for Automated Plasmid Mini-Preps

The Mini-Prep 24 is a fully automated bench-top instrument designed for purification of plasmid DNA directly from bacterial culture.

The instrument uses a revolutionary new method of nucleic acid purification based on modified agarose gel electrophoresis and subsequent recovery by electroelution.

The process utilizes premanufactured sample cassettes which allow for direct loading of up to 2 ml of culture.

Call now to learn how the New and Improved Mini-Prep 24 can provide you with great, high-quality DNA...while saving you a lot of time.



High Purity - sufficient for automated fluorescent and manual sequencing.

Easy Operation - begin prep with direct loading of bacteria culture - no centrifugation step saves you time.

Consistent Results - up to 6 µg of plasmid per ml.

Fast - up to 24 preps per hour, saving you time.

Quality - time and time again.

MacConnell
RESEARCH

1-800-466-7949


11339 Sorrento Valley Rd • San Diego, CA 92121 Phone: (619) 452-2603 Fax (619) 452-6753

www.macconnell.com


Circle No. 24 on Readers' Service Card

AmpliTaq Gold®

NEW
Reaction Buffer



Co-amplification of 8 human STR loci. 1,1, 2: male control DNA; 1,3, 4: female control DNA; L5: AmpliTaq negative control; L6, 7: male control DNA; L8, 9: female control DNA; L10: negative control.



Amplification of 10 copies HIV-1 DNA template in the presence of 1 µg human genomic DNA (from K-562 cells) using AmpliTaq DNA Polymerase with GeneAmp® Buffer II or AmpliTaq Gold DNA Polymerase using the NEW GeneAmp® PCR Gold Buffer.

**Now with GeneAmp®
10X PCR Gold Buffer**

- More robust amplification
- Reduced optimization
- Improves low copy detection

For PCR performance with higher yield, better specificity and more reliable results, discover AmpliTaq Gold® DNA Polymerase.

This new version of AmpliTaq® DNA Polymerase provides the specificity of Hot Start PCR, without all the extra steps. In most cases, you can substitute AmpliTaq Gold enzyme directly in existing PCR amplification protocols—without re-optimization.

You'll find AmpliTaq Gold DNA Polymerase saves time and money with dramatically lower drop-out rates, improved specificity, and easier multiplexing.

It also gives you consistently better PCR results. Because AmpliTaq Gold enzyme remains inactive until heated, conditions that lead to primer-dimer formation and mispriming are eliminated.

Circle No. 6 on Readers' Service Card

And of course, you have the continued assurance of knowing that AmpliTaq Gold DNA Polymerase is backed by PE Biosystems' exclusive PCR Performance Guarantee.

So discover AmpliTaq Gold DNA Polymerase, and discover high-performance PCR. To request information, call 1.800.345.5224. On the Internet, visit our home page at www.amplitaqgold.com, or e-mail pebio@perkin-elmer.com.

PE Biosystems

850 Lincoln Centre Drive
Foster City, CA 94404 USA
Tel. 800.345.5224
Tel. 650.638.5800
Fax. 650.638.5884
www.pebiosystems.com

Perkin-Elmer PCR reagents are developed and manufactured by Roche Molecular Systems, Inc., Branchburg, New Jersey, U.S.A.



AmpliTaq, AmpliTaq Gold and GeneAmp are registered trademarks of Roche Molecular Systems, Inc. The PCR process is covered by patents owned by Hoffmann-La Roche, Inc. and F. Hoffmann-La Roche Ltd. PE, PE Applied Biosystems and PE Biosystems are trademarks, and Perkin-Elmer is a registered trademark of The Perkin-Elmer Corporation. PE Applied Biosystems develops and produces its products in accordance with ISO 9000 quality system requirements.

Which best predicts **human** drug metabolism?

A or **B**

A.



Cryopreserved Human Hepatocytes



B.



Rat



Hint:

[It's not the rat.]

IVT OFFERS THE FOLLOWING PRODUCTS AND SERVICES:

Products (human and animal-derived):

• Cryopreserved hepatocytes • Plated hepatocytes • Microsomes & S9

Contract Research Services (in vitro):

• Human metabolite generation & identification • Drug-drug interaction
• CYP450 identification • CYP450 induction/inhibition
• Intestinal absorption (Caco-2) • Skin absorption • Toxicology
• Bioanalytical services (HPLC, LC/MS/MS)

Circle No. 11 on Readers' Service Card

Cryopreserved Human Hepatocytes

<http://www.hepatocytes.com>



IN VITRO
TECHNOLOGIES

1450 South Rolling Road
Baltimore, Maryland 21227

Phone: 410-455-1242


Fax: 410-455-1245

Internet: www.invitrotech.com

Email: info@invitrotech.com

BIO WHITTAKER

A **CAMBREX** Company



**Standardize
Your Research!
Call Your Clonetics®
Technical Specialist
Today.**

For Technical Information:
US & Canada (800) 852-5663
All other locations:
(301) 898-7025 ext. 7501

TO PLACE AN ORDER CALL:
(800) 344-6618
FAX: (301) 845-1008
E-MAIL: cs@biowhittaker.com
Internet:
www.biowhittaker.com, or
www.clonetics.com



©1999 BioWhittaker, Inc., A Cambrex Company.
All rights reserved.




Clonetics® Normal Human Neural Cell Systems!

• **Neural Progenitors** • **Astrocytes**


Neural Progenitor Features:

- ≥1,200,000 Frozen Cells Per Amp
- Serum Free Media, Growth/Differentiation Factors, Plating Reagents
- Express β Tubulin III and GFAP
- Differentiate into Neurons and Glia
- Test Negative for HIV-1, Hepatitis B & C (via PCR), Mycoplasma, Bacteria, Yeast and Fungi
- Quality Control Tested as a System of Cells and Medium
- Guaranteed Performance

**NEURAL SPHEROID PROLIFERATION IN
CLONETICS NEURAL PROGENITOR MEDIUM**



Note: Cells cease proliferation after differentiation into Neurons and Glia



▲ **Neural Progenitors one week
after plating**


**Normal Human Astrocytes at
40% confluence** ►

Applications:

- Gene Transfer
- Transplantation Studies
- Pharmacological Manipulations
- Brain Disease and Injury Studies
- Central Nervous System Physiology

Astrocyte Features:

- 500,000 Frozen Cells Per Amp
- Low Serum Media
- Available Proliferating in T-flasks or in Preseeded Plates
- Stain Positive for GFAP



Circle No. 10 on Readers' Service Card

Take Control of Expression!

Complete ControlTM

INDUCIBLE MAMMALIAN EXPRESSION SYSTEM

Now you can have precise, reliable and controlled mammalian gene expression with the Complete ControlTM inducible mammalian expression system.* This advanced ecdysone-inducible system combines innovative vector design with the tightest control available to give you a superior method for modulating expression in mammalian cells—without the leaky expression and pleiotropic effects from exogenous inducers that you find in other systems.

- Stable receptor cell lines available to get started quickly
- Optimal expression in a wide variety of tissues and cell lines
- Induction as high as 1,700-fold
- Specific and precise induction

Check out our web site at www.stratagene.com for all of the details on how you can take control of your experiments today!

STRATAGENE USA and CANADA

ORDER: (800) 424-5444 x3
TECHNICAL SERVICES: 800-394-1304

FAX: 619-535-0045

STRATAGENE EUROPE

Belgium, Germany, The Netherlands,
Switzerland, United Kingdom

ORDER: 00800 7000 7000

TECHNICAL SERVICES: 00800 7400 7400

FAX: 00800 7001 7001

Austria

ORDER: 0660 312 526


TECHNICAL SERVICES: 017 956 7036

FAX: 0660 312 527

INTERNET

eMAIL: techservices@stratagene.com

WEBSITE: www.stratagene.com



Double-stable CHO cells were created by first transfecting and selecting for the pERV receptor vector and then with a pEGSH expression construct containing the luciferase gene. These cells were induced with various amounts of ponasterone A and assayed for luciferase activity.

Circle No. 53 on Readers' Service Card



STRATAGENE[®]


www.stratagene.com

COMPLETE CONTROL[®] VECTOR KIT CAT.#217468

* Use of this system by commercial entities requires a commercial license.

Introducing Carl Zeiss Vision

The New Force In Imaging For Microscopy.



Come to Zeiss for a total solution to your imaging needs. The best microscopes with the best optics, the highest resolution cameras, and now the best imaging software.

The new **AxioVision** modular image acquisition and archiving system represents a new level of flexibility for digital microscopy. Two modes are available: the application mode for routine, and the expert mode for research. **AxioVision** is suitable for all biological and

materials applications. The system provides fully integrated control of the microscope and camera; efficient image enhancement; archiving of image, text, measuring data and all microscope parameters; and automatic report generation.

We also have specific programs for 3D deconvolution, confocal microscopy, materials analysis, and the **KS Elispot** for the automatic evaluation of T-lymphocyte immunoassays.

The Carl Zeiss Vision series covers the full range of counting, sizing, statistical evaluation, and image enhancement tools for all research and industrial applications.

For a free evaluation of your imaging needs, contact your Zeiss representative.

Call 800.233.2343 today.


Carl Zeiss Vision
total
solutions
for
imaging

Carl Zeiss, Inc.
Microscopy and Imaging Systems
One Zeiss Drive
Thornwood, NY 10594

800-233-2343
Fax 914-681-7446
micro@zeiss.com
www.zeiss.com
Circle No. 43 on Readers' Service Card



If you want
greater flexibility,
 we're on the same
wavelength.



Our SmartSpec™ 3000 spectrophotometer expands your options with a **wide wavelength range.**

- MOLECULAR CALCULATOR
- TOUCH-BUTTON OPERATION
- DNA/RNA, A₂₆₀/A₂₈₀
- PROTEIN ASSAYS
- WAVELENGTH SCANNING
- SIMPLE KINETICS
- O.D. 600

And that's just for starters. The new SmartSpec 3000 spectrophotometer lets you do a lot more for a lot less. It extends absorbance measurements from 200 to 800 nm, making it ideal for virtually any molecular biology or proteomic technique. It doubles as a molecular calculator for quick conversions to standard units. Better still, its intuitive, user-friendly interface lets you access any function with the touch of a button. See for yourself why SmartSpec 3000 is making waves. Contact your local representative for more information, or visit us on the web at discover.bio-rad.com and get a free lab timer!

BIO-RAD


Bio-Rad
Laboratories

Life Science
Group

Website www.bio-rad.com U.S. (800) 4BIORAD Australia 02 9914 2800 Austria (01)-877 89 01 Belgium 09-385 55 11 Canada (905) 712-2771
 China 86-10-62051850/51 Denmark 45 39 17 99 47 Finland 358 (0)9 804 2200 France 01 43 90 46 90 Germany 089 318 84-0
 Hong Kong 852-2789-3300 India (91-11) 461-0103 Israel 03 951 4127 Italy 39-02-216091 Japan 03-5811-6270 Korea 82-2-3473-4460
 Latin America 305-894-5950 Mexico 514-2210 The Netherlands 0318-540666 New Zealand 64-9-4152280 Russia 7-095-4585822
 Singapore 65-2729877 Spain 34-91-661-7085 Sweden 46 (0)8 627 50 00 Switzerland 01-809 55 55 United Kingdom 0800-181134 Sig 100798

Circle No. 29 on Readers' Service Card

SEQUENCHER™



The Best DNA Sequencing Software is NOW Available for Windows!

Now, Sequencer – the top software tool for DNA sequencing – can be used with whatever is on your desktop, WindowsNT, Windows95/98 or Macintosh. And you get the same high quality engineering and top-of-the-industry support you expect from Gene Codes Corporation.

Best known for its power and speed when handling shotgun and directed sequencing projects, Sequencer goes beyond simple fragment assembly and contig editing. It's the "complete software solution" that's become the industry standard – for major automated sequencing centers around the world, as well as for small labs doing short, manual sequencing. With Sequencer, your lab will increase productivity, reduce analysis time, and maximize the integrity and quality of your data.

See for yourself. Call for our demo program or visit our website and see what your colleagues have said about our software.

Gene Codes Corp., 640 Avis Dr, Ste 300, Ann Arbor, MI 48108
phone 734.769.7249 or 800.497.4939 • fax 734.769.7074
info@genecodes.com www.genecodes.com

- Fast, and accurate fragment assembly
- Heterozygote identification
- Vector and Transposon screening
- Translations, restriction analysis and ORF detection
- Scripting for simplifying repetitive tasks
- Compatibility with over a dozen sequence formats
- Exceptional interface to ABI, ALF, LI-COR, Vistra, Beckman and other auto-sequencers
- and a Scientist-friendly interface that makes it all easy to use!

Gene Codes Corporation

T C A G E N E
A G T C O D E S

Circle No. 40 on Readers' Service Card

DRUG DISCOVERY TOOLS



Contact us for our complete product listing or for
help with a novel reagent specific for your needs.

Visit www.chemicon.com • Phone: (909) 676-8080 • (800) 437-7500

Fax: (909) 676-9209 • E-Mail: custserv@chemicon.com

CHEMICON International, Ltd • (44) 181 8630 415
CHEMICON International, Germany • (49) 6192 2959 55

- Human Integrins Beta 1, 3, 5
- Functional Extracellular Matrix
- Antibodies to Integrins, MMPs and TIMPs
- Integrin Purification Gels

CHEMICON
International, Inc.

Circle No. 47 on Readers' Service Card

Human Tissues

normal and diseased

HYBRID-READY HUMAN TISSUES

Normal and diseased tissue sections mounted on slides,
ready for *in situ* hybridization or immunohistochemistry

International Distributors

Australia—Progen Industries Ltd. 7-3375-1888 • Austria—Bio-Trade 1-889-1819 • Belgium—Euro Biochem srl 10-43-70-60 • Denmark—Rie&Bernsen A-S 44-94-88-22 • Finland—Ya-Kemia Oy 9-350-3250 • France—TEBU 130-46-3900 • India—BioBusiness Development Agency 559-6820 • Israel—Megapharm 9-904596 • Italy—M-Medical 55-500-1871 • Japan—Takara Shuzo Co., Ltd. 77-543-7231 • Korea—BOHAN Biomedical 2-577-2002 • Netherlands—Omnilabo International BV 76-5795-795 • New Zealand—Intermed Scientific Ltd. 9-443-1284 • Norway—Bio-Test A/S 22-619-660 • Spain—bioNova científica, s.l. 91-551-5403 • Sweden—Labkemi AB 8-742-4200 • Switzerland—JURO Supply AG 41-226-20-60 • Taiwan—PROTECH 22-3810844 • UK—Cambridge Bioscience 1223-316855

For U. S. and Canada call Novagen at 800-526-7319 • For Germany call Calbiochem-NovaBiochem GmbH at 0800-6931-000
• For U.K. call Calbiochem-NovaBiochem U.K. at 0800-622935

Circle No. 37 on Readers' Service Card

 **Novagen**
800.526.7319

Immunohistochemical detection courtesy of Oncogene Research Products, Inc.

New Insights



LightcyclerTM

PCR perfected!

Get a new perspective on your PCR:

- ♦ High speed thermal cycling in glass capillaries (30 cycles in 20 minutes!)
- ♦ Real-time, on-line fluorescence monitoring
- ♦ Highly accurate quantification during log-linear PCR phase
- ♦ Mutation detection with continuous monitoring of hybridization kinetics

LightCycler



Internet: <http://biochem.roche.com>

Licensed from Idaho Technology, Inc.

LightCyclerTM is a trademark of and licensed from Idaho Technology, Inc.


ROCHE MOLECULAR BIOCHEMICALS

Austria (01) 277 870; Australia (02) 9899 7999; Belgium (02) 247 4930; Canada (514) 686 7050 (800) 361 2070; Denmark 45 13 82 32; France 04 76 76 30 86; Germany (0621) 759 856; Italy 039 247 4109-4181; Japan 03 5443 5284; Netherlands (036) 539 4911; Spain (93) 201 4411; Sweden (08) 404 8800; Switzerland 0 41/7 99 611; United Kingdom (0800) 521 578; USA (800) 428 5433

© 1998 Roche Diagnostics. All rights reserved



Circle No. 16 on Readers' Service Card



**Some see DNA
as a commodity.
We prefer a more
visionary approach.**

Call us idealistic, but we think your next custom oligo should be a thing of beauty.

Here's how we see it. Whether you need a couple of probes or primers, or a couple of thousand, you'll want impeccable quality and assured performance. That's why every Genosys oligo comes with its

1-800-853-3385
www.genosys.com

own certificate of Quality Assurance that includes a digitized PAGE analysis, quantitated yield, melting temperature and MW.

Of course, selection is important, too. So if your taste in oligos goes beyond plain vanilla primers, you'll be happy to know that our "standard flavors" include 6-FAM, HEX, TET, biotin, phosphorylation, fluorescein, amine labeling, Texas red and rhodamine. And if you're looking for S-oligos, 96-well plates, or even whole genes, look no further.

24-hour shipment? Routine for standard primers. Easy e-mail, web or fax ordering? Whichever you prefer. Knowledgeable, friendly tech support? Included free with every oligo.

Want to share the vision? Call today to

discuss your DNA requirements. For a really transcendental experience, let us handle your next order.

DNA micrographs are courtesy of Michael W. Davidson, director of the Optical Microscopy Division of the National High Magnetic Field Laboratory, a joint venture of The Florida State University, the University of Florida, and Los Alamos National Laboratory.

Circle No. 42 on Readers' Service Card

GENOSYS

Genosys Biotechnologies, Inc., The Woodlands, Texas, Phone: (281) 363-3693, E-Mail: info@genosys.com
Cambridge, U.K., Phone: (+44) (0)1223 839000, E-Mail: genosys@genosys.co.uk
Hokkaido, Japan, Phone: 81 133 73 5005, E-Mail: genosys@genosys.co.jp

OLIGOS DNA synthesis capacity



MASS OLIGOS™
 Deprotected, Desalted
 Adjusted Concentration
 Universal CPG
 (3' end secured)
 5 OD min.
 \$0.60 to \$1.00/base
 upon volume

*Order one or
 thousands of oligos.
 We'll produce them
 to the highest
 standards of quality
 and deliver them fast.*

With our proprietary
 UltraFast Parallel Synthesizer
 producing 24 oligos in about
 an hour, our post-synthesis

UltraFast Deprotector and robotic liquid handlers,
 our automated DNA synthesis plants can process more
 than 100,000 oligos per month. Bar code tracking of
 each oligo means complete control throughout produc-
 tion, quality control and packaging, and timely delivery.



A DIVISION OF GENSET

*Make GENSET OLIGOS
 YOUR
 GUARANTEED OLIGOS SOURCE.*

GENSET Corporation USA, CANADA, LATIN AMERICA
 ORDER- E-MAIL: oligos@gensetlj.com FAX: 1-800-551-5291
 PHONE: 1-800-995-0308

GENSET SA EUROPE, MIDDLE EAST, AFRICA
 ORDER- E-MAIL: oligos@genset.fr FAX: (33 1) 43 56 68 18
 PHONE: (33 1) 43 56 59 00

GENSET KK JAPAN
 ORDER- E-MAIL: oligos@po.ijnet.or.jp FAX: (81 3) 3585 5351
 PHONE: (81 3) 3585 5350

GENSET Singapore Biotech Pte Ltd. ASIA, OCEANIA
 ORDER- E-MAIL: oligos@genset.com.sg FAX: (65) 873 10 77
 PHONE: (65) 873 22 71

<http://www.genset.fr>

Circle No. 58 on Readers' Service Card

It's Here!



Our new **Mini-PROTEAN® 3** Mini Vertical Electrophoresis System is So Easy to Use, You Won't Want to Run With Anything Else.

For a run-through of what the Mini-PROTEAN 3 system can do for you, contact your local Bio-Rad representative or visit us on the web at discover.bio-rad.com.

BIO-RAD

**Bio-Rad
Laboratories**

Life Science
Group

Website www.bio-rad.com U.S. (800) 4BIORAD Australia 02 9914 2800 Austria (01)-877 89 01 Belgium 09-385 55 11 Canada (905) 712-2771
China 86-10-62051850/51 Denmark 45 39 17 99 47 Finland 358 (0)9 804 2200 France 01 43 90 46 90 Germany 089 318 84-0
Hong Kong 852-2789-3300 India (91-11) 461-0103 Israel 03 951 4127 Italy 39-02-216091 Japan 03-5811-6270 Korea 82-2-3473-4460
Latin America 305-894-5950 Mexico 514-2210 The Netherlands 0318-540666 New Zealand 64-9-4152280 Russia 7-095-4585822
Singapore 65-2729877 Spain 34-91-661-7085 Sweden 46 (0)8 627 50 00 Switzerland 01-809 55 55 United Kingdom 0800-181134 Sig 100798

Circle No. 28 on Readers' Service Card

Panorama™ Gene Arrays

EXPRESSION PROFILING

A Wider View of Expression Profiling


- Monitor the expression of hundreds to thousands of genes simultaneously
- Gene targets arranged in duplicate on nylon hybridization membranes
- Rapid analysis of differential gene expression
- Uses standard molecular biology techniques and equipment

S M A R T E R R E S O U R C E S F O R S M A R T E R R E S E A R C H .




E. coli ORFmers™

PCR* primers designed to amplify all E. coli ORFs. 4,290 primer pairs—representing each coding sequence—available in tubes or 96-well plates.



Panorama™ E. coli Gene Arrays
4,290 PCR-amplified open reading frames representing all E. coli protein coding sequences, arrayed in duplicate on nylon hybridization membranes.



****B. subtilis ORFmers™**

PCR primers designed to amplify all B. subtilis ORFs. 4,100 primer pairs—representing each coding sequence—available in tubes or 96-well plates.

GENOSYS

Genosys Biotechnologies, Inc. The Woodlands, Texas, Toll Free Phone: 877-260-0763 or (281) 363-3693, E-Mail: info@genosys.com

Nihon Genosys Biotechnologies, Inc., Hokkaido, Japan, Phone: 81 133 73 5005, E-Mail: genosys@genosys.co.jp

Genosys Biotechnologies (Europe) Ltd., Cambridge, U.K., Phone: (+44) (0)1223 839000, E-Mail: genosys@genosys.co.uk

Freephone (U.K.): 0800 269016 • Freephone (France): 0800 90 29 67 • Freephone (Germany): 0130 82 45 82

www.genosys.com

*The PCR process is covered by U.S. patents owned by Hoffman-LaRoche, Inc.

**Graphic reprinted by permission from Nature (230:252) ©1997 Macmillan Magazines Ltd.

Circle No. 41 on Readers' Service Card

QUESTIONS

Did the universe have a beginning? Is the universe designed? Are we alone?

In this conference eminent astronomers, cosmologists, philosophers, and religious scholars will review recent developments in astronomy and theoretical cosmology and discuss their significance for religious perspectives on the cosmos.

Cosmic Questions is made possible by support from the John Templeton Foundation, the Smithsonian Institution, the North American Montessori Teachers Association, *Science & Spirit Magazine*, and individual donors.

Cosmic Questions will feature more than 25 speakers. Here is a partial listing:

- **Steven Weinberg**, Nobel Laureate, University of Texas-Austin
- **Owen Gingerich**, Harvard-Smithsonian Center for Astrophysics
- **Joel Primack**, University of California-Santa Cruz
- **Rocky Kolb**, Fermi National Laboratory
- **Jill Tarter**, The SETI Institute
- **John Polkinghorne**, President Emeritus, Queens College, Cambridge

COSMIC

COSMIC QUESTIONS

When: **April 14-16, 1999**

Where: **National Museum of Natural History, Washington, DC**

For more information and a registration form, visit the Cosmic Questions Web Site at: www.aaas.org/spp/dspp/dbsr/events/cosmo or call **(202)326-6733**

Cosmic Questions Automated Registration Line
(Credit Card Only): **(202)326-7098**

Registration Fees		
	Non-Member	Member
Regular	\$150	\$135
Student*	\$90	\$80
* A limited number of student scholarships are available. For more information, call 202/326-6789.		



AMERICAN
ASSOCIATION FOR THE
ADVANCEMENT OF
SCIENCE


Registration Deadline: March 31

Circle No. 54 on Readers' Service Card

The
company
that
brought
you
Sequenase™
proudly
brings you

Thermo Sequenase™
The new standard
for manual cycle
sequencing.

... **Thermo**
Sequenase™
for "doubt-free"
³³P-Terminator
Sequencing.



Are your cycle sequencing results leaving you in doubt?
Strong stops, background sequence causing problems?

Thermo Sequenase Radiolabelled Terminator Cycle Sequencing Kit

The combination of USB's Thermo Sequenase Terminator Cycle Sequencing Kit with Amersham Pharmacia Biotech's radiolabelled terminators will give you sequencing results beyond a shadow of a doubt.

Thermo Sequenase - gives excellent performance at cycle sequencing temperatures. Helps eliminate secondary structure artifacts.

³³P-labelled terminators - yield detection of only truly terminated sequencing fragments. False terminations are not visible because these fragments are not labeled.

- Eliminate stops and BAFLs. (Bands Across all Four Lanes)
- Read the first base after the primer.
- Crisp, clean sequence.
- Detect heterozygotes and mutants more easily.
- Sequencing directly from bacterial colonies is possible.

No other combination for manual cycle sequencing gives more information per gel — long reads, accurate sequence. Doubt-free results.

Above: Elimination of bands across all four lanes (BAFLs) in plasmid template.

Fig. a) Using the Thermo Sequenase Radiolabelled Terminator Cycle Sequencing Kit.

Fig. b) Compared to standard ³²S-dATP cycle sequencing where a strong BAFL is observed.

To order within the United States contact USB directly at 1-800-321-9322 or e-mail us at customerserv@usbweb.com. For your authorized USB distributor outside the United States, please contact your local Amersham Pharmacia Biotech office.

USB Corporation, 26111 Miles Road
Cleveland, OH 4412
800-321-9322
www.usbweb.com



Rediscover the chemistry

Of course it leaves you speechless. It's designed for

visual

thinkers.

STARTING
@
\$3,395

For people whose world revolves around images not words, the new Silicon Graphics®

visual workstations offer the ultimate in graphics performance on Windows NT.® Featuring our Integrated Visual Computing (IVC) architecture with the Cobalt™ graphics chipset, these breakthrough machines move graphics data six times faster than other AGP 2X-based workstations*—a massive 3.2GB per second. Plus, I/O bandwidth is up to 12 times greater than other PC workstations.** So you can now get the outrageous 2D, 3D and digital media performance you expect from Silicon Graphics—for about the price of a mere clone. The Silicon Graphics 320™ workstation starts at just \$3,395. Or choose



the Silicon Graphics 540™ workstation†. Starting at \$5,995 and supporting up to four Intel® Pentium® II Xeon™ processors, it's the most scalable Windows NT workstation around. But enough words. After all, a picture is worth so much more.



Silicon Graphics 320 Visual Workstation

- Silicon Graphics Integrated Visual Computing architecture with Cobalt graphics chipset
- Supports up to two Intel® Pentium® II processors (up to 450MHz)
- Integrated 10/100 Ethernet, IEEE-1394†, USB and professional video I/O
- Ultra ATA or optional Ultra2 SCSI drives up to 28GB total capacity
- Microsoft® Windows NT® Workstation 4.0
- Add the SuperWide™ 17.3" Silicon Graphics 1600SW™ digital flat panel monitor (shown) for only \$2,495



SiliconGraphics

To get information, find a local reseller or to order, call 1 888 SGI-6315 or visit us at www.sgi.com/go/visual

*Bandwidth based on graphics to main memory. **I/O bandwidth as compared to 32-bit PCI. †For information on the 1999 availability of the Silicon Graphics 540 workstation, please call 1 888 SGI-6315. ‡Requires additional software under Windows NT Workstation 4.0. Prices quoted are for U.S. only.
© 1999 Silicon Graphics, Inc. All rights reserved. Silicon Graphics is a registered trademark, and the Silicon Graphics logo, Silicon Graphics 320, Silicon Graphics 540, SuperWide, Cobalt and Silicon Graphics 1600SW are trademarks, of Silicon Graphics, Inc. Intel, the Intel Inside logo and Pentium are registered trademarks, and Pentium II Xeon is a trademark, of Intel Corporation. Microsoft, Windows and Windows NT are registered trademarks of Microsoft Corporation. All other trademarks are property of their respective owners.

Circle No. 25 on Readers' Service Card