




Science

5 April 2002

Vol. 296 No. 5565
Pages 1-204 \$9



The
RICE
Genome



AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE


There's more to *PfuTurbo*[®] than high fidelity.

Robust PCR

PfuTurbo[®] DNA polymerase* consistently dominates the competition in high fidelity, PCR product yield and length. Competitors have introduced new PCR enzymes claiming high fidelity, but Stratagene's *PfuTurbo* DNA polymerase exhibits the highest fidelity with significantly better overall PCR performance. *PfuTurbo* DNA polymerase is a special formulation of high-fidelity cloned *Pfu* DNA polymerase* plus our patented ArchaeTect[™] polymerase enhancing factor*¹, which gives you greater yield over a broader range of target lengths and complexities.

Proven Performance

Pfu DNA polymerase has been helping scientists succeed for years. In fact, over 41% of researchers surveyed use *Pfu* DNA polymerase as their enzyme of choice.* Let Stratagene's expertise in high-fidelity PCR work for you. Choose *PfuTurbo* DNA polymerase for your high performance, high-fidelity PCR enzyme.



**PfuTurbo* DNA polymerase consistently dominates the competition in high fidelity, PCR product yield and length in side-by-side comparisons with competitor's DNA polymerases using manufacturer's recommended conditions.

Nothing beats *PfuTurbo*[®] DNA polymerase

- Highest fidelity
- Greater yields
- Longer targets
- Hotstart version available

STRATAGENE USA and CANADA
ORDER: (800) 424-5444 x3
TECHNICAL SERVICES: 800-894-1304

STRATAGENE EUROPE
Belgium, France, Germany,
The Netherlands, Switzerland,
United Kingdom

European Toll-Free Numbers
ORDER: 00800 7000 7000
TECHNICAL SERVICES: 00800 7400 7400
Austria
0800 312 526

**PfuTurbo* DNA polymerase: U.S. Patent Nos. 6,183,997 and 5,948,663 and 5,866,395 and 5,545,552 and patents pending; *Pfu* DNA polymerase: U.S. Patent Nos. 5,948,663, 5,866,395 and 5,545,552 and patents pending; ArchaeTect factor: U.S. Patent Nos. 6,183,997 and patents pending.

1. Hogrefe, H., et al. (2002) *Proc Natl. Acad. Sci. USA* 99, 596-601

‡ PhorTech International, DNA Amplification, 1999/2000 U.S. MSPPSA Series

www.stratagene.com

Stratagene invented high fidelity PCR.

PfuTurbo[®]
DNA POLYMERASE

PfuTurbo[®] DNA polymerase 100U **600250**
PfuTurbo[®] Hotstart DNA polymerase 100U **600320**



We reinvented the wheel... ...to accelerate your research




The PCR Workflow System

Automate your research with the **PCR Workflow System** from **Roche Applied Science**. Take advantage of our unique **LightCycler Sample Carousel** that integrates the **MagNA Pure LC** and **LightCycler Instruments**. Speed up your workflow to spend less time pipetting and more time advancing your research.

- Process virtually any sample material to isolate DNA, RNA, mRNA, and viral nucleic acids in less than two hours with the **MagNA Pure LC Instrument**.
- Perform real-time PCR analysis with the **LightCycler Instrument** in less than 40 minutes.
- Maximize flexibility in all applications (e.g., oncology research, infectious disease research, genotyping SNP).
- Generate reliable results with the highest reproducibility – day after day.

From Sample to Results in Three Hours


PCR Workflow System
Total: **185 min**
Hands on: **50 min**



Hands on ■
Automated ■

1. Extraction
2. PCR Setup
3. PCR Amplification and Detection


Manual Procedure
Total: **475 min**
Hands on: **220 min**



For detailed information, visit
www.roche-applied-science.com/pcr-workflowsystem
or contact your local sales representative

LightCycler and MagNA Pure are trademarks of a member of the Roche Group.
The technology used for the LightCycler is licensed from Idaho Technology, Inc.
© 2002 Roche Diagnostics Corporation. All rights reserved.

MagNA Pure LC Instrument



**LC Carousel
Centrifuge**



**LightCycler
Instrument**



Roche Diagnostics Corporation
Roche Applied Science
Indianapolis, IN

NEWS

NEWS OF THE WEEK

- 24 NIH DIRECTOR-DESIGNATE: Money, Mission, Management Top Zerhouni's Agenda
- 25 SMALLPOX VACCINES: New Cache Eases Shortage Worries
- ▼27
77
132 IMPACT HAZARD: Celestial Billiards Threaten Hit in 2880
- 27 SCIENCESCOPE
- 28 BIOETHICS: U.S. Questions Harvard Research in China
- 28 PROTEIN STRUCTURE: Harmless Proteins Twist Into Troublemakers

- ▼29
136 ANIMAL BEHAVIOR: Last Year's Food Guides This Year's Brood
- 31 PATENT LAW: Judge Casts Doubt on Scientist's Account
- 31 MIDDLE EAST: Science Foundation Sets Priority Areas

NEWS FOCUS

See The Rice Genome box on next page



- 39 MATHEMATICS: Erdős's Hard-to-Win Prizes Still Draw Bounty Hunters
- 41 ASTROPHYSICS: Closing In on the Cause of the Cosmos's Biggest Blasts
- 43 RANDOM SAMPLES



24

Zerhouni's agenda

SCIENCE'S COMPASS

LETTERS

Bowties and Brainiacs At the AAAS Meeting D. Helvarg. Regulating Natural Health Products J. C. Lashof, S. Margen, J. E. Swartzberg; I. Herskowitz. McClintock and Marriage S. B. McGrayne. Corrections and Clarifications

HISTORICAL ESSAY

- 49 PORTRAITS OF SCIENCE: "The Most Versatile Physicist of His Generation" M. J. Nye

BOOKS ET AL.

- 51 PHYSIOLOGY: *Lore of Running*, 4th ed. T. Noakes, reviewed by B. Heinrich
- 52 EDUCATION: *The Gold and the Blue* A Personal Memoir of the University of California, 1949–1967. Volume 1, Academic Triumphs C. Kerr, reviewed by R. L. Geiger

PERSPECTIVES

- ▼64
101 SOLAR PHYSICS: Order Amidst Turbulence J. Toomre
- ▼65
104
106 COLLOIDAL SYSTEMS: Playing Tricks with Designer "Atoms" D. Frenkel
- 66 POLYCRYSTALLINE MATERIALS: Grain Boundaries and Dislocations H. Van Swygenhoven
- ▼67
113 PALEOCLIMATE: Trans-Atlantic Climate Connections P. A. Baker
- REVIEW
- 69 BEHAVIORAL ECOLOGY: Breeding Together: Kin Selection and Mutualism in Cooperative Vertebrates T. Clutton-Brock
- 72 CONFLICT AND COOPERATION: Cooperation and Competition Between Relatives S. A. West, I. Pen, A. S. Griffin



69

Cooperation in animal societies

RESEARCH


BREVIA

- ▼77
27
132 Asteroid Hazard Mitigation Using the Yarkovsky Effect J. N. Spitale

REPORTS

- ▼101
64 Helioseismic Measurement of Solar Torsional Oscillations S. V. Vorontsov, J. Christensen-Dalsgaard, J. Schou, V. N. Strakhov, M. J. Thompson
- ▼104
65
106 Multiple Glassy States in a Simple Model System K. N. Pham, A. M. Puertas, J. Bergenholtz, S. U. Egelhaaf, A. Moussaïd, P. N. Pusey, A. B. Schofield, M. E. Cates, M. Fuchs, W. C. K. Poon

- ▼106
65
104 Layer-by-Layer Growth of Binary Colloidal Crystals K. P. Velikov, C. G. Christova, R. P. A. Dullens, A. van Blaaderen
- 109 Structure of Haloform Intercalated C₆₀ and Its Influence on Superconductive Properties R. E. Dinnebier, O. Gunnarsson, H. Brumm, E. Koch, P. W. Stephens, A. Huq, M. Jansen
- ▼113
67 A High-Resolution Paleoclimate Record Spanning the Past 25,000 Years in Southern East Africa T. C. Johnson, E. T. Brown, J. McManus, S. Barry, P. Barker, F. Gasse




101

As the Sun twists

RESEARCH (CONTINUED)

- ▼132 Asteroid 1950 DA's Encounter with Earth in 2880: Physical Limits of Collision Probability Prediction J. D. Giorgini, S. J. Ostro, L. A. M. Benner, P. W. Chodas, S. R. Chesley, R. S. Hudson, M. C. Nolan, A. R. Klemola, E. M. Standish, R. F. Jurgens, R. Rose, A. B. Chamberlin, D. K. Yeomans, J.-L. Margot
- ▼136 Evidence for the Effect of Learning on Timing of Reproduction in Blue Tits F. Grieco, A. J. van Noordwijk, M. E. Visser
- 138 Bile Acid Secreted by Male Sea Lamprey That Acts as a Sex Pheromone W. Li, A. P. Scott, M. J. Siefkes, H. Yan, Q. Liu, S.-S. Yun, D. A. Gage
- 141 Functional Annotation of a Full-Length *Arabidopsis* cDNA Collection M. Seki, M. Narusaka, A. Kamiya, J. Ishida, M. Satou, T. Sakurai, M. Nakajima, A. Enju, K. Akiyama, Y. Oono, M. Muramatsu, Y. Hayashizaki, J. Kawai, P. Carninci, M. Itoh, Y. Ishii, T. Arakawa, K. Shibata, A. Shinagawa, K. Shinozaki
- 145 Conserved Structure for Single-Stranded Telomeric DNA Recognition R. M. Mitton-Fry, E. M. Anderson, T. R. Hughes, V. Lundblad, D. S. Wuttke
- 148 The Interaction of Alba, a Conserved Archaeal Chromatin Protein, with Sir2 and Its Regulation by Acetylation S. D. Bell, C. H. Botting, B. N. Wardleworth, S. P. Jackson, M. F. White
- 151 Crystal Structure of the Extracellular Segment of Integrin $\alpha V\beta 3$ in Complex with an Arg-Gly-Asp Ligand J.-P. Xiong, T. Stehle, R. Zhang, A. Joachimiak, M. Frech, S. L. Goodman, M. A. Arnaout
- 155 Resolution of Lung Inflammation by CD44 P. Teder, R. W. Vandivier, D. Jiang, J. Liang, L. Cohn, E. Puré, P. M. Henson, P. W. Noble
- 158 Subnuclear Compartmentalization of Immunoglobulin Loci During Lymphocyte Development S. T. Kosak, J. A. Skok, K. L. Medina, R. Riblet, M. M. Le Beau, A. G. Fisher, H. Singh



THE RICE GENOME

EDITORIAL

- 13 DONALD KENNEDY
The Importance of Rice

NEWS FOCUS

- 32 RICE: Boiled Down to Bare Essentials
Rice Genomes Timeline
A Deal for the Rice Genome
- 36 BEIJING GENOMICS INSTITUTE: From Standing Start to Sequencing Superpower
- LETTERS

- 45 Why Finishing the Rice Genome Matters J. Leach, S. McCouch, T. Slezak, T. Sasaki, S. Wessler
Collaborating on the Rice Genome S. A. Goff

PERSPECTIVES

- 53 The Cereal of the World's Poor Takes Center Stage R. P. Cantrell and T. G. Reeves
- 54 World Poverty and Hunger—the Challenge for Science I. Serageldin
- 58 The Most Precious Things Are Not Jade and Pearls... P. Ronald and H. Leung
- 60 Opening the Door to Comparative Plant Biology J. Bennetzen

RESEARCH ARTICLES

- 79 A Draft Sequence of the Rice Genome (*Oryza sativa* L. ssp. *indica*) J. Yu, S. Hu, J. Wang, G. K.-S. Wong, S. Li, B. Liu, Y. Deng, L. Dai, Y. Zhou, X. Zhang, M. Cao, J. Liu, J. Sun, J. Tang, Y. Chen, X. Huang, W. Lin, C. Ye, W. Tong, L. Cong, J. Geng,

Y. Han, L. Li, W. Li, G. Hu, X. Huang, W. Li, J. Li, Z. Liu, L. Li, J. Liu, Q. Qi, J. Liu, L. Li, T. Li, X. Wang, H. Lu, T. Wu, M. Zhu, P. Ni, H. Han, W. Dong, X. Ren, X. Feng, P. Cui, X. Li, H. Wang, X. Xu, W. Zhai, Z. Xu, J. Zhang, S. He, J. Zhang, J. Xu, K. Zhang, X. Zheng, J. Dong, W. Zeng, L. Tao, J. Ye, J. Tan, X. Ren, X. Chen, J. He, D. Liu, W. Tian, C. Tian, H. Xia, Q. Bao, G. Li, H. Gao, T. Cao, J. Wang, W. Zhao, P. Li, W. Chen, X. Wang, Y. Zhang, J. Hu, J. Wang, S. Liu, J. Yang, G. Zhang, Y. Xiong, Z. Li, L. Mao, C. Zhou, Z. Zhu, R. Chen, B. Hao, W. Zheng, S. Chen, W. Guo, G. Li, S. Liu, M. Tao, J. Wang, L. Zhu, L. Yuan, H. Yang

92 A Draft Sequence of the Rice Genome (*Oryza sativa* L. ssp. *japonica*) S. A. Goff, D. Ricke, T.-H. Lan, G. Presting, R. Wang, M. Dunn, J. Glazebrook, A. Sessions, P. Oeller, H. Varma, D. Hadley, D. Hutchison, C. Martin, F. Katagiri, B. M. Lange, T. Moughamer, Y. Xia, P. Budworth, J. Zhong, T. Miguel, U. Paszkowski, S. Zhang, M. Colbert, W. Sun, L. Chen, B. Cooper, S. Park, T. C. Wood, L. Mao, P. Quail, R. Wing, R. Dean, Y. Yu, A. Zharkikh, R. Shen, S. Sahasrabudhe, A. Thomas, R. Cannings, A. Gutin, D. Pruss, J. Reid, S. Tavtigian, J. Mitchell, G. Eldredge, T. Scholl, R. M. Miller, S. Bhatnagar, N. Adey, T. Rubano, N. Tusneem, R. Robinson, J. Feldhaus, T. Macalma, A. Oliphant, S. Briggs



COVER 79, 92

Photograph of the Honghe Hani rice terraces in Yunnan Province, China. In this issue, two separate research groups report draft sequences of two strains of rice—*japonica* and *indica*. In addition, the Editorial, News Focus, Letters, and Perspectives highlight the significance of the rice genome to the world's population. [Image: Liwen Ma and Baoxing Qiu, Beijing Genomics Institute]

115

Pull-out section highlighting rice genomics

115 RICE GENOME POSTER

New on Science Express

Buzzing a comet: Glimpse of the nucleus



AMERICAN
ASSOCIATION FOR THE
ADVANCEMENT OF
SCIENCE

SCIENCE (ISSN 0036-8075) is published weekly on Friday, except the last week in December, by the American Association for the Advancement of Science, 1200 New York Avenue, NW, Washington, DC 20005. Periodicals Mail postage (publication No. 484460) paid at Washington, DC, and additional mailing offices. Copyright © 2002 by the American Association for the Advancement of Science. The title SCIENCE is a registered trademark of the AAAS. Domestic individual membership and subscription (51 issues): \$120 (\$66 allocated to subscription). Domestic institutional subscription (51 issues): \$390; Foreign postage extra: Mexico, Caribbean (surface mail) \$55; other countries (air assist delivery) \$85. First class, airmail, student, and emeritus rates on request. Canadian rates with GST available upon request, GST #1254 88122. Publications Mail Agreement Number 1069624. Printed in the U.S.A.

Change of address: allow 4 weeks, giving old and new addresses and 8-digit account number. Postmaster: Send change of address to SCIENCE, P.O. Box 1811, Danbury, CT 06813-1811. Single copy sales: \$9.00 per issue prepaid includes surface postage; bulk rates on request. Authorization to photocopy material for internal or personal use under circumstances not falling within the fair use provisions of the Copyright Act is granted by AAAS to libraries and other users registered with the Copyright Clearance Center (CCC) Transactional Reporting Service, provided that \$8.00 per article is paid directly to CCC, 222 Rosewood Drive, Danvers, MA 01923. The identification code for Science is 0036-8075/83 \$8.00. Science is indexed in the Reader's Guide to Periodical Literature and in several specialized indexes.

A N Y T H I N G I S P O S S I B L E

384

GIVES
YOU MORE

Greater throughput + greater efficiencies = MegaBACE 4000.

You'll get greater throughput with the new MegaBACE™ 4000 than you have ever achieved before. How much greater? Well, it processes 384 samples simultaneously. Put another way, it generates over 2.8 million bases per day. That's four times greater throughput than a 96-capillary instrument—which means you can sequence four times as many genomes in a week.

Receiving four times greater throughput from one system significantly boosts the efficiency of your operations. That means less labor, less space and less cost per base.

Time is money.

Go to **www.megabace.com** now for a special demonstration of the power of MegaBACE 4000.

Amersham Pharmacia Biotech
Life Sciences Division
Buckinghamshire, UK

 **Amersham**
Biosciences

Scienceonline

www.scienceonline.org

CONTENT HIGHLIGHTS AS OF 5 APRIL 2002

NOW!
Daily
coverage of
science and
science policy by
Science's news team
sciencenow
www.sciencenow.org

science magazine

www.sciencemag.org

SCIENCE EXPRESS

www.sciencexpress.org

Observations of Comet 19P/Borrelly by the Miniature Integrated Camera and Spectrometer Aboard Deep Space 1 L. A. Soderblom *et al.*

Images and spectra of the comet's nucleus, collected by the innovative, ion-propelled Deep Space 1 spacecraft, show that it is dark and dry, yet dynamic, spewing off fans and jets of dust and gas.

Genomic Instability in Mice Lacking Histone H2AX A. Celeste *et al.*

Absence of the DNA repair-associated protein, H2AX, leads to a poor DNA damage response and a range of genetic abnormalities.

Transcription Control by E1A and MAP Kinase Pathway via Sur2 Mediator Subunit J. L. Stevens *et al.*

The Mediator confers selectivity in its transcriptional activity through its various subunits that interact directly with specific transactivators.

TECHNICAL COMMENTS

The California Current, Devils Hole, and Pleistocene Climate

Herbert *et al.* (Research Articles, 6 July 2001, p. 71) found that sea surface temperatures (SSTs) in the area now dominated by the California Current warmed 10,000 to 15,000 years before deglaciation at the past five glacial maxima but did not rise before deglaciation south of the modern California Current front. Thus, they concluded, early warming in this area (including the trend recorded in vein calcite at Devils Hole, Nevada) likely represented not a global signal but a regional one tied to weakening of the California Current due to growth of the Laurentide ice sheet, and that the Devils Hole record "does not pose a fundamental challenge to the orbital ('Milankovich') theory of the ice ages." In a comment, Winograd, citing evidence of pre-deglaciation warming across a wide range of latitudes, asserts that the "intriguing notion linking California Current SSTs and the Laurentide ice sheet" is "open to question, as is [the] assertion" that the work resolves internal inconsistencies in the Milankovich theory. Herbert *et al.*, in response, offer several arguments in support of their "basic conclusion that regional climate processes are important and are expressed in the preserved paleoclimate record."

The full text of these comments can be seen at
www.sciencemag.org/cgi/content/full/296/5565/7a

science's next wave

www.nextwave.org

career resources for scientists

Global: Clinical Research Careers Edited by K. Cottingham

Throughout April: Find out how to break into the clinical research field, and read about the opportunities available for clinical scientists around the world.

US: GrantsNet's April News K. Cottingham

International research fellowships and grants from the NSF, a profile of an American Society for Microbiology Robert D. Watkins Minority Graduate Fellow, and the latest biomedical funding news.

Singapore: Clinical Research Odyssey in India S. Jayanthi, A. Kumar, N. Anand, V. Chhatwal

The clinical research scene in India emphasizes clinical trials and workforce development.

US: Raising Their Voices—The Emergence of Postdoc Grievance Procedures J. Gabriel

Postdocs are organizing to establish standardized grievance procedures, including postdoc representation on grievance committees.

KNOWLEDGE ENVIRONMENTS

science's sage ke

www.sageke.org

science of aging knowledge environment

The Werner Syndrome Helicase-Nuclease M. Fry

One protein, many mysteries.

A Jump-Start for Replication R. J. Davenport

When DNA-copying machinery conks out, Bloom's syndrome helicase brings repair proteins to the rescue.

science's stke

www.stke.org

signal transduction knowledge environment

Perspective: Calcium Release and Influx in Yeast—TRPC and VGCC Rule Another Kingdom M. Bonilla and K. W. Cunningham

Lessons in capacitative calcium entry from yeast.

Review: ADAP-ting TCR Signaling to Integrins E. K. Griffiths and J. M. Penninger

Inside-out signaling comes to the immune system.

Connections Map: Jak-STAT Pathway Pathway Authorities:

D. Aaronson and C. Horvath

A canonical view of cytokine signaling.

GrantsNet

www.grantsnet.org

RESEARCH FUNDING DATABASE

AIDScience

www.aidsience.com

HIV PREVENTION & VACCINE RESEARCH

Members Only!

www.AAASMember.org

AAAS ONLINE COMMUNITY

SCOPE

http://scope.educ.washington.edu

EXPLORING SCIENCE CONTROVERSIES

Functional Genomics

www.sciencegenomics.org

NEWS, RESEARCH, RESOURCES

ONLINE STAFF

SCIENCENOW EDITORS Martin Enserink, Laura Helmuth, Greg Miller, Erik Stokstad

SCIENCE'S NEXT WAVE EDITORIAL MANAGING EDITOR Crispin Taylor; EDITORS Eick von Ruschkowski (Germany), Kirstie Urquhart (UK), Ric Weibl (US); CONTRIBUTING EDITORS Lesley McCarney (Canada), Robert Metzke (Netherlands), Jennie Wong (Singapore); PROJECT EDITORS Jim Austin, Sibrina Collins, Katie Cottingham, Laure Haak; PRODUCTION ASSOCIATE Lily Han; MARKETING: MARKETING MANAGERS Karen Horting (Global), Hazel Crocker (Europe); PROGRAM DIRECTOR Lisa Kozlowski; MARKETING ASSOCIATE Angela Walker; PROGRAM ASSOCIATES Shajuan Martin, Tammy Minor

AIDSCIENCE SENIOR EDITOR Roberto Fernandez-Larsson; ASSOCIATE EDITOR Paula Werner

SCIENCE'S STKE EDITOR Bryan Ray; MANAGING EDITOR Nancy Gough; ASSOCIATE EDITORS Lisa Chong, John Nelson; PUBLICATIONS ASSISTANT Christopher Kenny

SCIENCE'S SAGE KE EDITORIAL DIRECTOR Kelly LaMarco; SENIOR NEWS EDITOR Evelyn Strauss; ASSOCIATE EDITOR R. John Davenport

ELECTRONIC MEDIA MANAGER Don Hemerway; INTERNET PRODUCTION MANAGER Betsy Hartman; ASSISTANT PRODUCTION MANAGER Wendy Stengel; SENIOR PRODUCTION ASSOCIATES Darcee O'Donnell, Lisa Stanford; ASSOCIATES Carla Cathey, Steve Kusek, Louis Williams; LEAD APPLICATIONS DEVELOPER Carl Saffell

PolyFect® Transfection Reagent

Fast and easy transfection, with optimized protocols for:

293
HeLa
COS-7
CHO
NIH/3T3

- Optimized protocols for COS-7, NIH/3T3, HeLa, 293, and CHO cells
- Fast procedure and easy handling
- No removal of complexes required



- Transfection in the presence of serum
- High efficiency
- Economical price

**Transfection has never been so simple —
order PolyFect Reagent today!**

Trademarks: QIAGEN®, PolyFect® (QIAGEN). © 2001 QIAGEN, all rights reserved.

www.qiagen.com

QIAGEN:

Australia
Tel. 03-9489-3666
Fax 03-9489-3888

Canada
Tel. 800-572-9613
Fax 800-713-5951

Germany
Tel. 02103-29-12400
Fax 02103-29-22022

Italy
Tel. 02-33430411
Fax 02-33430426

Switzerland
Tel. 061-319-3031
Fax 061-319-3033

UK
Tel. 01293-422-999
Fax 01293-422-922


France
Tel. 01-60-920-930
Fax 01-60-920-925

Japan
Tel. 03-5547-0811
Fax 03-5547-0818

USA
Tel. 800-426-8157
Fax 800-718-2056

Distributors:

Argentina Tecnolab S.A. (011) 4555 0010 Austria/Hungary/Slovenia R. u. P. MARGARITELLA Austria (01) 889 18 19 Belgium/Luxembourg Westlab b.v. 0800-1-9815 Brazil Uniscience do Brasil 011 3622 2320 China Gene Company Limited (852) 2896-6283 Cyprus Scientifica Ltd (02) 765 416 Czech Republic BIOCONSULT spol. s r.o. (02) 4447 1239 Denmark Merck EuroLab A/S 43 86 87 88 Egypt Clinicals 525 7212 Finland Merck EuroLab Oy (09) 804 551 Greece BioAnalytica S.A. (01) 640 03 18 India Genentix (011)-542 1714 or (011) 515 9346 Israel Westlab (Israel) Ltd. 08 66 50 813 or 1-800 20 22 20 Korea LRS Laboratories, Inc. (02) 924-86 97 Malaysia Research Biolabs Sdn. Bhd. (03) 7312099 Mexico Quimica Valaner S.A. de C.V. (5) 525 57 25 The Netherlands Westlab b.v. (033) 4950094 New Zealand Biolab Scientific Ltd. (09) 980 6700 or 0800 933 966 Norway Merck EuroLab AS 22 90 00 00 Poland Syngen Biotech Sp.z.o.o. (071) 351 41 06 or 0601 70 60 07 Portugal IZASA PORTUGAL, LDA (1) 424 73 64 Singapore Research Biolabs Pte Ltd 2731066 Slovakia Republic BIOCONSULT Slovakia spol. s r.o. (07) 50221 336 South Africa Southern Cross Biotechnology (Pty) Ltd (021) 671 5166 Spain IZASA, S.A. (93) 902.20.30.90 Sweden Merck EuroLab AB (08) 621 34 00 Taiwan TAIGEN Bioscience Corporation (02) 2880 2913 Thailand Theera Trading Co. Ltd. (02) 412-5672 In other countries contact: QIAGEN, Germany



TR052001SWW

THIS WEEK IN Science

edited by Phil Szuromi

Sun Twists

The Sun's 11-year cycle of magnetic activity is characterized not only by changes in sunspot concentration but also by latitudinal variations in the rotation rate of bands of plasma at the Sun's surface. Vorontsov *et al.* (p. 101; see the Perspective by Toomre) have used oscillations in the acoustic gravity waves in the Sun's interior collected by the Solar and Heliospheric Observatory (SOHO) to track these zonal bands over a period of almost 6 years into the Sun's interior. The entire convection zone of the Sun is involved in the differential rotation and latitudinal shifts of the zonal bands. Such measurements should help in connecting the observed surface features and the solar magnetic field with the deeper dynamo, which is thought to be located at the base of the convective zone.

Dodging Asteroids

Normally, it is difficult to predict the orbital path of asteroids beyond a few decades because of the nonlinear nature of the processes that perturb their orbits. Giorgini *et al.* (p. 132; see the news story by Kerr) have recognized that asteroid 29075 (1950DA) presents a special case. Extensive radar and optical data collected over decades, an inclined orbit that reduces perturbations, and a previously unrecognized gravitational resonance that modulates the uncertainty of its orbital path allow an accurate orbital path prediction that extends for hundreds of years—including a possible collision with Earth in March 2880. The largest uncertainty in the orbital path simulation is related to 29075's unknown surface characteristics because anisotropic thermal re-radiation of absorbed solar energy from the asteroid's surface (the Yarkovsky effect) can significantly alter the orbit over a long time. In Brevia, Spitale (p. 77) shows that we can potentially mitigate asteroid hazards by altering the asteroid's surface, and the Yarkovsky effect, enough to avoid a collision with Earth.

Into Africa

A record of biogenic silica accumulation from Lake Malawi assembled by Johnson *et al.* (p. 113; see the Perspective by Baker) suggests a connection between cold conditions in the Northern Hemisphere (between 11,000 and 25,000 years ago) and high diatom productivity caused by north winds over that part of equatorial

79
92

Going with the Grain


Two separate research groups report high-quality draft sequences of the rice genome that are expected to facilitate advances in the agriculture of this critical food grain. Yu *et al.* (p. 79) sequenced the indica variety of rice, and Goff *et al.* (p. 92) sequenced the japonica variety of rice. Comparisons will show how these two popular rice strains have diverged. The rice genome also provides a useful road map for investigating the larger genomes of related cereal grains such as wheat and maize [see also the report by Seki (p. 141) on *Arabidopsis* complementary DNA clones]. Related pieces include a pull-out wall chart that gives a summary of important aspects of rice genome research, the editorial discussing the agreements that govern accessibility to some of the data, and letters to the editor that call for continuing with the sequencing of rice to the point of a final, complete sequence. In the Perspectives, Bennetzen, as well as Ronald and Leung, discuss the implications for genomic and agricultural sciences, Cantrell and Reeves discuss the implications that sequencing the rice genome holds for promoting worldwide food security, and Serageldin discusses the interplay between worldwide economic development and food security. News stories discuss the background on what many of the various rice genome-sequencing groups have been doing, the Celera-type agreements that are governing release of the data for the Goff *et al.* paper, and a profile of the Chinese group that authored the Yu *et al.* paper.

Africa. More frequent north winds over the Malawi basin may have resulted from a more southerly migration of the Intertropical Convergence Zone. This interpretation is supported by the Nb-Ti ratio of the sediments. This ratio, which varies roughly in concert with the biogenic silica accumulation rate, is controlled by the input of windblown volcanic ash that originates from a region north of the lake.

Controlling Colloids

Insights into the formation of glasses and crystals by molecules can be gained by studying and manipulating similar processes in colloidal particles, which can be more readily imaged (see the Perspective by Frenkel). Different amorphous states of a material can be observed that have different density, as has been observed for water. Pham *et al.* (p. 104) confirm through a mixture of simulations, theory, and

experiments that two distinct glassy phases of a colloid can form that depend on the nature of the interparticle attractions. At high colloid concentrations, the particles jam against each other and produce a glass based on repulsive interactions. The addition of a noninteracting polymer first causes the glass to melt by causing additional repulsive interactions that open up the structure. Adding more polymer, however, produces a different glassy phase that forms through attractive interactions between the colloidal particles. In binary crystals of large (L) and small (S) colloids, crystal structures can be obtained that are analogous to the atomic structures, such as an LS_2 structure analogous to that of AlB_2 . Working with concentrated systems and a narrow range in the ratio of S:L sizes, Velikov *et al.* (p. 106) create colloidal crystals through a layer-by-layer deposition technique. They not only replicated the LS_2 and LS structures but also managed to obtain the LS_3 structure, which does not have an atomic analog. By building the crystal from two materials, they could selectively remove the small colloidal particles to leave behind a hexagonal, non-closed packed crystal.




On-target results. One easy process.

MagicMark™ Western Standard

Go ahead, take direct aim. With MagicMark™ Western Protein Standard from Invitrogen, you can bypass the maze of steps that conventional standards require. Visualize standard bands and your target protein simultaneously on western blots—and see superior results.

One-step accuracy. MagicMark™ lets you estimate molecular weight directly on your western blot. No aligning. No marking. No extra reagents. Simply load MagicMark™ onto your gel, then electrophorese and perform the immunoassay of your choice. You'll spot standard bands together with your target protein, minimizing errors.



**MagicMark™ detected
simultaneously
with target protein**

Superior results. Capture clear, sharp MagicMark™ bands with the same conjugated antibody and protocol you use to detect your target

protein. Choose chemiluminescent or colorimetric detection methods for your analysis. Because they are not chemically modified, MagicMark™ proteins offer superior sizing accuracy. Eight distinct bands let you easily identify your protein of interest on your blot (see figure).

Land on target. When you need a protein standard for western blots, get the accuracy and convenience that only MagicMark™ can deliver. To learn more, call Invitrogen today or visit us at www.invitrogen.com/magicmark.

Invitrogen
life technologies

Corporate Headquarters:
Invitrogen Corporation
1600 Faraday Avenue
Carlsbad, California 92008 U.S.A.
Tel: 1 760 603 7200
Tel (toll free): 1 800 955 6288
Fax: 1 760 603 7229
Email: tech_service@invitrogen.com

European Headquarters:
Invitrogen Ltd
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Tel: +44 (0) 141 814 6100
Fax: +44 (0) 141 814 6287
Email: eurotech@invitrogen.com


For an office near you go to:
www.invitrogen.com

Birds, Bile, and Breeding

Birds use a variety of physiological and environmental cues to determine when best to lay eggs, but learning also plays a role. Grieco *et al.* (p. 136; see the news story by Withgott) have discovered, in experiments with blue tits in a Netherlands forest, that a female bird's experience of food supply in 1 year can influence the date of egg-laying in the subsequent year. By increasing the availability of food during the nestling period, the birds were duped into laying eggs later in the following year. Female sea lampreys must travel long distances up streams to reach the nests that their male counterparts have already built. The females follow a chemical cue released by the males, but the known fish sex pheromones, which include steroids, only act at short ranges. Li *et al.* (p. 138) show that male lampreys release a bile acid that can signal both their location and reproductive status to females over such large migratory distances.

Cdc13 Joins the Fold

The ends of eukaryotic chromosomes are capped by specialized nucleoprotein complexes called telomeres. The telomeric single-stranded (ss) DNA binding protein used by the budding yeast *Saccharomyces cerevisiae*, Cdc13, is unrelated in amino acid sequence to the telomere end-binding proteins in other organisms. Mitton-Fry *et al.* (p. 145) determined the solution structure of a complex containing the Cdc13 DNA binding domain and telomeric ssDNA and found that despite its distinctive primary sequence, Cdc13 uses the same structural motif for DNA binding—the so-called "OB fold"—as the other proteins.



Tracking the Origins of Chromatin Acetylation

Modulation of chromatin structure plays key role in the regulation of gene expression in eukaryotes, whereas in archaea, gene regulation is thought to be effected by simple bacteria-like repressors. Sir2, which plays an important role in chromatin dynamics in yeast, is highly conserved across prokaryotes, archaea, and eukaryotes. Bell *et al.* (p. 148) show that Sir2 from the archaeon *Sulfolobus solfataricus* modulates (and increases) the DNA binding activity of the chromatin protein Alba by deacetylating two of the lysine residues in its amino terminus. This step in turn facilitates Alba's ability to repress transcription *in vitro*. These results suggest that regulation of chromatin structure by acetylation arose in the common ancestor of the archaea and eukaryotes.

The Outs and Ins of Signaling

Integrins, which are central to cell adhesion, are cell surface molecules comprised of α and β transmembrane subunits. Their extracellular domains bind to external ligands containing the Arg-Gly-Asp tripeptide motif, and this results in the canonical "outside-in" signaling common to many cell surface receptors: Internal signals transmitted to their intracellular domains alter the ligand-binding competency of the external regions. Xiong *et al.* (p. 151) present the structure of the complex between the extracellular portions of the integrin $\alpha V\beta 3$ and a cyclic peptide. The contacts made by the Arg-Gly-Asp residues lead to conformational changes in the relative orientations of the α and β subunits. **X**

Damage Control

Successful repair of injured tissue relies on removal of products generated in the turnover of the extracellular matrix. Teder *et al.* (p. 155) offer direct evidence that CD44, a receptor for the matrix protein hyaluronan, is pivotal in this process. In a model of lung injury that usually resolves within 2 weeks, mice lacking CD44 expression developed fatal pulmonary inflammation caused by persisting inflammatory cells and low-molecular-weight hyaluronan fragments at the damaged site. Normal removal of these agents was associated with transforming growth factor- β activation and could be mostly attributed to CD44 expression on hematopoietically derived cells, rather than parenchymal cells of the lung.

CREDIT: MITTON-FRY ET AL.

APOPTOSIS REVEALED

New Products

- Superoxide Dismutase Assay Kit
- VasoTACS™ *In Situ* Apoptosis Detection Kit

Visit us at The AACR
Booth #1220

TREVIGEN

has taken the guess work out of detecting apoptosis.

We have developed a product line that eliminates the confusion and difficulty in revealing apoptosis.


TREVIGEN...

remains committed to the highest level of product quality and expert technical support.

Contact us about
APOPTOSIS REVEALED



www.trevigen.com
1.800.TREVIGEN



Free trial download at www.phoretix.com

Phoretix 2D is specifically designed for the in-depth analysis of 2D electrophoresis gels. It provides the levels of accuracy, flexibility and control necessary for the most demanding and meticulous scientific research. Incorporating extensive tools for detailed profiling of protein expression changes, it remains the world leader in 2D gel analysis.

Visit the Phoretix web site to download your free trial software.

Nonlinear Dynamics - creators of *phoretix*™
tel: +44 (0)191 230 2121 | email: info@nonlinear.com

Nonlinear USA Inc
toll free: 1-866 GELS USA | tel: 919 474 5177 | email: inc@nonlinear.com


nonlinear
dynamics



DO MORE

The ABI PRISM® 3100 Genetic Analyzer.
Sixteen capillaries, endless possibilities.
Your research may now proceed.

We developed the ABI PRISM® 3100 Genetic Analyzer to give you more of what you need. With 16 capillaries and proven ABI PRISM® chemistries, the 3100 genetic analyzer gives you more throughput capacity, more sequencing and fragment analysis applications, and—best of all—more time. With the 3100 genetic analyzer, there are no gels to pour, load, or track. Just place your samples in the autosampler, set up the software, and start the run. The system will run multiple, unattended runs for up to 24 hours. You'll also get more peace of mind with time-tested Applied Biosystems reliability, and a customer support network that can help you take full advantage of everything the 3100 genetic analyzer has to offer. To find out more, visit www.appliedbiosystems.com/3100, or talk with an Applied Biosystems sales rep at **1.650.638.5800**.

AB Applied Biosystems
Data. Decision. Discovery.

Applied Biosystems is committed to providing the world's leading technology and information for life scientists. Applied Biosystems consists of the Applied Biosystems and Celera Genomics businesses. ABI Prism and its design and Applied Biosystems are registered trademarks and AB (Design) and Applied are trademarks of Applied Biosystems or its subsidiaries in the US and certain other countries. For Research Use Only. Not for use in diagnostic procedures. The ABI Prism® 3100 Genetic Analyzer includes patented technology licensed from Hitachi, Ltd., as a part of a strategic partnership between Applied Biosystems and Hitachi, Ltd., as well as patented technology of Applied Biosystems. © 2001 Applied Biosystems. All rights reserved.

Leaf Tissue to PCR in Under 15 Minutes

Extract-N-Amp™ for Plants

- No Freezing or Grinding
- Increased Throughput
- Increased Specificity

Visit sigma-aldrich.com/XNAPs for more information.

The **Extract-N-Amp Plant PCR Kits** contain reagents to rapidly extract genomic DNA from plant leaves and amplify targets of interest by PCR. This novel extraction system eliminates the need for freezing, mechanical disruption, organic extraction, centrifugation, column purification, filtration or alcohol precipitation. The kit also includes a PCR reaction mix, especially formulated for amplification directly from extract.



Features & Benefits:

- **Fast:** 15 minute prep allows greater speed and throughput
- **Convenient:** Simple procedure with no special equipment required
- **Specific:** Hot Start antibody for highly specific PCR amplification of genomic DNA
- **Flexible:** Tested on a wide variety of plant types

ORDERING INFORMATION

Product Code	Description	Size
XNAP2	Extract-N-Amp Plant PCR Kit	100 extractions/100 amplifications
XNAP	REExtract-N-Amp Plant PCR Kit (contains REDTaq)	100 extractions/100 amplifications

Inquire for bulk and high-throughput custom needs.



Extract-N-Amp Plant PCR Kit with 5 different plant specimens. Genomic DNA was extracted from 0.5 cm leaf disks that were cut using standard paper punch. DNA was extracted using the Extract-N-Amp Plant PCR Kit in less than 15 minutes. All samples were then amplified using the specially formulated Hot Start PCR mix. The products were generated from a 30-cycle duplex reaction containing primers specific to plant chloroplast DNA (upper band) and primers specific for *Cannabis sativa* DNA (lower band).
Data provided by Andy Hopwood, Forensic Science Service, Birmingham, England.

Contact Sigma for more information on
Extract-N-Amp for Plants at sigma-aldrich.com/XNAPs

availability, ordering & tracking
sigma-aldrich.com
scientific tools, products & information




LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI 63178 • USA



Purchase of these products is accompanied by a limited license to use them in the Polymerase Chain Reaction (PCR) process in conjunction with a thermal cycler whose use in the automated performance of the PCR process is covered by the up-front license fee, either by payment to Perkin-Elmer or as purchased, i.e., an authorized thermal cycler. These products are sold under licensing arrangements with F. Hoffmann-LaRoche Ltd., Roche Molecular System, Inc. and the Perkin-Elmer Corporation. Extract-N-Amp and REExtract-N-Amp are trademarks of Sigma-Aldrich Biotechnology, Inc.

We try not to let awards go to our heads.

(Unless of course that makes us smarter.)



Dr. Sergei Lokhov, Dr. Eugene Lukhtanov, Dr. Igor V. Kutyavin and Dr. Michael Reed are four Epoch Biosciences research scientists honored with the 2001 American Chemical Society's Northwest Industrial Innovation Award for their discovery of a new class of DNA probes for genetic analysis.

While others are making due with generic solutions, the scientists at Epoch Biosciences are actively researching and developing smarter technologies that improve the design, hybridization and detection of DNA.

Epoch Biosciences is not just another oligo house. Our team's expertise in nucleic acids has yielded novel oligonucleotide attachment chemistries, modified bases, quenchers, fluorophores, and the Minor Groove Binder (MGB™) that enhance the performance of today's leading genetic analysis platforms in use around the world.


But we aren't stopping there. Epoch Biosciences is now launching products of its own: MGB Eclipse™ Probe Systems for real-time PCR, DNAnchor™ for the production of microarrays, and custom oligonucleotides that incorporate our novel fluorescent dyes, modified bases, and the Eclipse™ Dark Quencher.

It's a new way of thinking. For faster, better results, call 800-562-5544 or visit www.epochbio.com



EPOCH
BIOSCIENCES

Smarter DNA
....



the incredible new precast gel system

Accelerate Your Research — Criterion Gels Have a Greater Capacity than Any Other Mini Gel.

Criterion sets the new standard for precast gel systems, exceeding all previous performance standards for precast gel electrophoresis. Easy-open cassettes and the patented leak-free design ensure unsurpassed quality and reproducibility, while highly visible wells simplify gel loading, further accelerating throughput. Have a question or want to set up a demo? Contact your Bio-Rad representative or sign up at criterion.bio-rad.com.

Visit us on the Web at discover.bio-rad.com.

Call toll free at 1-800-4BIORAD (1-800-424-6723); outside the US, contact your local sales office.

BIO-RAD


genomics

colony picking

microarraying

dna sequencing

dna library
management



The QArray^{mini} is our small bench-top microarrayer. This space saving, high specification arrayer is capable of arraying up to 54 slides from 5 source plates.

We can now offer you an arrayer to meet your own particular technical and budgetary requirements. All our arrayers are manufactured to the highest standards with the state of the art components, this with the full range of on-board accessories and a spotting resolution of 1 micron, gives you an excellent array quality time after time.

GENETIX

www.genetix.com

info@genetix.com

UK: +44 (0)1425 524700

US: +1 617 877 4400

Fax: +44 (0)1425 524700

US Fax: +1 (617) 522 7499

Genetix Ltd, Queensway, New Milton, BH25 5NN, UK

READER AND ADVERTISER RESOURCES

Science Careers:

SCIENCECAREERS.ORG

Science Careers Centre

Science is hosting a career event at the Drug Discovery Technology Europe meeting.

Science Careers Centre
DDT Europe
Stuttgart, Germany
15-18 April 2002

For more information, visit
www.drugdisc.com.

CAREERS & EVENTS IN DRUG DISCOVERY: EUROPE

AD SUPPLEMENT IN THIS ISSUE

This special feature addresses the key issues facing drug discovery in

Europe today and how they will impact job seekers and employers throughout Europe.

Look for it on
page 169.



Product News:

LAB TECHNOLOGY TRENDS

SPECIAL AD SECTION IN THE 12 APRIL ISSUE

Genomics & Bioinformatics

This section will cover the newest tools and techniques developed for the areas of genomics and bioinformatics.

Look for it in the 12 April issue.

WWW.SCIENCEPRODUCTLINK.ORG



Science Product Link

Science's online reader service system helps you find information on products advertised in the pages of *Science*. Search by product category, company name, keyword, or page number.

Find it all at
www.scienceproductlink.org.

Trade Shows & Conferences

LOOK FOR COPIES OF THE 3 MAY ISSUE OF *SCIENCE* AT THESE
UPCOMING TRADE SHOWS AND CONFERENCES:

IMPACT OF GENOMICS ON MEDICINE
13-14 May, Munich, Germany (CHI)

PHARMACOGENOMICS EUROPE
15-16 May, Munich, Germany (CHI)

CHI-Cambridge Healthtech Institute

For *Science* sales representatives' contact information, please see the previous page.



Science is published by The American Association for the Advancement of Science. To subscribe to *Science* and become a member of the AAAS, go to www.aaas.org.

The path to drug discovery has never been brighter.

Packard BioScience is now part of PerkinElmer Life Sciences.

Sometimes an idea appears that's so perfect it's as if it were written in the stars. So it is with the acquisition of Packard BioScience by PerkinElmer Life Sciences. This single strategic decision has created a global solutions provider to help you reach your research goals better, faster and with fewer obstacles than ever before.

Do more. Go further.

Our combined resources give you single source access to:

- *An even broader array of drug discovery tools, including:*
 - **Automated liquid handling and sample preparation**
 - **Advanced microarray products**
 - **Important drug screening and detection technologies**
 - **More extensive service network**
- *An impressive new product pipeline applying leading-edge technologies to advanced instrumentation and reagents tools*
- *An extensive, expert worldwide sales and support organization focused on enhancing your productivity*

As constant as the stars.

We're growing in size and strength, committed to bringing you the products and technologies you need to increase your productivity in research, drug discovery and genetic screening processes. With our worldwide sales and service capabilities, we can deliver products you need today and the service you'll depend on now and in the future.

See the world of difference we can make. Call 1-800-551-2121 or visit our website at www.perkinelmer.com/lifesciences.



PerkinElmer[™]
life sciences.

Worldwide Headquarters: PerkinElmer Life Sciences, Inc., 549 Albany Street, Boston, MA 02118-2512 USA (800) 551-2121

European Headquarters: PerkinElmer Life Sciences, Inc., Imperiastraat 8, BE-1930 Zaventem Belgium

© 2002 PerkinElmer Life Sciences, Inc.

www.perkinelmer.com/lifesciences



GenomeCanada

Genome Canada wishes to recognize the outstanding contribution of all scientists who participated in the International Peer Review Panel of Genome Canada's Competition II - 2002.

Thank you!

Mike Bevan
Norwich UK

Wout Boerjan
Gent Belgium

Philip E. Bourne
La Jolla CA USA

Steven D.M. Brown
Oxfordshire UK

Soren Brunak
Lyngby Denmark

Al Burlingame
San Francisco CA USA

C. Thomas Caskey, Chair
Houston TX USA

Ruth Chadwick
Lancaster UK

David Dyer
Oklahoma City OK USA

Jonathan Eisen
Rockville MD USA

Clair Francomano
Baltimore MD USA

Richard Gibbs
Houston TX USA

Trevor Hawkins
Walnut Creek CA USA

Stephen Kaattari
Gloucester Point VA USA

Mark Lathrop
Evry France

Stephen Lory
Boston MA USA

Francis Martin
Champenoux France

John Mattick
St. Lucia Australia

Joanne Meyer
Cambridge MA USA

Anthony Monaco
Headington Oxford UK

Steve Oliver
Manchester UK

Yoshiyuki Sakaki
Tokyo Japan

Hans-Martin Sass
Washington DC USA

W. Kelley Thomas
Durham New Hampshire USA

Paul B. Thompson
West Lafayette IN USA

James M. Tiedje
East Lansing MI USA

Gert-Jan van Ommen
Leiden The Netherlands

David Ward
New Haven CT USA

Barbara Weber
Philadelphia PA USA

David E. Wemmer
Berkeley CA USA


Rod Wing
Clemson SC USA

James Womack
College Station TX USA


Genome Canada is a not-for-profit corporation dedicated to developing and implementing a national strategy in genomics research for the benefit of all Canadians. Genome Canada has received \$300 million from the federal government to establish five research centres across the country. These Genome Centres (Atlantic, Québec, Ontario, Prairies and BC) work closely with other partners such as provincial governments, the private sector, the financial community, and national and international foundations to ensure that Canada becomes a world leader in genomics research in key selected areas such as agriculture, bioinformatics, environment, fisheries, forestry and health.

For more information, please visit Genome Canada's website at www.genomecanada.ca

CLONED?



The PSQ™ 96



The new PSQ™ 96MA

On the left is one of the most successful systems for SNP, mutation and sequence analysis. The one on the right is its successor. Identical? On the outside, yes.

On the inside, expanded capabilities mean that the new PSQ™ 96MA offers many new possibilities. Clone checking is just one of them.

You still get a simple to use solution, delivering unbeatable performance in terms of speed, accuracy and reliability. Only now, the possibilities are even greater:


<i>SNP analysis</i>	<i>bacterial typing</i>	<i>insertions/deletions</i>	<i>multiplexing</i>
<i>clone checking</i>	<i>forensic applications</i>	<i>allele quantification</i>	<i>genotyping</i>
<i>phage display</i>	<i>mtDNA analysis</i>	<i>tri/tetra-allelic analysis</i>	<i>haplotyping</i>
<i>viral typing</i>	<i>multiple SNPs</i>	<i>mutation analysis</i>	

To find out how Pyrosequencing can shed new light on your work, visit www.pyrosequencing.com



PYROSEQUENCING

A new light on DNA



Performance Enhancers

Power your cell culture system with the right nourishment.

Now you can minimize adaptation time, improve control, increase yield and simplify scale-up with the most trusted name in cell culture.

No matter what your cell type.

Hybridomas. CHO cells. HEK (293) cells. Insect cells. Neurons. Blood and bone marrow cells. Keratinocytes. ES cells. Hepatocytes. Endothelial cells. Mammalian cells for virus production.

You name it. We can make a serum-free, protein-free or chemically-defined GIBCO™ medium optimized for its culture.

What's more, we make it in a format to suit your system: dry powder, 1X liquid, concentrates, and new Advanced Granulation Technology™.

Use GIBCO™ specialty and serum-free nutrient media.

And pump up performance.

www.invitrogen.com



Better Cell Culture Through Innovation




Continuing to set the standard in nucleic acid markers.

For the highest quality, consistency, and reliability in nucleic acid markers, go with Invitrogen. We provide the Gibco BRL® DNA and RNA markers you've come to rely on for superior results.

Reach new heights. Invitrogen's selection of high-performance standards offers unsurpassed convenience and clarity. Markers ranging from 10 bp to 40 kb give you flexibility to size RNA as well as single-stranded, double-stranded and supercoiled DNA easily. Even increments and orientation bands remove obstacles for clear size estimations. Convenient Ready-Load™ formats store at room temperature and eliminate time-consuming

Ready-Load™ Formats



marker preparation steps, helping you to finish your analysis in record-breaking time.

Winning combination. Partner with Invitrogen and achieve your goals with confidence. To find out more about the entire line of nucleic acid markers, call Invitrogen today or visit us at www.invitrogen.com.

Invitrogen
life technologies

Corporate Headquarters:
Invitrogen Corporation
1600 Faraday Avenue
Carlsbad, California 92008 U.S.A.
Tel: 1 760 603 7200
Tel (Toll Free): 1 800 955 6288
Fax: 1 760 603 7229
Email: tech_service@invitrogen.com


European Headquarters:
Invitrogen Ltd
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Tel: +44 (0) 141 814 6100
Fax: +44 (0) 141 814 6287
Email: eurotech@invitrogen.com

For an office near you go to:
www.invitrogen.com

FLUOstar[★] & POLARstar[★] OPTIMA

*Multifunctional Microplate Reader
with Reagent Injectors*

- *Unique resolution with up to 50 readings per second*
- *Kinetic windows: collect the measurement points where it is important*
- *Up to 4 injection activities per well*
- *High-performance luminescence sensitivity and wider dynamic range (up to 2,000,000 counts per second)*



BMG Labtechnologies GmbH, Hanns-Martin-Schleyer-Str. 10, D-77656 Offenburg/Germany
Phone: +49 781 96968-0, Fax: +49 781 96968-67, germany@bmglabtech.com

AUSTRALIA
BMG Labtechnologies Pty. Ltd.
Phone: +61 3 59734744
Fax: +61 3 59734711
australia@bmglabtech.com

FRANCE
BMG Labtechnologies SARL
Phone: +33 1 48862020
Fax: +33 1 48864707
france@bmglabtech.com

UNITED KINGDOM
BMG Labtechnologies Ltd.
Phone: +44 1296 336650
Fax: +44 1296 336651
uk.irl@bmglabtech.com

UNITED STATES
BMG Labtechnologies Inc.
Phone: +1 919 8061735
Fax: +1 919 8068526
usa@bmglabtech.com

www.bmglabtech.com

bmG
BMG Labtechnologies



Economist Robert Shiller wonders why^{we} have

such faith in the utter rationality of markets when we ourselves can be kinda, you know, irrational. Bad judgment, lousy information, half-baked strategies—there are times when

money brings out the worst in people. That's why

Prof. Shiller's retirement dollars are invested with a company whose levelheaded thinking stands out in a world where impulse and intuition are bucking intelligence and insight.

Log on for ideas, advice, and results. TIAA-CREF.org or call 1.800.842.2776



*Managing money for people
with other things to think about.™*

RETIREMENT

INSURANCE

MUTUAL FUNDS

COLLEGE SAVINGS

TRUSTS

INVESTMENT MANAGEMENT

Robert Shiller became a participant in 1975. TIAA CREF Individual and Institutional Services, Inc., and Teachers Personal Investors Services, Inc., distribute securities products. ©2002 Teachers Insurance and Annuity Association - College Retirement Equities Fund (TIAA-CREF), New York, NY 10017. A charitable donation was made on behalf of Robert Shiller.

Fueling Innovation



It's in our genes.



Fueling innovation is at the very core of what we do for you, by supplying a wide range of quality products to the academic and industrial biosciences. Like ExoSAP-IT™. Our new Mut proteins. Sequenase™ T7 DNA Polymerase. And many other tools that fuel emerging technologies and applications. At USB, our commitment to innovation — along with our proven performance, integrity and responsiveness — have made us a leader in the supply of biochemicals and molecular biology reagents. Call USB today — Your Partner in Fueling Innovation™.

USB Corporation 26111 Miles Road, Cleveland OH 44128 800.321.9322 www.usbweb.com





2002

MI
Michigan

PFIZER DISCOVERS NEW ELEMENT FOR SUCCESS.

They call it Michigan.

Pfizer Inc. was looking to expand. So they did some research. And more research. The closer they looked, the more they discovered. In Michigan, they found a state that's committing a billion dollars over 20 years to the Life Sciences Corridor.

A billion-dollar commitment, a highly-skilled workforce and one of the nation's top ten life science industries are big advantages. And they're helping companies like Pfizer do big things.

How big? Well, Pfizer's decision to expand in Michigan is the biggest news since Boeing expanded in Chicago. Actually, it's bigger.

As a matter of fact, *Business Expansion* magazine has named the Pfizer investment the largest deal in the nation for 2001.

"This will bring the city of Ann Arbor and the state of Michigan economic growth and new jobs in the research and development sector, strengthening the Michigan Life Sciences corridor initiative", says David Canter, Pfizer senior vice president and director, Ann Arbor Laboratories.


For the full story, call 1.800.946.6829
or visit www.michigan.org.

www.michigan.org



MICHIGAN | GREAT LAKES. GREAT LOCATION.

Discovery. Transformed.



INTRODUCING FLA-5000

With Greater Versatility and Higher Resolution

Fujifilm's versatile new imaging system for protein electrophoresis and general array analysis allows more imaging methodologies than any other system. Greater versatility. Higher resolution. And designed to grow as your imaging needs change. Because the more you see, the more you know.

Largest Scanning Area – 40 x 46cm selectable scanning area allows 2D electrophoresis gel analysis.

Finest Pixel Size – User selectable at 25/50/100 or 200 microns.

Multi-format Scanning – Three standard lasers (SHG blue, SHG green and LD red) plus an optional fourth (LD) allow multi-format scanning with a single tool.

Dual Photo Multipliers – Allows concurrent or dual fluorescent detection. (Optional)

Modular Upgrades – Allows your system to grow as your needs grow.

 **FUJIFILM**
I&I – Imaging & Information

EXPRESSION PROFILING

Introducing Genome-Wide OLIGONUCLEOTIDE Libraries.



Designed by:
compugen

The Difference Is In The Design.

Accurate design is critical when performing high-throughput gene expression profiling. Sigma-Genosys and Compugen introduce three new Oligonucleotide Libraries, designed using Compugen's proprietary LEADS™ technology for transcriptome prediction and manufactured using Sigma-Genosys' patented Abacus™ oligo synthesis platform.

- Comprehensive collections of gene specific oligos
- Broader genome representation of mRNA's and EST's
- Alternative splice variant modeling
- Grouped by Gene Ontology™* assignment

Contact us at 1-877-710-1501 or go to www.labonweb.com/chips to find out more information about these new genome-wide oligo collections for high-throughput formats and applications.

* Gene Ontology is a trademark of Gene Ontology Consortium

Our New Oligonucleotide Libraries

- Human Library
(18,861 oligos)
- Mouse Library
(7,524 oligos)
- Rat Library
(4,854 oligos)

www.genosys.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE

1442 LAKE FRONT CIRCLE • THE WOODLANDS, TX 77380-3600
1-877-710-1501 • FAX: 281-363-2212 • EMAIL: GINFORMATION@SIAL.COM • A Member of the Sigma-Aldrich Family

SIGMA
GENOSYS

RICE

Science Coordinators: Pamela J. Hines and Barbara R. Jasny

Contributors: Gane Ka-Shu Wong and Jun Yu, Beijing Genomics Institute, Beijing 101300, China, Hangzhou Genomics Institute, Hangzhou 310007, China, and University of Washington Genome Center, Seattle, WA 98195, USA; Huanming Yang, Beijing Genomics Institute, Beijing 101300, China, Hangzhou Genomics Institute, Hangzhou 310007, China, Institute of Genetics, Chinese Academy of Sciences, Beijing 100101, China; Stephen A. Goff, Torrey Mesa Research Institute, Syngenta, San Diego, CA, 92121, USA

Illustrations: (Wall Chart) Terry E. Smith; Liwen Ma and Baoxing Qiu, Beijing Genomics Institute, Beijing 101300, China; Fossilized rice from the Hemudu Relics in Zhejiang Province (China), Sheng Zhao, Hangzhou Genomics Institute, Hangzhou 310007, China; Zhukuan Cheng and Jiming Jiang, University of Wisconsin, Madison, WI, USA; Rod A. Wing, Carol Soderlund, Mingsheng Chen, Clemson University Genomics Institute, Clemson, SC 29634, USA; M. Chen *et al.*, *Plant Cell* 14, 537 (2002); Santanu Dasgupta, Monsanto Research Centre, Bangalore, Karnataka, India, and Gerard Barry, Monsanto, St. Louis, MO, USA.; X. Ye, *et al.*, *Science* 287, 303 (2000); W.Y. Song *et al.*, *Science* 270, 1804 (1995); Tien-Hung Lan and Molly Dunn, Torrey Mesa Research Institute, Syngenta, San Diego, CA, 92121, USA; M. D. Gale and K. M. Devos, *Science* 282, 656 (1998); Gladys Lucille Smith, California Academy of Sciences; Grant Heilman. (This Page) Rice Field, Sonia Guimil and Molly Dunn, Torrey Mesa Research Institute; Calligraphy, Yongxiang Lu, Chinese Academy of Sciences.

Reviewer: Chris Somerville, Carnegie Institution, Stanford, CA 94305, USA.

Design and Art Direction: C. Faber Smith; Production Assistance: Debra Morgenegg, Alan Stonebraker; Copyediting: Harry Jach.

©2002 Science, a publication of the American Association for the Advancement of Science.

sponsored by
 Amersham
 Biosciences

民以食為天

TO THE
 PEOPLE,
 FOOD IS
 HEAVEN.

—ancient Chinese
 proverb

A N Y T H I N G I S P O S S I B L E



GIVES YOU MORE

Greater throughput + greater efficiencies = MegaBACE 4000.

You'll get greater throughput with the new MegaBACE™ 4000 than you have ever achieved before. How much greater? Well, it processes 384 samples simultaneously. Put another way, it generates over 2.8 million bases per day. That's four times greater throughput than a 96-capillary instrument—which means you can sequence four times as many genomes in a week.

Receiving four times greater throughput from one system significantly boosts the efficiency of your operations. That means less labor, less space and less cost per base.

Time is money.

Go to www.megabace.com now for a special demonstration of the power of MegaBACE 4000.

Amersham Pharmacia Biotech UK Limited,
Amersham Place, Little Chalfont,
Buckinghamshire, England HP7 9NA

 **Amersham**
Biosciences



**Yes, this is an advertisement from
Amersham Pharmacia Biotech.
But you may have noticed that
something's different.**

We're now Amersham Biosciences. A name a little shorter than our previous one. And a name that embraces our broad range of expertise - systems for the biosciences, from genes to drugs. With it we have a new logo, one that expresses our role as a provider of integrated products and services. But that's all that's changed. We will still be providing the high quality products and services you have come to expect from us. We will continue to have an open and collaborative business approach. And we'll still have excellent science, that we'll continue to use to create new and innovative solutions.


2001 www.wolf-schins.com



Amersham Pharmacia Biotech UK Limited,
Amersham Place, Little Chalfont,
Buckinghamshire, England HP7 9NA



www.amershambiosciences.com

pennies a day...



For about the cost of a week of daily newspaper deliveries, you can have the latest-breaking news in the scientific world delivered to you in **Science** each week. Join **AAAS** today and get a full selection of headline-making member benefits as well as all the news, policy and research you need to keep you on top of science worldwide. Plus, a classified section with the best employment listings available each week. Visit promo.aaas.org/getscience for more information.

AAAS and **Science**—all the news and research you need each week.
(sorry, no crossword puzzle)



promo.aaas.org/getscience



AMERICAN ASSOCIATION FOR THE
ADVANCEMENT OF SCIENCE

Science
Discover the Future