

Now available...

Adenosine Triphosphate (ATP); Amygdalin; Amylase; Animal Lecithin; Ascorbic Acid Oxidase; Bacitracin; BAL; Biotin cryst.; Caffeic Acid; Carotene cryst.; Catalase; Cellulase; Chlorogenic Acid; Chorionic Gonadotropin; Circulatory Hormone; Colchicine; Cytidylic Acid; Cytochrome C; Dehydroascorbic Acid; Dihydroxyacetone; Dihydroxyphenylalanine (DOPA); Dipyriddy; Edestin; Emulsin; Erythritol; N-Ethylpiperidine; Fibroin; Folic Acid; Galacturonic Acid; Gentisic Acid; Girard Reagents P and T; Gliadin; Glucose Pentaacetate; Glucuronic Acid; Glyceraldehyde; Glyceric Acid; Heparin; Hordenine; Hyaluronidase; 2-Hydroxyadipaldehyde; Humulon; Isoascorbic Acid; Isopropylarterenol; Kojic Acid; Kynurenic Acid; Lanthionine; Lipase; Lysozyme; Lyxose; Malononitrile; Maltase; Melezitose; Mesobilirubinogen; Muscle Adenylic Acid; Myanesis; p-Nitrophenylphosphate; Nucleoprotein; Orcinol; Pancreatin; Pangestin; Pantothenyl Alcohol; Penicillinase; Peroxidase; Phenazine; Phenylpyruvic Acid; Phloridzin; Phosphorylase; Piperin; Porphyrindene; Protamines; Pyridoxal; Pyridoxamine; Pyrocatechuic Acid; Pyruvic Aldehyde; Ribonuclease; Saccharic Acid; Salmine; Serine Phosphoric Acid; Spermidine; Spermine; Thioacetic Acid; Thiocytosine; Thyroxine; Trigonellin; Triphenyltetrazolium Chloride; Tripyridyl; Trypsinogen; Tyrosinase; Tyrothricin; Urease; Uricase; Uridine; Vitellin; Xanthosine.

Ask us for others!

DELTA CHEMICAL WORKS

23 West 60th St. New York 23, N.Y.

Telephone PLaza 7-6317

PHOTOVOLT Line-Operated Electronic MULTIPLIER-PHOTOMETER

Mod. 520-M



For the exact measurement
of extremely low light values
down to 1/10,000 microlumen.

Write for Bulletin #360 to

PHOTOVOLT CORP.

95 Madison Ave.

New York 16, N. Y.

*Now . . . get a firm, lasting grip on small objects
. . . with this new **AUTOMATIC PICK-UP TWEEZER***


You'll appreciate the many uses of this truly practical automatic pick-up tweezer. Picks up small or odd-shaped objects difficult to grip otherwise, and small sterilized objects where finger handling is objectionable. It is ideal for holding fractional weights, crystals, small vessels, etc. Its automatic closing gives you a firm, lasting grip **without continued pressure** making it useful for lengthy gripping of an object, such as holding wire in difficult or awkward positions for soldering. You'll easily perform many operations that are difficult with an ordinary tweezer.

36-601

Automatic tweezer
of stainless steel

f.o.b. New York, N. Y.

\$2.25



*can be
steam
sterilized
and boiled.*

The Automatic Pick-up is the size of an automatic pencil and is provided with a pencil clip for pocket storage. Depressing the top cap causes three strong hooked spring steel wires to extend from the tip permitting them to flare out uniformly. Releasing the pressure on the cap causes them to withdraw toward the tip and close together. Any object within their grasp is thus **firmly** held until released by again depressing the cap. The wire fingers are extendable up to $\frac{3}{4}$ inch between tips.



E. MACHLETT & SON

ESTABLISHED 1897

Laboratory APPARATUS • SUPPLIES • CHEMICALS

220 EAST 23rd STREET • NEW YORK 10, N. Y.