

Reg. U. S.

GBI

Pat. Off.

AMINO ACIDS

hydroxyglutamic acid
 hydroxy-L-proline
 D isoleucine
 DL-isoleucine
 L-isoleucine
 DL-lanthionine
 D leucine
 DL-leucine
 L-leucine
 DL leucyl-glycine
 L leucyl-glycine
 DL leucyl-glycyl-glycine
 L leucyl-glycyl-glycine
 D leucyl-L-tyrosine
 L leucyl-L-tyrosine
 DL-lysine Hydrochloride
 L-lysine Hydrochloride
 D-methionine
 DL-methionine
 L-methionine
 DL-norleucine
 DL-norvaline
 DL-ornithine Hydrochloride
 L ornithine dihydrochloride
 D-phenylalanine
 DL-phenylalanine
 L-phenylalanine
 L-proline
 sarcosine
 D-serine
 DL-serine
 L-serine
 Sodium glutamate
 taurine
 DL-threonine
 L-threonine
 DL-thyroxine
 D-tryptophane
 DL-tryptophane
 L-tryptophane
 DL-tyrosine
 L-tyrosine
 DL-valine
 L-valine

See *Science* June 26 issue
for additional listings


A READY REFERENCE
That will Save Time for You

Use this catalog as a "one stop" source of Amino Acids, Vitamins, Carbohydrates, Adenylates, Nucleates, Purines, Pyrimidines, Tetrazolium Salts, Enzymes, Microbiological and Bacteriological Media, Complete Animal Test Diets and Ingredients for investigational use.



WRITE FOR
YOUR COPY
TODAY

GENERAL BIOCHEMICALS, INC.
60 LABORATORY PARK • CHAGRIN FALLS, OHIO

*The Original***VOLUMETRIC
FRACTIONATOR**

Constant volume (1 to 25 c.c.) per sample—not possible with drop counter, as drop size changes. An illuminated glass cylinder is filled to the level of a phototube which is adjustable in height. The phototube actuates a mechanism which opens a glass ball valve and empties the cylinder.

Column and stainless steel volumetric apparatus stationary.

A simple mechanism moves the spun stainless steel test tube rack in a combined rotation and straight line motion to produce a smooth continuous spiral arrangement of 200 tubes.

Special models available.

Please send for further data.

GME

714 Market Place
MADISON 3, WIS.