


8 October 1982 • Vol. 218 • No. 4568

\$2.50

SCIENCE

AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE



You get real RPM only from Beckman: Reliability Plus Memory.

Start with a choice of speeds up to 80,000 rpm—forces to 600,000 g. Add the exclusive Ultra-Smooth frequency-controlled induction drive and a 16-billion revolution warranty.

Then add *true* microprocessor control with the exclusive, convenient Memory-Pac™ capability, and you have Reliability *Plus* Memory.

The Beckman L8 preparative ultracentrifuges deliver extraordinarily smooth, fast acceleration and energy-efficient high performance. You have ease of operation with touch-sensitive controls and Memory-Pac™ solid state modules to store your exact instructions for repeat runs.

From diagnostic indicators to ω^2t , the true microprocessor L8 gives you

informed control during every phase of the run.

Choose from the largest variety of rotors available, plus tubes, accessories and chemicals, all from a single source—Beckman.

You can expect service to be as prompt as it is expert because there are nearly 500 trained Beckman service representatives in Canada and the U.S.A.

Combined with our 35 years of technical and applications experience, it's no wonder we're called the world leader in centrifuge technology.

Get real RPM. Call your Beckman representative or write: Beckman Instruments, Inc., Spinco Division, P.O. Box 10200, Palo Alto, California 94304.



BECKMAN

Circle No. 119 on Readers' Service Card

Monoclonal Antibodies: Quality to match the quantity

To give you choices, we are building the most comprehensive line of monoclonal antibodies and immunological reagents. To give you performance, we are bringing 25 years of quality control experience to bear on every product in the line.

Monoclonal Anti-Mouse

Anti-mouse Thy 1.2 (IgM)
Anti-mouse Thy 1.2, fluorescein labeled
Anti-mouse Thy 1.2, biotin labeled
Anti-mouse Thy 1.1
Anti-mouse Lyt 1.1 (IgG_{2b})
Anti-mouse Lyt 1.2 (IgG_{2b})
Anti-mouse Lyt 2.1 (IgG_{2b})
Anti-mouse Lyt 2.2 (IgM)
Anti-mouse Lyt 3.2 (IgM)
Anti-mouse Lyb 2.1 (IgG₂)
Anti-mouse Lyb 2.3 (IgG₂)
Anti-mouse Ly 5.1 (IgG_{2a})
Anti-mouse Ly 5.2 (IgG_{2a})
Anti-mouse TL (IgG_{2a})
Anti-mouse β_2 -Microglobulin (IgG_{2a})

Monoclonal Anti-Human

Anti-human Ia (IgG_{2b})
Anti-human Lyt 1 (IgG_{2a})
Anti-human Lyt 2 (IgG_{2a})
Anti-human Lyt 3 (IgG_{2b})
"Pan T Cell" antibody
Anti-human IgE (IgG₁)

Monoclonal Rat Anti-Mouse

Rat anti-mouse (IgG_{2a})
Rat anti-mouse (IgG₂)
Rat anti-mouse (IgA)
Rat anti-mouse (IgG_{2b})

Screening Reagents and Systems T and B Cell Typing System

FIA

Sheep F(ab')₂ anti-mouse Ig, Fluorescein Conj.

Protein A Binding Assays

Protein A₁[¹²⁵I]-
Protein A, HRP Conj.

Isotope Release Assays

Ultra Pure Chromium-51
Indium Oxine, [¹¹¹In]-

RIA

Goat anti-mouse Ig₁[¹²⁵I]-
Goat anti-rabbit IgG₁[¹²⁵I]-
Goat F(ab')₂ anti-rabbit IgG, [¹²⁵I]-
Goat anti-human IgG₁[¹²⁵I]-
Sheep anti-mouse Ig₁[¹²⁵I]-
Sheep F(ab')₂ anti-mouse Ig₁[¹²⁵I]-
Rabbit anti-mouse IgG₁[¹²⁵I]-

Avidin

Avidin₁[¹²⁵I]-

ELISA

Available as components or systems.

Anti-HRP
HRP Conj. Goat anti-mouse Ig
HRP Conj. Sheep F(ab')₂ anti-mouse Ig
Alk. Phos. Conj. Goat anti-mouse Ig
Alk. Phos. Conj. Sheep anti-mouse Ig
Alk. Phos. Conj. Sheep F(ab')₂ anti-mouse Ig
 β -Gal. Conj. Goat anti-mouse Ig
 β -Gal. Conj. Sheep F(ab')₂ anti-mouse Ig

Not for use in humans or clinical diagnosis

New England Nuclear, 549 Albany Street, Boston, MA 02118
Call toll free: 800-225-1572, Telex: 94-0996
Mass. and Internat'l: 617-482-9595

Europe: NEN Chemicals GmbH, D-6072, W. Germany
Postfach 401240, Tel. (06103) 803-0, Telex 4-17993 NEN D

NEN Canada: 2453 46th Avenue, Lachine, Que. H8T 3C9
Tel. 514-636-4971, Telex 05-821808 © 1982 NEN

NEN New England Nuclear®
a Du Pont company

Circle No. 131 on Readers' Service Card

LETTERS	DNA Database: <i>E. Jordan</i> and <i>C. Carrico</i> ; <i>W. A. Beyer</i> ; Elephant Grass: <i>L. V. Crowder</i> ; The Inevitability of Cancer: <i>C. H. Fox</i> ; Empiricism or Diversity?: <i>R. B. Glassman</i> ; Rabies Treatment: <i>R. D. Lele</i> ; <i>Index Medicus</i> : Essential Tool: <i>A. D. Nevill</i> ; Engineering Programs in Arizona: <i>R. H. Gallagher</i>	108
EDITORIAL	Engineering Education and National Policy: <i>J. G. Truxal</i> and <i>M. Visich, Jr.</i>	113
ARTICLES	Some Emerging Applications of Lasers: <i>C. P. Christensen</i>	115
	Sedimentation Field Flow Fractionation: Applications: <i>J. J. Kirkland</i> and <i>W. W. Yau</i>	121
	Reexamination of Acoustic Evidence in the Kennedy Assassination: <i>Committee on Ballistic Acoustics, National Research Council</i>	127
NEWS AND COMMENT	Los Alamos: The Winds of Mutiny	134
	The Perils of Clinch River	137
	<i>Briefing</i> : ICSU Accepts China and Psychology; Soviets Clamp Down on Dissident Groups; Shultz Names Kennedy Nonproliferation Chief	139
	Who's Who in Physics and Math	140
RESEARCH NEWS	Planetary Rings Explained and Unexplained	141
	How Thin Is Thin?	143
	Coagulation as a Common Thread in Disease	145
AAAS NEWS	Social Issues Are Themes for Seminars: <i>J. Wrather</i> ; AAAS Summer Fellows; NSF Chautauqua Courses for 1982-1983; Nomination of AAAS Fellows Invited; Major Scientific Meetings	147

BOARD OF DIRECTORS

D. ALLAN BROMLEY Retiring President, Chairman	E. MARGARET BURBIDGE President	ANNA J. HARRISON President-Elect	LAWRENCE BOGORAD EDWARD E. DAVID, JR.	NANCIE L. GONZALEZ DAVID A. HAMBURG
--	-----------------------------------	-------------------------------------	--	--

CHAIRMEN AND SECRETARIES OF AAAS SECTIONS

MATHEMATICS (A) Felix E. Browder Lynn Arthur Steen	PHYSICS (B) Donald N. Langenberg Rolf M. Sinclair	CHEMISTRY (C) Charles G. Overberger William L. Jolly	ASTRONOMY (D) Irwin I. Shapiro Donat G. Wentzel
PSYCHOLOGY (J) Eleanor J. Gibson Bert F. Green	SOCIAL, ECONOMIC, AND POLITICAL SCIENCES (K) Thomas C. Schelling David L. Sills	HISTORY AND PHILOSOPHY OF SCIENCE (L) Erwin N. Hiebert David L. Hull	ENGINEERING (M) Robert W. Dunlap W. Edward Lear
EDUCATION (Q) Elaine W. Ledbetter Roger G. Olstad	DENTISTRY (R) Paul Goldhaber Harold M. Fullmer	PHARMACEUTICAL SCIENCES (S) Louis A. Luzzi Robert A. Wiley	INFORMATION, COMPUTING, AND COMMUNICATION (T) Marilyn C. Bracken Madeline M. Henderson

DIVISIONS

ARCTIC DIVISION	PACIFIC DIVISION	SOUTHWESTERN AND ROCKY MOUNTAIN DIVISION
Arthur M. Pearson President	Herbert Baker President	Klaus D. Timmerhaus President
Gunter E. Weller Executive Secretary	Alan E. Leviton Executive Director	M. Michelle Bal Executive Office

SCIENCE is published weekly on Friday, except the last week in December, by the American Association for the Advancement of Science, 1515 Massachusetts Avenue, NW, Washington, D.C. 20005. Second-class postage (publication No. 484460) paid at Washington, D.C., and at an additional entry. Now combined with **The Science Monthly**. Copyright © 1982 by the American Association for the Advancement of Science. Domestic individual membership and subscription (51 issues): \$48. Domestic institutional subscription (51 issues): \$85. Foreign postage extra: Canada \$24, other (surface mail) \$27, air-surface via Amsterdam \$55. First class, airmail, school-year, and student rates on request. Single copies \$2.50 (\$3 by mail); back issues \$3 (\$3.50 by mail); classroom rates on request. **Change of address:** allow 6 weeks, giving old and new addresses and seven-digit account number. **Postmaster:** Send Form 3579 to *Science*, 1515 Massachusetts Avenue, NW, Washington, D.C. 20005. *Science* is indexed in the *Reader's Guide to Periodical Literature* and in several specialized indexes.

AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE

BOOK REVIEWS	Comets, reviewed by <i>D. K. Yeomans</i> ; Uranus and the Outer Planets, <i>D. Hunten</i> ; Large White Butterfly, <i>M. D. Bowers</i> ; A Companion to Animal Physiology, <i>W. R. Dawson</i> ; Books Received	149
REPORTS	Orthopyroxene Exsolution in Augite: A Two-Step, Diffusion-Transformation Process: <i>A. M. Isaacs</i> and <i>D. R. Peacor</i>	152
	Are Ions Involved in the Gating of Calcium Channels?: <i>Y. Saimi</i> and <i>C. Kung</i>	153
	Bilirubin-Induced Modulation of Cerebral Protein Phosphorylation in Neonate Rabbits in vivo: <i>L. Morphis</i> et al.	156
	Calanoid Copepods, Feeding Currents, and the Role of Gravity: <i>J. R. Strickler</i>	158
	Suppression of Ovulation in the Rat by an Orally Active Antagonist of Luteinizing Hormone-Releasing Hormone: <i>M. V. Nekola</i> et al.	160
	Complete Amino Acid Sequence of Urotensin I, a Hypotensive and Corticotropin-Releasing Neuropeptide from <i>Catostomus</i> : <i>K. Lederis</i> et al.	162
	Synergistic Action of Potassium Chloride and Magnesium Sulfate on Parasitoid Wasp Oviposition: <i>W. C. Nettles, Jr.</i> et al.	164
	Reversal of Morphine Disruption of Maternal Behavior by Concurrent Treatment with the Opiate Antagonist Naloxone: <i>R. S. Bridges</i> and <i>C. T. Grimm</i>	166
	Male Lek Formation and Female Calling in a Population of the Arctiid Moth <i>Estigmene acrea</i> : <i>M. A. Willis</i> and <i>M. C. Birch</i>	168
	Ecological Correlates of Paternal Investment of Urates in a Tropical Cockroach: <i>C. Schal</i> and <i>W. J. Bell</i>	170
	Serial Position Curve in Rats: Role of the Dorsal Hippocampus: <i>R. P. Kesner</i> and <i>J. M. Novak</i>	173
	Postural Asymmetry and Movement Disorder After Unilateral Microinjection of Adrenocorticotropin 1-24 in Rat Brainstem: <i>Y. F. Jacquet</i> and <i>G. M. Abrams</i>	175
	Focal Cortical Seizures Cause Distant Thalamic Lesions: <i>R. C. Collins</i> and <i>J. W. Olney</i>	177
	Discrimination and Imitation of Facial Expressions by Neonates: <i>T. M. Field</i> et al.	179
MEETINGS	Gordon Research Conferences: 1983 Winter Schedule: <i>A. M. Cruickshank</i>	182


R. E. MASSEY E. SAWYER	SHEILA E. WIDNALL HARRIET ZUCKERMAN	WILLIAM T. GOLDEN Treasurer	WILLIAM D. CAREY Executive Officer
AGRICULTURE AND GEOGRAPHY (E) A. Socolow Mas Dutro, Jr.	BIOLOGICAL SCIENCES (G) Carl Gans Walter Chavin	ANTHROPOLOGY (H) John W. Bennett Priscilla Reining	
PHYSICAL SCIENCES (N) A. Tepperman I. Lowenstein	AGRICULTURE (O) Duane Acker Coyt T. Wilson	INDUSTRIAL SCIENCE (P) Ward J. Haas Robert L. Stern	
PSYCHOLOGY (U) Rosenblatt Laser	ATMOSPHERIC AND HYDROSPHERIC Frederic Sanders Glenn R. Hilst	GENERAL (X) Daniel Alpert S. Fred Singer	

COVER

Female calanoid copepod (*Eucalanus attenuatus*) from the Coral Sea, Australia. This herbivorous planktonic microcrustacean lives in a nutritionally dilute environment. Note the sensor arrays on the outstretched first antennae and the elaborate mouthparts. These animals create feeding currents that help them detect algae. See page 158. [Paul Dixon, Australian Institute of Marine Science, Townsville, Australia]

SmithKline

SmithKline Beckman Corporation/World Headquarters



Eye Care. Improved visual health is one likely outcome of the unification of two world leaders in technology—SmithKline and Beckman Instruments—under the name SmithKline Beckman Corporation.

We entered ophthalmic research equipped with extensive chemical technology and unique capabilities in H_2 receptor technology. We are


now exploring the potential for these technologies in a number of diagnostic and therapeutic applications.

Ongoing research is aimed at finding new pharmaceutical agents for controlling intraocular pressure—a principal factor in glaucoma. Above is an eye treated with sodium fluorescein for study of the flow dynamics of aqueous humor.

Converging

Beckman

One Franklin Plaza/Philadelphia, PA 19101



Physiological Monitoring. Precision instruments and associated products by SmithKline Beckman perform a remarkable variety of tasks in health care and research, speeding and simplifying procedures in hospitals and laboratories.

Many of these systems deliver instantaneous assessments of critical physiological parameters and supply indispensable information for

medical diagnosis and therapy.


A Beckman rate nephelometer is the first high-speed clinical system for quantitation of body-fluid proteins *and* for therapeutic drug monitoring. The mobile, programmable system above measures metabolic and respiratory variables.


Look into SmithKline Beckman. And look into the future.

Technologies

**Guinness Book of World Records, please note:
This is probably the largest selection of
plastic laboratory animal cages in the world.**

29 cages that accommodate mice.
22 cages that accommodate rats.
18 cages that accommodate hamsters.
12 cages that accommodate guinea pigs.





Plus: Many types of covers, filter caps, bedding, racks, accessories, watering systems, metabolism cages, etc., etc., etc.

Plus: Our housing systems and accessories adhere to all NAS, NIH, and Public Law 89-544 requirements.

Plus: We manufacture all our cages and covers and virtually everything else listed above. We control the quality of everything we sell.

Please write or call for details.
 Lab Products Inc., 255 West Spring Valley
 Avenue, Maywood, New Jersey 07607
 Phone: 201/843-4600

lab products inc
 a **LabMedic** company

Circle No. 85 on Readers' Service Card

79-1432

Daily nucleotides [^{32}P]: You pick the day.

We're pleased to announce that a new and better production method enables us to make our dNTP, [α - ^{32}P]- and ATP, [γ - ^{32}P]- available from stock every day at catalog specifications.

Now you can order any day for shipment within 24 hours. Or order Thursday for Monday delivery. You'll get the purest and most stable ^{32}P nucleotides available.

M	T	W	T	F
2	3	4	5	6
9	10	11	12	13
16	17	18	19	20
23	24	25	26	27
30	31			

HPLC makes the difference

By introducing preparative HPLC technology into our manufacturing process, we not only improved production efficiency, but purity and stability as well. In fact, extensive testing convinces us that no other manufacturer can match the reproducibility and

reliability of these nucleotides.

Send for our complete ^{32}P nucleotides listing.

Ready when you are

ATP, [γ - ^{32}P]-

>7000Ci/mmol	Tricine
1000-3000Ci/mmol	Tricine
1000-3000Ci/mmol	Aqueous
10-50Ci/mmol	Aqueous
2-10Ci/mmol	Aqueous

dNTP, [α - ^{32}P]-

>600Ci/mmol	Ethanol:Tricine
>600Ci/mmol	Tricine
~3000Ci/mmol	Tricine

Circle No. 134 on Readers' Service Card

Not for use in humans or clinical diagnosis

New England Nuclear
549 Albany Street, Boston, MA 02118
Call toll free: 800-225-1572, Telex: 94-0996
Mass. and Internat'l: 617-482-9595
Europe: NEN Chemicals GmbH, D-6072, W. Germany
Postfach 401240, Tel. (06103) 803-0, Telex 4-17993 NEN D
NEN Canada: 2453 46th Avenue, Lachine, Que H8T 3C9
Tel. 514-636-4971, Telex 05-821808 © 1982 NEN

NEN New England Nuclear®
a Du Pont company

LEADER IN
BIOTECHNOLOGY

VWR presents:



pH meters that measure up.

Only the new Beckman/Altex family of seven pH_i meters have an exclusive labyrinth sealed case and spill-proof front panel that keep harmful vapors and liquids away from sensitive electronics.

But the new Beckman pH_i series of meters offer you more than sealed-in reliability. They represent a significant advancement in pH meter technology thanks to microcomputer instrumentation and easy-to-use controls. They automatically standardize, calibrate, compute, recognize 4, 7, and 10 buffers and compensate readings for temperature (with

optional ATC probe). The meters also feature a large area LCD to clearly display all essential data, system status, errors and weak or defective electrode messages. Add user-serviceability, and computer/printer compatibility and you have the most reliable, cost-effective family of pH_i meters available.

To find out more about Beckman's family of unique pH_i meters, send for your free copy of "The Beckman Handbook of Applied Electrochemistry." This 128-page book will give you a complete analysis of your

analytical requirements, to help you select the right analyzer, electrodes, and accessories for your particular needs. Just call or write one of VWR's 29 sales offices located nationwide today.

VWR delivers



VWR Scientific Inc.
subsidiary of **Univar**
P.O. Box 7900, San Francisco, CA 94120

Science Magazine in cooperation with
Scherago Associates announces
The First Annual Conference on
COMPUTERS IN SCIENCE
December 6-9, 1982
Conrad Hilton, Chicago, Illinois

Chairman: Dr. Dennis Smith, Group Leader, Discovery Research Systems, Lederle Laboratories

Vice Chairman: Dr. Peter Friedland, Research Associate, Dept. of Computer Science, Stanford University


The conference's primary objective is to provide information on how changing computational technologies will influence future scientific research. Sessions will examine state of the art applications of existing technologies to specific areas of research as well as the impact of new developments on scientific research in the decades ahead.

Each conference day is devoted to a specific topic with lectures and presentations providing in-depth coverage of the subject:

- **Tuesday, December 7th — Products of the Technological Revolution: Building Blocks of Future Computer Systems.** A look at the underlying technology of future computer systems, emphasizing hardware developments and their role in providing powerful, distributed systems.
- **Wednesday, December 8th — Computational Systems: Man/Machine Synergism and the Conduct of Scientific Research.** Discussion on building useful research and development systems, focusing on the interaction of scientists with computers.
- **Thursday, December 9th — Scientific Communication and Collaboration: Conducting Research in the New Computational Environment.** An examination of the influence of computers, on how research is conducted, covering scientific collaboration, communication, resource sharing and the sociology of research in a new computational environment.

On Monday, December 6th, the day immediately preceding the main conference three Tutorials will be presented on the subjects of hardware, software, and communication technology. The Tutorials will provide intensive instruction on these subjects to provide a background for the main conference sessions to follow.

In addition, comprehensive exposition of equipment and software used in scientific computer applications will run concurrently with the conference. Company technical personnel will be on hand to assist attendees and to give demonstrations. Contributed papers in the form of posters or demonstrations are invited. Please check registration form to receive Tutorial information and/or Abstract forms where indicated.



Tuesday, December 7th — Morning lectures:

Microelectronics — Dr. Gordon Moore, Chairman of the Board/CEO, Intel Corp.

Personal Computers — Dr. Adele Goldberg, Principal Scientist/Area Manager, Software Concepts Group, Xerox Palo Alto Research Center

Local Area Networks — Dr. Robert M. Metcalfe, Chairman of the Board, 3COM Corp.

Shared Resources — Dr. C. Gordon Bell, Vice President of Engineering, Digital Equipment Corp.

Wednesday, December 8th — Morning Lectures:

Symbols and Software for Science — Prof. Edward A. Feigenbaum, Dept. of Computer Science, Stanford University

Methodology of Programming — Dr. Ira Goldstein, Manager, Application Technology Dept., Computer Research Center, Hewlett-Packard

Expert Systems/Artificial Intelligence — Prof. Bruce G. Buchanan, Dept. of Computer Science, Stanford University

Computer Graphics — Prof. Robert Langridge, Dept. of Pharmaceutical Chemistry, University of California, San Francisco

Thursday, December 9th — Morning Lectures:

General Perspectives — Dr. Ralph E. Gomory, IBM Vice President and Director of Research

Evolution of Computer Networks — Dr. Robert Kahn, Director, Information Processing Techniques, Defense Advanced Research Projects Agency

Scientific Collaborations — Lynn Conway, Research Fellow/Manager, VLSI System Design Area, Xerox Palo Alto Research Center

Sociological Impact of Computing — Prof. J.C.R. Licklider, Laboratory for Computer Science, MIT

CONFERENCE REGISTRATION FORM

☐ Conference registration fee: \$200

☐ Send Information on Tutorials

☐ Send Abstract Forms

Attendance will be limited.

Please reserve _____ space(s). Registration fees must be enclosed. Please make checks payable to "Computers in Science".

☐ Please send exhibit space prospectus.

Name _____

Affiliation _____

Street _____ City _____

State _____ Zip _____ Telephone _____

Return to: Computers in Science; Scherago Associates, Inc.; 1515 Broadway; New York, New York 10036.

**A new generation of microscopes with
advanced optical engineering and ergonomic design:
Olympus BH-2 System Series.**



This could well be a most significant contribution to the science of microscopy and photomicrography.

The Olympus BH-2.

Fulfillment of operational, illumination and brightness control goals.

Unsurpassed totally new LB optics provide vastly superior image quality.

New operational superiority improves viewing comfort and offers a natural eyepoint position.

Low position focusing and stage movement knobs provide natural

hand positioning.

Improved Koehler illumination is developed for ultra-low magnifications (1x to 4x); and standard magnifications (2x to 100x).

Add precision L.E.D. brightness controls and a host of accessories (phase contrast, Nomarski interference contrast, Hoffman modulation contrast, reflected and transmitted light fluorescence, polarized light and more) and you have today's state-of-the-art development

inviting your examination.

Contact the Olympus Authorized Dealer listed in your Yellow Pages, or Precision Instrument Division, Olympus Corporation of America, 4 Nevada Drive, New Hyde Park, NY 11042.

In Canada, W. Carsen Co., Ltd., Ontario

OLYMPUS®
The Science Company

Circle No. 37 on Readers' Service Card

Shopped for la

Look at a few of the new items



New Corning Mega-Pure™ DF system
Deionization/filtration of Type I water. 0.2 micron filter. Digital resistivity meter continuously shows purity levels. Easy to remove cartridges. No housings to clean. Wall or bench mount.

Circle Reader Service Card #351

New Corning X-EL™ electrodes

Premium performance. New AgCl internal with exclusive ion exchange barrier eliminates silver chloride precipitation, clogging and contamination for long, consistent life. Quick, accurate response, low drift. -5°C to 100°C. Measures accurately in low ionic strength solutions.


Choice of six, reference or combination.

Circle Reader Service Card #352

New Corning Cell Wells™ plates

Sterile 24 and 96 wells for hybridoma/monoclonal antibodies, EIA, ELISA, titrations, and more. Lid rings and raised wells prevent contamination. Rank and file markings for easy identification. Uniform bottoms provide accurate microscope scanning.

Circle Reader Service Card #353



FREE X-EL™ ELECTRODE

when you buy a Corning Tall Meter. 10% off any new X-EL™ electrode you purchase separately. Clip and send with literature coupon for details.

SE-10-82-TB-A

FREE PLATE

Clip and send with literature coupon for free Cell Wells™ □ 24-well or □ 96-well flat bottom plate. Check preference. One per customer.

SE-10-82-TB-B

TWO FREE CARTRIDGES

when you buy the new Corning Mega-Pure™ DF system. Clip and send with literature coupon for details.

SE-10-82-TB-C

bware lately?

in our 43rd catalog.





New Pyrex® disposable shorty pipets

New negative graduations increase capacity or eliminate blowout. Crisp, easy to read markings. Ideal for use under hoods. Withstand heat and chemicals that plastics can't. 1, 5, 10, and 25 ml sizes.

Circle Reader Service Card #354

New 1 and 2 liter Pyrex® funnels. Rotaflo® stopcock won't freeze or leak. You can see the seal set. Only glass and PTFE touch your solutions. The Pyrex® stopper is strong and light. Stands on its head, keeps your hands free.

Circle Reader Service Card #355



Circle Reader Service Card #356

Shop The Book

Send today for details on how to receive *The Lab Book*, the new 43rd edition Corning catalog, with over 4,000 trusted tools for your lab.

Ask for literature on:

- | | |
|---|---|
| <input type="checkbox"/> The Lab Book | <input type="checkbox"/> New Cell Wells™ plates |
| <input type="checkbox"/> New Corning X-EL™ electrodes | <input type="checkbox"/> New Pyrex® shorty pipets |
| <input type="checkbox"/> New Pyrex® funnels | <input type="checkbox"/> New Corning Mega-Pure™ DF system |

Attach coupon(s) to your letterhead with name and address.

Corning Glass Works, P.O. Box 1050, Elmira, NY 14902.

GIANT SWEEPSTAKES!

Win a rare First Edition of one of the giant names in science—Darwin's *On The Origin of Species*, Faraday's *Experimental Researches in Chemistry and Physics*, and more! Full details with catalog. No purchase necessary. Sweepstakes ends 4/11/83.

**FREE
PACK**

Clip and send with literature coupon for free PYREX shorty pipets. One pack per customer.

SE-10-82-TB-D

The Most Trusted Tools of Science

CORNING

SE-10-82-TB-E

NEW, HPLC & SDS ELECTROPHORESIS

MOLECULAR WEIGHT STANDARD KITS

Product No. 30175


MW (SDS) MARKER PROTEIN KIT

A single vial mixture of five highly purified, colored proteins for the determination of molecular weights by SDS-polyacrylamide disc gel electrophoresis. The kit consists of cytochrome c from horse heart and various cross-linked cytochrome c's. Each vial contains 100 micrograms of each protein and approximately 10 mg of sucrose.

1 kit 45.00
3 kits 120.00

Protein	M.W.
cytochrome c	12,400*
cytochrome c dimer	24,800
cytochrome c trimer	37,200
cytochrome c tetramer	49,600
cytochrome c hexamer	74,400

*Margoliash, E., Smith, E.L., Kreil, G. and Tuppy, H. *Nature* 192, 1125 (1961).



Densitometer scan of the SDS—polyacrylamide disc gel electrophoresis pattern of a MW (SDS) Marker Protein Kit.

Product No. 30180


MW (HPLC) MARKER PROTEIN KIT

A single vial mixture of five highly purified proteins for the determination of molecular weights by gel filtration chromatography using HPLC. Each vial contains a total of approximately 200 micrograms of protein.

1 kit 87.50
5 kits 400.00

Column: TSK-Gel G3000 SW
(0.75 x 50 cm)
Detector: UV 280 nm, 0.08 AUFS
Chart speed: 0.5 cm/min
Mobile phase: 0.1 M K-PO₄ +
0.2 M NaCl (pH 7)
Flow rate: 0.7 ml/min.
Sample: MW (HPLC) Marker
reconstituted to 0.1 ml
in mobile phase.

Peak	Protein	M.W.
1	Glutamate dehydrogenase	290,000
2	Lactate dehydrogenase	140,000
3	Enolase	67,000
4	Adenylate kinase	32,000
5	Cytochrome c	12,400



For research chemical use only. Not for drug, food or human use.

HPLC pattern of a MW (HPLC) Marker Protein Kit.



USB

We're part of the solution.

United States Biochemical Corporation

P.O. Box 22400 • Cleveland, Ohio 44122 • 216-663-0330 • Telex: 980718

Write for your free USB Catalog containing over 5,000 biochemicals.



Hellma—the largest assortment of highest precision glass and quartz cells.
Standard · Flow-through · Constant-temperature
Anaerobic · Special Designs
Also available—ULTRAVIOLET LIGHT SOURCES
Deuterium Lamps · Mercury Vapor Lamps
Hollow Cathode Lamps · Power Supplies

HELLMA
CELLS, INC.

Write for literature
Box 544
Borough Hall Station
Jamaica, New York 11424
Phone (212) 544-9534

Circle No. 102 on Readers' Service Card

ULTRA CLEAN COAL

- less than 1% each ash and sulfur
- assay attractive for research seeking lower limit combustion properties
- ideally suited for preparing:

**COAL WATER SLURRIES
COAL OIL MIXTURES
GASIF./LIQUIF. FEEDSTOCK**



OTISCA
industries limited

PO BOX 127
SALINA STATION
SYRACUSE NY
13208

Circle No. 130 on Readers' Service Card

Catecholamines and Metabolites or other neurochemically significant compounds

WITHIN

Brain Tissue
Adrenal Tissue
Body Fluids
Dosage Forms
Plant Matter
Insect Tissue

WITHOUT

Derivatives
Radioisotopes
Fluorescence
Gas Chromatography
Mass Spectrometry
A Rich Uncle

BAS Analyzers based on liquid chromatography/electrochemistry (LCEC) have replaced radiochemical and fluorescence procedures in all of the above areas. Detailed procedures, sample preparation materials, and support from our own analytical laboratory compliment each complete instrument package. These systems are capable of monitoring many phenolic natural products and drugs. All of the catecholamine metabolites and many tryptophan metabolites can be determined. In addition, the activity of many neurologically important enzymes can be measured. Let us know your specific requirements.


DON'T CIRCLE THE BINGO CARD! Call or write today for a detailed reply!

**DOPAMINE • NOREPINEPHRINE • EPINEPHRINE • VMA • HVA • DOPAC • MHPG • M • NM • 3-MT • 5-HT • 5-HIAA • TH • DBH
COMT • NMT • MORPHINE • PHENOTHIAZINES • TRICYCLIC ANTIDEPRESSANTS • GLUTATHIONE • GABA**



bioanalytical systems

1205 Kent Ave • West Lafayette • IN 47906 USA • (317) 463-2505
JAPAN: Tokai Irika SWITZERLAND: Paul Bucher ITALY: Analytical Control





Protein Transfers with S&S BA85. New uses for the proven standard.

Since it was named as the reference standard in hybridization studies almost 20 years ago, researchers have continually found new uses for S&S BA85 nitrocellulose. Now this versatile medium is used to help investigators transfer and immobilize proteins for identification.

For example, proteins separated using two dimensional electrophoresis can be transferred with BA85 by capillary blotting or by electrophoretic means. Staining reaction with appropriate antibodies and other analytical methods can then be employed for identification.

These exciting new methodologies promise to impact on virtually every area of the life sciences. Researchers are currently analyzing body fluids with the ambitious goal of mapping every protein of human origin. Others are using protein studies in the development of monoclonal antibodies for immunologic work. In basic research, protein analysis helps to identify the translation products

of organisms that have undergone genetic manipulation.

In all these areas from pure research through commercial and eventual clinical applications, BA85 makes a vital contribution. And, as with all S&S products, the company's commitment only begins with the product. Full technical support, topical reference literature and a proven willingness to work side by side with scientists to answer special needs are an integral part of the way S&S does business.

Call or write for our protein transfer methods—capillary, contact diffusion and electroblotting—all of which make exciting new use of BA85, the proven standard.

Schleicher & Schuell

Schleicher & Schuell, Inc.
Keene, New Hampshire 03431 (603) 352-3810

Schleicher & Schüll GmbH, D-3354
Dassel, West Germany

Schleicher & Schüll AG, CH-8714
Feldbach ZH, Switzerland

Schleicher & Schüll Nederland BV
s'-Hertogenbosch, The Netherlands

Shown—Preparation for the electrophoretic transfer of multiple Coomassie blue stained 2-D gel separations of basic human and platelet proteins onto S&S nitrocellulose using DALT electrophoresis system. (Photo courtesy of N.L. Anderson, Molecular Anatomy Program, Division of Biological and Medical Research, Argonne National Laboratory.)

AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE

Science serves its readers as a forum for the presentation and discussion of important issues related to the advancement of science, including the presentation of minority or conflicting points of view, rather than by publishing only material on which a consensus has been reached. Accordingly, all articles published in *Science*—including editorials, news and comment, and book reviews—are signed and reflect the individual views of the authors and not official points of view adopted by the AAAS or the institutions with which the authors are affiliated.

Editorial Board

1982: WILLIAM ESTES, CLEMENT L. MARKERT, JOHN R. PIERCE, BRYANT W. ROSSITER, VERA C. RUBIN, MAXINE F. SINGER, PAUL E. WAGGONER, ALEXANDER ZUCKER

1983: FREDERICK R. BLATTNER, BERNARD F. BURKE, CHARLES L. DRAKE, ARTHUR F. FINDEIS, E. PETER GEIDUSCHEK, GLYNN ISAAC, MILTON RUSSELL, WILLIAM P. SLICHTER, JOHN WOOD

Publisher

WILLIAM D. CAREY

Associate Publisher: ROBERT V. ORMES

Editor

PHILIP H. ABELSON

Editorial Staff

Assistant Managing Editor: JOHN E. RINGLE

Production Editor: ELLEN E. MURPHY

Business Manager: HANS NUSSBAUM

News Editor: BARBARA J. CULLITON

News and Comment: COLIN NORMAN (deputy editor), WILLIAM J. BROAD, LUTHER J. CARTER, CONSTANCE HOLDEN, ELIOT MARSHALL, R. JEFFREY SMITH, MARJORIE SUN, JOHN WALSH

European Correspondent: DAVID DICKSON

Research News: ROGER LEWIN (deputy editor), RICHARD A. KERR, GINA KOLATA, JEAN L. MARX, THOMAS H. MAUGH II, ARTHUR L. ROBINSON, M. MITCHELL WALDROP

Administrative Assistant, News: SCHERRAINE MACK; Editorial Assistant, News: FANNIE GROOM

Senior Editors: ELEANORE BUTZ, MARY DORFMAN, RUTH KULSTAD

Associate Editors: SYLVIA EBERHART, CAITILIN GORDON, LOIS SCHMITT

Assistant Editors: MARTHA COLLINS, STEPHEN KEPPEL, EDITH MEYERS

Book Reviews: KATHERINE LIVINGSTON, Editor; LINDA HEISERMAN, JANET KEGG

Letters: CHRISTINE GILBERT

Copy Editor: ISABELLA BOULDIN

Production: NANCY HARTNAGEL, JOHN BAKER; ROSE LOWERY; HOLLY BISHOP, ELEANOR WARNER; BEVERLY DURHAM, JEAN ROCKWOOD, LEAH RYAN, SHARON RYAN

Covers, Reprints, and Permissions: GRAYCE FINGER, Editor; GERALDINE CRUMP, CORRINE HARRIS

Guide to Scientific Instruments: RICHARD G. SOMMER
Assistants to the Editors: SUSAN ELLIOTT, DIANE HOLLAND

Membership Recruitment: GWENDOLYN HUDDLE

Member and Subscription Records: ANN RAGLAND
EDITORIAL CORRESPONDENCE: 1515 Massachusetts Ave., NW, Washington, D.C. 20005. Area code 202. General Editorial Office, 467-4350; Book Reviews, 467-4367; Guide to Scientific Instruments, 467-4480; News and Comment, 467-4430; Reprints and Permissions, 467-4483; Research News, 467-4321. Cable: Advancesci, Washington. For "Information for Contributors," write to the editorial office or see page xi, *Science*, 25 June 1982.

BUSINESS CORRESPONDENCE: Area Code 202. Membership and Subscriptions: 467-4417.

Advertising Representatives

Director: EARL J. SCHERAGO

Production Manager: GINA REILLY

Advertising Sales Manager: RICHARD L. CHARLES

Marketing Manager: HERBERT L. BURKLUND

Sales: NEW YORK, N.Y. 10036: Steve Hamburger, 1515 Broadway (212-730-1050); SCOTCH PLAINS, N.J. 07076: C. Richard Callis, 12 Unami Lane (201-889-4873); CHICAGO, ILL. 60611: Jack Ryan, Room 2107, 919 N. Michigan Ave. (312-337-4973); BEVERLY HILLS, CALIF. 90211: Winn Nance, 111 N. La Cienega Blvd. (213-657-2772); DORSET, VT. 05251: Fred W. Dieffenbach, Kent Hill Rd. (802-867-5581).

ADVERTISING CORRESPONDENCE: Tenth floor, 1515 Broadway, New York, N.Y. 10036. Phone: 212-730-1050.

Engineering Education and National Policy

We are witnessing a significant transfer of sources of support from the federal government to private industry. But the objectives of industry may differ radically from our previous national policy. Industry is often primarily interested in quick practical applications, and accordingly focuses resources on schools where the research quality is very high and the existing research parallels industrial interests. We need only note the sequence of industrial commitments to major universities such as Massachusetts Institute of Technology or Stanford to envision a narrow concentration of support for excellent engineering education on a few schools by the end of this decade.

The growth of industrial support for engineering education occurs just after the abolition of the National Science Foundation's (NSF's) Directorate for Science and Engineering Education. Since the mission-oriented federal agencies have responsibilities analogous to those of industry, the NSF change means that we cannot anticipate future programs giving support to average or just-above-average institutions.

Possibly this change in national policy is desirable. It may well be true that, with limited resources, we must concentrate on a few institutions and educate a small number of excellent students to the fullest degree possible: the country's economic and social problems may demand relatively few, truly outstanding engineers in the future. But the policy change may lead to a serious national shortage of solidly educated engineers. There may also be less obvious impacts: second-echelon colleges have frequently been the strongest contributors to improving precollege mathematics and science education, to opening engineering careers to minorities and women, and to contributing technological capability to local and state governments and citizen groups. Moreover, the second-echelon institutions frequently carry on research that fills critical gaps. For example, during the 1950's and 1960's, electrical energy research and education remained in a low priority position at the most prestigious schools, where Department of Defense and NASA interests were strong driving forces. When energy problems became a national priority in the early 1970's, a history of strong emphasis existed only at a few schools and state universities.

If the trend toward concentration of excellence in a few schools continues, we anticipate the emergence by the 1990's of distinct groups of engineering schools: perhaps 25 superb schools, 150 schools offering weak programs, and 75 schools which really should not be offering engineering at all. We also anticipate that there will be a large difference between the superb schools and the weaker ones. What are the conditions for future excellence?

Three factors guide us in predicting which engineering colleges will be outstanding: (i) Quality and visibility of present research. These criteria depend on the level of research activity, available facilities, quality of graduate students, and, especially, quality of faculty. To some extent, areas of research are important simply because of strong national interest. (ii) Sources of support. Certain institutions have unusually strong sources of support. Some are located in cities with exceptional civic pride and have developed rapidly with help from local industrial and financial communities. Others are able to capitalize on remarkable state political commitments to the importance of quality engineering education. (iii) Leadership. Occasionally an engineering college has such exceptional leadership that it is rapidly improving research quality and obtaining long-term sources of support. The leadership understands national priorities, has contacts with industrial and political decision-makers, and is able to allocate resources.

There has been no open discussion of the existence, let alone the desirability, of the changes that are occurring. Are the national trends reversible? Certainly, but only if there is open discussion and formulation of an overt national policy on engineering education.—J. G. TRUXAL and M. VISICH, JR., *College of Engineering and Applied Sciences, State University of New York, Stony Brook 11794*

Image analysis with a human touch

Extract information at a touch

The Zeiss Videoplan is a compact system combining the capabilities of a powerful microcomputer with the discrimination faculties of the operator's own eyes and intelligence. Analyze directly through a microscope, or from photographs, negative material, projections, drawings, X-rays, even directly from video cameras mounted on electron or light microscopes.

Display information at a touch

A single touch provides direct access to all instructions and subroutines of several programs, including Basic Measurement and Evaluation, High-level Statistics, Distribution, and Stereology. A 64K computer memory coupled with dual floppy disc drives is designed for maximum storage, speed, and flexibility.

Nationwide service.

The great name in optics

ZEISS

West Germany

Carl Zeiss, Inc., One Zeiss Drive, Thornwood, NY 10594 (914) 747-1800. Branches: Atlanta, Boston, Chicago, Houston, Los Angeles, San Francisco, Washington, D.C. In Canada: 45 Valleybrook Drive, Don Mills, Ontario, M3B 2S6. Or call (416) 449-4660.



For literature circle reader service number 91
For a demonstration circle reader service number 92