


AMERICAN
ASSOCIATION FOR THE
ADVANCEMENT OF
SCIENCE

SCIENCE

21 JUNE 1991

\$6.00

VOL. 252 ■ PAGES 1585-1756



Frontiers in
Biotechnology

AS LIGHT AS A FEATHER



A Lightweight, Compact, Power Supply for High Voltage Vertical Gel Electrophoresis

The Stratagene FeatherVolt™ 2000 power supply is a compact, lightweight power supply designed specifically for high voltage vertical electrophoresis of sequencing gels. The power supply may be operated in any of three modes: constant power, constant voltage, or constant current, with automatic crossover among all three.

When Bench Space Is At A Premium-- Weighing in at 2.5 Kg., The FeatherVolt™ 2000 Power Supply offers an alternative to bulky, heavyweight high voltage power supplies and is easily portable within your lab setting. External dimensions are 20.3cm W X 24.8cm D X 6.4cm H.

Designed With Your Safety In Mind-- The FeatherVolt™ 2000 incorporates state of the art short circuit detection technology to ensure complete safety at your workstation.

The FeatherVolt 500™ is also available and has been designed for: Laemmli (SDS), O'Farrell 2nd dimension (SDS) and Weber Osborn (native continuous) protein gels; vertical and horizontal DNA and RNA gels, electro-elution, and pulsed field electrophoresis.



Figure Legend : Autoradiograph of sequencing gel using the BaseAce Sequencer and the FeatherVolt 2000 power supply. 1.0 µg M13 ssDNA was labeled with S-35 ATP and run for 3 hrs at 65 watts on a 6% acrylamide (7M Urea, 1X TBE) gel.



Customized Power Supplies are also available. Please contact Stratagene for more information on Volume Pricing.

Catalog numbers:

FeatherVolt™ 2000 (110V)	# 400610
FeatherVolt™ 2000 (220V)	# 400612
FeatherVolt™ 500 (110V)	# 400620
FeatherVolt™ 500 (220V)	# 400622
BaseAce™ Sequencer	# 400120

For more information on our complete line of products contact Stratagene at 1-800-548-1113.



Corporate Headquarters:
11099 North Torrey Pines Rd.
La Jolla, CA 92037
Ordering: 800-424-5444
Tech. Services: 800-548-1113
FAX: 619-535-5430
TELEX: 9103809841

Stratagene GmbH
Postfach 105466
D-6900 Heidelberg
(Federal Republic of Germany)
Telephone: (06221) 40 06 34
Telefax: (06221) 40 06 39
Technical Services: (06221) 40 60 37

Stratagene Ltd.
Cambridge Innovation Center
Cambridge Science Park
Milton Road
Cambridge CB4 4GF
Telephone: (0223) 420955
Telefax: (0223) 420234
Telex: (81417) INNEN G

Circle No. 253 on Readers' Service Card

Primer Specificity Breakthrough.

Perfect Match™ Polymerase Enhancer

Dramatically Reduces False Priming • Enhances Specificity of PCR Amplification

Perfect Match™ polymerase enhancer* is an additive formulated to greatly increase the specificity of primer extension reactions (1,2). Primer-template complexes that contain mismatched nucleotides near the 3' terminus of the primer are destabilized by Perfect Match polymerase enhancer (3). However, perfect and near perfect template-primer matches or primers with regions of non-homology at the 5' end are not affected. Thus, the frequency of false priming events is dramatically decreased. Perfect Match polymerase enhancer may be used to minimize false priming events in such applications as cDNA synthesis and primer extension reactions including those involving PCR amplification (4).

1. Nielson, K. and Mathur, E.J. (1990) Strategies 3:17-19.

2. Nielson, K. and Mathur, E.J. Manuscript in preparation.

3. Nielson, K. and Mathur E.J. (1989) U.S. patents filed.

4. Mullis, K.B., and Faloona, F.A. (1987) Meth. Enzymol. 155:335-350.

Figure 1 shows two examples of *in vitro* amplification reactions that are significantly enhanced by the addition of Perfect Match polymerase enhancer to the polymerase preparation. Note that in lanes 1 and 2, the desired PCR product cannot be detected unless Perfect Match polymerase enhancer is added to the amplification reaction.

In lanes 3 and 4, Perfect Match polymerase enhancer not only increases the intensity of the desired amplification products, but dramatically reduces the background artifacts generated by non-specific priming events.

Size	Catalog #
100 units	600129
200 units	600130




STRATAGENE®

Corporate Headquarters:
11099 North Torrey Pines Road
La Jolla, CA 92037
Ordering: 800-424-5444
Tech. Services: 800-548-1113
FAX: 619-535-5430
TELEX: 9103809841

Stratagene GmbH
Postfach 105466
D-6900 Heidelberg
(Federal Republic of Germany)
Telephone: (06221) 40 06 34
Telefax: (06221) 40 06 39

Stratagene Ltd.
Cambridge Innovation Centre
Cambridge Science Park
Milton Road
Cambridge CB4 4GF
Telephone: (0223) 420955
Telefax: (0223) 420234
Telex: 81417 INNEN G

Circle No. 254 on Readers' Service Card



*Patents Pending

†GeneAmp is a trademark of Perkin-Elmer Cetus

1591 This Week in *Science*

Editorial

1593 Frontiers in Biotechnology

Letters

1599 Citation Counting: T. C. MORRILL; R. M. BATEMAN ■ Satellite Oceanography in the 1980s: R. E. CHENEY AND B. C. DOUGLAS ■ Regulation of Biotechnology: H. I. MILLER, R. H. BURRIS, A. K. VIDAVER ■ The Trabi: Not a Problem: K. KONOPKA ■ *Journal of Biological Chemistry* and Protein Crystallization Papers: H. TABOR ■ Antinoise and Energy Expenditure: J. T. DURKIN

ScienceScope

1603 CEBAF schedule uncertainty; feds formalizing DNA fingerprinting rules; etc.

News & Comment

1604 A "Chilly Climate" for Women?
1606 Scientists Get Mad at OSI: Richards Panel Out of the Loop? ■ Did Imanishi-Kari Get a Fair "Trial"?
1608 French Ban Immunotherapy Treatment Steps Toward a Cooler Greenhouse
1609 To Stop Kuwait's Fires, First Clear the Mines
1610 Ad Hoc Team Revives SSC Competition
1611 Gene Mapping Japan's Number One Crop
1612 *Briefings*: Back to the Future ■ A Tasty Tufted Tuber ■ Women's Health: A World Crisis ■ Researcher's Legal Battle Ends ■ The AIDS-Seal Connection ■ A Nuclear Cure for Nuclear Waste ■ When Do You Send Blue Roses?

Research News

1614 A Genetic Survey of Vanishing Peoples ■ Scientific Split Over Sampling Strategy
1616 Yanomami People Threatened
1617 A Core-Mantle Link?
1618 Gambling on a Shortcut to Genome Sequencing
1619 A Well-Rounded Worm: 220 Kilobases and Counting ■ The Education of *C. elegans* ■ Roundworm as Canary?
1620 Can Earth's Internal Heat Drive Ocean Circulation?

Articles

Frontiers in Biotechnology

1643 Recent Advances in the Polymerase Chain Reaction: H. A. ERLICH, D. GELFAND, J. J. SNINSKY
1651 Complementary DNA Sequencing: Expressed Sequence Tags and Human Genome Project: M. D. ADAMS, J. M. KELLEY, J. D. GOCAYNE, M. DUBNICK, M. H. POLYMERPOULOS, H. XIAO, C. R. MERRIL, A. WU, B. OLDE *et al.*
1657 Monoclonal Antibodies in Diagnosis and Therapy: T. A. WALDMANN
1662 Vaccinia Virus: A Tool for Research and Vaccine Development: B. MOSS
1668 Toward a Science of Metabolic Engineering: J. E. BAILEY

■ **SCIENCE (ISSN 0036-8075)** is published weekly on Friday, except the last week in December, by the American Association for the Advancement of Science, 1333 H Street, NW, Washington, DC 20005. Second-class postage (publication No. 484460) paid at Washington, DC, and additional mailing offices. Copyright © 1991 by the American Association for the Advancement of Science. The title SCIENCE is a registered trademark of the AAAS. Domestic individual membership and subscription (51 issues): \$82 (\$50 allocated to subscription). Domestic institutional subscription (51 issues): \$150. Foreign postage extra: Canada \$46, other (surface mail) \$46, air freight \$90. First class, airmail, school-year, and student rates on request. **Change of address:** allow 6 weeks, giving old and new addresses and 11-digit account number. **Postmaster:** Send change of address to *Science*, P.O. Box 2033, Marion, OH 43305-2003. **Single copy sales:** \$6.00 per issue prepaid includes surface postage; Guide to Biotechnology Products and Instruments, \$20. Bulk rates on request. **Authorization to photocopy** material for internal or personal use under circumstances not falling within the fair use provisions of the Copyright Act is granted by AAAS to libraries and other users registered with the Copyright Clearance Center (CCC) Transactional Reporting Service, provided that the base fee of \$1 per copy plus \$0.10 per page is paid directly to CCC, 27 Congress Street, Salem, Massachusetts 01970. The identification code for *Science* is 0036-8075/91 \$1 + .10. *Science* is indexed in the *Reader's Guide to Periodical Literature* and in several specialized indexes.

■ The American Association for the Advancement of Science was founded in 1848 and incorporated in 1874. Its objectives are to further the work of scientists, to facilitate cooperation among them, to foster scientific freedom and responsibility, to improve the effectiveness of science in the promotion of human welfare, to advance education in science, and to increase public understanding and appreciation of the importance and promise of the methods of science in human progress.



COVER Inspired by Leonardo da Vinci's "Proportions of the Human Figure According to Vitruvius," the cover image symbolizes that the creativity of scientists is at the center of technological advances. Scientists' curiosity about biological processes has often led to discoveries that proved useful for biotechnology. One-half of the figure's appendages are robot-like, showing that technology is an extension of its creator. The figure holds a flask of microorganisms, an antibody, and oligonucleotide primers, representing some of the tools of biotechnology discussed in this issue; the fistful of dollars signifies the entrepreneurial aspect of biotechnology. See pages 1643 to 1681. [Illustration by Julie Cherry]

Research Article

- 1675 Network Rigidity and Metabolic Engineering in Metabolite Overproduction: G. STEPHANOPOULOS AND J. J. VALLINO

Reports

- 1682 Class II Aminoacyl Transfer RNA Synthetases: Crystal Structure of Yeast Aspartyl-tRNA Synthetase Complexed with tRNA^{Asp}: M. RUFF, S. KRISHNASWAMY, M. BOEGLIN, A. POTERSZMAN, A. MITSCHLER *et al.*
- 1690 Impacts, Tsunamis, and the Haitian Cretaceous-Tertiary Boundary Layer: F. J.-M. R. MAURASSE AND G. SEN
- 1693 Identification of Widespread Pollution in the Southern Hemisphere Deduced from Satellite Analyses: J. FISHMAN, K. FAKHRUZZAMAN, B. CROS, D. NGANGA
- 1696 Identity Elements for Specific Aminoacylation of Yeast tRNA^{Asp} by Cognate Aspartyl-tRNA Synthetase: J. PÜTZ, J. D. PUGLISI, C. FLORENTZ, R. GIEGÉ
- 1700 Prevention of Xenograft Rejection by Masking Donor HLA Class I Antigens: D. FAUSTMAN AND C. COE
- 1703 Selective Elimination of HIV-1-Infected Cells with an Interleukin-2 Receptor-Specific Cytotoxin: R. W. FINBERG, S. M. WAHL, J. B. ALLEN, G. SOMAN, T. B. STROM, J. R. MURPHY, J. C. NICHOLS
- 1705 Requirement of Heparan Sulfate for bFGF-Mediated Fibroblast Growth and Myoblast Differentiation: A. C. RAPRAEGER, A. KRUFKA, B. B. OLWIN
- 1708 Identification of p53 as a Sequence-Specific DNA-Binding Protein: S. E. KERN, K. W. KINZLER, A. BRUSKIN, D. JAROSZ, P. FRIEDMAN, C. PRIVES, B. VOGELSTEIN
- 1711 Mapping of DNA Instability at the Fragile X to a Trinucleotide Repeat Sequence p(CCG)_n: E. J. KREMER, M. PRITCHARD, M. LYNCH, S. YU, K. HOLMAN, E. BAKER, S. T. WARREN, D. SCHLESSINGER, G. R. SUTHERLAND *et al.*
- 1715 Structural Determinants of Ion Flow Through Recombinant Glutamate Receptor Channels: T. A. VERDOORN, N. BURNASHEV, H. MONYER, P. H. SEEBURG, B. SAKMANN
- 1718 Facilitation of the Induction of Long-Term Potentiation by GABA_B Receptors: D. D. MOTT AND D. V. LEWIS

Book Reviews

- 1727 The Invisible Hand, reviewed by L. G. NEUBERG ■ The Creation of the Anglo-Australian Observatory, R. W. SMITH ■ The Physics of Musical Instruments, D. E. HALL ■ Combustion, O. I. SMITH ■ Books Received

Products & Materials

- 1731 Transgenic Mice ■ Organics Removal and Concentration Kit ■ Configurations Adapter for 3-D Imaging ■ Primer Design Software ■ DNA and Protein Sequence Analysis Software Updated ■ Dot Blotter ■ Visual Comparison of Chromatograms ■ Literature

Board of Directors

Donald N. Langenberg
*Retiring President,
Chairman*

Leon M. Lederman
President

F. Sherwood Rowland
President-elect

Mary Ellen Avery
Francisco J. Ayala
Eugene H. Cota-Robles
Robert A. Frosch
Joseph G. Gavin, Jr.
Florence P. Haseltine
Jean'ne M. Shreeve
Warren M. Washington

William T. Golden
Treasurer

Richard S. Nicholson
Executive Officer

Editorial Board

Charles J. Arntzen
Elizabeth E. Bailey
David Baltimore
William F. Brinkman
E. Margaret Burbidge
Pierre-Gilles de Gennes
Joseph L. Goldstein
Mary L. Good
Harry B. Gray
John J. Hopfield
F. Clark Howell
Paul A. Marks
Yasutomi Nishizuka
Helen M. Ranney
Robert M. Solow
Edward C. Stone
James D. Watson


Board of Reviewing Editors

John Abelson
Frederick W. Alt
Don L. Anderson
Stephen J. Benkovic
Floyd E. Bloom
Henry R. Bourne
James J. Bull
Kathryn Calame
Charles R. Cantor
Ralph J. Cicerone
John M. Coffin
Robert Dorfman
Bruce F. Eldridge
Paul T. Englund
Fredric S. Fay


Douglas T. Fearon
Harry A. Fozzard
Theodore H. Geballe
Roger I. M. Glass
Stephen P. Goff
Corey S. Goodman
Stephen J. Gould
Eric F. Johnson
Stephen M. Kosslyn
Konrad B. Krauskopf
Charles S. Levings III
Richard Losick
Anthony R. Means
Mortimer Mishkin
Roger A. Nicoll
William H. Orme-Johnson III
Yeshayau Pocker

Dennis A. Powers
Erkki Ruoslahti
Thomas W. Schoener
Ronald H. Schwartz
Terrence J. Sejnowski
Thomas A. Steitz
Robert T. N. Tjian
Emil R. Unanue
Geerat J. Vermeij
Bert Vogelstein
Harold Weintraub
Zena Werb
George M. Whitesides
Owen N. Witte
William B. Wood
Keith Yamamoto

Brinkmann: Quality products today with value for tomorrow.



Eppendorf® Micro Centrifuges Quiet, safe, reliable operation



Model 5415C delivers durable and convenient operation for lab or cold-room use. Refrigerated Model 5402 provides precise temperature control to -9°C.

- Variable speed to 14,000 rpm; 16,000 x g
- Covered rotors reduce noise and help control aerosols
- Quick-release interchangeable rotors for 1.5 mL and 500, 400, and 250 µL tubes

New refrigerated model brings the cold room to your benchtop.

Eppendorf is a registered trademark of Eppendorf-Netheler-Hinz GmbH.
Circle No. 158 on Reader's Service Card


Call 800-645-3050 for more information. In New York, call 516-334-7500. Or write Brinkmann Instruments, Inc., One Cantiague Road, P.O. Box 1019, Westbury, NY 11590-0207. (In Canada: 416-675-7911; 50 Galaxy Blvd., Rexdale, Ont. M9W 4Y5.)

New Eppendorf® Autoclavable Pipettes Now calibratable and serviceable

Adjustable and fixed-volume models from 0.5 to 2500 µL. Fully autoclavable for sterilization and decontamination; can be calibrated and serviced in the lab.

- UV and chemical resistant
- Fast, quiet volume adjustment won't snag safety gloves
- Top-view digital display

Special Ultra Micro models for 0.5 to 10 µL.




Eppendorf is a registered trademark of Eppendorf-Netheler-Hinz GmbH.
Circle No. 157 on Reader's Service Card

New Tuttnauer/Brinkmann Autoclaves Convenient tabletop steam sterilization

Seven mechanical and electronic models with four chamber diameters from 7" to 15". No special plumbing or venting required.

- 316 L stainless steel chambers for longer life
- Fast and slow exhaust cycles
- Double-locking door for safety

Convenient, safe, and reliable steam sterilization.



Circle No. 159 on Reader's Service Card

BRINKMANN
Quality products for research and control.

This Week in SCIENCE

The big chill

IT seems to be the rule rather than the exception that, professionally, women astronomers and physicists are treated as second class citizens. At meetings and in the workplace, they are subjected to a "chilly climate;" they are often excluded from informal social events, where, as all scientists are aware, much of the inside information of science is passed along. Flam reports on two recent surveys—one conducted by an astronomer and the other performed under the auspices of the American Astronomical Society—that indicate a high level of discontent among women astronomers and physicists (page 1604). In general the complaint is not one of overt discrimination but of forms of discrimination that are more subtle and more difficult to document. In a companion piece, Selvin describes the legal battles of Jenny Harrison, a mathematician who was denied tenure at the University of California at Berkeley in 1986. She is currently suing the university, charging overt sex discrimination in the math department's decision to deny her tenure (page 1606).

Workhorse virus

IN 1980, the World Health Organization declared that smallpox had been eradicated from the earth; the last naturally occurring case had developed in Somalia in 1977. Successful control of the disease was achieved with a vaccine made from a virus (vaccinia) related to but not the same as the smallpox virus (variola). Even though smallpox vaccines are no longer administered, the vaccinia virus is not being consigned to retirement. In one of the 6 articles in this *Frontiers in Biotechnology* issue (see pages 1643 to 1681 and the editorial on page 1593), Moss describes some of the research projects and clinical trials that are now underway or in the planning stages with this versatile virus (page 1662). One of the most interesting is the proposed use of vaccinia in vaccines for diseases that are unrelated to smallpox. The virus has a

large genome and can easily accommodate a dozen or more average-size genes. Thus vaccinia viruses could serve as live vectors, multiplying in hosts and carrying with them the genes for immunogenic substances of one or more pathogens. Perhaps a cocktail of many recombinant vaccinia viruses might one day be used in a multi-pathogen vaccine. Risks from the vaccine were acceptable for smallpox because smallpox was such a serious disease; for future uses, some modifications of vaccinia viruses may be necessary to further reduce side effects of these viral agents.

Southern Hemisphere atmosphere pollution

BIOMASS burning in Africa and in southeastern Brazil has caused pollution not only locally in the troposphere but throughout the Southern Hemisphere (page 1693). Evidence of widespread pollution was obtained with satellite data and ozonesonde measurements made at Brazzaville, Congo, and at Ascension Island in the Atlantic Ocean. The ozone that is produced over southern Africa during the dry season (August to October) is transported by low-troposphere easterlies to the eastern South Atlantic Ocean. Some ozone then rises into the upper troposphere where westerlies can transport it long distances, even as far as Australia. In addition to ozone, tropospheric abundances of two other gases—carbon monoxide and methane—vary seasonally. Fishman *et al.* note that these three pollutants typically have been associated with human activities in industrialized countries of the Northern Hemisphere; in the Southern Hemisphere their seasonal variabilities are best explained by the seasonal differences in the amount of biomass that gets burned.

Accepting xenografts


IT is not easy to graft tissue between species (xenografts), because tissue that is so unlike self typically elicits a powerful reaction from the immune

system of the graft recipient. One way around this problem is to suppress the immune system of the host, but this leaves the host vulnerable to infections. A different approach is to try to mask the foreignness of the graft such that the immune system is not activated. Faustman and Coe have had a remarkable success with the second approach (page 1700). They report that when human pancreatic islet cells (the cells that were to be grafted) were pretreated with specific antibody fragments, these cells could then be successfully transplanted into mice; survival exceeded 200 days. The antibody fragments were those directed against the cells' surface HLA class I antigens, which are among the target antigens against which the host's immune response is usually directed. This result is surprising and exciting; if it is broadly applicable, it could simplify tissue grafting across many types of transplantation barriers.

Fusion toxin for HIV-1-infected cells

A chimeric molecule has been found to halt the effects of the human immunodeficiency virus type 1 (HIV-1) on lymphoid cells in vitro (page 1703). The fusion toxin consists of the receptor-binding domain of interleukin-2 (IL-2) and portions of the diphtheria toxin. Receptors for IL-2 are expressed on HIV-1-infected cells after the cells are activated. The fusion toxin binds selectively to activated but not quiescent cells and then specifically kills the cells. Finberg *et al.* report that production of viral proteins was blocked by this treatment and thus the proliferation of infectious viruses was inhibited. Because the IL-2 receptors appear on cell surfaces only after the cells have been activated, normal T lymphoid cells are not affected by the treatment. This form of chemotherapy might be valuable not only for the treatment of the early stages of infection with HIV-1 but for other diseases in which the activation of T lymphoid cells also marks the beginning of the disease process.

■ RUTH LEVY GUYER



It always comes out the same.
Superior.

Sure Blot® Blue. The True Blue.

It is easy to see why Sure Blot® Blue for non-isotopic DNA/RNA hybridization and detection consistently yields superior results.

We couple our positively charged nylon membrane, the industry standard, to our matched hybridization and non-isotopic detection reagents to create a single tailor-made product for Southern and Northern hybridization experiments.

The total system approach plus rigorous quality control testing of every lot makes all the following come true: same consistent level of reproducibility, extremely low background noise, and sensitivity down to at least 0.2 picograms of target nucleic acids...even 0.05 picograms is possible. Unmatched sensitivity and consistency, 5 years of experience, and Oncor® service and support make Sure Blot® Blue the True Blue.

Your future holds superior results.

Just call 1-800-77-ONCOR today and find out more about Sure Blot® Blue.

The Molecular Pathology Company

oncor

209 Perry Parkway Gaithersburg, MD 20877

©1991, Oncor, Inc. Oncor & Sure Blot are registered trademarks of Oncor, Inc.

Circle No. 27 on Readers' Service Card

Pure mRNA in Minutes...

...Directly from Small or Large Samples of Cells or Tissue.

FastTrack™ and MicroFastTrack™ set the industry standard in high quality mRNA isolation.

MicroFastTrack™*: 20 Reactions

- Ideal for PCR, Northern and cDNA synthesis
- Isolation from samples ranging in size from 10^3 - 10^6 cells or 10-250mg of tissue.
- Reproducible yields of high quality mRNA.

FastTrack™*: 6 Reactions

- mRNA isolation for Northern, cDNA, library construction, PCR, microinjection, RNA protection studies and *in vitro* translation.
- Isolation from samples ranging in size from 10^7 - 10^8 cells or 0.4-1.0 gram of tissue.
- Fast, efficient recovery of large amounts of polyA+ RNA from a variety of sources.

Both systems offer:

- High yields of intact mRNA with low ribosomal contamination.
- Eliminate the need for total RNA isolation or the use of toxic chemicals.
- The most cost effective means of generating high quality mRNA.
- Consistency, convenience and the fastest isolation time.

For the very best in direct mRNA isolation FastTrack™ and MicroFastTrack™ are the choice of thousands of research labs worldwide. When the quality of your mRNA is important, turn to the original source for purity, reliability and convenience; turn to Invitrogen.

Toll Free 1-800-955-6288




3985 • B Sorrento Valley Blvd. • San Diego, CA 92121
(619) 597-6200 Phone • (619) 597-6201 Fax

 BRITISH BIOTECHNOLOGY LTD, UK - TEL: 44-235529449 • AMS BIOTECHNOLOGY UK LTD, UK - TEL: 44-993822786 • BDH INC., CANADA - TEL: 800-268-0310 • BIO-TRADE, AUSTRIA - TEL: 43-2228284694 • CELBIO, ITALY - TEL: 39-24048646 • FUNAKOSHI PHARMACEUTICALS, JAPAN - TEL: 81-356841622 • ITC BIOTECHNOLOGY GMBH, GERMANY - TEL: 06221-303907 • KEBO LABS AB, SWEDEN - TEL: 46-86213400 • MEDOS COMPANY PTY LTD, AUSTRALIA - TEL: 61-38089077

*patent pending. mRNA model courtesy of BIOSYM

Circle No. 241 on Readers' Service Card



4 Out of 5 Electroporation References Specify the Bio-Rad Gene Pulser® System



Gene Pulser References



Check our references. You'll find that more scientists use the Gene Pulser for bacterial electro-transformation than all other electroporation units combined.* And our track record for electroporation of eukaryotic cells is just as impressive.

Why? Because we've pioneered a practically designed transfection system that gives you the feedback you need to standardize conditions to those of other labs and to reproduce them from day to day. And because we've established




References to All Other Electroporation Units Combined

the most extensive applications support available—we'll even supply you with otherwise unpublished protocols.

No matter what cells type you work with, from *Acholeplasma* to *Yersinia*, alder tree protoplasts to primary human lymphocytes, Bio-Rad will get you up to speed quickly. Request your free copy of our comprehensive applications binder by calling 1-800-4BIORAD.

* The Gene Pulser is cited in 86% of all published bacterial electroporation articles (gathered from on-line data bases).



BIO-RAD

**Chemical
Division**

U.S. (800) 4BIORAD • California • Ph. (415) 232-7000 • Fx. (415) 232-4257; New York • Ph. (516) 756-2575 • Fx. (516) 756-2594; Canada • Ph. (416) 624-0713 • Fx. (416) 624-3019; Australia • Ph. 61-2-805-5000 • Fx. 61-2-805-1920; Austria • Ph. 43/222/82 89 010 • Fx. 43/222/82 85 629; Belgium • Ph. 32/2/91 85 55 11 • Fx. 32/2/91 82 65 54; France • Ph. 33/1/49 60 68 34 • Fx. 33/1/46 71 24 67; Germany • Ph. 49/89/31 88 40 • Fx. 49/89/31 88 41 00; Kowloon • Ph. 852/789/3300 • Fx. 852/789/1257; Italy • Ph. 39/2/213 87 51 • Fx. 39/2/213 90 32; Spain • Ph. 34/1/661 7085 • Fx. 34/1/661 9698; Japan • Ph. 81-3-534-7240 • Fx. 81-3-534-8037; The Netherlands • Ph. 31/8385-40666 • Fx. 31/8385-42216; Switzerland • Ph. 41/1/810 16 77 • Fx. 41/1/810 19 33; England • Ph. 44/442/23 25 52 • Fx. 44/442/59118; New Zealand • Ph. 64/9/443/3099 • Fx. 64/9/443/3097

Circle No. 195 on Readers' Service Card

THE WHOLE WORLD
IS WAITING
FOR YOUR NEXT
DISCOVERY.

(CAN YOU AFFORD TO WAIT FOR YOUR COMPUTER?)


After doing all the research, data acquisition and analysis, you're about to make the biggest breakthrough in your field. But first, you have to wait. And wait. And wait.

What you're waiting for is time on your department's supercomputer - that is, unless you've discovered VAX 6000 and VAX 9000 vector systems.

VAX vectors can prove to be a major boost for you because they perform analyses up to 400 times faster than a general processing system. And they can be added to an existing VAX 6000 or VAX 9000, or purchased as a complete system. Giving you the supercomputing power you need - up to 500 megaflops - at a price that's easy to justify.

VAX vectors also deliver proven performance in two other areas - accessibility and application availability. You get high-speed access to the VAX vectors from your visualization workstation.

VAX SUPER- COMPUTERS AND NAS OPEN A WORLD OF VECTOR PROCESSING POSSIBILITIES.



You can work with your favorite world-class vectorized applications, including CHARMm, DISCOVER, FIDAP, FLUENT, GAUSSIAN 90, IMSL Libraries, MOPAC, NAG Libraries, PHOENICS, and others. And, what's more, you can work with over 10,000

existing VAX applications.

You get breakthrough performance with VAX vectors since they work with Digital's Network Application Support (NAS), the industry's most open computing environment. An environment that's based not only on industry standards,

but a commitment to letting different applications on different computing systems from different companies all share information and work together.

With all that VAX vectors offer, it's easy to see why Digital has the lead in supercomputing market share after just one year.

Right now, scientists in research centers all around the world are using VAX vectors for their latest and greatest breakthroughs. So call 1-800-332-4636, extension 289 for more information on VAX vectors. Because you can't afford to wait any longer.

digitalTM

THE OPEN ADVANTAGE.

FIDAP IS A TRADEMARK OF FLUID DYNAMICS INTERNATIONAL. FLUENT IS A REGISTERED TRADEMARK OF creare.x. GAUSSIAN 90 IS A TRADEMARK OF GAUSSIAN, INC. IMSL IS A REGISTERED TRADEMARK OF IMSL, INC. NAG IS A TRADEMARK OF NUMERICAL ALGORITHMS GROUP, LTD. PHOENICS IS A TRADEMARK OF CHAM LTD.
Circle No. 208 on Readers' Service Card

OPTICAL FUSION

The storage of important scientific data should never be left to chance. Which is why more researchers are using optical disks for their secure data storage. With the fusion of WORM and erasable functionality into a 1 GByte disk drive, Maximum Storage's new *DUETTE*™ stands alone as the premier optical disk subsystem.

DUETTE makes optical storage easy. Just connect the drive to your computer, load the operating software and you're in business. *DUETTE* works just like a standard disk drive, so you already know how to use it.


Need archival permanence? Insert a WORM disk and write away. Want reusable storage? Use an erasable disk instead. With 1 GByte of capacity, these 5 1/4-inch removable disks are ideal for data-intensive scientific applications.

What's more, *DUETTE* is available for both PC and Sun SPARCstation platforms. Data interchange between these two popular computing environments is fully supported, so disks written by a PC can be read by a SPARCstation and vice versa. And support for other popular computing platforms will be available later this year.

Since 1987, Maximum Storage optical disk drives have been the choice of the research, medical and scientific community. With *DUETTE*, this choice is made even easier. And with our direct sales hot line, you can be assured of prices that will not send your budget into hyperspace.

If you're tired of risking your valuable data to bad sectors, stray magnetic fields and head crashes, then call or write for our informative data pack.

Or better yet, pick up the phone and order your *DUETTE* today.



**AVAILABLE
NOW!**

**To Order, Call Our
Direct Sales
Hot Line**

1-800-THE-MAXX

 **MAXIMUM**
Optical Solutions With a Future

All trademarks registered by their respective companies.

Maximum Storage, Inc. 5025 Centennial Blvd., Colorado Springs, CO 80919 (719) 531-6888

Circle No. 235 on Readers' Service Card

AMPLITAG DNA POLYMERASE

THE MOST VALUABLE PCR ENZYME IN THE WORLD.

Not because AmpliTaq® DNA Polymerase is the only recombinant *Taq* DNA Polymerase and therefore the most consistent in performance. Not because there are hundreds of published PCR protocols for its use. Not because AmpliTaq DNA Polymerase is backed by the unique PCR Performance Guarantee, but...BECAUSE YOUR RESEARCH IS PRICELESS.



AmpliTaq DNA Polymerase.
Backed by our PCR Performance Guarantee.

In the U.S., call PE XPRESS at 1-800-762-4002 to order.

Or call 1-800-762-4001 for technical information.

Outside the U.S., contact your local Perkin-Elmer representative.

PERKIN ELMER CETUS

Europe: Walpersdorf, Germany Tel: 49-8106-381-112 Fax: 49-8106-6697
Canada: Montreal, Canada Tel: 514-737-7575 Fax: 514-737-9726
Far East: Melbourne, Australia Tel: 61-3-560-4566 Fax: 61-3-560-3231
Latin America: Mexico City, Mexico Tel: 52-5-651-7077 Fax: 52-5-593-6223
AmpliTaq is a registered trademark of Cetus Corporation. The PCR process is covered by U.S. patents issued to Cetus Corporation.

Circle No. 218 on Readers' Service Card



We can help you get your absorbable

If your plans to explore absorbable polymers aren't moving fast enough, do what scores of companies—including 8 of the world's top 10 pharmaceutical firms—have already done. Team up with Medisorb Technologies International.

MTI is a partnership of Stolle Medical Products and The DuPont Company.

So we have the resources it takes to get controlled-release drugs and bioabsorbable devices to market fast. In fact, we can assist you every step of the way: from



polymer program off to a faster start.

concept, design, and development, right through to FDA approvals and production. Small wonder that many of our projects are already in advanced clinicals. And that four others are expected to gain FDA approval later this year. To get

a program of your own started, call us at 1-800-772-5091, or send us a fax at 513-489-7267. But please don't delay. Because the only thing less fortunate than missing out on absorbables would be getting in too late.



This may be the last inverted microscope you'll ever buy.

Designed to expand to the newest technologies, the Nikon Diaphot makes certain you're on the cutting edge—this year, every year.

That's why Nikon Diaphot has become the standard in research-level inverted microscopy. Diaphot doesn't just handle what's happening today. Its design is virtually timeless, letting it adapt easily to even the most innovative technologies. In fact, a Nikon Diaphot purchased 10 years ago can be configured to meet any of today's applications in about five minutes.

Expanding easily from

an economical cell culture microscope into integrated workstations for analytical fluorescence.



NT-88 Micromanipulation Workstation



investigations, micromanipulation, calcium ion ratioing and image analysis, Diaphot is a versatile and effective addition to any laboratory.


Which is why top researchers are choosing Nikon Diaphot above all others to interface with quantitative analysis instrumentation from Photon Technology International Inc., Spex Industries and Tracor Northern.

And, in addition to its versatile flexibility, Diaphot is renowned for delivering incredibly crisp, clear images.

That's the Nikon Diaphot. Designed and engineered with one thing in mind: the needs of researchers everywhere...today and tomorrow. And that's something well worth looking into.

For more information, contact your local dealer or Nikon Inc., Instrument Group,

1300 Walt Whitman Road
Melville, New York 11747
(516) 547-8500.



CA++ Ion Fluorescence Ratioing Workstation



Circle No. 154 on Readers' Service Card

Realize Your Potential.

In just 18 months, more than 2,000 scientists worldwide have discovered how to achieve a research performance level at a QC price: the Perkin-Elmer Lambda 2 Series of UV/VIS spectrometers.

Lambda 2 is the answer for all your basic method needs for routine analysis, including wavelength programming, quantitative analysis, and routine scanning. But that's just the beginning. Its new features and applications software make it capable of performing autosampling, routine kinetics, and automatic dissolution testing. And its simple design makes it reliable, easy to use, and affordable.

And with PE XPRESS, all in-stock accessories and supplies are shipped in 24 hours.

More productivity, greater flexibility. Two vital advantages. The Lambda 2 Series.

For more information on our complete product line, contact your local Perkin-Elmer office. For product literature call 1-800-762-4000 or FAX us at 203-762-6037.

PERKIN ELMER

The Perkin-Elmer Corporation
Norwalk, CT 06859-0012

233 Information Only
234 Sales Call



The Lambda 2 is PC compatible, enhancing archiving and data handling.



The Lambda 2 Routine Kinetics System, for determination of enzyme and substrate activity in biochemical and body fluid samples.



The Lambda 2 Automatic Dissolution System, combining intelligent software with accessory versatility.




The Lambda 2 Autosampling System, featuring a new multipurpose auto-sampler holding up to 137 test tubes.



Perkin-Elmer combines experience, expertise, and worldwide resources to provide unmatched products, services, and support for analytical laboratories.

Monoclonal Antibodies to Rat Leukocytes and Over 2000 Other Reagents.



A vital part of your research™

It's vital that we have a model to study the immune system for identification, delineation, and function of the various T and B cells. If your research involves use of the rat model for allografts, rejections, lymphocyte subsets, and other immunoregulation studies, let Bioproducts for Science help. We have the most complete line of Rat Leukocyte cell surface markers available today.

This extensive line of SEROTEC LTD. monoclonal antibodies was developed by researchers, for researchers, with economy and versatility in mind. These antibodies are well characterized, extensively evaluated, and cited in many journal articles. Many are conjugated either to FITC or PE for analysis by either the use of fluorescent microscopy or by flow cytometry.

In addition to Rat Leukocyte antibodies, we have over 2,000 other immunological reagents ready for quick shipment.


Let us join your research team. Call our technical specialists at (317) 894-7536; or write Bioproducts for Science, Inc., P.O. Box 29176, Indianapolis, Indiana 46229.



**BIOPRODUCTS
FOR SCIENCE, INC.**

A Vital Part of Your Research™

No pushover



Only Nalge could introduce a new low-profile, high-flow rate 115-mL filter unit...and roll back prices too!

Just because our old NALGENE™ 115-mL filter unit has been the vacuum filter unit of choice ever since 1965 (and is still available) doesn't mean it's perfect. And you told us so. We listened, and we made improvements.

Greater stability

The low-profile design of the new NALGENE 115-mL filter unit means it won't fall down on the job.

Higher throughput

We've boosted the flow rate for increased filtration efficiency and productivity by eliminating the fibrous support pad and increasing the effective diameter of the membrane. We've also increased spout diameter for faster, easier pouring.


Improved biosafety

New NALGENE 115-mL filter units are certified radiation sterilized, non-cytotoxic and non-pyrogenic...another Nalge exclusive.

Wider membrane selection—including exclusive SFCA

We've kept the standard cellulose nitrate membranes and added low-extractable nylon and low-protein-binding, surfactant-free cellulose acetate, which is a Nalge exclusive.

Environmentally considerate

Like many NALGENE labware products the new unit carries SPI  coding for correct plastics recycling. And it's packaged in easy-to-open, recyclable shipping materials, too.

All your favorite features

Of course, you also told us you liked most of the old unit's features—such as its graduations on both upper chamber and receiver, color-coded membrane identification, and a convenient single pour/vent spout with quick disconnect, multiple-size tubing adaptor.

And we were very careful not to change anything that would require you to requalify the new unit for research protocols or manufacturing SOPs.

All this at '86 prices

One of our major design criteria was cost reduction through better, more efficient manufacturability to actually bring you '90s filtration at '86 prices. This led us to eliminate several expensive, performance-limiting components. Our success means that, now, you actually pay less for more. Available from Authorized NALGENE™ Labware Dealers worldwide.

Ask for your FREE sample

Call or write for a FREE sample of the new NALGENE 115-mL filter unit (0.2- μ m pore size only) in your choice of cellulose nitrate, nylon or surfactant-free cellulose acetate.

Call Nalge Company at 1-800-TRY-A-115.

Nalge Company, A Subsidiary of Sybron Corporation, Department SFCA, P.O. Box 20365, Rochester, N.Y. 14602-0365

Nalge Company


NALGENE® Brand Products

NALGENE™ Filterware


Circle No. 163 on Readers' Service Card

Reason #7 “Cutting-Edge” Technology for Clonal Selection


ACAS 570 Interactive Laser Cytometer... The Cytometer For All Reasons



Two-dimensional, pseudo-color fluorescence image of FDG-labeled, *lacZ*⁺, NIH3T3 cells.



Histogram display of β -galactosidase positive, NIH3T3 cells.



An example of ACAS cell sorting using the “Cookie-Cutter”™ technique to isolate single cells in culture while preserving their normal anchorage dependent growth state.

The ACAS is *the* pioneer instrument for Interactive Laser Cytometry, giving you new techniques for quantitative fluorescence analysis and adherent cell sorting.

It is the only instrument you can use for the automated selection and cloning of adherent cells in culture. With the ACAS you can:

- Screen cells,
- Sort cells in culture based on physical and biochemical measurements,
- Select mutant cells at frequencies of one in a million,
- Clone cells without altering morphology, phenotype, or viability,
- Isolate sub-populations of cells for quantitative fluorescence analysis and,
- Automatically store cell locations for repeated scans of the same cells over time.

Meridian Instruments Inc., introduces a major breakthrough in cloning transfected cells with the ACAS 570. The ACAS was used for the screening, selection, and cloning of transfected mammalian cells carrying the *E. coli lacZ* gene.¹ The product of expression of the *lacZ* gene, β -galactosidase, was measured using the fluorescent probe fluorescein di- β -D-galactopyranoside (FDG).² Selected cells were easily isolated and cloned using the patented ACAS “Cookie-Cutter”™ technique. The ACAS provided quantitative criteria to select *lacZ*⁺ and *lacZ*⁻ cells, and to distinguish between high and low expressers for cloning.

You can use this powerful interactive

technique for the automated sorting and cloning of:

- Mutant cells
- Transfected cells
- Hybridomas
- Rare cells

In addition, you get *all* these reasons that make the ACAS 570 *the* cytometer for cell and molecular biology research:

- #1 – Interactive Laser Cytometry
- #2 – Quantitative Fluorescence Measurements
- #3 – Sensitive Fluorescence Measurements
- #4 – Comprehensive Interactive Software
- #5 – Versatile Integrated Instrumentation
- #6 – Quantitative Confocal Microscopy

Only the ACAS 570 brings the methods of cell sorting and analysis to your cells in culture. These are a few of the many reasons why the ACAS 570 is **The Cytometer For All Reasons**. Call us now at 800/274-7757, or in Michigan 517/349-7200, to discover your reasons.



Meridian Instruments Europe, Inc.
Industriepark — West 75
B-2700 St. Niklaas BELGIUM
Phone: 011-32-3-7801760
Fax: 011-32-3-7781727


Japan Scientific Instrument Co., Ltd.
6-7-5 Higashi Kasai
Edogawa ku Tokyo, 134
Phone: 011-81-3-877-8111
Fax: 011-81-3-877-9236

MERIDIAN[™]
Instruments, Inc.
*Bringing Analytical Power to
the Cell in Culture*

2310 Science Parkway
Okemos, Michigan 48864
Phone: 800-274-7757
517-349-7200
Fax: 517-349-5967

¹Courtesy of James M. Wilson, Howard Hughes Medical Institute, University of Michigan, and Asmina H. Jiwa, Meridian Instruments, Inc.

²Molecular Probes, Inc.



HIGHER
INTELLIGENCE
IN
HIGH MASS
ANALYSIS.

HOW TO ACHIEVE ENHANCED CHARACTERIZATION OF BIOMOLECULES.

The Electrospray System from Finnigan MAT simplifies tedious sequencing processes, and lets you produce accurate and intelligent data in a fraction of the time.

Picomole and femtomole sensitivity in molecular weight determination, coupled with structural elucidation achieved in hours—not days or weeks—makes the Electrospray System a powerful tool.

The Electrospray System combines electrospray ionization (ESI) with our high-performance TSQ[®] 700 mass

spectrometer to provide molecular weight determination of biomolecules, such as peptides and proteins with mass accuracy of 0.01%.

And the innovative Finnigan MAT data processing software extracts meaningful information and presents it in a format tailored for the biochemist, letting you spend more time on science and less time crunching numbers.

To seek higher intelligence in high mass analysis, call a Finnigan MAT office listed below or FAX (408) 433-4823.



A subsidiary of Thermo Instrument Systems, Inc.

California (408) 433-4800 • Georgia (404) 424-7880 • Ohio (513) 891-1255 • Illinois (708) 310-0140 • New Jersey (201) 740-9177 • Maryland (301) 698-9760
Germany 421-54931 • UK 442-233555 • France 1-6941-9800 • Italy 6-601-1742 • Netherlands 838-527266 • Sweden 08-680-0101

Circle No. 186 on Readers' Service Card

- ✓ For microcentrifuge tubes
- ✓ For 96-well plates
- ✓ User-friendly software. Select from 5 languages
- ✓ Large 40 character LCD screen
- ✓ Monitor system variables during experiments
- ✓ Store and print at least 63 programs
- ✓ 4 models in stock
- ✓ Order from your local dealer

TECHNE PROGRAMMABLE HEAT/COOL DRI-BLOCK UNITS

*Temperature cycling
We've made it simple!*



Techne, Dri-Block Reg™ Techne Corp.

For catalog and prices write to Techne Incorporated, 3700 Brunswick Pike, Princeton, NJ 08540 or call toll-free 1-800-225-9243, or fax 1-609-987-8177.

Circle No. 37 on Readers' Service Card

Fast, easy quantitation of fluorescence-based assays

Get fast, easy quantitation of soluble or cell-associated fluorescence in your multiwell plate with the CytoFluor™ 2300 Fluorescence Measurement System. This computer-driven scanner allows you to perform more assays in less time than with traditional techniques.

- This automatic system can read a 96-well plate in under 1 minute. Time course readings are possible. And wavelengths can be mixed and matched for dual or triple-labelling experiments.
- The system can scan 6-, 12-, 24-, 48-, or 96-well standard tissue culture plates, and can be programmed for non-standard plates. Because the system scans plates from the bottom, the cover can be left on to maintain sterility of the sample.
- It can be used in a variety of applications including cytotoxicity testing, fluorescent ELISA, calcium measurement, nucleic acid quantitation, proteinase assays, cell adhesion assays, *in vitro* toxicity assays, gene expression, and intracellular oxidative burst.



For more details and a copy of our descriptive videotape, call toll-free **800-225-1380** (in Massachusetts, dial 617-275-9200).



IF ELECTROPHORESIS IS PART OF YOUR RESEARCH, A TRULY FLEXIBLE POWER SUPPLY CAN SAVE YOU VALUABLE TIME.

Buchler Programmable Power Supplies for electrophoresis bring new flexibility, speed, and convenience to scientific methodology.

Now, you can achieve consistently reproducible results faster by eliminating guesswork from complex separations, focusing, and sequencing.

Each Buchler Programmable Power Supply has its own microprocessor capable of handling up to eight different programs, each with as many as six parameters. Plus, it easily interfaces

with desktop computers.

Available in 1000V, 3000V, and 5000V models. We also offer 250V and 500V units without microprocessors, as well as vertical and horizontal cells, gel dryers, reagents, and accessories.

If electrophoresis is part of your research, make Buchler Instruments part of your lab.

Available through laboratory supply dealers worldwide or call BUCHLERMAN toll free at 800-732-0031.



8811 Prospect • Kansas City, Missouri 64132
Phone 816-333-8811 • FAX 816-363-0130

Circle No. 118 on Readers' Service Card

Enhance your vision.

The precision optics of the Olympus BH-2 Series for diversified microscopy. From distortion-free image quality to reduced eye fatigue, they give you performance that is anything but routine.

Then consider the ergonomic design, the exceptional versatility and operational simplicity, the unquestioned dependability. And it's easy to see why the Olympus BHTU model, for example, has become the preferred laboratory microscope of its type in the market.

The Olympus BH-2 Series microscopes — they'll change the way you look at things.


For a demonstration, call toll-free:
1-800-446-5967. Or write for literature to

Olympus Corporation,
Precision Instrument
Division, 4 Nevada
Drive, Lake Success,
New York 11042-1179.
In Canada: W. Carsten Co., Ltd.,
Ontario.

OLYMPUS
The Image of Quality



For information circle reader service number 155
For a demonstration circle reader service number 156



Catch the next wave of balance technology. And hang 10 μ g.

Surf's Up! And so is productivity! If you're working with samples that require readabilities as fine as 0.1 or 0.01 mg, the RC balance from Sartorius will be welcome news. No more limited sample access or cramped weighing chambers—the tall, cylindrical RC weighing chamber can be rotated a full 360° for practically all-around access—even from above.

No more plastic doors to slide open and shut either—just a touch of a button (or signal from a robot) is all it takes to open the motorized, glass doors wide enough to accommodate any human hand or robotic arm. Quickly, conveniently.

But the RC balance is more than convenient. When you also consider its accuracy and range of on-board electronic features, you'll agree it's positively AWESOME.

It's all in the board (and chip) you choose. Sartorius is the only balance company to use microcomputer technology. With a 40 MHz microchip, the new RC balance is as powerful as most desktop computers! That translates into faster vibration filtration, faster temperature compensation, and faster final stabilization than any other balance on the planet!

By combining advanced electronics and superior ergonomics, the new RC balance eliminates inaccuracies and inconvenience.

For more information or a demonstration of how you can ride the next wave of balance technology for faster, more accurate results, call or circle the number below: Sartorius Instruments, 1430 Waukegan Rd., McGaw Park, IL 60085
Phone: 800•544•3409 FAX: 708•689•2038

sartorius

Smart. State-of-the-Art. Sartorius.

Circle No. 183 on Readers' Service Card

VOLUME 436 • NUMBER 27

64 pages
35 cents

The Boston Digest

FRIDAY, MARCH 15, 1991

WEATHER REPORT

Today: Clouding up, mid 40's
Saturday: Partly Sunny, 45-50
High Tide: 7:24 am 1:51 pm
Full report: Page 54

NEW ENGLAND BIOLABS DISCOVERS *Asc* I, A NEW 8-BASE CUTTER

Wetland Yields New Endonuclease

**GG[▼]CGCGCC
CCGCGC[▲]GG**

BOSTON - An exciting new restriction endonuclease has been found in an *Arthrobacter* species isolated from a Midwestern wetland. The researcher, Rick Morgan of New England Biolabs, reports that the bacteria was isolated from a small pond, dug nine years ago behind his brother's house. A recent family visit gave Morgan the opportunity to collect a water

sample from the mature pond. In Boston, NEB scientists isolated several bacterial strains from the "swamp" water. A number of endonucleases were identified, including isoschizomers for *Sna*B I, *Stu* I and



Scientists confirm global warming trend

By Judy Kramer
DIGEST STAFF

1990 set a record as the warmest year for the earth as a whole according to three out of four research teams monitoring global temperature trends. More significant, however, is the marked warming trend over the past decade which is leading some researchers to conclude that rising temperatures are caused in part by the greenhouse effect.

Infrequent Cutter of Mammalian DNA.

The presence of CpGs in the recognition site makes *Asc* I an infrequent cutter of mammalian DNA.

Produces Ends Compatible with *Bss*H II Cloning Site.

Asc I digestion produces a CGCG 5' extension compatible with *Bss*H II produced ends.

Linkers Available.

Asc I linkers (phosphorylated and non-phosphorylated) are available in three reading frames.



Certified for Digestion of Agarose-Embedded Chromosomal DNA for PFG Electrophoresis.

Asc I digests agarose-embedded chromosomal DNA without any nonspecific degradation.

Stable in Overnight Digestions.

Asc I retains >90% activity at 37°C for 16 hours, making it a preferred choice when cutting agarose-embedded DNAs.

New England Biolabs now offers four 8-Base Cutters.

Asc I (GG[▼]CGCGCC)
Not I (GC[▼]GGCCGC)
Pac I (TTAAT[▼]TAA)
Sfi I (GGCCN4[▼]NGGCC)



NEW ENGLAND

Biolabs INC.

32 Tozer Road, Beverly, MA 01915 USA


800-NEB LABS (US and MA)
Tel. (508) 927-5054 Fax (508) 921-1350

New England Biolabs Ltd., Canada
Tel. (800) 387-1095 (416) 672-3370
Fax (416) 672-3414

New England Biolabs GmbH, Germany
Tel. 49 (06196) 3031 Fax (06196) 83639

Asc I Digestion of Chromosomal DNAs

G + C Content 68% 50% 34%



Agarose-embedded chromosomal DNAs of varying G+C content were digested with *Asc* I and then electrophoresed on a PFG for 29 hrs. Low range PFG marker (lane 1); *Asc* I digested *R. sphaeroides* DNA, 68% G+C (lane 2); *Asc* I digested *E. coli* DNA, 50% G+C (lane 3); *Asc* I digested *S. aureus* DNA 34% G+C (lane 4); and Lambda ladder PFG Marker (lane 5).

DISTRIBUTORS: AUSTRALIA GENESEARCH Tel. (075) 37 5499 / FINLAND, SWEDEN, DENMARK, USSR FINNZYMES (Finland) Tel. (0) 437-5312 / FRANCE OZYME Tel. (1) 30 57 0025 / INDIA BIOTECH INDIA Tel. (542) 311473 / ITALY C.A.M.Bio S.R.L. Tel. (02) 487 06070 / ISRAEL GAMIDOR LTD. Tel. (03) 535-1205 / JAPAN DAIICHI PURE CHEMICALS CO. LTD. Tel. (03) 3272-0671 / KOREA KORAM BIOTECH CORP. Tel. (02) 556-0311 / THE NETHERLANDS WESTBURG B.V. Tel. (033) 95 00 94 / NEW ZEALAND BIOLAB SCIENTIFIC Tel. (09) 418-3039 / NORWAY ING. F. HEIDENREICH Tel. (02) 22 04 11 / PEOPLE'S REPUBLIC OF CHINA CHINA UNITED BIOTECH. CORP. Tel. (1) 256 1627 / PORTUGAL ISODER LDA. Tel. (01) 363-8788 / SPAIN LANDERDIAGNOSTICO S.A. Tel. (01) 594 08 06, (03) 256 9706 / SWITZERLAND FLOW LABORATORIES AG Tel. (061) 634713 / TAIWAN LONG CHAIN INTERNATIONAL CORP. Tel. (02) 565-2605 / UK CP LABORATORIES Tel. (0279) 758200


Circle No. 258 on Readers' Service Card

The resolution to succeed

Epithelium.
Objectiv Achroplan 40X/0.65.



The routine microscope with some "not-so-routine" features.



The Standard 25 microscope is flexible, dependable, surprisingly affordable and anything but routine! Actually, it combines features which give you new-found ease and comfort for your specific applications.

The high performance objectives and the integrated 20 W halogen illumination guarantee high-contrast, bright images. This means reliable assessment of your results and reliable diagnosis. Everything you expect from Carl Zeiss.

A selection of focusing condensers for brightfield, phase contrast 1, 2, 3 and darkfield, and a quintuple nosepiece provide precisely the right choice of equipment for all laboratory techniques.

And the new ergonomic pyramid design ensures that you can work for hours on end without discomfort.

Get closer to the unique story of the Standard 25 — Circle the Reader Service Number. Because you deserve a "not-so-routine" microscope.

Carl Zeiss, Inc.
One Zeiss Drive
Thornwood, NY 10594

The Geometry
for Microscopes

ZEISS

Standard 25
Laboratory Microscope

Circle No. 179 on Readers' Service Card

Carl Zeiss

DNA by Operon.

**Right Price.
Right Now.** **\$3.60**
per base

Now the world's leading supplier of synthetic DNA is also the price leader. Operon's custom DNA is now \$3.60 per base with a \$20.00 set-up fee per sequence, and free domestic delivery. Same outstanding customer service. Same high product quality. New low price. Call for your free researcher kit.

1-800-688-2248

OPERON

OPERON TECHNOLOGIES, INC.

1000 Atlantic Ave., Suite 108 · Alameda CA 94501
Tel. (415) 865-8644 · Fax. (415) 865-5255—NIHBPA 263-00033233

WORLD'S LEADING SUPPLIER OF SYNTHETIC DNA.

Circle No. 15 on Readers' Service Card

TAHITI!

1991



Explore **Tahiti** on board *S/V Wind Song*,
September 27–October 7. Enjoy Papeete
• Raiatea • Bora Bora • & Moorea!
\$2,195 (plus air fare).

(800) 252-4910

Travels with AAAS by **BETCHART**

Recombinant Collagenase

Eliminates Lot To Lot Variation



BioSource International, a leading supplier of recombinant proteins for biomedical research, is pleased to announce the availability of recombinant collagenase.




BIOSOURCE
INTERNATIONAL
950 Flynn Road, Unit A
Camarillo, CA 93120
1-(800) 242-0607

- High specificity
- No contaminating protease activities
- Long shelf life
- High purity
- Bulk quantities available

Developed by Trigen, Inc.
CAUTION: FOR RESEARCH
USE ONLY IN LABORATORY
ANIMALS OR IN VITRO
TESTING. NOT FOR HUMAN
OR DIAGNOSTIC USE.

Circle No. 269 on Readers' Service Card



Advancing the Leading Edge for GMP Peptide Synthesis

**Setting the standard for purity,
consistency and value.**

For over 20 years Bachem California has pioneered peptide technology to consistently offer you products of unequalled quality for your research or clinical applications.*

And we're still setting the standard.

For clinical applications, our unique capabilities in large-scale synthesis and purification give you ultra-pure peptides in quantities ranging from milligrams to kilograms.

We have over 30 drug master files already on file with the FDA and other regulatory agencies throughout the world. Our experience gained from the synthesis of literally hundreds of GMP-grade peptides in our FDA-approved facilities means you can count on Bachem California for the documentation and regulatory support your clinical studies demand.

From GMP-grade peptides with documented therapeutic applications such as Calcitonin, Somatostatin, GRF, CRF and LHRH to custom sequences, only Bachem has the technology and experience to meet your most demanding requirements.

Take advantage of the Bachem edge. Call us today to discuss your individual requirements.



Bachem California

3132 Kashiwa Street • Torrance, CA 90505
Tel: (213) 539-4171 • TWX: 910-347-7320 BACHEM TRNC
Telefax: (213) 530-1571

*GMP-grade products are for manufacturing, processing or repacking in the preparation of a new drug limited by Federal law to investigational use or as a finished raw material in pharmaceutical production. Products are not sterile filtered nor tested for pyrogenicity.

Bachem (UK), Ltd.

A subsidiary of Bachem California
69 High Street • Saffron Walden
Essex CB10 1AA, England
Tel: (0799) 26465 • TLX: 817205 BACHEM G
Telefax: (0799) 26351

Bissendorf Biochemicals, GmbH

A subsidiary of Bachem California
Feodor-Lynen-Strasse 14 • 3000 Hannover 61
Germany
Tel: (0511) 530080 • Telefax: (0511) 5300888

PFANSTIEHL

Octyl Glucosides: This is no soap story

If you're doing protein isolation, you know the choice of detergent is critical. You want complete solubilization without denaturation, as well as ease of removal. Octyl glucosides may be just what you've been looking for.

We now offer both Octyl- β -D-glucopyranoside and 1-S-Octyl- β -D-thiogluco-pyranoside at attractive prices. These mild, nonionic detergents are easily removed by dialysis and offer superior solubilization of even membrane-bound proteins. In addition, the thioglycoside offers increased pH stability and is not degraded by beta-glucosidases. Write or call for more information and pricing.

PFANSTIEHL LABORATORIES, INC.

The source for carbohydrate chemistry

1219 Glen Rock Avenue/Waukegan, IL 60085-0439
Tel.: 1-708/623-0370/Toll Free: 1-800/383-0126
FAX: 708/623-9173/Telex 25-3672 Pfanlab
80-W



MAKING THINGS WORK

Circle No. 170 on Readers' Service Card

Glaxo Inc Research Institute and the University of North Carolina at Chapel Hill present

The Third Annual Frontiers in Chemistry and Medicine Symposium

November 3-5, 1991

To Be Held At: The University of North Carolina at Chapel Hill

"Enzymes and Abzymes"
Stephen J. Benkovic, Penn State

"Towards an Understanding of the
Black Box of Signal Transmission"
Stuart Schreiber, Harvard

"Chemistry of the Ene-diyne Antibiotics"
K. C. Nicolaou, UC-San Diego

"Fiction and Fact"
Carl Djerassi, Stanford

"Catalytic Antibodies"
Richard A. Lerner,
Scripps Research Institute

"A Chemical Approach to the Specific
Recognition of Double Helical DNA"
Peter B. Dervan, Cal Tech

"Transition Metal Complexes—
Pseudoenzymes of the Synthesis
of Bioactive Molecules"
Barry Trost, Stanford


For Registration and additional information call 919-990-5041 (Glaxo/Wanda Holman) or 919-962-2172 (UNC/Becky Smith) or write to: Ms. Becky Smith, Symposium Coordinator/UNC, CB #3290, University of North Carolina, Chapel Hill, 27599-3290.

Make

THE SOURCE

Your
**Source for
STABLE ISOTOPES**

Isotec Inc. offers the scientific community one of the world's largest selections of stable isotopes and labelled compounds



You Get More From

THE SOURCE

- > 99% chemical purity
- High isotopic enrichment
- Competitive prices
- Fast delivery

Plus . . .

metal stable isotopes, helium-3 and other noble gas isotopes, other enriched stable isotopes, multiply-labelled compounds and custom synthesis


Call or write for your **NEW** complete price list today!

ISOTEC INC.
A Matheson, USA Company

3858 Benner Road, Miamisburg, OH 45342
(513) 859-1808 (800) 448-9760
Telex: 288278 Fax: (513) 859-4878

Circle No. 35 on Readers' Service Card

Characterization of the vascularization status with monoclonal antibodies from BMA



BMA BIOMEDICALS AG

Rheinstrasse 28-32
CH-4302 Augst
Switzerland
Tel. (061) 811 62 22
Fax. (061) 811 60 06

Fermentation Control Has Come A Long Way...

The Proteus 2000[®] Benchtop Fermentation System

We can't promise you Chateau LaFite Rothschild but, the Wheaton Proteus 2000 Fermentation Control System can advance the control capability of any fermentor.

Totally automatic with a redundant control structure for maximum reliability and flexibility, the Wheaton Proteus 2000 is so versatile that it can be retrofitted to virtually any bench or production scale fermentor on the market.

For more information call toll-free 1-800-225-1437.




WHEATON

Manufacturers Since 1888


1501 N. Tenth Street
Millville, NJ 08332, USA

Call Toll-Free: 1-800-225-1437, Ext. 2972
TLX: 55-1295 (WHEATON US)
FAX: 1-609-825-1368



For demonstration circle service number 225
For literature circle reader service number 226

HELLMA[®] Dye Laser Cells.



A complete line of stock and special custom cells for all standard and special lasers. Also inquire about our light sources.

For over 20 years...

HELLMA[®] Quality you can trust.

Box 544, Borough Hall Station, Jamaica, N.Y. 11424
Telephone (718) 544-9534 or (718) 544-9166.

Circle No. 216 on Readers' Service Card

PERSONAL BIBLIOGRAPHIC DATABASES

THE FIRST GENERATION:

SCI-MATE • DMS 4 CITE • PRO-CITE
REFERENCE MANAGER • REF-11

THE SECOND GENERATION:

PAPYRUSTM VERSION 6.0

- Manages up to 2 million reference citations. Stores up to 16,000 characters and 100 keywords per reference.
- Fast, powerful **search** capabilities.
- Dozens of **predefined output formats**, plus the ability to easily design your own additional formats.
- 100% compatible with **WordPerfect, Microsoft Word, PC-Write, WordStar 2000, WordStar Pro, XyWrite.**
- Can also be used with virtually all other word processors.
- Able to **import** references from national databases, other bibliography programs, or almost any other database or text file.

Full Registered System **\$99⁰⁰**

Full money-back guarantee on purchase of Registered System.

Demo System **\$25⁰⁰**

Demo price credited toward subsequent Registered System purchase.

Outside North America, add \$10 shipping charge - U.S. funds, on a U.S. bank.

IBM PC and compatibles VAX-VMS

(Macintosh version planned for mid-1991)

Research Software Design

2718 SW Kelly Street, Suite 181
Portland, OR 97201

(503) 796-1368 FAX: 503-241-4260


Circle No. 22 on Readers' Service Card

Vibration Isolation

Micro-g® Floor Platforms



- Better resolution
- Patented Gimbal Piston® Isolators
- Reduced operator fatigue
- Cleanroom compatible
- More than 60 standard models for all commonly used TEM's and SEM's



TMC


Technical Manufacturing Corporation
15 Centennial Drive, Peabody, MA 01960
1-800-542-9725 • FAX 508-531-8682

Send for New Catalog

TMC-35-R

Circle No. 172 on Readers' Service Card

ELECTROPHORESIS ENGINEERING EXCELLENCE



Serva introduces a complete line of electrophoresis chambers and power supplies. Offered are the latest in vertical and mini-vertical slab chambers, submerged agarose gel units, semi-dry blotter, DNA-sequencer, IEF apparatus and power supplies up to 3500 volts. Innovative design and sturdy construction make these units user-friendly and durable.

SERVA

BIOCHEMICALS • SEPARATION TECHNOLOGY

50 A & S Drive • Paramus, New Jersey 07652
Phone: (201) 967-5900 and 800-645-3412
Telex #275658 Serv • Telefax # (201) 967-8858

Circle No. 256 on Readers' Service Card

Don't Risk Endotoxin Contamination. Introducing PyroBind.™


New PyroBind endotoxin removal products give you selective endotoxin binding—even from concentrated protein solutions. Laboratory to process

scale, 1 ml to hundreds of liters.

For complete information, check the capacity you want and send back the coupon. Or call (508) 481-6700.

- ☐ PyroBind-ST 5,000 EU syringe tip for up to 10 ml
- ☐ PyroBind-LS 200,000 EU pump or syringe for up to 2 liters

- ☐ PyroBind Model 300 1 million EU operation up to 20 liters with Trio™ Bioprocessing System
- ☐ Larger industrial scale modules



Name _____

Company _____

Street _____

City/State/Zip _____

Phone _____



SEPRACOR INC.

Mail to: Sepracor, Inc., 33 Locke Dr., Marlborough, MA 01752

Circle No. 222 on Readers' Service Card

CALBIOCHEM

Expands Range Of Products For Second Messenger Researchers

Elucidate the role of protein phosphorylation and signal transduction with:

- Protein Kinase Inhibitors, H7, H8, H9, H89, HA1004, K-252a, Staurosporine
- CaM Antagonists
W5, W7, W12, Calmidazolium
CaM Binding Domain
- Protein Kinase Activators
PMA, SC-9, SC-10
- *NEW* Caged Calcium
DM-Nitrophen, NITR 5, NITR 5/AM,
NITR 7
- Enzymes/Proteins
Calmodulin, Protein Kinase C
sn-1,2-DG Kinase
- Caged Compounds
InsP₃, cAMP, cGMP, ATP, GTP
- *NEW* cGMP Modulator
LY 83583
- Diacylglycerol Kinase Inhibitor
R59022, R59949
- Tyrosine Kinase Inhibitors
Erbstatin, Genstein,
Lavendustin A, Tyrphostin

CALL OR FAX TODAY:
CALBIOCHEM Corporation
10933 North Torrey Pines Rd.
La Jolla, CA 92037
Phone: (800) 854-3417
FAX: (619) 453-3552
Circle No. 247 on Readers' Service Card



Conformations and Forces in Protein Folding

Barry T. Nall and Ken A. Dill, editors

Protein folding, the self-directed transition from disorganized chains to highly ordered and functional biological structures, is of increasing practical concern for the biotechnology industry and for interpreting DNA sequences. In the biological sciences folding is of major importance in the "self-assembly" process that produces the protein catalysts that facilitate and regulate cellular chemistry. Folding plays a role in such diverse cellular processes as macromolecular transport and assembly, targeting of proteins to intra- or extracellular locations, and in vivo stability of proteins.


Several aspects of folding addressed include forces and interactions important to protein stability and function, methods for determining proteins, studies of alterations in structure in mutant proteins, mechanistic investigations of the folding process, and analyses of auxiliary factors that modify or catalyze protein folding.

1991; ca. 272 pp.; indexed and illustrated; #91-05S — softcover; \$34.95 (members \$27.95); ISBN 0-81768-394-6


CONTENTS

- Preface—G. Rose
- I.** Compact States, Electrostatics, and Folding — K. A. Dill
1. Calculation of the Total Electrostatic Energy of a Macromolecular System: Solvation Energies, Binding Energies, and Conformational Analysis — M. K. Gilson and B. Honig
 2. Electrostatic Effects and Allosteric Regulations in E. Coli Aspartate Transcarbamylase — M. P. Glackin, J. B. Matthew, N. M. Allewell
 3. Charge Effects on Folded and Unfolded Proteins — D. Stigter and K. A. Dill
 4. Compact Polymers — H. S. Chan and K. A. Dill
- II.** Relation of Amino Acid Sequence to Structure and Folding — L. M. Gierasch
5. Can Molecular Evolution Provide Clues to the Folding Code? — S. C. Hardies and L. D. Garvin
 6. Folding and Activity of Hybrid Sequence, Disulfide-Stabilized Peptides — J. H. B. Pease, R. W. Storrs, D. E. Wemmer
 7. ¹H NMR Assignments and Three-Dimensional Structure of Ala14/Ala38 Bovine Pancreatic Trypsin Inhibitor Based on Two-Dimensional NMR and Distance Geometry — H. M. Naderi, J. F. Thomason, B. A. Borgias, S. Anderson, T. L. James, I. D. Kuntz
 8. The Tryptophan Synthase $\alpha_2\beta_2$ Multienzyme Complex: Relationship of the Amino Acid Sequence and Folding Domains to the Three-Dimensional Structure — E. W. Miles
- III.** Folding Mechanisms — C. N. Pace
9. Protein Structure Determination in Solution by Two-Dimensional and Three-Dimensional NMR — A. M. Gronenborn and G. M. Clore
 10. The Mechanism of Protein Folding — O. B. Ptitsyn, G. V. Semisotnov
 11. Conformation States in Acid-Denatured Proteins — A. L. Fink, L. J. Calciano, Y. Gotto, D. Palleros
 12. Characterization of Unfolded and Partially Folded States of Proteins by NMR Spectroscopy — C. M. Dobson, C. Hanley, S. E. Radford, J. Baum, P. A. Evans
- IV.** Auxiliary Factors and Folding: Membranes and Catalysis — B. T. Nall
13. Teaching Proteins to Fold — P. M. Horowitz
 14. Alternate Folding Motifs for Gramicidin: Crystallographic and Spectroscopic Analyses of Polymorphism — B. A. Wallace
 15. Prolyl Isomerase: Role in Protein Folding and Speculation on Its Function in the Cell — F. X. Schmid, K. Lang, T. Kiefhaber, S. Mayer, E. R. Schonbrunner
 16. Protein-Disulfide Isomerase: An Enzyme That Catalyzes Protein Folding in the Test Tube and in the Cell — R. B. Freedman

Conformations and Forces in Protein Folding



Edited by
Barry T. Nall and Ken A. Dill
AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE



Protein Folding: Deciphering the Second Half of the Genetic Code

Lila M. Gierasch and Jonathan King, editors

1990; 342 pp.; indexed and illustrated; color plates #89-18S — softcover; \$39.95 (members \$31.95); ISBN 0-81768-353-9

"Readers of the book will get a vivid impression of the diversity of systems under study and of the impressive variety of techniques employed in folding research."

—Trends in Biochemical Sciences

Order from: AAAS Books, Dept. A40, PO Box 753, Waldorf, MD 20604 (FAX: 301-843-0159). To order by phone (VISA/MasterCard only), call 301-645-5643 (9am-4pm ET) and ask for AAAS. Please specify item number. Individuals must prepay or use VISA/MC. Add \$4.00 postage/handling. For shipments to CA, add applicable sales tax.

American Association for the Advancement of Science

AMINCO•Bowman® Series 2

Luminescence Spectrometer

Exclusive dual light source design:

- flash lamp for phosphorescence and fluorescence
- continuous lamp for high sensitivity fluorescence

Innovative software with a graphical user interface for multitasking, automation, and third party programs


T-Optics™ sample compartment for a wide range of accessories:

- automated sample sipper and sample changer
- stopped flow accessory
- automated polarizers
- and the entire family of T-Optics accessories

Automatic intracellular analyses using fluorescent probes

Aggressive pricing for maximum value in its class!

For those of you who know the name, you can expect the same performance breakthroughs for the 1990s that marked the introduction of the original AMINCO-Bowman®.



For those of you who don't know the name, we would be happy to introduce you.


Call us.

SIMON AMINCO®

810 West Anthony Drive, Urbana, IL 61801 USA
(217) 384-7730 (800) 637-7689 Telex: 206079 FAX: (217) 384-7744

Circle No. 196 on Readers' Service Card

How to become a master manipulator



Locked-in stability

Drift can really put the squeeze on cells, especially in lengthy patch clamp and intracellular recordings. The Newport MX300 is the only Huxley-style micromanipulator with a locking Z-axis that completely eliminates the largest source of drift, gravity-induced vertical "creep." To minimize drift even more, exclusive ball-bearing contact surfaces in the drivetrain linkage prevent the build-up of residual mechanical stresses. Together, these two exclusive refinements deliver unsurpassed overall stability for long-term experimentation.

- Ultra-low drift
- True submicron resolution
- Adjustable probe angle
- Right- and left-handed configurations
- Compatible with Newport motorized drives



Another first: True submicron control

Superfine, 80-pitch threads yield twice the basic resolution of any other Huxley-style manipulator—less than 0.1 μ m per degree of fine rotation.

With the MX300, it's easy to be a master manipulator. Call us for details at **714/965-4631** or fax a message to **714/964-3071** for faster service.



Newport.
Bio-Instruments

Corporate Headquarters
18235 Mt. Baldy Circle
Fountain Valley, CA 92708
Ph: 714/965-4631

Canada: Ph. 416/567-0390
Japan: Ph. 06/359-0270

European Headquarters
Ger.: Ph. 06151-1540
U.K.: Ph. 0582-769995
Switz.: Ph. 01-740-2283
Neth.: Ph. 03402-50588

©1991 03912
Newport Corporation

Circle No. 219 on Readers' Service Card

New Journals in the LIFE SCIENCES

Write or Call for Free Sample Copies

New in 1991

Biological Control

Theory and Application in Pest Management

EDITORS

Raghavan Charudattan, University of Florida, Gainesville

Harry K. Kaya, University of California, Davis

W. Joe Lewis and **Charlie E. Rogers**, USDA-ARS, Insect Biology and Population Management Research Lab, Tifton, Georgia

Biological Control focuses on the reduction or mitigation of pests and pest effects through the use of natural enemies. Research articles and reviews are featured.

Volume 1 (1991), 4 issues
In the U.S.A. and Canada: \$146.00

ISSN 1049-9644
All other countries: \$166.00

Biomedical Science and Technology

*Research, Engineering, and Manufacturing for
Medical Devices and Instruments*

EDITOR

Robert D. Gold, University of Wisconsin, Madison

Biomedical Science and Technology covers the development and production of devices and instruments used in medical practice and publishes peer-reviewed, technical articles and other information of interest and applicability to a broad range of those working to create, develop, and manufacture the materials and products used by the medical community.

Volume 1 (1991), 4 issues
In the U.S.A. and Canada: \$95.00

ISSN 1051-2020
All other countries: \$112.00

New Journals to Appear in 1992

Consciousness and Cognition

EDITORS: **Bernard J. Baars**, The Wright Institute, Berkeley, California
William P. Banks, Pomona College and the Graduate School, Claremont, California

Providing a forum for a natural science approach to the issues of *consciousness, voluntary control, and self*
ISSN 1053-8100

NeuroImage

EDITOR-IN-CHIEF: **Arthur W. Toga**, University of California School of Medicine, Los Angeles

Focusing on the visualization of all neuroscientific data, the journal will publish articles describing neurobiological insights revealed through the use of imaging.
ISSN 1053-8119

Molecular Phylogenetics and Evolution

EDITOR-IN-CHIEF: **Morris Goodman**, Wayne State University School of Medicine, Detroit, Michigan

Establishing a forum for molecular studies that advance our understanding of phylogeny and evolution



Sample copies and privileged personal rates are available on request.

For more information, please write or call:

ACADEMIC PRESS, INC., Journal Promotion Dept.
1250 Sixth Avenue, San Diego, CA 92101, U.S.A.
(619) 699-6742

All prices are in U.S. dollars and are subject to change without notice.

S1277

Circle No. 48 on Readers' Service Card

Drummond Auto/Oocyte Injector

**MICROPROCESSOR
CONTROLLED
NANOLITER
VOLUMES**

- Remote control and non-rotating plunger assure minimum vibration.
- Easily attached to most micro-manipulators, weight only 3.5 oz.
- Model 203-X injects one volume, 46.0 nanoliters, other volumes possible.
- Model 203-XV volume selectable from 4.6 to 73.6 nanoliters, in 16 steps.
- Total capacity .5 μ l
- Supplied with sufficient glass to pull 300 tips.

**2
MODELS
NOW
AVAILABLE**

Catalog #3-00-203-X
Catalog #3-00-203-XV

Price \$495.00
Price \$695.00



**DRUMMOND
SCIENTIFIC COMPANY**

500 Parkway, Box 700, Broomall, PA 19008
1-800/523-7480, In PA 1-215/353-0200
FAX: 1-215/353-6204
MADE IN U.S.A.

Circle No. 175 on Readers' Service Card

BURN YOUR REFERENCE CARDS!

REF-11 TM

Computerizes your REFERENCES
and prepares your BIBLIOGRAPHIES

- ☐ Maintains a data base of references
- ☐ Searches for any combination of authors, years of publication, reference title, publication title, keywords or abstract
- ☐ Formats bibliographies exactly as you want them
- ☐ Reads your paper, inserts citations into the paper, and prepares a bibliography of the references cited (optional)
- ☐ Downloads references from any Online data base including NLM, BRS, DIALOG and Laserdisc formats (optional)

IBM PC/XT/AT, MS-DOS, CP/M 80 ... **\$195.⁰⁰**

RT-11, TSX-Plus, RSX-11, P/OS **\$250.⁰⁰**

VAX/VMS (native mode) **\$650.⁰⁰**






ANY MANUAL & DEMO \$20.⁰⁰

322 Prospect Ave., Hartford, CT 06106
(203) 247-8500

Connecticut residents add 8% sales tax.

Circle No. 31 on Readers' Service Card

Three who mastered science for better medicine



For their achievements, they receive the *Edgar M. Queeny Award*...and mankind receives an innovative new drug.

For the seventh time since 1982, scientific research resulting in a commercial product from Monsanto has

earned the *Edgar M. Queeny Award*. The award was named for the son of Monsanto's founder.

This year's recipients, all from Monsanto's pharmaceutical subsidiary, **Searle**, are: Dr. Paul Collins, section head for gastrointestinal diseases; Dr. Esam Z. Dajani, director of clinical research; and Dr. Raphael Pappo, a retired senior research fellow.

The three developed a valuable preventive for NSAID-induced gastric ulcers. Their work exemplifies Monsanto's belief that science, imaginatively applied, benefits everyone.

Monsanto

ACCUGENE™: THE PERFECT SOLUTIONS FOR MOLECULAR BIOLOGY


AccuGENE™ Solutions Turn Preparation Time Into Productive Time.

How much valuable time does your laboratory spend screening raw materials and preparing solutions for molecular biology applications? Do you sometimes suspect that inconsistencies and impurities in those solutions may influence your results?

Whittaker Bioproducts' complete line of ready-to-use, certified-pure AccuGENE solutions eliminate time-consuming preparation and bring new levels of performance consistency to molecular biology applications. With AccuGENE, preparation time is converted into productive time, and solution variability is replaced with the constant of Whittaker Bioproducts quality.

AccuGENE solutions are produced in meticulous compliance with FDA Guidelines and Good Manufacturing Practices, using only molecular biology-grade reagents and ultrapure water. Each AccuGENE lot is screened and certified to be free of detectable DNase, RNase and protease activity.

There simply isn't an easier, more convenient, or more reliable solution for molecular biology than Whittaker Bioproducts' AccuGENE. Call 1-800-638-8174 or 800-654-4452 for orders and information. Outside U.S. FAX: 301-845-8291. TELEX: 6919009 MABIO UW.



Whittaker
BIOPRODUCTS

Circle No. 270 on Readers' Service Card