




AMERICAN  
ASSOCIATION FOR THE  
ADVANCEMENT OF  
SCIENCE


# SCIENCE

14 AUGUST 1992  
VOL. 257 • PAGES 849-1012

\$6.00



# How are you measuring inflammation?



You should be measuring Serum Amyloid A (SAA) in:

- cholesterol metabolism
- rheumatoid arthritis
- sepsis
- renal allograft rejection
- cytokine research
- bone marrow transplantation
- early infection
- chronic inflammation
- cystic fibrosis
- lung inflammation

## Serum Amyloid A ELISA Kit Exclusively from BioSource International

Serum amyloid A (SAA) is a cysteine-free single polypeptide of 11,685 daltons that is produced by the liver during an acute phase response. SAA is almost entirely associated with high density lipoprotein and can reach a concentration of up to 1.0 mg/ml (100-1000 times the normal physiological level) in the sera of acutely ill patients. Experiments using whole animals, liver sections, cultured hepatocytes, or hepatocyte-derived tumor cell lines have revealed that changes in SAA biosynthesis are caused by alterations in gene transcription regulated primarily by IL-6, often in conjunction with IL-1 and tumor necrosis factor (TNF). Measurements of SAA are useful in studies examining the synthesis and regulation of

inflammatory cytokines, the interaction of glucocorticoids in the enhancement and regulation of inflammation, as well as the overall assessment of the inflammatory response.

### Other ELISA Kits available from BioSource

- hTNF- $\alpha$
- hTNF- $\beta$
- sTNF-R
- hIL-2
- hIL-6
- hIFN- $\alpha$
- hIL-8/NAP-1
- hIFN- $\gamma$
- hICAM-1

If you don't see the kit you're interested in, call us - it's probably under development

*Sensitive*

*Specific*

*Simple*

*Accurate*

*Rapid*

*Economical*

- 5 ng/ml sensitivity
- 5-300 ng/ml dynamic range
- No cross reactivity with other proteins has been observed
- Pre-coated strip-well plates
- Combined step procedure
- C.V. inter-assay 7.8%
- C.V. intra-assay 4.9%
- Only 2 hours required for incubation steps
- Two plates and reagents for 192 determinations for about the same cost as other one-plate kits

It's no wonder that more researchers are switching to BioSource International for their immunology research products.

**(800) 242-0607**




**BIOSOURCE**  
INTERNATIONAL

950 Flynn Rd • Camarillo, CA 93012


CAUTION: FOR RESEARCH USE ONLY IN LABORATORY ANIMALS OR IN-VITRO TESTING. NOT FOR HUMAN OR DIAGNOSTIC USE.

Circle No. 33 on Readers' Service Card



## MILESTONES IN PCR PERFORMANCE.

**Quantitation.** Set new standards in quantitative PCR. The advanced technology of the GeneAmp™ PCR System 9600 with thin-walled MicroAmp™



142-bp HIV-1 PCR products (lanes 1, 2) generated after 40 cycles on the GeneAmp PCR System 9600 using SK145/431 primers.

Reaction Tubes translates into unsurpassed accuracy and uniformity in temperature cycling with precise system control. **Specificity.** Achieve a clearly detectable specific PCR product with your DNA and RNA samples.

**Proven PCR tech-**

nology, combined with the GeneAmp PCR System 9600 and GeneAmp® PCR Reagents, provides the most



sensitive and specific amplification every time. **Simplicity.** Oil-free operation, unparalleled speed and low volume PCR amplification in a convenient, easy-to-program system. Backed by our PCR Performance Guarantee. To order in the U.S., please contact your local Perkin-Elmer sales representative. For PCR technical information in the U.S., call 1-800-



762-4001. Or call 1-800-762-4000 for literature. Outside the U.S., contact your local Perkin-Elmer sales representative.

HIV-1 sample was provided by SRA Technologies Inc., Rockville, Maryland. GeneAmp is a registered trademark and GeneAmp and MicroAmp are trademarks of Hoffmann-La Roche Inc. The PCR process is covered by U.S. patents owned by Hoffmann-La Roche Inc.


### PERKIN ELMER

Europe Vaterstetten, Germany Tel: 49-8106-381-115 Fax: 49-8106-6697  
Canada Montreal, Canada Tel: 514-737-7575 Fax: 514-737-9726  
Far East Melbourne, Australia Tel: 61-3-560-4566 Fax: 61-3-560-3231  
Latin America Mexico City, Mexico Tel: 52-5-651-7077 Fax: 52-5-593-6223

Perkin-Elmer PCR reagents are developed and manufactured by Roche Molecular Systems, Inc., Branchburg, New Jersey, U.S.A.



Circle No. 15 on Readers' Service Card



Number 1 in a Series

# CENTRIFUGATION

*Why do 7 out of 10 choose  
L/XL Optima™ Series Ultracentrifuges?*

**D**r. John Noti, Ph.D. of the Guthrie Foundation in Sayre, Pennsylvania, USA, puts it succinctly: *"The XL-90 is in a class by itself. The short run times we're achieving give us more time to perform actual experiments. Programming is simple and fast. This instrument is very reliable and service is, as always, excellent. A remarkable machine!"*

## RUN TIMES SHORTENED BY AS MUCH AS 90%!

Beckman's ESP™ Efficient Sedimentation Program lets you use high performance rotors to full advantage without danger of precipitating gradient-forming salts. ESP designs protocols for optimum efficiency by simulating run conditions at the rotor's maximum speed. The software determines how the salts will redistribute during the run and automatically adjusts speed to avoid precipitation. By applying significantly higher forces at the early stages of the separation process, dramatic reductions in run times are achieved.

To simplify the set-up of subsequent runs, optimized protocols designed by ESP can be stored for future recall.

## ENVIRONMENTALLY SENSITIVE - THE ONLY ULTRACENTRIFUGE WITH NO CFC'S

As the leading supplier of ultracentrifuges, Beckman's concern about the environmental impact of these instruments led to the development of the L and XL Optima Series. These instruments

are designed with no chlorofluorocarbons. Instead they have an air-cooled drive and solid-state thermoelectric cooling and heating modules.


Power consumption is reduced to 3100 BTU/hr with the L and XL — less than a third that of other ultracentrifuges. And the benchtop Optima TLX generates a low 2000 BTU/hr.

## THE WORLD'S FIRST AND ONLY IMBALANCE-TOLERANT DRIVE

The driving force behind Optima Series instruments is Beckman's proven, vacuum-encased induction drive, now brought to a new level of stress tolerance. This drive is so resilient it eliminates the need for tedious sample weighing. Tubes can be conveniently balanced by eye to within  $\pm 10\%$  or 5 mL, whichever is greater.

## MORE USEFUL INNOVATIONS FROM BECKMAN

Beckman pushes the limits of centrifuge technology in a commitment to increase the efficiency of this essential laboratory process. Still the most widely used separation technique, continual Beckman innovation allows its forces to be applied in new ways.



*ESP-designed protocol with VT1 90  
Rotor reduced run times from 36 hours  
to just 2.5 hours.*

**For more information contact: Beckman Instruments, Inc., 2500 Harbor Blvd., Fullerton, CA, 92634-3100. (800) 742-2345. Service and distribution offices worldwide.**

# BECKMAN



**873**

Tree of modern Human  
origins

## NEWS & COMMENT

Fetal Transplants Show Promise **868**  
The Tissue Bank's Shaky Underpinnings

Chemistry Curricula Edge Toward **871**  
a New World

NSF Holds Plan Close to Its Vest **872**

Mitochondrial Eve: Wounded, But Not **873**  
Dead Yet  
Reading the Bones for Modern Human Origins

## RESEARCH NEWS

Huge Impact Tied to Mass Extinction **878**

Getting Some "Backbone": **880**  
How MHC Binds Peptides

Microgravity Materials Science Strives to **882**  
Stay in Orbit

Possible Evolutionary Role Explored for **884**  
"Jumping Genes"

Science Innovation '92: **885**  
The San Francisco Sequel

## PERSPECTIVES

Intermolecular Interactions **887**  
W. Klemperer

Splicing Takes a Holiday **888**  
J. A. Steitz

## ARTICLES

The Nature of the Metal-Metal Bond **897**  
in Bimetallic Surfaces  
J. A. Rodriguez and D. W. Goodman

## PATENT POLICY

Genes, Patents, and Product **903**  
Development  
R. S. Eisenberg

Genome Research: Fulfilling the **908**  
Public's Expectations for Knowledge and  
Commercialization  
R. G. Adler

Patents on Random Complementary **915**  
DNA Fragments?  
T. D. Kiley



## DEPARTMENTS

**THIS WEEK IN SCIENCE** **857**

**EDITORIAL** **859**  
Science Is Great, But Scientists Are Still People:  
A. Kornberg

**LETTERS** **860**  
Drugs from Third World Plants: I. S. Johnson •  
Computing in Science: J. Ross • Condom Use: F.  
L. Sonenstein • NSF and Duplicate Grant Submis-  
sions: D. Jaron • Combating Epidemic Cholera:  
R. A. Finkelstein • Nucleoside Diphosphate  
Kinase: Conclusions Withdrawn: P. A. Randazzo,  
R. A. Kahn, J. K. Northrup

**SCIENCESCOPE** **867**


**RANDOM SAMPLES** **876**  
Bush Administration Weighs in on AIDS •  
Societies Sound Alarm on Biodiversity, etc.

**BOOK REVIEWS** **981**  
*The Code of Codes*, reviewed by B. D. Davis •  
*Everyday Cognition in Adulthood and Late Life*, T. A.  
Salthouse • *Andean Magmatism and Its Tectonic*  
*Setting*, B. S. Singer • Vignettes: Tales of Invention  
• Books Received

**PRODUCTS & MATERIALS** **987**

**878 & 954**

Dating the K-T impact crater



## AAAS Board of Directors

Leon M. Lederman  
*Retiring President,*  
*Chairman*  
F. Sherwood Rowland  
*President*  
Eloise E. Clark  
*President-elect*  
Mary Ellen Avery  
Francisco J. Ayala  
Robert A. Froesch

Florence P. Haseltine  
Alan Schriesheim  
Jean'ne M. Shreeve  
Chang-Lin Tien  
Warren M. Washington

William T. Golden  
*Treasurer*  
Richard S. Nicholson  
*Executive Officer*

John Abelson  
Frederick W. Alt  
Don L. Anderson  
Stephen J. Benkovic  
David E. Bloom  
Floyd E. Bloom  
Henry R. Bourne  
James J. Bull  
Kathryn Calame  
C. Thomas Caskey  
Dennis W. Choi

John M. Coffin  
Bruce F. Eldridge  
Paul T. Englund  
Richard G. Fairbanks  
Douglas T. Fearon  
Harry A. Fozzard  
Victor R. Fuchs  
Theodore H. Geballe  
Margaret J. Geller  
John C. Gerhart  
Roger I. M. Glass

Stephen P. Goff  
Corey S. Goodman  
Stephen J. Gould  
Ira Herskowitz  
Eric F. Johnson  
Stephen M. Kosslyn  
Michael LaBarbera  
Charles S. Levings III  
Harvey F. Lodish  
Richard Losick  
Anthony R. Means

Mortimer Mishkin  
Roger A. Nicoll  
William H. Orme-Johnson III  
Stuart L. Pimm  
Yeshayau Pocker  
Dennis A. Powers  
Ralph S. Quatrano  
V. Ramanathan  
Erkki Ruoslahti  
Ronald H. Schwartz  
Terrence J. Sejnowski

Thomas A. Steitz  
Richard F. Thompson  
Robert T. N. Tjian  
Emil R. Unanue  
Geerat J. Vermeij  
Bert Vogelstein  
Harold Weintraub  
Zena Werb  
George M. Whitesides  
Owen N. Witte  
Keith Yamamoto

## Board of Reviewing Editors


## COVER


Computer simulation of a molecular beam of the benzene-water dimer. The sizable penetration of the van der Waals surfaces, depicted by the dotted surfaces in the bottom right corner (light blue for water and pink for benzene), illustrates the hydrogen bond interaction of

water with the benzene  $\pi$  electron cloud. See page 942 and the Perspective on page 887. [Computer graphics: Geoffrey A. Blake and Siddharth Dasgupta. Production: Mario Blanco and Carolyn Sherby, Molecular Simulations Inc.]



## RESEARCH ARTICLES


- Crystal Structures of Two Viral Peptides**  **919**  
in Complex with Murine MHC Class I H-2K<sup>b</sup>  
D. H. Fremont, M. Matsumura, E. A. Stura,  
P. A. Peterson, I. A. Wilson

- Emerging Principles for the Recognition**  **927**  
of Peptide Antigens by MHC Class I Molecules  
M. Matsumura, D. H. Fremont, P. A. Peterson,  
I. A. Wilson

## REPORTS

- A Lunar Occultation of the Dust-Scattering** **935**  
**Halo Around GX 5-1 Observed with ROSAT**  
P. Predehl, J. H. M. M. Schmitt, S. L. Snowden,  
J. Trümper


- High-Velocity Pulsars in the Galactic Halo** **937**  
D. Eichler and J. Silk

- Benzene Forms Hydrogen Bonds**  **942**  
with Water  
S. Suzuki, P. G. Green, R. E. Bumgarner,  
S. Dasgupta, W. A. Goddard III, G. A. Blake

- Synthesis of Fluoropolymers in** **945**  
**Supercritical Carbon Dioxide**  
J. M. DeSimone, Z. Guan, C. S. Elsbernd


- Internal Stark Effect Measurement of the** **947**  
**Electric Field at the Amino Terminus of an  $\alpha$  Helix**  
D. J. Lockhart and P. S. Kim

- The Reversal and Splitting of Waves in an** **951**  
**Excitable Medium Caused by an Electrical Field**  
H. Sevcíková, M. Marek, S. C. Müller

- Coeval <sup>40</sup>Ar/<sup>39</sup>Ar Ages of 65.0 Million**  **954**  
**Years Ago from Chicxulub Crater Melt Rock**  
**and Cretaceous-Tertiary Boundary Tektites**  
C. C. Swisher III, J. M. Grajales-Nishimura,  
A. Montanari, S. V. Margolis, P. Claeys, W. Alvarez,  
P. Renne, E. Cedillo-Pardo, F. J.-M. R. Maurrasse,  
G. H. Curtis, J. Smit, M. O. McWilliams

- Pseudo-Half-Knot Formation with RNA** **958**  
D. J. Ecker, T. A. Vickers, T. W. Bruice,  
S. M. Freier, R. D. Jenison, M. Manoharan,  
M. Zounes

- Combining Experimental Information from** **961**  
**Crystal and Solution Studies: Joint X-ray**  
**and NMR Refinement**  
B. Shaanan, A. M. Gronenborn, G. H. Cohen,  
G. L. Gilliland, B. Veerapandian, D. R. Davies,  
G. M. Clore

- A Critical Role for Conserved Residues**  **964**  
**in the Cleft of HLA-A2 in Presentation**  
**of a Nonapeptide to T Cells**  
F. Latron, L. Pazmany, J. Morrison, R. Moots,  
M. A. Saper, A. McMichael, J. L. Strominger

- Differential Display of Eukaryotic** **967**  
**Messenger RNA by Means of the**  
**Polymerase Chain Reaction**  
P. Liang and A. B. Pardee

- Biosynthesis of Human Papillomavirus** **971**  
**from a Continuous Cell Line Upon**  
**Epithelial Differentiation**  
C. Meyers, M. G. Frattini, J. B. Hudson,  
L. A. Laimins

- Rapamycin-Induced Inhibition of the** **973**  
**70-Kilodalton S6 Protein Kinase**  
D. J. Price, J. R. Grove, V. Calvo, J. Avruch,  
B. E. Bierer

- A Point Mutation of the  $\alpha_2$ -Adrenoceptor** **977**  
**That Blocks Coupling to Potassium But Not**  
**Calcium Currents**  
A. Surprenant, D. A. Horstman, H. Akbarali,  
L. E. Limbird




**880 & 927**  
Peptide-in-a-pocket

 Indicates accompanying feature

■ **SCIENCE** (ISSN 0036-8075) is published weekly on Friday, except the last week in December, by the American Association for the Advancement of Science, 1333 H Street, NW, Washington, DC 20005. Second-class postage (publication No. 484460) paid at Washington, DC, and additional mailing offices. Copyright © 1992 by the American Association for the Advancement of Science. The title **SCIENCE** is a registered trademark of the AAAS. Domestic individual membership and subscription (51 issues): \$87 (\$47 allocated to subscription). Domestic institutional subscription (51 issues): \$195. Foreign postage extra: Mexico, Caribbean (surface mail) \$50; other countries (air assist delivery) \$95. First class, airmail, student and emeritus rates on request. Canadian rates with GST available upon request, GST #1254 88122. **Change of address:** allow 6 weeks, giving old and new addresses and 11-digit account number. **Postmaster:** Send change of address to *Science*, P.O. Box 2033, Marion, OH 43305-2033. **Single copy sales:** \$6.00 per issue prepaid includes surface postage; Guide to Biotechnology Products and Instruments, \$20.

Bulk rates on request. **Authorization to photocopy** material for internal or personal use under circumstances not falling within the fair use provisions of the Copyright Act is granted by AAAS to libraries and other users registered with the Copyright Clearance Center (CCC) Transactional Reporting Service, provided that the base fee of \$1 per copy plus \$0.10 per page is paid directly to CCC, 27 Congress Street, Salem, MA 01970. The identification code for *Science* is 0036-8075/83 \$1 + .10. *Science* is indexed in the *Reader's Guide to Periodical Literature* and in several specialized indexes.

■ The American Association for the Advancement of Science was founded in 1848 and incorporated in 1874. Its objectives are to further the work of scientists, to facilitate cooperation among them, to foster scientific freedom and responsibility, to improve the effectiveness of science in the promotion of human welfare, to advance education in science, and to increase public understanding and appreciation of the importance and promise of the methods of science in human progress.



# HIV detection just got more sensitive.

## Introducing the newest AIDS research products from Du Pont.

HIV research depends on sensitive detection. That's why HIV researchers depend on innovative AIDS research products from Du Pont.



### New RT-Detect™

Reproducible non-radiometric detection of reverse transcriptase. Only Du Pont's RT-Detect makes it a practical reality.

This new microplate-based ELOSA (Enzyme-Linked Oligonucleotide Sorbent Assay) kit offers high sample throughput – with minimum hands-on time and excellent sensitivity. And its microplate format is compatible with your existing lab equipment.



### HIV-1 p24 Antigen ELISA Enhancements

Improve on an already good thing? It's easy with two new enhancements to our HIV-1 p24 Antigen ELISA.

Our easy-to-use HIV-1 p24 ADD (Acid Disruption Difference) ICD reagent kit increases

detection of p24 through acid disruption of immune complexes in a convenient one-day format.

Or, boost sensitivity up to fifteen-fold with our ELAST™ ELISA Amplification System.

### And more

We also offer an exceptionally wide range of monoclonals for viral detection – including gp120-neutralizing monoclonal antibodies and recombinant soluble CD4.\*

And, for high-throughput screening of potential anti-AIDS pharmaceuticals, our CD4/gp 120 NENQUEST™ Drug Discovery System redefines the state of the art.

Why not detect the Du Pont difference for yourself? For information by fax, 24 hours a day, 7 days a week, call Du Pont FaxBack® at 1-800-666-6527 (or 302-892-0616) and request document #3000. We'll treat you with maximum sensitivity.

\* CD4 IS SUPPLIED BY AGREEMENT WITH PROGENICS PHARMACEUTICALS, INC.

Products are for research use only. Not for use in diagnostic procedures.

United States 1-800-551-2121 • Canada 1-800-387-8391 • Australia +61 (008) 257149 • Belgium +32 (02) 7242717  
Denmark +45 31506610 • France +33 (01) 69825450 • Germany +49 (06172) 87-2600 • Italy +39 (055) 2478044  
Japan +81 (03) 3224-8763 • Latin America/Asia Pacific Fax +1 (617) 426-2464 • Netherlands +31 (073) 206550  
Sweden +46 (08) 7503700 • Switzerland +41 (01) 8410330 • United Kingdom +44 (0438) 734865.



## Patently obvious?

Patent applications filed by the National Institutes of Health (NIH) for thousands of partial complementary DNA sequences of unknown function have led to a broad discussion of whether this "interim policy" should be pursued; three articles this week explore this issue. Eisenberg (p. 903) presents an overview of the impact of NIH taking a lead role on patent protection rather than private industry and how this may affect product development. Adler (p. 908) discusses the advantages of filing the gene patent and puts this issue in the historical perspective of patent law. Kiley (p. 915) argues that these patent applications do not meet the original intention of the usefulness criterion and proposes a number of corrective legal measures.

## Stardust scattering

Much like the halos around streetlamps on a foggy night, rings of scattered light occur

around stars because of dust grains in the interstellar medium. Predehl *et al.* (p. 935) report measurements of interstellar dust halos at x-ray wavelengths taken with the Röntgen satellite. The bright x-ray source GX 5-1, positioned in the bulge of the Milky Way, was observed as the moon passed between it and the x-ray detectors. This lunar occultation revealed extended emission as GX 5-1 was blocked by the moon. From the details of the scattered radiation, the authors can deduce the dust grain characteristics of these halos.

## Gamma rays and pulsars

More than a decade ago, astronomers discovered an unusual class of objects that emit bursts of gamma rays; the nature of these "busters" is still hotly debated. Eichler and Silk (p. 937) propose that these gamma-ray sources are related to another astronomical puzzle, a population of pulsars (pulsating neutron stars) outside but mov-

ing toward the Milky Way with high velocity. They suggest that the pulsars are formed by the coalescence of compact stellar remnants that may form part of the dark matter thought to be present in galaxies. The mysterious gamma-ray bursts might be produced when these coalesced neutron stars release some of their stored energy.

## Greener chemistry

Fluoropolymers and chlorofluoropolymers, such as Teflon, are usually synthesized in chlorofluorocarbon (CFC) solvents (freons) to maintain solubility. DeSimone *et al.* (p. 945) synthesized these polymers through free-radical chemistry using supercritical CO<sub>2</sub> as the solvent. The advantage of the very high solubility that many compounds have in supercritical CO<sub>2</sub> has normally been exploited in chromatography. The extension to synthesis may allow a reduction in the use of CFCs, which harm the ozone layer when released to the atmosphere.

## Stark results

A direct measure of the electric field at the amino terminus of an  $\alpha$ -helical peptide has been made by Lockhart and Kim (p. 947). They took advantage of the Stark effect; an electric field shifts the absorption spectrum of a covalently attached probe molecule between its ground and excited states. The field at the amino terminus of the peptide was about an order of magnitude greater than what would be expected from the dielectric properties of bulk water. The electric field at the amino terminus did not increase with peptide length.

## A well-timed event

The 180-kilometer-diameter Chicxulub structure buried beneath sediments on the Yucatán Peninsula has emerged as the most likely candidate for the Cretaceous-Tertiary (K-T) boundary impact site. Its age, however, has been poorly constrained, and thus the identification has been uncertain. Swisher *et al.* (p. 954) used <sup>40</sup>Ar/<sup>39</sup>Ar incremental heating to date impact melt rock obtained from drill cores into the structure and to perform an intralaboratory comparison of tektites (glassy droplets likely produced by an impact) from nearby K-T boundary sites. The best ages obtained, 65.0 million years, are indistinguishable from ages for the tektites. Kerr (p. 878) discusses the age and the 15-year search for the impact site.

## Cultural advance

Molecular analysis of the life cycle of viruses requires the propagation of viruses in tissue culture rather than in animals; Meyers *et al.* (p. 971) have developed an in vitro system for the production of human papillomaviruses (HPVs), which are associated with a number of different anogenital tumors. Organotypic (or raft) cell cultures, which re-create aspects of the cellular environment, have not been able to duplicate all of the stages of epithelial cell differentiation linked to virion production. The HPV's life cycle proceeds through the formation of keratinocytes, the epidermal skin cells that contain fibrous proteins such as keratin. Addition of 12-O-tetradecanoyl phorbol-13-acetate to the culture induced morphological differentiation and expression of keratin 10, and also induced virion production.

## Hydrogen bonding and antigen binding

The major histocompatibility complex (MHC) molecules recognize and present the foreign peptides to the cells that mount the immune system's attack. Two studies address how class I MHC molecules recognize foreign peptides. Fremont *et al.* (p. 919) present x-ray structures of two complexes of the mouse MHC class I molecule H-2K<sup>b</sup> with peptides, one with an eight-residue viral peptide and the other with a nine-residue peptide. Analysis of these structures by Matsumura *et al.* (p. 927) suggests that deep, highly conserved pockets at the ends of the peptide-binding groove form hydrogen bonds with the peptide that dictate its orientation. A deep polymorphic pocket in the middle of the groove plays a major role in binding and specificity. In a separate study, Latrone *et al.* (p. 964) studied the effect of mutating conserved residues in a human MHC class I molecule of known structure, HLA-A2, on presentation of influenza virus peptides to T cells. Mutations that appear to disrupt hydrogen bonding in the A pocket, which binds the amino terminal of the peptides, greatly reduced recognition (see news story by Barinaga, p. 880).




# AN EFFICIENT CONSUMER

*Consider the hummingbird: Rapid, efficient liquid handling. Extremely flexible operation. Economical in its use of effort and fluids. The same can be said of the Eppendorf ECOSYN™ D-300 DNA/RNA Synthesizer. Reagent use is minimized through syringe-metered delivery and optical sensors. Operation is completely under the user's control via the flexible program design. Fast cycle times are achieved using zero dead volume valves and short liquid paths. ECOSYN D-300 is ideal for all applications requiring synthetic oligonucleotides.*



ECOSYN™ D-300

**eppendorf®**  
**1-800-421-9988**



Eppendorf North America, Inc. • 545 Science Drive • Madison, WI 53711 • (608) 231-1188 • Fax (608) 231-1339  
 Eppendorf - Netheler - Hinz GmbH • P.O. Box 65 06 70 • D-2000 Hamburg 65 • Tel. (40) 5 38 01-0 • Fax (40) 5 38 01 556 • Teletex-no. 40 30 61=EGHAM

Circle No. 34 on Readers' Service Card

# Pure mRNA in Minutes...

**...Directly from  
Small or Large Samples  
of Cells or Tissue.**

FastTrack™ and MicroFastTrack™ set the industry standard in high quality mRNA isolation.

**MicroFastTrack™\*: 20 Reactions**

- Ideal for PCR, Northern and cDNA synthesis
- Isolation from samples ranging in size from  $10\text{-}3 \times 10^6$  cells or 10-250mg of tissue.
- Reproducible yields of high quality mRNA.

**FastTrack™\*: 6 Reactions**

- mRNA isolation for Northern, cDNA, library construction, PCR, microinjection, RNA protection studies and *in vitro* translation.
- Isolation from samples ranging in size from  $10^7\text{-}10^8$  cells or 0.4-1.0 gram of tissue.
- Fast, efficient recovery of large amounts of polyA+ RNA from a variety of sources.

**Both systems offer:**

- High yields of intact mRNA with low ribosomal contamination.
- Eliminate the need for total RNA isolation or the use of toxic chemicals.
- The most cost effective means of generating high quality mRNA.
- Consistency, convenience and the fastest isolation time.

For the very best in direct mRNA isolation FastTrack™ and MicroFastTrack™ are the choice of thousands of research labs worldwide. When the quality of your mRNA is important, turn to the original source for purity, reliability and convenience; turn to Invitrogen.

**Toll Free 1-800-955-6288**

**Invitrogen**  
CORPORATION

3985 • B Sorrento Valley Blvd. • San Diego, CA 92121  
(619) 597-6200 Phone • (619) 597-6201 Fax

  
BRITISH BIO-TECHNOLOGY  
PRODUCTS LTD  
UK - Tel: (0865) 781045 Fax: (0235) 533420  
France - Numéro Vert 05 90 72 49  
Sweden - 020-Linjen 020 793149  
Norway - Ring Grønt Nummer 050 11033  
Denmark - Grønt Nummer 80 01 85 92  
Belgium - Numéro Vert/Groen Nummer 78 11 04 68

\*patent pending. mRNA model courtesy of BIOSYM


  
CEL BIO  
Italy  
Tel: 39-238103171  
Fax: 39-2381014651

  
ITC  
Germany - Tel: 62-21303907  
Fax: 62-21303511  
Switzerland - Tel: 155-5044

  
FUNAKOSHI CO., LTD.  
Japan  
Tel: 81-356841622  
Fax: 81-356841633

BDH INC., CANADA - TEL: 800-565-7933 • BIO-TRADE, AUSTRIA - TEL: 43-1-8891819 • FINNZYMES, FINLAND - TEL: 35804208077 • MEDOS, AUSTRALIA - TEL: 61-38089077 • SANBIO BV, NETHERLANDS - TEL: 31-413251115 • TAL RON, ISRAEL - TEL: 972-8-472563 • TDI, SPAIN - TEL: 34-14091251 • ECOGEN, SPAIN - TEL: 34-934560607 • UNITED RESEARCH/GOODMAN BIOTECHNOLOGIES, INDIA - TEL: 59 1107

Circle No. 25 on Readers' Service Card



# FTS takes the guesswork out of freeze-drying

FTS gives you full microprocessor control for automatic freeze-drying operation. This affordable system means you won't waste valuable time designing freeze-drying protocol, monitoring start-up and shut-down and watching for overload conditions.

**Check these features against any other system:**

- Five modes of operation including automatic, semi-automatic and programmable
- Shelf temperature stability of  $\pm 0.5^{\circ}\text{C}$
- Reliable electro-mechanical stoppering — without platens or air bladders
- Vacuum, temperature and time can be digitally displayed at any point in the cycle
- Titanium condensers provide superior corrosion resistance

Learn how microprocessor control can improve freeze-drying performance in your lab. Call TOLL-FREE 800-251-1531. In NY, 914-687-0071.

**Systems Engineered for Results**



FTS Systems, Inc.  
PO Box 158  
Stone Ridge, NY 12484  
TWX 510-247-0547  
FAX 914-687-7481

Circle No. 11 on Readers' Service Card





## Now Yields Up to 80% Recombinants

Invitrogen's new version of the TA Cloning™ kit, the most advanced method of direct cloning of PCR amplified nucleic acids, now offers improved cloning efficiencies.

■ The TA Cloning™ kit allows direct cloning of PCR products without modification of primers, restriction digestion or the use of modifying enzymes and new improvements yield the highest cloning efficiencies available.

■ The new pCR™ II vector provides an f1 origin for rescue of single stranded DNA, T7 and Sp6 RNA promoters for sense and anti-sense transcription, the *lacZ* gene for blue/white selection, and a choice of ampicillin or kanamycin resistance to eliminate background when amplifying from plasmids.

■ Invitrogen's TA Cloning™ kit eliminates the need for traditional PCR cloning procedures which are time-consuming, expensive and inefficient and provides the fastest, most reliable method for direct cloning of PCR products.

### Polylinker Region of the New TA Cloning™ Vector, pCR™ II



**Toll Free 1-800-955-6288**  
(619) 597-6200 Phone • (619) 597-6201 Fax

**Invitrogen**  
CORPORATION

3985•B Sorrento Valley Blvd., San Diego, CA 92121

**BRITISH BIO-TECHNOLOGY**  
PRODUCTS LTD.  
UK- Tel: (0865) 781045 Fax: (0235) 533420  
France - Numéro Vert 05 90 72 49  
Sweden - 020-Linjen 020 793149  
Norway - Ring Grant Nummer 050 11033  
Denmark - Grant Nummer 80 01 85 92  
Belgium - Numéro Vert/Groen Nummer 78 11 04 68

**CEL BIO**  
Italy  
Tel: 39-238103171  
Fax: 39-2381014651

**ITC**  
Germany/Switzerland  
Tel: 62-21303907  
Fax: 62-21303511


**FUNAKOSHI CO., LTD.**  
Japan  
Tel: 81-356841622  
Fax: 81-356841633

BDH INC., CANADA - TEL: 800-565-7933 • BIO-TRADE, AUSTRIA - TEL: 43-1-8891819 • FINNZYMES, FINLAND - TEL: 35804208077 • MEDOS, AUSTRALIA - TEL: 61-38089077 • SANBIO BV, NETHERLANDS - TEL: 31-413251115 • TAL RON, ISRAEL - TEL: 972-8-472563 • TDI, SPAIN - TEL: 34-14091251 • ECOGEN, SPAIN - TEL: 34-934560607 • UNITED RESEARCH/GOODMAN BIOTECHNOLOGIES, INDIA - TEL: 59 1107


Circle No. 26 on Readers' Service Card

# Vent<sup>TM</sup> DNA Polymerases

## Thermostability, Fidelity & Versatility From The Ocean Depths



Heat stability of various thermal stable DNA polymerases. All were incubated at 95°C under standard assay conditions.



Reversion frequency reflects error rate in DNA synthesis and was measured by the opal codon reversion assay. Kunkel et al. (1987) Proc. Natl. Acad. Sci. USA 84, 4865-4869. Mattila P. et al. (1991) Nucleic Acids Res. 19, 4967-4973.

For high fidelity, thermostability up to 100°C, primer extensions consistently in the 5-10 Kb range, thermal cycle DNA sequencing and high temperature dideoxy DNA sequencing, discover New England Biolabs' complete range of cloned thermostable Vent<sup>TM</sup> DNA polymerases:

- Vent<sup>TM</sup> DNA Polymerase
- Deep Vent<sup>TM</sup> DNA Polymerase
- Vent<sup>TM</sup> (exo<sup>-</sup>) DNA Polymerase

For thermal cycle DNA sequencing, choose our new easy-to-use system:

- CircumVent<sup>TM</sup> DNA Sequencing Kit



Tube worms at 2010 meters around a hydrothermal vent where the native organism containing Deep Vent<sup>TM</sup> DNA Polymerase was isolated.

### Vent<sup>TM</sup> DNA Polymerase

High fidelity thermostable DNA polymerase with fidelity levels up to 20 fold higher than Taq DNA polymerase due to a 3'→5' proofreading exonuclease activity.

- Discovered and cloned at NEB
- High fidelity
- Primer extensions up to 13 Kb
- Stable up to 100°C
- Half life of 6.7 hours at 95°C

### Deep Vent<sup>TM</sup> DNA Polymerase

NEB's second thermostable DNA polymerase with a 3'→5' proofreading exonuclease activity. Originally isolated from an ocean submarine thermal vent at 2010 meters, Deep Vent<sup>TM</sup> DNA Polymerase is even more thermostable than Vent<sup>TM</sup> DNA polymerase at temperatures of 95°C to 100°C.

- Discovered and cloned at NEB
- Superior stability up to 100°C
- Half life of 23 hours at 95°C

### Vent<sup>TM</sup> (exo<sup>-</sup>) DNA Polymerase

Genetically modified form (exo<sup>-</sup>) of the cloned Vent<sup>TM</sup> DNA polymerase with fidelity levels up to 2-fold higher than Taq DNA polymerase. Preferred form for high-temperature dideoxy DNA sequencing and thermal cycle DNA sequencing.

- Discovered and cloned at NEB
- High temperature dideoxy sequencing
- Primer extensions up to 13 Kb
- Stable up to 100°C

### CircumVent<sup>TM</sup> DNA Sequencing Kit

The only thermal cycle sequencing kit that uses the Vent<sup>TM</sup> (exo<sup>-</sup>) DNA polymerase.

- Easier than conventional methods
- Fast - no need to collapse double stranded plasmids, eliminates centrifugation steps and independent priming steps. You save valuable research time.
- Requires only nanograms of template (only femtomoles when using <sup>32</sup>P end labelled primers)
- Allows direct sequencing from colonies, plaques, cosmids, or DNA fragments isolated from agarose gels
- Diminished secondary structure effects due to high temperature reaction
- Compatible with radiolabel or chemiluminescent detection

These products do not carry a license to perform the Polymerase Chain Reaction (PCR) process. PCR is covered by U.S. patents issued to the Cetus Corporation.

Circle No. 27 on Readers' Service Card

- New England Biolabs Inc. 32 Tozer Road, Beverly, MA 01915 USA 1-800-NEB-LABS (US and MA) Tel. (508) 927-5054 Fax (508) 921-1350
- New England Biolabs Ltd., Canada Tel. (800) 387-1095 (416) 672-3370 Fax (416) 672-3414
- New England Biolabs GmbH, Federal Republic of Germany Tel. (06196) 3031 Fax (06196) 83639

**DISTRIBUTORS:** AUSTRALIA GENESEARCH Tel. (075) 94 0299 FINLAND, SWEDEN, DENMARK, FINNZYMES (Finland) Tel. (0) 420-8077 / FRANCE OZYME Tel. (1) 30 57 0025 INDIA BIOTECH INDIA Tel. (542) 311473 ISRAEL GAMIDOR Tel. (03) 535-1205 ITALY C.A.M.Bio Tel. (02) 38103171 / JAPAN DAIICHI PURE CHEMICALS CO. LTD. Tel. (03) 3272-0671 KOREA KORAM BIOTECH Tel. (02) 556-0311 THE NETHERLANDS WESTBURG Tel. (033) 95 00 94 NORWAY ING. F. HEIDENREICH Tel. (02) 22 04 11 PEOPLE'S REPUBLIC OF CHINA CUBC Tel. (1) 256-1627 PORTUGAL QUIMIGRANEL Tel. (1) 858 15 64 SPAIN LANDERDIAGNOSTICO Tel. (01) 594 08 06 SWITZERLAND FLOW LABORATORIES Tel. (061) 4814713 / TAIWAN LONG CHAIN INTERNATIONAL Tel. (02) 552-2605 UK CP LABORATORIES Tel. (0279) 758200



## BURN YOUR REFERENCE CARDS!

NEW VERSION

# REF-11™

SINCE 1983

Computerizes your REFERENCES  
and prepares your BIBLIOGRAPHIES

- ☐ Maintains a data base of references
- ☐ Searches for any combination of authors, years of publication, reference title, publication title, keywords or abstract
- ☐ Formats bibliographies exactly as you want them
- ☐ Reads your paper, inserts citations into the paper, and prepares a bibliography of the references cited (optional)
- ☐ Downloads references from any Online data base including NLM, BRS, DIALOG and Laserdisc formats (optional)

IBM PC/XT/AT, MS-DOS, CP/M 80 ... **\$195.00**

RT-11, TSX-Plus, RSX-11, P/OS ..... **\$250.00**

VAX/VMS (native mode) ..... **\$650.00**



ANY MANUAL & DEMO ..... **\$20.00**

322 Prospect Ave., Hartford, CT 06106  
(203) 247-8500

Connecticut residents add 8% sales tax.

Circle No. 3 on Readers' Service Card

## American Type Culture Collection

### BIOLOGICAL CULTURES

"NOW AVAILABLE"  
ATCC CULTURE CATALOGUE ON DISKETTE

#### SECTIONS


- ALGAE
- ANIMAL VIRUSES
- BACTERIA
- CELL LINES
- DNA MATERIAL
- FUNGI
- PLANT VIRUSES
- PROTOZOA
- VIRUS ANTISERA
- YEASTS

- IBM COMPATIBLE- 3 1/2" OR 5 1/4" DISKETTES

- BUILT-IN SEARCH SOFTWARE

- \$20-\$25 PER SECTION

- 20% DISCOUNT ON MULTIPLE SECTIONS




#### ORDERING INFO:

CALL  
1 (800) 638-6597  
OR (301) 881-2600  
FAX (301) 231-5826

ATCC, 12301 Parklawn Dr.  
Rockville, MD. 20852 USA

Circle No. 37 on Readers' Service Card

# Show off your results!



"If you need a tool for creating good-looking technical presentations, SlideWrite Plus makes an excellent choice."

PC Magazine, 1991

It's easy with SlideWrite Plus – the PC presentation software that supports all popular printers, plotters, and film recorders. SlideWrite Plus gives you technical graphing, word charting, science figure packs, scalable fonts, plus freehand drawing. No wonder it's the perfect presentation tool for the scientific and technical community! If you want to turn your data into powerful presentations that get noticed and remembered, you need SlideWrite Plus!

Order SlideWrite Plus for just \$445, satisfaction guaranteed, or get a FREE working demo.

**(619) 931-1919**

FAX (619) 931-9313

Advanced Graphics Software, Inc.

5825 Avenida Encinas, Suite 105, Carlsbad, CA 92008

NOW  
AVAILABLE IN  
WINDOWS™  
OR DOS

**SlideWrite Plus™**  
We make your data look terrific.

Circle No. 35 on Readers' Service Card



# The Future of Peptide Technology is in the Palm of Your Hands

## Look Into the Future of Peptide Technology

Imagine peptides and reagents of the highest possible quality and value. Imagine new product breakthroughs in all aspects of peptide research and chemistry. Imagine a company so advanced that it handles everything from custom synthesis at the milligram level, to large-scale peptide production of kilogram quantities quickly, confidentially and at the most competitive prices.

## The Future Is Now

Bachem California's state-of-the-science technology and 20 years' experience brings you the most comprehensive, up-to-date selection of peptides, offering unequalled quality and value. For basic research or clinical applications, our advanced technology, combined with our unique scale-up capability, is your assurance of quality and consistency, regardless of the complexity of the peptide or the quantity you require.

**When it's Peptide Quality and Value You're Looking for Count on  
Bachem California and Take Control of Your Future**



**Bachem California**  
3132 Kashiwa Street • Torrance, CA 90505  
Tel: (310) 539-4171  
Telefax: (310) 530-1571

**Bachem (UK), Ltd.**  
A subsidiary of Bachem California  
69 High Street • Saffron Walden  
Essex CB10 1AA, England  
Tel: (0799) 526465  
Telefax: (0799) 526351

**Saxon Biochemicals, GmbH**  
A subsidiary of Bachem California  
Feodor-Lynen-Strasse 14  
3000 Hannover 61, Germany  
Tel: (0511) 530080 • Telefax: (0511) 5300888

Circle No. 13 on Readers' Service Card

# REAL-TIME WINDOWS

**Continuous data acquisition at 1 KHz while Windows multitasks in real-time!**

With LABTECH for Windows, you get a real-time environment for Microsoft Windows, and an icon-driven application package that lets you start collecting data and controlling your system immediately.

With support for over 50 vendors of data acquisition hardware, you'll never be locked

in to anyone's proprietary software. We support hardware from ADAC, Advantech, Analog Devices, Burr-Brown, Computer Boards, Cyber Research, Data Translation, Keithley/Metrabyte, Hewlett-Packard, Industrial Computer Source, Fluke, National Instruments, and Omega Engineering, as well as GPIB and RS-232 devices.


For help setting up a real-time system for Windows, call us.

 **1-800-TRY-LABTECH**

1-800-879-5228 or 508-657-5400

Circle No. 30 on Readers' Service Card

# DNA SEQUENCING



Call  
**1-800-5410-DNA**  
to discuss how  
we can work as an extension  
of your research laboratory.  
(Outside U.S. call:  
**+1-301-330-0111**)

**Sequencing strategies  
custom-tailored for your  
projects and budget.**

- Manual sequencing (Sanger Method)
- Automated (ABI 373A)
- Rapid turn-around time
- Sequences guaranteed
- Complete confidentiality assured

Services to support your Research Projects:

- Sub-Cloning
- In-vitro Mutagenesis
- Phage & Plasmid DNA Production
- Probe Labeling
- Oligonucleotide Synthesis



**Lofstrand Labs Limited**

*"Highest quality sequencing at the lowest price"*

Circle No. 36 on Readers' Service Card

# The Meeting Becomes...

Cardiovascular Biology  
Cell Injury  
Epithelial Cell Biology  
Growth, Development  
and Aging  
Inflammation  
Mechanisms of Molecular  
Regulation  
Metabolism  
Neuroregulation

Experimental  
Biology

93

The American Physiological  
Society  
American Society for  
Pharmacology and Experimental  
Therapeutics  
American Association of  
Pathologists  
American Institute of Nutrition and  
The American Society for  
Clinical Nutrition

The Biomedical Engineering  
Society  
Society for Experimental Biology  
and Medicine  
International Society for  
Bioelectricity  
North American Society for  
Biorheology

March 28 – April 1, 1993

New Orleans, Louisiana

Detailed announcements, with abstract, registration and housing forms will be mailed to members of all the participating Societies listed above. Members of the American Society for Biochemistry and Molecular Biology, the American Association of Immunologists and the American Society for Cell Biology may submit abstracts for the meeting. **ASBMB**, **AAI** and **ASCB** members will *not* automatically be sent the announcement. Those members who wish to receive details and forms should complete and return the Request for Forms. Nonmembers should also complete and return the Request for Forms.

Detach and mail in  
envelope addressed to:

## FASEB

Office of Scientific Meetings  
Room 3200  
9650 Rockville Pike  
Bethesda, MD 20814-3998  
USA

## Request for Forms – Experimental Biology 93

Deadline for receipt of abstracts: Tuesday, November 17, 1992

### Please indicate Society affiliation:

- |                                |                                |
|--------------------------------|--------------------------------|
| <input type="checkbox"/> APS   | <input type="checkbox"/> ASBMB |
| <input type="checkbox"/> ASPET | <input type="checkbox"/> AAP   |
| <input type="checkbox"/> AIN   | <input type="checkbox"/> AAI   |
| <input type="checkbox"/> ASCB  | <input type="checkbox"/> BMES  |
| <input type="checkbox"/> SEBM  | <input type="checkbox"/> ISB   |
| <input type="checkbox"/> NASB  |                                |

### Please indicate Forms required:

- |   | Quantity |
|---|----------|
| <input type="checkbox"/> Abstract booklet containing all forms and instructions | _____    |
| <input type="checkbox"/> Abstract forms only                                    | _____    |
| <input type="checkbox"/> Advance registration card only                         | _____    |
| <input type="checkbox"/> Hotel application form only                            | _____    |

Not included in booklet, unless requested – Placement Service Information:  
Candidate \_\_\_\_\_ or Employer \_\_\_\_\_


### Please print your name and complete mailing address with Zip Code:

NAME \_\_\_\_\_  
DEPARTMENT \_\_\_\_\_  
INSTITUTION \_\_\_\_\_  
ADDRESS \_\_\_\_\_  
CITY/STATE/ZIP \_\_\_\_\_  
COUNTRY \_\_\_\_\_

Journals come in various shapes, styles and sizes.

# Journal

Ours you plug in.



## THE ONLINE JOURNAL OF CURRENT CLINICAL TRIALS

Editorial excellence, technological know-how

1992 subscription price is \$95 plus telecommunication charges  
For subscription information call (202) 326-6446



The Online Journal of

**Launched  
July 1**

# CURRENT CLINICAL TRIALS

## July Contents

### Articles

Intensive therapy with Cisplatin, Interleukin-2 and Interferon-Alpha-2A in patients with metastatic melanoma: A phase II study, **Sznol M, Steis RG, Smith JW, et al.**, 1992 Jul 1 **Doc. No. 9**

Does episiotomy prevent perineal trauma and pelvic floor relaxation?, **Klein MC, Gauthier RJ, Jorgensen SH, et al.**, 1992 Jul 1 **Doc. No. 10**

A clinical trials database as a research tool in health care, **Morris RD, Lau J, Arena NJ, et al.**, 1992 Jul 17 **Doc. No. 14**

### Editorials

Is the medical world ready for electronic journals?, **Huth EJ**, 1992 Jul 1 **Doc. No. 7**

The publisher's perspective, **Nicholson RN**, 1992 Jul 1 **Doc. No. 8**









Will publication bias vanish in the age of online journals?, **Berlin JA**, 1992 Jul 8 **Doc. No. 12**

Searching the Online Journal of Current Clinical Trials for "old" documents, **Huth EJ**, 1992 Jul 14 **Doc. No. 13**

For information on  
manuscript submission, contact the  
Managing Editor at (202) 326-6735,  
fax (202) 842-2868.

A joint venture of the American Association for  
the Advancement of Science (AAAS) and  
OCLC Online Computer Library Center, Inc.,  
1333 H Street, NW, Washington, DC 20005


# Improved Chemical Substance Searching from CAS



Circle No. 1 on Readers' Service Card

# A P E X

## THE MOST INTELLIGENT APPROACH TO YOUR NEXT SYNTHESIS



### Taking You Way Beyond Traditional QSAR

#### What is being first to market with your new product worth to you?

BIOSYM brings computer - based lead generation into your lab. Apex is a much more powerful tool than presently available conventional programs. The Apex activity prediction system can speed up your search for the next compound for synthesis. Apex goes beyond traditional QSAR methods to function as an expert system - a system that accumulates and extracts knowledge on structure - activity relationships to predict activity and design novel bioactive compounds.

And Apex is backed by the world's largest and preeminent supplier of advanced molecular modeling software. BIOSYM has 10 offices worldwide and 80 PhDs on a staff of over 200 people to support you.

- **Advanced Knowledge - Based System**
- **Dramatically Shortening Your Time to Synthesis**
- **Less Time Means Lower Costs**
- **Seamless, User - Friendly Interface**
- **Completely Integrated With *Insight***
- **Adds Our Intuition To Yours**

When it comes to modeling software go with the leader. For a free demo videotape or more information call our hotline at 619-625-4555 (San Diego) or +89-688-2201 (Munich).



Corporate Address:  
9685 Scranton Road  
San Diego, CA 92121  
Fax: 619-458-0136

Apex has been developed by  
DCL Systems International.