# **VACUUM GAUGE HEAD , BAYARD-ALPERT TYPE**

Glass envelope, ultra-high vacuum gauge head of the Bayard-Alpert type. Measuring range  $10^{-3}$  torr to  $10^{-10}$  torr; sensitivity approx. 12 per torr.

The gauge head is provided with an electrically conductive layer deposited on the inside of the glass envelope. By applying a fixed potential to the layer, excess primary electrons are attracted directly to the envelope rather than oscillating around the collector thereby leading to very stable measurements of low pressure.


Moreover the gauge head features a low thermal inertia and a low filament power


Moreover the gauge head features a low thermal inertia and a low filament power consumption.

### **CHARACTERISTICS**

Pressure range	10-3 to 10-10	torr
Sensitivity (for nitrogen)	approx. 12 per	torr
Emission current range $1~\mu\mathrm{A}$		mA
Filament characteristics	see page 3	
Insulation resistance		
Collector to other electrodes	min. $10^{14}$	Ω
Grid to other electrodes	min. $10^{12}$	Ω

### TYPICAL OPERATING CONDITIONS





Emission current (see also page 3)

measurement

 $100~\mu\text{A}$ , 1~mA or 10~mA

outgassing

75 mA

### LIMITING VALUES

Pressure (filament litt)
Filament voltage
Emission current
Grid input power
Bulb temperature during operation

 $max. 10^{-3}$ torr 8 V max. 75 max. mΑ W max. 40  $^{\rm o}C$ max. 100 oС 450 max.


Dimensions in mm


### MECHANICAL DATA


Bake-out temperature

Material tubulation G28 glass

Filament Tungsten











## **IOG12**

page	sheet	date
1	1	1969.01
2	2	1970.01
3	3	1969.01
4	FP	2001.05.17