

cleanser keeps your calibrated glassware intact.

AVAILABLE IN

 Etching, corroding and weakening of the glass are eliminated. Alconox is ideal for washing all types of laboratory equipment, filters, glass, hypo needles, slides, pipettes, test tubes and other instruments and equipment.

SAVE TIME—SAVE LABOR with ALCONOX. Simply wash and rinse. Alconox leaves no streaks or film. It is kind to the skin. Economical too, a spoonful makes a gallon.


If your dealer cannot supply you, write for literature and samples.

ALCONOX cleans better cleans quicker

Dept. S1 ALCONOX, INC. 61 Cornelison Ave., Jersey City, N. J.



GLASS ABSORPTION CELLS made KLETT



Klett-Summerson Photoelectric Colorimeters— Colorimeters— Nephelometers— Fluorimeters— Bio-Colorimeters— Comparators— Glass Standards—Klett Reagents.

Klett Manufacturing Co. 179 East 87 Street, New York, New York