


*there
must be
a reason*



Here are eight distinct reasons why **Haemo-Sol** is the preferred cleaner for laboratory glassware

- Completely Soluble
- Leaves No Residue
- Fully Haemo-Digestive
- Lower Surface Tension
- Scientifically Formulated
- Controlled pH
- Multiple Laboratory Uses
- Greater Economy

*for literature and samples
write*

MEINECKE & COMPANY, INC.

225 Varick St. • New York 14, N. Y.

Now available...


Adenosine Triphosphate (ATP); Amygdalin; Amylase; Animal Lecithin; Ascorbic Acid Oxidase; Bacitracin; BAL; Biotin cryst.; Caffeic Acid; Carotene cryst.; Catalase; Cellulase; Chlorogenic Acid; Chorionic Gonadotropin; Circulatory Hormone; Colchicine; Cytidylic Acid; Cytochrome C; Dehydroascorbic Acid; Dihydroxyacetone; Dihydroxyphenylalanine (DOPA); Dipyrldyl; Edestin; Emulsin; Erythritol; N-Ethyl-piperidine; Fibroin; Folic Acid; Galacturonic Acid; Gentisic Acid; Girard Reagents P and T; Gliadin; Glucose Pentaacetate; Glucuronic Acid; Glyceraldehyde; Glyceric Acid; Heparin; Hordenine; Hyaluronidase; 2-Hydroxyadipaldehyde; Humulon; Isoascorbic Acid; Isopropylarterenol; Kojic Acid; Kynurenic Acid; Lanthionine; Lipase; Lysozyme; Lyxose; Malononitrile; Maltase; Melezitose; Mesobilirubinogen; Muscle Adenylic Acid; Myanasin; p-Nitrophenylphosphate; Nucleoprotein; Orcinol; Pancreatin; Pangestin; Pantothenyl Alcohol; Penicillinase; Peroxidase; Phenazine; Phenylpyruvic Acid; Phloridzin; Phosphorylase; Piperin; Porphyrindene; Protamines; Pyridoxal; Pyridoxamine; Pyrocatechuic Acid; Pyruvic Aldehyde; Ribonuclease; Saccharic Acid; Salmine; Serine Phosphoric Acid; Spermidine; Spermine; Thioacetic Acid; Thiocytosine; Thyroxine; Trigonellin; Triphenyltetrazolium Chloride; Tripyridyl; Trypsinogen; Tyrosinase; Tyrothricin; Urease; Uricase; Uridine; Vitellin; Xanthosine.

Ask us for others!


DELTA CHEMICAL WORKS

23 West 60th St. New York 23, N. Y.

Telephone PLaza 7-6317



Give point to your lectures with the
Ednalite **PROJECTION
POINTER**



The perfect precision pointer for the lecturer who requires a sharply defined arrow-image which clearly indicates the specific point of discussion, at greater projection distances.

- Uses standard, readily-available, easily replaceable 32 candle-power, low-voltage automotive lamp.
- Operates on ordinary house current, 115 Volt AC, with transformer built into DELUXE, permanent carrying case. Complete unit Underwriters' approved.
- Made from precision, lathe-machined, light-weight metals, completely anodized. Perfectly hand-balanced with feather-touch momentary switch.
- Aspheric, highest-efficiency optical system, projecting a crisp, sharp, intensely brilliant arrow on the screen for extra-long projection throw.
- Helical, micro-focusing ground and polished objective.
- Supplied with full 25 feet of detachable line cord—all contained in solid wood-framed leatherette-lined fitted instrument case.



\$49.50

LIST PRICE

complete with line
cord, transformer,
and carrying case.

Manufactured by *Ednalite* OPTICAL CO., INC.

200 North Water St., Peekskill, N. Y.