

Now available...

Adenosine Triphosphate (ATP); Amygdalin; Amylase; Animal Lecithin; Ascorbic Acid Oxidase; Bacitracin; BAL; Biotin, cryst.; Caffeic Acid; Carotene, cryst.; Catalase; Cellulase Chlorogenic Acid; Chorionic Gonadotropin; Circulatory Hormone; Colchicine; Cytidylic Acid; Cytochrome C; Dehydroascorbic Acid; Diacetonamine; Dihydroxyacetone; Dihydroxyphenylalanine (DOPA); Dipyrityl; Edestin; Emulsin; Erythritol; N-Ethylpiperidine; Fibroin; Folic Acid; Galacturonic Acid; Gentisic Acid; Girard Reagents P and T; Gliadin; Glucose Pentaacetate; Glucuronic Acid; Glyceraldehyde; Glyceric Acid; Heparin; Horde-nine; Hyaluronidase; Hydrindene; 2-Hydroxyadipaldehyde; Humulon; Indan; Isoascorbic Acid; Isopropylarterenol; Kojic Acid; Kynurenic Acid; Lanthionine; Lipase; Lysozyme; Lyxose; Malononitrile; Maltase; Melezitose; Mesobilirubinogen; Muscle Adenylic Acid; p-Nitrophenylphosphate; Nucleoprotein; Orcinol; Pancreatrin; Pantothenyl Alcohol; Penicillinase; Peroxidase; Phenazine; Phenylpyruvic Acid; Phloridzin; Phosphorylase; Piperin; Porphyrindine; Protamines; Protoporphyrin; Pyridoxal; Pyridoxamine; Pyrocatechuic Acid; Pyruvic Aldehyde; Ribonuclease; Saccharic Acid; Salmine; Serine Phosphoric Acid; Spermidine; Spermine; Thioacetic Acid; Thiocytosine; Thyroxine; Trigonelline; Triphenyltetrazolium Chloride; Tripyridyl; Trypsinogen; Tyrosinase; Tyrothricin; Urease; Uricase; Uridine; Vitellin; Xanthosine.

Ask us for others!

DELTA CHEMICAL WORKS

23 West 60th St. New York 23, N.Y.
Telephone PLaza 7-6317

LaMOTTE

COMBINATION SWIMMING POOL SET

for small pools and those who do not care to install separate units for chlorine and alkalinity.



This unit is applicable to the ranges of chlorine and alkalinity actually required for swimming pool water, and will enable the operator to keep his alkalinity and chlorine values between the limits required by State Health Department specifications.


They are complete with 3 color standards for each test representing the upper and lower limits of chlorine (0.2, 0.5, 0.7) and alkalinity (pH 7.2, 7.6, & 8.0) together with reagents and necessary glassware.

These units are contained in a plastic comparator block, the top of which serves for making tests. The plastic blocks are moisture resistant and ideal for use around pools.

LaMotte Chemical Products Company

Dept. "H"

Towson, Baltimore 4, Md.



it's uniform...
**insist on PURINA
Laboratory Chow**

Controlled conditions are a "must" for exacting results in research using laboratory animals.

One factor you can count on is a uniform diet when you feed Purina Laboratory Chow. It's the same formula day after day, month after month.

Buy Laboratory Chow Checkers from your local Purina Dealer. Pressed in a convenient pellet form, it's clean and easy to feed.

RALSTON PURINA COMPANY • St. Louis 2, Missouri