




Get excited

Transfected and functional neurons

- > Primary neurons and PC-12 cells
- > Non-viral transfection
- > Up to 60% efficiency
- > Maintenance of functionality
- > Long-term survival
- > Delivery of DNA straight into the nucleus

Convince yourself.

Try the Nucleofector™ technology.

Check www.amaxa.com/neurobiology




Figure shows primary dissociated mouse neurons of mixed glial cells nucleofected with eGFP. 48 hours post nucleofection, the cells were analyzed by light (A) and fluorescence microscopy (B). (Courtesy of Dr. A. Dityatev, Center for Molecular Neurobiology, Hamburg, Germany).

amaxa Europe / Export
techservice@amaxa.com
+49 (0) 221 99199-400

amaxa USA
techservice.US@amaxa.com
240-632-9110

a m a x a

gene transfer begins here

visit us

www.amaxa.com

AWARDS & PRIZES

GENOMICS

WEST COAST

Submit Search

Call my agent

Science @
CAREERS

Meetings and Announcements
@ www.sciencemeetings.org

SCIENCE'S COMPASS

Oct., p. 356) highlights the inequitable situation of women in the natural sciences. However, unlike in sports, women in science do not need separate teams to be established to enable them to succeed. They will achieve parity once they are provided with a "level playing field" upon which they can compete directly with men. My observation is that the scientific contributions of female scientists tend to be undervalued and their opinions discounted by men. This unconscious bias among many male scientists allows the men to believe that there is no discrimination against women and that they are merely acting to "maintain academic standards." The bar is set higher for women scientists, and it continues to be raised throughout their careers. Independent objective assessment of the achievements of men and women scientists and their contributions within their institutions could help to alleviate this bias.

MARGARET L. SMITHER-KOPPERL

Smither-Kopperl Consulting, 19024 SW 191 Avenue, Williston, FL 32696, USA. E-mail: MSKConsultant@aol.com

Atacamite in Jaws and Printed Wiring Boards

IN THEIR REPORT "HIGH ABRASION RESISTANCE with sparse mineralization: copper biomineral in worm jaws" (11 Oct., p. 389), H. C. Lichtenegger *et al.* identify atacamite in the jaws of *Glycera dibranchiata*, list four polymorphs, and mention its formation in very arid climates and in seawater. Atacamite also forms under very specific electrochemical conditions as conductive anodic filament (CAF) in printed wiring boards (1). These filaments occur between conductors held at a potential difference and result from an electrochemical migration process that initiates at the anode and proceeds along separated glass fiber/epoxy interfaces. Dendritic growth differs from CAF formation. In dendritic growth, the metal ions go into solution at the anode but plate out at the cathode, forming needles or tree-like formations. In contrast, CAF forms initially at the anode and the filament extends with time to the cathode, producing an electrical short circuit. One wonders whether the conditions of this inorganic route to atacamite shed light on the biological processes used in forming atacamite as a biomineral.

In printed wiring boards, temperatures of 75°C, 85°C, or 95°C combined with 85% relative humidity, 150 or 200 V potential difference between copper conductors (0.5- or 0.75-mm spacing), and different solder flux formulations produce CAF in printed wiring boards. A necessary con-